

*The first address for
Precision measurement*

Laboratory Products



SCHOTT
Instruments

Welcome to SCHOTT Instruments



With this new catalogue we are presenting you our complete laboratory products: The catalogue covers the product areas of electrochemical meters and electrodes, titrators, spectrophotometers, laboratory hotplates and stirrers as well as our extensive range of viscometry products – including capillary glass viscometers and viscosity measurement systems.

Electrochemical measuring methods and capillary viscometry are two areas of measuring technology which have become increasingly important in fields such as general science, research and production monitoring. We have been involved in these areas right from the very beginning and have repeatedly succeeded in generating innovative impulses in the form of new products and technology. Since 2006 spectrophotometers complete our traditional programme. With this catalogue we would like to draw your special attention to our latest product developments, which you will come across in all product sections and hope you enjoy discovering these novelties.

Our customers, to whom we would like to express our heartfelt thanks, have made an enormous contribution to our success. Your analytical requirements, thoughts and experience have encouraged us to rise to the challenge again and again. The result of this fruitful dialogue can be seen in this catalogue.

We at SCHOTT Instruments look forward to a continuation of this close relationship with our customers and hope that our new catalogue will help us to support your needs.

SCHOTT Instruments GmbH

Dr. Robert Reining
Managing Director

Contents

News		Page 2
Meters and electrodes for electrochemistry	Laboratory meters	Page 6
	Portable meters for field analysis	Page 28
	Laboratory electrodes	Page 40
	IoLine pH electrodes	Page 56
	ScienceLine electrodes	Page 60
	BlueLine electrodes	Page 80
	Connection cables	Page 84
	Solutions	Page 86
	Electrolyte bridges, other accessories	Page 92
	Notes for successful measurement	Page 94
	Index laboratory electrodes	Page 96
	ProcessLine electrodes for analytics in process	Page 98
Titration	Titration	Page 100
Titration software	TitroLine and TITRONIC®	Page 102
	TitroLine alpha <i>plus</i>	Page 118
	TW alpha <i>plus</i> sample changer	Page 128
	TitriSoft 2.6	Page 130
	TitriSoft 2.6 P	Page 134
Spectrophotometry	Spectrophotometers	Page 140
Capillary viscometry	Automatic viscosity measurement systems	Page 168
	Thermostats and accessories	Page 194
	Capillary viscometers	Page 198
Laboratory equipment	Hotplates and stirrers for the laboratory	Page 216
SCHOTT Instruments	www.schottinstruments.com – your online information platform	Page 220
	Innovative electrochemistry, innovative viscometry	Page 222
	Company's history	Page 224


Something to be measured against: the first pH meter with integrated PC

- ProLab 3000
2-channel pH / ISE
- ProLab 4000
2-channel pH / ISE / Conductivity




Operate the new ProLabs like a PC

- You have a 300 MHz processor and a storage volume of 64 MB.
- On a brilliant QVGA colour graphic display the operation is handled via a Windows menu structure. The display can be configured in the same manner as on a PC.
- Operate by using the numeric key pad or using an external keyboard or optionally use the mouse, which is part of the delivery scope.



Electrode and user identification

- When operating SCHOTT Instruments ID electrodes their individual data is transmitted to the meter for the calculation of the measurement values. Hence, optimum measurement accuracy is achieved.
- The automatic user recognition via transponder pendant plus password protection enables access control and allocation of measurement values, calibrations etc. for each user to meet GLP requirements.



More information on pages 16 to 21.



loLine, the pH electrode for the most challenging measuring tasks

- ▶ Iodine/iodide reference system with patented three-chambers system
- ▶ 100% free of metal ions

loLine

The new loLine product family from SCHOTT Instruments with iodine/iodide reference system for demanding measurements in the laboratory is 100% metal ion free. This makes the electrode perfectly suitable for measurements not allowing any metal ion contamination of the samples. Besides other convincing features, the loLine electrodes give the right answer to the demand for a suitable electrode for pH measurements in Tris buffers. loLine electrodes are the right choice for all pH measurements with highest quality requirements in the area of research and quality control of pharmacy, biotechnology or food industry.

loLine electrodes persuade with the following advantages compared to the conventional Ag/AgCl reference system:

- ▶ Patented three-chambers system with iodine storage reservoir in the reference electrode for prolonging the product lifecycle.
- ▶ Fast response times and most precise measurements with remarkably stable measuring values, especially with changing temperatures or a temperature deviating from the calibration through the metal ion free iodine/iodide system with a low temperature dependency of the reference potential.
- ▶ Free choice of bridge electrolyte and consequently optimal adaptation to each sample thanks to the double electrolyte system.

By using different membrane glass shapes and types with different glass shaft lengths and connection opportunities you always have the perfectly matching electrode for your corresponding application at hand. The comprehensive technical data of the loLine series are as follows:

pH application area	0 to 14
Temperature range	-5 up to 100 °C
Reference system	Iodine/iodide
Diaphragm	Platinum or ceramic
Glass membrane	Type: A or H glass; Shape: ball, spear or micro (cylinder)
Application area	Utmost precise and fast pH measurements in manifold media in the research and quality control areas of pharmacy, biotechnology or food industry


For further information please refer to pages 46/47 and the section starting at page 56.

TitroLine KF *trace*

The new coulometric KF-titrator from SCHOTT Instruments

- ▶ TitroLine KF *trace* – coulometric KF titrator for determining even smallest water contents.
- ▶ Easiest handling: Determination of concentration (titrant) of the titration solution is obsolete.

- ▶ Use without delay through pre-programmed methods.
- ▶ Large graphic display with online curve and clear instructions for navigation.
- ▶ GLP compliant documentation incl. curve print out.
- ▶ With optional password safety to avoid unwanted altering of methods.
- ▶ USB port for PC



Sample ID: 18
PPM
dy=50
Water:
Mean01/01:
SD:
RSD:
Weight:
Date:
Time:
Blank:
Meas.drift
Meas.time
Device ad
Serial n
Method



More information on pages 108 to 116.

UviLine 9100 and 9400

The new spectrophotometers from SCHOTT Instruments

Two completely new high-performance single beam photometers

- ▶ UviLine 9100 for measurements in the VIS range of 320 – 1100 nm
- ▶ UviLine 9400 for measurements in the UV-VIS range of 190 – 1100 nm



- ▶ Absorption and transmission measurements
- ▶ Concentration measurements of up to 8 standards
- ▶ Multi wavelength measurements
- ▶ Spectra adding with online graphic
- ▶ Kinetics measurements
- ▶ Extensive evaluation tools
- ▶ Large storage capacity for data and spectra
- ▶ USB-A and USB-B interfaces for communication

More information on pages 144 to 151.

Laboratory Meters of the Lab- and ProLab Series:

The new standard for electrochemical measurements

Electrodes and meters by SCHOTT Instruments

Electrodes from SCHOTT Instruments are used by professionals in laboratories all over the world. This comes as no surprise as we have been involved in the development, optimization and production for electrodes for more than 70 years. A know-how all around glass from which our customers can benefit. What once began with the patent on the pH electrode, has grown into an extensive range of products including several hundred electrodes designed to meet standard and special applications. Whether it is wastewater, varnish or wine - we have the appropriate electrode to meet our customers' requirements. Only the optimal interaction between electrode, calibration and meter allows precise measurements. Consequently the challenge was to develop meters perfectly matching our electrodes and buffer solutions.

The result

The new pH, ISE, conductometer and multi parameter measuring instruments from SCHOTT Instruments. These instruments again are setting standards for electrochemical measuring technology. Along with our successful electrodes and buffer solutions these instruments guarantee an ideal measurement result – fast, convenient and precise.



ProLab 2000



ProLab 3000




Lab 860

Performance Overview	Lab 850	Lab 860	Lab 870	Lab 960	Lab 970	ProLab 1000	ProLab 2000	ProLab 3000	ProLab 4000
Measuring Parameters & Special Functions page	8	8	8	10	10	12	14	16	20
pH	■	■	■			■	■	■	■
Dead Stop function						■			
2-channel pH measuring (galvanically separated)								■	■
16 pre-programmed pH buffer sets	■	■	■						
22 pre-programmed pH buffer sets						■	■	■	■
Automatic buffer recognition and display	■	■	■			■	■	■	■
pH calibration points max.	3	3	3			5	5	5	5
VariCal: manual calibration with selected buffer								■	■
mV	■	■	■			■	■	■	■
mV differential measurement								■	■
2-channel mV- (galvanically separated)								■	■
ISE							■	■	■
Conductivity				■	■		■		■
D.O.							■		
Temperature - simultaneous display	■	■	■	■	■	■	■	■	■
GLP and User Convenience									
Automatic recognition of ID sensors			■		■	■	■	■	■
Automatic user recognition with electronic ID card						■	■	■	■
Additional password entry with user recognition								■	■
CalClock - at a glance: sensor evaluation and calibration timer	■	■	■	■	■	■	■	■	■
Selectable calibration interval	■	■	■	■	■	■	■	■	■
Display of actual calibration data incl. date/time	■	■	■	■	■	■	■	■	■
Display of calibration history (10 data sets) incl. date/time								■	■
Measuring with stability control	■	■	■	■	■	■	■	■	■
Display accuracy of measuring value adjustable	■	■	■	■	■	■	■	■	■
Display	LCD	LCD	LCD	LCD	LCD	black/white graphic	black/white graphic	QVGA colour graphic	QVGA colour graphic
Windows user surface with optional mouse operation								■	■
Recorder function (display of measuring sequences over display)								■	■
Tactile response as well as optical and acoustical signal	■	■	■	■	■	■	■	■	■
Numerical keypad with 12 keys								■	■
Data storage		■		■		■	■	■	■
USB (slave) and RS232 interface		■	■	■	■	■	■	■	■
USB host interface: plug and play connection of USB hub, USB printer, USB memory, keyboard, mouse, USB stick								■	■
Quality and Service									
Complete scope of delivery: - instrument with cover, universal power supply unit and stand - set additionally with electrode and buffer	■	■	■	■	■	■	■	■	■
IQ and OQ documents available	■	■	■	■	■	■	■	■	■
Warranty 3 years	■	■	■	■	■	■	■	■	■

Measuring pH with the Lab 8xx series

pH Meter Lab 850

The new standard in the beginners' class



Quality saves time

All meters from the Lab 8xx series feature our advanced measuring algorithm specially adapted to SCHOTT Instruments electrodes. The result is a calibration and measurement in usual precision – but in much shorter time!

Keeping reliability in view: The CalClock...

This unique combination of sensor evaluation and calibration timer controls all calibration relevant settings. The sensor evaluation shows the electrode quality in increments. The quality criteria are slope, zero point and reaction time. The preset calibration interval is displayed as a countdown in 6 steps. When the time interval has elapsed, an optical and an acoustical signal is activated. The CalClock combines both functions at a glance.



CalClock

pH Meter Lab 860

Perfect communication for fast documentation

Direct data transmission

All instruments with data transmission function are provided with a conventional RS 232 plus a modern USB PC interface. No adapter is required.

USB and RS 232-
interfaces



pH Meter Lab 870

Laboratory practice improved – GLP with SCHOTT Instruments

- ▶ Sensor is automatically recognized

The quality controller's dream: A measuring instrument that automatically recognizes the electrode and the user! The new sensors with distinct identification send their specific data to the measuring instrument wirelessly. Wrong measurements are excluded.

Wireless sensor recognition



Special stand S4D



- ▶ **Reliable and precise measuring values**
due to perfect matching of the measuring instruments with our electrodes and buffer solutions.
- ▶ **Keeping reliability in view with „CalClock“**
sensor recognition and calibration timer combined.
- ▶ **Easy documentation complying with GLP**
through perfect communication via USB (slave) and RS 232 (except Lab 850).
- ▶ **Highest safety for measuring and calibrating**
through wireless sensor recognition (Lab 870):
ID electrodes and meters with automatic identification and data exchange.
- ▶ **Scope of delivery, set**
 - Measuring instrument
 - Electrode with integrated temperature sensor
 - Buffer solutions
 - Stand
 - Universal power supply unit
 - Cover

Advantages
Lab 8xx

Determining Conductivity: Lab 9xx

Conductivity Meter Lab 960

GLP documentation on highest level: the interface all-rounder

■ Dual interface included

Whether conventional via RS 232 or modern via USB – both interfaces are integrated in the Lab 960. No adapter is required.

■ Reliable documentation

All GLP compliant calibration protocols and up to 800 stored data sets with time and identification can be recalled.

```
14.03.2008 08:53:54
Lab 960 02320025

CALIBRATION COND
Cal Time : 14.03.2008 08:22:14
Cal Interval: 180d
Cal Std.: 0.01 mol/l KCL
          40.0 °C
Conduct./Tref25: 1413µS/cm
Cell Const : 0.650 1/cm
Probe : +++
```

Calibration protocol

USB and RS 232
interfaces



Scope of delivery, set

Conductivity Meter Lab 970

Laboratory practice at its best – GLP with SCHOTT Instruments

Electrode is automatically recognized

Ideally suited for measurements according to USP 28. A measuring instrument that automatically recognizes the electrode! The latest microtechnology enables storing of the calibration data in the sensor. For each measurement the sensor logs in with its ID (type and serial no.) and by sending its calibration data assures that only the sensor specific data are used to calculate the measured value.

Wireless sensor recognition



Direct data transmission

The Lab 970 sends the data directly via USB or RS 232 to PC (no storing function). The instrument is therefore very easy to handle and can be integrated into automated systems such as i.e. LIMS. Extensive administrator levels are not necessary.

Service taken serious

For the Lab 970 also the complete IQ and OQ documentation is available. The service from SCHOTT Instruments is pleased to assist you.

Reliable and precise measuring values

due to perfect matching of the Lab instruments with our conductivity cells and conductivity test solutions.

Fast and easy documentation complying with GLP

through perfect communication via USB (slave) and RS 232

Complete support

including IQ and OQ documentation and qualified service.

Highest safety for measuring and calibrating

through wireless sensor recognition (Lab 970): ID electrodes and meters with automatic identification and data exchange.

Scope of delivery, set

- Measuring instrument
- Conductivity cell
- Conductivity testing solution
- Stand
- Universal power supply unit
- Cover

Advantages
Lab 9xx

ProLab 1000: pH for the professional

pH Meter ProLab 1000:

Superior technique in a classy metal body

Precise measurements

With a measuring range of -2.000 up to +20.000 pH and an accuracy of 0.003 pH the ProLab is especially suited for demanding measurements in research and process control.

Reproducible results

The automatic calibration at up to 5 points and the automatic buffer recognition from more than 16 selectable buffer sets guarantee trustworthy measurement results on highest level.

User recognition with electronic ID card



Utmost reliability

The wireless sensor recognition and user identification allows to assign a user name to each measured value. ProLab 1000 also complies with further requirements, e.g. user levels allowing to allocate administrator rights or user identification via an electronic ID card.

Glass meets metal


Housed in a solid metal body, the big and brilliant graphic display is protected by a one-piece glass slide. The ProLab 1000's weight of approx. 2.5 kg guarantees a firm standing. The S4D stand can be attached without tools either to the left or the right side of the instrument to make up an immovable unit.



Direct data transmission

Whether conventional via RS 232 or modern via USB – both interfaces to the PC are integrated in the ProLab 1000. No adapter is required.

The integrated dead stop function increases the operation area also for manual titration tasks. (wine industry, food industry).




Wireless sensor recognition



Documentation complying with GLP

The calibration protocols comply with GLP as well as date and time and identification of up to 1500 data entries can be recalled. The storage is according to 21 CFR Part 11 in user levels with user identifications.



Highest safety level with measurements and calibration through:

- **automized user identification**
user recognition with electronic ID card through transponder technology.
- **wireless sensor recognition**
ID electrodes and measuring instrument with automatic identification and data exchange.

High-precision and reliable measuring values through perfect matching with the measuring system.

Fast and easy documentation complying with GLP through perfect communication via USB (slave) and RS 232

Scope of delivery, set

- Measuring instrument
- Electrode with integrated temperature sensor
- Buffer solutions
- Stand
- Universal power supply unit
- Cover

Advantages
ProLab 1000

ProLab 2000: the multi talented ...

Multi-parameter instrument ProLab 2000:


One for all

Versatile

The ProLab 2000 is a true all-rounder. pH, ISE, conductivity and D.O. measuring is integrated in one single instrument. Up to 4 parameters can be measured simultaneously and displayed in the graphic display. The parameters are selected via menu.

Professional

With a measuring range from -2.000 to +20.000 pH and an accuracy of 0.003 pH, the ProLab 2000 is especially suited for manifold applications in research and development. The automatic pH and ISE calibration at up to 3 points guarantees trustworthy measurements at the highest stage. The galvanic D.O. sensor and multi-pole conductivity measuring cells widen the measuring functions for most applications.




User recognition with electronic ID card



▲ Innovative - Sensors are recognized automatically

The latest microtechnology enables storing of calibration data in the individual sensor. Upon connection the sensor logs in with its ID (type and serial no.), and by transmitting its specific data assures a correct measurement. The sensors are identified unmistakably via a shielded, short-range radio connection. No input is needed. Even two sensors side by side can be identified unequivocally.



*pH combination electrode
A 161 ... ID*



*Combined conductivity cell
and D.O. sensor
LFOX 1400 ID*



▲ Highest safety level with measurements and calibration through:

- **automized user identification**
with electronic ID card through transponder technology.
- **wireless sensor recognition**
ID electrodes and measuring instrument with automatic identification and data exchange.

▲ High-precision and reliable measurement of pH, ISE, conductivity and oxygen with sensors from SCHOTT Instruments perfectly matching the instrument

▲ Safety at a glance with „CalClock“

Sensor evaluation and calibration timer combined.

▲ Scope of delivery, set

- Measuring instrument
- Electrodes with integrated temperature sensor
- Buffer solutions
- Stand
- Universal power supply unit
- Cover

Advantages
ProLab 2000

ProLab 3000: measure pH as at a PC





High end pH/ION meter ProLab 3000:

The first pH meter with integrated PC

Operate it like „Windows“

The menu structure which is similar to windows can be operated just like a PC via the menu keys on the instrument or by using the mouse, which is part of the delivery scope. For example the entry of text and numbers for inserting file names can be carried out either by using the numeric keypad or an external keyboard. Mouse and keyboard can be used at the same time, as the USB (host) interface can be extended by a hub.



Flexible display of measurements on a brilliant colour graphic display

The colour intensive QVGA colour graphic display (320 x 240 pixels) is backlit for high resolution and gives good readability even from a side distance of 2 to 3 meters. The measurement display can easily be switched to full display of a measuring channel or multiple display of different measuring parameters (pH, mV, ISE). Another option is the completely flexible setting of the "recorder" display regarding the measurement sequences of all parameters against time.

Plug and Play


Easy connection of peripheral devices with the automatic recognition – no configuration required. The integrated USB host, USB slave and RS232 interfaces enable the instrument to communicate with mouse, printer and barcode reader – when interconnecting a hub to the USB host even a parallel operation is possible.



ProLab 3000: measure pH as at a PC

▲ Differential measurement made easy

Two galvanically separated pH entries allow simultaneous measurement of both pH sensors, even in the same vessel without causing any perturbation – with SCHOTT Instruments' ID electrodes providing automatic sensor recognition even unequivocally. On every pH channel a differential measurement against a reference value is possible.



GLP thought out

▶ User recognition with electronic ID card incl. password entry

The automatic user recognition via transponder pendant enables access control and allocation of measurement values, calibrations etc. for each user. Hence it is automatically documented which user did what at a certain time. The security guaranteed through the access control is furthermore increased by the password protection.

▶ Sensor recognition – wireless and automatic

The new ID electrodes send their specific data wireless to the ProLab 3000 and 4000. It is therefore guaranteed that each electrode always uses the correct calibration and erroneous measurements are excluded – unequivocal even with 2-channel measurements!

The ProLab user specific display automatically adapts to the ID sensor and activates only the necessary and admissible operation structures - this guarantees a higher usage convenience through increased transparency.



▶ Highest safety level with measurements and calibration through:

- **automized user identification**
with electronic ID card through transponder technology with password entry
- **wireless sensor recognition**
ID electrodes and measuring instrument with automatic identification and data exchange

▶ Measuring pH, mV, ISE – high precision with many special functions:

- 2-channel pH/mV measurement (galvanically separated)
- differential measurement
- professional ISE measuring with various addition/subtraction procedures

▶ Operation via mouse or keyboard as at a PC familiar menu structures and clear menu navigation

▶ Plug and play by state-of-the-art technology

▶ Scope of delivery, set

- Measuring instrument (incl. mouse)
- Electrode with integrated temperature sensor
- Buffer solutions
- Stand
- Universal power supply unit
- Cover

Advantages
ProLab 3000

ProLab 4000: and Cond on top...

ProLab 4000: High end pH/ION/Conductivity meter

pH, ISE and conductivity measuring on highest level

Equivalent to ProLab 3000 – added with another feature ...

The ProLab 4000 offers the high-quality measurement technology known with ProLab 3000 and is furthermore completed with conductivity measurement for highest standards: A measuring range from 0.000 $\mu\text{S}/\text{cm}$... 2000 mS/cm , TDS and salinity measurement as well as various functions for temperature compensation and setting of the cell constant form the standard for a measuring instrument of this performance class. ProLab 4000 offers even more ...

Special functions for determining the dependency of conductivity regarding temperature or concentration



The conductivity of aqueous solutions is subject to temperature and concentration of dissolved substances. For comparing measurement values, which were determined at different temperatures, all values have to be recalculated onto the same reference temperature. Both common procedures linear and

non-linear (acc. to EN 27888) compensation should however only be applied with diluted measuring solutions, as otherwise the concentration dependency goes unnoticed. Furthermore, the measurement temperature range should lay within $\pm 10\text{K}$ of the reference temperature. The ProLab 4000 eliminates these restrictions thanks to special compensation procedures. These methods are featured for example by:

- Usage of two pre-programmed temperature coefficients for HCl, NaOH, NaCl and KCl for a temperature range of 0...40°C
- Possibility of entering literature values for two temperature coefficients for additional solutions.
- Determination of temperature coefficients by:
 - Setting the temperature range and intervals
 - Measurement of one or multiple solutions of known and even unknown concentrations (through equidistant dilution).

ProLab 4000 is therefore qualified for operation in science and for monitoring industrial processes. The instrument enables high precision conductivity measurements in a wide temperature and concentration range.





- ▶ **pH, mV, ISE measuring – same as ProLab 3000 ...**
- ▶ **Measuring conductivity on highest level:**
Determining the dependency of temperature and concentration through
 - **using a temperature coefficient**
(stored value or entry of own literature values)
 - **Self-determination of temperature coefficients** in various standards of known or unknown concentrations in a self-defined temperature range and span
- ▶ **Highest safety level with measurements and calibration through:**
 - **automized user identification**
with electronic ID card through transponder technology with password entry
 - **wireless sensor recognition**
ID electrodes and measuring instrument with automatic identification and data exchange
- ▶ **Operation via mouse or keyboard as at a PC**
familiar menu structures and clear menu navigation
- ▶ **Plug and play by state-of-the-art technology**
- ▶ **Scope of delivery, set**
 - Measuring instrument (incl. mouse)
 - Electrodes with integrated temperature sensor
 - Buffer- and conductivity test solutions
 - Stand
 - Universal power supply unit
 - Cover

Advantages
ProLab 4000

Performance black on white ...

Measuring technology in detail...	Lab 850	Lab 860	Lab 870	Lab 960
<i>page</i>	<i>p. 8/9</i>	<i>p. 8/9</i>	<i>p. 8/9</i>	<i>p. 10/11</i>
pH measurement	■	■	■	
Range / Accuracy	-2.000 ... +19.999 pH -2.00 ... +19.99 pH	-2.000 ... +19.999 pH -2.00 ... +19.99 pH	-2.000 ... +19.999 pH -2.00 ... +19.99 pH	
Accuracy (for each measuring area) (±1 digit)	± 0.005 pH ± 0.01 pH	± 0.005 pH ± 0.01 pH	± 0.005 pH ± 0.01 pH	
Calibration: pre-programmed pH buffer sets	16	16	16	
Automatic buffer recognition and display	■	■	■	
pH calibration points max.	3	3	3	
VariCal: manual calibration with selectable buffers	–	–	–	
Dead stop function	–	–	–	
2-channel pH measurement (galvanically separated)	–	–	–	
mV measurement	■	■	■	
Range / Accuracy	-999.9 ... +999.9 mV -1999 ... +1999 mV	-999.9 ... +999.9 mV -1999 ... +1999 mV	-999.9 ... +999.9 mV -1999 ... +1999 mV	
Accuracy (for each measuring area) (±1 digit)	± 0.3 mV ± 1 mV	± 0.3 mV ± 1 mV	± 0.3 mV ± 1 mV	
AutoRange function (can be switched off)	■	■	■	
mV differential measurement	–	–	–	
2-channel mV measurement (galvanically separated)	–	–	–	
ISE measurement				
Range / Accuracy				
Display in %, ppm, mg/kg, mol/l				
Two separate ISE channels (with dedicated separate temperature channels)				
Methods				
ISE calibration points				
Standard concentrations				
Conductivity measurement				■
Range / Accuracy				0.000 µS/cm...500 mS/cm
TDS measurement with factor 0.4 ... 1.0				■
Salinity measurement acc. to Natural Sea Water Scale (UNESCO 1966b)				■
Accuracy in % from measuring value (± 1 digit)				0.5
Calibrated cell constant 0.450 ... 0.500 cm ⁻¹ ; 0.585 ... 0.715 cm ⁻¹ ; 0.800 ... 1.200 cm ⁻¹ (calibration with control standard) d 0.01 mol KCl:				■
Adjustable cell constant 0.250 ... 2.500 cm ⁻¹ and 0.090 ... 0.110 cm ⁻¹				■
Fixed cell constant 0.010 cm ⁻¹				■
Temperature compensation nLF / Lin (0.001 ... 3.000 %/K) / selectable				■
Temperature compensation purity water				■
Pre-programmed temperature coefficients for HCl, NaOH, NaCl and KCl				
Determination of temperature coefficients for one or multiple standards and known or unknown concentrations at various temperatures				
Reference temperature 20°C or 25 °C selectable				■

... the technical data

Lab 970	ProLab 1000	ProLab 2000	ProLab 3000	ProLab 4000
p. 10/11	p. 12/13	p. 14/15	p. 16 – 19	p. 20/21
■	■	■	■	■
	-2.000 ... +20.000 pH -2.00 ... +20.00 pH -2.0 ... +20.0 pH	-2.000 ... +20.000 pH -2.00 ... +20.00 pH -2.0 ... +20.0 pH	-2.000 ... +20.000 pH -2.00 ... +20.00 pH -2.0 ... +20.0 pH	-2.000 ... +20.000 pH -2.00 ... +20.00 pH -2.0 ... +20.0 pH
	± 0.003 pH ± 0.01 pH	± 0.003 pH ± 0.01 pH	± 0.002 pH ± 0.01 pH	± 0.002 pH ± 0.01 pH
	22	22	22	22
■	■	■	■	■
	5	5	5	5
	–	–	■	■
■	–	–	–	–
–	–	–	■	■
■	■	■	■	■
	-1999.9 ... +1999.9 mV -1999 ... +1999 mV	-1999.9 ... +1999.9 mV -1999 ... +1999 mV	-2200.0 ... +2200.0 mV -2200 ... + 2200 mV	-2200.0 ... +2200.0 mV -2200 ... + 2200 mV
	± 0.2 mV ± 1 mV	± 0.2 mV ± 1 mV	± 0.1 mV ± 1 mV	± 0.1 mV ± 1 mV
■	■	■	–	–
–	–	–	■	■
–	–	–	■	■
		■	■	■
		0.000 ... 10000 mg/l	1.0E-40 ... 9.9E39 mg/l	1.0E-40 ... 9.9E39 mg/l
		–	■	■
		–	■	■
		–	Std. Add., Double Std. Add., Std. Sub., Sample Add., Sample Sub., Blank Add., Blank Corr., Ref. Measur.	Std. Add., Double Std. Add., Std. Sub., Sample Add., Sample Sub., Blank Add., Blank Corr., Ref. Measur.
		2 ... 3	2 ... 9	2 ... 9
		0.01 ... 10 000 mg/l in 19 selectable concentrations	1.00E-30 ... 1.00E30 mg/l can be inserted	1.00E-30 ... 1.00E30 mg/l can be inserted
■		■		■
0.000 µS/cm...500 mS/cm		0.000 µS/cm...2000 mS/cm		0.000 µS/cm...2000 mS/cm
■		■		■
■		■		■
0.5		0.5		0.5
■		■		■
■		■		■
■		■		■
■		■		■
■		■		■
				■
				■
■		■		■

The technical data continued...

Measuring technology in detail...	Lab 850	Lab 860	Lab 870	Lab 960
<i>page</i>	<i>p. 8/9</i>	<i>p. 8/9</i>	<i>p. 8/9</i>	<i>p. 10/11</i>
D.O. measurement (O₂ dissolved)				
Range / Accuracy:				
O ₂ concentration				
O ₂ saturation				
O ₂ partial pressure				
Accuracy in % from measuring value (±1 digit) at an ambient temperature of 5...30 °C				
Salinity correction				
Calibration in calibration vessel with water vapor-saturated air				
Temperature measurement	■	■	■	■
Range / Accuracy	-5.0 ... +120.0 °C	-5.0 ... +120.0 °C	-5.0 ... +120.0 °C	-5.0 ... +120.0 °C
Accuracy (±1 digit)	± 0.1 °C	± 0.1 °C	± 0.1 °C	± 0.1 °C
Two separated temperature channels	–	–	–	–
Selectable °C / °F (Fahrenheit)	■	■	■	■
Automatic switch-over to manual temperature input when no temperature sensor is connected	■	■	■	■
Design & Quality				
Display	LCD 75 x 60 mm	LCD 75 x 60 mm	LCD 75 x 60 mm	LCD 75 x 60 mm
Contrast adjustment via menu	–	–	–	–
Glass vision panel	–	–	–	–
Vision panel integrated in keyboard plastic foil	■	■	■	■
Measuring value storage (manual/automatic)	–	800 data sets, storage intervals from 5 s ... 60 min	–	800 data sets, storage intervals from 5 s ... 60 min
USB (slave) and RS232 interface		■	■	■
USB host interface: plug and play connection of USB hub, USB printer, USB memory, keyboard, mouse, USB stick				
Lower casing	plastic	plastic	plastic	plastic
Plastic foil keypad (polyester) with tactile response	■	■	■	■
Power supply: external universal power supply unit (medical approval) with country specific primary adapters, (primary: 100-240V, 50/60 Hz, secondary: 9V=1,5A)	■	■	■	■
built-in real-time clock (processor solution) battery buffered, exchangeable battery	■	■	■	■
Battery operation possible (4 mignon)	■	■	■	■
Battery switch-off automatic (adjustable 10 min ... 24 h, default 1 h, cannot be switched off)	■	■	■	■
Dimensions (W x H x D mm)	240 x 190 x 80	240 x 190 x 80	240 x 190 x 80	240 x 190 x 80
Weight	approx. 1.0 kg	approx. 1.0 kg	approx. 1.0 kg	approx. 1.0 kg
Compliance	CE, cETLus	CE, cETLus	CE, cETLus	CE, cETLus
Safety	protection class III, EG guidelines 73/23, EN 61010-1: 2001	protection class III, EG guidelines 73/23, EN 61010-1: 2001	protection class III, EG guidelines 73/23, EN 61010-1: 2001	protection class III, EG guidelines 73/23, EN 61010-1: 2001
Climate class	2 (VDI/VDE 3540)	2 (VDI/VDE 3540)	2 (VDI/VDE 3540)	2 (VDI/VDE 3540)
Complete delivery scope: - Instrument with cover, power supply unit and stand - set additionally with electrode and buffer	■	■	■	■
IQ and OQ documents available	■	■	■	■
Warranty 3 years	■	■	■	■

Lab 970	ProLab 1000	ProLab 2000	ProLab 3000	ProLab 4000
p. 10/11	p. 12/13	p. 14/15	p. 16 – 19	p. 20/21
		■		
		0...20.00 mg/l / 0.01 mg/l 0...90.0 mg/l / 0.1		
		0...200.0 % / 0.1 % 0...600 % / 1 %		
		0...200.0 mbar / 0.1 mbar 0...1250 mbar / 1 mbar		
		0.5		
		■		
		■		
■	■	■	■	■
-5.0 ... +120.0 °C	-10.0 ... +120.0 °C	-10.0 ... +120.0 °C	-35.0 ... +150.0 °C	-35.0 ... +150.0 °C
± 0.1 °C	± 0.1 °C	± 0.1 °C	± 0.1 °C	± 0.1 °C
–	–	–	■	■
■	■	■	■	■
■	■	■	■	■
LCD 75 x 60 mm	black & white graphic 120 x 90 mm with lighting	black & white graphic 120 x 90 mm with lighting	QVGA colour graphic display 120 x 90 mm with lighting	QVGA colour graphic display 120 x 90 mm with lighting
–	■	■	–	–
–	■	■	■	■
■	–	–	–	–
–	1500 data sets, storage intervals from 1 s ... 60 min	1500 data sets, storage intervals from 1 s ... 60 min	> 10000 data sets, storage intervals from 1 s ... 60 min	> 10000 data sets, storage intervals from 1 s ... 60 min
■	■	■	■	■
			■	■
plastic	metal diecast	metal diecast	metal diecast	metal diecast
■	■	■	■	■
■	■	■	■	■
■	■	■	■	■
■	–	–	–	–
■	–	–	–	–
240 x 190 x 80	240 x 280 x 70	240 x 280 x 70	240 x 280 x 70	240 x 280 x 70
approx. 1.0 kg	approx. 2.5 kg	approx. 2.5 kg	approx. 2.5 kg	approx. 2.5 kg
CE, cETLus	CE, cETLus	CE, cETLus	CE, cETLus	CE, cETLus
protection class III, EG guidelines 73/23, EN 61010-1: 2001	protection class III, EG guidelines 73/23, EN 61010-1: 2001	protection class III, EG guidelines 73/23, EN 61010-1: 2001	protection class III, EG guidelines 73/23, EN 61010-1: 2001	protection class III, EG guidelines 73/23, EN 61010-1: 2001
2 (VDI/VDE 3540)	2 (VDI/VDE 3540)	2 (VDI/VDE 3540)	2 (VDI/VDE 3540)	2 (VDI/VDE 3540)
■	■	■	■	■
■	■	■	■	■
■	■	■	■	■

Order information

Type no.	Order no.	Product	Description
Lab series			
Lab 850	285201300	Laboratory pH Meter	Measuring parameters pH, mV, temp., microprocessor, DIN 19262 connection. Including cover Z 880, stand S4D Z 865 and power supply Z 850.
Lab 850 Set	285201310	Laboratory pH Meter	Measuring parameters pH, mV, temp., microprocessor, DIN 19262 connection. Including cover Z 880, stand S4D Z 865, power supply Z 850, pH-temp. combination electrode BlueLine 14 pH, calibration solutions (DIN).
Lab 850 BNC	285201360	Laboratory pH Meter	Measuring parameters pH, mV, temp., microprocessor, BNC connection. Including cover Z 880, stand S4D Z 865 and power supply Z 850.
Lab 850 BNC Set	285201370	Laboratory pH Meter	Measuring parameters pH, mV, temp., microprocessor, BNC connection. Including cover Z 880, stand S4D Z 865, power supply Z 850, pH-temp. combination electrode BlueLine 15 pH, calibration solutions (DIN).
Lab 860	285201320	Laboratory pH Meter	Measuring parameters pH, mV, temp., microprocessor, RS 232 C and USB (slave) interface, data storage for 800 data sets, GLP conform, DIN 19262 connection. Including cover Z 880, stand S4D Z 865 and power supply Z 850.
Lab 860 Set	285201330	Laboratory pH Meter	Measuring parameters pH, mV, temp., microprocessor, RS 232 C and USB (slave) interface, storage for 800 data sets, GLP conform, DIN 19262 connection. Incl. cover Z 880, stand S4D Z 865, power supply Z 850, pH-temp. combination electrode BlueLine 14 pH, calibration solutions (DIN).
Lab 860 BNC	285201380	Laboratory pH Meter	Measuring parameters pH, mV, temp., microprocessor, RS 232 C and USB (slave) interface, data storage for 800 data sets, GLP conform, BNC connection. Including cover Z 880, stand S4D Z 865 and power supply Z 850.
Lab 860 BNC Set	285201390	Laboratory pH Meter	Measuring parameters pH, mV, temp., microprocessor, RS 232 C and USB (slave) interface, storage for 800 data sets, GLP conform, BNC connection. Including cover Z 880, stand S4D Z 865, power supply Z 850, pH-temp. combination electrode BlueLine 15 pH, calibration solutions (DIN).
Lab 870	285201340	Laboratory pH Meter	Electrode recognition. Measuring parameters pH, mV, temp., microprocessor, RS 232 C and USB (slave) interface, GLP conform, DIN 19262 connection. Including cover Z 880, stand S4D Z 865 and power supply Z 850.
Lab 870 Set	285201350	Laboratory pH Meter	Electrode recognition. Measuring parameters pH, mV, temp., microprocessor, RS 232 C and USB (slave) interface, GLP conform, DIN 19262 connection. Including cover Z 880, stand S4D Z 865, power supply Z 850, pH-temp. combination electrode BlueLine 14 pH ID, calibration solutions (DIN).
Lab 870 BNC	285201400	Laboratory pH Meter	Electrode recognition. Measuring parameters pH, mV, temp., microprocessor, RS 232 C and USB (slave) interface, GLP conform, BNC connection. Including cover Z 880, stand S4D Z 865 and power supply Z 850.
Lab 870 BNC Set	285201410	Laboratory pH Meter	Electrode recognition. Measuring parameters pH, mV, temp., microprocessor, RS 232 C and USB (slave) interface, GLP conform, BNC connection. Including cover Z 880, stand S4D Z 865, power supply Z 850, pH-temp. combination electrode BlueLine 15 pH ID, calibration solutions (DIN).
Lab 960	285201420	Laboratory Conductivity Meter	Measuring ranges 0.000 µS/cm...500 mS/cm, salinity, total dissolved solids (TDS), temp., RS 232 C and USB (slave) interface, microprocessor, data storage for 800 data sets, GLP conform. Including cover Z 880, stand S4D Z 865 and power supply Z 850.
Lab 960 Set	285201430	Laboratory Conductivity Meter	Measuring ranges 0.000 µS/cm...500 mS/cm, salinity, total dissolved solids (TDS), temp., RS 232 C and USB (slave) interface, microprocessor, data storage for 800 data sets, GLP conform. Including cover Z 880, stand S4D Z 865, power supply Z 850, conductivity cell LF 413 T and conductivity testing solution.
Lab 970	285201440	Laboratory Conductivity Meter	Sensor recognition. Measuring ranges 0.000 µS/cm...500 mS/cm, salinity, total dissolved solids (TDS), temp., RS 232 C and USB (slave) interface, microprocessor, GLP conform. Including cover Z 880, stand S4D Z 865 and power supply Z 850.
Lab 970 Set	285201450	Laboratory Conductivity Meter	Sensor recognition. Measuring ranges 0.000 µS/cm...500 mS/cm, salinity, total dissolved solids (TDS), temp., RS 232 C and USB (slave) interface, microprocessor, GLP conform. Including cover Z 880, stand S4D Z 865, power supply Z 850, conductivity cell LF 413 T ID and conductivity testing solution.
ProLab series			
ProLab 1000	285201700	Digital laboratory pH Meter	Electrode recognition and user identification. Measuring parameters pH, mV, temp., microprocessor, RS 232 C and USB (slave) interface, GLP conform, DIN 19262 connection. Including cover Z 881, stand S4D Z 865 and power supply Z 850.
ProLab 1000 Set	285201710	Digital laboratory pH Meter	Electrode recognition and user identification. Measuring parameters pH, mV, temp., microprocessor, RS 232 C and USB (slave) interface, GLP conform, DIN 19262 connection. Including cover Z 881, stand S4D Z 865, power supply Z 850, pH-temp. combination electrode A 161 1M-DIN-ID, calibration solutions (DIN).
ProLab 1000 BNC	285201720	Digital laboratory pH Meter	Electrode recognition and user identification. Measuring parameters pH, mV, temp., microprocessor, RS 232 C and USB (slave) interface, GLP conform, BNC connection. Including cover Z 881, stand S4D Z 865 and power supply Z 850.
ProLab 1000 BNC Set	285201730	Digital laboratory pH Meter	Electrode recognition and user identification. Measuring parameters pH, mV, temp., microprocessor, RS 232 C and USB (slave) interface, GLP conform, BNC connection. Including cover Z 881, stand S4D Z 865, power supply Z 850, pH-temp. combination electrode A 161 1M-BNC-ID, calibration solutions (DIN).
ProLab 2000	285201740	Digital laboratory multi Meter	Electrode recognition and user identification. Measuring parameters pH, mV, ISE, conductivity, D.O. and temp., microprocessor, RS 232 C and USB (slave) interface, GLP conform, DIN 19262 connection. Including cover Z 881, stand S4D Z 865 and power supply Z 850.
ProLab 2000 Set	285201750	Digital laboratory multi Meter	Electrode recognition and user identification. Measuring parameters pH, mV, ISE, conductivity, D.O. and temp., microprocessor, RS 232 C and USB (slave) interface, GLP conform, DIN 19262 connection. Including cover Z 881, stand S4D Z 865, power supply Z 850, pH-temp. combination electrode A 161 1M-DIN-ID, combined conductivity and D.O. sensor LFOX 1400 ID, calibration solutions (DIN), conductivity testing solutions.

Type no.	Order no.	Product	Description
ProLab 2000 BNC	285201760	Digital laboratory multi Meter	Electrode recognition and user identification. Measuring parameters pH, mV, ISE, conductivity, D.O. and temp., microprocessor, RS 232 C and USB (slave) interface, GLP conform, BNC connection. Including cover Z 881, stand S4D Z 865 and power supply Z 850.
ProLab 2000 BNC Set	285201770	Digital laboratory multi Meter	Electrode recognition and user identification. Measuring parameters pH, mV, ISE, conductivity, D.O. and temp., microprocessor, RS 232 C and USB (slave) interface, GLP conform, BNC connection. Including cover Z 881, stand S4D Z 865, power supply Z 850, pH-temp. combination electrode A 161 1M-BNC-ID, combined conductivity and D.O. sensor LFOX 1400 ID, calibration solutions (DIN), conductivity testing solutions.
ProLab 3000	285203600	Digital laboratory pH Meter	Electrode recognition and user identification. QVGA colour display. Menu based operation. Recorder function. Measuring parameters: double-pH, mV, Temp, ISE. RS232, USB host + USB slave interfaces. DIN connection. Incl. Z880, Z865 + Z850.
ProLab 3000 Set	285203610	Digital laboratory pH Meter	Electrode recognition and user identification. QVGA colour display. Menu based operation. Recorder function. Measuring parameters: conductivity + double-pH, mV, Temp, ISE. RS232, USB host + USB slave interfaces. DIN connection. Incl. Z880, Z865 + Z850. Z880, Z865, Z850, IL-pHT-A170MF-DIN-N, DIN buffers.
ProLab 3000 BNC	285203620	Digital laboratory pH Meter	Electrode recognition and user identification. QVGA colour display. Menu based operation. Recorder function. Measuring parameters: double-pH, mV, Temp, ISE. RS232, USB host + USB slave interfaces. BNC connection. Incl. Z880, Z865 + Z850..
ProLab 3000 BNC Set	285203630	Digital laboratory pH Meter	Electrode recognition and user identification. QVGA colour display. Menu based operation. Recorder function. Measuring parameters: double-pH, mV, Temp, ISE. RS232, USB host + USB slave interfaces. BNC connection. Incl. Z880, Z865, Z850, IL-pHT-A170MF-BNC-N, DIN buffers.
ProLab 4000	285203640	Digital laboratory multi Meter	Electrode recognition and user identification. QVGA colour display. Menu based operation. Recorder function. Measuring parameters: conductivity + double-pH, mV, Temp, ISE. RS232, USB host + USB slave interfaces. DIN connection. Incl. Z880, Z865 + Z850.
ProLab 4000 Set	285203650	Digital laboratory multi Meter	Electrode recognition and user identification. QVGA colour display. Menu based operation. Recorder function. Measuring parameters: conductivity + double-pH, mV, Temp, ISE. RS232, USB host + USB slave interfaces. DIN connection. Incl. Z880, Z865, Z850, IL-pHT-A170MF-DIN-N, LF413TID, DIN buffer, conductivity testing solution.
ProLab 4000 BNC	285203660	Digital laboratory multi Meter	Electrode recognition and user identification. QVGA colour display. Menu based operation. Recorder function. Measuring parameters: conductivity + double-pH, mV, Temp, ISE. RS232, USB host + USB slave interfaces. BNC connection. Incl. Z880, Z865 + Z850.
ProLab 4000 BNC Set	285203670	Digital laboratory multi Meter	Electrode recognition and user identification. QVGA colour display. Menu based operation. Recorder function. Measuring parameters: conductivity + double-pH, mV, Temp, ISE. RS232, USB host + USB slave interfaces. BNC connection. Incl. Z880, Z865, Z850, IL-pHT-A170MF-BNC-N, LF413TID, DIN buffer, conductivity testing solution.
Accessories			
Logbook Lab 850	285201800	Logbook	for Lab 850 (DIN and BNC) incl. review by SCHOTT Instruments after resending the filled in documents.
Logbook Lab 860	285201810	Logbook	for Lab 860 (DIN and BNC) incl. review by SCHOTT Instruments after resending the filled in documents.
Logbook Lab 870	285201820	Logbook	for Lab 870 (DIN and BNC) incl. review by SCHOTT Instruments after resending the filled in documents.
Logbook Lab 960	285201840	Logbook	for Lab 960 incl. review by SCHOTT Instruments after resending the filled in documents.
Logbook Lab 970	285201850	Logbook	for Lab 970 incl. review by SCHOTT Instruments after resending the filled in documents.
Logbook ProLab 1000	285201830	Logbook	for ProLab 1000 (DIN and BNC) incl. review by SCHOTT Instruments after resending the filled in documents.
Logbook ProLab 2000	285201860	Logbook	for ProLab 2000 (DIN and BNC) incl. review by SCHOTT Instruments after resending the filled in documents.
Logbook ProLab 3000	285203680	Logbook	for ProLab 3000 (DIN and BNC) incl. review by SCHOTT Instruments after resending the filled in documents.
Logbook ProLab 4000	285203690	Logbook	for ProLab 4000 (DIN and BNC) incl. review by SCHOTT Instruments after resending the filled in documents.
Z 390	285201560	Cable for connection to PC	RS 232 6pole cable for connection to PC for Lab 860, Lab 870, Lab 960 and Lab 970 as well as for all instruments being part the ProLab series
Z 396	285201580	Software	Software for documentation for Lab 860, Lab 870, Lab 960, Lab 970, handylab 12 as well as for all instruments being part of the ProLab series
Z 850	285204889	Power supply	Universal power supply unit, 230 and 120 V for the Lab- and ProLab-meter family
Z 865	285201520	Stand set S4D	Stand set S4D, including arm and electrode holder for the Lab- and ProLab-meter family
Z 875	285201540	USB cable	for Lab 860, Lab 870, Lab 960 and Lab 970 as well as for all instruments of the ProLab-meter series with USB (slave)
Z 876	285201890	Transponder	User recognition transponder for ProLab instruments
Z 880	285201550	Cover	for the Lab-meter family
Z 881	285201880	Cover	for the ProLab-meter family
Z 890	285203700	Universal paper printer	Star SP-712 (9-matrix printer). Easy paper load. Serial interface. Dimensions: 160 (width) x 245 (depth) x 152 (height) mm. Weight 2.96 kg. Integrated power supply..
Z 891	285203710	ink ribbon (black)	for printer Z 890. Product life cycle: 3 million charachters.
Z 892	285203720	printer paper role	for printer Z 890, 1 piece. Universal paper. Width 76 mm. External diameter 80 mm, inner core 12 mm.
Z 893	285203730	Connecting cable	for connection of printer Z 890 to the ProLab meters (except Lab 850) and ProLab-meter family.
	285209081	Manufacturer certificate	for pH-Meter, conductivity meter and pH/mV simulators from SCHOTT Instruments

handylab – pocket-size all-rounders ...

handylab – the mobile, multi-functional mini-laboratories

The latest handylab generation is available in seven different models – all with a new modern look and expanded features. These compact, battery powered, pocket-size meters were specially designed for field work.

They are available as a set in a practical carrying case with the respective combination electrode and all of the requisite accessories, thus providing the user with a high performance mini-laboratory.

The handylab pH/LF12 multi-parameter portable meter can be used to determine pH values, redox potential, conductivity and temperature. The handylab multi12, which is a real all-rounder, can also be used to measure oxygen concentrations.

The measurement parameters pH, mV and °C make the handylab pH11 and pH12 pH meters suitable for a variety of uses. The fully automatic one to three point calibration including preprogrammed DIN or technical buffers makes practical work substantially easier.

The handylab OX 12 oxygen meter automatically takes influence variables such as temperature and air pressure into account during measurements. The influence of a higher salinity level on the oxygen measurement can also be corrected by entering the salinity value, after determining it using a conductivity meter.

All of the handylab 12 models have a data memory, which means that measurements can be recorded manually or automatically using a timer control, and can be evaluated at a later time. In addition, they have a serial interface, and an optional power pack is available for utilization at stationary measurement locations.

Features and applicability of the portable handylab pH meters and conductivity meters

handylab	pH 11	pH 12	LF 11	LF 12	OX12	pH/LF 12	Multi 12
pH	■	■	-	-	-	■	■
ORP	■	■	-	-	-	■	■
Temperature	■	■	■	■	■	■	■
Conductivity	-	-	■	■	-	■	■
Dissolved Oxygen	-	-	-	-	■	-	■
AutoRead	■	■	■	■	■	■	■
Battery operation	■	■	■	■	■	■	■
Mains power connection (power pack optional)	-	■	-	■	■	■	■
Data memory	-	-	-	■	■	■	■
RS 232	-	■	-	■	■	■	■

All of the handylab pH meters and conductivity meters are also available as a complete, cost-effective set in a carrying case.

... for measuring pH values and redox potential,
conductivity and dissolved oxygen

Contents

Portable pH meters with GLP functions handylab pH 11 and pH 12	Page 30
Technical data handylab pH 11, handylab pH 12	Page 31
Portable conductivity meters with GLP functions handylab LF 11 and LF 12	Page 32
Technical data handylab LF 11, handylab LF 12	Page 33
Portable oxygen meter with GLP functions handylab OX12	Page 34
Technical data 9009/61 O ₂ sensor	Page 34
Technical data handylab OX12	Page 35
Multi-parameter portable meters with GLP functions handylab pH/LF 12 and handylab multi 12	Page 36
Technical data handylab pH/LF 12, handylab multi 12	Page 37
Order overview handylab pH meters, conductivity meters and oxygen meters	Page 38
Order overview handylab multi-parameter portable meters	Page 39



Portable pH meters with GLP functions

handylab pH 11 and handylab pH 12

These pocket-size meters in shock-proof, water-tight casings are ideally suited for field work.

Measurement parameters

The pH, mV and °C measurement parameters mean that the pocket pH meters from SCHOTT Instruments have a variety of uses.

Measurement memory and interface

In comparison with the handylab pH 11, the handylab pH 12 additionally has a data memory, which makes it possible to save measurements manually or automatically using a timer control, and then evaluate them at a later time. Furthermore, this pH meter has a configurable interface with a recognition function (RS 232) so that it can be connected to a computer (bidirectional) or a recorder.

Measurement reliability

The special AutoRead function, which can be additionally activated, serves to monitor the drift of the combination electrode. The measured value is only released when the stability criteria are fulfilled. This ensures the reproducibility of measurement results.

Temperature compensation

Measurements can be performed with and without the temperature sensor. Temperature compensation of pH measurements can be effected automatically or manually. The type of temperature sensor that is connected (Pt1000 or NTC 30) is recognized automatically.

Calibration

The first option is a fully automatic one to three point calibration using buffers that have already been programmed into the meter in accordance with DIN or technical buffers from SCHOTT Instruments. The meter recognizes the buffer solutions automatically. Alternatively, conventional calibration with optionally selectable buffers can also be used. A sensor symbol indicates the status of the pH combination electrode after the automatic calibration. The adjustable calibration timer of the handylab pH 12 can be set to remind the user of any calibration that is due to be performed.

Power supply

The battery-powered meters allow the user to work independently of a mains power supply for thousands of hours. When the batteries are changed, all of the calibration data are retained in the memory. A message appears on the display in good time to remind the user to replace the batteries.

An optional power pack is also available for the handylab pH 12.



Sensors

We have a comprehensive product range of precision pH electrodes and can offer you the right sensor for every type of application. We would be pleased to advise you about your specific application.

Separately or as a set

The handylab pH 11 and handylab pH 12 pocket-size pH meters are not only available separately, but can also be purchased in a complete cost-effective set that includes a combination electrode, buffer solutions and measuring beakers in a carrying case. With this set, you can get to work right away.

Technical data

handylab pH 11, handylab pH 12

ph meters		handylab pH 11	handylab pH 12
Measuring ranges			
pH	range	-2.000...+19.999 pH	-2.000...+19.999 pH
	resolution max.	0.001 pH	0.001 pH
	accuracy	+0.005/±0.01 pH	+0.005/±0.01 pH
mV	range max.	-1999...+1999 mV	-1999...+1999 mV
	resolution max.	0.1 mV	0.1 mV
	accuracy	+0.3/±1 mV	+0.3/±1 mV
temperature	range	-5.0...+105.0 °C	-5.0...+105.0 °C
	resolution	0.1 K	0.1 K
	accuracy (with NTC 30)	±0.1 K	±0.1 K
	manual adjustment	-20...+130 °C	-20...+130 °C
drift control	can be switched off	yes	yes
slope matching		85...105 %	85...105 %
zero point matching		± 30 mV	± 30 mV
sensor evaluation	via symbol in the display	yes	yes
input resistance		>10 ¹² Ω	>10 ¹² Ω
offset current		< 10 ¹² A	< 10 ¹² A
Calibration			
buffer sets	DIN (1.68/4.01/6.87/9.18)	1-/2-/3 point	1-/2-/3 point
	technical (2.00/4.00/7.00/10.01)*	1-/2-/3 point	1-/2-/3 point
	selectable buffers	1-/2 point	1-/2 point
calibration interval control		-	1...999 days
saving calibration data in memory		-	yes
Real time clock	integrated with time/date	-	yes
Data storage			
storage by depression of key		-	800 data records
time controlled storage	in 7 intervals (5 sec. ...60 min)	-	800 data records
Connections			
electrode (socket in accord. with DIN 19 262)		yes	yes
temperature sensor (NTC 30/Pt 1000, 2 x 4 mm banana plug)		yes	yes
Interface			
for analogue recorder cable Z 394		-	4 poles socket
for RS-232 cable Z 395, bi-directional		-	4 poles socket
Ambient temperature			
operating temperature		-10...+55 °C	-10...+55 °C
relative humidity (annual average)		< 90 %	< 90 %
Power supply			
battery operated (type AA)		4 x 1.5 V mignon cells	4 x 1.5 V mignon cells
battery life time (data is saved even if batteries are changed)		approx. 5.000 h	approx. 5.000 h
power supply (no akku)		-	optionally
automatic switch-off during battery operation		60 min	60 min
Housing			
dimensions (H x W x D)		172 mm x 80 mm x 37 mm	172 mm x 80 mm x 37 mm
weight		approx. 0.3 kg	approx. 0.3 kg
Display			
LCD multi-function display		60 mm x 45 mm	60 mm x 45 mm
Instrument safety	protection class	3, EN 61010-1 A2	3, EN 61010-1 A2
	protection type	IP 66, EN 60529	IP 66, EN 60529
approvals/marks of conformity		cETLus, CE	cETLus, CE
instrument warranty		3 years	3 years

* SCHOTT Instruments

Portable conductivity meters with GLP functions handylab LF 11 and LF 12

The handylab LF 11 and LF 12 portable conductivity meters in shock-proof, water-tight casings are ideally suited for field work.

Measurement parameters

The versatile conductivity meters can be used to measure electrical conductivity, total dissolved solids (TDS), salinity and temperature.

Measurement memory and interface

In comparison with the handylab LF 11, the handylab LF 12 additionally has a data memory, which makes it possible to save measurements manually or automatically using a timer control, and then evaluate them at a later time. Furthermore, this conductivity meter has a configurable interface with a recognition function (RS-232) so that it can be connected to a computer (bidirectional) or a recorder.

Measurement reliability

The special AutoRead function, which can be additionally activated, serves to monitor the drift of the combination electrode. The measure value is only released when the stability criteria are fulfilled. This ensures the reproducibility of measurement results.

Temperature compensation

The automatic temperature compensation works in various selectable modes:

- with an adjustable linear temperature coefficient,
- with a fixed non-linear temperature coefficient or
- with the temperature compensation deactivated.

A reference temperature of 20 °C or 25 °C can be selected.

Calibration


The cell constant can be set within a wide range. In addition, there is a fixed constant of 0.01. Combination cells with a cell constant of 0.475 or 1 can also be calibrated automatically. The integrated adjustable calibration timer in the handylab LF 12 can be set to remind the user of any calibration that is due to be performed.

Power supply

The conductivity meters can be used for about 2,500 hours without a mains power supply using four conventional batteries. A reminder is shown on the display when the batteries have to be replaced. When the batteries are changed, the calibration data are retained in the memory. The handylab LF12 can also be operated with the optional powerpack.

Sensors

Either type LF 513 T electrodes (two pole technology) or type LF 613 T electrodes (four pole technology) can be utilized alternatively. Both types have an integrated temperature sensor. We would be pleased to advise you about your specific application.



Included in the set

The LF 11 and LF 12 conductivity meters are also available as part of a cost-effective set in a carrying case, which includes a combination electrode, calibration solutions and a measuring beaker. With this set, you can get to work right away.

Technical data

handylab LF 11, handylab LF 12

Parameter		handylab LF 11	handylab LF 12
Measuring ranges			
conductivity	in 5 ranges or AutoRange	0.0 µS/cm...500 mS/cm	0.0 µS/cm...500 mS/cm
	at k = 0.1 and k = 0.01	0.00 µS/cm...19.99 µS/cm	0.00 µS/cm...19.99 µS/cm
	at k = 0.01	0.000 µS/cm...1.999 µS/cm	0.000 µS/cm...1.999 µS/cm
specific resistance		0.000...1999 MΩ·cm	0.000...1999 MΩ·cm
salinity	acc. to IOT table	0.0...70.0	0.0...70.0
TDS	factor adjustable 0.40...1.00	0...1999 mg/l	0...1999 mg/l
temperature	automatic, 3 modes selectable	-5.0...+105.0 °C	-5.0...+105.0 °C
	resolution	0.1 K	0.1 K
	manual adjustment	-5...+100 °C	-5...+100 °C
Cell constants	adjustable	0.01; 0.090...0.110;	0.01; 0.090...0.110;
		0.250...2.500	0.250...2.500
	calibrate	0.450...0.500 ; 0.800...1.200	0.450...0.500 ; 0.800...1.200
	calibration interval control	-	1...999 days
Accuracy			
	conductivity	± 0.5 % of measured value	± 0.5 % of measured value
	salinity	± 0.2	± 0.2
	TDS	± 2 %	± 2 %
	temperature (NTC 30)	±0.1 K	±0.1 K
Reference temperature	selectable	20 °C or 25 °C	20 °C or 25 °C
temperature compensation mode			
	non-linear function natural water	acc. to EN 27 888 (DIN 38 404)	yes
	linear compensation	0.001...3.000 %/K	0.001...3.000 %/K
	no compensation	yes	yes
real time clock	integrated with time/date	-	yes
Data storage			
storage by depression of key		-	800 data records
time controlled storage	in 7 intervals (5 sec ... 60 min)	-	800 data records
Connections			
for 2 and 4 pole cells			
with/without temperature sensor (NTC 30)		8 poles socket	8 poles socket
Interface			
for analogue recorder cable Z 394		-	4 poles socket
for RS-232 cable Z 395, bi-directional		-	4 poles socket
Ambient temperature			
operating temperature		-10...+55 °C	-10...+55 °C
relative humidity (annual average)		< 90 %	< 90 %
Power supply			
battery operated		4 x 1.5 V mignon cells	4 x 1.5 V mignon cells
battery life time (data remain when changing batteries)		approx. 2,500 h	approx. 2,500 h
power supply (no akku)		-	optionally
automatic switch-off at operation		60 min	60 min
Housing			
dimensions (H x W x D)		ABS, water-tight key pad 172 mm x 80 mm x 37 mm	ABS, water-tight key pad 172 mm x 80 mm x 37 mm
weight		approx. 0.3 kg	approx. 0.3 kg
Display			
LCD multi-function display		60 mm x 45 mm	60 mm x 45 mm
Instrument safety			
	protection class	3, EN 61010-1 A2	3, EN 61010-1 A2
	protection type	IP 66, EN 60529	IP 66, EN 60529
approvals/marks of conformity		cETLus, CE	cETLus, CE
instrument warranty		3 years	3 years

Portable oxygen meter with GLP functions handylab OX 12

The handylab OX12 portable oxygen meter in a shock-proof, water-tight casing is also ideally suited for on-site oxygen measurements in rivers, lakes or effluent, as well as for BOD measurements.

Measurement parameters

The oxygen concentration, saturation index and temperature measurement parameters mean that the SCHOTT Instruments handylab OX 12 has a variety of uses.

Measurement memory and interface

The meter has a data memory, whereby measurements can be saved manually or automatically using a timer control, and then evaluated later on. In addition, the oxygen meter has a configurable interface with a recognition function (RS-232) so that it can be connected to a computer (bi-directional) or a recorder.

Measurement reliability

The special AutoRead function, which can be additionally activated, serves to monitor the drift of the combination electrode. The measured value is only released when the stability criteria are fulfilled. This ensures the reproducibility of measurement results.

Measurements

During the measurement process, influence variables such as temperature and air pressure are automatically taken into account and compensated. Even the influence of a higher salinity level on oxygen determination can be corrected by entering the salinity that has been determined using a conductivity meter.

Calibration

Calibration of the handylab OX12 can be performed easily on-site using the air calibration vessel. The vessel ensures a defined humidity and therefore ideal calibration conditions. After automatic calibration, a sensor symbol indicates the status of the oxygen combination electro-

de. The adjustable calibration timer can remind the user when the next calibration is due to be performed.

Power supply

The handylab OX12 can be operated for at least 2,000 hours independently of a mains power supply using four conventional batteries. A reminder for the user to replace the batteries appears on the display in good time. The calibration data are retained when the batteries are changed.

Sensor

The modern, zero current free, galvanic sensor 9009/61, which is included with the meter and can be used immediately for measuring purposes, ensures precise, reliable and rapid measurement of oxygen concentrations.

As a set

The handylab OX 12 portable oxygen meter is available as a complete set in a carrying case together with the 9009/61 sensor, the OX 925 maintenance set and the OxiCal®-SL calibration vessel.



Technical data	9009/61 O₂ sensor
Measuring principle	membrane covered galvanic sensor
temperature compensation	IMT
measurement range	0...50 mg/l O ₂
temperature range	0...50 °C
max. over-pressure	6 bar
immersion depth	min. 6 cm max. 20 m water depth
Material	membrane head and shaft: POM membrane FEP thermistor housing VA steel (1.4571)
Dimensions	shaft length: 145 mm diameter: 15.25 mm membrane thickness: 13 µm
Cable connection	fixed cable length: 1.5 m (standard); max. length: 20 m
Approach velocity	> 3 cm/s at 10 % measuring accuracy 10 cm/s at 5 % measuring accuracy 18 cm/s at 1 % measuring accuracy
Specifications of sensor when new	
zero signal	< 0.1 % of saturation value
reaction time at 20 °C	t ₉₀ (90 % of final value) after < 10 sec. t ₉₅ (95 % of final value) after < 16 sec. t ₉₉ (99 % of final value) after < 60 sec.
internal consumption	0.008 µg/h
drift	approx. 3 % per month under operating conditions
service life	min. 6 months per electrolyte filling
polarization time	not required; sensor can be used immediately

Technical data handylab OX 12

Oxygen meter		handylab OX 12
Measuring ranges		
O ₂ concentration	ranges	0.00...19.99 mg/l / 0...90.0 mg/l
	resolution max.	0.01
	accuracy	± 0.5 % of measured value
O ₂ saturation index	ranges	0.0...199.9 % / 0...600 %
	resolution max.	0.1 %
	accuracy	± 0.5 % of measured value
O ₂ partial pressure	ranges	0.0...199.9 mbar / 0...1250 mbar
temperature	range	0...+50.0 °C
	resolution	0.1 K
	accuracy	± 0.1 K
drift control	can be switched off	yes
Correction functions		
air pressure	automatic (built-in pressure sensor)	500...1100 hPa
temperature	automatic (IMT)	0...+40 °C
salinity	using setting keys	0.0...70.0
Calibration		
procedure		air calibration procedure
slope range		0.60...1.25
calibration interval control		1...999 days
calibration data storage		yes
sensor evaluation	via symbol on display	yes
real time clock	integrated with time/date	yes
Data storage		
storage by depression of key		800 data records
time controlled storage	in 7 intervals (5 sec...60 min)	800 data records
Connections		
oxygen sensor		8 pole socket
Interface		
for analogue recorder cable Z 394		4 pole socket
for RS-232 cable Z 395, bi-directional		4 pole socket
Ambient temperature		
operating temperature		-10...+55 °C
relative humidity (annual average)		< 90 %
Power supply		
battery operated		4 x 1.5 V mignon cells (type AA)
battery life time (data remain when changing batteries)		approx. 2,000 h
power supply (no akku)		optionally
automatic switch-off at operation		60 min
Housing		
		ABS, water-tight key pad
dimensions (H x W x D)		172 mm x 80 mm x 37 mm
weight		approx. 0.3 kg
Display		
LCD multi-function display		60 mm x 45 mm
Instrument safety	protection class	3, EN 61010-1 A2
	protection type	IP 66, EN 60529
approvals/marks of conformity		cETLus, CE
instrument warranty		3 years

Multi-parameter portable meters with GLP functions

handylab pH/LF 12 and handylab multi 12

The multi-parameter pocket-size meters – handylab pH/LF 12 and handylab multi 12 – in shock-proof, water-tight casings are ideally suited for field work.

Measurement parameters

The measurement parameters pH, redox potential, temperature and conductivity means that the handylab pH/LF 12 has a variety of uses. The handylab multi 12 can also be used to measure the oxygen concentration in solutions.

Measurement memory and interface

The meters have a data memory, which means that measurements can be saved manually or automatically by using a timer control, and evaluated later on. Both meters have a serial RS-232 interface (bi-directional) for data transfer purposes.

Measurement reliability

The special AutoRead function, which can be additionally activated, serves to monitor the drift of the combination electrode. The measured value is only released when the stability criteria are fulfilled. This ensures the reproducibility of measurement results.

Calibration

For calibration of the pH measurement, there is a one or two point calibration with technical buffers. For calibration of the conductivity sensor and the oxygen sensor, if needed, there is an automatic

calibration function. After automatic calibration, a sensor symbol indicates the status of the calibrated sensor. The adjustable calibration timer can remind the user that calibration is due to be performed.

Power supply


The multi-parameter portable meters can be operated for approx. 2,500 hours independently of a mains power supply using four conventional batteries, or with the optional power pack. The display shows a reminder to replace the batteries in good time. The calibration data are retained when the batteries are changed.

Sensors

We supply the handylab pH/LF12 complete with a suitable pH combination electrode and conductivity sensor. The handylab multi 12 additionally includes an oxygen sensor.

As a set

The handylab pH/LF 12 and handylab multi 12 multi-parameter meters are available as a complete set in a carrying case with all of the requisite sensors, calibration and maintenance accessories. With the set, you can get to work immediately.



handylab

Technical data

handylab pH/LF 12, handylab multi 12

Parameters		handylab pH/LF 12	handylab multi 12
Measuring ranges			
pH/mV	pH range/resolution	-2.00...+19.99 pH	-2.00...+19.99 pH
	accuracy (±1 digit)	±0.01 pH	±0.01 pH
	mV range/resolution	-1999...+1999 mV	-1999...+1999 mV
	accuracy (±1 digit)	±1 mV	±1 mV
temperature	measuring range	-5.0...+105.0 °C	-5.0...+105.0 °C
	manual setting	-20...+130 °C	-20...+130 °C
oxygen	concentration: ranges/resolution	-	0.00...19.99 mg/l/0...90.0 mg/l
	saturation: ranges/resolution	-	0.00...199.9 %/0.0...600 %
	accuracy (±1 digit)	-	±0.5 % of measured value
	temp. compensation, automatically	-	0.0...50.0 °C
conductivity	4 ranges/Auto range	1 µS/cm...500 mS/cm	1 µS/cm...500 mS/cm
	salinity acc. to IOT table	0.0...70.0	0.0...70.0
	accuracy (±1 digit)	±0.5 % of measured value	±0.5 % of measured value
	temperature compensation modes	linear, nonlinear, no compensation	nonlinear
	cell constant, calibrate	0.450...0.500	0.450...0.500
drift control	deactivatable	yes	yes
sensor evaluation	by symbol on display	yes	yes
Calibration			
pH	technical (2.00/4.00/7.00/10.01) ^{*)}	1-/2-point	1-/2-point
	DIN (1.68/4.01/6.87/9.18)	1-/2-point	-
oxygen	automatic calibration	-	yes
conductivity	automatic calibration	yes	yes
calibration interval control		1...999 days	1...999 days
calibration data store		yes	yes
real time clock	integrated with time/date	yes	yes
Serial interface			
type		RS 232, bi-directional	RS 232, bi-directional
baud rate		adjustable	adjustable
Data storage			
storage by depression of key		500 data records	500 data records
time controlled storage	in 7 levels (5 sec...60 min)	500 data records	500 data records
Input			
pH/redox-electrode (opt. with temperature sensor)		socket acc. to DIN 19262 + socket 4 mm	socket acc. to DIN 19262 + socket 4 mm
conductivity/oxygen sensor		8 pole socket	8 pole socket
Output			
for RS-232 cable Z 395, bi-directional		4 pole socket	4 pole socket
for recorder cable Z 394		4 pole socket	4 pole socket
Ambient temperature			
operating temperature		-10...+55 °C	-10...+55 °C
relative humidity (annual average)		< 90 %	< 90 %
Power supply			
battery operated		4 x 1.5 V mignon cells (type AA)	4 x 1.5 V mignon cells (type AA)
battery life time (data remain when changing batteries)		approx. 2,500 h	approx. 2,500 h
automatic switch off at battery operation		60 min	60 min
power supply (no akku)		optionally	optionally
Housing			
dimensions (H x W x D)		172 mm x 80 mm x 37 mm	172 mm x 80 mm x 37 mm
weight		approx. 0.3 kg	approx. 0.3 kg
Display			
LCD multi-function display		60 mm x 45 mm	60 mm x 45 mm
Instrument safety			
	protection class	3, EN 61010-1	3, EN 61010-1
	protection type	IP 66, EN 60529	IP 66, EN 60529
approvals/marks of conformity		cETLus, CE	cETLus, CE
instrument warranty		3 years	3 years

^{*)} SCHOTT Instruments

Order overview handylab pH meters, conductivity meters and oxygen meters

pH meters	Type no.	Order no.
handylab pH 11, individual meter	handylab pH 11	28 520 2871
handylab pH11, individual meter with carrying case	handylab pH 11/K	28 520 2863
handylab pH 11, case set, complete, ready to use with pH combination electrode BlueLine 23 pH, calibration solutions and plastic beakers	handylab pH 11/23 pH	28 520 2917
handylab pH 11, case set, complete, ready to use with pH combination electrode BlueLine 24 pH, calibration solutions and plastic beakers	handylab pH 11/24 pH	28 520 2982
handylab pH 11, case set, complete, ready to use with pH combination electrode BlueLine 14 pH, calibration solutions and plastic beakers	handylab pH 11/14 pH	28 520 2999
handylab pH 12, individual meter	handylab pH 12	28 520 2896
handylab pH12, individual meter with carrying case	handylab pH 12/K	28 520 2888
handylab pH 12, case set, complete, ready to use with pH combination electrode BlueLine 24 pH, calibration solutions and plastic beakers	handylab pH 12/24 pH	28 520 3054
handylab pH 12, case set, complete, ready to use with pH combination electrode BlueLine 14 pH, calibration solutions and plastic beakers	handylab pH 12/14 pH	28 520 3062
Conductivity meters		
handylab LF 11, individual meter	handylab LF 11	28 520 3292
handylab LF11, individual meter with carrying case	handylab LF 11/K	28 520 3276
handylab LF11, case set, complete, ready to use with 4-pole conductivity cell LF 413 T, calibration solutions and plastic beaker	handylab LF 11/413 T	28 520 3310
handylab LF 11, case set, complete, ready to use with 2-pole conductivity cell LF 513 T, calibration solutions and plastic beaker	handylab LF 11/513 T	28 520 3321
handylab LF 11, case set, complete, ready to use with 4-pole conductivity cell LF 613 T, calibration solutions and plastic beaker	handylab LF 11/613 T	28 520 3346
handylab LF 12, individual meter	handylab LF 12	28 520 3362
handylab LF12, individual meter with carrying case	handylab LF 12/K	28 520 3354
handylab LF12, case set, complete, ready to use with 4-pole conductivity cell LF 413 T, calibration solutions and plastic beaker	handylab LF 12/413 T	28 520 3330
handylab LF 12, case set, complete, ready to use with 4-pole conductivity cell LF 613 T, calibration solutions and plastic beaker	handylab LF 12/613 T	28 520 3379
Oxygen meter		
handylab OX12, individual meter with carrying case	handylab LF 12/K	106 3835
handylab OX 12, case set, complete, ready to use with oxygen sensor 9009/61, calibration and maintenance accessories	handylab OX12-Set	28 520 2793

Order overview handylab multi-parameter portable meters

Multi-Parameter meters	Type no.	Order no.
handylab pH/LF 12, individual meter	handylab pH/LF 12	28 520 3465
handylab pH/LF 12, case set, complete, ready to use with pH combination electrode BlueLine 24-3 pH, 4-pole conductivity cell LF 413-3 T, calibration and maintenance accessories	handylab pH/LF 12-Set	28 520 3473
handylab multi 12, individual meter	handylab multi 12	28 520 3502
handylab multi 12, case set, complete, ready to use with pH combination electrode BlueLine 24-3 pH, 4-pole conductivity cell LF 413-3 T, oxygen sensor 9009/63, calibration and maintenance accessories	handylab multi 12-Set	28 520 3519
Accessories		
Redox combination electrode with plug head	BlueLine 31 Rx	28 512 9311
Plug cable combination e.g. for BlueLine 31 Rx, 1 m cable, DIN plug	LB 1 A	28 512 2653
Electrolyte solution KCl 3 mol/l, 1000 ml DURAN® bottle	L 300	28 513 8554
Technical buffer solutions pH 4.00 / 7.00, 2 x 30 ampoules	L 4690	28 513 8398
Redox test solution 180, 430, 600 mV Pt/calomel; 220, 470, 640 mV Pt/Ag/AgCl, 3 x 20 ampoule	L 4648	28 513 8784
Conductivity test solutions KCl 0.01 / 0.1 / 1 mol/l (1.41 mS/cm / 12.9 mS/cm / 112 mS/cm), 3 x 6 ampoules	LF 995	28 512 6293
Field armouring with holder carrying handle and shoulder strap, for handylab pH meters	Z 384	28 520 4848
Protective armouring with holder and carrying handle, for handylab pH meters	Z 385	28 520 4856
Holder set for protective armouring, for handylab OX12 meters	Z 386	28 520 4864
Rubberized elastic protective armouring with handle support, for all handylab meters	Z 387	28 520 4872
Universal mains power supply unit, 100...240 V for all handylab 12 models	Z 850	28 520 4889
Connecting cable for analogue recorder, for handylab pH 12, LF 12, OX12	Z 394	28 520 4942
Connecting cable for PC, for all handylab 12 models (software included)	Z 395	28 520 4959

Subject to technical changes.

DURAN® is a trademark of the SCHOTT AG, Mainz, Germany.

The ideal set for reliable measuring results:

Electrodes and meters

from SCHOTT Instruments



Contents laboratory electrodes

Survey of product families BlueLine, ScienceLine and IoLine	Page 42
Selection table electrodes/applications	Page 48
ID technology and overview of ID electrodes	Page 52
IoLine pH combination electrodes	Page 56
IoLine pH combination electrodes with temperature sensor	Page 58
ScienceLine pH combination electrodes	Page 60
ScienceLine pH combination electrodes with temperature sensor	Page 62
ScienceLine micro, spear tip and surface pH combination electrodes	Page 64
ScienceLine metal combination electrodes	Page 66
ScienceLine single electrodes: pH glass electrodes, metal electrodes	Page 68
ScienceLine single electrodes: reference electrodes	Page 70
ScienceLine conductivity measuring cells	Page 72
ScienceLine sensors for ammonia, sodium, oxygen, ion-selective electrodes	Page 76
ScienceLine resistance thermometers	Page 78
BlueLine pH combination electrodes	Page 80
BlueLine, special sensors	Page 82
Connection cables	Page 84
Solutions	Page 86
Electrolyte bridges/other accessories	Page 92
Tips and notes for successful measurement	Page 94
Index	Page 96
ProcessLine electrodes:	
Insight into our extensive programme for analytics in process	Page 98

Laboratory electrodes from SCHOTT Instruments: application orientated and perfectly matching


The standards for modern pH measurement are very high regarding precision, reproducibility, speed, comfortable handling and reliability. Every measurement is different. Different compositions, temperatures, conductivities and viscosities of samples and different measuring conditions make for a million of different applications. Only application orientated and perfectly matching systems of electrodes, meters and buffer solutions can meet these standards. Our concern at SCHOTT Instruments is to supply such systems.

The pH electrode is a very important part of the system as it comes in direct contact with the sample and delivers the measurement signal. For more than 70 years our focus has been set on the electrode and we have dedicated ourselves to the development and manufacturing of glass electrodes. For a remarkable long time our electrodes are being used for the most demanding tasks in labs throughout the world where quality really matters, and our customers benefit from this know-how.



The glassmaker's skill is as important as ever: Our first instruction booklet appeared in 1938. In those days the electrochemical pH measuring and the potentiometric titration still needed to be explained.

It all started with a patent on pH electrodes – today it is a range of several hundred different sensors. Our electrode programme comprising the three product families BlueLine, ScienceLine and IoLine is as manifold as your applications. Whether for ultrapure water, jam, wine, creme or drinking water, SCHOTT Instruments offers the right electrode for every conceivable use.



Even today glass blowing talent is still indispensable.

BlueLine

Attractive shape, reliable function

The compact BlueLine range is a basic series comprising electrodes for the usual laboratory applications which makes the user's choice easier.

An accurate and fast measurement is ensured by precision manufacturing and employment of high quality material, e.g. low-resistance A-type membrane glass or the unique platinum diaphragm.

The BlueLine family comprises the robust electrodes with gel electrolyte and plastic shaft for general use, the liquid-electrolyte sensors for more critical measurements and special sensors. The gel and the liquid electrolyte sensors are available with different connections (e.g. plug head or fixed cable with either DIN or BNC type plug) and with different cable lengths and optionally with built-in temperature sensors (NTC 30 k Ω or Pt 1000). The special electrodes range includes pH electrodes for surface measurements, for small sample amounts, for ultrapure water and emulsions or measurements in semi-solid samples (insert measurements).



- ▶ **Compact basic series**
for the usual applications, made from the universal A-type membrane glass
- ▶ **Comfortable handling, stylish design and reliable function**
- ▶ **Gel electrolyte, liquid electrolyte and special sensors**
- ▶ **Liquid electrolyte electrodes with unique platinum diaphragm and refill port slider**
- ▶ **Each electrode with individual serial number**

Advantages
BlueLine

ScienceLine

The proven high end laboratory electrodes

In research and development, manufacturing and quality control, our ScienceLine electrodes have long since become standard for the most demanding measuring tasks. Each electrode has an individual serial number and pH- and metal combination electrodes are supplied with a quality certificate, thus making documentation easier and better traceable. We have kept on improving the glass membrane shapes and types to make the electrodes even more

robust, durable and easier to clean. Furthermore, they achieve stable measurement values even faster.

The ScienceLine electrodes from SCHOTT Instruments guarantee not only high measurement accuracy and stability and a long service life, but are highly adaptable to your measuring tasks. Today we can offer you a range of electrodes unmatched in versatility and quality.





Here are some examples:

- pH electrodes with a length of up to 600 mm for measurements in very deep vessels
- On the other hand, the N 6003 electrode allows pH measurements even in NMR tubes or other small sample vessels. Likewise, the A 157 is a micro electrode with an integrated temperature sensor with only 5 mm in diameter.
- For more demanding media, you can choose among different diaphragms and membrane glasses. For measurements in samples low on ions, you have the choice between the N 64 and A 164. Both types feature a ground joint diaphragm, and the A 164 additionally even a temperature sensor.
- A wide selection of separate reference and glass electrodes completes the programme.

The faster and more stable display of the measured value under critical circumstances with Science Line electrodes, as well as their longer lifespan are due to their Silamid reference system. In contrast to a chlorinated silver wire, as it is used in the silver/silver chlorine reference system of the BlueLine series, the ScienceLine employs discharge cartridges. The silver coating of the inner tube makes for a five times larger silver surface compared to the silver wire. Hence, the stability of the potential is much higher.

- ▶ **Proven high-end electrodes for demanding measuring tasks**
- ▶ **Silamid reference system for fast and stable acquiring of measured values and for longer electrode lifespan.**
- ▶ **Utmost versatility** is achieved by a large selection of diaphragms, membrane glass types and shapes, shaft lengths and diameters, ground joints, plug connections and integrated temperature sensors.
- ▶ **Each pH and metal combination electrode with individual serial number and quality certificate.**
- ▶ **Large selection** of separate glass and reference electrodes, metal combination electrodes, conductivity sensors, ion selective indicator electrodes and ammonia, sodium and oxygen sensors.


Advantages
ScienceLine

IoLine pH electrodes for the most demanding measuring tasks

Patented three-chambers system with Iodine reservoir in the iodine/ iodide reference electrode

The reference system is a very important part of the pH electrode. The standard hydrogen electrode has proven too difficult in practical use to gain more than a mere theoretical importance. The Ag/AgCl system, which is nowadays almost exclusively used, can cause measuring instabilities originating from potential variations with changing temperatures or reactions between the silver ions and the measuring solution in the area of the diaphragm.

IoLine electrodes, in contrast, have the advantage of a much lower temperature sensitivity and a metal ion free reference system. The reference system is based on the following reaction:



The ORP is thereby described by the Nernstian equation:

$$\text{EH} = \text{E}^\circ + \text{RT}/\text{zF} * \ln ([\text{I}_3^-] / [\text{I}^-]^3)$$


Whereby

$\text{E}^\circ = 0.536 \text{ V}$, $\text{R} = 8.314472 \text{ J/(K} \cdot \text{mol)}$,
 T in K, $\text{z} = 2$ und $\text{F} = 96485.34 \text{ C/mol}$.

The stability of the reference system potential even at changing temperatures is the key to the IoLine electrodes' superior response speed and measurement stability and, additionally, higher accuracy compared to conventional Ag/AgCl electrodes.


A further advantage is the fact that the component determining the potential, i.e. the Iodine, is continuously being resupplied from the patented three-chambers reservoir system. The first chamber contains a supply which is used to fresh up the reference electrolyte and the reference system of the second chamber with iodine in order to maintain the necessary I_3^-/I^- concentration and thereby a stable reference potential. The electrolyte link between the second and the third chamber containing the exchangeable bridge electrolyte is provided by an interior diaphragm. The bridge electrolyte also provides the contact with the sample via the diaphragm.

The interior diaphragm is designed to allow only a minimal diffusion of I_3^-/I^- (Tri-iodide/Iodide) into the bridge electrolyte. Hence, the iodine consumption in the reference electrolyte is very low and needs only very little re-supplying from the first chamber. The supply is practically unlimited and guarantees a high stability of the potential and a long lifespan of the IoLine electrodes.



Platinum diaphragm for fast response and high stability

Further responsible for the high stability and fast response of the iodine/iodide reference system is the platinum diaphragm, which has been developed by SCHOTT. The platinum diaphragm makes for remarkably constant and reproducible measuring characteristics of the electrode. It contains twisted platinum wires being fused into the glass shaft of the electrode. The defined spaces of the platinum wires guarantee a continuous and steady electrolyte flow and high stability of the reference system in all media and at changing temperatures.



- ▶ **Unique iodine/iodide reference system**
in connection with the patented three-chambers reservoir system including the large iodine reservoir offering unbeatable stability, fast response times and high accuracy at a higher speed compared to the electrodes with the usual Ag/AgCl reference system. Furthermore independent from sample composition and temperature.
- ▶ **100% metal ion free reference system**
avoids the measuring media of being contaminated by metal ions i.e. optimal for measurements in Tris buffer.
- ▶ **Exchangeable bridge electrolyte**
enables a matching of electrolyte solution and measuring sample.
- ▶ **Wide application area**
Ideal for most precise pH measurements in manifold media for research and quality control i.e. in pharmacy, biotechnology or food industry.
- ▶ **Electrode head with integrated refill port**
enables an easy refilling of the reference system in connection with the refill port slider.
- ▶ **Manifold selection:**
Many variants regarding the connection, the membrane glass types and shapes as well as diaphragms.
- ▶ **Extensive delivery scope:**
Liquid vessel with bayonet connector to avoid drying out and for comfortable storage of the electrode and certificate are included in the delivery scope.

Advantages
IoLine

The corresponding sensor for any application:

Application recommendations for pH and ORP electrodes

The table provides an orientation to the large variety of our electrodes. The listed electrodes give an example for similar measuring models, i.e. variation only regarding the connecting system or the integrated temperature sensor. The electrode BlueLine 11 pH for instance also represents the versions 12 pH, 14 pH, 15 pH, 17 pH, 18 pH and 19 pH. Regarding the ScienceLine and loLine pH electrodes there is special attention to be drawn to the versions N 62 and H 62 as well as to IL-pH-A120MF and IL-pH-H120MF; these models are also available with larger shaft lengths.

An extension of length under the same application conditions delivers faster and more stable measuring results; additionally it enables a longer lifespan of the electrode. The higher electrolyte stand along with the increased electrolyte outflow reduces unwanted diffusion potentials on the diaphragm and rinses it free.

Some applications may require other electrode recommendations due to certain operation conditions, as identical applications can differ fundamentally with varying concentrations and temperatures. Please also note the material resistance of the sensor towards the measuring media. The recommended and additional sensors with the corresponding technical data are stated on the following pages of our catalogue. Please contact us via telephone, fax or email whenever you cannot find your application or have any queries regarding certain operation conditions.



... and conductivity cells

Electrode series		IoLine		pH measuring ScienceLine														BlueLine					ORP ScienceLine				BL*		Conductivity ScienceLine											
Application area	Sensor example	IL-pH-A120MF	IL-pH-H120MF	IL-Micro-pH-A	IL-SP-pH-A	A 157	A 7780	H 62	H 64	N 1048 A	L 32	L 39	L 6880	L 8280	N 62	N 64	N 6000 A	N 6003	11 pH	22 pH	13 pH	16 pH	21 pH	27 pH	Ag 6280	Pt 62	Pt 6140	Pt 8280	Pt 5900 A	31 RX	32 RX	LF 213 T	LF 313 T	LF 313 T NFTC	LF 413 T	LF 613 T	LF 713 T			
	Application																																							
Chemistry	Etching and degreasing baths	■	■					■	■						■	■			■	■					■					■								■		
	Bleach and dyeing solutions	■	■					■	■						■	■			■	■					■					■								■		
	Cutting oil emulsions	■												■	■	■	■		■	■					■					■								■		
	Cyanide detoxification	■	■					■	■						■	■			■	■					■					■								■		
	Dispersion paint	■	■					■	■						■	■			■	■					■					■								■		
	Emulsions, water-based	■	■					■	■					■	■	■			■	■					■					■								■		
	Emulsions, partly water-based	■													■	■			■	■					■					■							■			
	Paint/varnish, water-soluble	■	■					■	■						■	■			■	■					■					■								■		
	Fixing bath	■	■					■	■						■	■			■	■					■					■								■		
	Varnish, water-based	■	■					■	■						■	■			■	■					■					■								■		
	Varnish, partly water-based	■													■	■			■	■					■					■								■		
	Lye, extreme		■					■	■						■	■			■	■					■					■								■		
	Oil/water-emulsions	■													■	■			■	■					■					■								■		
	Organic percentile high	■													■	■			■	■					■					■								■		
	Paper extract	■	■					■	■						■	■			■	■					■					■									■	
	Acid, extreme	■	■					■	■						■	■			■	■					■					■									■	
	Sulphide containing liquid	■	■												■	■			■	■					■					■									■	
	Suspension, water-based	■	■					■	■					■	■	■			■	■					■					■									■	
	Ink	■	■					■	■						■	■			■	■					■					■									■	
	Viscose samples	■								■						■			■	■					■					■										■
Field measurements	Beck	■					■					■		■	■	■			■	■					■			■		■	■								■	
	Ground water	■					■					■		■	■	■			■	■					■			■		■	■								■	
	Lake water	■					■					■		■	■	■			■	■					■			■		■	■								■	
	Seawater	■					■					■		■	■	■			■	■					■			■		■	■								■	
	Rain water	■					■					■		■	■	■			■	■					■			■		■	■								■	
Drinks production	Beer	■					■					■		■	■	■			■	■					■			■		■	■									■
	Fruit juice	■					■					■		■	■	■			■	■					■			■		■	■								■	
	Vegetable juice	■					■					■		■	■	■			■	■					■			■		■	■								■	
	Lemonades/soda	■					■					■		■	■	■			■	■					■			■		■	■								■	
	Mineral water	■					■					■		■	■	■			■	■					■			■		■	■								■	
	Juice	■					■					■		■	■	■			■	■					■			■		■	■								■	
	Spirits	■					■					■		■	■	■			■	■					■			■		■	■								■	
	Wine	■					■					■		■	■	■			■	■					■			■		■	■								■	

* BL = BlueLine

* BL = BlueLine

Further application recommendations for pH

Electrode series		IoLine		pH measuring ScienceLine										BlueLine			ORP ScienceLine			BL*	Conductivity ScienceLine																	
Application area	Sensor example	IL-pH-A120MF	IL-pH-H120MF	IL-Micro-pH-A	IL-SP-pH-A	A 157	A 7780	H 62	H 64	N 1048 A	L 32	L 39	L 6880	L 8280	N 62	N 64	N 6000 A	N 6003	11 pH	22 pH	13 pH	16 pH	21 pH	27 pH	Ag 6280	Pt 62	Pt 6140	Pt 8280	Pt 5900 A	31 RX	32 RX	LF 213 T	LF 313 T	LF 313 T NFTC	LF 413 T	LF 613 T	LF 713 T	
	Application																																					
Cosmetics	Creme	■			■					■		■	■		■	■			■	■	■	■	■	■		■				■				■	■			
	Hair dye	■													■	■			■		■					■				■				■	■			
	Hair gel	■			■					■		■	■			■	■			■				■	■	■				■				■				
	Hair mousse	■													■	■			■	■						■				■				■	■			
	Lotions	■													■	■			■		■	■				■				■				■	■			
	Make-up	■													■	■			■	■						■				■					■	■		
	Mouth wash	■													■	■			■	■						■				■				■	■			
	Shaving foam	■													■	■			■	■						■				■				■	■			
	Sun lotion	■													■	■			■	■						■				■				■	■			
	Tooth paste	■			■						■		■	■			■			■	■					■				■				■	■			
Agriculture	Ground (extract/slug)	■					■				■			■	■	■			■	■	■					■		■		■	■					■		
	Fertilizer solution	■					■					■		■	■	■			■	■						■		■		■	■					■		
	Vegetables	■			■					■		■	■						■	■						■		■		■	■							
	Liquid manure	■					■				■			■	■	■			■	■						■		■		■	■					■		
	Fruit	■			■					■			■											■			■				■	■						
Food production	Bread/dough/pastry	■			■					■			■											■			■											
	Vinegar	■					■								■	■			■	■	■					■		■		■					■		■	
	Grease	■													■	■			■	■						■				■					■			
	Fish	■			■					■			■											■			■											
	Meat	■			■					■			■											■			■											
	Honey	■			■					■			■													■					■					■		
	Margarine	■													■	■			■	■						■				■						■	■	
	Coffee extract	■													■	■			■	■						■		■		■					■	■		
	Jam/marmelade	■													■	■			■	■						■				■						■	■	
	Mayonnaise	■													■	■			■	■						■		■		■						■	■	
	Sausage	■			■					■		■	■											■	■		■		■									
Dairy	Butter	■			■										■	■			■	■						■				■						■		
	Yoghurt	■					■				■				■	■			■	■						■				■						■		
	Cheese	■			■					■		■	■		■	■								■	■	■	■		■							■		
	Milk	■					■				■				■	■			■	■						■				■					■			
	Cream	■													■	■			■	■						■				■					■			
Surface	Skin											■												■														
	Leather											■												■														
	Paper											■												■														
	Textiles											■												■														

* BL = BlueLine

* BL = BlueLine

... ORP electrodes and conductivity cells

Electrode series		IoLine		pH measuring ScienceLine										BlueLine					ORP ScienceLine				BL*		Conductivity ScienceLine														
Application area	Sensor example	IL-pH-A120MF	IL-pH-H120MF	IL-Micro-pH-A	IL-SP-pH-A	A 157	A 7780	H 62	H 64	N 1048 A	L 32	L 39	L 6880	L 8280	N 62	N 64	N 6000 A	N 6003	11 pH	22 pH	13 pH	16 pH	21 pH	27 pH	Ag 6280	Pt 62	Pt 6140	Pt 8280	Pt 5900 A	31 RX	32 RX	LF 213 T	LF 313 T	LF 313 T NFTC	LF 413 T	LF 613 T	LF 713 T		
	Application																																						
Pharmacy, biology, biotechnology, medicine, microbiology	Agar-agar gel				■					■		■	■										■	■			■												
	Enzyme solution	■	■	■																■		■				■		■			■	■						■	
	Infusion solutions	■	■	■											■	■				■		■				■				■	■							■	
	Small vessels/sample quantity			■		■												■	■				■						■										
	Bacteria cultures	■	■	■		■	■								■	■	■	■	■	■		■	■	■			■				■	■					■		
	Gastric juice		■												■	■				■		■				■					■	■					■		
	NMR tubes																	■									■												
	Precision measurement	■	■	■												■	■									■					■	■							
	Protein containing liquid	■	■	■												■	■									■					■	■							
	Serum	■	■	■		■										■	■	■	■				■			■				■	■							■	
	Tris puffer	■		■																							■											■	
	Urine	■														■	■			■		■				■											■		
	Vials			■		■												■	■				■							■									
Technical	Cooling water	■					■								■	■			■		■					■					■			■					
	Lye, hot		■					■	■																	■					■						■		
	Acid, hot		■					■	■																	■					■						■		
Washing agents	Detergents	■													■	■			■		■	■				■				■				■	■	■	■	■	
	Disinfectant	■													■	■			■		■	■				■				■				■	■	■	■		
	Cleaning agent	■													■	■			■		■	■				■				■				■	■	■	■		
	Soap solution	■													■	■			■		■	■				■				■				■	■	■	■		
	Dishwashing liquid	■													■	■			■		■	■				■				■				■	■	■	■		
	Tenside solution	■													■	■			■		■	■				■				■				■	■	■	■		
Water	Waste water, general	■					■	■			■			■	■	■			■	■	■	■				■	■		■		■	■				■	■	■	
	Aquarium water	■					■				■				■	■			■	■	■	■				■	■		■		■	■			■	■	■		
	Demineralization/ion exchanger	■													■	■			■		■	■				■				■			■	■					
	pH values, extreme		■					■	■													■				■					■						■		
	Media containing low ions	■					■								■	■			■		■	■				■				■			■	■					
	Boiler feed water	■					■								■	■			■		■	■				■				■		■	■	■					
	Condensate	■					■								■	■			■		■	■				■				■		■	■	■					
	Purity water	■															■		■		■	■				■				■		■							
	Salt solution		■				■	■			■				■					■		■					■		■		■	■			■	■			
	Drinking water	■					■									■	■			■		■				■		■		■	■		■	■	■	■			
	Drops											■												■															

* BL = BlueLine

* BL = BlueLine

ID electrodes – reliable and precise pH measurements through automatic electrode recognition


Measurements of utmost accuracy with matching systems from SCHOTT Instruments

The demand for accuracy, reproducibility and primarily the stability of the pH measurement is exceedingly high. It has become even more important to have an application focussed measuring system at hand, consisting of electrodes, measuring instrument and buffer solutions, as each measurement is unique. SCHOTT Instruments has taken this into account and hence offers premium components ideally corresponding to all your applications. Only a perfectly harmonising system will enable measurements of utmost accuracy.

Automatic electrode recognition guarantees the accuracy of the measurement

Basis for the accuracy of pH measurement is the calibration. So far the main efforts for optimizing the calibration focussed on instrument-based help functions, such as for example the automatic buffer recognition. Still the uncertainty regarding the electrode remained, not knowing for sure if the connected electrode was the one matching the calibration data stored with the instrument. For assuring that the electrode's slope and zero point corresponded to the instrument stored data for calculating the pH value it was necessary to calibrate once more.

The Lab 870 and 970 and the instruments belonging to the ProLab family can recognize the ID electrodes automatically – making the wishes of all quality control reality. The electrode recognition works with especially for this purpose developed ID sensors. Individual sensor data is now transferred to the instrument automatically and wireless thanks to a tiny transponder integrated in the cable plug. The sensor data with pH electrodes is for example the slope and zero point, data of last calibration, buffer used as well as sensor type and serial number. The measuring instrument always uses this specific data for each corresponding ID electrode in order to calculate the pH value out of the measured mV, regardless whether using the ID electrode with multiple instruments or various ID electrodes on one instrument.





Continuous updating of sensor data

When the ID electrode is calibrated afresh with the Lab 870/970 or the ProLabs, the data on the electrode are accordingly updated. The instrument will recognize and access the data with the next measurement. This process runs in the background without any interaction of the user. The result shows stable measurements and no necessity for repeated calibrations when changing the sensor.

Records include sensor type and serial number

The instruments also fulfil the increasing documentation requirements. The entire data, such as sensor type and serial number along with measuring values, data and time are part of the measuring record, which can be transferred to the PC via interfaces.

- ▶ **Each ID sensor with unique identification.**
- ▶ **Qualified quality**
- ▶ **Highest comfort** – data exchange between sensor and instrument runs fully automatic in the background.
- ▶ Several ID sensors can interact with one instrument and one ID sensor can interact with multiple instruments via recognition, **without having to calibrate at every changing process.**
- ▶ **Accurate and safe measurements** through continuous usage of sensor specified data.
- ▶ **GLP at its best:**
Automatic and complete documentation of calibrations and measurements including the used electrode (model and serial number) with date, time and measuring values.

Advantages
ID

ID electrodes for highest safety

... with fixed cable and integrated electrode recognition

ID electrodes for pH measurement

Shaft material: glass
Zero point: pH = 7.0 ± 0.3
pH range: 0 ... 14
Reference system¹⁾: iodine/iodid,
 Silamid®,
 Ag/AgCl
Reference electrolyte: KCl 3 mol/l
Fixed cable: 1 m long,
 with DIN or
 BNC plug and
 banana plug
 with the versions
 including an
 integrated
 temperature
 sensor

¹⁾ please view the following pages for
 the technical data of each electrode



IL-pHT-
A120-DIN-N

IL-pHT-A120-
BNC-N

IL-pHT-
A170-DIN-N

IL-pHT-A170-
BNC-N

A 7780
1M-DIN-ID

A 7780
1M-BNC-ID

A 161
1M-DIN-ID

A 161
1M-BNC-ID

A 164
1M-DIN-ID

A 164
1M-BNC-ID

BlueLine 14
pH ID

BlueLine 15
pH ID

ID electrodes for conductivity measuring with temperature sensor

Temperature sensor: NTC 30 kΩ
Fixed cable: 1 m long,
 8-pole plug

* LFOX 1400 ID additionally with oxy-
 gen measurement



LF 213 T-ID

LF 313 T-ID

LF 413 T-ID

LF 913 T-ID

LFOX 1400 ID*

A selection of the SCHOTT Instruments ID electrodes programme

ID electrodes for pH measuring

Micro, spear tip and surface combination electrodes

Shaft material: glass
(except
BlueLine 21:
plastic shaft)

Zero point: pH = 7.0 ± 0.3

pH range: 0 ... 14
(except
BlueLine 21
and 27:
1 ... 13 pH)

Reference system¹⁾: iodine/iodid,
Silamid®,
Ag/AgCl

Reference electrolyte: KCl 3 mol/l,
gel or Referid®

Fixed cable: 1 m long,
with DIN or
BNC plug and
banana plug
with the versions
including an
integrated
temperature
sensor



IL-Micro- pHT-A- DIN-N	L 6880 1M-DIN-ID	N 1048 1M-DIN-ID	L 39 1M-DIN-ID	N 6000 1M-DIN-ID	BlueLine 21 pH 1M-DIN-ID	BlueLine 27 pH 1M-DIN-ID
IL-Micro-pHT- A-BNC-N	L 6880 1M-BNC-ID	N 1048 1M-BNC-ID	L 39 1M-BNC-ID	N 6000 1M-BNC-ID	BlueLine 21 pH 1M-BNC-ID	BlueLine 27 pH 1M-BNC-ID

¹⁾ please view the following pages for
the technical data of each electrode

electrodes

IoLine pH combination electrodes

pH combination electrodes

Reference system: iodine/iodide
Zero point: pH = 7,00 ± 0.25
pH range: 0 ... 14
Temp. range: -5 ... 100 °C
Shaft material: glass

A IL-pH-A120-MF
 IL-pH-A120

B IL-pH-A170-MF
 IL-pH-A170

C IL-pH-A120-MF-DIN
 IL-pH-A120-DIN
 IL-pH-A120-MF-BNC
 IL-pH-A120-BNC


D IL-pH-A170-MF-DIN
 IL-pH-A170-DIN
 IL-pH-A170-MF-BNC
 IL-pH-A170-BNC

E IL-Micro-pH-A
 IL-Micro-pH-A-DIN
 IL-Micro-pH-A-BNC

F IL-SP-pH-A
 IL-SP-pH-A-DIN
 IL-SP-pH-A-BNC

Glossary

IL	IoLine
Micro	micro pH electrode for measuring in small sample vessels
SP	spear tip pH electrode for measuring in solid and semi-solid samples
pH	pH combination electrodes
pHT	pH combination electrodes with temperature sensor
A	A-type membrane glass
H	H-type membrane glass
120	120 mm overall length
170	170 mm overall length
MF	platinum diaphragm (multi-flow)
DIN	DIN instrument plug
BNC	BNC instrument plug
N	4 mm banana plug
CI	cinch plug



Type No.	Order No.	Length L [mm]	Ø [mm]	Dia- phragm	Membrane glass	Membrane glass resistance	Membrane shape	Connection	Appli- cation	Shape
IL-pH-A120MF	285114140	120	12	platinum	A	200 MΩ	sphere	Screw plug head S7	■	A
IL-pH-A120	285114150	120	12	ceramic	A	200 MΩ	sphere	Screw plug head S7	■	A
IL-pH-A170MF	285114180	170	12	platinum	A	200 MΩ	sphere	Screw plug head S7	■	B
IL-pH-A170	285114190	170	12	ceramic	A	200 MΩ	sphere	Screw plug head S7	■	B
IL-pH-A120MF-DIN	285113810	120	12	platinum	A	200 MΩ	sphere	DIN ¹⁾	■	C
IL-pH-A120-DIN	285113820	120	12	ceramic	A	200 MΩ	sphere	DIN ¹⁾	■	C
IL-pH-A120MF-BNC	285114160	120	12	platinum	A	200 MΩ	sphere	BNC ¹⁾	■	C
IL-pH-A120-BNC	285114170	120	12	ceramic	A	200 MΩ	sphere	BNC ¹⁾	■	C
IL-pH-A120MF-R	285114410	120	12	platinum	A	200 MΩ	sphere	Metrohm plug ¹⁾	■	C
IL-pH-A170MF-DIN	285113830	170	12	platinum	A	200 MΩ	sphere	DIN ¹⁾	■	D
IL-pH-A170-DIN	285113840	170	12	ceramic	A	200 MΩ	sphere	DIN ¹⁾	■	D
IL-pH-A170MF-BNC	285114340	170	12	platinum	A	200 MΩ	sphere	BNC ¹⁾	■	D
IL-pH-A170-BNC	285114350	170	12	ceramic	A	200 MΩ	sphere	BNC ¹⁾	■	D
IL-pH-A170MF-R	285114420	170	12	platinum	A	200 MΩ	sphere	Metrohm plug ¹⁾	■	D
IL-MICRO-pH-A	285114280	170 (130/40)	12/6	platinum	A	400 MΩ	cylindrical	Screw plug head S7	■	E
IL-MICRO-pH-A-DIN	285113930	170 (130/40)	12/6	platinum	A	400 MΩ	cylindrical	DIN ¹⁾	■	E
IL-MICRO-pH-A-BNC	285114290	170 (130/40)	12/6	platinum	A	400 MΩ	cylindrical	BNC ¹⁾	■	E
IL-SP-pH-A	285114320	120 (50/70)	12/8	ceramic	A	400 MΩ	spear	Screw plug head S7	■	F
IL-SP-pH-A-DIN	285113940	120 (50/70)	12/8	ceramic	A	400 MΩ	spear	Screw plug head S7	■	F
IL-SP-pH-A-BNC	285114330	120 (50/70)	12/8	ceramic	A	400 MΩ	spear	Screw plug head S7	■	F

- general applications,
low ion media
- small sample quantities
- insert measurement

¹⁾ with 1 m fixed cable

IoLine pH combination electrodes with temperature sensor

pH combination electrodes with temperature sensor

Reference system: iodine/iodide
 Zero point: pH = 7,00 ±0.25
 pH range: 0 ... 14
 Temp. range: -5 ... 100 °C
 Shaft material: glass



- A
 - IL-pHT-A120MF-DIN-N
 - IL-pHT-A120-DIN-N
 - IL-pHT-A120MF-BNC-N
 - IL-pHT-A120-BNC-N
 - IL-pHT-H120MF-DIN-N
 - IL-pHT-H120-DIN-N
 - IL-pHT-H120MF-BNC-N
 - IL-pHT-H120-BNC-N
- B
 - IL-pHT-A170MF-DIN-N
 - IL-pHT-A170-DIN-N
 - IL-pHT-A170MF-BNC-N
 - IL-pHT-A170-BNC-N
 - IL-pHT-H170MF-DIN-N
 - IL-pHT-H170-DIN-N
 - IL-pHT-H170MF-BNC-N
 - IL-pHT-H170-BNC-N
- C
 - IL-MICRO-pHT-A-DIN-N
 - IL-MICRO-pHT-A-BNC-N

Glossary	
IL	IoLine
Micro	micro pH electrode for measuring in small sample vessels
SP	spear tip pH electrode for measuring in solid and semi-solid samples
pH	pH combination electrodes
pHT	pH combination electrodes with temperature sensor
A	A-type membrane glass
H	H-type membrane glass
120	120 mm overall length
170	170 mm overall length
MF	platinum diaphragm (multi-flow)
DIN	DIN instrument plug
BNC	BNC instrument plug
N	4 mm banana plug
CI	cinch plug

Type No.	Order No.	Length L [mm]	Ø [mm]	Dia- phragm	Membrane glass	Membrane glass resistance	Membrane shape	Temp. sensor	ID Function	Connection with 1 m fixed cable	Appli- cation	Shape
IL-pHT-A120MF-DIN-N	285113890	120	12	platinum	A	200 MΩ	sphere	Pt 1000	yes	DIN + banana plug	■	A
IL-pHT-A120-DIN-N	285113900	120	12	ceramic	A	200 MΩ	sphere	Pt 1000	yes	DIN + banana plug	■	A
IL-pHT-A120MF-BNC-N	285113850	120	12	platinum	A	200 MΩ	sphere	Pt 1000	yes	BNC + banana plug	■	A
IL-pHT-A120-BNC-N	285113860	120	12	ceramic	A	200 MΩ	sphere	Pt 1000	yes	BNC + banana plug	■	A
IL-pHT-H120MF-DIN-N	285113870	120	12	platinum	H	300 MΩ	sphere	Pt 1000	yes	DIN + banana plug	■	A
IL-pHT-H120-DIN-N	285113880	120	12	ceramic	H	300 MΩ	sphere	Pt 1000	yes	DIN + banana plug	■	A
IL-pHT-H120MF-BNC-N	285114200	120	12	platinum	H	300 MΩ	sphere	Pt 1000	yes	BNC + banana plug	■	A
IL-pHT-H120-BNC-N	285114210	120	12	ceramic	H	300 MΩ	sphere	Pt 1000	yes	BNC + banana plug	■	A
IL-pHT-A120MF-BNC-CI	285114370	120	12	platinum	A	200 MΩ	sphere	NTC 30 kΩ		BNC + Cinch	■	A
IL-pHT-A120MF-R-NN	285114390	120	12	platinum	A	200 MΩ	sphere	Pt 1000		Methrom plug + 2 banana plug	■	A
IL-pHT-A170MF-DIN-N	285113910	170	12	platinum	A	200 Mohm	sphere	Pt 1000	yes	DIN + banana plug	■	B
IL-pHT-A170-DIN-N	285113920	170	12	ceramic	A	200 MΩ	sphere	Pt 1000	yes	DIN + banana plug	■	B
IL-pHT-A170MF-BNC-N	285114220	170	12	platinum	A	200 MΩ	sphere	Pt 1000	yes	BNC + banana plug	■	B
IL-pHT-A170-BNC-N	285114230	170	12	ceramic	A	200 MΩ	sphere	Pt 1000	yes	BNC + banana plug	■	B
IL-pHT-H170MF-DIN-N	285114240	170	12	platinum	H	300 MΩ	sphere	Pt 1000	yes	DIN + banana plug	■	B
IL-pHT-H170-DIN-N	285114250	170	12	ceramic	H	300 MΩ	sphere	Pt 1000	yes	DIN + banana plug	■	B
IL-pHT-H170MF-BNC-N	285114260	170	12	platinum	H	300 MΩ	sphere	Pt 1000	yes	BNC + banana plug	■	B
IL-pHT-H170-BNC-N	285114270	170	12	ceramic	H	300 MΩ	sphere	Pt 1000	yes	BNC + banana plug	■	B
IL-pHT-A170MF-BNC-CI	285114380	170	12	platinum	A	200 MΩ	sphere	NTC 30 kΩ		BNC + Cinch	■	B
IL-pHT-A170MF-R-NN	285114400	170	12	platinum	A	200 MΩ	sphere	Pt 1000		Methrom plug + 2 banana plug	■	B
IL-MICRO-pHT-A-DIN-N	285114300	170 (130/40)	12/6	platinum	A	400 MΩ	cylindrical	Pt 1000	yes	DIN + banana plug	■	C
IL-MICRO-pHT-A-BNC-N	285114310	170 (130/40)	12/6	platinum	A	400 MΩ	cylindrical	Pt 1000	yes	BNC + banana plug	■	C

- general applications, low ion media
- small sample quantities
- high temperatures, optimized for alkaline area

ScienceLine pH combination electrodes

pH combination electrodes with plug head and fixed cable

Reference system: Silamid®
Shaft material: glass
Zero point: pH = 7.0 ± 0.3
Electrolyte: KCl 3 mol/l
 (except N 6250: KCl 4.2 mol/l
 A 7780 and L 7780: gel electrolyte
 L 8280: Referid® electrolyte)
Membrane shape: sphere
pH range: 0 .. 14
Connection cable for plug head: e.g. L 1 A
 (See also page with connection cables)
fixed cable: 1 m long, with plug A acc. to DIN 19262 or with BNC plug



H 61
 H 62
 H 63
 N 61
 N 62
 H 6180
 H 6280
 H 6380
 N 6180
 N 6250
 N 6280
 N 42 A
 N 42 BNC
 N 50 A
 N 52 A
 N 52 BNC
 N 61 eis
 H 61-500
 H 61-600

H 64
 H 64 1M-DIN-ID
 H 64 1M-BNC-ID
 N 64
 N 6480 eis
 N 6480 eth

N 65
 H 65
 H 6580
 N 6580

L 32

A 7780
 L 7780

N 6980

L 8280

Order No.	Type No.	Length L [mm]	Ø [mm]	Dia- phragm	pH- glass	Temp. range [°C]	Connection	Remarks
285101260	A 7780	120	12	3 x ceramic	A	-5 ... +80	plug head	gel electrolyte
285100207	H 61	170	12	platinum	H	+10 ... +100	plug head	
285092583	H 61-500	500	12	platinum	H	0...+100	plug head	
285092591	H 61-600	600	12	platinum	H	0...+100	plug head	
285102524	H 6180	170	12	ceramic	H	+10 ... +100	plug head	
285100215	H 62	120	12	platinum	H	+10 ... +100	plug head	
285102532	H 6280	120	12	ceramic	H	+10 ... +100	plug head	
285100223	H 63	320	12	platinum	H	+10 ... +100	plug head	
285102549	H 6380	320	12	ceramic	H	+10 ... +100	plug head	
285100231	H 64	170	12	ground joint	H	+10 ... +100	plug head	
285130220	H 64 1M-DIN-ID	170	12	ground joint	H	+10 ... +100	DIN plug ²⁾	ID function
285130230	H 64 1M-BNC-ID	170	12	ground joint	H	+10 ... +100	BNC plug ²⁾	ID function
285100248	H 65	103 ¹⁾	10	platinum	H	+10 ... +100	plug head	standard taper NS 14.5
285102565	H 6580	103 ¹⁾	10	ceramic	H	+10 ... +100	plug head	standard taper NS 14.5
1061093	L 32	120	12	fibre	A	-5 ... +50	plug head	plastic shaft
285101252	L 7780	120	12	ceramic	L	-5 ... +80	plug head	gel electrolyte
285101277	L 8280	120	12	KPG®	L	-5 ... +80	plug head	Referid® electrolyte
285100437	N 42 A	120	12	ceramic	A	-5 ... +100	DIN plug ²⁾	
285101544	N 42 BNC	120	12	ceramic	A	-5 ... +100	BNC plug ²⁾	
285100453	N 50 A	108	12	ceramic	A	-5 ... +100	DIN plug ²⁾	for portable Knick pH meters
285100494	N 52 A	120	12	platinum	A	-5 ... +100	DIN plug ²⁾	
285105451	N 52 BNC	120	12	platinum	A	-5 ... +100	BNC plug ²⁾	
285100001	N 61	170	12	platinum	A	-5 ... +100	plug head	
285100018	N 6180	170	12	ceramic	A	-5 ... +100	plug head	
285100034	N 62	120	12	platinum	A	-5 ... +100	plug head	
285100112	N 6250	120	12	ceramic	A	+15 ... +40	plug head	calomel ref., for TRIS buffers
285100042	N 6280	120	12	ceramic	A	-5 ... +100	plug head	
285100059	N 64	170	12	ground joint	A	-5 ... +100	plug head	
285100067	N 65	103 ¹⁾	10	platinum	A	-5 ... +100	plug head	standard taper NS 14.5
285102516	N 6580	103 ¹⁾	10	ceramic	A	-5 ... +100	plug head	standard taper NS 14.5
285101709	N 6980	103 ¹⁾	10	ground joint	A	-5 ... +100	plug head	standard taper NS 14.5
285092661	N 61 eis	170	12	3 x platinum	A	+10 ... +40	plug head	electrolyte L 5014, Ag/AgCl ref.
285092337	N 6480 eis	170	12	ground joint	A	+10 ... +40	plug head	electrolyte L 5014, Ag/AgCl ref.
285092329	N 6480 eth	170	12	ground joint	A	0 ... +40	plug head	electrolyte L 5014, Ag/AgCl ref.

¹⁾ Length from upper end of standard taper

²⁾ with 1 m fixed cable

ScienceLine pH combination electrodes with temperature sensor

pH combination electrodes with temperature sensor

Reference system: Silamid®
Shaft material: glass
Diameter: 12 mm
Zero point: pH = 7.0 ± 0.3
Electrolyte: KCl 3 mol/l
Temperature sensor: Pt 1000
Membrane shape: sphere
pH range: 0 ... 14
Connection cable:
for SMEK-plug head: e.g. LS 1 ANN
 (See also page with connection cables)
fixed cable: 1 m long, with plug A acc. to DIN 19262 or with BNC plug, as well as plug for temperature sensor



N 1042 A
 N 1041 A
 N 1041 BNC
 N 1042 BNC
 N 1050 A
 N 1051 A
 N 1051 BNC
 N 1052 A
 N 1052 BNC
 N 2041 A
 N 2042 A
 N 1041 A - 600
 N 1043 A

A 162
 A 161
 H 161
 H 162
 A 161 1M DIN ID
 A 161 1M BNC ID
 H 161 1M DIN ID
 H 161 1M BNC ID

A 164
 A 164 1M DIN ID
 A 164 1M BNC ID

A 7780 1M DIN ID
 A 7780 1M BNC ID

Sc

Order No.	Type No.	Length L [mm]	Dia- phragm	pH- glass	Temp.- range [°C]	Connection	Remarks
285129517	A 161	170	platinum	A	-5 ... +100	SMEK plug head	
285130240	A 161 1M-DIN-ID	170	platinum	A	-5 ... +100	DIN ¹⁾ - + 4-mm plug	ID function
285130250	A 161 1M-BNC-ID	170	platinum	A	-5 ... +100	BNC ¹⁾ - + 4-mm plug	ID function
285129525	A 162	120	platinum	A	-5 ... +100	SMEK plug head	
285129600	A 164	170	ground joint	A	-5 ... +100	SMEK plug head	
285130280	A 164 1M-DIN-ID	170	ground joint	A	-5 ... +100	DIN ¹⁾ - + 4-mm plug	ID function
285130290	A 164 1M-BNC-ID	170	ground joint	A	-5 ... +100	BNC ¹⁾ - + 4-mm plug	ID function
285130200	A 7780 1M-DIN-ID	120	3 x ceramic	A	-5 ... +80	DIN ¹⁾ + 4-mm plug	ID function
285130210	A 7780 1M-BNC-ID	120	3 x ceramic	A	-5 ... +80	BNC ¹⁾ + 4-mm plug	ID function
285129590	H 161	170	platinum	H	+10 ... +100	SMEK plug head	
285130260	H 161 1M-DIN-ID	170	platinum	H	+10 ... +100	DIN ¹⁾ - + 4-mm plug	ID function
285130270	H 161 1M-BNC-ID	170	platinum	H	+10 ... +100	BNC ¹⁾ - + 4-mm plug	ID function
285129580	H 162	120	platinum	H	+10 ... +100	SMEK plug head	
285100486	N 1041 A	170	ceramic	A	-5 ... +100	DIN ¹⁾ + 4-mm plug	
285093111	N 1041 A-600	600	ceramic	A	-5...+100	DIN ¹⁾ + 4-mm plug	Ag/AgCl ref.
285100531	N 1041 BNC	170	ceramic	A	-5 ... +100	BNC ¹⁾ + 4-mm plug	
285104541	N 1042 A	120	ceramic	A	-5 ... +100	DIN ¹⁾ + 4-mm plug	
285105476	N 1042 BNC	120	ceramic	A	-5 ... +100	BNC ¹⁾ + 4-mm plug	
285093009	N 1043 A	320	ceramic	A	-5...+100	DIN ¹⁾ + 4-mm plug	
285100375	N 1050 A	108	ceramic	A	-5 ... +100	DIN ¹⁾ + 4-mm plug	for portable Knick pH Meter
285100510	N 1051 A	170	platinum	A	-5 ... +100	DIN ¹⁾ + 4-mm plug	
285100500	N 1051 BNC	170	platinum	A	-5 ... +100	BNC ¹⁾ + 4-mm plug	
1054512	N 1052 A	120	platinum	A	-5 ... +100	DIN ¹⁾ + 4-mm plug	
285100380	N 1052 BNC	120	platinum	A	-5 ... +100	BNC ¹⁾ + 4-mm plug	
285100342	N 2041 A	170	ceramic	A	-5 ... +100	DIN ¹⁾ + 2-mm plug	
285100359	N 2042 A	120	ceramic	A	-5 ... +100	DIN ¹⁾ + 2-mm plug	

scienceLine

¹⁾ with 1 m fixed cable

ScienceLine micro, spear tip and surface pH combination electrodes

Micro, spear tip and surface pH combination electrodes

Reference system: Silamid®
Shaft material: glass
 (except L 39:
 plastic shaft)
Zero point: pH = 7.0 ± 0.3
Electrolyte: KCl 3 mol/l
 (except L8880:
 Referid®)
**Type of
 membrane glass:** A
**Connection cable:
 for SMEK plug head:** e.g. LS 1 ANN
 (See also page
 with connec-
 tion cables)
**for plug head
 versions:** e.g. L 1 A
 (See also page
 with connec-
 tion cables)
fixed cable: 1 m long, with
 plug A acc.
 to DIN 19262
 or with
 BNC plug, as
 well as plug
 for temperature
 sensor



A 157 1M
BNC ID
 A 157
 A 157 1M
 DIN ID

N 5800 A
 N 5800 BNC
 N 5900 A

N 6000 1M
DIN ID
 N 6000 1M
 BNC ID
 N 6000 A
 N 6000 BNC

N 6003

L 6880
 L 6880 1M-
 DIN-ID
 L 6880 1M-
 BNC-ID
 L 8880

N 1048 A
 N 1048 1M
 DIN ID
 N 1048 1M
 DIN ID
 N 48 A
 N 48 BNC

L 39
 L 39 1M
 DIN ID
 L 39 1M
 DIN ID

Order No.	Type No.	Length L [mm]	Ø [mm]	Dia- phragm	pH- glass	Membrane shape	Temp.- range [°C]	Range [pH]	Connection	Remarks
Micro										
285129610	A 157 ¹⁾	40/130	12/5	platinum	A	cylindrical	-5 ... +100	0 ... 14	SMEK plug head	
285130160	A 157 1M-DIN-ID ¹⁾	40/130	12/5	platinum	A	cylindrical	-5 ... +100	0 ... 14	DIN plug ³⁾	ID function
285130170	A 157 1M-BNC-ID ¹⁾	40/130	12/5	platinum	A	cylindrical	-5 ... +100	0 ... 14	BNC plug ³⁾	ID function
285105127	N 5800 A	96 ²⁾	5	3 x platinum	A	spear	-5 ... +100	0 ... 14	DIN plug ³⁾	Ag/AgCl ref.
285105579	N 5800 BNC	96 ²⁾	5	3 x platinum	A	spear	-5 ... +100	0 ... 14	BNC plug ³⁾	Ag/AgCl ref.
285105135	N 5900 A	96 ²⁾	5	platinum	A	sphere	-5 ... +100	0 ... 14	DIN plug ³⁾	Ag/AgCl ref.
285105151	N 6000 A	96 ²⁾	3	platinum	A	cylindrical	-5 ... +100	0 ... 14	DIN plug ³⁾	Ag/AgCl ref.
285105632	N 6000 BNC	96 ²⁾	3	platinum	A	cylindrical	-5 ... +100	0 ... 14	BNC plug ³⁾	Ag/AgCl ref.
285130180	N 6000 1M-DIN-ID	96 ²⁾	3	platinum	A	cylindrical	-5 ... +100	0 ... 14	DIN plug ³⁾	ID function, Ag/AgCl ref.
285130190	N 6000 1M-BNC-ID	96 ²⁾	3	platinum	A	cylindrical	-5 ... +100	0 ... 14	BNC plug ³⁾	ID function, Ag/AgCl ref.
285105176	N 6003	180	3	ceramic	A	cylindrical	-5 ... +100	0 ... 14	plug head	Ag/AgCl ref.
Spear tip										
285101211	L 6880	70/50	12/8	3 x ceramic	A	spear	-5 ... +100	0 ... 14	plug head	
285130100	L 6880 1M-DIN-ID	70/50	12/8	3 x ceramic	A	spear	-5 ... +100	0 ... 14	DIN plug ³⁾	ID function
285130110	L 6880 1M-BNC-ID	70/50	12/8	3 x ceramic	A	spear	-5 ... +100	0 ... 14	BNC plug ³⁾	ID function
285101285	L 8880	70/50	12/8	hole	A	spear	-5 ... +80	2 ... 13	plug head	
285104611	N 1048 A ¹⁾	120	12	ceramic	A	spear	-5 ... +100	0 ... 14	DIN- ³⁾ + 4-mm plug	
285130120	N 1048 1M-DIN-ID ¹⁾	120	12	ceramic	A	spear	-5 ... +100	0 ... 14	DIN- ³⁾ + 4-mm plug	ID function
285130130	N 1048 1M-BNC-ID ¹⁾	120	12	ceramic	A	spear	-5 ... +100	0 ... 14	BNC- ³⁾ + 4-mm plug	ID function
285100445	N 48 A	120	12	ceramic	A	spear	-5 ... +100	0 ... 14	DIN plug ³⁾	
285101569	N 48 BNC	120	12	ceramic	A	spear	-5 ... +100	0 ... 14	BNC plug ³⁾	
Surface										
1061094	L 39	120	12	fibre	A	flat	-5 ... +50	1 ... 13	plug head	
285130140	L 39 1M-DIN-ID	120	12	fibre	A	flat	-5 ... +50	1 ... 13	DIN plug ³⁾	ID function
285130150	L 39 1M-BNC-ID	120	12	fibre	A	flat	-5 ... +50	1 ... 13	BNC plug ³⁾	ID function

¹⁾ with integrated temperature sensor Pt 1000

²⁾ Length from upper end of standard taper (Standard taper NS 7.5)

³⁾ with 1 m fixed cable

ScienceLine metal combination electrodes

Metal combination electrodes with plug head and connection cable

Temperature range: -5 ... +100 °C
(except Pt 6140:
+10 ... +40 °C)

Reference system: Silamid®

Shaft material: glass

Electrolyte: KCl 3 mol/l
(See also remarks)

Connection cable:
for plug head: e.g. L 1 A
(See also page
with connection
cables)

fixed cable: 1 m long, with
plug A acc. to
DIN 19262 or
with BNC plug



AgCl 62
AgCl 65
Ag 42 A
Ag 6180
Ag 6280
Ag 6580
AgCl 6280
Au 6280

Pt 61
Pt 62
Pt 6180
Pt 6280
Pt 6580
Pt 42 A

Pt 6880
Pt 6980
Pt 48 A

Pt 6140

Pt 8280

Pt 5900 A
Pt 5900 BNC
Pt 5901

Order No.	Type No.	Length L [mm]	Dia- phragm	Ø [mm]	Sensor Metal, shape	Connection	Remarks
285102051	Ag 42 A	120	ceramic	12	Ag, cap, 5 mm Ø	DIN plug ⁴⁾	electrolyte L 2114, Ag/AgCl ref.
285102208	Ag 6180	170	ceramic	12	Ag, cap, 5 mm Ø	plug head	electrolyte L 2114, Ag/AgCl ref.
285102343	Ag 6280	120	ceramic	12	Ag, cap, 5 mm Ø	plug head	electrolyte L 2114, Ag/AgCl ref.
285102216	Ag 6580	103 ¹⁾	ceramic	10	Ag, cap, 5 mm Ø	plug head	electrolyte L 2114, Ag/AgCl ref.
285102351	AgCl 6280 ³⁾	120	ceramic	12	Ag, cap, 5 mm Ø	plug head	electrolyte L 2114, Ag/AgCl ref.
285102413	AgCl 62 ³⁾	120	platinum	12	Ag, cap, 5 mm Ø	plug head	electrolyte L 2114, Ag/AgCl ref.
1061051	AgCl 65 ³⁾	103 ¹⁾	platinum	12	Ag, cap, 5 mm Ø	plug head	electrolyte L 2114, Ag/AgCl ref.
285102121	Au 6280	120	ceramic	12	Au, pole, 2 mm Ø	plug head	
285102302	Pt 42 A	120	ceramic	12	Pt, pole, 1 mm Ø	DIN plug ⁴⁾	
285102224	Pt 48 A	120	ceramic	12	Pt, ring, 6 mm Ø	DIN plug ⁴⁾	Ag/AgCl ref.
285105192	Pt 5900 A	96 ²⁾	platinum	5	Pt, pole, 1 mm Ø	DIN plug ⁴⁾	Ag/AgCl ref.
285105702	Pt 5900 BNC	96 ²⁾	platinum	5	Pt, pole, 1 mm Ø	BNC plug ⁴⁾	Ag/AgCl ref.
285105065	Pt 5901	160 ²⁾	platinum	5	Pt, pole, 1 mm Ø	plug head	
285102002	Pt 61	170	platinum	12	Pt, pole, 1 mm Ø	plug head	
285102019	Pt 62	120	platinum	12	Pt, pole, 1 mm Ø	plug head	
285097162	Pt 6140	150/20	platinum	12/5	Pt, pole, 1 mm Ø	plug head	for spear tip, electrolyte L420
285102232	Pt 6180	170	ceramic	12	Pt, pole, 1 mm Ø	plug head	
285102249	Pt 6280	120	ceramic	12	Pt, pole, 1 mm Ø	plug head	
285102257	Pt 6580	103 ¹⁾	ceramic	10	Pt, pole, 1 mm Ø	plug head	
285100075	Pt 6880	120	ceramic	12	Pt, ring, 6 mm Ø	plug head	
285102265	Pt 6980	170	ceramic	12	Pt, ring, 6 mm Ø	plug head	
285102281	Pt 8280	120	KPC®	12	Pt, round, 6 mm Ø	plug head	electrolyte Referid®

¹⁾ Length from upper end of standard taper; standard taper NS 14.5

²⁾ Length from upper end of standard taper; standard taper NS 7.5

³⁾ Sensor coated with AgCl

⁴⁾ with 1 m fixed cable

ScienceLine single electrodes: pH glass electrodes and metal electrodes

ScienceLine single electrodes

pH glass electrodes

Reference system: Silamid®
Shaft material: glass, 12 mm Ø
Zero point: pH = 7.0 ± 0.3
Membrane shape: sphere
Connection cable: e.g. L 1 A

Metal electrodes

Shaft material: glass, 12 mm Ø
(See remarks)



A 1180
H 1180

Ag 1100

KF 1100

Pt 1400
Pt 1200

Pt 1800

Sc

Order No.	Type No.	Length L [mm]	pH Glass	Range [pH]	Temp.- range [°C]	Remarks
1057997	A 1180 ¹⁾	120	H	0 ... 14	0 ... +80	plug head
285103212	H 1180	120	H	0 ... 14	10 ... +100	plug head

Order No.	Type No.	Length L [mm]	Sensor Metal	Sensor shape	Temp. range [°C]	Remarks
285103607	Ag 1100	120	Ag	cap, 4 mm Ø	-5 ... +100	plug head, cable e.g. L 1 A
285102030	KF 1100	96 ¹⁾	Pt ²⁾	2 pole, 1 mm Ø	-30 ... +135	shaft 5 mm Ø, standard taper NS 7.5, fixed cable, 2x 4-mm plug
285103512	Pt 1200	120	Pt ²⁾	2 pole, 1 mm Ø	-30 ... +135	plug head, cable e.g. L 1 NN
285103537	Pt 1400	103 ¹⁾	Pt ²⁾	2 pole, 1 mm Ø	-30 ... +135	shaft 10 mm Ø, standard taper NS 14.5, cable e.g. L 1 NN
285103553	Pt 1800	120	Pt	ring, 6 mm Ø	-30 ... +135	plug head, cable e.g. L 1 A

ScienceLine

¹⁾ Length from upper end of standard taper

²⁾ Double platinum electrode

ScienceLine single electrodes: Reference electrodes

Reference electrodes

Shaft material: glass

**Electrolyte depending on
reference system:**

Ag/AgCl: KCl 3 mol/l,
e.g. L 300

Calomel: KCl 4.2 mol/l,
e.g. L 420

Hg/Hg₂SO₄: K₂SO₄ 0.6 mol/l,
e.g. L 1254

pH range: 0 ... 14

Connection cable: e.g. L 1 N



B 2220+



B 2420+



B 2810+
B 2820+
B 2910+
B 2920+



B 3420+
B 3410+
B 3510+
B 3520+
B 3610+



B 3920+

Sc

Order No.	Type No.	Length L [mm]	Ø [mm]	Temp. range [°C]	Dia- phragm	Reference system	Remarks
1069994	B 2220+	120	12	-5 ... +100	platinum	Ag/AgCl	double electrolyte system
1070028	B 2420+	120	12	-5 ... +100	ground joint	Ag/AgCl	
1070029	B 2810+	120	12	+15 ... +40	ceramic	Calomel	
1070044	B 2820+	120	12	-5 ... +100	ceramic	Ag/AgCl	
1070077	B 2910+	120	12	+15 ... +40	platinum	Calomel	
1070046	B 2920+	120	12	-5 ... +100	platinum	Ag/AgCl	
1070048	B 3410+	103 ¹⁾	10	+15 ... +40	ceramic	Calomel	standard taper NS 14.5
1070070	B 3420+	103 ¹⁾	10	-5 ... +100	ceramic	Ag/AgCl	standard taper NS 14.5
1070100	B 3510+	103 ¹⁾	10	+15 ... +40	platinum	Calomel	standard taper NS 14.5
1070073	B 3520+	103 ¹⁾	10	-5 ... +100	platinum	Ag/AgCl	standard taper NS 14.5
1070074	B 3610+	103 ¹⁾	10	+15 ... +40	ceramic	Hg/Hg ₂ SO ₄	standard taper NS 14.5
1070075	B 3920+	103 ¹⁾	10	-5 ... +100	ground joint	Ag/AgCl	double electrolyte system, standard taper NS 14.5

ScienceLine

¹⁾ Length from upper end of standard taper

ScienceLine conductivity measuring cells with fixed cable

Conductivity measuring cells with fixed cable and 8-pole plug

Shaft: 12 mm Ø
(except LF 413
T-3 and LF 413 T:
15.3 mm)
Temperature sensor: NTC 30 kΩ



LF 213 T
LF 213 T ID

LF 313 T NTEC
LF 313 T
LF 313 T ID

LF 413 T-3
LF 413 T
LF 413 T ID

LF 513 T
LF 613 T
LF 813 T

LF 713 T
LF 713 T-250

LF 913 T
LF 913 T ID

LFOX 1400
LFOX 1400 ID

Sc

Order No.	Type No.	Length L [mm]	Ø [mm]	Sensor	Cell const. approx. [cm ⁻¹]	Temp. range [°C]	Meas. range ¹⁾ [µS/cm] . . [mS/cm]	Remarks
285106150	LF 213 T	120	12	Stainless steel	0.01	0 ... +100	0 ... 0.03	Trace conductivity cell with integrated flow-through vessel, stainless steel, 1.5 m cable
285106160	LF 213 T ID	120	12	Stainless steel	0.01	0 ... +100	0 ... 0.03	Trace conductivity cell with integrated flow-through vessel, stainless steel, 1.5 m cable, ID function
285414360	LF 313 T	120	12	Stainless steel	0.1	0 ... +100	0 ... 0.2	Ultrapure water conductivity cell with flow-through vessel, stainless steel shaft, fixed cable 1.5 m
285130300	LF 313 T-ID	120	12	Stainless steel	0.1	0...+100	0 ... 0.2	Ultrapure water conductivity cell with flow-through vessel, stainless steel shaft, cable 1.5 m, ID function
285414351	LF 313 T NFTC	120	12	Stainless steel	0.1	0 ... +100	0 ... 0.2	Ultrapure water conductivity cell without flow-through vessel, stainless steel shaft, fixed cable 1.5 m
285106172	LF 413 T	120	15.3	4 x Graphite	0.475	-5 ... +80	1 ... 2000	Plastic shaft, 1.5 m cable
285130310	LF 413 T-ID	120	15.3	4 x Graphite	0.475	-5 ... +80	0 ... 2000	Plastic shaft, 1.5 m cable, ID function
285106148	LF 413 T-3	120	15.3	4 x Graphite	0.475	-5 ... +80	1 ... 2000	Plastic shaft, fixed cable 3 m
285106037	LF 513 T	120	12	2 Pt rings	1.0	-5 ... +80	1 ... 200	Plastic shaft, 1 m cable
285106131	LF 613 T	120	12	4 Pt rings	1.0	-5 ... +80	1 ... 2000	Plastic shaft, 1 m cable
285106189	LF 713 T	120	12	4 Pt rings	1.0	-30 ... +135	1 ... 2000	Glass shaft, 1 m cable
285106190	LF 713 T-250	250	12	4 Pt rings	1.0	-30 ... +135	1 ... 2000	Glass shaft, 1 m cable
285106250	LF 813 T	120	12	5 Pt rings	0.650	-5 ... +80	1 ... 2000	Plastic shaft, 1 m cable
285106260	LF 913 T	120	12	5 Pt rings	0.650	-30 ... +135	1 ... 2000	Glass shaft, 1 m cable
285130320	LF 913 T-ID	120	12	5 Pt rings	0.650	-30 ... +135	1 ... 2000	5-pole cell, glass shaft, 1 m cable, ID function
285104630	LFOX 1400	145	15.3	Graphite	0.475	0 ... +50	1 ... 2000	Combined 4-pole conductivity cell and galvanic D.O. sensor LFOX 1400 ID, plastic shaft, fixed cable 3 m
285130330	LFOX 1400 ID	145	15.3	Graphite	0.475	0...+50	1 ... 2000	Combined 4-pole conductivity cell and galvanic D.O. sensor LFOX 1400 ID, plastic shaft, fixed cable 3 m, ID function

¹⁾ Outside the recommended ranges measuring errors >10% can occur with these conductivity measuring cells.

ScienceLine conductivity measuring cells with plug head

Conductivity measuring cells
with plug head

Shaft: 12 mm Ø



Order No.	Type No.	Length L [mm]	Ø [mm]	Sensor	Cell const. approx. [cm ⁻¹]	Temp. range [°C]	Meas. range ¹⁾ [μS/cm]. . .[mS/cm]	Remarks
1069976	LF 1100+	120	12	2 Pt plates	1.0	-30 ... 135	0 ... 200	SMEK plug head
1069977	LF 1100T+	120	12	2 Pt plates	1.0	-30 ... 135	0 ... 200	SMEK plug head
1069978	LF 4100+	120	12	2 Pt plates	1.0	-30 ... 135	0 ... 200	SMEK plug head, flow-through cell
1069979	LF 5100+	120	12	2 Pt rings	1.0	-5 ... 80	0 ... 200	SMEK plug head, plastic shaft
1069990	LF 5100T+	120	12	2 Pt rings	1.0	-5 ... 80	0 ... 200	SMEK plug head, plastic shaft

ScienceLine

¹⁾ Outside the recommended ranges measuring errors >10% can occur with these conductivity measuring cells.

ScienceLine sensors for ammonia, sodium, oxygen, ion-selective indicator electrodes

Ammonia combination electrode with plug head

Shaft material: plastic, 12 mm Ø
Connection cable: e.g. L 1 A

Sodium combination electrode with plug head

Reference system: Silamid®
Shaft material: glass, 12 mm Ø
Zero point: pNa = 2.0
Membrane shape: sphere
Connection cable: e.g. L 1 A

Oxygen electrodes

Shaft material: plastic (POM)

ISE measuring cells

Shaft material: plastic
Length: 120 mm
Fixed cable: 1 m long,
with DIN plug

ISE combination electrodes with plug head

Shaft material: plastic
Length: 120 mm



NH 1100

Na 61

OX 1100+

9009/61

Cu 1100 A

F 60

Ca 1100 A

Cl 60

F 1100 A

NO 60

Pb 1100 A

K 60

CA 60

CN 60

AG-S 60

I 60

BR 60

CU 60

PB 60

ScienceLine

Order No.	Type No.	Length L [mm]	Temp. range [°C]	Meas. range [mg/l]	Remarks
285102808	NH 1100	120	0 ... +50	0.1 ... 1,000	membrane module replaceable

Order No.	Type No.	Length L [mm]	Diaphragm	Membrane Glass	Temp. range [°C]	Meas. range [pNa]	Remarks
285100026	Na 61	170	platinum	Na	-10 ... +80	0 ... 6	electrolyte KCl 3 mol/l, aqueous solution NaCl 0.1 mol/l

Order No.	Type No.	Length L [mm]	Temp. range [°C]	Meas. range [mg/l]	Remarks
1069975	OX 1100+	120	0 ... +45	0 ... 60	galvanic sensor, Pt cathode, Ag anode, SMEK plug head, temperature compensated (NTC 100kΩ), shaft 12 mm Ø, measuring current at saturation approx. 100 nA, minimum flow rate 10 cm/s, connection cable e.g. LS 1 ST4 OX (for CG 867)
285111664	9009/61	145	0 ... +50	0 ... 50	amperometric sensor, Au cathode, Pb anode, fixed cable 1.5 m ¹⁾ with 8-pole plug, IMT temperature compensation, shaft 15.25 mm Ø, membrane FEP, 13 µm thick, accuracy 1% at 18 cm/s flow rate.

Order No.	Type No.	Parameter	Temp. range [°C]	pH-range	Measuring range [mg/l]
285216314	Ca 1100 A	Calcium	0 ... +40	2.5 ... 11	0.02 ... 40,000
285216312	Cu 1100 A	Copper	0 ... +80	2 ... 6	0.0006 ... 6,400
285216313	F 1100 A	Fluoride	0 ... +80	5 ... 7	0.02 ... saturated
285216315	Pb 1100 A	Lead	0 ... +80	4 ... 7	0.1 ... 20,000

Order No.	Type No.	Parameter	Temp. range [°C]	pH-range	Measuring range [mg/l]
285130340	F 60	Fluoride	0 ... +80	5 ... 7	0.02 ... saturated
285130350	Cl 60	Chloride	0 ... +80	2 ... 12	2 ... 35,000
285130360	NO 60	Nitrate	0 ... +40	2.5 ... 11	0.4 ... 62,000
285130370	K 60	Potassium	0 ... +40	2 ... 12	0.04 ... 39,000
285130380	CA 60	Calcium	0 ... +40	2.5 ... 11	0.02 ... 40,000
285130390	CN 60	Cyanide	0 ... +80	0 ... 14	0.2 ... 260
285130400	AG-S 60	Sulfide	0 ... +80	2 ... 12	0.003 ... 32,000
285130410	I 60	Iodide	0 ... +80	0 ... 14	0.006 ... 127,000
285130420	BR 60	Bromide	0 ... +80	1 ... 12	0.4 ... 79,000
285130430	CU 60	Copper	0 ... +80	2 ... 6	0.0006 ... 6400
285130440	PB 60	Lead	0 ... +80	4 ... 7	0.2 ... 20,000

¹⁾ Other cable lengths available on request

Resistance thermometers

Resistance thermometers
with SMEK plug head

Resistance thermometers
with 1 m fixed cable

Resistance thermometer
with coaxial plug head



W 2030+
W 2130+

W 5780 NN

W 5790 NN
W 5790 PP
W 5791 NN

W 5980 NN

W 2180-KOAX

Science

Resistance thermometers with SMEK plug head

Order No.	Type No.	Length L [mm]	Ø [mm]	Sensor	Temp. range [°C]	Shaft material	Connection cable e.g.
1069991	W 2030+	120	12	Pt 100	-30 ... +135	glass	LS 1 N6
1069992	W 2130+	120	12	Pt 1000	-30 ... +135	glass	LS 1 N6

Resistance thermometers with 1 m fixed cable

Order No.	Type No.	Length L [mm]	Ø [mm]	Sensor	Temp. range [°C]	Shaft material	Connection plug
285105221	W 5780 NN	120	6	Pt 1000	-30 ... +135	glass	2 x 4 mm Ø
285105254	W 5790 NN	120	4	Pt 1000	-30 ... +135	stainless steel	2 x 4 mm Ø
285105776	W 5790 PP	120	4	Pt 1000	-30 ... +135	stainless steel	2 x 2 mm Ø
285105262	W 5791 NN	170	4	Pt 1000	-30 ... +135	stainless steel	2 x 4 mm Ø
285105287	W 5980 NN	96 ¹⁾	5 NS 7.5	Pt 1000	-30 ... +135	glass	2 x 4 mm Ø

Resistance thermometer with coaxial plug head

Order No.	Type No.	Length L [mm]	Ø [mm]	Sensor	Temp. range [°C]	Shaft material
285119030	W 2180-KOAX	120	12	Pt 1000	-30 ... +135	glass

eLine

¹⁾ length from upper end of standard taper

BlueLine pH combination electrodes

The robust electrodes for general applications

pH range	0...14
Temperature range	-5...+80 °C
Shaft	Noryl, 12 mm Ø
Shaft length L	120 mm
Zero point	pH = 7.0 ± 0.3
Diaphragm	fibre
Reference system	Ag/AgCl
Reference electrolyte	gel (KCl), low maintenance, not refillable
Shape of glass membrane	cylindrical
Resistance of glass membrane (25 °C)	400 MΩ
Type of membrane glass	A

The liquid electrolyte electrodes for demanding measurements

pH range	0...14
Temperature range	-5...+100 °C
Shaft	glass, 12 mm Ø
Shaft length L L	120 mm
Zero point	pH = 7.0 ± 0.3
Diaphragm	platinum
Reference system	Ag/AgCl
Reference electrolyte	KCl 3 mol/l
Shape of glass membrane	conical
Resistance of glass membrane (25 °C)	200 MΩ
Type of membrane glass	A



BlueLine
28 pH
22 pH
23 pH
23-2 pH
23-5 pH-S
24 pH
24-3 pH
25 pH
25-2 pH
25-5 pH
26 pH
26 pH-Cinch
28 pH-P
28-5 pH
29 pH
29 pH-P

BlueLine
18 pH
11 pH
12 pH
14 pH
14 pH ID
15 pH
15 pH ID
15 pH Cinch
17 pH
17 pH-R
19 pH

BlueLine

Order No.	BlueLine Type No.	Temperature sensor integrated	Connection
285129225	22 pH	no	plug head, recommended cable: e.g. LB1A
285129233	23 pH	no	1 m fixed cable with DIN plug 19 262
1063462	23-2 pH	no	2 m fixed cable with DIN plug
1066411	23-5 pH-S	no	5 m fixed cable with S plug
285129241	24 pH	NTC 30 k Ω	1 m fixed cable with DIN plug 19 262 + banana plug
285129533	24-3 pH	NTC 30 k Ω	3 m fixed cable with DIN plug 19 262 + banana plug
285129258	25 pH	no	1 m fixed cable with BNC plug
1063461	25-2 pH	no	2 m fixed cable with BNC plug
285129540	25-5 pH	no	5 m fixed cable with BNC plug
285129266	26 pH	NTC 30 k Ω	1 m fixed cable with BNC plug + banana plug
285095712	26 pH-Cinch	NTC 30 k Ω	1 m fixed cable with BNC plug + cinch plug
285129282	28 pH	Pt 1000	1 m fixed cable with DIN plug 19 262 + banana plug
1065896	28 pH-P	Pt 1000	1 m fixed cable with DIN plug 19 262 + 2-mm pole plug
285129570	28-5 pH	Pt 1000	5 m fixed cable with DIN plug 19 262 + banana plug
1065895	29 pH	Pt 1000	1 m fixed cable with BNC plug + banana plug
1065894	29 pH-P	Pt 1000	1 m fixed cable with BNC plug + 2-mm pole plug

Order No.	BlueLine Type No.	Temperature sensor integrated	Connection
285129114	11 pH	no	plug head, recommended cable: e.g. LB1A
285129122	12 pH	no	1 m fixed cable with DIN plug 19 262
285129147	14 pH	NTC 30 k Ω	1 m fixed cable with DIN plug 19 262 + banana plug
285129440	14 pH ID	NTC 30 k Ω	1 m fixed cable with DIN plug + 4-mm banana plug, ID function
285129155	15 pH	NTC 30 k Ω	1 m fixed cable with BNC plug + banana plug
285129450	15 pH ID	NTC 30 k Ω	1 m fixed cable with BNC plug + 4-mm banana plug, ID function
285095730	15 pH Cinch	NTC 30 k Ω	1 m fixed cable with BNC plug + cinch plug
285129171	17 pH	no	1 m fixed cable with BNC plug
1064746	17 pH-R	no	1 m fixed cable with Metrohm plug
285129188	18 pH	Pt 1000	1 m fixed cable with DIN plug 19 262 + banana plug
285129190	19 pH	Pt 1000	1 m fixed cable with BNC plug + banana plug

ueLine

BlueLine

Special sensors

The specialists
for special applications

Zero point of
pH electrodes

pH = 7.0 ± 0.3

Connection cable
for pH/Redox
electrodes

e.g. LB 1 A



BlueLine
13 pH

BlueLine
16 pH

BlueLine
21 pH

21 pH 1M
DIN ID

21 pH 1M
BNC ID

BlueLine
27 pH

27 pH 1M
DIN ID

27 pH 1M
BNC ID

BlueLine
31 Rx

BlueLine
32 Rx

BlueLine
48 LF

Precision electrode BlueLine 13 pH	Glass shaft, screw ground joint diaphragm, electrolyte KCl 3 mol/l, Ag/AgCl reference system, spherical membrane, A-glass, plug head, length 170 mm, 12 mm Ø, -5 ... +100 °C, 0 ... 14pH, Order No. 285129139
Micro electrode BlueLine 16 pH	Glass shaft, platinum diaphragm, electrolyte KCl 3 mol/l, Ag/AgCl reference system, spherical membrane, L-glass, plug head, length 40/80 mm, 12/5 mm Ø, -5 ... +100 °C, 0 ... 14pH, Order No. 285129163
Spear tip electrode BlueLine 21 pH	Glass shaft, hole diaphragm, Referid® electrolyte, Ag/AgCl reference system, Spear membrane, L-glass, plug head, length 65/25 mm, 12/5 mm Ø, -5 ... +80 °C, 2 ... 13pH, Order No. 285129217
Spear tip electrode with sensor recognition BlueLine 21 pH 1M-DIN-ID	Glass shaft, hole diaphragm, Referid® electrolyte, Ag/AgCl reference system, Spear membrane, L-glass, 1 m fixed cable with DIN plug and sensor recognition, length 65/25 mm, 12/5 mm Ø, -5...+80 °C, 2...13 pH Order No. 285129930
Spear tip electrode with sensor recognition BlueLine 21 pH 1M-BNC-ID	Glass shaft, hole diaphragm, Referid® electrolyte, Ag/AgCl reference system, Spear membrane, L-glass, 1 m fixed cable with BNC plug and sensor recognition, length 65/25 mm, 12/5 mm Ø, -5...+80 °C, 2...13 pH Order No. 285129940
Surface electrode BlueLine 27 pH	Glass shaft, KPG® annular gap diaphragm, Referid® electrolyte, Ag/AgCl reference system, flat membrane, L-glass, plug head, length 120 mm, 12 mm Ø, -5 ... +50 °C, 2 ... 13pH, Order No. 285129274
Surface electrode with sensor recognition BlueLine 27 pH 1M-DIN-ID	Glass shaft, KPG® annular gap diaphragm, Referid® electrolyte, Ag/AgCl reference system, flat membrane, L-glass, 1 m fixed cable with DIN plug and sensor recognition, length 120 mm, 12 mm Ø, -5...+50 °C, 2...13 pH Order No. 285129950
Surface electrode with sensor recognition BlueLine 27 pH 1M-BNC-ID	Glass shaft, KPG® annular gap diaphragm, Referid® electrolyte, Ag/AgCl reference system, flat membrane, L-glass, 1 m fixed cable with BNC plug and sensor recognition, length 120 mm, 12 mm Ø, -5...+50 °C, 2...13 pH Order No. 285129960
Redox electrode BlueLine 31 Rx	Glass shaft, ceramic diaphragm, electrolyte KCl 3 mol/l, Ag/AgCl reference system, sensor platinum disk 4 mm Ø, plug head, length 120 mm, 12 mm Ø, -5 ... +100 °C, Order No. 285129311
Redox electrode BlueLine 32 Rx	Plastic shaft, fibre diaphragm, gel electrolyte, Ag/AgCl reference system, sensor platinum pin 1 mm Ø, plug head, length 120 mm, 12 mm Ø, -5 ... +80 °C, Order No. 285129320
Conductivity cell for low ionic media BlueLine 48 LF	Stainless steel shaft, 2-pin cell, 1 m fixed cable with 8-pole plug, sensor stainless steel, cell constant 0.1 cm ⁻¹ , temperature sensor NTC 30 kΩ, length 120 mm, 12 mm Ø, -5 ... +80 °C, measuring range 0 ... 300 μS/cm, Order No. 285129488

Connection cables



① Electrode socket/plug

Coaxial plug for pH, redox, ammonia and sodium combination electrodes, pH and redox single electrodes as well as reference electrodes in *Plus* series

plug L



SMEK plug for pH combination electrodes with temperature sensor as well as conductivity measuring cells, resistance thermometers and oxygen sensors from *Plus* series

plug LS



Electrode plug for reference electrodes from the predecessor series, i.e. "non-*Plus*" versions

plug B



Plug for resistance thermometers in conductivity measuring cells without temperature sensor, for older models

plug 9907/00



Plug for conductivity measuring cells with temperature sensor and oxygen cells, for older models

plug 9909/00



② Instrument connector/plug

A (DIN)



BNC



EE (Radiometer)



R (Metrohm)



S (UK socket without extension)



N (4-mm banana plug)



P (2-mm pole plug)



8-pole (for Handylab and Lab and ProLab conductometer)



9910/00



Not illustrated:

X (without instrument plug, free cable end)

The connecting cables are available in various combinations, existing of electrode plug, instrument connection and cable length. Should you i.e. require a coaxial cable for connecting a pH electrode to a meter, please select a cable type L 1 A or L 2 A. The „L“ as part of the type description stands for the coaxial plug (plug L, please refer to page 82) of the electrode, the middle number stands for the cable length and the „A“ for the instrument connection (in this example for a DIN connection).

In case you do not find your desired cable combination listed below, please contact us.

Order No.	Type No.	① Electrode socket/plug	② Instrument connector/plug	Cable length and type
285122904	A 1 A	DIN instrument plug (A)	DIN instrument plug (A)	1 m coax. cable
285123793	A 1 BNC	DIN instrument plug (A)	BNC instrument plug	1 m coax. cable
285121916	B 1 N	reference electrode plug (B)	4 mm banana plug (N)	1 m single conductor cable
285122012	B 1 P	reference electrode plug (B)	2 mm Pole plug (P)	1 m single conductor cable
285121813	B 1 X	reference electrode plug (B)	free end (X)	1 m single conductor cable
285122456	L 1 A	electrode plug (L)	DIN instrument plug (A)	1 m coax. cable
285122497	L 1 BNC	electrode plug (L)	BNC instrument plug	1 m coax. cable
285122501	L 1 EE	electrode plug (L)	Radiometer instrument plug (EE)	1 m coax. cable
285122457	L 1 N	electrode plug (L)	4 mm banana plug (N)	1 m coax. cable
285122489	L 1 NN	electrode plug (L)	2 x 4 mm banana plug (N)	1 m coax. cable
285122534	L 1 R	electrode plug (L)	Metrohm instrument plug (R)	1 m coax. cable
285122407	L 1 X	electrode plug (L)	free end (X)	1 m coax. cable
285122464	L 2 A	electrode plug (L)	DIN instrument plug (A)	2 m coax. cable
285122448	L 2 NN	electrode plug (L)	2 x 4 mm banana plug (N)	2 m coax. cable
285122653	LB 1 A	electrode plug (LB)	DIN instrument plug (A)	1 m coax. cable
285122661	LB 1 BNC	electrode plug (LB)	BNC instrument plug	1 m coax. cable
285122678	LB 3 A	electrode plug (LB)	DIN instrument plug (A)	3 m coax. cable
285122707	LS 1 ANN	SMEK electrode plug	DIN (A) + 2 x 4 mm banana plug (N)	1 m cable KA19
285122715	LS 3 ANN	SMEK electrode plug	DIN (A) + 2 x 4 mm banana plug (N)	3 m cable KA19
285122723	LS 1 BNCNN	SMEK electrode plug	BNC + 2 x 4 mm banana plug (N)	1 m cable KA19
285122731	LS 3 BNCNN	SMEK electrode plug	BNC + 2 x 4 mm banana plug (N)	3 m cable KA19
1066726	LS 1 D8	SMEK electrode plug	8-pole instrument plug	1 m cable
1066728	LS 1 N6	SMEK electrode plug	6 x 4 mm banana plug (N)	1 m cable KA09
285122756	LS 1 RNN	SMEK electrode plug	Metrohm (R) + 2 x 4 mm banana plug (N)	1 m cable KA19
1069104	LS 1 ST4LF	SMEK electrode plug	4-pole incremental plug	1 m cable
1066727	LS 1 ST4OX	SMEK electrode plug	4-pole incremental plug	1 m cable KA10
285124716	9907/21	electrode plug (9907/00)	2 x 4-mm plug (N) for LF cells	1 m two-conductor cable
285125618	9909/31	electrode plug (9907/00)	2 x 4-mm plug (N)	1 m two-conductor cable
285125515	9910/11	electrode plug (9909/00)	9910	1 m four-conductor cable
285125215	9910/21	electrode plug (9909/00)	9910	1 m four-conductor cable, shielded
285125523	9919/21	electrode plug (9907/00)	8-pole instrument plug	1 m two-conductor cable
285125548	9919/41	electrode plug (9909/00)	8-pole instrument plug	1 m four-conductor cable

Other plug/cable combinations available on request


Solutions

Buffer solutions in the unique double-end ampoules offer a particularly high degree of reliability and measuring accuracy.

The exactness of the pH measurement is mainly dependent on the accuracy of calibration. This again highly depends on the reliability of the buffer.

Hermetically sealed in the glass ampoule and sterilized with hot steam, same as a pharmaceutical product, the buffer solutions free of preservation agent have an extremely long shelf life and guarantee continuously error-free characteristics.

The ampoules can be easily opened at the breaking point. Tools are not required. Since refilling is not possible, you are always ensured of maximum calibration reliability.




Standard buffer solutions according to DIN 19 266

Hot steam sterilized for longer stability, no preservation agents used.

Order No.	Type No.	pH value at 25 °C	Contents
285137977	L 4791	1.68	60 FIOLAX® ampoules à 20 ml*, with manufacturer's certificate
285138246	L 4794	4.01	60 FIOLAX® ampoules à 20 ml*, with manufacturer's certificate
285138254	L 4796	6.87	60 FIOLAX® ampoules à 20 ml*, with manufacturer's certificate
285138262	L 4799	9.18	60 FIOLAX® ampoules à 20 ml*, with manufacturer's certificate
285138402	L 4790	4.01/6.87	2 x 30 FIOLAX® ampoules à 20 ml*, with manufacturer's certificate
285137985	L 4797	1.68/6.87/9.18	3 x 20 FIOLAX® ampoules à 20 ml*, with manufacturer's certificate
285138238	L 4798	4.01/6.87/9.18	3 x 20 FIOLAX® ampoules à 20 ml*, with manufacturer's certificate
285138279	L 4893/Set	4.01/6.87	2 x 9 FIOLAX® ampoules à 20 ml*, with manufacturer's certificate, with electrolyte solution L 3008

Order No.	Type No.	pH value at 25 °C	Contents
285137841	L 168	1.68	1000 ml in DURAN® glass bottle, with manufacturer's certificate
285137677	L 1684	1.68	250 ml in DURAN® glass bottle, with manufacturer's certificate
285138098	L 401	4.01	1000 ml in DURAN® glass bottle, with manufacturer's certificate
285138008	L 4014	4.01	250 ml in DURAN® glass bottle, with manufacturer's certificate
285138102	L 687	6.87	1000 ml in DURAN® glass bottle, with manufacturer's certificate
285138016	L 6874	6.87	250 ml in DURAN® glass bottle, with manufacturer's certificate
285138119	L 918	9.18	1000 ml in DURAN® glass bottle, with manufacturer's certificate
285138024	L 9184	9.18	250 ml in DURAN® glass bottle, with manufacturer's certificate



- ▲ Reliability and measuring safety
- ▲ Extremely long storage times, thanks to hot-steam sterilization
- ▲ Without preservative agent
- ▲ A maximum of calibration safety

Advantages
FIOLAX®

Technical buffer solutions

Hot steam sterilized for longer stability, no preservation agents used.

Order No.	Type No.	pH value at 25 °C	Contents
285138213	L 4694	4.00	60 FIOLEX® ampoules à 20 ml*
285138221	L 4697	7.00	60 FIOLEX® ampoules à 20 ml*
285138205	L 4691	10.01	60 FIOLEX® ampoules à 20 ml*
285138398	L 4690	4.00/7.00	2 x 30 FIOLEX® ampoules à 20 ml*
285138192	L 4698	4.00/7.00/10.01	3 x 20 FIOLEX® ampoules à 20 ml*
285138632	L 4895/Set	4.00/7.00	2 x 9 FIOLEX® ampoules à 20 ml* with electrolyte solution L 3008

Order No.	Type No.	pH value at 25 °C	Contents
285138727	L 400	4.00	1000 ml in DURAN® glass bottle
285138032	L 4004	4.00	250 ml in DURAN® glass bottle
285138735	L 700	7.00	1000 ml in DURAN® glass bottle
285138049	L 7004	7.00	250 ml in DURAN® glass bottle
285138719	L 100	10.01	1000 ml in DURAN® glass bottle
285138057	L 1004	10.01	250 ml in DURAN® glass bottle

* 20 ml volume = approx. 17 ml content

Solutions

Color-coded technical buffer solutions in plastic bottles

Order No.	Type No.	pH value at 25 °C	Contents
285139156	LC 4004 K	4.01	250 ml in PE bottle
285139189	LC 7004 K	7.00	250 ml in PE bottle
285139218	LC 1004 K	10.01	250 ml in PE bottle



**Electrolyte solutions, aqueous
for reference electrodes and as electrolyte bridges**

Order No.	Type No.	Description	Contents
285136956	L 101	potassium chloride solution 1 mol/l	1000 ml in DURAN® glass bottle, sterilized
285138649	L 1254	potassium sulfate solution 0.6 mol/l	250 ml in DURAN® glass bottle
285138151	L 200	low temperature electrolyte (-30 °C)	1000 ml in DURAN® glass bottle
285138365	L 2004	low temperature electrolyte (-30 °C)	250 ml in DURAN® glass bottle
285138349	L 2114	2 mol/l KNO ₃ + 0.001 mol/l KCl for Ag combination electrodes	250 ml in DURAN® glass bottle
285136923	L 2214	2 mol/l KNO ₃ + 0.001 mol/l KCl for Ag combination electrodes, thickened	250 ml in DURAN® glass bottle
285138332	L 2224	potassium chloride solution 2 mol/l	250 ml in DURAN® glass bottle
285138554	L 300	potassium chloride solution 3 mol/l	1000 ml in DURAN® glass bottle, sterilized
285138427	L 3004	potassium chloride solution 3 mol/l	250 ml in DURAN® glass bottle, sterilized
285138505	L 3008	potassium chloride solution 3 mol/l	50 ml in PE bottle
285138419	L 3014	potassium chloride solution 3 mol/l, Ag/AgCl saturated	250 ml in DURAN® glass bottle
285138468	L 310	potassium chloride solution 2 mol/l, gel for sterilizable electrodes	1000 ml in DURAN® glass bottle
285138484	L 3104	potassium chloride solution 2 mol/l, gel for sterilizable electrodes	250 ml in DURAN® glass bottle
285138702	L 320 K	potassium chloride solution 2 mol/l, gel for Ag ₂ S electrodes	1000 ml in DURAN® glass bottle
285138143	L 350	potassium chloride solution 3.5 mol/l	1000 ml in DURAN® glass bottle, sterilized
285138127	L 3504	potassium chloride solution 3.5 mol/l	250 ml in DURAN® glass bottle, sterilized
285138587	L 420	potassium chloride solution 4.2 mol/l	1000 ml in DURAN® glass bottle
285138608	L 4204	potassium chloride solution 4.2 mol/l	250 ml in DURAN® glass bottle
285138590	L 911	storage electrolyte solution, sterilized	1000 ml in DURAN® glass bottle
285138560	L 9114	storage electrolyte solution, sterilized	250 ml in DURAN® glass bottle



Solutions

Electrolyte solutions, organic

for measurements in organic solutions for reference electrodes and as electrolyte bridges

Order No.	Type No.	Description	Contents
285138324	L 5014	LiCl saturated in glacial acetic acid	250 ml in DURAN® glass bottle
285138308	L 5034	LiCl saturated in ethanol	250 ml in DURAN® glass bottle

Solutions for oxygen measurements

Order No.	Type No.	Description	Contents
285138513	L 6708	electrolyte for oxygen electrodes OX 1100/OX 1100+/OX 1101	50 ml in PE bottle
285126606	OX 920	electrolyte for oxygen electrodes 9009/61	50 ml in PE bottle
285126614	OX 921	cleaning solution for oxygen electrodes 9009/61	30 ml in PE bottle
285138287	OX 060	zero point solution for oxygen electrodes OX 1100/OX 1100+	60 FIOLAX® ampoules à 20 ml

Solutions for ammonia measurements

Order No.	Type No.	Description	Contents
285137344	L 6408	electrolyte for ammonia combination electrodes	50 ml in PE bottle



Solutions and accessories for conductivity measurements

Order No.	Type No.	Description	Contents
285126503	LF 990	test solution KCl 0.001 mol/l (147 μ S/cm)	3 x 6 FIOLAX® ampoules à 20 ml
285126511	LF 991	test solution KCl 0.01 mol/l (1.41 mS/cm)	3 x 6 FIOLAX® ampoules à 20 ml
285126528	LF 992	test solution KCl 0.1 mol/l (12.9 mS/cm)	3 x 6 FIOLAX® ampoules à 20 ml
285126293	LF 995	test solutions KCl 0.01/0.1/1 mol/l (1.41/12.9/112 mS/cm)	3 x 6 FIOLAX® ampoules à 20 ml
285126166	LF 1000/Set	same as LF 999/set, in addition platinizing vessel and cable B 1 N	3 x 6 FIOLAX® ampoules à 20 ml
285136907	LF 1024	test solution KCl 0.01 mol/l (1.41 mS/cm)	250 ml in PE bottle
285126530	LF CSKC13	test solution KCl 1.3 μ S/cm	250 ml in PE bottle
285126540	LF CSKC5	test solution KCl 5.0 μ S/cm	500 ml in PE bottle

Redox test solutions

Order No.	Type No.	Redox voltage		Contents
		Pt/Calomel (KCl sat.)	Pt/Ag/AgCl (KCl 3 mol/l)	
285138373	L 4619	180 mV	220 mV	60 FIOLAX® ampoules à 20 ml, acc. to DIN 38 404-C6
285138357	L 4643	430 mV	470 mV	60 FIOLAX® ampoules à 20 ml,
285138381	L 4660	600 mV	640 mV	60 FIOLAX® ampoules à 20 ml
285138784	L 4648	180, 430, 600 mV	220, 470, 640 mV	3 x 20 FIOLAX® ampoules à 20 ml
285138184	L 430	430 mV	470 mV	1000 ml in DURAN® glass bottle
285138168	L 4304	430 mV	470 mV	250 ml in DURAN® glass bottle

Cleaning solutions for combination electrodes and reference electrodes

Order No.	Type No.	Description	Contents
285138538	L 510	pepsin/hydrochloric acid solution	1000 ml in DURAN® glass bottle
285138295	L 5104	pepsin/hydrochloric acid solution	250 ml in DURAN® glass bottle

Electrolyte bridges, other accessories

Electrolyte bridges

Shaft: glass, 12 mm Ø

Electrode vessel

Shaft: plastic



B 511



Z 451
Z 461



Z 462



Z 453

Order No.	Type No.	Length L [mm]	Dia- phragm	Remarks
285104209	B 511	103 ¹⁾	ceramic	standard taper NS 14.5 and sleeve NS 14.5 for electrode installation
285104217	B 521	120	ceramic	plastic sleeve and sleeve NS 14.5 for electrode installation
285104225	B 522	120	Pt lateral	plastic sleeve and sleeve NS 14.5 for electrode installation
285104233	B 524	120	ground joint	plastic sleeve and sleeve NS 14.5 for electrode installation

¹⁾ Length from upper end of standard taper

Order No.	Type No.	Description
285123806	BXX	plug for reference electrodes, single pole
285123703	KXX	coaxial plug for combination electrodes and indicator electrodes
285126482	NH 928	electrolyte for ammonia electrodes in 50 ml plastic bottle, 3 membrane modules
285126499	NH 995	membrane module set: 3 membrane modules, 3 caps
285126639	OX 923	3 spare membrane heads for oxygen electrodes 9009/61
285126655	OX 925	maintenance set (OX 920, OX 921, OX 923 and SF 300) for oxygen electrodes 9009/61
285126277	OX 929	5 spare membrane heads for oxygen electrodes OX 1100/OX 1100+/OX 1101
285126647	OxiCal® SL	calibrating vessel for oxygen electrodes 9009/61
285126622	SF 300	grinding foil for oxygen electrodes 9009/61
285123728	SXX	coaxial plug for extension cable and for UK socket
285215229	TZ 1520	taper adapter NS 14.5 of PTFE for electrodes with Ø 12 mm shaft
285123103	Z 341	stainless steel clamp for NS 7.5/16
285123136	Z 451	measuring and storage vessel with sleeve NS 7.5/16
285123170	Z 453	electrode vessel for storing electrodes with Ø 12 mm shaft
285123152	Z 461	measuring and storage vessel with sleeve NS 14.5/23
285123169	Z 462	flow-through measuring vessel with sleeve NS 14.5/23
285123185	Z 472	watering cap for electrodes with Ø 12 mm shaft
285122961	Z 50	Knick electrode adapter
285123193	Z 501	O-Ring seal 10.5/1.5 for electrode plug head
285123214	Z 506	plug head sealing cap with male thread for KXX and BXX plugs
285129509	Z 512	plug head sealing cap with female thread for BlueLine electrodes

Notes for the successful measuring with pH and Redox electrodes:

Decisive for a most extensive lifetime of the electrode and therefore especially solid and reproducible measurements is the maintenance and care, besides the appropriate electrode for the according application. The following points are of high importance:

Preparation and general

The electrode is delivered with a mounted watering cap. Prior to measuring and calibrating/adjusting remove this cap. It contains an aqueous solution. In most cases a 3 mol/l KCl solution (L 911) is used. For individual cases please check the technical data of the electrode. You are also welcome to contact us for further assistance.

Measuring

Open the refill opening on refillable electrodes before processing measurements. Immerse the sensor into the solution to be measured at least up to the diaphragm. When using refillable sensors please pay attention that the electrolyte fill level is 5 cm above the level of the medium to be measured (see refilling). Rinse the sensor with distilled water between measurements, however do not wipe off. Carefully dab off excess drops.

Calibration and setting

In order to minimize false measuring results with pH electrodes being used under extreme conditions or at the limit of the specified application areas, the electrodes should be calibrated more frequently and when necessary adjusted also for guaranteeing the appropriate condition for measuring safety. The measurement can only be as precise as the exactness of the last setting and its

actuality. For determining whether the pH electrode has changed since the last setting, the electrode should be calibrated i.e. checked using the according buffer solution. When deviations are discovered, a setting i.e. calibration of the electrode data (slope and zero point) on the pH meter is necessary. If buffer solutions are taken from a bottle, ensure that the bottle is reclosed immediately after removal. Never refill the used buffer solution into the bottle, always throw it away. For an accurate calibration we recommend to operate with our certified buffer ampoules, which are sterilized in superheated steam acc. to DIN 19 266. These buffer ampoules correspond to national standards. Only use fresh buffer solutions and pay attention to the temperature stability. When using electrodes without an integrated temperature sensor, the pH meter must be set to the actual temperature of the buffer solution.



Redox sensors are not calibrated. They can be checked for proper function using appropriate test solutions.

Refilling of liquid electrolyte electrodes

The missing solution in the electrolyte area of the reference system must be refilled regularly, in order to guarantee that the fill level of the electrolyte solution is at least 5 cm above the level of the medium being measured. BlueLine electrodes can be refilled simply by pumping electrolyte solution into them with a small dispensing bottle (cf. fig.).



The electrolyte solution must be completely renewed, when electrodes are used rarely and therefore the fill level varies only faintly. Any crystals in the electrolyte area can be disintegrated through warming of the electrode in a water base. The electrolyte solution should thereafter be renewed, whereby it should be rinsed a couple of times with fresh electrolyte solution prior to filling the new remaining solution into the electrolyte area.

Storage and maintenance

The electrodes should be stored between 0 and 40°C in aqueous solution (L 911). They must not be stored in distilled water. Depending on the storage conditions (temperature and humidity) the aqueous solution in the cap can run dry preterm. In this case the electrode must be rinsed for a minimum of 24 hours (L 911). The usability of the electrode must then be checked.

Cleaning

Contaminations on the pH glass membrane/Redox sensor (measuring electrode) and the diaphragm cause measuring deviations. Before heading to further cleaning steps, the glass membrane should be cleaned with an ethanol tinctured cloth. For the case that further deviations should appear with the calibration, the cleaning process should be continued by the following steps. Depending on the degree of contamination, submerge only the measuring electrode or the complete electrode i.e. including also the diaphragm into the cleaning solution:

- Coverings can be removed with diluted inorganic acid (i.e. hydrochloric acid 0,1 mol/l or NaOH 0,1 mol/l).
- Organic contamination can be dissolved with suitable solutions.
- Fat is removed with tenside solution.
- Protein removed with hydrochloric pepsin solution (cleaning solution L 510).

When cleaning please note:

- Ensure that any cleaning agent, which has leaked into the electrode, does not come into contact with the reference system; if necessary rinse out the reference electrode with electrolyte solution.
- After cleaning the electrode please rinse it with distilled water, do not rub dry.
- Any blocked ceramic diaphragms can be made functional through careful rubbing with sandpaper or using a diamond file.
- The pH glass membrane must not be scratched!
- Platinum diaphragms must not be maintained mechanically. A chemical cleaning (i.e. with diluted hydrochloric acid) will rinse and reopen (i.e. extraction with vacuum).



The platinum diaphragm developed by SCHOTT gives electrodes particularly constant and reproducible measuring characteristics.

It consists of twisted platinum wires potted into the glass shaft of the electrode. The defined intermediate space between the platinum wires ensures a continuously uniform electrolyte flow rate in all mediums and at all temperatures, which remains constant over the entire service life of the electrode.

Index

electrodes

Type No.	Order No.	Page	Type No.	Order No.	Page	Type No.	Order No.	Page
9009/61	285111664	77	BlueLine 27 pH			IL-pHT-A170MF-		
A 1180	1057997	69	1M-DIN-ID	285129950	83	BNC-CI	285114380	59
A 157 1M-BNC-ID	285130170	65	BlueLine 28 pH	285129282	81	IL-pHT-A170MF-		
A 157 1M-DIN-ID	285130160	65	BlueLine 28 pH-P	1065896	81	BNC-N	285114220	59
A 157	285129610	65	BlueLine 28-5 pH	285129570	81	IL-pHT-A170MF-DIN-N	285113910	59
A 161	285129517	63	BlueLine 29 pH	1065895	81	IL-pHT-A170MF-R-NN	285114400	59
A 161 1M-BNC-ID	285130250	63	BlueLine 29 pH-P	1065894	81	IL-pHT-H120-BNC-N	285114210	59
A 161 1M-DIN-ID	285130240	63	BlueLine 31 Rx	285129311	83	IL-pHT-H120-DIN-N	285113880	59
A 162	285129525	63	BlueLine 32 Rx	285129320	83	IL-pHT-H120MF-		
A 164	285129600	63	BlueLine 48 LF	285129488	83	BNC-N	285114200	59
A 164 1M-BNC-ID	285130290	63	BR 60	285130420	77	IL-pHT-H120MF-DIN-N	285113870	59
A 164 1M-DIN-ID	285130280	63	Ca 1100 A	285216314	77	IL-pHT-H170-BNC-N	285114270	59
A 7780	285101260	61	CA 60	285130380	77	IL-pHT-H170-DIN-N	285114250	59
A 7780 1M-BNC-ID	285130210	63	CI 60	285130350	77	IL-pHT-H170MF-		
A 7780 1M-DIN-ID	285130200	63	CN 60	285130390	77	BNC-N	285114260	59
Ag 1100	285103607	69	Cu 1100 A	285216312	77	IL-pHT-H170MF-DIN-N	285114240	59
Ag 42 A	285102051	67	CU 60	285130430	77	IL-SP-pH-A	285114320	57
Ag 6180	285102208	67	F 1100 A	285216313	77	IL-SP-pH-A-BNC	285114330	57
Ag 6280	285102343	67	F 60	285130340	77	IL-SP-pH-A-DIN	285113940	57
Ag 6580	285102216	67	H 1180	285103212	69	K 60	285130370	77
AG-S 60	285130400	77	H 161	285129590	63	KF 1100	285102030	69
AgCl 62	285102413	67	H 161 1M-BNC-ID	285130270	63	L 32	1061093	61
AgCl 6280	285102351	67	H 161 1M-DIN-ID	285130260	63	L 39	1061094	65
AgCl 65	1061051	67	H 162	285129580	63	L 39 1M-BNC-ID	285130150	65
Au 6280	285102121	67	H 61	285100207	61	L 39 1M-DIN-ID	285130140	65
B 2220+	1069994	71	H 61-500	285092583	61	L 6880	285101211	65
B 2420+	1070028	71	H 61-600	285092591	61	L 6880 1M-BNC-ID	285130110	65
B 2810+	1070029	71	H 6180	285102524	61	L 6880 1M-DIN-ID	285130100	65
B 2820+	1070044	71	H 62	285100215	61	L 7780	285101252	61
B 2910+	1070077	71	H 6280	285102532	61	L 8280	285101277	61
B 2920+	1070046	71	H 63	285100223	61	L 8880	285101285	65
B 3410+	1070048	71	H 6380	285102549	61	LF 1100+	1069976	75
B 3420+	1070070	71	H 64	285100231	61	LF 1100T+	1069977	75
B 3510+	1070100	71	H 64 1M-BNC-ID	285130230	61	LF 213 T	285106150	73
B 3520+	1070073	71	H 64 1M-DIN-ID	285130220	61	LF 213 T-ID	285106160	73
B 3610+	1070074	71	H 65	285100248	61	LF 313 T	285414360	73
B 3920+	1070075	71	H 6580	285102565	61	LF 313 T NFTC	285414351	73
BlueLine 11 pH	285129114	81	I 60	285130410	77	LF 313 T-ID	285130300	73
BlueLine 12 pH	285129122	81	IL-MICRO-pH-A	285114280	57	LF 4100+	1069978	75
BlueLine 13 pH	285129139	83	IL-MICRO-pH-A-BNC	285114290	57	LF 413 T	285106172	73
BlueLine 14 pH	285129147	81	IL-MICRO-pH-A-DIN	285113930	57	LF 413 T-3	285106148	73
BlueLine 14 pH ID	285129440	81	IL-MICRO-pHT-			LF 413 T-ID	285130310	73
BlueLine 15 pH	285129155	81	A-BNC-N	285114310	59	LF 5100+	1069979	75
BlueLine 15 pH Cinch	285095730	81	IL-MICRO-pHT-			LF 5100T+	1069990	75
BlueLine 15 pH ID	285129450	81	A-DIN-N	285114300	59	LF 513 T	285106037	73
BlueLine 16 pH	285129163	83	IL-pH-A120	285114150	57	LF 613 T	285106131	73
BlueLine 17 pH	285129171	81	IL-pH-A120-BNC	285114170	57	LF 713 T	285106189	73
BlueLine 17 pH-R	1064746	81	IL-pH-A120-DIN	285113820	57	LF 713 T-250	285106190	73
BlueLine 18 pH	285129188	81	IL-pH-A120MF	285114140	57	LF 813 T	285106250	73
BlueLine 19 pH	285129190	81	IL-pH-A120MF-BNC	285114160	57	LF 913 T	285106260	73
BlueLine 21 pH	285129217	83	IL-pH-A120MF-DIN	285113810	57	LF 913 T-ID	285130320	73
BlueLine 21 pH			IL-pH-A120MF-R	285114410	57	LFOX 1400	285104630	73
1M-BNC-ID	285129940	83	IL-pH-A170	285114190	57	LFOX 1400 ID	285130330	73
BlueLine 21 pH			IL-pH-A170-BNC	285114350	57	N 1041 A	285100486	63
1M-DIN-ID	285129930	83	IL-pH-A170-DIN	285113840	57	N 1041 A -600	285093111	63
BlueLine 22 pH	285129225	81	IL-pH-A170MF	285114180	57	N 1041 BNC	285100531	63
BlueLine 23 pH	285129233	81	IL-pH-A170MF-BNC	285114340	57	N 1042 A	285104541	63
BlueLine 23-2 pH	1063462	81	IL-pH-A170MF-DIN	285113830	57	N 1042 BNC	285105476	63
BlueLine 23-5 pH-S	1066411	81	IL-pH-A170MF-R	285114420	57	N 1043 A	285093009	63
BlueLine 24 pH	285129241	81	IL-pHT-A120-BNC-N	285113860	59	N 1048 1M-BNC-ID	285130130	65
BlueLine 24-3 pH	285129533	81	IL-pHT-A120-DIN-N	285113900	59	N 1048 1M-DIN-ID	285130120	65
BlueLine 25 pH	285129258	81	IL-pHT-A120MF-			N 1048 A	285104611	65
BlueLine 25-2 pH	1063461	81	BNC-CI	285114370	59	N 1050 A	285100375	63
BlueLine 25-5 pH	285129540	81	IL-pHT-A120MF-			N 1051 A	285100510	63
BlueLine 26 pH	285129266	81	BNC-N	285113850	59	N 1051 BNC	285100500	63
BlueLine 26 pH-Cinch	285095712	81	IL-pHT-A120MF-DIN-N	285113890	59	N 1052 A	1054512	63
BlueLine 27 pH	285129274	83	IL-pHT-A120MF-R-NN	285114390	59	N 1052 BNC	285100380	63
BlueLine 27 pH			IL-pHT-A170-BNC-N	285114230	59	N 2041 A	285100342	63
1M-BNC-ID	285129960	83	IL-pHT-A170-DIN-N	285113920	59	N 2042 A	285100359	63

electrodes

Type No.	Order No.	Page
N 42 A	285100437	61
N 42 BNC	285101544	61
N 48 A	285100445	65
N 48 BNC	285101569	65
N 50 A	285100453	61
N 52 A	285100494	61
N 52 BNC	285105451	61
N 5800 A	285105127	65
N 5800 BNC	285105579	65
N 5900 A	285105135	65
N 6000 1M-BNC-ID	285130190	65
N 6000 1M-DIN-ID	285130180	65
N 6000 A	285105151	65
N 6000 BNC	285105632	65
N 6003	285105176	65
N 61	285100001	61
N 6180	285100018	61
N 61eis	285092661	61
N 62	285100034	61
N 6250	285100112	61
N 6280	285100042	61
N 64	285100059	61
N 6480 eis	285092337	61
N 6480 eth	285092329	61
N 65	285100067	61
N 6580	285102516	61
N 6980	285101709	61
Na 61	285100026	77
NH 1100	285102808	77
NO 60	285130360	77
OX 1100+	1069975	77
Pb 1100 A	285216315	77
PB 60	285130440	77
Pt 1200	285103512	69
Pt 1400	285103537	69
Pt 1800	285103553	69
Pt 42 A	285102302	67
Pt 48 A	285102224	67
Pt 5900 A	285105192	67
Pt 5900 BNC	285105702	67
Pt 5901	285105065	67
Pt 61	285102002	67
Pt 6140	285097162	67
Pt 6180	285102232	67
Pt 62	285102019	67
Pt 6280	285102249	67
Pt 6580	285102257	67
Pt 6880	285100075	67
Pt 6980	285102265	67
Pt 8280	285102281	67
W 2030+	1069991	79
W 2130+	1069992	79
W 2180-KOAX	285119030	79
W 5780 NN	285105221	79
W 5790 NN	285105254	79
W 5790 PP	285105776	79
W 5791 NN	285105262	79
W 5980 NN	285105287	79

accessories

Type No.	Order No.	Page	Type No.	Order No.	Page
9907/21	285124716	85	L 4798	285138238	86
9909/31	285125618	85	L 4799	285138262	86
9910/11	285125515	85	L 4893/Set	285138279	86
9910/21	285125215	85	L 4895/Set	285138632	87
9919/21	285125523	85	L 5014	285138324	90
9919/41	285125548	85	L 5034	285138308	90
A 1 A	285122904	85	L 510	285138538	91
A 1 BNC	285123793	85	L 5104	285138295	91
B 1 N	285121916	85	L 6408	285137344	90
B 1 P	285122012	85	L 6708	285138513	90
B 1X	285121813	85	L 687	285138102	86
B 511	285104209	93	L 6874	285138016	86
B 521	285104217	93	L 700	285138735	87
B 522	285104225	93	L 7004	285138049	87
B 524	285104233	93	L 911	285138590	89
BXX	285123806	93	L 9114	285138560	89
KXX	285123703	93	L 918	285138119	86
L 1 A	285122456	85	L 9184	285138024	86
L 1 BNC	285122497	85	LB 1 A	285122653	85
L 1 EE	285122501	85	LB 1 BNC	285122661	85
L 1 N	285122457	85	LB 3 A	285122678	85
L 1 NN	285122489	85	LC 1004 K	285139218	88
L 1 R	285122534	85	LC 4004 K	285139156	88
L 1 X	285122407	85	LC 7004 K	285139189	88
L 100	285138719	87	LF 1000/Set	285126166	91
L 1004	285138057	87	LF 1024	285136907	91
L 101	285136956	89	LF 990	285126503	91
L 1254	285138649	89	LF 991	285126511	91
L 168	285137841	86	LF 992	285126528	91
L 1684	285137677	86	LF 995	285126293	91
L 2 A	285122464	85	LF CSKC13	285126530	91
L 2 NN	285122448	85	LF CSKC5	285126540	91
L 200	285138151	89	LS 1 ANN	285122707	85
L 2004	285138365	89	LS 1 BNCNN	285122723	85
L 2114	285138349	89	LS 1 D8	1066726	85
L 2214	285136923	89	LS 1 N6	1066728	85
L 2224	285138332	89	LS 1 RNN	285122756	85
L 300	285138554	89	LS 1 ST4LF	1069104	85
L 3004	285138427	89	LS 1 ST4OX	1066727	85
L 3008	285138505	89	LS 3 ANN	285122715	85
L 3014	285138419	89	LS 3 BNCNN	285122731	85
L 310	285138468	89	NH 928	285126482	93
L 3104	285138484	89	NH 995	285126499	93
L 320 K	285138702	89	OX 060	285138287	90
L 350	285138143	89	OX 920	285126606	90
L 3504	285138127	89	OX 921	285126614	90
L 400	285138727	87	OX 923	285126639	93
L 4004	285138032	87	OX 925	285126655	93
L 401	285138098	86	OX 929	285126277	93
L 4014	285138008	86	OxiCal® SL	285126647	93
L 420	285138587	89	SF 300	285126622	93
L 4204	285138608	89	SXX	285123728	93
L 430	285138184	91	TZ 1520	285215229	93
L 4304	285138168	91	Z 341	285123103	93
L 4619	285138373	91	Z 451	285123136	93
L 4643	285138357	91	Z 453	285123170	93
L 4648	285138784	91	Z 461	285123152	93
L 4660	285138381	91	Z 462	285123169	93
L 4690	285138398	87	Z 472	285123185	93
L 4691	285138205	87	Z 50	285122961	93
L 4694	285138213	87	Z 501	285123193	93
L 4697	285138221	87	Z 506	285123214	93
L 4698	285138192	87	Z 512	285129509	93
L 4790	285138402	86			
L 4791	285137977	86			
L 4794	285138246	86			
L 4796	285138254	86			
L 4797	285137985	86			

ProcessLine – Process electrodes

for measuring pH, temperature and redox potentials

One for all applications – for highest demands

ProcessLine electrodes are low maintenance sensors for heavy duty process applications, as they are especially present in the chemical industry.

They are ideally suitable for measuring media with extreme ionic strength – whether boiler feed water or brine – also in strongly oxidizing acid and alkali containing media.

The ProcessLine electrodes' special design with regard to accuracy, stability, rapidness and durability is very close to the one of liquid electrolyte electrodes, although the ProcessLine does not require refilling the electrolytes and its complex pressure sequence regulation. Therefore the ProcessLine electrodes require only low-maintenance, including calibration and adjusting efforts, hence offering a high potential for cost savings.

Duralid solid electrolyte with high content of KCl and special formulation

The solid reference electrolyte Duralid does not require a special diaphragm – the reference system holds a direct contact to the measuring media via the two open connections. This minimizes the risk of contamination/blockage of the diaphragm – the main source for measuring failures and even outfall – and guarantees long durability and high accuracy.

The long lifetime and small liquid junction potential resulting in high accuracy measurements of the ProcessLine electrodes is based on the special formula and fabrication of the Duralid electrolyte:

- High content of the conductivity salt potassium chloride in polymer and therefore high electrolyte output into the measuring media, reduces the interferences of the measurement through diffusion potentials between the junction of reference electrode and measuring media.
- The special distribution of the potassium chloride in the Duralid polymer counteracts positively against a reduced durability of the reference system, by releasing a high quantity of electrolyte.

This special attribute given by the Duralid, does not only improve the durability and the response characteristic, it also enables stable measuring values – even under most difficult conditions such as changing flow rate/rotational frequency of stirrer or with measurements in organic solvents.

Cushion for pressure compensation in the reference electrode


Pressure and temperature fluctuations can easily be managed by the ProcessLine electrodes due to the integrated pressure compensation cushion inside the reference electrode.

Reliable H membrane glass

Besides the reference electrode, also the measuring electrode is of great importance regarding rapidness and accuracy of the measurement. The glass electrode of the ProcessLine series features a H membrane glass, a very high-quality and approved special glass. It excels by its high-temperature application range and very low alkaline errors.

The special ball shape enables an optimal membrane resistance of 300 MOhm and ensures an easy cleaning.





- ▲ **Low maintenance**, i. e. no refilling of electrolyte or installation of complicated pressure sequence regulations.
- ▲ **Hole diaphragms**, therefore no contamination or blockage of the reference electrode.
- ▲ **Duralid electrolyte with high proportion of KCl and special consistency**
Long durability as well as fast and stable measuring values. Furthermore, no substances of animal origin.
- ▲ **Buffer in the reference electrode** for compensation of pressure and temperature fluctuation.
- ▲ **Approved H membrane glass** with very low alkaline error and optimized ball shape.
- ▲ **Wide application range for media** with extreme ionic strength, strong oxidation character, high alkaline or acid components and also organic solvents.
- ▲ **Certificate for temperature and pressure resistance** of 10 bar at 0 to 130 °C.
- ▲ **Shaft length** from 120, 225, 325, 360 and 425 mm suitable for all assembling conditions.
- ▲ **Versions with Pt 100 and Pt 1000 temperature sensor** with SMEK- and also VP plug head for high flexibility.

Advantages
ProcessLine

Process Electrodes

The ProcessLine electrodes take up a small part of our extensive process electrode programme:

Further information is available in our specialized catalogue "Process Electrodes" displayed and downloadable on our website; it can also be delivered personally on request.

Contents titration

Selection table titration	Page 101
TITRONIC® <i>basic</i>	Page 103
TITRONIC® <i>universal</i>	Page 104
TitroLine <i>easy</i>	Page 106
KF Titration: Selection guide for Coulometry and Volumetry	Page 108
TitroLine <i>KF</i> and Titroline <i>KF trace</i>	Page 110
Order information TITRONIC®, TitroLine	Page 117
TitroLine <i>alpha plus</i>	Page 118
Selection table for titration electrodes	Page 120
Water analysis according to Karl Fischer with the TitroLine <i>alpha plus</i>	Page 124
TITRONIC® 110 – the piston burette with the <i>plus</i>	Page 126
Technical data TitroLine <i>alpha plus</i> and TITRONIC® 110 <i>plus</i>	Page 127
TW <i>alpha plus</i> sample changer	Page 128
TitriSoft 2.6 titration software	Page 130
TitriSoft 2.6 P titration software	Page 134
Order information TitroLine <i>alpha plus</i>	Page 138
Order information TW <i>alpha plus</i>	Page 139

Selection table titration

Overview table piston burettes

Application	TITRONIC® <i>basic</i>	TITRONIC® <i>universal</i>	TITRONIC® 110 <i>plus</i>
Manual titration	■	■	■
Automatic titration ⁽¹⁾		■	■
Dosing of pre-selected volumes		■	■
Pre-titration		■	■
Variable dosing and filling speed		■	■
20 ml dosing unit	■	■	■
50 ml dosing unit		■	■
1, 5, 10, 20 and 50 ml changer units			■
Results output via RS 232	■	■	■
Remote control via RS 232		■	■

¹⁾ if piston burette is connected to a TitroLine alpha *plus* or TitriSoft

Overview table of titrators

Application	TitroLine <i>easy</i>	TitroLine <i>KF</i>	TitroLine <i>KF trace</i>	TitroLine alpha <i>plus</i>
pH/mV titration aqueous (acidity, hydrochloric acid, citric acid, "Kjeldahl", ammonia ...)	■			■
pH/mV titration non-aqueous (TAN/TBN, FFA, titrations with perchloric acid ...)				■
Redox titrations (iodometry, permanganometry ...) ⁽²⁾	■			■
Halogenide titrations (chloride, "salt", bromide ...)	■			■
H ₂ S and mercaptan				■
pH-stat applications (enzyme kinetics, soil samples, biotechnology)				■
Water analysis according to KF Volumetric method (10 ppm – 100 %)				■
Water analysis according to KF Volumetric method (0.01 % – 100 %)		■		■
Water analysis according to KF Coulometric method (1 ppm – 5 %)			■	
Bromine number				■
Titration with <u>more</u> than one end point or equivalence point (phosphoric acid ...)				■
Applications with several piston burettes				■
Applications with sample changer				■
Applications with TitriSoft				■

⁽²⁾ except COD and sulphuric acid (SO₂)

Dosing, titrating and water analysis according to Karl Fischer can be so easy

Innovative electrochemistry – right from the start

With the development of the glass electrode more than 70 years ago, we created the basis for the success of electrochemical measurement. Since then, with our range of efficient pH glasses, innovative electrodes and instruments such as pH meters, conductometers, oxygen measuring instruments, we have turned the electrochemical measurement into an indispensable, trouble-free and reliable procedure that is now being used throughout the world.

Building on this know-how, we have also developed a range of reliable laboratory instruments for dosing, titrating and Karl Fischer water analysis. The coulometric KF titrator TitroLine KF *trace* is the most recent instrument of this series.

KF titrators from SCHOTT Instruments combine ease of use with maximum accuracy, and the robustness required for daily operation in the laboratory. Benefits, which far outweigh the cost of these instruments.

For complex application such as difficult, nonaqueous titrations and for automatic measuring stations, the TitroLine alpha *plus* titration system is also available.

Just what you need to make your routine daily work simpler and even better

Like the TitroLine *easy* and TitroLine *KF* titrators, the TITRONIC® *basic* and TITRONIC® *universal* piston burettes are robust tools for dosing and titrating, which were designed specifically for daily routine use in the laboratory. Despite their robustness, these are high precision instruments. Even the simplest burette is equipped with an UV-protected precision glass cylinder made of DURAN® and a motor-controlled 3/2-way valve made of extremely resistant PTFE/ETFE. But we have also focussed our attention on the importance of easy, trouble-free operation – so a manual is something you may never need to use.



TITRONIC® and DURAN® are registered trademarks.
Subject to technical modifications.

TITRONIC® *basic*

The burette with the ›Mouse‹

Anyone with a TITRONIC® *basic* in the laboratory will leave bottle-top burettes and classic glass burettes on the shelf. Manual titration can be performed more reliably and accurately with the TITRONIC® *basic*, and the results can be documented when required.

Operation is so easy

The titration process is carried out by pressing a button on the ›mouse‹ – the handheld device TZ 3680. You can monitor the dosed quantity conveniently on the large display. The TITRONIC® *basic* is equipped with an RS232-C serial interface so that you can document your results. Here, for example, you can also connect our small, practical TZ 3460 rollpaper printer or any other printer with a serial RS-232-C interface. Needless to say, you can also connect your PC to the TITRONIC® *basic*.

Precision is integrated

The accuracy of the TITRONIC® *basic* is guaranteed by the precision glass cylinder made of DURAN® borosilicate glass with its measurement deviation of less than 0.1 %. And the motor-driven, chemically resistant 3/2-way valve also provides its contribution for precise, reproducible values: It enables unpresurised drawing and dosing and therefore effectively prevents outgassing of liquids and vapour formation due to excessive vacuum pressure.

The magnetic stirrer is available as an accessory

The TM 96 magnetic stirrer is available as an accessory. It is connected directly to the burette, which provides the necessary power.



The complete workplace: precision analysis with no shaky compromises. With 8000-step resolution, precision glass cylinder with UV-protection, motor-controlled 3/2-way valve of extremely resistant PTFE/ETFE and an interface for documentation of the results. It's better to be on the safe side! (The bottle set must be ordered separately as an accessory.)

Technical data

Hand control element	miniature 4 pole round socket, conforming to DIN standards
RS-232-C	for connecting a printer with serial interface or PC for documentation
Display	four digit LCD display, 20 x 48 mm, height of digits: 12.7 mm
Volume display	0.01 ... 999.9 ml
Resolution	0.01 ml
Cylinder	20 ml DURAN® borosilicate glass cylinder with UV protection sleeve
Tubing	FEP with UV protection systematic error 0.1 %,
Dosing accuracy	random error 0.05 %, determined according to EN ISO 8655-6
Valve	3/2-port directional control valve made of PTFE/ETFE
Housing material	polypropylene and polyflam RPP 371 NT, 20 % talcum
Front foil	polyester 135 x 310 x 205 mm (W x H x D),
Dimensions	including dosing unit, without stirrer
Weight	approx. 2.1 kg
Ambient temperature	+10 ... +40 °C (for operation and storage)
Power supply	230 V~; 50/60 Hz or 115 V~; 50/60 Hz
Appliance safety	corresponds to Protection Class II in accordance with DIN EN 61010, part 1
Conformity	EN ISO 8655-3

TITRONIC® *universal*

Titration manually, dosing perfectly

The TITRONIC® *universal* is a perfect motor-driven burette for manual titration and an extremely precise dosing instrument for dosable liquids, solvents and titrating agents. However, the TITRONIC® *universal* not only first-rate as a stand-alone instrument – it also thrives as the heart of a computer-controlled dosing or titrating system.

Adjusting easily, dosing precisely

With the TITRONIC® *universal* you can preselect any dosing volume from 0.01 ml to 999.99 ml easily with the keypad and you can adjust the dosing speed to a continuously controllable setting. Furthermore, with the TITRONIC® *universal* you can define the waiting time between the volume steps, a useful tool for incremental dosing tasks. The dosing process is carried out precisely upon being called. This, by the way, is also extremely practical in the case of manual titration with the hand-held device: Using a precisely adjusting pre-titrating volume, which can be called up at the press of a button before each titration, you can reduce your titration times considerably.

Documenting results reliably

To document your results, simply connect our small, practical TZ 3460 rollpaper printer or any other printer with a serial RS-232-C interface.

The TITRONIC® *universal* gets on quite well with the PC

We have equipped the TITRONIC® *universal* with **two** serial RS-232-C interfaces. This allows you to not only connect a printer in order to document data in the stand-alone mode but also extends the available range of use of the TITRONIC® *universal* quite considerably. For instance, you can use a PC to control all functions of the TITRONIC® *universal* via one of the two serial interfaces. The

address is set automatically or manually. But the TITRONIC® *universal* can do a lot more: For complex dosing and titrating processes, up to 16 burettes can be connected in series whenever required. The burettes are connected to one another via the RS-232-C interfaces according to the daisy chain principle. Accordingly, each instrument can then be addressed separately and provides independent feedback data – without an additional data cable.



Technical data

Designed for maximum precision

All components of the TITRONIC® *universal* are designed for maximum precision. This begins with the dosing attachments, which are available in 20 ml or 50 ml volumes. The glass cylinders made of DURAN® borosilicate glass are precisely calibrated and provided with an UV protective coating. The dosing piston is driven by a step motor with a resolution of 8,000 steps. The motor-controlled 3/2-way valve is made of extremely resistant PTFE/ETFE. This 3/2-way valve enables unpressurised drawing and dosing so that outgassing of liquids is prevented as well as vapour formation due to excessive vacuum pressure.

Made for robust laboratory operation

All parts of the TITRONIC® *universal* that come into contact with liquids are made of chemically resistant materials. A polyester front foil protects the keypad and display, and the tubing is in FEP with UV protection.

The magnetic stirrer is available as an accessory

The TM 96 magnetic stirrer is available as an accessory. It is connected directly to the burette, which provides the necessary power.

Keyboard connection	miniature 4 pole round socket, conforming to DIN standards
Stirrer connection	plug-and-socket connection with integrated low-voltage power supply (15 V DC) for the TM 96 magnetic stirrer
RS-232-C-1	for connecting a printer with serial interface or a PC to document consumption in ml or for data backup
RS-232-C-2	connection of additional piston burettes TITRONIC® <i>universal</i> ('Daisy Chain')
Configuration of the RS-232-C interface	connection: miniature 4 pole round socket preset: 1 stop bit adjustable: baud rate: 1200, 2400, 4800 or 9600 baud word length: 7 or 8; parity: no, even or odd
Display	8-line LCD display, 69 x 39 mm, 64 x 128 pixel, with background illumination and contrast adjustment
Volume display	00.00 ... 999.9 ml
Display resolution	0.01 ml
Dosing volume	0.0 ... 999.99 ml
Dosing speed	0.1 ... 40 ml/min (with 20 ml dosing unit) 0.1 ... 100 ml/min (with 50 ml dosing unit)
Filling time	30 s to 999 s adjustable (100% in relation to the cylinder volume)
Pre-titrating volume	0.1 ml to 99.99 ml
Increment volume	0.01 ... 999.99 ml
Waiting time between the increments	0.1 ... 999.9 s
Cylinder	20 ml or 50 ml DURAN® borosilicate glass cylinder with UV protection sleeve
Dosing accuracy	systematic error 0.1 %, random error 0.05 %, determined according to EN ISO 8655-6
Valve	3/2-port directional control valve made of PTFE / ETFE
Tubing	FEP with UV protection
Housing material	polypropylen and polyflamm RPP371 NT, 20 % talcum
Front foil	polyester
Dimensions	135 x 310 x 205 mm (W x H x D), including dosing unit, without stirrer
Weight	approx. 2.1 kg
Ambient temperature	+10 ... +40 °C (for operation and storage)
Power supply	230 V~; 50/60 Hz or 115 V~; 50/60 Hz
Power consumption	18 VA
Conformity	EN ISO 8655-3


With the TITRONIC® *universal* in the stand-alone mode, you can use the keypad to input all settings conveniently on the instrument. The practical hand-held device can be used for manual titrations or to start and stop a dosing task. (The bottle set must be ordered separately as an accessory.)

TitroLine *easy*

The intelligent titrator for your routine daily work

Quick and easy as its name suggests

The TitroLine *easy* is the ideal titrator for your routine daily work. This instrument provides you with the perfect combination of a piston burette, a pH/mV meter and integrated intelligence. Ten titration methods for various applications are preinstalled and can be called up easily as required. The methods are pre-parameterised. You only need to select your titration procedure: with a self-searching end point, with a pre-set end point, or manual titration with the »mouse«. The titration process begins as soon as you press the start button. This saves you time and money.



With the TitroLine Chart software (option), the titration curve can be displayed on the monitor of a connected PC and the titration data can be processed.



Practical and compact: A complete measuring unit. The magnetic stirrer is included. It is connected to the TitroLine *easy*. The bottle set must be ordered separately as an accessory.

Suitable applications for the TitroLine *easy* include:

- salt content in foods
(cheese, soya sauce, ketchup)
- total acid in wine and beverages
- nitrogen according to Kjeldahl



Technical data

easv

Measuring amplifier	measuring input pH/mV electrode: pH-input with 12-bit converter for highly accurate resolution of the measuring signal during titration measuring range pH: 0.00 ... 14.00 measuring range mV: -1400 ... +1400 electrode socket according to DIN 19262 or BNC-socket and reference electrode 1 x 4 mm measuring input temperature sensor Pt 1000, measuring range: -30 ... +115 °C connection socket 2 x 4 mm and 1 x 2 mm
Keyboard connection	miniature 4 pole round socket, conforming to DIN standards
Stirrer connection	plug-and-socket connection with integrated low-voltage power supply (15 V DC) for the magnetic stirrer TM 96
RS-232-C interface	for connecting a printer with serial interface or PC for documentation
Configuration of the RS-232-C interface	preset: 4800 baud, 7-bit word length, 2 stop bits, no parity
Display	matrix-LCD display 69 x 39 mm, 64 x 128 Pixel background illumination and contrast adjustment
Volume display	00.00 ... 999.9 ml
Display resolution	0.01 ml
Cylinder	20 ml DURAN® borosilicate glass cylinder with UV protection
Dosing accuracy	systematic error 0.1 %, random error 0.05 % determined according to EN ISO 8655-6
Calibration	two-point calibration, selection of eight stored buffer solutions in conformity with DIN 19 266 and NBS
Valve	3/2-port directional control valve made of PTFE / ETFE
Tubing	FEP with UV protection
Housing material	polypropylene and polyflamm RPP 371 NT, 20 % talcum
Front foil	polyester
Dimensions	135 x 310 x 205 mm (W x H x D), including dosing unit, without stirrer
Weight	approx. 2.4 kg
Ambient temperature	+10 ... +40 °C (for operation and storage)
Power supply	230 V~; 50/60 Hz or 115 V~; 50/60 Hz
Power consumption	24 VA
Appliance safety	corresponds to Protection Class II in accordance with DIN EN 61 010, Part 1
Conformity	EN ISO 8655, part 3

The sensors – from SCHOTT Instruments

Suitable sensors include pH combination electrodes with or without integrated temperature sensors (Pt 1000), redox combination electrodes, Ag combination electrodes or separate measuring or reference electrodes.

Stored data: the buffer solutions

Data for 2.00 / 4.00 / 4.01 / 6.87 / 7.00 / 9.18 / 10.01 / 12.45 buffers, including temperature coefficients are already stored in the TitroLine *easy*.

Maximum precision for reproducible results

All components of the TitroLine *easy* are designed for maximum accuracy. The glass cylinders made of DURAN® borosilicate glass are precisely calibrated and provided with an UV protective coating. The motor-controlled 3/2-way valve is made of extremely resistant PTFE/ETFE. This 3/2-way valve enables unpressurised drawing and dosing so that outgassing of liquids is prevented as well as vapour formation due to excessive vacuum pressure.

As robust as required for laboratory operation

All parts of the TitroLine *easy* that come into contact with liquids are made of chemically resistant materials. A polyester front foil protects the keypad and display, and the tubing is in FEP with UV protection.

Karl-Fischer Titration – the method for determining water

Quite some experienced analyst will be unpleasantly reminded by the pyridine smell, when hearing the name Karl Fischer. However, modern reagents and most user friendly analysing instruments have definitely cleared that picture. Nowadays all applications can be handled and processed very easily, fast and accurate by using the **coulometric** and **volumetric** Karl Fischer titration instruments. Thanks to its selectivity and precision, the Karl Fischer titration has achieved to be established as the most important method for determining water and humidity.


In the following, we would like to ease your decision for either the coulometric (TitroLine KF trace) or the volumetric (TitroLine KF) KF titrator a little.

The basic principle of the water determination according to Karl Fischer (short: **KF**) is a reaction of iodine with water in an alcoholic solution with presence of sulfuric acid and a base.

With the **volumetric** method the iodine can be accurately added through a piston burette or **coulometric** directly produced in the reaction vessel. The

difference between the volumetry and coulometry mainly exists in the manner of dosing the iodine for the titration.

The illustration shows the different types of dosing:




TitroLine KF




TitroLine KF trace

In practice small differences occur between the two methods which are displayed in the table. The advantages of the volumetry lie in the different types of sample adding and solvent variations, offering more flexible operation potentials. Where on the other hand the coulometry can handle yet lower detection limits and score with the even simpler handling. The compared work flow with coulometry and volumetry are shown with the following illustration. The clearly shorter and easier sequence is noticable with the coulometry.

Coulometric KF titration



Volumetric KF titration



Comparison: Coulometric and volumetric Karl-Fischer-titration		
Property	Coulometry	Volumetry
Water amount and sample amount	Small water amount Small sample amounts	Medium and large water amounts Adapted sample amount
Sample types	Liquid Gaseous (i.e. KF oven) Solid samples with oven	Solid Liquid
Sample addition and preparation	Direct with syringe Gas inlet with oven External extraction Solid samples are evaporated with an oven	Solid samples are added directly Sample preparation with homogenisator Working at higher temperature Direct with syringe
Working method	Very fast Very simple	Fast Simple
Working range	µg range 10 µg up to 5 mg water	mg range 200 µg up to 50 mg water
Trueness	Pretty good for small water amounts > 400 µg Wasser (± 0,5%)	Pretty good for water amounts > 5 mg water (± 0,5%, standardization required!)
Reproducibility	Typical RSD of appr. 1% for water > 400 µg	Typical RSD of appr. 1% for water > 5 mg

TitroLine KF *trace*

In dialogue the coulometric Karl Fischer Titration is quite easy!

Karl Fischer titrations made easy

With the new TitroLine KF *trace* the coulometric water determination according to Karl Fischer you cannot go wrong:

The large display shows every work step ahead in a dialogue structure. The pre-parameterized methods are easily recalled and enhance the total work process. Also the versatility makes the both KF titrators a trouble-free KF measuring place for nearly all areas in the industry, such as for pharmaceutical, chemical and petroleum industry.

The coulometric Karl Fischer titrator TitroLine KF *trace* is the dedicated instrument for determining even smallest water content in your samples. As the coulometric determination of water does not require a standardisation of a titrant, the handling is easier compared with the volumetric titration: Once the instrument is installed, the reagents are inserted into the titration cell and the instrument is switched-on. The TitroLine KF *trace* starts to operate immediately. The conditioning is triggered in the background and automatically determines the drift. Only few minutes later the TitroLine KF *trace* is ready for the first samples.




Conveniently with methods

The TitroLine KF is already programmed with the following methods: sample titration, titre water, titre liquid standard, titre tartrate dihydrate, blank value open and blank value solvent. The methods for titre determination cease to apply for the TitroLine KF *trace*. All methods are already set with the commonly used parameters. However, should it become necessary, the parameters can of course be changed.

Parameterization – just in case it will be necessary

The large display gives a clear overview of the next steps in process. Parameterization using the arrow key and the enter/F1 and ESC/F4 key is fairly easy. Taking a look into the operating manual is almost not necessary.



Live titration process

The TitroLine KF *trace* allows watching over the titration by displaying a real-time measuring curve. With just a keystroke you have the option to switch between the standard or graphic display.

Documentation – exactly the way you require it



You can decide if you want to print the results in either a short form, as standard form with a curve (only TitroLine KF *trace*) or as complete GLP printout including all method parameters. Of course all results are also indicate with the mean value and the drift on the display.

Automatic selection of the correct calculation formula

Two different may be used to calculate the results. When choosing the method, the correct formula is automatically selected and pre-assigned with the corresponding values. Measurement units for the result can be selected from %, ppm, mg, mg/l, mg/pc (pc= piece) μg (TitroLine KF *trace*) or ml (TitroLine KF). The blank value is calculated in ml or μg and automatically subtracted from the sample titration result.

Statistics

For the statistical evaluation of the analysis the mean value, standard deviation or the relative standard deviation can be determined. The mean value of the titre at the TitroLine KF is the automatic reference for the calculation of the sample result.



Titration stand and titration vessel: Accessories made to match

Titrated samples can be extracted through pushing a button on the titration stand TM KF (standard with TitroLine KF and KF *trace* Module 2 + 4). A further keystroke provides a new reagent. An integrated magnetic stirrer in TM KF takes care of the balanced distribution of reagents and sample.



The titration vessels are hermetically sealed and avoid the penetration of moisture extensively (low drift!). The removable glass vessel from the TitroLine KF is available in two sizes and it is easy to clean. For the TitroLine KF *trace* are to be used two different glass vessels with 3 and 5 openings. Both have a very low drift.

Connection of analytical balances, printers, PC KF oven ...


The two RS-232C interfaces and one USB-port* allow you to connect a balance and for automatic transfer of the weighing data and a printer at the same time. It is also possible to connect a PC via the additional USB-port*.

* only TitroLine KF trace

PC control

Both titrators can be easily connected at a PC. The software "KF-Soft"* allows an easy transfer of the data and titration curve on a PC and store it in the integrated database. *

* only TitroLine KF



Technical Specifications TitroLine KF / TitroLine KF *trace*

Hardware

	TitroLine KF	TitroLine KF <i>trace</i>
Display	High contrast 8-lines LCD with 64 x 128 pixel and background illumination ; contrast adjustable	High contrast 8-lines LCD with 64 x 128 pixel and background illumination ; contrast adjustable
Interfaces	2 x RS232 for PC or printer, balance and further devices („daisy chain“)	2 x RS232 for PC or printer, balance and further devices („daisy chain“) 1 x USB („slave“) for PC
Indicator electrode	Dual platinum electrode Connection 2 x mm socket	Dual platinum electrode Connection 2 x mm socket
Generator electrode		Generator electrode Connection 2 x mm socket
Keyboard	5-pole DIN-socket for PC-keyboard	For PC-keyboard with PS/2-plug such as TZ 2835
Stirrer/pump	Stirrer TM 135 respectively stirrer and pump of the titration stand TM KF	Stirrer TM 135 respectively stirrer and pump of the titration stand TM KF
Cylinder	20 ml made out of DURAN®	
Valve	Motor driven 3/2 way valve made out of PTFE/ETFE	
Dimensions	310 x 265 x 205 mm (h x w x d), height with titration stand and titration vessel	200 x 265 x 205 mm (h x w x d) with titration stand TM 135/TM KF 310 x 265 x 205 mm (h x w x d), height with titration vessel
Weight	2,1 kg for basic unit; appr. 3,2 kg for complete unit with titration stand TMKF	Appr. 1,4 kg for basic unit; appr. 2,5 kg for complete unit with titration stand TMKF (module 2 and 4)
Casing	Polypropylene	Polypropylene
Front foil	Polyester	Polyester
Temperature	Ambient temperature: + 10 ... + 40 C for operation and storage	Ambient temperature: + 10 ... + 40 C for operation and storage
Power supply	Mains 230 V, 50/60 Hz or 115 V ; 50/ 60 Hz power drain: 30 VA	Universal power adapter 100-140 V; 50/60 Hz, power drain: 30 VA

Software

	TitroLine KF	TitroLine KF <i>trace</i>
Measuring range	100 ppm – 100 %	10 µg – 100 mg / 1 ppm – 5 % (recommended)
Number of methods	8 (3 x sample, 3 x titre, 2 x blank value)	10 (9 x sample, 1 x blank value)
Conditioning	On start, automatic drift correction	Automatic after switch on, drift correction
End criteria	End point delay, drift	Drift, drift stop tolerance
Autostart after sample addition	Only after conformation of sample weight	■
Statistic	Mean value, standard deviation and rel. relative standard deviation	Mean value, standard deviation and rel. relative standard deviation
Recalculation	After new weighted sample or sample volume has been entered Erasing of one result from a series of measurement	After new weighted sample or sample volume has been entered Erasing of one result from a series of measurement
Online curve		■
Documentation	GLP	GLP + curve print out
Result output	%, ppm, mg, mg/l, mg/pc (pc= piece), ml	%, ppm, mg, mg/l, mg/pc (pc= piece, µg
Password		■
Update software	EPROM-change	Update via RS232 and USB



Ordering information TitroLine KF and TitroLine KF *trace*

TitroLine KF and TitroLine KF <i>trace</i>		Order no.
TitroLine KF <i>trace</i> M1 Complete module for coulometric Karl-Fischer titration	Scope of delivery: TitroLine KF <i>trace</i> basic unit, magnetic stirrer TM 135, generating electrode TZ 1752 without diaphragm, titration vessel TZ 1751, micro dual platinum-electrode KF 1150, connection cable generating electrode	285212258
TitroLine KF <i>trace</i> M2 Complete module for coulometric Karl-Fischer titration	Scope of delivery: TitroLine KF <i>trace</i> basic unit, titration stand with pump TM KF, generating electrode TZ 1752 without diaphragm, titration vessel TZ 1754, micro dual platinum-electrode KF 1150, connection cable generating electrode	285212268
TitroLine KF <i>trace</i> M3 Complete module for coulometric Karl-Fischer titration	Scope of delivery: TitroLine KF <i>trace</i> basic unit, magnetic stirrer TM 135, generating electrode TZ 1753 with diaphragm, titration vessel TZ 1751, micro dual platinum-electrode KF 1150, connection cable generating electrode	285212278
TitroLine KF <i>trace</i> M4 Complete module for coulometric Karl-Fischer titration	Scope of delivery: TitroLine KF <i>trace</i> basic unit, titration stand with pump TM KF, generating electrode TZ 1753 with diaphragm, titration vessel TZ 1754, micro dual platinum-electrode KF 1150, connection cable generating electrode	285212288
TitroLine KF-230 V Volumetric KF-Titrator	Scope of delivery: titration unit, titration stand with integrated stirrer and pump TM KF, titration vessel TZ 1770, micro dual platinum-electrode KF 1100 and starter kit	285212248
TitroLine KF-115 V Volumetric KF-Titrator	Scope of delivery: titration unit, titration stand with integrated stirrer and pump TM KF, titration vessel TZ 1770, micro dual platinum-electrode KF 1100 and starter kit	285212231

Accessories for TitroLine KF and TitroLine KF *trace*

TZ 2835	PC-keyboard (with PS2/DIN-adapter for TitroLine KF)	1007852
TZ 1052	Evaporation oven for water determination according to Karl-Fischer	285214721
TZ 1060	Accessory for Evaporation oven TZ 1052	285218115
TZ 2073	KF-Soft for TitroLine KF	285221733
TZ 3460	RS-232-C printer for TitroLine KF, incl. connection cable TZ 3090, 230 V	285225608
TZ 3461	RS-232-C printer for TitroLine KF <i>trace</i> , incl. connection cable TZ 3090, 230 V	285225610
TZ 3465	RS-232-C printer for TitroLine KF, incl. connection cable TZ 3090, 115 V	285225657
TZ 3466	RS-232-C printer for TitroLine KF <i>trace</i> , incl. connection cable TZ 3090, 115 V	285225660

Ordering information TITRONIC[®], TitroLine

TITRONIC[®] basic and TITRONIC[®] universal	Order no.
TITRONIC [®] basic module 1, (230 V)	285212572
TITRONIC [®] basic module 2, same as module 1, with magnetic stirrer TM 96, (230 V)	285212823
TITRONIC [®] universal 20 ml module 1, (230 V)	285212429
TITRONIC [®] universal 20 ml module 2, same as module 1, with magnetic stirrer TM 96, (230 V)	285212437
TITRONIC [®] universal 50 ml module 1, (230 V)	285212445
TITRONIC [®] universal 50 ml module 2, same as module 1, with magnetic stirrer TM 96, (230 V)	285212494
TITRONIC [®] basic module 1, (115 V)	285212564
TITRONIC [®] basic module 2, same as module 1, with magnetic stirrer TM 96, (115 V)	285212815
TITRONIC [®] universal 20 ml module 1, (115 V)	285211921
TITRONIC [®] universal 20 ml module 2, same as module 1, with magnetic stirrer TM 96, (115 V)	285211962
TITRONIC [®] universal 50 ml module 1, (115 V)	285211979
TITRONIC [®] universal 50 ml module 2, same as module 1, with magnetic stirrer TM 96, (115 V)	285211987
TitroLine easy	
TitroLine easy module 1 without electrode, (230 V)	285212597
TitroLine easy module 2 for pH titration, same as module 1, with one pH-electrode and buffer set, (230 V)	285212848
TitroLine easy module 3 for halogenide titration, same as module 1, with one silver combination electrode, (230 V)	285212864
TitroLine easy module 2 for pH titration, same as module 1, with one pH-electrode and buffer set, (115 V)	285212831
TitroLine easy module 3 for halogenide titration, same as module 1, with one silver combination electrode, (115 V)	285212856
Accessories for TITRONIC[®] basic, TITRONIC[®] universal and TitroLine easy	
TZ 2005, bottle top adapter, GL 45	285221055
TZ 2008, bottle top adapter, S 40	285221088
TZ 2004, bottle top adapter GL 45, with 1 L reagent bottle, brown	285221047
TZ 3460, RS-232 printer including data cable, (230 V)	285225608
TZ 2074, TitroLine Chart for TitroLine easy	1015738

Easier Titration at any level of complexity:

TitroLine alpha *plus*

Innovative electrochemistry – from the very beginning

By developing the glass electrode more than 70 years ago, we laid the foundation for the success of electrochemical measurement. With high-performance pH glasses, innovative electrodes and electrochemical measuring instruments such as pH meters, conductometers, oxygen measuring instruments, piston burettes and titrators we have since made sure that electrochemical measurement today is an indispensable, trouble-free and reliable procedure all over the world.

Based on this know-how we have developed the automatic titrator **TitroLine alpha *plus*** which combines the ease-of-use of its predecessor TitroLine alpha with the robustness of the TITRONIC® 110 and TITRONIC® 200 precision piston burettes and exceeds the performance of the now almost legendary TPC 2000 titration system.

The right choice for simple and complex titrations

The **TitroLine alpha *plus*** is compact, flexible, very robust and universally applicable. Its capabilities range from simple end-point titrations (EP), such as the determination of the total acid in wine, to complex and difficult, non-aqueous titrations such as the determination of the acid and base numbers in oils (TAN/TBN). Of course, the automatic titrator **TitroLine alpha *plus*** is also the ideal choice for pH-stat applications such as the determination of the enzyme activity or for “dead-stop” titrations such as water determination according to Karl Fischer (KF).



Whatever your titration needs are, it will be worth your while to take a closer look at the **TitroLine alpha plus**, especially if your tasks include one of the following applications.

Environmental and water analytics

- Chloride in tap water and sewage water
- Calcium and magnesium hardness
- pH values
- Alkalinity (“p and m values”)
- Permanganate index
- COD

Foodstuff and beverages

- Salt content (NaCl) in soya sauce, cheese, ketchup, spices and other foodstuff
- Peroxide number, saponification number, iodine and acid numbers in fats and oils
- Formol number in fruit juices
- Calcium in milk products
- Ascorbic acid (Vitamin C)
- Alpha acids in hop



The addition of up to five piston burettes for dosing and titrating transforms the stand alone instrument into a team player.

Galvanics

- Determination of copper, zinc, nickel and aluminium with Cu-selective electrode
- Boric acid and chloride in nickel baths
- Alkali in degreasing baths

Petrochemistry

- Acid and base number (TAN and TBN)
- Bromine index
- Water determination according to Karl Fischer (KF)

Pharmaceutics

- Content determination of pharmaceutical products with perchloric acid in pure acetic acid
- Chloride
- Water determination according to Karl Fischer (KF)

General chemistry and plastics

- Titration of strong acids and alkaline solutions
- Epoxy number, isocyanates, acid number, hydroxyl number and saponification number
- Amino end groups
- Carboxyl end groups

Paper industry

- “White, green and black liquor”

plus

The right electrode for your titration application

The applicable electrode for the titration application is a decisive factor for the accuracy and reproducibility of the results. In order to support you with selecting the appropriate electrode, we have summarized the according electrodes for the most important applications in the following.



Application	Electrode (w.o. temp.-sensor)	Electrode with integrated. temp.-sensor
Acid-base-titrations		
Aqueous, general strong acid and bases	A 7780	–
Kjeldahl	A 7780	–
Alkalinity	N 62, N 61	N 1052 A, N 1051 A
Aqueous, difficult applications	IL-pH-A120MF IL-pH-A170MF	IL-pHT-A120MF-DIN-N IL-pHT-A170-DIN-N
Low ionic liquids	IL-pH-A120MF IL-pH-A170MF	IL-pHT-A120MF-DIN-N IL-pHT-A170-DIN-N
Small sample amounts	N 5900 A	A 157 IL-MICRO-pHT-A-DIN-N
Titration with sample changer (100 – 250 ml vessels)	N 65	N 1051 A IL-pHT-A170-DIN-N
Titration with sample changer (50 ml vessels, micro)	N 5900 A	–
Non aqueous acid base-titrations		
TAN (ASTM 664)	N 6480 eth	–
OH-No, NCO-No, FFA saponification No. ...	N 6480 eth	–
TBN (ISO 3771/ASTM 2896)	N 6480 eis	–
Epoxy value	N 6480 eis	–
Titration with perchloric acid/acetic acid	N 6480 eis	–
Precipitation titrations		
Halogenides (chloride, "salt")	AgCl 62	–
Halogenides, sample changer	AgCl 65	–
Pseudo halogenides (cyanide ...)	Ag 6280	–
Detergents	TEN 1100*	–
Redox titrations		
General, iodometric, permanganometric, cerimetric	Pt 62 Pt 6280	–
Iodine number, peroxid number	Pt 61	–
COD	Pt 61	–
Sample changer, general	Pt 6580	–
Sample changer, COD	Pt 5901	–
Dead stop (SO ₂ bromine no. ...) general	Pt 1200	–
Dead stop (SO ₂ bromine no. ...) sample changer, general and titration vessels	Pt 1400	–
Dead stop (SO ₂ bromine no. ...) sample changer micro	KF 1100	–
KF-titrations	KF 1100	–
Complexometric titrations		
Water hardness (Ca/Mg separated)	Ca 1100 A*	–
Water hardness, total	Cu 1100 A*	–
Copper, zinc, nickel, alumina ...	Cu 1100 A*	–


* An applicable reference electrode is required: B 2920+ respectively. B 3520+

TitroLine alpha *plus*: So adaptable ...

Working with the TitroLine alpha *plus* is so easy

Take a closer look at the large graphical display to see how easy it is to work with the TitroLine alpha *plus*. Everything you need to know is visible in clear text. Just press a few buttons to select the desired function: the method, the log you want, the kind of output ...

Two arrow buttons are enough to navigate you through the self explanatory menu. Use the enter button to confirm your selection and press ESC to leave a menu item. After set up, start your titration with the separate start/stop button. Parameterize your method from a connected keyboard.



Use the arrow buttons (centre keys) to navigate up and down in the menu, and confirm your selection with Enter (lower key). Use ESC (upper key) to leave a menu item.



During the titration, you can watch the entire procedure in real time at the titration curve shown on the large display. In this way, you are always in control and don't have to wait for the curve printout.

The TitroLine alpha *plus* adapts itself to your applications

For optimal adaptation of the titration to your application, the TitroLine alpha *plus* provides an extensive database with the most important titration methods pre-programmed. From this database of 100 methods you can download up to 50 methods into the free method memory and modify each method as required to meet your own specific needs. Needless to say you can also enter and save your own tried and tested methods. Our application database in the internet is a useful source of methods and informations which are available for free download.

The right titration control for any method

Reagent can be added after a fixed waiting period or drift-controlled, in linear titration steps or with dynamic adaptation to the curve slope. Additionally, you can select other forms of control for end-point titrations for pH, mV and μA , and for KF and pH-stat titrations.

Up to five equivalence points can be preselected for **turning-point titrations**, and up to three end points for **end-point titrations**.

plus

... as precise and robust as you need it.

Correct results – good documentation

To calculate the results, you can choose from eight preset formulas. Additionally, the formula editor allows the creation of your own formulas. 50 variables are available, for example to store blank values, titre and means for other calculations and applications.

You can generate your own logs to document the measured results:

The **Brief Log**

contains the result, originally weighed-in quantity, sample name, date and time.

The **Standard Log**

in addition to the above also includes the titration curve with first derivation.

The **Detailed Log**

in addition to the above also includes the calculation formula, calibration data, dates of preparation and change of the method.

The **GLP Log**

includes all titration parameters of the method.

Method link to solve complex tasks

For complex tasks, the TitroLine alpha *plus* facilitates easy combination of methods. For example, in a first method you can determine the alkalinity ("m" value) with an end-point titration to pH 4.3 with HNO_3 . Following this, a second method can be automatically started ("linked") to determine the chloride content with silver nitrate.



Unlock ...



... take off ...



... everything under control!

Most precise and robust – the exchangeable dosing units.

The TitroLine alpha *plus* is available with a choice of five exchangeable units for the reagents, with volumes of 1, 5, 10, 20 and 50 ml. The dosing cylinders of the exchangeable units are made of high-precision calibrated Duran®. This is a speciality which enables you to dose your reagents with the highest accuracy. As only highly resistant materials (PTFE/PCTFE, FEP and FPA) are used for all other wetted parts, you can use practically any measurable liquids (except HF).

Changing the reagents on the TitroLine alpha *plus* is really child's play: Simply press the unlock button on the left side of the unit, and you can remove the unit with a flick of the wrist. Thanks to the robust design, you'll always have a firm grip even on well-filled bottles.

Fitting a new unit is just as easy. Not only does the unit lock itself automatically but the logical encoding corresponding to the volume is also automatically transferred to the titrator or piston burette. There's no need to adjust the titrator. And, by the way, the units of the TITRONIC® 100, TITRONIC® 110 and TITRONIC® 200 piston burettes are compatible with each other and can be used when working with the TitroLine alpha *plus*.

Water determination according to Karl Fischer – starting at 10 ppm with the TitroLine alpha *plus* KF

Just add a few accessories and your TitroLine alpha *plus* becomes a precise KF Titrator

The TMKF KF titration stand featuring solvent addition at the push of a button and automatic discharge of titrated samples, the TZ 1770 KF titration vessel and the KF 1100 double platinum electrode will transform your TitroLine alpha *plus* into a fully functional, most powerful volumetric titrator for water determination according to Karl Fischer (KF).

KF titration parameters – exactly as required

All parameters required to optimally adapt the method to your sample are available for your KF titration: extraction time, drift stop, endpoint delay, stop current (μA), adjustable pole voltage, maximum and minimum titration time. The drift determined can be automatically corrected.

Versatile and very precise

The TitroLine alpha *plus* KF is a perfect choice for all volumetric KF applications in the fields of pharmaceuticals, chemistry, petrochemistry, foodstuff and plastics industry. The very high precision of its 5 ml and 10 ml dosing units allow determination of water contents from as low as 10 ppm with excellent reproducibility. The upper limit is 100 %.



The KF drying oven extends the range of applications

Using the TZ 1052 drying oven allows you to analyse samples which cannot be titrated directly, e.g. samples of plastic or oil containing additives.


The TitroLine alpha *plus* gets along with everyone

The TitroLine alpha *plus* is well connected

The TitroLine alpha *plus* is top as a stand-alone device but it rises above itself a team member with the support of two RS 232 interfaces, combined with the possibility of concatenation (Daisy chain) of titrator, piston burettes and sample changer and a Centronics interface.

For example, the two RS-232-C interfaces allow simultaneous connection of a PC and a balance to automatically enter the weight of the sample. On the second RS-232-C interface you can connect additional TITRONIC® type piston burettes (except TITRONIC® *basic*), a sample changer and a balance.

Your Epson and HP-compatible printer (see Technical data) can be connected to the Centronics interface.



We are happy to support you with your applications

The staff in our Application Laboratory will be glad to assist you and impart their many years of practical experience.

You can also find much of this practical experience in our application database on the internet:

www.schottinstruments.com

Support for device qualification

In connection with quality management systems, more and more importance is being given to the traceability of analysis. We support your needs with a logbook that provides you with forms for **IQ** (Installation Qualification), **OQ** (Operational Qualification) and **PQ** (Performance Qualification). Using these instruments you can effectively document commissioning, routine work and inspections of the TitroLine alpha *plus*.

TITRONIC® 110 – the piston burette with the *plus*

Titrating and dosing

The TITRONIC® 110 *plus* is the piston burette for your precise dosing and titration. It can be used as stand-alone device, in combination with a titrator and additional piston burettes or in connection with a PC.

The TR 160 manual controller allows manual titrations to visual end point or in combination with a pH meter.

Very precise and robust

The resolution of 10,000 steps, the high-precision calibrated DURAN® glass cylinders – one of our specialties – the quick and easy to change units and the high-quality workmanship make the TITRONIC® 110 *plus* a piston burette unrivalled in accuracy and robustness.

Exchangeable units compatible

By the way, the exchangeable dosing units are compatible with those of the TitroLine alpha *plus* titrator and the TITRONIC® 100, TITRONIC® 110 and TITRONIC® 200 piston burettes.

As dosing and titration burette with TitroLine alpha *plus* and TitriSoft

You can also use the TITRONIC® 110 *plus* as a dosing burette for exact dispensing of reagents, as a titration burette in combination with the TitroLine alpha *plus*, or as a dosing and titration burette within the TitriSoft titration system.



PC control and concatenation (Daisy Chain)

All functions of the TITRONIC® 110 *plus* can be controlled via a PC serial interface, so the TITRONIC® 110 *plus* can be used within systems of other manufacturers, e.g. as dosing and titration burette. For complex applications, concatenation (Daisy chain) of up to 16 devices is possible. The devices are simply connected via the second serial interface. In this way, each device can be addressed separately and reply on its own without the need for an additional data line to the PC.

Special dosing applications

With the use of PC keyboard, dosing tasks can be performed at the push of a button. You can optimize dosing and filling speed for precise measurement even with very viscous liquids such as concentrated sulphuric acid, making the TITRONIC® 110 *plus* most suitable for sample preparation in viscometry.

Technical data

TitroLine alpha *plus* and TITRONIC® 110 *plus*

Conformity	ISO 8655, mark of conformity
CE sign	
Valve	motor-driven 3/2-way valve made of PTFE/ETFE
Hoses	FEP with UV protection
Keyboard	PS2 socket for connection of a PC keyboard. Connection TZ 2825 possible with adapter
RS-232-1	PC, input for concatenation of several devices (Daisy chain)
RS-232-2	piston burettes types TITRONIC® 110, TITRONIC® 110 <i>plus</i> , TITRONIC® 200 and TITRONIC® <i>universal</i> sample changer types TW 280, TW alpha und TW alpha <i>plus</i> TitroLine alpha <i>plus</i> : balances (Mettler, Sartorius, Kern, Ohaus, others on request)
Power supply	mains: 230 V~, 50 / 60 Hz; or 115 V~, 50 / 60 Hz, power consumption: 43 VA
Housing	polypropylene
Front foil	polyester
Housing dimensions	145 x 260 x 270 mm (W x H x D), only exchangeable unit 145 x 360 x 295 mm (W x H x D) height inclusive of exchangeable unit
Weight	basic device approx. 4.1 kg, complete device with exchangeable unit approx. 5.1 kg
Climate	ambient temperature: + 10 ... + 40 °C for operation and storage
Units	1, 5, 10, 20 and 50 ml units with calibrated glass cylinder made of DURAN® (borosilicate glass) size coding allows automatic detection of unit
Burette resolution	1/10,000, smallest step 0.1 µl with 1 ml burette size
Dosing accuracy	trueness: 0.1 ... 0.3 %, referred to nominal volume (in dependence on burette size) precision: 0.05 ... 0.1 % (in dependence on burette size)

Achievable accuracies in the entire system with exchangeable unit

Exchangeable unit	Volume	Tolerances of inside diameter of the glass cylinder	Dosing error referred to 100 % volume	Reproduce- ability
TA 01	1.00 ml	± 0.003 mm	± 0.3 %	0.10 %
TA 05 <i>plus</i>	5.00 ml	± 0.003 mm	± 0.15 %	0.07 %
TA 10 <i>plus</i>	10.00 ml	± 0.003 mm	± 0.1 %	0.05 %
TA 20 <i>plus</i>	20.00 ml	± 0.003 mm	± 0.1 %	0.05 %
TA 50 <i>plus</i>	50.00 ml	± 0.003 mm	± 0.1 %	0.05 %

TitroLine alpha *plus* only

Display	matrix LCD display, 69 x 69 mm, with background illumination, contrast adjustable via keyboard
Measuring input A	pH/mV input with electrode socket in accordance with DIN 19 262/or BNC
Measuring input B	pH/mV input with electrode socket in accordance with DIN 19 262/or BNC, galvanic separated
Measuring input KF/µA	Karl-Fischer (dead-stop) connection for double-platinum electrode (connection sockets: 2 x 4 mm), polarization voltage adjustable
Measuring input Pt 1000	temperature sensor connection of resistance thermometer Pt 1000 (connection sockets: 2 x 4 mm)
Printer connection	centronics interface for connection of an Epson (ESC/P2 and Raster) and HP (PCI 3) -compatible printers

TITRONIC® 110 *plus* only

Display	LCD display, 4-digit with floating point
I/O multifunction port	15-pole sub D-socket for connection of the TR 160 manual controller for manual titration Special applications on request
Volume display	00.00 ... 9.999 ml
Indication resolution	0.000 ... 9.999 ml
Dosing volume	0.01 ... 9.999 ml
Dosing speed	0.01 ml/h ... 100 ml/min (in dependence on burette size)
Filling speed	30 ... 999 s, freely selectable

TITRONIC® and DURAN® are registered trademarks. Subject to technical changes without notice.

TW alpha *plus* sample changer – automatic titration in series

The number of samples to be processed is growing constantly while at the same time the demands on reliability are increasing in accordance with GLP and ISO 900X standards. The TW alpha *plus* sample changer by SCHOTT Instruments helps you meet these increased requirements and relieve qualified employees from routine work.

Control by titrator or by PC

You can control the sample changer from the TitroLine alpha *plus* titrator or from a PC with the Titriflex software.

Higher flexibility due to exchangeable sample racks

With four sample racks for up to 24 samples and titration head fittings for a variety of beaker and titrator vessels you get the flexibility your lab needs. A mere flick of the wrist is sufficient to change the sample racks and titrator heads. The size of the rack can be selected in the TitroLine alpha *plus* or in the ›Titration Center‹ of the Titriflex software.

Stirring from “above” or “below”

As standard, the TW alpha *plus* comes with an integrated magnetic stirrer to stir the samples from “below”. Alternatively, you can use a rod stirrer which enables stirring from “above”.



plus



Washing the electrode and the titration tip

To ensure accuracy of the results, the electrodes and the titration tips are rinsed after each titration. This can, for example, be done by immersing the electrodes and titration tips into a wash-ing solution. The number of rinsing positions to be used (up to a maximum of three) and the rinsing time are set in the method. Direct and fast rinsing of the electrodes and titration tips can be ensured by using the MP 25 washing unit that rinses directly after the titration. In addition to this, a waiting position may also be used for example to immerse the pH electrodes into a KCl solution.

Up to 24 samples in 50 ml glass beakers or 16 samples in 250 ml glass beakers will fit in the rotating sample tray. A sample tray for 24 COD containers is also available.

TitriSoft 2.6 – convincingly simple ...

The TitriSoft 2.6 titration software is the optimum solution for your titration tasks. The software can be used with WINDOWS 98/ME and WINDOWS 2000/XP/Vista and supports your daily work procedures during sample preparation, titration and evaluation of the results. The software has been developed to be clear, logical and user-friendly.


Connection possibilities

Using TitriSoft 2.6 you can control the following devices from a PC:

- **Titrators** (TitroLine alpha *plus*, TitroLine alpha, TR 250)
- **Sample changers** (TW alpha *plus*, TW alpha, TW 280)
- **Piston burettes** (TITRONIC® 110 *plus* and TITRONIC® universal, TITRONIC® 110, TITRONIC® 200)
- **Balances**

You can connect the titration hardware to any of your PC's available serial interfaces. Each of the serial interfaces allows different combinations of devices (configurations). To automate a titration procedure the software may be used to control the TitroLine alpha *plus* in connection with the TW alpha *plus* sample changer. For more complex titration tasks with sample preparation you can dose with piston burettes followed by titration with a TitroLine alpha *plus*. Of course, you can also use the software for dosing only.

The image below shows possible device configurations.



System requirements

For optimal and fast working with the TitriSoft 2.6 software your system should be equipped as shown below:

Interface: 1 free serial RS-232-C interface per configuration

Computer: Pentium II or higher

Operating system: WINDOWS 98/ME, WINDOWS 2000/XP or Vista

RAM: minimum 256 MB

Hard disk: minimum free storage place 100 MB

Graphics card: minimum resolution 1024 x 768

... strong benefits ...



›Navigator‹, the main menu

The different software tasks are assigned to four different centers:


- the Maintenance Center,
- the Revision Center,
- the Analysis Center and
- the Titration Center.

The centers can be accessed from the main menu, the Navigator.


›Maintenance Center‹, the system configuration

In the Maintenance Center, the software is set up for operation prior to running the first application, i.e. a configuration is set up with the connected hardware. The configuration of the attached hardware is automatically detected in a hardware scan. Each of these hardware configurations allows any number of "methods" and "work lists". Different configurations can work in parallel (see Connection Possibilities).

All TitriSoft users can be listed by their names. TitriSoft supports three user types. The Administrator has access to all configuration and software operation options. The "Administrator" has access to all configuration and software operation options. The "User" or "Advanced User" has the same rights as the Administrator but is not allowed to delete results, methods and worklists. Users are restricted to operation of the Titration Center which very much simplifies matters.



... clearly structured ...




»Analysis Center», your method center

This is where you set up and save your titration methods. Even complex methods can be installed with a few mouse clicks. Adjustment of the titration parameters is facilitated by the use of symbolic slide controls. Functions such as waiting time, IF loops, repetition, do-sings and measurements in addition to the titration parameters and calculation formulas provide virtually unlimited options for method procedures.


»Revision Center», your database

Titration curves, results, measured values and used methods of all titrations ever carried out are stored in the database. These data can be selected by sample name, date, user and method and loaded in a few seconds.

Information on titrations performed can be displayed in the form of a diagram, results list or measured value list. You can optimize stored titration information in accordance with your requirements, e.g. add and store subsequent calculations or analyze titration curves and print it out together. Additionally, an export of the data to Excel and ASCII is also available.



... highly productive: TitrISOft 2.6



›Titration Center‹, your clearly structured workplace

The ›Titration Center‹ is the place where you carry out your daily jobs, i.e. select methods, enter sample names and origin weighed-in quantities, start the work list and display (and print if desired) the results at the end of a titration. The work list shows the individual samples with the associated methods and their characteristics such as sample name, number, status, date, time, results and events and other freely configurable sample data, e.g. density.

During the titration you can observe the titration process in an on-line diagram. You can, however, simply allow the samples to be processed in the background and use the PC for other tasks or start an additional titration with another configuration in parallel.

When working with the TW alpha *plus* sample changer, you can adjust various settings such as skip empty items, rinse and waiting options.

Documentation, which is in accordance with GLP and ISO 9000 directives, can be produced in a number of different forms; tables, lists, curves or individual print-outs with curves. In addition results can be saved in ASCII or CSV format, external documentation programs may be accessed and results transferred directly, e.g. into a LIMS.

TitriSoft 2.6 P – simply safe ...

In this case, the “P” does not mean “professional”, but rather “pharmaceutical”. Nevertheless, the performance of the “P” version is of course just as good and as professional as that of the standard version “TitriSoft 2.6”. Additionally to the standard version, the TitriSoft 2.6 P fully meets all requirements of the FDA 21 CFR Part 11 regulation regarding „Electronic Records“, „Electronic Signature“ and „Audit Trail“.

The FDA (i. e. Food and Drug Administration of the USA) 21 CFR Part 11 regulations describe how to deal with electronically stored data („Electronic Records“) and how to prepare electronic signatures („Electronic Signature“). These regulations are binding for all companies offering medical, pharmaceutical or food products and services in the USA.

System requirements

For optimal and fast working with the TitriSoft 2.6 P software your system should be equipped as shown below:

Interface: 1 free serial RS-232-C interface per configuration

Computer: Pentium II or higher

Operating system: WINDOWS 2000 and XP Pro/Vista

RAM: minimum 256 MB

Hard disk:
minimum free storage place 100 MB

Graphics card:
minimum resolution 1024 x 768

Comparison between TitriSoft 2.6 and 2.6 P


Functions	TitriSoft 2.6	TitriSoft 2.6 P
Electronic Record		■
Electronic Signature		■
Audit Trail		■
Controlled Access		■
Copies of Records		■
Manual with forms for SOP's, IQ, OQ, PQ and validation reports		■
Straightforward procedure	■	■
All types of titrations	■	■
Comfortable worklists	■	■
Online titration curves	■	■
Clear documentation	■	■
Perfect titration control by PC	■	■

TitriSoft

Controlled Access

The controlled access guarantees that only authorized individuals have access to the software functions, according to your company's security policy and the FDA requirements.

TitriSoft 2.6 P has 5 different access levels: The "Operator" level does only allow to carry out the routine titrations, whereas the "Advanced User" level is entitled to approve the methods. The highest level, the "Administrator" may set up the users and assign them the user rights. He even has the permission to delete records, but only after a copy of the database has been generated.



Titration controller: TitrSoft 2.60 - [Maintenance Center]

Maintenance Center

Service Help Logout Change password

Navigation: Configurations Users System traces Audit trail

Audit trail Built trail to pdf


ID	User	Date	Time	Timezone	Requested action	Description	Comment
720	ad: Administrator	18.03.2008	18:45:15	GMT +1	Added	Analysis done and measuring time added to configuration 1	new method
712	ad: Administrator	18.03.2008	20:00:27	GMT +1	Added	Analysis time measure added to configuration 1	new method
688	ad: Administrator	13.03.2008	21:48:51	GMT +1	Added	Analysis demo method Pitron 2006 added to configuration 1	new method British
279	Christoph Kapteina	18.03.2008	11:21:13	GMT +1	Added	Analysis pit-Trend 20 added to configuration 1	new method
701	ad: Administrator	15.03.2008	15:47:28	GMT +1	Edited	An analysis with name pit-Trend has been updated in configuration 1	no changes
801	Stefan Kaus	20.03.2008	08:45:35	GMT +1	Edited	An analysis with name time measure has been updated in configuration 1	review from Kaus
443	Christoph Kapteina	18.03.2008	12:59:34	GMT +1	Edited	An analysis with name Acid Base with NaOH has been updated in configuration 1	save
406	Christoph Kapteina	18.03.2008	13:31:54	GMT +1	Edited	An analysis with name Acid Base with NaOH has been updated in configuration 1	save
481	Christoph Kapteina	13.03.2008	11:35:08	GMT +1	Edited	An analysis with name Redox with Thiosulfate has been updated in	save
513	Christoph Kapteina	13.03.2008	14:22:37	GMT +1	Edited	An analysis with name Ca and Mg with EDTA has been updated in configuration 1	save
461	Christoph Kapteina	13.03.2008	09:22:53	GMT +1	Edited	An analysis with name Calibration fluoride has been updated in	save new conc.
518	Christoph Kapteina	13.03.2008	14:27:58	GMT +1	Edited	An analysis with name Ca and Mg with EDTA has been updated in configuration 1	use end value
442	Christoph Kapteina	18.03.2008	12:58:37	GMT +1	Edited	An analysis with name Acid Base with NaOH has been updated in configuration 1	indexing

• Table : Configuration
 • Table : Electrode
 • Table : Measurement
 • Table : prog_user
 • Table : property
 • Table : worklist

ad: Administrator

Audit Trail

The 21 CFR Part 11 prescribes that each creating, saving oder modifying of records (e.g. creating methods, modifying passwords or saving results), generates an entry in the Audit Trail. TitrSoft 2.6 P automatically generates an entry in the Audit Trail table as soon as an access to the database has taken place. The local time and the GMT are automatically stored together with this entry in the Audit Trail. Each entry also asks for a comment. The Audit trail or parts of it can be printed out, or a "human" readable digital copy of it, e.g. a PDF file can be generated.



Electronic Records


The 21 CFR Part 11 prescribes how to safeguard and store the generated results over a long timespan. Besides regularly making backup copies of the complete database, is it possible to generate readable digital copies of the results, methods, worklists, the Audit Trail, the user administration and the configuration(s). For that purpose, a PDF writer is already integrated in the software. The purchase of expensive third-party software for generating PDF files ist not necessary.

Of course the database is password protected against unauthorized access.

Tit

Electronic Signature

Digital analysis results have to be as reliable as classical, manually checked results with a handwritten signature. A digital signature, which is as safe as a handwritten one, can be placed to approve all electronic records. The approver has to enter the name and an additional password. The electronic signature is stored together with the signer's function, the reason of signing and the date and time.



Ordering information TitroLine alpha *plus*

TitroLine alpha <i>plus</i>		Order no.
TitroLine alpha <i>plus</i>	TitroLine alpha <i>plus</i> basic unit <u>without</u> exchange unit, 230 V	285216952
TitroLine alpha <i>plus</i>	TitroLine alpha <i>plus</i> basic unit <u>without</u> exchange unit, 115 V	285216969

Scope of delivery: TitroLine alpha *plus* incl. stand rod with holder, titration clamp, PC keyboard TZ 2835.

TitroLine alpha <i>plus</i>	TitroLine alpha 05 <i>plus</i> with 5 ml exchange unit, (230 V)	285212934
TitroLine alpha <i>plus</i>	TitroLine alpha 10 <i>plus</i> with 10 ml exchange unit, (230 V)	285216944
TitroLine alpha <i>plus</i>	TitroLine alpha 20 <i>plus</i> with 20 ml exchange unit, (230 V)	285216977
TitroLine alpha <i>plus</i>	TitroLine alpha 50 <i>plus</i> with 50 ml exchange unit, (230 V)	285212983
TitroLine alpha <i>plus</i>	TitroLine alpha 05 <i>plus</i> with 5 ml exchange unit, (115 V)	285215467
TitroLine alpha <i>plus</i>	TitroLine alpha 10 <i>plus</i> with 10 ml exchange unit, (115 V)	285215475
TitroLine alpha <i>plus</i>	TitroLine alpha 20 <i>plus</i> with 20 ml exchange unit, (115 V)	285215631
TitroLine alpha <i>plus</i>	TitroLine alpha 50 <i>plus</i> with 50 ml exchange unit, (115 V)	285215648

Scope of delivery: As TitroLine alpha *plus* basic unit with 5, 10, 20 or 50 ml exchange unit, incl. brown glass bottle for titrant, GL 45 bottle adapter, hoses, drip glass and titration tip.

TitroLine alpha KF *plus*

TitroLine alpha <i>plus</i>	TitroLine alpha KF 05 <i>plus</i> with 5 ml exchange unit, (230 V)	285212991
TitroLine alpha <i>plus</i>	TitroLine alpha KF 10 <i>plus</i> with 10 ml exchange unit, (230 V)	285213109
TitroLine alpha <i>plus</i>	TitroLine alpha KF 05 <i>plus</i> with 5 ml exchange unit, (115 V)	285215656
TitroLine alpha <i>plus</i>	TitroLine alpha KF 10 <i>plus</i> with 10 ml exchange unit, (115 V)	285215664

Scope of delivery: As TitroLine alpha *plus* basic unit with 5 or 10 ml exchange unit, incl. brown glass bottle for titrant, GL 45 bottle adapter, hoses, drip glass and titration tip, titration stand TMKF incl. supply and waste bottle, micro-double-platinum electrode KF 1100, titration vessel TZ 1770.

TITRONIC® 110 *plus*

TITRONIC® 110 <i>plus</i>	TITRONIC® 110 <i>plus</i> basic unit <u>without</u> exchange unit, 230 V	1007302
TITRONIC® 110 <i>plus</i>	TITRONIC® 110 <i>plus</i> basic unit <u>without</u> exchange unit, 115 V	1007303

Scope of delivery: TITRONIC® 110 *plus* incl. stand rod with holder and titration clamp

Exchange units for TitroLine alpha *plus* and TITRONIC® 110 *plus*

TA 01	Exchange unit with 1 ml glass cylinder incl. reagent bottle	285211313
TA 05 <i>plus</i>	Exchange unit with 5 ml glass cylinder incl. reagent bottle	285211038
TA 10 <i>plus</i>	Exchange unit with 10 ml glass cylinder incl. reagent bottle	285211046
TA 20 <i>plus</i>	Exchange unit with 20 ml glass cylinder incl. reagent bottle	285211054
TA 50 <i>plus</i>	Exchange unit with 50 ml glass cylinder incl. reagent bottle	285211062

Software TitriSoft

TitriSoft 2.6	Titration software for TitroLine alpha <i>plus</i> and TitroLine alpha	285221717
TitriSoft 2.6 P	Titration software according to CFR 21 Part 11	285221720

Accessories for TitroLine alpha *plus* and TITRONIC® 110 *plus*

TZ 2835	PC keyboard for TitroLine alpha <i>plus</i> and TITRONIC® 110 <i>plus</i>	1007852
TM 135	Magnetic stirrer	285211013
TM 128	Titration clamp/rod stirrer combination	285215167
TMKF	Titration stand Karl-Fischer with suction pump and stirrer incl. supply and waste bottle	285216611
TZ 1770	KF titration vessel 30 - 150 ml	285216677
TZ 1772	KF titration vessel 80 - 200 ml	285216693
TZ 1052	Drying oven for water determination according to Karl-Fischer, 230 V	285214721
TZ 1060	Accessory for drying oven TZ 1052	285218115
Z 303	Titration clamp for TL alpha <i>plus</i> and TITRONIC® 110 <i>plus</i>	1007304

Ordering information TW alpha *plus*

Sample changer TW alpha <i>plus</i>		Order no.
TW alpha <i>plus</i>	Sample changer basic unit, 230 V	1007290
TW alpha <i>plus</i>	Sample changer basic unit, 115 V	1007291

Scope of delivery: Sample changer basic unit TW alpha *plus* with integrated magnetic stirrer and connection cable TZ 1581 for rod stirrer.

TW alpha <i>plus</i> 12	TW alpha <i>plus</i> basic unit with sample rack TZ 1452 for 12 samples, incl. titration head TZ 1463, connection cable and 20 beakers 250 ml, 230 V	1007292
TW alpha <i>plus</i> 16	TW alpha <i>plus</i> basic unit with sample rack TZ 1459 for 16 samples, incl. titration head TZ 1463, connection cable and 20 beakers 150 ml, 230 V	1007294
TW alpha <i>plus</i> 24	TW alpha <i>plus</i> basic unit with sample rack TZ 1454 for 24 samples, incl. titration head TZ 1462 and 30 beakers 50 ml, 230 V	1007296
TW alpha <i>plus</i> COD	TW alpha <i>plus</i> basic unit with sample rack TZ 1444 for 24 COD vessels in accordance with DIN, incl. titration head TZ 1461, rod stirrer TZ 1846, redox electrode Pt 5901, titration tip TZ 1648 and connection cable, 230 V	1007298
TW alpha <i>plus</i> MP	TW alpha <i>plus</i> basic unit with sample rack TZ 1459 for 16 samples, incl. titration head TZ 1467, membrane pump MP 25, connection cable and 20 beakers 150 ml, 230 V	1007305
TW alpha <i>plus</i> 12	TW alpha <i>plus</i> basic unit with sample rack TZ 1452 for 12 samples, incl. titration head TZ 1463, connection cable and 20 beakers 250 ml, 115 V	1007293
TW alpha <i>plus</i> 16	TW alpha <i>plus</i> basic unit with sample rack TZ 1459 for 16 samples, incl. titration head TZ 1463, connection cable and 20 beakers 150 ml, 115 V	1007295
TW alpha <i>plus</i> 24	TW alpha <i>plus</i> basic unit with sample rack TZ 1454 for 24 samples, incl. titration head TZ 1462 and 30 beakers 50 ml, 115 V	1007297
TW alpha <i>plus</i> COD	TW alpha <i>plus</i> basic unit with sample rack TZ 1444 for 24 COD vessels in accordance with DIN, incl. titration head TZ 1461, rod stirrer TZ 1846, redox electrode Pt 5901, titration tip TZ 1648 and connection cable, 115 V	1007299
TW alpha <i>plus</i> MP	TW alpha <i>plus</i> basic unit with sample rack TZ 1459 for 16 samples, incl. titration head TZ 1467, membrane pump MP 25, connection cable and 20 beakers 150 ml, 115 V	1007306

Accessories for TW alpha *plus*

TZ 1444	Sample rack for 24 COD vessels in accordance with DIN 38 409	285213836
TZ 1452	Sample rack for 12 sample vessels, incl. 20 beakers 250 ml	285214927
TZ 1454	Sample rack for 24 sample vessels, incl. 30 beakers 50 ml	285213844
TZ 1459	Sample rack for 16 sample vessels, incl. 20 beakers 150 ml	285213166
TZ 1461	Titration head for COD sample rack TZ 1444	285213621
TZ 1462	Titration head for 24-sample rack TZ 1454	285213639
TZ 1463	Titration head for 12- (TZ 1452) and 16-sample rack (TZ 1459)	285213647
TZ 1467	Titration head for 12- (TZ 1452) and 16-sample rack (TZ 1459) incl. splash shield with membrane pump MP 25	285213671
MP 25	Membrane pump MP 25 with accessories for rinsing, 230 V	285216005
TZ 1847	Glass stirrer rod for 12-, 16- and 24-sample rack	285215175
TZ 1846	Glass stirrer rod for COD sample rack	285215134
TZ 1545	Magnetic stirring rods (10 pcs.)	285214232

Data cable

TZ 3088	Data cable TitroLine alpha <i>plus</i> , TW alpha <i>plus</i> or TITRONIC® 110 <i>plus</i> ↔ PC, 5 m	1007972
TZ 3089	Data cable TitroLine alpha <i>plus</i> , TW alpha <i>plus</i> or TITRONIC® 110 <i>plus</i> ↔ PC, 10 m	1007973
TZ 3084	Data cable TitroLine alpha <i>plus</i> , TW alpha <i>plus</i> , TITRONIC® 110 <i>plus</i> ↔ TitroLine alpha <i>plus</i> , TW alpha <i>plus</i> , TITRONIC® 110 <i>plus</i> , 1.5 m	1007974
TZ 3086	Data cable TitroLine alpha <i>plus</i> , TW alpha <i>plus</i> , TITRONIC® 110 <i>plus</i> ↔ TitroLine alpha, TW alpha, TITRONIC® 110/ TITRONIC® 200, 1.5 m	1007975
TZ 3087	Data cable TitroLine alpha <i>plus</i> , TW alpha <i>plus</i> , or TITRONIC® 110 <i>plus</i> ↔ TITRONIC® universal, 1.5 m	1007976
TZ 3082	Data cable TitroLine alpha <i>plus</i> , TW alpha <i>plus</i> , or TITRONIC® 110 <i>plus</i> ↔ Sartorius balances, 5 m	1007977
TZ 3083	Data cable TitroLine alpha <i>plus</i> , TW alpha <i>plus</i> , or TITRONIC® 110 <i>plus</i> ↔ Mettler AT, PM balances, 5 m	1007978
TZ 3081	Data cable TitroLine alpha <i>plus</i> , TW alpha <i>plus</i> , or TITRONIC® 110 <i>plus</i> ↔ Mettler AB-S, PG balances, 5 m	1007979

Informations on electrodes for titration and other sensors see section "Laboratory electrodes"

Content Spectrophotometers

Selection table Spectrophotometers	Page 141
PRIM	Page 142
UviLine	Page 144
Uvi Light	Page 152
UVIKON	Page 156



Selection table Spectrophotometers

Technical Specifications	PRIM Light/ PRIM Advanced	Uvi Light XS/ Uvi Light XT	UviLine 9100/ UviLine 9400	UVIKON XS/ UVIKON XL
Wavelength range	VIS	UV-VIS	VIS/ UV-VIS	UV-VIS
Technique	single beam	single beam	single beam	double beam
Display	2 lines	2-lines VGA internal/external and PC	QVGA	PC
Data storage		■ (only with XT and PC)	■	■
Methods	50 methods	50 methods XS 250 for XT and unlimited for PC	100 methods; unlimited with external data storage	unlimited
Interfaces	RS232	RS232 and Centronics for XT	2 x USB and RS232	RS232 for PC
Cell changer		■ (XT and PC only)	■	■
Sipper	■ (manual)	■	■	■
Absorbance / Transmission	■	■	■	■
Concentration	■	■	■	■
Spectra scanning	■ (Advanced only)	■	■	■
Kinetics	■ (Advanced only)	■	■	■
Multi wavelength	■ (Advanced only)	■	■	■



Compact Visible Spectrophotometers

PRIM Light and PRIM Advanced

PRIM Light & PRIM Advanced spectrophotometers combine a high level of performance with a simple and intuitive user interface. Compact and light, these new spectrophotometers are ideal for standard applications in education or in the laboratory.



Internal Applications

All applications, basic or standard, are included as standard and are immediately available on each spectrophotometer.

50 User Methods

All PRIM spectrophotometers can store up to 50 methods in memory. Stored applications can be recalled at any time without re-programming the method parameters. Prior to running a method, it is possible to check the stored parameters using only one button from the navigator keys.

Real-time Display

The display of absorbance and transmission is processed in real-time with any measurement mode.

1 Key = 1 Function

The keyboard is clearly arranged with each button corresponding to one specific single function. Operation is therefore simpler and faster.

Safety – Low Voltage Power Supply

Mandatory for the area of education and safeness of trainees and of course also applicable for any laboratory.

Compact

Light at only 2.5 kg and compact A4 size footprint, PRIM can be handled very easily.

Internal Calibration Filter

In order to certify accurate and repeatable results, the spectrophotometer calibrates itself automatically at each start-up using the internal didymium filter; a complete report is automatically printed on the external printer if connected.

Advantages
PRIM

Choose between 2 PRIM units:

PRIM Light:

Basic internal software including standard spectrophotometric applications in absorbance, transmission and single standard concentration.

PRIM Advanced:

Advanced applications for absorbance, transmission, multi-standard concentration, kinetics, multi-wavelength and spectrum scanning.

A wide choice of measuring modes

Kinetics

- Analysis of absorbance variation as a function of the time.
- Program to include lag time and reaction time.
- Automatic calculation of absorbance variation during each segment of time or the total time.

Spectrum Scan*)

- A curve of absorbance values as a function of the wavelength with detection of absorbance maxima and minima.
- User-definable scan mode using the entire or partial visible range: 330 to 900 nm, in 1 nm steps, with automatic recording of the baseline.

Multi-Wavelength Mode

- Measure of ratio and difference of absorbance at 2 wavelengths.
- Simultaneous display of calculation results as well as individual absorbance values at each wavelength

*) on optional external printer.

Accessories

A complete range of accessories:

Printer, tube-holder, manual aspiration system, thermostatted cuvette holder are available on option.

Technical Specifications	PRIM Light	PRIM Advanced
Spectral range	330 – 990 nm	
Bandwidth	10 nm	
Accuracy	± 1.5 %	
Precision	± 1 nm	
Photometric range	–0.3 ... 2.5 Abs, 0 ... 200 %T	
Accuracy	± 2 %	
Drift	< 0.003 A/h @ 500 nm	
Stray Light	0.5 % T @ 340 & 400 nm	
Display	Alphanumeric, LCD, back-lit, 2 lines, height 8 mm, 16 characters	
Zero	Automatic	
Light source	Halogen	
Detector	Silicon diode	
Interface	Serial RS232C	
Cell holder	1 cuvette 10 mm	
Power	115/230V ... 50/60Hz	
H x W x D, Weight	180 x 280 x 220 mm, 2.5 kg	
Software equipment		
Absorbance	■	■
% Transmission	■	■
Concentration with factor	■	■
Concentration with 1 standard	■	■
Concentration with 1 to 8 standards		■
Kinetics		■
Multi-wavelengths		■
Spectrum scanning		■
Peaks and valleys detection		■
Multi-language	■	■
Automatic stand-by	■	■

Type no.	Order no.	Product	Description
PRIM Light	285600190	Spectrophotometer	VIS Spectrophotometer with 10 mm cell holder
PRIM Advanced	285600200	Spectrophotometer	VIS Spectrophotometer with 10 mm cell holder

UviLine 9100/9400:

The new series of Spectrophotometers from SCHOTT Instruments

Two complete new high-performance single beam photometer:

UviLine 9100

for measurements at VIS range from 320 – 1100 nm

UviLine 9400

for measurements at UV-VIS-range from 190 – 1100 nm

The new spectrophotometers from SCHOTT Instruments convince with an exceptional wide performance range.

- Measurement of absorbance and transmission
- Concentration measurements up to 8 standards
- Multi-wavelength
- Spektra scan with online graphic
- Kinetics

And additionally:

Extensive evaluation functions

The option of storing more than 100 own methods and up to 1000 measuring values/saving approx. 30 spectra (extendable via USB!)



UviLine



“Hard” facts:

Large backlit graphic display

The large display enables a clearly structured navigation and enhances i.e. the graphic evaluation using arrow keys and zoom function.


Clear and well-arranged keyboard

With four function keys F1 – F4 for fast and direct function triggering, such as settings, tools, zoom etc.

Alphanumeric keys for entering the wavelengths, sample ID's and other parameters.

Specially separated cursor section with ESC and START/ENTER keys for fast navigation.

Also 5 additional keys with direct functions such as i.e. print, zero (reference), save etc.



Modern interfaces

Besides the USB slave interface (USB-B) for connecting a PC, both UviLines also feature a USB master interface (USB-A). This USB master interface enables the connection of various instruments via USB:

- USB memory sticks and USB hard disks. This enables an easy extension towards the internal memory volume and also facilitates the data exchange (bi-directional).
- Standard printers with USB interface
- External PC keyboard
- A RS232 interface is also included

Powerful optics

The UviLine comprises an excellent optical design in this range:

- 1200 lines/mm holographic, concave grating for minimal stray light
- wide wavelength range from 190 – 1100 nm
- high optical resolution of 4 nm
- automatic compensation of ambient light
- automatic wave length calibration

Wide measuring range and extensive accessories

The wide measuring range enables the usage of cell holders with a length of max 100 mm. This also enables measurements with very low concentrations. There are 4 different cell holders/sample changers available, which can be replaced very easily:


- single cell holder, 10 mm (is included in the scope of delivery), thermostatable
- single cell holder 100 mm
- single cell holder, 10 mm thermostatable (peltier)
- automatic sample changer, 5+1



“Soft” facts

Besides the standard applications such as absorption, transmission and concentration measurement, both UviLine models comprise the complete functionality for the spectra, kinetics and multi wavelength analytics.

Home (admin)	02/25/08 15:09
<div>Concentration</div> <div>Absorbance / % Transmission</div> <div>Multi wavelengths</div> <div>Spectrum</div> <div>Kinetics</div>	
General setup	Logout AQA Info




With the **multi wavelength analytics** up to 4 absorption values with different wavelengths can be measured and saved. For each application you have the possibility to programme an own optimized method. By using own programmed methods, the evaluation of the measuring results becomes easier.

The determination of the reaction time values is enabled through the **kinetics** programme section. The kinetics measurement can be stored and in the following evaluated using the cursor.


Edit method (5 of 6)	02/25/08 15:09
Number: 2002 Name: PROT Version: 1.0 Citationform: Protein Unit: mg/ml Resolution: 0.1 Cell: 10 mm	
Back	Next

Edit method (6 of 6)	02/25/08 15:09
Function: $1.550 * A(280 \text{ nm}) - 0.757 * A(260 \text{ nm})$ R= 1.000	
Back	Next



Whenever the optimal wavelength for the determination of the concentration is unknown or the purity of a solution is to be determined, the adding of a **spectrum** becomes necessary. By using the zoom function and the cursor it is quite easy to edit the spectra.

Furthermore, there is a whole range of additional evaluation functions such as min/max recognition, add spectrum and subtract spectrum, derivative calculation (1-3. derivative), peak surface calculation, as well as multiply fixed value.



Technical Specifications

Software

Concentration	From 0 to 10 standards with interpolations. Graphic calibration curve management.
Kinetics	Dynamic graphic curve display, graphic management: zoom, slope calculation, current Abs
Spectrum scanning	Dynamic graphic curve display, graphic management: zoom, derivative, current Abs, Maximum and Minimum values
Multi Wavelength	Up to 10 wavelengths – flexible calculation of results
GLP compliant	User login with 3 levels
Storage capacity	Internal: 100 methods/30 graphics/1000 Data – with USB stick: limited to the size of the USB Stick
Methods	> 100 methods

Spectrophotometer

Technical Specifications	UviLine 9100	UviLine 9400
Wavelength range	320 – 1100 nm	190 – 1100 nm
Light source	Halogen	Xenon
Technique	single beam	
Optical resolution	4 nm	
Wavelength accuracy	± 1 nm	
Wavelength repeatability	< ± 0,2 nm	
Photometric range	-3.3 to 3.3 A	
Photometric accuracy	0.3 % or ± 0.003 A (from 0 – 0,6 A)	
Photometric linearity	< 1% at 2 Abs between 340 – 900 nm	
Stray light	< 0.1 % at 340 and 400 nm	< 0.1 % at 220, 340 and 400 nm
Display	graphic display with backlit, 320x240 pixel	
Update	via internet and USB Stick	
Interfaces	1 x USB-A, 1 x USB-B, 1 x RS232C	
Power supply	110 – 220 V, 50/60 Hz	
Temperature range	Use: 10 °C to 35 °C / Storage -25 °C to 65 °C	
Dimensions	404 x 197 x 314 mm (W x H x D)	
Weight	4 kg	
Accessories	Sipper, 5+1 cell changer, cell holder (thermostatable), cell holder 100 mm	

Ordering Information

Type No.	Order No.	Product	Description
UviLine 9100-EU	285700100	Spectrophotometer	UviLine 9100, single beam spectro photometer, 4 nm with a measuring range from 320 to 1100 nm and EU-power cord
UviLine 9100-US	285700110	Spectrophotometer	UviLine 9100, single beam spectro photometer, 4 nm with a measuring range from 320 to 1100 nm and US-power cord
UviLine 9400-EU	285700120	Spectrophotometer	UviLine 9400, single beam spectro photometer, 4 nm with a measuring range from 190 to 1100 nm and EU-power cord
UviLine 9400-US	285700130	Spectrophotometer	UviLine 9400, single beam spectro photometer, 4 nm with a measuring range from 190 to 1100 nm and US-power cord
SZ 2100	285700200	Automatic cell changer	5+1 cell changer, controlled from UviLine
SZ 2110	285700210	Cell holder	Single, 10 – 100 mm
SZ 2120	285700220	Cell holder	Single, 20 mm
SZ 2130	285700230	Cell holder	Single, 10 mm, (included in delivery of all UviLine)
SZ 2140	285700240	Cell holder	Single 10 mm, thermostatable
SZ 2150	285700250	Sipper	controlled from UviLiner
SZ 2160	285700260	Cell holder	Single 10 mm, thermostatable (Peltier)



Uvi Light –

The universal single beam instruments

▲ A flexible range of instruments.

Any application can find a solution with one of the 7 UviLight models:


- Powerful software for any use level
- Huge cell compartment for any kind of sample
- Large range of accessories

▲ The best performance

Accurate results thanks to SCHOTT Instruments' best technological choices.

▲ Designed for most applications:

- Agro-food
- Biotechnology
- Chemistry & petro-chemistry
- Control labs
- Clinical biology
- Research
- Teaching



The best technological choices:

Smart electronic design

- New electronic board design
- Only one control card
- Fewer contactors, fewer connectors
- Less background noise
- Higher energy output

Large Cell Compartment

- More space for accessories settings
- Position for automatic 9-cells changer
- Drain for liquid spillage

Compact & High Perf. Monochromator

- New design with small footprint
- Holographic concave with 1200 lines/mm
- Sealed against dust

Halogen and Deuterium Lamps

- Preadjusted lamps, no tools required for replacement or adjustment
- Very continuous signal over wavelength for spectrum scanning efficiency

Seiya Namioka Monochromator

- 1200 l/mm holographic concave grating:
- High energy
- Low stray light

Light Beam Chopper


- Gives fixed light frequency
- Eliminates stray light from other sources

Two bandwidths


5 nm for the best energy, 2 nm for high-resolution applications

For the best performances


Linearity 0.99 % up to 3 Abs @ 250 nm and 540 nm




Wavelength accuracy $< \pm 1.0$ nm



Base line reproducibility ± 0.002



Stray light better than 0.05 % @ 200 and 340



A wide range of models and accessories

A choice of 3 different software/interfaces:



XS: ready to use software for routine applications

- Absorbance/Transmittance
- Concentration up to 8 standards
- Spectrum scanning
- Kinetics
- 50 methods programmable



XTD: Fully programmable advanced applications

- Absorbance/Transmittance, permanent reading
- Multi-standard concentration with graphic functions
- Spectrum scanning with calculation & graphic

- Kinetics with calculation & graphic
- Unlimited methods programmable with hard disk & floppy
- Macro-commands & programmable software
- Methods and data export



PC: Enter the Windows world

- Absorbance, permanent reading
- Multi-standard concentration with graphic functions
- Spectrum scanning with calculation & graphic
- Kinetics with calculation & graphic
- Unlimited methods programmable with hard disk
- Methods and data transfer

A large range of accessories

SZ 1520 9 positions 10 mm cell holder

Accepts up to 9 standard 10 mm optical path cells. Requires automatic cell changer SZ 1550.

SZ 1550 Automatic cell changer for 9 positions cell holder

For use with 9-cell holder SZ 1520 – quick and simple installation with 1/4 turn clips.

SZ 1580 Aspiration system with 80 µl quartz cuvette

Automatic sample aspiration system delivered with peristaltic pump, 80 µl quartz flow-through cell and tubing kit.

SZ 1010 LabPower software

The multiplatform Windows NT/2000 or XP compatible software package, LabPower provides a convenient and easy-to-use well proven software environment. Basic routine applications including control of UviLight accessories.

SZ 1910 Serial Black & white KYOLINE thermal printer

40 columns with cable for Uvi Light XS.

SZ 1915 KYOLINE thermal printer

40 columns with cable for Uvi Light XTD.

Technical Specifications	
Band width	2 or 5 nm
Spectral range	190 – 900 nm
Accuracy	± 1 nm
Reproducibility	± 0.1 nm
Display resolution	1 nm (0.1 nm in spectral mode)
Photometric range	-0.1 to 3,000 A, 0.1 to 110 % T
Accuracy	± 0.005 A
Drift	< 0.0003 A/h at 500 nm
Stray light	< 0.4 % at 200 nm < 0.005 % at 220 and 340 nm
Base line stability	± 0.002 A
Memories	50 methods on XS, unlimited on XT/XTD and PC
Light source	Visible tungsten-halogen pre-adjusted lamp UV deuterium pre-adjusted lamp
Monochromator	Concave holographic grating
Resolution	1200 lines/mm
Detector	Silicon diode
Interfaces	RS232C & Parallel
Power	115 V/60 Hz to 230 V/50 Hz
Dimensions	(H x L x D) 295 x 500 x 415 mm (11,6 x 19,7 x 16,3 inches)
Weight	7.5 kg (15 lbs)

Uvi Light Models	XS	XTD/ XT Color	PC
Range	UV-VIS	UV-VIS	UV-VIS
Optics	Single beam	Single beam	Single beam
WL, nm	190 – 900	190 – 900	190 – 900
BW, nm	5 or 2	5 or 2	5 or 2
Display	2 lines	XTD: external XT: internal	PC
Software	Abs/Trans, Concentration, Wavelength scan, Kinetics, Multi-wavelength (2), Storage 50 methods	Abs/Trans, Concentration, Wavelength scan with graphic, Kinetics with graphic, Multi-wavelength (10), Storage 250 methods, Storage of results, Macro commands	Abs/Trans, Concentration, Wavelength scan with graphic, Kinetics with graphic, Multi-wavelength (7), Storage of methods, Storage of results, Windows environment

Type no.	Order no.	Product	Description
Uvi Light XS5	285600130	Spectrophotometer	Uvi Light XS single beam spectrophotometer 5 nm, with 10 – 100 mm standard cell holder
Uvi Light XS2	285600140	Spectrophotometer	Uvi Light XS single beam spectrophotometer 2 nm, with 10 – 100 mm standard cell holder
Uvi Light XTD5	285600170	Spectrophotometer	Uvi Light XTD single beam spectrophotometer 5 nm, PS/2-mouse, keyboard and 15" TFT screen, 10 – 100 mm cell holder
Uvi Light XTD2	285600180	Spectrophotometer	Uvi Light XTD single beam spectrophotometer 2 nm, PS/2-mouse, keyboard and 15" TFT screen, 10 – 100 mm cell holder
Uvi Light XT5	285600175	Spectrophotometer	Uvi Light XT single beam spectrophotometer 5 nm, with integrated LCD colour screen and 10 – 100 mm cell holder
Uvi Light XT2	285600185	Spectrophotometer	Uvi Light XT single beam spectrophotometer 2 nm, with integrated LCD colour screen and 10 – 100 mm cell holder
Uvi Light PC2	285600160	Spectrophotometer	Uvi Light PC single beam spectrophotometer 2 nm, with 10 – 100 mm standard cell holder, LabPower software and PC-connection cable
SZ 1010	285600300	Software	LabPower for Uvi Light and Uvikon (delivered with UVIKON XS) with PC-connection cable SZ 1430
SZ 1520	285600810	Cell holder	9 positions, 10 mm, for Uvi Light (requires cell changer SZ 1550), the Uvi Light XS requires software LabPower for control
SZ 1530	285600820	Cell holder	9 positions, 10 – 100 mm, for Uvi Light (requires cell changer SZ 1550), the Uvi Light XS requires software LabPower for control
SZ 1540	285600830	Cell holder	9 positions, thermostatable 10 mm, for Uvi Light (requires cell changer SZ 1550), Uvi Light XS requires software LabPower for control
SZ 1550	285600840	Automatic cell changer	For all 9 positions cell holder
SZ 1560	285600850	Sipper	Contents peltier thermostated single cell holder, 80 µl flow-through quartz cell, peristaltic pump and tubing kit
SZ 1580	285600870	Sipper	Contents 80 µl flow-through quartz cell, peristaltic pump and tubing kit
SZ 1910	285601200	Printer	KYOLINE b&w thermal printer 40 columns with connection cable
SZ 1915	285601205	Printer	Black & white KYOLINE thermal (for Uvi Light XT and XTD only)

UVIKON – A true double beam instrument based on long-term experience

Experience

Since more than 70 years SCHOTT Instruments provides innovative solutions for the laboratory analysis. Now SCHOTT Instruments introduces the fourth-generation UVIKON UV-VIS Spectrophotometers.

The UVIKON gives you confidence in your data. Its superior and highly efficient true double beam optical design is matched with the latest digital signal processing technology (ADSP, Advanced Digital Signal Processing). ADSP guarantees reliable data acquisition, outstanding sensitivity and extremely low signal noise for consistent linear measurement results – even for solutions with high ground absorption. These unique advantages give you the flexibility to accomplish the analyses you would like to process.

The right instrument

Two UVIKON models – the UVIKON XS and XL – are available to match your budget and analytical needs. Both models include full software for simple reads, standard curves, wavelength scans, wavelength programming, time drive and kinetics. Its software is built to make it quick and easy for all to use – regardless of the analytical method – without sacrificing capabilities.

The UVIKON XS is provided with the software LabPower. The UVIKON XL is provided with the software UVVISION. The UVVISION software is also available in a 21 CFR part 11 compliant version (please refer also to the technical specs of the software).

Make your next UV-VIS an instrument built to give you confidence in your data.

Let us prove value on how the UVIKON is the instrument of choice to meet your present and future UV-VIS needs.



*UVIKON XS/XL
Double beam UV-VIS Spectrophotometer*

Ease of operation

UVIKON XS

The UVIKON XS offers double beam capability at a very competitive price. It sets the standard for both value and performance. The UVIKON XS is quality-engineered with optimized simplicity to provide more productivity to every laboratory. Its advanced double-beam technology, high accuracy and precision, and quick usability assure top-quality results every time.

UVIKON XL

The UVIKON XL sets the standard for excellence in optics, software and performance. With its low stray light, superior photometric linearity, and high resolution, the UVIKON XL easily meets your lab's advanced UV-VIS requirements. It gives every user the power to accomplish both routine and demanding applications with ease.

The UVIKON XL's automation and data handling capabilities also provide maximum versatility and productivity for your laboratory.

The UVIKON's status display always shows the operating status of the system



The UVIKON's sample compartment is engineered to save you time. There is plenty of room to insert samples easily, and it lets you quickly exchange accessories whenever you want. And of course it is designed for easy maintenance.

Superior Optics

Superior precision and accuracy begin with a superior optical design. SCHOTT Instruments sets the standard for uncompromised optical performance with the UVIKON's state-of-the-art optical design.

Measurement stability

- Reproducibility with every analysis
- Accurate data for long kinetics
- Double-beam optics with dark current correction
- Thermally insulated optical bench
- High-sensitivity detector

Precise data

- Accurate peak wavelength values
- Superb resolution for fine molecular structure data
- State-of-the-art grating drive system
- Precise data uncompromised by noise

Measurements even with open sample compartment lid

- Easy use of extra-large accessories
- Convenient access for larger cuvettes
- Optimised design of quartz windows in sample compartment

Direct usability with microcuvettes


- Microcuvettes fit directly into the standard cell holder without adapters or beam masks to enable sensitive small-volume measurements.

High energy throughput

- Extended photometric range
- Accurate measurements even at very high concentrations
- Modern computer-optimized optics
- Compact optics design
- Proven, efficient ZeissTM monochromator

Consistent performance

- Consistent performance now and in the future.
- Optics protected from lab environment
- Reproducible and reliable signal processing
- True symmetrical optics and electronics
- Automatic wavelength calibration at switch on








Lamps compartment thermally isolated from optical bench

Deuterium und Halogen lamp

High-efficiency Zeiss® grating designed for maximum optical throughput and low stray light




Chopper system corrects for any change in lamp output and dark current 100 times per second (120 at 60Hz)




Multifunctional Accessories

SCHOTT Instruments offers a complete range of accessories for the UVIKON spectrophotometers. Designed for reliability and built to exact standards, each accessory is pre-aligned. All accessory components can be easily replaced in the sample shaft in its original operation position, not requiring any tools.



- Peristaltic sipper with sample return and thermostatable flowcell holders
- Temperature sensors
- Magnetic stirrer
- Long path rectangular cell holder (up to 100 mm)
- Long path cylindrical cell holder (up to 100 mm)
- Filter and solid sample holder



The UVIKON's QuickLock feature allows installation and exchange of accessories without tools in a matter of seconds.

Gilson auto sampler for laboratory automation





Thermopack – external Peltier thermostating system




Performance validation kit with filters and software. Every Uvikon comes fully tested and includes a factory Certificate of Performance. All test data is archived at the factory for full traceability



Automatic water-thermostated sample changer for simultaneous analysis of up to 10 enzyme reactions



Thermosystem – Software controlled, Peltier thermostated accessory for nucleic acid thermal melts



Software UV VISION

Simple compilation of methods and automation of complex calculations:

All instrument parameters used for measuring, also including the formulas for calculating the results, are defined in a method and stored in the database. The registered measuring data is stored with the corresponding method and moreover the method version. As such, it is always traceable to see which method was used for generating the measuring data. In the event that a method is altered, a new version of the method will be released automatically.

The integrated formula editor enables to accomplish even complex calculations by automizing them. Regardless whether calculations of multiple samples and/or the calculation of one result of various wavelengths – every calculation is fully automated by the formula editor.

Documentation and storage of measuring data:

Upon completion of a single measurement or a measurement sequence, an automatic report can be produced. The report records the curves, the original measurement data, the calculated results, the method with all measurement parameters, the current user, and the date and time of the measurement. Further, specific client data such as the company name, department and company logo can be included.

The measurement data, the corresponding results and methods used, as well as the report are stored in a save database. To guarantee further data safety, the database can be automatically stored by backup on two separate paths. This makes a loss of data practically impossible.


Data export:

All raw data, graphics and results can be directly exported to Excel and as such be made available to other MSOffice applications.


UVVISION Pharma:

The version UVVISION Pharma fully complies with the requirements of the FDA Regulation 21 CFR Part 11 in respect to:


- **User Management:**
Three separate user groups with varying rights. The administrator can limit the creation of methods to a small group, for example.
- **Electronic Records:**
All measurement data, reports, results and the corresponding methods are stored in the database. Double-assured data security.
- **Electronic Signature:**
The reports can be given a digital signature (login name and password).
- **Audit Trail:**
Every alteration of method or report will be automatically registered in the audit trail with date/time/time zone process, user name and a detailed description.




The UV VISION software supports all measurement methods in photometry such as:



- *Fixed Wavelength:*
Readings at a fixed wavelength, with calibration curve
- *Wavelength Program:*
Readings at numerous wavelengths



- *Wavelength Scan:*
Spectra recording



- *Time Drive:*
Sequential readings with rapid kinetics
- *Auto Rate Assay:*
Parallel measurements from up to 10 samples at slower assay of activity

Software Lab Power

The Lab Power software supports all basic routine applications in photometry such as:

Fixed Wavelength:

Measurements at a single wavelength including calibration curves up to 33 standards.

Wavelength Program:

Allows data acquisition at 2–14 wavelengths)

Wavelength Scan:

Adding of spectra with real time calculation and post-run calculations of already existing spectra. Maximum, minimum and derivative evaluations as well as curve overlay maximum are possible.

Time Drive:

Allows sequential measurements of kinetics data at very high speed (up to 2000 data points/min) or over a long period (15 hours).

Auto Rate Assay:

Provides parallel kinetics measurements up to 10 samples using the cell changer

Further evaluation functions:


- Enzyme activity calculations by linear regression
- Calibration function via linear or quadratic regression as well as linear interpolation
- Curve addition, subtraction, multiplication and division
- Wavelength ratio, difference and base point

Report Capability:

Automatic printing and storage of individual analytical reports. All measuring data are stored with specific comments and file extension.

Data Export:

Simple export into Excel or ASCII.



Specifications

Optical Design

- Well-proven true symmetric double beam optics with only 6 reflecting surfaces per beam.
- High energy diffraction Zeiss® monochromator design using a holographic blazed grating with 1300 lines per mm.
- Very sensitive photodiode (XS)/photo-multiplier (XL) with a rotating chopper providing 100% energy in sample and reference beams plus dark current correction.
- External source bench providing thermal insulation with optical bench and isolation of the whole photometric assembly against contamination by ozone emitted by deuterium lamp.
- Optical bench base manufactured to a tolerance of 0.01 mm in 1 meter.
- Complete optical system sealed against atmospheric contamination from dust and volatile compounds.

Sample Compartment

- 140 mm depth, 332 mm width, 162 mm height, with 120 mm beam separation.
- Total and easy access: front top and back covers are quickly removable.
- Three-point positioning of accessories for high reproducibility of results.
- Quartz windows at the beam entrance and exit to the sample compartment prevent the optics from being contaminated.

Accessories

- Thermostatable standard cell holder directly useable with micro and ultramicro cells for small samples.
- Large family of optional Quick Lock accessories including a variety of cell holders, stirrer, thermostatable automatic cell changer, temperature sensor, thermostatable sipper, Thermopack and Thermosystem (PeltierTechnology).

Ambient temperature

- Operating temperature: 15 to 30°C.

Configuration requirements

- Recommended: PC with Windows XP Pro, 1 free RS 232 port, 128 MB RAM, 50 MB free disk space.
- Optional printer: all standard Windows printers.

Technical Specification UVIKON XS

Technical Data UVIKON XS

Wavelength range	190 – 1100 nm
Wavelength steps	0.05 – 10 nm
Scan speed	5 – 2000 nm/min
Transfer speed	5000 nm/min
Spectral bandwidth	1.8 nm (fixed)
Wavelength accuracy (Holmium oxide filter at 10 nm/1 nm slit)	± 0.3 nm
Wavelength reproducibility (with Holmium oxide filter at 44.4 nm and second wavelength at 536.7 nm; positioning speed 7,000 nm/min)	± 0.03 nm
Baseline flatness, 200 – 800 nm (1 nm Step, 200 nm/min scanning speed))	< ±1 mAbs
Stray light (NaI (10g/l), 220 nm, 1.8 nm, 1 sec)	< 0.03 %T
Linearity (250 nm, 1.8 nm, 1 sec. $r^2 > 0.999$)	3.3 Abs
Photometric range	±3.5 Abs
Photometric accuracy (0.85 Abs, certified filters, 1 sec., 590 nm)	±3 mAbs
Photometric precision (0.85 Abs, certified filters, 1 sec., 590 nm)	± 0.5 mAbs
Noise (RMS), (measured at 580 nm for 5 min, 1 sec, ASTM E685)	< 30 µAbs
Drift (580 nm, 1 sec, after 2 hrs. of warm-up)	< 0.1 mAbs/h
Response time	0.02 – 5 sec
Light sources	Deuterium and Tungsten-Halogen
Monochromator	High energy, low stray light diffraction grating with 1300 lines/mm, 175 mm focal length
Dimensions (W x D x H)	680 mm x 565 mm x 275 mm
Weight	35 kg net
Mains Power	100/240 VAC ± 10%, 50/60 Hz
Power Consumption	Max. 200 VA

... UVIKON XL and Ordering Information

Technical Data UVIKON XL

Wavelength range	180 – 900 nm
Wavelength steps	0.05 – 10 nm
Scan speed	5 – 2000 nm/min
Transfer speed	7000 nm/min
Spectral bandwidth	0.2 – 0.5 – 1 – 2 – 4 – 6 nm; 0.5 – 1 – 2 – 4 nm, reduced height
Wavelength accuracy (Holmium oxide filter at 10 nm/1 nm slit)	± 0.25 nm
Wavelength reproducibility (with Holmium oxide filter at 44.4 nm and second wavelength at 536.7 nm; positioning speed 7.000 nm/min)	± 0.025 nm
Baseline flatness, 200 – 800 nm (1 nm Step, 200 nm/min scanning speed)	< ±1 mAbs
Stray light (NaI, 220 nm, 1 nm, 1 sec)	< 0.015 %T
Linearity (250 nm, 2 nm, 1 sec. $r^2 > 0.999$)	3.3 Abs
Photometric range	± 5 Abs
Photometric accuracy (0.85 Abs, certified filters, 1 sec., 590 nm)	± 3 mAbs
Photometric precision (0.85 Abs, certified filters, 1 sec., 590 nm)	± 0.5 mAbs
Noise (RMS), (measured at 580 nm for 5 min, 1 sec, ASTM E685)	40 µ Abs
Drift (580 nm, 1 sec, after 2 hrs. of warm-up)	< 0.1 mAbs/h
Response time	0.02 – 10 sec
Light sources	Deuterium and Tungsten-Halogen
Monochromator	High energy, low stray light diffraction grating with 1300 lines/mm, 175 mm focal length
Dimensions (W x D x H)	680 mm x 565 mm x 275 mm
Weight	35 kg net
Mains Power	100/240 VAC ± 10% 50/60 Hz
Power Consumption	Max. 200VA

Type no.	Order no.	Product
Uvikon XS (CE) (US/CSA)	285600110	Includes LabPower software package. The packages includes: LabPower software, DNA/RNA and validation applications
Uvikon XS-UVVision	285600100	Delivered with Software UV VISION
Uvikon XL (CE) (US/CSA)	285600120	Includes UVVISION Software Package (standard version) All UVIKON models include standard thermostatable cell holder compatible with micro and ultra micro cells, accessory controller to allow connection of all UVIKON accessories.
SZ 1015	285600305	UVVISION software
SZ 1016	285600306	UVVISION Pharma, software according to 21 CFR, part 11
SZ 1017	285600307	Upgrade to UVVISION Pharma
SZ 1035	285600325	Thermo system (6x6 Auto cuvette changer, Peltier thermostated)
SZ 1036	285600326	Thermo pack (external water bath with Peltier technology)
SZ 1110	285600400	6+6 automatic cell changer, water thermostatable
SZ 1140	285600430	Sipper system

Capillary viscometry from SCHOTT Instruments – know how from the very beginning

Innovative capillary viscometry – from the outset

The viscosity of Newtonian fluids can be most precisely determined using capillary viscometers. This method of measurement, measures the time taken for a defined quantity of fluid to flow through a capillary with a known diameter and known length. With the industrial production of such precisely calibrated capillary viscometers, we have created the conditions to enable this measuring method to establish itself worldwide as a reliable procedure.

With the development of the first automatic measuring systems, we replaced the stopwatch with automatic registration of the fluid at the start of the 1970's. Since then, subjective measuring errors have been a thing of the past.

Further developments and improvements of viscometers, measuring instruments and accessories led to a range of products whose excellent performance is universally recognized. It is therefore no wonder that viscosity measurement systems from SCHOTT Instruments have become indispensable production control and quality insurance tools worldwide, whether in the mineral oil industry, for polymer manufacturers and processors, in the pharmaceutical or food industry.

SCHOTT Instruments capillary viscometers are the worldwide basis for precise viscosity measurements of Newtonian fluids.



Contents viscometry

Viscometers within quality assurance systems	Page 170
Polymer applications for the AVS measurement systems	Page 171
ViscoClock	Page 172
ViscoSystem® AVS 370	Page 174
Order information ViscoSystem® AVS 370	Page 184
ViscoSystem® AVS 470	Page 180
Order information ViscoSystem® AVS 470	Page 185
ViscoSystem® AVS 270	Page 186
Automatic sampler AVSPro II	Page 188
Technical data AVSPro II	Page 192
Transparent thermostats	Page 194
Accessories	Page 197
Viscometers and their range of use	Page 198
Ubbelohde viscometers, normal form	Page 199
Ubbelohde viscometers, normal form (ASTM)	Page 200
Ubbelohde viscometers, with additional tube and threads	Page 201
Ubbelohde viscometers with TC sensors	Page 202
Micro-Ubbelohde viscometers, viscometers for dilution series	Page 205
Cannon-Fenske viscometers	Page 206
Ostwald viscometers	Page 208
Accessories	Page 209
AVS measuring stands and tube sets	Page 214
AVS measuring stands	Page 215









































































Viscometers within quality assurance systems

Business sector	Product	Example
Automotive engineering	motor oil (fresh and used) high polymer plastics	bumpers
Brewery	original wort hop-wort	beer beer
Electrical engineering and electronics	high polymer plastics of all types	chips, casings
Power supply	turbine oil transformer oil	generators
Film	gelatine as pigment-bearing agent carrier film for film material	color films
Plastics manufacturers	high polymer plastics of all types	
Plastics processors	high polymer plastics of all types	injection molding
Food industry	starch gelatine packaging materials milk products fruit and fruit juice concentrates gelatinizing agents	instant flour thickeners jelly bears yoghurt containers yoghurt drink pectin
Aviation	high polymer plastics of all types fuels hydraulic fluids	kerosene horizontal stabilizers and undercarriages
Mechanical engineering	motor oil hardening emulsions hydraulic fluids	mill trains stamp shops
Medicine	body fluids injection solutions tinctures and drops blood substitute materials	blood, bile insulin nose, eyes blood plasma
Mineral oil	light motor oil turbine oil liquid fuels of all types	gasoline, diesel fuel
Textile	high polymer plastics of all types cotton	for mixed fibers
Entertainment	high polymer plastics	CDs, videotapes

The table on the right illustrates the extensive area of high polymer plastics and the large variety of testing methods.

Polymer applications for the AVS measurement systems

Polymers, their applications and utilization of automatic systems from SCHOTT Instruments GmbH

Type	Abbr.	Solvent	Capillary	Operating temperature	Standards	Suitability of the AVS measurement systems			
						VC*	370	470	Pro
Cellulose EWNN	C I	Cuen Couxam	0c CAN CM I Micro	20 °C	SNV 195 598S 15:88				
Cellulose acetate	CA	Dimethyl-chloride/ methanol	0c I I Micro	25 °C	DIN 53 728/1				
Polyamide	PA	Sulphuric acid (96%)	II IIc	25 °C	DIN 53 727 ISO 307				
Polyamide	PA	Formic acid (90%)	I Ic	25 °C	DIN 53 727 ISO 307				
Polyamide	PA	m-cresol	II IIc	25 °C	DIN 53 727 ISO 307				
Polybutylene terephthalate	PBT	Phenol/dichloro benzene (50:50)	Ic II	25 °C	DIN 53 728/3 ISO 1628-4				
Polycarbonate	PC	Dichloromethane	0c I	25 °C	DIN 7744/2 ISO 1628-4				
Polyethylene	PE	Decahydro-naphthalene	I Ic	135 °C	DIN 53 728/5 ISO 1191ASTM D 1601				
Polyethylene terephthalate	PET	m-cresol	II IIc IIc Micro	25 °C	DIN 53 728/3 ISO 1628-5ASTM D 4603				
Polyethylene terephthalate	PET	Phenol/dichloro benzene (50:50)	Ic II	25 °C	DIN 53 728/3 ISO 1628-5ASTM D 4603				
Polyethylene terephthalate	PET	Dichloroacetic acid	II IIc Micro	25 °C					
Polymethyl methacrylate	PMMA	Chloroform	0c I	25 °C	DIN 7745/2 ISO 1628-6				
Polymethyl methacrylate	PMMA	Acetophenone	0c I	25 °C	DIN 7745/2 ISO 1628-6				
Polypropylene	PP	Decahydro-naphthalene	IIc	135 °C	DIN 53 728/4 ISO 1191				
Polyphenyl sulphide	PPS	Ortho dichloro naphthalene	IIc	230 °C					
Polystyrene	PS	Toluene	I Ic	25 °C					
Polysulphone	PSU	Chloroform	IIc	25 °C					
Polyvinyl chloride	PVC	Cyclohexanone	IIc	25 °C	DIN 53 726 ISO 1628-2ASTM D 1243				
Styrene-acrylonitrile copolymer	SAN	Ethyl methyl ketone	0c I	25 °C					
Styrene-butadiene copolymer	SB	Toluene	0c I	25 °C					

VC* = ViscoClock

 excellent suitability;  can be used;  limited suitability for application related reasons.

This table makes no claim to completeness.

ViscoClock.

If you need more accuracy:

The ViscoClock is the economically priced introductory model in the field of automatic viscosity measurements. Manual measurements with a stopwatch and a trained eye is therefore something of the past because time is money.

The ViscoClock

The ViscoClock is an electronic time-measuring unit used to determine absolute and relative viscosity. It consists of a stand which is used to mount a viscometer and the electronic measuring unit. The two measuring levels are integrated in the stand made of high-quality PPA synthetic material, and the electronic measuring unit is included in a PP casing. The large LCD display allows the measured values to be read off easily.

Range of use

The ViscoClock is designed for the use of an Ubbelohde viscometer, a Micro-Ubbelohde viscometer or a Micro-Ostwald viscometer made by SCHOTT Instruments. The ViscoClock automatically measures the flow-through time of temperature-stabilized liquids through the capillaries of the viscometer at temperatures ranging from $-40\text{ }^{\circ}\text{C}$ to $150\text{ }^{\circ}\text{C}$.



For temperature stabilization in the thermostatic bath, the following tempering liquids are suitable: water, alcohol water (e.g. ethanol, methanol), paraffin oil, and silicone oil. Liquids can be measured that qualify for use with the viscometer being used in each instance.

Accuracy

The most precise method used to determine the viscosity of liquids is their measurement in capillary viscometers; the ViscoClock functions according to this method. The operating time is indicated with a resolution of $1/100\text{ sec.}$ with quartz precision. The accuracy of $0.01\text{ }\%$ of the measured time used to calculate the absolute and relative viscosity is indicated as measuring uncertainty with a confidence level of $95\text{ }\%$.

Absolute viscosity

Only the calibrated viscometers made by SCHOTT Instruments are suitable for the calculation of absolute viscosity in the temperature-stabilized, transparent thermostatic bath.

Relative viscosity

For the measurement and calculation of relative viscosity, all Ubbelohde viscometers, uncalibrated and calibrated, can be used for manual or automatic measurements.

The ViscoClock can be used with any SCHOTT Instruments thermostats bath. The viscometer stand is included.

Technical data ViscoClock

Measuring range - time	up to 999.99 s; resolution 0.01 s
Accuracy of time measurement	$\pm 0.01 \text{ s} / \pm 1 \text{ digit}$; however no more precise than 0.1 %; indicated as measuring uncertainty with a confidence level of 95 %
Measuring range - viscosity	0.35 ... 10,000 mm ² /s (cSt) the absolute, kinematic viscosity is additionally dependent on the uncertainty of the numerical value of the viscometer constant and on the measuring conditions, in particular the measuring temperature
Display	5-digit LCD display, 20 x 48 mm (H x W), digit height 12.7 mm, seconds indication with 2 decimal digits after the decimal point, resolution 0.01 s
Voltage supply	low voltage U: 9 V
Plug-in connection	socket for low voltage connection: jack plug, internal contact $\varnothing = 2.1 \text{ mm}$, plus pole at pin contact, for connection of a TZ 1848 or TZ 1859 power supply unit
Power supply	in accordance with class of protection III. degree of protection for dust and humidity IP 50 in accordance with DIN 40 050 power supply unit 230 V, 50-60 Hz (TZ 1848) power supply unit 115 V, 50-60 Hz (TZ 1859), with US-plug not suitable for use in areas subject to explosion hazards
RS-232-C interface	for connection of a printer with serial interface or of a computer (PC) for documentation of the data
Plug-in connections	4 pole circular plug, mini, DIN
Configuration of RS-232-C interface, permanently set	4800 baud, 7 bit word length, 2 stop bits, no parity; after each measurement, the measured value is transmitted automatically. the string of digits consists of 4 digits before the decimal point, 2 digits after the decimal point, and the terminating characters CR and LF.
Ambient temperature	+10 ... +40 °C for storage and transport
Operating temperature	stand: -40 ... +150 °C electronic measuring unit: +10 ... +40 °C
Air moisture	in accordance with EN 61 010, Part 1; max. relative humidity 80 % for temperatures up to 31 °C, decreasing linearly to 50 % of relative humidity at a temperature of 40 °C
Materials	stand: polyphthalamide (PPA) casing*: polypropylene (PP) sealing membrane: silicone
Dimensions	approx. 490 x 95 x 50 mm (H x W x D)
Weight	approx. 450 g (without viscometer) power supply unit: approx. 220 g
Country of origin	Federal Republic of Germany
CE symbol	in accordance with Guideline 89/336/EEG (electromagnetic compatibility EMC): emitted interference in accordance with Standard EN 50 081, Part 1 interference immunity in accordance with Standard EN 50 082, Part 2, in accordance with Guideline 93/23/EEG (low voltage guideline), last altered by Guideline 93/68/EEG: Testing basis EN 61 010, Part 1
Viscometer types	Ubbelohde (DIN; ISO; ASTM; Micro), Micro-Ostwald
Transparent thermostatic baths	the ViscoClock can be used in every transparent thermostatic bath made by SCHOTT Instruments.

* Use in heat carrier liquids can result in discoloration of the synthetic material. The discoloration does not, however, have any effect on the function and quality of the ViscoClock.
DURAN® is a registered trademark of SCHOTT Glaswerke Mainz, Germany. Subject to technical changes.

ViscoSystem® AVS 370 makes maximum precision ...

Well equipped for every viscosity determination

With the ViscoSystem® AVS 370 we have created a measuring device, which not only measures as precisely and consistently as you expect from SCHOTT Instruments, but also offers you maximum flexibility and possibilities for future extensions. Furthermore, it also saves valuable space on the laboratory bench.

Now possible for the first time ever: "suction" and "pressure" measurement – with one device

The ViscoSystem® AVS 370 is the first viscosity measuring device, which can be used for both "suction" and "pressure" measurement. This enables simple adjustment of the method of measurement to each sample. This significantly reduces investment costs for measuring stations at which pressure and suction methods are to be used. In most cases, using the ViscoSystem® AVS 370 also achieves noticeable savings in setting up time.

Easy with a modular concept for future expansion

The ViscoSystem® AVS 370 has a modular design. The basic version is available with one ViscoPump II module in optical or in TC version. Up to 3 other ViscoPump II modules can be inserted in the compact 19" housing. This means a

measuring station can be adapted to increasing requirements at any time. The modular concept also significantly reduces the space required, and measuring instrument set-ups can be more easily and clearly arranged, for example for parallel and comparison measurements.



... easier and more flexible, with provision for future extension!

Can be extended from an affordable single measuring station up to an 8-sample station

The basic version of the ViscoSystem® AVS 370 is an affordable starter model, which can be used to measure high or low viscosity liquids. In the version for

TC viscometers, for example, it is ideal for measuring opaque and black fluids. If necessary, each single measuring station can be extended to form a multiple measuring station with PC-controlled-multitasking. The WinVisco 370 software included in the standard equipment

enables parallel operation of two fully equipped AVS 370, with a total of eight ViscoPump II modules. Each module can measure a different sample using a different method. All the results can be quickly and easily evaluated and documented independently of each other. It could hardly be more flexible!

Compatible with existing accessories

Existing accessories (thermostats, stands, flow through cooler, etc.) can continue to be used with the ViscoSystem® AVS 370. Also, virtually all customary capillary viscometers can be used.



The ViscoSystem® AVS 370 from SCHOTT Instruments. Up to 4 ViscoPump II modules can be integrated in the compact 19" housing. With a PC and the WinVisco 370 software, all kinds of different samples can be measured, evaluated and documented in parallel, independently of each other.

ViscoSystem® AVS 370 – the right solution for all situations

Anyone working with the ViscoSystem® AVS 370 is perfectly equipped for all tasks involved in determining viscosity using capillary viscometers.

How to automatically achieve the right results

PC-controlled, the ViscoSystem® AVS 370 determines the time which the liquid to be examined requires to flow through the measuring distance in the capillary viscometer with quartz precision. The time is displayed with a resolution of 0.01 s (1 digit).

Measurement of the flow time of the liquid's meniscus can be scanned optoelectronically or with TC sensors. (During optoelectronic scanning the meniscus is detected by glass light fibres, with TC sensors the sensor detects the different thermal conductivity of the sample and air.) Therefore the ViscoSystem® AVS 370 offers an extraordinary broad field of uses, which range from viscosity measurement of clear fluids through to black or fully opaque liquids.

New: Two working principles with the same device.

For the first time ever, with the ViscoSystem® AVS 370 you can use the same device to work with "pressure" or "suction". This gives you more flexibility and better adjustment to the liquids to be examined.

In the "pressure" method of working an overpressure of up to 0.1 bar is applied to the liquid in the capillary, this is particularly advantageous for fluids with a low boiling point. In the "suction" principle the sample is sucked up into the capillary by a vacuum. A greater reproducibility of results is achieved using the "suction" method for higher viscosity samples. A further advantage of the "suction" principle is that it guarantees formation of the "hanging ball level" in Ubbelohde viscometers, even for the minimum required sample quantities.

Working with the ViscoSystem® AVS 370 is easy

The ViscoSystem® AVS 370 is very easy to use. The whole measuring procedure takes place automatically, subjective measuring errors are reliably precluded. The PC starts the measurement. After the set pretempering period has expired the entered number of measurements are carried out and the measured values saved.

The system can be protected against accidental overpumping or oversuction by means of a capacitive sensor. This prevents the sample to be measured from getting into the vessel containing the tempering liquid or inside the device.

"Suction" or "pressure"? A comparison of preferred applications

	"pressure"	"suction"
Highly viscous samples e.g. oils, polymers	■	■
Solvent (examples):		
highly volatile	■	–
Dichlormethane	■	–
Chloroform	■	–
Sulfuric acid	–	■
Dichloroethanoic acid	–	■
Toluene	■	■
Hexafluorisopropanol	–	■
m-cresol	–	■
Formic acid	–	■
Phenol-dichlorobenzene	–	■
Phenol-tetrachloroethane	–	■



Technical data

Unique flexibility

In the PC-controlled multiple measuring station, the ViscoSystem® AVS 370 offers you unique flexibility while working in a very small space: Up to eight modules, which equates to two fully equipped ViscoSystem® AVS 370, can be run parallel with the WinVisco 370 software. Each module can measure the same or different samples using "pressure" or "suction", fully independently of each other. In this way, series of measurements can be prepared extremely quickly and immediately evaluated and documented in the computer. This significantly reduces the time required to carry out viscosity measurements, especially for in process controls and quality assurance.



ViscoSystem® AVS 370	
Measuring range (time)	up to 9,999.99 s; resolution 0.01 s
Measuring range (viscosity)	pressure: 0.35 ... 1,800 mm ² /s (cSt) suction: 0.35 ... approx. 5,000 mm ² /s (cSt)
Measured parameter	flow through time [s]
Accuracy of the time measurement	± 0.01 %
Measured value display	via PC
Display accuracy	± 1 digit (0.01 s)
Pump pressure	automatically controlled
Preselectable tempering period	0 ... 20 min
Preselectable number of measurements	up to 10
Connections	
Pneumatic connections	threaded connections for viscometers
Electrical connections	circular connector with bayonet lock for measuring stands and TC viscometers
RS-232-C interface	9-pin
Mains connections	plug in accordance with EN 60320
Pump connection	socket outlet in accordance with EN 60320
Data Input/Output	serial to EIA RS-232-C
Ambient conditions	
Ambient temperature	+10 ... +40 °C
Air humidity	max. 85 % rel.
Housing	
Material	coated aluminum plate
Dimensions (for 1 ... 4 modules)	(W x H x D) approx. 255 x 205 x 320 mm
Weight (incl. 1 module)	approx. 5.4 kg
Power supply	90 ... 240 V ~, 50 ... 60 Hz
Equipment safety	EMC-Compatibility according to the Directive 89/336/EEC of the Council; low-voltage directive according to the Directive 73/23/EEC of the Council, as amended by the Directive 93/68/EEC of the Council
Multi-tasking	for 1 ... 8 ViscoPump II modules, with WinVisco 370 software

The following viscometers can be used with the ViscoSystem® AVS 370: Ubbelohde viscometer to DIN, Ubbelohde viscometer to ASTM, micro Ubbelohde viscometer to DIN, micro Ostwald viscometer, Cannon-Fenske routine viscometer, TC-Ubbelohde viscometer, TC-micro Ubbelohde viscometer.

We reserve the right to make technical changes.
ViscoSystem® is a registered trademark.

Compact, space-saving viscosity measuring station with the ViscoSystem® AVS 370. The measuring device is attached to a support table (recommended accessory). All the connections are within view and easily controlled. If necessary the sample can be automatically sucked away and the viscometer flushed e.g. with the TITRONIC® universal or the TITRONIC® 110 plus burette.

Real multitasking for up to 8 measurements in parallel mode ...


Easy to understand, proven in practice: The WinVisco 370 software

WinVisco 370 is the ideal software for the ViscoSystem® AVS 370*). It is supplied as part of the standard equipment. WinVisco 370 is easy to understand and can be quickly learned. Up to eight viscosity measurement modules can be controlled with only a few operating steps. The device parameters are easy to enter: Constants, t_0 flow time, number of measurements, pretempering period, type of viscometer, date and sample labeling for each measuring station.



WinVisco 370 works in real multitasking mode. This makes it possible for each measurement to be processed independently from the others. It also means that time-consuming measurements can be carried out from the same PC, without hindering the course of other, faster measurements. During the measurements you can change the monitor displays, start or stop other measurements, print out or save measured values. All data provided by the software can be passed on to an LIMS system.

WinVisco 370 supports three groups of users. For simple use, access is limited to: select viscometer, measure, load and save methods as well as enter parameters. In the highest level, users with administrator status can access all the software's facilities. Each user is given a user ID, an access level and a password.


*) The language (English or German) can be chosen after installation over the programme menu.



All the important parameters required for the measurement are displayed on the "Methods/Results" page. If necessary, the parameter editor can be called up using "Add Parameter", in order to enter non-standard or customer specific formulae.




All the measurements currently running can be monitored in parallel in the overview.




The viscometer data required for the evaluation can be stored in a table. This guarantees perfect allocation of e.g. the t_0 runtime, viscometer constants, the series number, etc. for each individual viscometer being used.


... with the practically proven WinVisco 370 software



The password protection prevents unwanted or confusing changes to the important measurement parameters.



The parameters can be individually adjusted to the measurement for each measuring position.




Each rinsing/dry step can be individually preselected. Even the application dependent quantity of solvent and the drying time can be separately determined.

With the ViscoSystem® AVS 370 and WinVisco 370 you will even quickly find the right connection for rinsing

With the daisy chain link of the ViscoSystem® AVS 370, further devices can be integrated in the system and controlled using the WinVisco 370 software. For example, when working in suction mode the viscometers can be rinsed using the TITRONIC® universal and TITRONIC® 110 plus burettes. The TITRONIC® universal is preferably used for light solvents, the TITRONIC® 110 plus for solvents with a viscosity $>3 \text{ mm}^2/\text{s}$. A special interchangeable unit (TA 50V) is available for highly aggressive solvents.

A vacuum pump (accessory) integrated in the system is used to conveniently suck away samples and solvents.



Transparent thermostat CT 54 with four measuring stands

Two basic concepts are available for the rinsing:

- A ViscoSystem® AVS 370 with four ViscoPump II modules (four measuring positions) and eight burettes, which enable each viscometer to be rinsed with two solvents. Time-consuming removal of the transparent thermostat for external rinsing of the viscometer is no longer necessary.
- Two ViscoSystem® AVS 370 complete with four ViscoPump II modules each (eight measuring positions), which enables semi-automatic rinsing of the viscometer with the next sample or solvent.

Precise Capillary Viscometry – Easy, Flexible and Independent: ViscoSystem® AVS 470



That's New: "Suction" and "Pressure" Measurements With Just One Instrument, No need for a PC

The ViscoSystem® AVS 470 is the first viscosity measuring device that allows "suction" and "pressure" measurements completely independent of a PC. This makes for maximum independence and

flexibility, allowing you to set up a measuring station that meets highest requirements even under difficult conditions, e.g. to monitor production or control quality in the polymers and mineral oil industry.

Perfectly Equipped For Fully Automatic Viscosity Measurements

The ViscoSystem® AVS 470 is a measuring system that includes almost everything you need to take precise and reproducible measurements. All common types of viscosity calculation are already integrated into the device, a small PS2 keyboard is all you need to



The ViscoSystem® AVS 470 needs no PC and therefore requires just a little more space than a sheet of paper.

Keyboard and printer are available as accessories.

enter additional data. A serial printer can be used to conveniently document your measuring results.

So, in a minimum of space, you can set up a measuring station equal in every way to complex measuring installations in terms of precision and reproducibility.

Simple and updateable Modular Concept

The ViscoSystem® AVS 470 is of a modular design and an optional optical or TC version ViscoPump II module can be used to adapt your measuring station to new requirements at any time. It is possible to connect the new

cleaning system AVS 270. You can use your existing accessories such as thermostats, stands, flow-through coolers or automatic cleaners e.g. AVS 26 with the ViscoSystem® AVS 470. Also, virtually all customary capillary viscometers can be used.


Precise and Reliable – The ViscoSystem® AVS 470

Working With the ViscoSystem® AVS 470 Is Easy

The ViscoSystem® AVS 470 is very easy to handle. The desired measuring method can be preselected and started on the device. The entire measurement is taken automatically to rule out subjective measurement errors. Once the set pre-heating time is reached, the desired number of measurements are taken and the viscometer automatically cleaned if required. The status of the measurements is continuously indicated on the LC.

If required, individual parameters may be input via a PS 2 key board (in scope of delivery). A serial printer can be used to print measurement logs.

The connections are on the front panel of the device for easy control. Over-pumping and oversuction are prevented by means of a capacitive sensor (optional).




Right figure:
The print-out shows everything you need for reliable documentation of your test. Simultaneously it demonstrates the unique performance of the ViscoSystem® AVS 470.

“Suction” or “Pressure”? Preferred applications in comparison

		“Pressure”	“Suction”
highly viscous samples e.g. oils, polymers		■	■
Solvents (examples):	highly volatile	■	–
	Dichloromethane	■	–
	Chloroform	■	–
	Sulfuric acid	–	■
	Dichloroacetic acid	–	■
	Toluene	■	■
	Hexafluoro-isopropanol	–	■
	m-cresol	–	■
	Formic acid	–	■
	Phenol-dichlorobenzene	–	■
	Phenol-Tetrachloroethane	–	■


Technical data

Clear user guidance, clear status – even without PC!




AVS 470

After switching on the AVS a self test is run and then an entry prompt appears.



AVS 470

The parameters can be set in the test mode. The t_0 value is determined automatically.



AVS 470

All setup parameters can be preset conveniently, e.g. pressure/suction, velocity, waiting time between two tests, etc.



AVS 470

The readings can be read off conveniently on the display regardless of whether or not a printer is connected.

ViscoSystem® AVS 470	
Measuring range (time)	up to 9,999.99 s; resolution 0.01 s
Measuring range (viscosity)	pressure: 0.35 ... 1,800 mm ² /s (cSt) suction: 0.35 ... approx. 5,000 mm ² /s (cSt)
Measured parameter	flow-through time [s]
Time measuring accuracy	± 0.01 %
Measured value display	LC-Display
Display accuracy	± 0.01 s, ± 1 Digit, but not exceeding 0.01%
Pumping pressure	fully automatically controlled suction up to approx. -160 mbar, pressure up to approx. +160 mbar
Preselectable tempering period	0 ... 20 min
Preselectable no. of measurements	1 to 99 for each sample
Connections	
Pneumatic connections	threaded connections for viscometers
Electrical connections	circular connector with bayonet lock for viscometer 4-pin DIN socket for TC viscometer 4-pin circular connector for capacitive sensor 7-pin circular connector for AVS 26, with bayonet lock
RS-232-C interface	9-pin for serial printer
Mains connection	connector in acc. with EN 60320
Pump connection	socket outlet in accordance with EN 60320
Ambient Conditions	
Ambient temperature	+10 ... +40 °C for operation and storage
Air humidity	max. 80 % in acc. with EN 61010, Part 1
Housing	
Material	steel aluminium housing; with chemically resistant 2-component coating
Dimensions	(W x H x D) Approx. 255 x 205 x 320 mm
Weight (incl. pump module)	approx. 5.4 kg
Power supply	
	90 ... 240 V ~, 50 ... 60 Hz
Equipment safety	
	EMC in acc. with Council Directive 89/336/EWG; low-voltage directive

The following viscometers can be used with the ViscoSystem® AVS 470:

Ubbelohde viscometer to DIN, Ubbelohde viscometer to ASTM, micro Ubbelohde viscometer to DIN, micro Ostwald viscometer, Cannon-Fenske routine viscometer, TC Ubbelohde viscometer, TC micro Ubbelohde viscometer.

We reserve the right to make technical changes.

ViscoSystem® is a registered trademark of SCHOTT Instruments.

Ordering information ViscoSystem® AVS 370



ViscoSystem® AVS 370

The AVS 370 viscosity measuring system is made up of individual components, which have to be individually ordered. Please always ask for a detailed offer.

Description	Order No.
ViscoSystem® AVS 370 basic unit, housing incl. one ViscoPump II module and WinVisco 370 software, for optoelectronic scanning	1056509
ViscoSystem® AVS 370 basic unit, housing incl. one ViscoPump II module and WinVisco 370 software, for TC scanning	1056515
ViscoPump II module for optical scanning, VZ 8511	1054306
ViscoPump II module for TC scanning, VZ 8512	1054304
Accessories	
Support table	1057903
Vacuum pump (230 V)	1057901
Vacuum pump (115 V)	1057902
Measuring stand AVS/S	28 541 0502
Measuring stand AVS/SK	28 541 0876
Measuring stand AVS/SK-CF	28 541 0892
Measuring stand AVS/SK-V	28 541 0905
Fixing frame	28 540 5043
Holder VZ 7191	28 542 1968

Ordering information ViscoSystem® AVS 470



ViscoSystem® AVS 470

The AVS 470 viscosity test station is composed of individual components, which must also be ordered separately. Please always request a detailed offer.

Description	Order No.
ViscoSystem® AVS 470 basic unit, housing incl. one ViscoPump II module for opto-electronic sensing, Version: 95 V to 230 V/50-60 Hz	28 541 5709
ViscoSystem® AVS 470 basic unit, housing incl. one ViscoPump II module for TC sensing, Version: 95 V to 230 V/50-60 Hz	28 541 5708
ViscoPump II module for optical sensing, VZ 8511	1054306
ViscoPump II module for TC sensing, VZ 8512	1054304
Accessories	
Support table	1057903
RS-232-C Data printer (230 V), TZ 3460 R	28 522 5608
Vacuum pump (230 V)	1057901
Vacuum pump (115 V)	1057902
Measuring stand AVS/S	28 541 0502
Measuring stand AVS/SK	28 541 0876
Measuring stand AVS/SK-CF	28 541 0892
Measuring stand AVS/SK-V	28 541 0905
Fixing frame	28 540 5043
Holder VZ 7191	28 542 1968

ViscoSystem® AVS 270 – the automatic cleaning system

The automatic rinsing system AVS 270 replaces the proven rinsing automats AVS 24/26.

The used viscometer can be rinsed immediately after completion of the series of measurements or on demand. This has the advantage that the measured values can be examined first and the rinsing process started afterwards.

The modular construction of the AVS 270 allows the fast and easy exchange of system components.

The AVS 270 allows the usage of both detection systems (optical and thermal).



Technical data AVS 270

Country of origin	Germany
Display	7-segment LED for status message
Mains connection	Connector in accordance with EN 60320
Pump connection	Socket outlet in accordance with EN 60320
Power supply	90 ... 240 V~, 50 ... 60 Hz complies with protection class 1 acc. to DIN 57 411, part 1 / VDE 0411, part 1
Protection class	IP 20 in accordance with EN 60 529
Power consumption	At 230 V [VA]: 100; at 115 V [VA]: 200
Data transmission	Bidirectional serial interface in accordance with EIA RS-232-C Baud rate: 4.800 Parity: no Word length: 7 Bit Stop bit: 2
Pneumatic connections	Threaded connections for viscometers and solvent
Dimensions	205 x 255 x 350 mm (H x W x D)
Weight	6 kg
Housing Material	Steel aluminium housing with chemically resistant 2-component coating, stackable
Keypad	Multi-color polyester foil
Ambient conditions	Ambient temperature: +10 ... 40°C (for operation und storage) Air humidity in acc. with DIN EN 61 010, part 1 Max. relative humidity 80 % for temperatures up to 31°C, declining linearly to 50 % relative humidity at a temperature of 40°C
Equipment safety	EMC in acc. with Council Directive 89/336/EWG; Low-voltage directive in acc. with Council Directive 73/23/EWG; as amended by Council Directive 93/68/EWG Examination based on EN 61 010

Ordering information ViscoSystem® AVS 270

The following viscometers can be employed:

- DIN Ubbelohde
- ASTM Ubbelohde
- TC Ubbelohde without rinsing tube
- Micro TC Ubbelohde without rinsing tube
- DIN Micro Ubbelohde
- Cannon Fenske routine
- Micro Ostwald

Safety:

The vacuum operation guarantees optimum safety. The materials in contact with the permitted solvents are according to the highest requirements of a modern lab. The screwing caps for instance are made from PTFE, PCTFE, PP and PPS.

Parameter:

No PC is required when an AVS 470 is used. The rinsing parameters can be set easily and comfortably via the PS2 keyboard of the AVS 470 or via the foil keypad of the AVS 270.

Flexibility:

Standard-Viscometer according to DIN, ASTM and ISO can be employed. No viscometers of special forms (i.e. with threads or special connections) are required.

ViscoSystem® AVS 270

The automatic rinsing system AVS 270 is composed of individual components, which must also be ordered separately. Please always request a detailed offer.

Description	Order No.
ViscoSystem® AVS 270 basic unit for connection to viscosity measuring device	28 541 5697
ViscoSystem® AVS 470 (Scope of delivery: All necessary connection tubes/cables and solvent bottles)	
Accessories	
Vacuum pump VZ 8630, 230V/50Hz	105 7901
Vacuum pump VZ 8631, 115V/60Hz	105 7902

Selecting the solvents:

Proven solvents are, e.g., alcohols, ketons, esters, chlorinated hydrocarbons, aliphatic and aromatic hydrocarbons. Make sure the boiling point of the selected solvent is appropriate to the temperature of the thermostatic bath. The first solvent should have the maximum possible solving power (detergence) and the boiling point should be at least 50° above the bath temperature.

The second solvent is used mainly to eliminate the first, high boiling point solvent and its temperature should be slightly higher than that of the thermostatic bath.

Exclusively organic solvents or aqueous solvents are used with non-corrosive characteristics in respect to glass, brass and stainless steel. Only such solvents

without solid components must be used. Otherwise there is a risk of contaminating the viscometer capillaries and valves.

When the according purity of the solvent cannot be guaranteed, it is recommended to use filters. A solvent that is now used will fulfil both requirements. It must be considered that the solvent temperature must not exceed the 80 °C.

Due to safety reasons it is not allowed to use any solvents with an ignition temperature of under 250°C. Regarding the solvent storage or disposal vessels it is authorized to use the standard laboratory bottles or customized vessels.

Automatic viscosity measurement has been improved ...

The AVSPro II automatic sampler is a fully automated measuring instrument for determining the viscosity of Newtonian fluids with capillary viscometers. In spite of the high sample throughput, the AVSPro II provides maximum accuracy and reproducibility. Furthermore, working with the AVSPro II is easy and even allows unsupervised 24-hour operation.

Particularly with time consuming measurement runs, the AVSPro II helps to substantially reduce the burden on qualified employees. An additional advantage is the increased level of safety when handling aggressive media, e.g. sulphuric acid, which is achieved through the fully automatic measurement procedure.

The ProClean system and the micro dosing make routine operation safer. The filtration of solutions, which occasionally may be harmful, can thereby be omitted.

The capacitive sensors in the suction pipe effectively prevent any damage of the measuring system.

Because of its large throughput capacity and its functional reliability, which it has demonstrated in the course of continuous practical operations, the AVSPro II has proven itself to be an indispensable instrument for day-to-day utilization, particularly in the petroleum and plastics industries.



The AVSPro II automatic sampler works with the capillary method, which is the most precise method for determining the viscosity of Newtonian liquids in terms of physical chemistry. Using this method, measurements with an accuracy of more than 0.1% can be achieved. The great versatility offered by viscometers with optical and TC sensing systems opens up an extremely wide range of applications. This includes measurements of clear liquids as well as opaque petroleum products.

The viscosity measurement requirements of the polymer and petroleum industries in particular have been incorporated into the design of the AVSPro II. The main feature of the automatic unit is the three-axis positioning mechanism of the sample dosing system. The X-Y-Z positioning mechanism allows operation of up to four Micro TC viscometers in two thermostatic baths, which can be set at two different measurement temperatures. This method is used in the oil industry in order to determine the viscosity index.

... with the AVSPro II Automatic Sampler:



Two different sample racks are supplied:
a) one rack with 56 positions for 20 ml sample bottles for micro-viscometer applications



b) one rack with 16 positions to accommodate 100 ml sample bottles for normal volume applications



The electric sample lift ensures positioning of the samples in the rack at a convenient and easily monitored working height.


The AVSPro II allows the operator to select optionally the sample sequence and which sample is to be filled into which viscometer. The dosing system is available in either normal or micro construction and operates without a valve. It is thus suitable for nearly any type of sample.

The AVSPro II is equipped with opto-electronic and TC scanning (TC = thermal conductivity method) functions for the meniscus passage in the capillary viscometer. The samples are positioned in the sample rack, which is easy to load using the electric motorized lifting mechanism. If needed, the rack can be temperature-regulated.

The operator interface and control logistics are spatially and logically separate. This ensures a high degree of flexibility with regard to the installation location, and serves to reduce the environmental influences on the measurement results.

AVSPro II

Working with the AVSPro II is ...



Operating the AVSPro II is extremely easy. The operator controls the process at a PC connected via the RS-232-C interface. The intuitive user interface of the operation software guides the user clearly through the program. All data inputs are made using the computer keyboard and mouse.

A faulty operating status is indicated by acoustic or optical signals such as arrows, icons and other status messages or request messages. During the entire work sequence, the respective status of the AVSPro II is documented on the computer screen. Furthermore, status indicators can be selected for each

individual measuring position, which provide additional information on the operating status.

For the respective type of measurement, pre-parameterized sets of parameters depending on the viscometers, temperature and other measurement criteria are already provided. In addition, all parameters can be individually adjusted to special requirements at a special menu level. All of the standard calculation methods are available.

Screenshot: 16 sample rack

The AVSPro II allows individual allocation between the characteristics of the sample and the viscometers that are currently in operation.

In practice, this means that it is not only possible to simultaneously test the characteristics of samples with greatly differing viscosity, but also to perform measurements in various different capillary sizes and types of viscometers. This even applies to a combination of optical and thermal scanning. Therefore, preliminary sorting of the samples with regard to viscosity and the size of capillary required for the testing process is no longer necessary.

It is possible to "individually" allocate each sample to a capillary viscometer that is currently being used by means of the conventional MS-Windows® "drag and drop" method. This procedure makes it possible to increase the sample throughput.

The allocation between the sample and the viscometer is shown on the status display.

... easy, reliable and safe

The proved and tested AVSPro II software also makes it possible to prepare additional individually selected calculations, such as:

- mean value,
- standard deviation,
- outlier test (A %),
- Hagenbach correction,
- absolute viscosity, dynamic viscosity (density value required),
- viscosity index (measurement at two temperatures required),
- SUS and SFS,
- relative viscosity,
- specific viscosity,
- reduced viscosity (viscosity number),
- inherent viscosity and
- K-value.

During the entire process, all of the parameters (depending on the menu level) and the respective status of the individual measuring positions, the temperature regulation system and the sample transfer system are either visible or can be selected.

The operator interface of the AVSPro II is available in German and English. Commercially available printers for which Windows drivers are available are suitable for documentation purposes.

Screenshot: selection of methods


This mode is used to specify the number of measurements, the preliminary temperature regulation pe-riod, the allowable standard deviation, the maximum allowable temperature tolerance, the rinsing type and method of the viscometer.

Screenshot: options

This mode is used to specify what monitoring parameters are to be activated, e.g. if the temperature control of the thermostats is supposed to be handled via the PC.

Screenshot: dosing parameters

This mode is used to specify the filling quantity of the viscometer, the dosing speed depending on the viscosity and the type of rinsing.



Precision, reproducibility and comparability are in compliance with the DIN 51 562-1(1995-08), ASTM D 445 and ISO/DIS 3105 standards.

The AVSPro II is built in accordance with international equipment safety standards: CE symbol (equipment safety, low voltage safety, emitted interference and interference immunity).

The AVSPro II is produced by a manufacturer that is certified in accordance with DIN/ISO9001.

If requested, the AVSPro II automatic sampler can be supplied with a manufacturer's inspection certificate based on direct comparison with normal viscometers of the first order in accordance with DIN 51 532 - 4: 1995-08.

Technical data AVSPro II

AV

Sampling system	
Sample bottles	100 ml screw-type and bottles with standard ground joint (16 pcs per rack) 20 ml round bottom glass pieces (56 pcs. per rack)
Sample rack	for 100 ml screw-type and bottles with standard ground joint for 100 ml screw-type and bottles with standard ground joint (temperature controlled up to 135°C) for 20 ml round bottom glass pieces
Measured value recording	
Method	meniscus scanning by means of opto-electronic system or thermal conductivity (TC)
Measuring parameter	
	throughput time in seconds [s] temperature in degrees Celsius [°C]
Calculated parameters	
	mean value, standard deviation, outlier test (A %), Hagenbach correction, absolute viscosity, dynamic viscosity (knowledge of density required), viscosity index (measurement at two temperatures required) SUS and SFS, relative viscosity, specific viscosity, reduced viscosity (viscosity number), inherent viscosity, K-value
Selection parameters	
	by means of PC keyboard, mean value, standard deviation, outlier test (A %), Hagenbach correction, absolute viscosity, dynamic viscosity (knowledge of density required), viscosity index (measurement at two temperatures required) SUS and SFS, relative viscosity, specific viscosity, reduced viscosity (viscosity number), inherent viscosity, K-value, rack position, date/time, temperature regulation period, number of measurements, number of rinsing operations, start, stop/reset
Number of measurements	1 ... 99
Temperature regulation period	0 ... 99 min., selectable in increments of 1 min.
Number of Viscometer tests	0 ... 9 with next sample (observe sample quantity) or with preselected rack position
Data memory	by means of PC
Viscosity measurement range	
Time	0.35 to 1,200 mm ² /s (at room temperature of samples) up to 9999.99 s, resolution = 0.01 s
Vacuum pressure	automatically controlled
Viscometers available for use	
	Ubbelohde viscometer in accordance with DIN standards Ubbelohde viscometer in accordance with ASTM standards Micro-Ubbelohde viscometer in accordance with DIN standards Micro-Ostwald viscometer Cannon-Fenske-Routine viscometer Cannon-Fenske-Routine viscometer TC Ubbelohde viscometer TC Micro-Ubbelohde viscometer

/S Pro II

Measuring accuracy	$\pm 0.01 \text{ s} \pm 1 \text{ digit}$, but not more precise than 0.01% The measuring uncertainty for measurements of absolute kinematic viscosity is also dependent on the uncertainty of the numeric value for the viscometer constant and on the measuring conditions, especially the measuring temperature.
Evaluations / results	
Correction	Hagenbach correction (HC for Ubbelohde, Cannon-Fenske-Routine, Micro-Ubbelohde and Micro-Ostwald viscometers)
Statistical evaluation	standard deviation, outlier search
Ambient conditions	
Ambient temperature	10 ... + 40 °C
Air humidity	max. 85 % relative humidity
Equipment safety	
CE-symbol	in accordance with Guideline 89/336/EEC of the Council (EMC compatibility) in accordance with Standard EN 50 081, Part 1; interference immunity in accordance with Standard EN 50 082, Part 2; in accordance with Guideline 73/23/EEC of the Council (low-voltage guideline)
Housing	plastic/stainless steel / aluminium casing with chemically resistant two-component coating of the plastic pieces
Dimensions	w = 1,300 mm, h = 1,100 mm, d = 610 mm (approx. 51" x 43" x 24")
Weight	dependent on the number of measuring positions approx. 70 kg
Connections	
Pneumatic connections	screw-type connections for viscometer
Electric connections	circular connectors with bayonet lock for measuring stand and TC viscometer
Viscometers	up to 8 viscometers connected by individual control units
Temperature	via serial interface RS-232-C of suspended thermostat, type: 1 pc, CT 1654 or up to 2 pcs. CT 53 made by SCHOTT Instruments
Interfaces	control system using PC with 2 x RS-232-C interfaces
Safety	overfilling safety device or waste bottle
Mains connection	European built-in plug DIN 49 457 6 with fuse
Data transmission	
Interface internal	bidirectional serial interface in accordance with EIA RS-232-C (daisy chain concept)
Interface external	via PC, bidirectional serial interface in accordance with EIA RS-232-C
Power supply	
Mains voltage	230 V (AC) or 115 V (AC), 50 ... 60 Hz (AC)

Transparent thermostats – CT series

High temperature constancy and visual observation

Transparent thermostats manufactured by SCHOTT Instruments have been specially designed to measure the viscosity of Newtonian liquids in capillary viscometers. They can be used for both manual measurements and, when used in connection with viscosity measuring equipment, for automatic measurements. The most important characteristics of the transparent thermostats are their ability to maintain a constant temperature and capability of visual observation of the flow of the fluid in the viscometer.

The transparent thermostats of the series CT 53 and CT 54 are suitable for viscosity measurements in compliance with DIN 51 562 (Part 1) and ASTM D 445. They are comprised of a stainless steel bath with insulating glass, a coated steel casing and a bath thermostat. In addition, the CT 54 has an integrated discharge outlet on one side to drain the bath. An RS-232-C interface enables it to be connected to a PC.



The viscosity of Newtonian liquids is extremely dependent on temperature. Depending on the measurement medium, a deviation of 0.5 to 2% can be expected for a temperature deviation of 0.1 K. For this reason, the significant influence of temperature on the viscosity of a fluid must be taken into consideration when selecting a thermostat. All SCHOTT Instruments thermostats have the possible temperature stability of 0.01 K (see Technical Specifications) under optimum ambient conditions.

CT 52

The transparent thermostat CT 52 is made of acrylic glass and it is able to take up to two automatic measurement positions or brackets for manual measurements. Due to its design its ability to keep the temperature constant is not quite as great and it can only be used up to temperatures of +60 °C. If temperature constancy is not a top priority, the CT 52 is a cost-effective alternative.

CT 52



CT 53

This thermostat is for use with temperatures between +5 °C and +102 °C. Between +5 °C and +50 °C cooling is recommended to maintain the temperature constancy. Either a flowthrough cooler (e.g. CK 300, see Accessories) or simple cooling with circulated water can be used.

CT 53 HT

The high temperature version of the thermostats is used to measure viscosity at temperatures above +80 °C (see Technical specification).

2 or 4 measuring points

All CT 53 models enable the user to position 2 measurement stands or brackets into the thermostats. Up to 4 micro-TC viscometers can be positioned in the thermostats using the special VZ 7191 holder.

CT 53 TT

These thermostats have been specifically designed for use at temperatures well below room temperature. Viscosity measurements can be measured between -40 °C and +102 °C. For measurements under +5 °C, a cryostat is necessary to reach the low temperatures.

Transparent thermostats – CT series

CT 54

The main technical features of this thermostat are identical to those of the CT53. The main differences are the number of measuring positions (the number doubles) and the additional discharge outlet to drain the bath. Up to 8 micro-TC-viscometers can be positioned in the thermostats if 2 special VZ 7191 holders are used.



Recommended temperature equalization fluids

Fluid	Alcohol	Water	Paraffin oil	Silicon oil
Temperature range	-40 °C ...+10 °C	+5 °C ...+ 80 °C	+40 °C ...+150 °C	+80 °C ...+150 °C

The following applies to all temperature equalization fluids:

The viscosity of the temperature equalization fluid should be a max. 10 mm²/s (cSt) at 25 °C.

Technical specifications	CT 52	CT 53 TT**	CT 53	CT 53 HT	CT 54
Operating temperature	+10 ...+60 °C	-40 ...+102 °C	+5 ...+102 °C	+5 ...+150 °C	+5 ...+102 °C
Measurement points for AVS/S	2	2	2	2	4
Measurement points for TC	2	2	2	2	4
Measurement points Micro-TC	2	4*	4*	4*	8*
Temperature constancy in compliance with DIN 58 966 at 25 °C	± 0.02 K	± 0.01 K	± 0.01 K	± 0.01 K	± 0.01 K
Dimensions (W x H x D in mm)	355 x 370 x 250	355 x 370 x 250	355 x 370 x 250	355 x 370 x 250	605 x 370 x 250
Volume	18 l	15 l	15 l	15 l	27 l
Weight (empty)	approx. 5 kg	approx. 13.5 kg	approx. 13.5 kg	approx. 13.7 kg	approx. 28 kg

When using in the normal temperature range it is necessary to cool the system (+5 °C to approx. +50 °C) to maintain the temperature constancy. This can be achieved either by using circulated water or a flow-through cooler (e.g. CK300).

* When using with 4 or 8 micro-TC-viscometers, a special holder (Type: VZ 7191) is required for two of the existing measurement positions.

** When using at very low temperatures (well below room temperature) a cryostat is required. Cryostats can be included in delivery (at the list price of the manufacturer).

Accessories



Flow-through cooler CK300

The fluor-hydro-carbon gas free flow-through cooler serves to mechanically cool the bathfluid and is filled with environment compatible coolant R 134a. It works at ambient temperatures between +5 and +50 °C. Refrigeration power is 300 W at 20 °C. The CK300 is very compact (200 x 430 x 300 mm, W x H x D) and very stable (approx. 25 kg).

The flow-through cooler is available for all standard international voltage and frequency ranges (230 V, 50 Hz; 115 V, 60 Hz). Other voltage and frequency combinations are only available upon request as special productions.



AVS Measurement stands

For the use of capillary viscometers for opto-electronic measurement recording.

Type no.: AVS/S	Metal measurement stand, preferably for non-aqueous bath fluids	For use with temp. range: -80 ...+80 °C
Type no.: AVS/SK	PVDF measurement stand, corrosion free, suitable for both aqueous and non-aqueous bath fluids	For use with temp. range: 0 ...+80 °C
Type no.: AVS/SK-CF	PVDF measurement stand specially for use of Cannon-Fenske-Routine viscometers	For use with temp. range: 0 ...+80 °C
Type no.: AVS/SK-V	PVDF measurement stand specially for the use of dilution viscometers	For use with temp. range: 0 ...+80 °C
Dimensions	90 x 90 x 460 mm (W x D x H)	
Weight	approx. 1.1 kg	

We reserve the right to make technical changes.

Viscometers and their range of use

Measurement substance property	Viscometer type							
	Ubbelohde	Micro Ubbelohde	TC Ubbelohde	Ostwald	Micro Ostwald	Cannon-Fenske	Cannon-Fenske-Routine	BS/IP-U tube reverse flow
Transparent liquids manual measurement	++	++	-	+	+	+	o	o
Transparent liquids automatic measurement	++	++	+	-	+	+	-	-
Opaque liquids manual measurement	-	-	-	-	-	-	+	+ ²⁾
Opaque liquids automatic measurement	-	-	++ ¹⁾	-	-	-	-	-
Foaming liquids	o	o	o	+	+	+	o	o
Liquid mixture with highly volatile components	o	o	o	+	+	+	o	o
Minimum measurement substance and/or rinsing agent quantities	-	++	-	-	+	-	-	-
High-temperature or low- temperature measurements	+	+	+	o	o	o	o	o

Selection of glass capillary viscometers

++ use preferably
+ highly suitable
o less suitable
- unsuitable


¹⁾ up to 30,000 mm²/s

²⁾ above 30,000 mm²/s

Ubbelohde viscometers, normal form

Viscometers with suspended ball level for determination of absolute and relative kinematic viscosity of liquids with Newtonian flow behavior. The calibrated viscometers are delivered with manufacturer's certificate in accordance with DIN 55 350, Part 18.

All viscometers are provided with ring marks. This ensures that viscometers for automatic measurements can also be checked by means of manual measurements. The recommended minimum flowthrough time is 200 s.



Ubbelohde-Viskosimeter (DIN)

- in accordance with DIN 51 562 Part 1, ISO/DIS 3105 (BS-IP-SL)
- filling quantity: 15 ... 20 ml
- overall length: approx. 290 mm

calibrated,
with constant,
for manual measurements

calibrated
with constant,
for manual measurements;
automatic measurement with
stand AVS/SK-HV

$$v = K \cdot t$$

$$K = \frac{v}{t}$$

$$t = \frac{v}{K}$$

v = kinematic viscosity in mm^2/s
 K = constant [mm^2/s]
 t = flow-through time in s

Type No.	Order No.	Type No.	Order No.	Capillary No. acc. DIN	acc ISO	Capillary $\varnothing \pm 0,01$ [mm]	Constant K (approx.)	Measuring range [mm^2/s] (approx.)
501 00	285400004	-	-	0	-	0.36	0.001	0.3 ... 1
501 03	285400012	-	-	0c	-	0.47	0.003	0.5 ... 3
501 01	285400029	-	-	0a	-	0.53	0.005	0.8 ... 5
501 10	285400037	-	-	I	I	0.63	0.01	1.2 ... 10
501 13	285400045	-	-	Ic	Ia	0.84	0.03	3 ... 30
501 11	285400053	-	-	Ia	-	0.95	0.05	5 ... 50
501 20	285400061	-	-	II	II	1.13	0.1	10 ... 100
501 23	285400078	-	-	IIc	Ila	1.50	0.3	30 ... 300
501 21	285400086	-	-	Ila	-	1.69	0.5	50 ... 500
501 30	285400094	-	-	III	III	2.01	1	100 ... 1000
501 33	285400107	-	-	IIIC	IIla	2.65	3	300 ... 3000
501 31	285400115	-	-	IIla	-	3.00	5	500 ... 5000
501 40	285400123	-	-	IV	IV	3.60	10	1000 ... 10000
-	-	502 43	285400131	IVc	IVa	4.70	30	3000 ... 30000
-	-	502 41	285400148	IVa	-	5.34	50	6000 ... 30000
-	-	502 50	285400156	-	V	6.30	100	> 10000

not calibrated,
without constant;
for determination of rela-
tive viscosity

calibrated,
with constant for auto-
matic measurements

$$v = K \cdot t$$


$$K = \frac{v}{t}$$

$$t = \frac{v}{K}$$

v = kinematic viscosity in mm^2/s
 K = constant [mm^2/s]
 t = flow-through time in s

Type No.	Order No.	Type No.	Order No.	Capillary No. acc. DIN	acc ISO	Capillary $\varnothing \pm 0,01$ [mm]	Constant K (approx.)	Measuring range [mm^2/s] (approx.)
-	-	532 00	285400164	0	-	0.36	0.001	0.3 ... 1
530 03	285400304	532 03	285400201	0c	-	0.47	0.003	0.5 ... 3
530 01	285400312	532 01	285400218	0a	-	0.53	0.005	0.8 ... 5
530 10	285400329	532 10	285400226	I	I	0.63	0.01	1.2 ... 10
530 13	285400337	532 13	285400234	Ic	Ia	0.84	0.03	3 ... 30
-	-	532 11	285400172	Ia	-	0.95	0.05	5 ... 50
530 20	285400345	532 20	285400242	II	II	1.13	0.1	10 ... 100
530 23	285400353	532 23	285400259	IIc	Ila	1.50	0.3	30 ... 300
-	-	532 21	285400189	Ila	-	1.69	0.5	50 ... 500
530 30	285400361	532 30	285400267	III	III	2.01	1	100 ... 1000
530 33	285400378	532 33	285400275	IIIC	IIla	2.65	3	300 ... 3000
-	-	532 31	285400197	IIla	-	3.00	5	500 ... 5000
530 40	285400386	532 40	285400283	IV	IV	3.60	10	1000 ... 10000

Ubbelohde viscometers, normal form (ASTM)



Ubbelohde Viscometer (ASTM)

- in accordance with ISO/DIS 3105, ASTM D 2515, ASTM D 446
- filling quantity: 15 ... 20 ml
- overall length: approx. 285 mm

calibrated,
with constant for manual
measurements

not calibrated,
without constant for
determination of relative
Viscosity


calibrated,
with constant for automa-
tic measurements

Type No.	Order No.	Type No.	Order No.	Type No.	Order No.	Capillary No.	Capillary Ø i ± 0,01 [mm]	Constant K (approx.)	Measuring range [mm ² /s] (approx.)	
525 00	285400501	526 00	285400707	527 00	285401255	0	0.24	0.001	0.35 ...	1
525 03	285400518	526 03	285400715	527 03	285401271	0c	0.36	0.003	0.6 ...	3
525 01	285400526	526 01	285400723	527 01	285401263	0b	0.46	0.005	1 ...	5
525 10	285400534	526 10	285400731	527 10	285401152	I	0.58	0.01	2 ...	10
525 13	285400542	526 13	285400748	527 13	285401169	Ic	0.78	0.03	6 ...	30
525 20	285400559	526 20	285400756	527 20	285401177	II	1.03	0.1	20 ...	100
525 23	285400567	526 23	285400764	527 23	285401185	IIc	1.36	0.3	60 ...	300
525 30	285400575	526 30	285400772	527 30	285401193	III	1.83	1	200 ...	1000
525 33	285400583	526 33	285400789	527 33	285401288	IIIc	2.43	3	600 ...	3000
525 40	285400591	526 40	285400797	527 40	285401296	IV	3.27	10	2000 ...	10000
525 43	285400604	526 43	285400801	527 43	285401309	IVc	4.32	30	6000 ...	30000

Ubbelohde viscometers, with additional tube and threads

Viscometers with suspended ball level for determination of absolute or relative kinematic viscosity. These viscometers are preferably used for automatic measurements when an AVS 24, AVS 26 or AVS 270 automatic viscometer cleaner is used simultaneously. The additional filling and cleaning tube and the glass

thread ensure safe operational use. The calibrated viscometers are delivered with manufacturer's certificate in accordance with DIN 55 350, Part 18. The ring marks that are also present serve as auxiliary marks in case the viscometers must be checked by means of manual measurements.




Ubbelohde viscometer (DIN)

- in accordance with ISO/DIS 3105, DIN 51 562, Part 1, BS 133, NFT 60-100
- filling quantity: 18 ... 22 ml
- overall length: approx. 290 mm

calibrated,
with constant for automatic measurements

Type No.	Order No.	Capillary No. acc. DIN	acc. ISO	Capillary Ø i [mm]	Constant K (approx.)	Measuring range [mm ² /s] (approx.)
541 03	285401925	0c	-	0.47	0.003	0.5 ... 3
541 01	285401917	0a	-	0.53	0.005	0.8 ... 5
541 10	285401933	I	I	0.63	0.01	1.2 ... 10
541 13	285401941	Ic	Ia	0.84	0.03	3 ... 30
541 20	285401958	II	II	1.13	0.1	10 ... 100
541 23	285401966	IIc	IIa	1.50	0.3	30 ... 300
541 30	285401974	III	III	2.01	1	100 ... 1000
541 33	285401982	IIIc	IIIa	2.65	3	300 ... 3000
541 40	285401999	IV	IV	3.60	10	1000 ... 6000




Ubbelohde viscometer (ASTM)

- the technical measurement characteristics are in accordance with ISO/DIS 5105, ASTM D 2515, ASTM D 446
- filling quantity: 15 ... 22 ml
- overall length: approx. 290 mm

calibrated,
with constant for automatic measurements

Type No.	Order No.	Capillary No. acc. DIN	Capillary Ø i [mm]	Constant K (approx.)	Measuring range [mm ² /s] (approx.)
545 00	285402005	0	0.24	0.001	0.35 ... 1
545 03	285402021	0c	0.36	0.003	0.6 ... 3
545 01	285402013	0b	0.46	0.005	1 ... 5
545 10	285402038	I	0.58	0.01	2 ... 10
545 13	285402046	Ic	0.78	0.03	6 ... 30
545 20	285402054	II	1.03	0.1	20 ... 100
545 23	285402062	IIc	1.36	0.3	60 ... 300
545 30	285402079	III	1.83	1	200 ... 1000
545 33	285402087	IIIc	2.43	3	600 ... 3000
545 40	285402095	IV	3.27	10	2000 ... 10000
545 43	285402108	IVc	4.32	30	6000 ... 30000

Ubbelohde viscometers with TC sensors



Viscometers with suspended ball level for determination of absolute and relative kinematic viscosity of liquids with Newtonian flow behaviour. The measuring levels are marked by TC sensors. The meniscus passage is detected due to the different conductivity of the liquid phase and the gas phase. A measurement stand of the type series AVS/S is not required. TC viscometers can be used to determine the kinematic viscosity of all liquids with Newtonian flow behaviour.

They are especially suitable for liquids that cannot be detected with other systems: untransparent and/or black and/or electric conductive measuring samples.

TC viscometers are manufactured from technical glass types with an expansion coefficient of $\alpha = \text{ca. } 9 \cdot 10^{-6}$. Due to the electric properties of TC sensors, it is important to make sure that a type is selected that is suitable for the required application temperature.

TC viscometers with additional filling and cleaning tube and with glass thread

- the technical measurement characteristics are in accordance with DIN 51 562, part 1, ISO/DIS 3105 (BS-IP-SL)
- for use in combination with an automatic viscosity measuring instrument and an AVS 24, AVS 26 or AVS 270 automatic viscometer cleaner
- filling quantity: 18 ... 22 ml
- overall length: approx. 355 mm

calibrated,
with constant for automatic measurements


Type No.	Order No.	Type No.	Order No.	Type No.	Order No.	Capillary No.	Capillary Ø i [mm]	Constant K (approx.)	Measuring range [mm ² /s] (approx.)	
+10 ... +80 °C		-40 ... +30 °C		+70 ... +150 °C						
562 03	285423120	–	–	–	–	0c	0.47	0.003	0.5	... 3
562 10	285423130	563 10	285423240	564 10	285423330	I	0.54	0.01	1,2	... 10
562 13	285423140	563 13	285423250	564 13	285423340	Ic	0.84	0.03	3	... 30
562 20	285423150	563 20	285423260	564 20	285423350	II	1.15	0.1	10	... 100
562 23	285423170	563 23	285423270	564 23	285423360	IIc	1.51	0.3	30	... 300
562 21	285423160	–	–	–	–	IIa	1.69	0.5	50	... 500
562 30	285423180	563 30	285423280	564 30	285423370	III	2.05	1	100	... 1000
562 33	285423200	563 33	285423290	564 33	285423380	IIIc	2.7	3	300	... 3000
562 31	285423190	–	–	–	–	IIIa	3.0	5	500	... 5000
562 40	285423210	563 40	285423300	564 40	285423390	IV	3.7	10	1000	... 10000
562 43	285423230	563 43	285423320	564 43	285423400	IVc	4.9	30	3000	... 20000
562 41	285423220	563 41	285423310	–	–	IVa	5.3	50	5000	... 30000

Ubbelohde viscometers with TC sensors

Viscometers with suspended ball level for determination of absolute and relative kinematic viscosity of liquids with Newtonian flow behaviour. The measuring levels are marked by TC sensors. The meniscus passage is detected due to the different conductivity of the liquid phase and the gas phase. A measurement stand of the type series AVS/S is not required. TC viscometers can be used to determine the kinematic viscosity of all liquids with Newtonian flow behaviour.

They are especially suitable for liquids that cannot be detected with other systems: untransparent and/or black and/or electric conductive measuring samples.

TC viscometers are manufactured from technical glass types with an expansion coefficient of $\alpha = \text{ca. } 9 \cdot 10^{-6}$. Due to the electric properties of TC sensors, it is important to make sure that a type is selected that is suitable for the required application temperature.




TC viscometers

- the technical measurement characteristics are in accordance with DIN 51 562, part 1, ISO/DIS 3105 (BS-IP-SL)
- for use in combination with an automatic viscosity measuring instrument and an AVS 24, AVS 26 or AVS 270 automatic viscometer cleaner
- filling quantity: 18 ... 22 ml
- overall length: ca. 355 mm

calibrated,
with constant for automatic measurements

Type No.	Order No.	Type No.	Order No.	Type No.	Order No.	Capillary No.	Capillary Ø i [mm]	Constant K (approx.)	Measuring range [mm ² /s] (approx.)
+10 ... +80 °C		-40 ... +30 °C		+70 ... +150 °C					
567 03	285423420	–	–	–	–	0c	0.47	0.003	0.5 ... 3
567 10	285423430	568 10	285423540	569 10	285423630	I	0.64	0.01	1.2 ... 10
567 13	285423440	568 13	285423550	569 13	285423640	Ic	0.84	0.03	3 ... 30
567 20	285423450	568 20	285423560	569 20	285423650	II	1.15	0.1	10 ... 100
567 23	285423470	568 23	285423570	569 23	285423660	IIc	1.51	0.3	30 ... 300
567 21	285423460	–	–	–	–	IIa	1.69	0.5	50 ... 500
567 30	285423480	568 30	285423580	569 30	285423670	III	2.05	1	100 ... 1000
567 33	285423500	568 33	285423590	569 33	285423680	IIIc	2.7	3	300 ... 3000
567 31	285423490	–	–	–	–	IIIa	3.0	5	500 ... 5000
567 40	285423510	568 40	285423600	569 40	285423690	IV	3.7	10	1000 ... 10000
567 43	285423530	568 43	285423620	569 43	285423700	IVc	4.9	30	3000 ... 20000
567 41	285423520	568 41	285423610	–	–	IVa	5.3	50	5000 ... 30000

Micro-Ubbelohde viscometers with TC sensors



Viscometers with suspended ball level for determination of absolute and relative kinematic viscosity of liquids with Newtonian flow behaviour. The measuring levels are marked by TC sensors. The meniscus passage is detected due to the different conductivity of the liquid phase and the gas phase. A measurement stand of the type series AVS/S is not required. TC viscometers can be used to determine the kinematic viscosity of all liquids with Newtonian flow behaviour.

They are especially suitable for liquids that cannot be detected with other systems: untransparent and/or black and/or electric conductive measuring samples.

TC viscometers are manufactured from technical glass types with an expansion coefficient of $\alpha = \text{ca. } 9 \cdot 10^{-6}$. Due to the electric properties of TC sensors, it is important to make sure that a type is selected that is suitable for the required application temperature.

Micro TC viscometers

- the technical measurement characteristics are in accordance with DIN 51 562, Part 2
- for use in combination with an automatic viscosity measuring instrument
- filling quantity: 3 ... 4 ml
- overall length: approx. 350 mm

calibrated,
with constant for automatic measurements

Type No.	Order No.	Type No.	Order No.	Type No.	Order No.	Capillary No.	Capillary Ø i [mm]	Constant K (approx.)	Measuring range [mm²/s] (approx.)		
+10 ... +80 °C		-40 ... +30 °C		+70 ... +150 °C							
572 10	285423710	573 10	285423780	574 10	285423850	M I	0.40	0.01	0.4	...	6
572 13	285423720	573 13	285423790	574 13	285423860	M Ic	0.53	0.03	1.2	...	18
572 20	285423730	573 20	285423800	574 20	285423870	M II	0.70	0.1	4	...	60
572 23	285423740	573 23	285423810	574 23	285423880	M IIc	0.95	0.3	12	...	180
572 30	285423750	573 30	285423820	574 30	285423890	M III	1.26	1	40	...	800

Micro-Ubbelohde viscometers


Viscometers for dilution viscometry

Viscometers with suspended ball level for determination of absolute and relative kinematic viscosity of liquids with Newtonian flow behavior. Due to their design, these viscometers are especially suitable for measurement of small liquid quantities and for particularly short running times. All viscometers are provided with ring marks. This ensures that viscometers for automatic measurements can also be checked by means of manual measurements.

Micro-Ubbelohde viscometers (DIN)

The calibrated viscometers are delivered with manufacturer's certificate in accordance with DIN 55 350, Part 18. For measurements with automatic viscosity measuring instruments, another constant is valid. This constant is determined by multiplication of the constant K with the correction factor F.

- in accordance with DIN 51562, Part 2
- filling quantity: 3 ... 4 ml
- overall length: approx. 290 mm



calibrated,
with constant for manual
measurement

calibrated,
with constant for
automatic measurement

not calibrated,
without constant;
for determination of
relative viscosity


Type No.	Order No.	Type No.	Order No.	Type No.	Order No.	Capillary No.	Capillary Ø i [mm]	Constant K (approx.)	Measuring range [mm²/s] (approx.)
536 10	285401009	537 10	285401103	538 10	285401206	M I	0.40	0.01	0.4 ... 6
536 13	285401017	537 13	285401111	538 13	285401214	M Ic	0.53	0.03	1.2 ... 18
536 20	285401025	537 20	285401128	538 20	285401222	M II	0.70	0.1	4 ... 60
536 23	285401033	537 23	285401136	538 23	285401239	M IIc	0.95	0.3	12 ... 180
536 30	285401041	537 30	285401144	538 30	285401247	M III	1.26	1	40 ... 800

Viscometers for dilution viscometry

Viscometers with suspended ball level designed according to the principle of the Ubbelohde viscometers for determination of the limit viscosity number of polymers. The limit viscosity number is determined automatically in combi-

nation with viscosity measuring instruments and piston burette, type AVS 20, made by SCHOTT Instruments.


- filling quantity: 15 ... 75 ml
- overall length: approx. 290 mm



not calibrated,
for automatic measurements

Type No.	Order No.	Capillary No.	Capillary Ø i [mm]	Constant K (approx.)	Measuring range [mm²/s] (approx.)
531 00	285401403	0	0.36	0.001	0.35 ... 0.6
531 03	285401428	0c	0.47	0.003	0.5 ... 2
531 01	285401411	0a	0.53	0.005	0.8 ... 3
531 10	285401436	I	0.64	0.01	1.2 ... 6
531 13	285401444	Ic	0.84	0.03	3 ... 20
531 20	285401452	II	1.15	0.1	10 ... 60

Cannon-Fenske viscometers



Cannon-Fenske routine viscometers comply with standards ISO/DIS 3105, ASTM D 2515, BS 188 with respect to technical measuring specifications.


- are suitable for all Newtonian liquids with a viscosity of 0.35...20,000 mm²/s

- the present design has, as a supplement to the standard, a deepening in the lower bend. Accordingly, these viscometers can also be used for automatic measurements.
- filling quantity: approx. 7 ... 10 ml
- overall length: approx. 245 mm

calibrated,
with ring mark,
for manual measurements

with constant
for automatic measurements

Type No.	Order No.	Type No.	Order No.	Capillary No.	Capillary Ø i [mm]	Constant K (approx.)	Measuring range [mm ² /s] (approx.)
513 00	285403507	520 00	285403704	25	0.30	0.002	0.4 ... 1.6
513 03	285403515	520 03	285403712	50	0.44	0.004	0.8 ... 3.2
513 01	285403523	520 01	285403729	75	0.54	0.008	1.6 ... 6.4
513 10	285403531	520 10	285403737	100	0.63	0.015	3 ... 15
513 13	285403548	520 13	285403745	150	0.78	0.035	7 ... 35
513 20	285403556	520 20	285403753	200	1.01	0.1	20 ... 100
513 23	285403564	520 23	285403761	300	1.27	0.25	50 ... 200
513 21	285403572	520 21	285403778	350	1.52	0.5	100 ... 500
513 30	285403589	520 30	285403786	400	1.92	1.2	240 ... 1200
513 33	285403597	520 33	285403794	450	2.35	2.5	500 ... 2500
513 40	285403601	520 40	285403807	500	3.20	8	1600 ... 8000
513 43	285403618	520 43	285403815	600	4.20	20	4000 ... 20000



Cannon-Fenske reverse flow viscometers

- Comply with standards ISO/DIS 3105, ASTM D 2515, ASTM D 446,

NF T 60 - 100 with respect to technical measuring specifications.

- filling quantity: approx. 12 ml
- overall length: approx. 295 mm


calibrated,
with 3 ring marks,
with 2 constants,
only for manual measurement

Type No.	Order No.	Capillary No.	Capillary Ø i [mm]	Constant K (approx.)	Measuring range [mm ² /s] (approx.)
511 00	285403001	25	0,31	0.002	0.4 ... 1.6
511 03	285403018	50	0,42	0.004	0.8 ... 3.2
511 01	285403026	75	0,54	0.008	1.6 ... 6.4
511 10	285403034	100	0,63	0.015	3 ... 15
511 13	285403042	150	0,78	0.035	7 ... 35
511 20	285403059	200	1,02	0.1	20 ... 100
511 23	285403067	300	1,26	0.25	50 ... 200
511 21	285403075	350	1,48	0.5	100 ... 500
511 30	285403083	400	1,88	1.2	240 ... 1200
511 33	285403091	450	2,20	2.5	500 ... 2500
511 40	285403104	500	3,10	8	1600 ... 8000
511 43	285403112	600	4,00	20	4000 ... 20000

Cannon-Fenske routine viscometers

comply with standards ISO/DIS 3105, ASTM D 2515, BS 188 with respect to technical measuring specifications. These viscometers are preferably used for automatic measurements when an AVS 24, AVS 26 or AVS 270 automatic viscometer cleaner is used simultaneously. The additional filling and cleaning tube and the glass thread ensure safe operational use. The calibrated viscometers are delivered with manufacturer's certificate in accordance with DIN 55 350, Part 18.


- are suitable for all Newtonian liquids with a viscosity of 0.35 ... 20,000 mm²/s.
- filling quantity: approx. 7 ... 12 ml
- overall length: approx. 245 mm



calibrated,
with ring marks,
with constant for automatic measurements

Type No.	Order No.	Capillary No.	Capillary Ø i [mm]	Constant K (approx.)	Measuring range [mm ² /s] (approx.)
546 00	285402116	25	0.30	0.002	0.4 ... 1.6
546 03	285402132	50	0.44	0.004	0.8 ... 3.2
546 01	285402124	75	0.54	0.008	1.6 ... 6.4
546 10	285402149	100	0.63	0.015	3 ... 15
546 13	285402157	150	0.78	0.035	7 ... 35
546 20	285402165	200	1.01	0.1	20 ... 100
546 23	285402181	300	1.27	0.25	50 ... 200
546 21	285402173	350	1.52	0.5	100 ... 500
546 30	285402198	400	1.92	1.2	240 ... 1200
546 33	285402202	450	2.35	2.5	500 ... 2500
546 40	285402219	500	3.20	8	1600 ... 8000
546 43	285402227	600	4.20	20	4000 ... 20000

Ostwald viscometers




Ostwald viscometers

- filling quantity: 3 ml
- overall length: approx. 220 mm

with ring marks,
without constant,
for manual measurements

Type No.	Order No.	Capillary Ø i [mm]	Transit time for water approx. [s]	Constant K (approx.)	for use from [mm ² /s] (approx.)
509 03	285404006	0.3	250	0.004	0.3
509 04	285404014	0.4	75	0.01	1
509 05	285404022	0.5	30	0.03	2.5
509 06	285404039	0.6	15	0.07	5.5
509 07	285404047	0.7	10	0.1	10



Micro-Ostwald viscometers

- are suitable for measurements of small liquid quantities even with extreme formation of foam.
- filling quantity: 2 ml
- overall length: approx. 290 mm

calibrated,
with ring marks,
with constant,
for manual measurements

calibrated,
with ring marks,
with constant,
for automatic measurements

Type No.	Order No.	Type No.	Order No.	Capillary No.	Capillary Ø i [mm]	Constant K (approx.)	Measuring range [mm ² /s] (approx.)		
516 10	285404203	517 10	285404306	I	0.43	0.01	0.4	...	6
516 13	285404211	517 13	285404314	Ic	0.60	0.03	1.2	...	18
516 20	285404228	517 20	285404322	II	0.77	0.1	4	...	60
516 23	285404236	517 23	285404339	IIc	1.00	0.3	12	...	180
516 30	285404244	517 30	285404347	III	1.36	1	40	...	800

Accessories

Brackets and stands

All brackets and stands are designed to ensure that the viscometers are held vertically. They also protect the viscometers from breakage. The maximum deviation is $< 1^\circ$. In application in conjunction with SCHOTT Instruments and other commercially available see-

through thermostats the viscometers can only be used with the appropriate stand or bracket.

For DIN Ubbelohde viscometers that are used as reference measuring standard, specifically modified bracket (VZ 5840) must be used.

Brackets made of stainless steel

suitable for use with all Ubbelohde viscometers

for manual and automatic measurements

Type No.	Order No.
053 92	285405043
VZ 5840 (accessory for reference measuring standard)	285417201

suitable for use with Ubbelohde viscometers with TC sensors

Type No.	Order No.
053 93	285405035

suitable for use with all reverse flow viscometers


(Cannon-Fenske and BS/IP U-tube viscometers) for manual and automatic measurements (not illustrated)

Type No.	Order No.
053 96	285405019

suitable for use with Micro-Ostwald viscometers

for manual and automatic measurements

Type No.	Order No.
053 97	285405027



Accessories



Stands made of PTFE

suitable for use with Cannon-Fenske routine viscometers
for automatic measurements only (not illustrated)

Type No.	Order No.
065 99	285405113

Brackets for reference measuring standard

DIN Ubbelohde viscometers which are used as testing standard should be stored in a specially modified viscometer bracket (053 92) according to official inspection / calibration authorities. The

extension set for the test standard (VZ 5840) guarantees vertical slope with a maximum deviation of $< 1^\circ$ and the centered positioning of the capillaries.

Type No.	Order No.
VZ 5840	285417201

Control thermometers

Type No.	Order No.	Measuring range °C			Graduation °C
VZ 2801	285415763	- 5	...	+ 38	1/10
VZ 2802	285415771	+ 33	...	+ 67	1/10
VZ 2803	285415788	+ 66	...	+ 102	1/10
VZ 2804	285415796	+ 95	...	+ 152	1/10
VZ 2901	285415809	+ 20	...	+ 25	1/100
VZ 2907	285417078	+ 22	...	+ 27	1/100
VZ 2908	285415825	+ 37	...	+ 42	1/100
VZ 2905	285415841	+ 45	...	+ 50	1/100
VZ 2906	285415858	+ 97	...	+ 101	1/100
VZ 2909	285417094	+ 132	...	+ 137	1/100

Accessories

LabPump

The LabPump VZ 5655 (not illustrated) used in manual measurements and semi-automatic measurements to suck and pump up solutions:

- filling of viscometers
- rinsing with the next sample
- sucking up between manual measurements
- emptying of viscometers without removing them from the thermostatic bath

Since the materials used and the connections of the LabPump VZ 5655 are made of PTFE or stainless steel, the pump is suitable for use with aggressive mediums.

The range of use for semi-automatic processing of samples, e.g. with a viscosity measuring instrument AVS 360, AVS 370 or AVS 470, is possible up to a viscosity of 30,000 mm²/s. For semi-automatic processing work, the PTFE tube combination with stand (see illustration) and the waste bottle, type no. VZ 5624, are used.

Type No.	Order No.
VZ 5655	1040755

Polyamide bracket



for use with Cannon-Fenske routine viscometers,
Cannon-Fenske reverse flow viscometers and all Ostwald viscometers
for manual measurements only

Type No.	Order No.
064 99	285405105


PTFE bracket

for use with Cannon-Fenske routine viscometers,
for automatic measurements only (not illustrated)

Type No.	Order No.
065 99	285405113



Accessories



Temperature stabilization jackets

In the absence of a see-through thermostat the temperature of capillary viscometers can be stabilized in this type of jacket using circulation thermostats in the temperature range 0 to 180 °C. The shape of the jacket and the number of holes in the support plate depend

upon the type of viscometer being used. The support plate has been designed to facilitate changing the viscometer when required. An additional hole is provided in the support plate so that a control thermometer can be fitted. A quick-action seal simplifies changing viscometers.

Temperature stabilization jacket with support plate for Ubbelohde viscometers

Type No.	Order No.	Item No.	Comment
577 00	285405508		complete, without viscometer
Component parts			
577 01	285405516	1	temperature stabilization jacket, straight
238 00	285405524	2	support plate with 4 silicone rings (d = 4, 6, 8 and 10 mm)
225 34	285405532	3	silicone O-ring, ND 60
072 34	285405549	4	quick-action seal, NW 60

Accessories


Temperature stabilization jacket with support plate for Cannon-Fenske reverse flow viscometers and Ostwald viscometers

Component parts


Type No.	Order No.	Item No.	Comment
577 01	285405516	1	temperature stabilization jacket, straight
225 34	285405532	3	silicone O-ring, ND 60
072 34	285405549	4	quick-action seal, NW 60

Silicone rings

Type No.	Order No.	d mm	D mm	h mm
228 11	285405808	4	10	5
228 14	285405816	6	16	5
228 16	285405824	8	16	5
228 17	285405832	10	16	5



AVS measuring stands and tube sets



AVS measuring stands

Measuring stands of the type series AVS/S can be used to measure the flow-through time in viscometers automatically.

The measuring stands can be connected to all measuring instruments made by SCHOTT Instruments for automatic measurement of viscosity and operate with all standard viscometers for repetitive measurements.

Automatic measurements have the following advantages:

- the repetitive standard deviation is less than for manual measurements
- the measurement is free from subjective factors of influence
- the results can be printed out and/or be automatically documented on data memory system
- automatic processing of sample series is available.

The use of different materials ensures unproblematic adaptation to existing measurement temperatures and applications.

The measuring stands or brackets can be exchanged at random.

The distance between the levels of the automatic optoelectronic unloading system is $40.00 \text{ mm} \pm 0.03 \text{ mm}$. When measuring stands are replaced, this results in a standard deviation of VK = 0.05 % for Ubbelohde viscometers.

For repetitive measurements with viscosity measuring instruments and Ubbelohde viscometers with measuring stands, the standard deviation VK = 0.03 %.

Manually calibrated Ubbelohde viscometers can also be used in AVS measuring stands. If the automatic sensing levels do not correspond to the ring marks, the superimposed meniscus detection system will provide a higher constant. The difference amounts to 0.1 % per millimeter of height offset.

Required tube/cable combinations

Required tube/cable combinations					
Viscometer type					
	517 ...	540 ...	542 ...	547 ...	531 ... ⁽³⁾
	520 ...	541 ...	543 ...	548 ...	
	530 ...	545 ...	544 ...	549 ...	
	532 ...	546 ...		552 ...	
	537 ...			553 ...	
				554 ...	
Instrument	Tube/cable combinations				
AVS 300 and	VZ 5505 ⁽¹⁾ or	VZ 5621 ⁽¹⁾ and	–	–	VZ 5857 ⁽¹⁾
AVS 310	VZ 5501 ⁽²⁾	VZ 5505 ⁽¹⁾			
AVS 350	VZ 5505 ⁽¹⁾ or	VZ 5623 ⁽²⁾	VZ 5606 ⁽¹⁾	VZ 5505 ⁽¹⁾ and	VZ 5857 ⁽¹⁾
	VZ 5501 ⁽²⁾			VZ 6226	
AVS 360 and	VZ 5104 ⁽¹⁾ or	VZ 5623 ⁽²⁾	VZ 5623 ⁽²⁾	VZ 5104 ⁽¹⁾ or	VZ 5104 ⁽¹⁾ or
AVS 361	VZ 5622 ⁽²⁾			VZ 5622 ⁽²⁾	VZ 5622 ⁽²⁾
AVS 400 and	VZ 5505 ⁽¹⁾ or	VZ 5621 ⁽¹⁾ and	–	–	VZ 5857 ⁽¹⁾
AVS 410	VZ 5501 ⁽²⁾	VZ 5505 ⁽¹⁾			
AVS 440 and	VZ 5505 ⁽¹⁾ or	VZ 5621 ⁽¹⁾ and	VZ 5606 ⁽¹⁾	VZ 5505 ⁽¹⁾ and	VZ 5857 ⁽¹⁾
AVS 450	VZ 5501 ⁽²⁾	VZ 5505 ⁽¹⁾		VZ 6226	

(1) Silicon tube

(2) PTFE tube (aggressive mediums)


(3) The necessary connection tube TZ 1607 (l = 1.5 m) is included at the piston burette ViscoDoser AVS 20.

AVS measuring stands

	Measuring stands AVS/S	AVS/S-HT	AVS/SK	AVS/S-CF	AVS/SK-V
Available viscometers	Ubbelohde viscometers in accordance with DIN, ASTM, ISO 3105, Micro-Ubbelohde viscometers, Micro-Ostwald viscometers			Cannon-Fenske-routine viscometer	Ubbelohde-dilution viscometer
Temperature range	-80 ...+100 °C	-80 ...+200 °C	0 ...+60 °C	-80 ...+100 °C	0 ...+60 °C other temperature ranges available on request
Suitable for use with the measuring units	AVS 310, AVS 350, AVS 360, AVS 361, AVS 370, AVS 410, AVS 450, AVS 470, AVS 500, AVSPro				
Suitable for use with the thermo-static baths	CT 52, CT 53, CT 53 HT, CT 53 TT, CT 54				
Suitable brackets (type no.)	05392 05397			no bracket required	
Electrical connection	Cable VZ 6225 for all measuring stands to all instruments (is included in hose sets VZ 5505, VZ 5622 and VZ 5857), control lamp as function display				
Distance between measuring levels	40.00 mm ± 0.03 mm at 25 °C				
Signal transmission	Optically using optical fibres from the measuring level in the stand head, converted into analogue signal from stand to measuring instrument				
Material	Aluminium, TiO ₂ -anodized		PVDF, stainless steel	Aluminium, TiO ₂ -anodized	PVDF, stainless steel
Dimensions (W x H x D) mm	90 x 447 x 90	90 x 496 x 90	90 x 447 x 90	90 x 447 x 90	90 x 447 x 90
Weight (kg) appr.	1.0	1.25	0.8	1.0	0.8
Accessories included in scope of delivery	Bracket Type No. 05392 for Ubbelohde viscometers, tube/cable combination VZ 5505			tube/cable combination VZ 5505	tube/cable combination VZ 5857, magnetic stirring rods, fastening springs for viscometer

Note:

When TC viscometers are being used, a bracket type no. 05393, with the necessary tube set is required only. A measuring stand is not required.



Hotplates and stirrers:

Attractive design for the laboratory with glass ceramic heating surface

Perfect function, exclusive design

The laboratory hotplates and stirrers from SCHOTT Instruments demonstrate the kind of creative design that results when practical laboratory experience is paired with the most advanced materials know-how. The perfectly nonporous surface of the unique Ceran glass ceramic material makes the surface nearly indestructible. And the touch control panel gives the laboratory hotplates a truly exclusive touch. But you will hardly even notice that anymore once you have experienced how eminently practical they are.

Our know-how – your benefit

Both laboratory hotplates SLK 1 and SLK 2 enable speedy heating up. The laboratory stirrer (SLR) is the optimum solution for a careful to intense mixing of liquids. It can also be used for speedy heating up or controlled temperature adaptation.

Both product families have the benefits of the glass ceramic heating surface which has proven in millions of households. Chemical resistance, a high surface quality, and a resistance to temperature shocks of more than 700 °C provide the user with maximum benefits compared to conventional heating surface materials. The poreless smooth surface enables even most stubborn dirt to be easily removed. The high infrared permeability ensures that the heating energy is transferred quickly and with a low loss rate, i.e. it heats liquids faster than other heating surface materials, and thus saves time and energy.


Quality and safety

As a matter of course, our laboratory hotplates and laboratory stirrers bear the CE sign and are developed and produced according to high international quality standards.

A residual-heat display protects the user from the hazard of injuries (burning).

The corrosion-resistant, solid, casing with a hermetically sealed, non-inflammable top made of duroplastic in the case of the laboratory hotplates and of coated die-cast aluminium in the case of the laboratory stirrer ensure a long and trouble-free use of the devices.

If the laboratory hotplates are to be used in an aggressive environment, we recommend the use of the compressed-air connector which is available as an option. The stand rod holder on the back panel of the device can be used to connect accessories such as a temperature sensor conveniently and securely.



Overview of features

Hotplate / stirrer	SLK1	SLK2	SLR
Glass ceramic	■	■	■
Heating	■	■	■
Number of heating plates	1	1	1
Temperature control	-	-	■
Stirring	-	-	■
Touch control	■	■	-
Control knobs	-	-	■

Fast heating using the SLK 1 and SLK 2

To heat up liquids, the SLK 1 and SLK 2 laboratory hotplates from SCHOTT Instruments are the optimum solution. The heating power of the infrared radiation heating element can be adjusted in nine steps, with an average heating output of 1.2 kW or 1.8 kW, respectively, on step 9. During the heating process, the temperature distribution across the hotplates surface is almost homogenous.

Technical data	SLK 1	SLK 2
Heating function		
heating power (kW)	1.2	1.8
heated zone (mm)	Ø 165	Ø 200
max. hot plates temperature (°C)	approx. 600	approx. 600
min. time to boiling point 1 l H ₂ O* (min)	approx. 10	approx. 7
hot plate material	glass ceramic	glass ceramic
hot plates area (mm)	280 x 280	280 x 280
General data		
dimensions (L x W x H in mm)	395 x 295 x 110	395 x 295 x 110
weight (kg)	approx. 3.6	approx. 3.6
max. load (kg)	25	25
admissible ambient temperature (°C)	10 - 40	10 - 40
admissible air humidity (%)	85	85
protection type	IP 20	IP 20
protection class	1	1
housing material	SMC	SMC
thread for stand attachment	M 8	M 8
cable connector socket	cold appliances	cold appliances
mains connection (V/Hz)	230 V, 50/60 or 115 V, 50/60	230 V, 50/60 or 115 V, 50/60
Order No. 230 V	28 541 6316	28 541 6324
Order No. 115 V***	28 541 6213	28 541 6221

* measured in 3 l glass beaker at 25 °C ambient temperature and 1 bar air pressure

** depending on liquid quantity, heat capacity, and ambient temperature

*** the above measurement values may deviate slightly when using the 115 V version

• CE sign

Council Directive 89/336/EMC (electromagnetic compatibility)

Council Directive 73/23/EMC (low-voltage directive), last modified by the Council Directive 93/68/EMC



Stirring, heating and controlling using the SLR

In addition to fast or temperature-controlled heating, a lot of applications also require liquids to be stirred. Using our laboratory stirrer SLR with heating, the process of mixing liquids can be selected from careful to intense, and the device can also be used for speedy heating up or controlled temperature adaptation.

All functions can be viewed and monitored on the large and clear LCD display. The stirrer and heating are controlled separately by convenient turning knobs.


The rpm range of the stirrer stretches from 100 to 1100 min⁻¹ and can be set in steps of 10 min⁻¹. The mean rpm is also indicated in the form of a bar graph.

ph. Even in the lower rpm range, the smooth-running properties of the device are excellent.

The heating power can be set in 24 steps and reaches an average heating output of 0.9 kW at step 24. If a Pt 1000 temperature sensor (optional) is connected, temperature-controlled work with fluctuations of $\pm 2\text{ }^{\circ}\text{C}$ to $5\text{ }^{\circ}\text{C}$ as a function of liquid volume, heat capacity, and ambient temperature is possible between $25\text{ }^{\circ}\text{C}$ and $200\text{ }^{\circ}\text{C}$. The display will inform the user at intervals of 5 s alternately about the set and the actual temperature. In this case, too, the mean rpm is indicated in the form of a bar graph.

Technical data	SLR
Heating function	
heating power (kW)	0.9
heated zone (mm)	Ø 155
max. hot plates temperature (°C)	approx. 550
min. time to boiling point 1 l H ₂ O* (min)	approx. 15
temperature sensor connector	yes, Pt 1000
setting accuracy with temperature sensor (°C)	1
controlling accuracy with temperature sensor** (°C)	$\pm 2 \dots 5$
hot plate material	glass ceramic
hot plates area (mm)	235 x 235
digital set/actual temperature display (temperature sensor connector)	yes
Stirring function	
max. rpm (min ⁻¹)	100 - 1100
setting accuracy rpm (min ⁻¹)	10
max. stirring volume (l H ₂ O)	10
digital set/actual rpm display	yes
General data	
dimensions (L x W x H in mm)	370 x 240 x 85
weight (kg)	approx. 3.8
max. load (kg)	25
admissible ambient temperature (°C)	10 - 40
admissible air humidity (%)	85
protection type	IP 20
protection class	1
housing material	die-cast
cable connector	fixed cable
mains connection (V/Hz)	230 V, 50/60 or 115 V, 50/60
Order No. 230 V	28 541 6373
Order No. 115 V***	28 541 6279

- * measured in 3 l glass beaker at $25\text{ }^{\circ}\text{C}$ ambient temperature and 1 bar air pressure
- ** depending on liquid quantity, heat capacity, and ambient temperature
- *** the above measurement values may deviate slightly when using the 115 V version
- CE sign
Council Directive 89/336/EMC (electromagnetic compatibility)
Council Directive 73/23/EMC (low-voltage directive), last modified by the Council Directive 93/68/EMC



Accessories

Description	Type No.	Order No.
Temperature sensor stainless steel shaft (V4A), Pt 1000 sensor, 1 m fixed cable with 2 x 4 mm banana plug, length 170 mm, Ø 4 mm, -30 ...+ 200 °C	W5791NNHT	28 510 5308
Temperature sensor glass shaft, Sensor Pt 1000, 1 m fixed cable with 2 x 4 mm banana plug, length 250 mm, Ø 6 mm, -30 ...+ 200 °C	W5780NNHT	28 510 5238
Stand rod with fixing nut (M8) stainless steel, Ø 10 mm, length 450 mm	Z 601	28 541 6492
Temperature sensor holder clamp with extension rod made of stainless steel, connector	Z 602	28 541 6505
Magnetic stirrer rod set for standard applications AlNiCo5, circular cross-section, PTFE coated, 1 piece 15, 20, 30, 40, 50, 60, 70, 80 mm each	Z 603	28 541 6554
Magnetic stirrer ripe for medium volumes SmCo, circular cross-section, PTFE coated, 5 pieces 9 x 15 mm each	Z 604	28 541 6562
Magnetic stirrer for large volumes SmCo, elliptic cross-section, PTFE coated, 1 piece 19 x 75 mm each	Z 605	28 541 6579
Compressed-air connector (only SLK) for use in an aggressive environment (subsequent installation only by manufacturer)	Z 607	28 541 6595
Compressed-air connector set (only SLK)	Z 608	28 541 6608



www.schottinstruments.com

The first address for electrochemical measurements

SCHOTT Instruments

ProLab series Lab series Automatic recognition Stand Electrodes Buffers Home

The first address for pH meters

Crystal-clear benefits:

Electrodes and meters from SCHOTT Instruments

Electrodes from SCHOTT Instruments are used by professionals in laboratories all over the world. This comes as no surprise as we have been involved in the production and development of electrodes for almost 70 years: a vast know how in glass which our customers benefit directly from. What once began with the patent on the pH electrode, has grown into an extensive range of products including several hundred electrodes designed to meet standard- and special applications. From water to wine – we have the appropriate electrode to meet our customers' requirements.

Product innovations in detail!

SCHOTT Instruments present new products, further developments and and much more ...

up-to-date

and detailed

ProLab series

- ProLab 1000 - highlights
- ProLab 1000 - special details
- ProLab 2000 - highlights
- ProLab 2000 - special details
- ProLab 3000 - Highlights
- ProLab 3000 - Special details
- ProLab 3000 - Recorder features
- ProLab 4000 - Highlights

You would like to have more information?

Crystal clear Benefits: High-End-pH/ION Meter ProLab 3000

Performance black on white - technical data

Technical Data: Overview	
pH measurement	Range: -2.000 pH up to +20.000 pH -2.00 pH up to +20.00 pH -2.0 pH up to +20.0 pH Accuracy (±1 digit)
mV measurement	Range: -2200.0 mV up to +2200.0 mV -2200 mV up to +2200 mV Accuracy (±1 digit) in mV 2-channel-mV-measurement (galvanically separated)
ISE measurement	1.0E-40 ... 9.9E39 mg/l
Temperature Measurement	Measuring range -35.0°C up to + 150.0 °C Accuracy in °C (±1 digit)
Interfaces	USB (slave) and RS-232-C interface USB host interface Plug and play connection of USB hub, USB printer


SCHOTT Instruments

A further click displays all technical product details and ordering data:

Stored as PDF file

All essential information summarized on one sheet of paper – ready for printout.

Information around the clock.



Come and visit our download area!




Here you can find a lot of useful information, such as e.g.:

- ▶ Catalogues and product brochures
- ▶ Operating manuals
- ▶ Application reports and published articles

A further download service is the detailed section „Software“

- ▶ Just one look at the software introduction page will quickly indicate, whether a new software update is available for your instrument.
- ▶ All software updates can be downloaded for free at any time – the installation directions will explain all necessary steps for your update.



Innovative electrochemistry, innovative viscometry – from the very beginning.



Thorough quality control checking is an essential element from start to finish. The zero point and response time of **every** electrode is checked – as shown here with a batch of BlueLine laboratory electrodes.



Even today glass blowing talent is still indispensable.

The success story of electrochemical measurement began more than 70 years ago with the development of the glass electrode at SCHOTT.

It is hard to imagine – but in 1936 SCHOTT revolutionized the area of chemical measurements with a glass electrode that looks more like a glass bulb to our modern eyes. Glass electrodes made of the newly created, electrically

conductive pH glass were developed at the *Jena^{er} Glaswerk SCHOTT & Gen.*, thus making it possible to achieve “sufficient accuracy” for pH measurements “with conventional pointer-type galvanometers”.

In 1938 our first small brochure described how this was achieved. The development was based on the experience which we had gained from close co-operation with pioneer users in the industry.

Glass know-how was also the idea behind another measurement process that we pioneered in 1940: capillary viscometry. With precisely calibrated glass capillaries it was possible to determine the viscosity of Newtonian liquids more accurately than with any other method known at that time. You only need to measure the time required for a liquid sample to flow through a calibrated capillary with a defined constant. Then the required time was measured with a stop-watch. We have changed that.

The success story of our meters began in the seventies

In addition to our pH electrodes and viscometers, advances in the field of microelectronics in the seventies paved

the way for the development of our first measurement instruments.

To enable us to react faster and more flexible to our customer needs SCHOTT Geräte GmbH was established in 1973 as a separate company. Our newly developed electronic instruments – such as laboratory meters, pH meters or the automatic viscosity measurement system AVS – caused a sensation and rapidly conquered their respective markets.

Our first microprocessor-controlled piston burette and our titration systems were a small sensation in the industry. Our inexpensive, portable, pocket-size pH meters and conductivity meters were an instant success. In 1988, SCHOTT Geräte presented the first PC-supported titration system.

In recent years we have repeatedly demonstrated our remarkable power of innovation. With the SILAMID reference system, the multi-functional electrodes, the SMEK plug system, the complete BlueLine laboratory electrode program, the new Type A pH glass, the SteamLine electrodes for CIP and SIP applications, the ScienceLine range of products and the new Ioline product family with



The new measurement method had to be explained: in 1938 we published our first instructions for electrochemical pH measurement and potentiometric titration.



Our buffer solutions are hermetically sealed in double-pointed ampoules and sterilized with superheated steam. You can rest assured that you will always have a reliable buffer solution on hand for calibration.

the iodine/iodide reference system, we have provided new impulses for electrochemistry.

Our innovative AVS products have also made life much easier for our customers who need to measure viscosity. Examples include the practical ViscoClock, the AVSPro II, an apparatus for automatic viscosity measurements that is top-of-the-line worldwide, or the latest modular measurement systems AVS 370 and AVS 470. This is also true of our most recent developments in the area of titration equipment, where we have set new standards for top-of-the-line equipment with the TitroLine alpha *plus* titrator, the T110 *plus* piston burette and the TW alpha *plus* sample changer, the titration software TitriSoft 2.6 as well as the new coulometric KF titrator, TitroLine KF *trace*.



With the processed calibration, the viscometers are provided with an ID number and a certificate documenting the specific characteristics.

The success story continues – SCHOTT Geräte has become SCHOTT Instruments and is now even more international.

In October 2003 SCHOTT Geräte GmbH transferred to SCHOTT Instruments GmbH and became part of the Nova Analytics group of companies. This was a logical step for us – because although our historical roots lie in the area of glass manufacturing, our core competence has long since been measurement technology for the laboratory and production industry. Our relationship with the affiliated companies within the Nova Analytics Group will result in synergy effects that will benefit our customers and distributors worldwide. In particular, it will be possible to improve service and logistics in America and Asia even further.



Everything O.K. A batch of TITRONIC® basic piston burettes after final inspection

More than 70 years of research and development and a long-standing tradition

The list of our innovations is long: today our electrodes are smaller, more precise, faster and more stable; our measurement equipment offers incomparably higher performance. Over the years the electrochemical measurement methods and the viscometry that we initiated have established themselves as problem-free and reliable methods throughout the world, and they have become indispensable for an incredibly varied range of applications. Nevertheless, since our pioneering days, one thing has invariably remained the same – our tradition of working very closely together with those who use our products in order to create something new. We intend to remain true to this tradition in the future.



A customer satisfaction center. In our application laboratory, for example, new methods can be developed for our customers or the suitability of existing methods for new applications can be tested.

Just a brief excerpt from our company's history

Since 1936 – new impulses from research and development time and time again

- 1936** Beginning of development and production of pH glass electrodes at *Jena^{er} Glaswerk SCHOTT & Gen.* in Jena.
- 1940** Beginning of viscometer production using capillaries that were manufactured in accordance with the calibrated precision glass method that SCHOTT had developed.
- 1952** Development and production of the first gel-filled, low-maintenance reference electrodes.
- 1962** The unique platinum diaphragm makes substantially faster response times possible, among other things.
- 1964** Double electrolyte system for reference electrodes.
- 1970** Introduction of semiconductor preamplifiers for pH measurement technology.
- 1972** Buffer solutions in double-pointed ampules sterilized with superheated steam guarantee reliable calibration – even after several years in storage.
- 1972** Plug system from SCHOTT, copied time and again.
- 1973** SCHOTT Geräte GmbH established as an independent company.
- 1973** Beginning of viscometer calibration using PTB tested reference measurement standards. (German Physical Technical Institute).
- 1974** Development and production of electronic laboratory pH meters.
- 1975** Market launch of the first automatic viscosity measurement apparatus for aggressive and corrosive solvents (AVS/G and AVS/PA).
- 1977** Development and production of portable electronic pH meters.
- 1978** Production start for the first titration control unit TR 155 and the T 100 piston burette with exchange unit.

- 1982** The first microprocessor-controlled viscosity measurement apparatus (AVS 300).
- 1983** Development of the new Type S pH glass for hot alkaline solutions with extraordinarily high reliability and useful life, and Type H pH glass, robust and small alkali error.
- 1984** Combination measurement and reference pH electrode with integrated Pt 1000 as standard.
 - SCHOTT Geräte presents the first thermal scanning method for viscosity measurement.
 - The first stand-alone viscosity measurement apparatus with integrated computing function (AVS 400 and AVS 440) are introduced on the market.
 - Compact T 80/T 90 piston burettes and simple control unit TR 85.
- 1988** Presentation of the first PC-controlled titration system TPC 2000 at the Achema 1988.
- 1989** With the AVS 500, the tradition of successful automatic samplers for determination of the viscosity of aggressive polymer solutions was continued.
- 1990** REFERID[®] electrodes with polymer electrolyte, low-maintenance.
- 1991** Low-impedance Type L pH glass for low temperatures and ultrapure water.
 - Automatic sampler TW 280.
- 1992** TT electrodes, capable of withstanding up to -60 °C.
 - T 200 and T 110 piston burettes and universal titration control unit TC 1200.
- 1993** Combination pH electrodes with temperature sensor and plastic shaft.
- 1994** Compact TitroLine alpha titrator.

View over Mainz with premises of
SCHOTT AG/SCHOTT Instruments

Source of image: SCHOTT AG

1995 SILAMID®, potential-stable reference system.

- First Windows titration software TitrSoft 1.0 (WIN 3.1).

1996 New SMEK 6-pin plug system, shielded.

1997 New BlueLine range of laboratory electrodes and VP plug system.

- Electrodes with certified pressure and temperature range.
- Market launch of the Visco-Clock for capillary viscosity measurement.

1998 Development of TitrSoft 2.0 software (as of WIN 95).

1999 Range of industrial electrodes up to 10 bar and 135 °C, SMEK plug system in IP 68 version.

- New Type A pH glass, rapid reaction in drinking water.
- Market launch of the fully automatic AVSPro viscosity measurement system for high sample throughput.

2000 Introduction of a completely new series of compact, simple piston burettes and titrators: TITRONIC® *basic*, TITRONIC® *universal* and TitroLine *easy*.

- Introduction of the Karl Fischer titration system TitroLine *KF*.

2001 Development and production of SteamLine process electrodes for CIP and SIP applications in the pharmaceutical, food and chemical sectors.

2002 Sales launch of newly developed “*plus*” product line: TitroLine *alpha plus*, T 110 *plus*.

- Introduction of TW *alpha plus* sample changer.
- Market launch of TitrSoft 2.5 software.

2003 The compact and highly flexible AVS 370 viscosity measurement system is presented to the market.

- Change of company name to SCHOTT Instruments GmbH, Mainz, integration into the internationally active Nova Analytics Group.

2004 Amalgamation and further development of the laboratory electrode product range for the most exacting requirements in the “ScienceLine” product line.

- The new generation of automatic viscosity measurement systems is rounded off with the AVS 470.

2005 The Lab meters family is introduced:

A wireless sensor recognition guarantees the optimum interplay between electrode and meter.

2006 Introduction of the ProLab family of instruments: Multi-functional measuring instruments with integrated user recognition guarantee utmost flexibility and reliability of measurement.

2007 ProLab 3000 and 4000 high-end laboratory instruments signify the cutting-edge standard for pH/ionic and conductivity measuring and for the first time combine highest measuring quality with functionality, also providing a user-friendly navigation comparable to a Windows PC.

- The new Karl Fischer titrator, TitroLine *KF trace* from SCHOTT Instruments, also offers a coulometric technique for determining even smallest water content.

2008 The new IoLine electrodes with their patented iodine/iodid three-chambers reference system represent the perfect solution for accomplishing the ultimate challenging measuring tasks in i.e. pharmacy, biotechnology and food industry.






*Precision measurement
in a complete system*

SCHOTT Instruments

SCHOTT Instruments GmbH
Hattenbergstraße 10 · 55122 Mainz · Germany
Phone +49 (0)6131665111 · Fax +49 (0)6131665001
support@schottinstruments.com

A Nova Analytics company  NOVA
ANALYTICS