

## Thermo Scientific Value Laboratory Refrigerator and Freezers

Thermo Scientific Value laboratory refrigerators, freezers and combination units provide an economical solution for safe, sample storage.

Our range of models are particularly suited for routine laboratory applications with sizes from small benchtop to full-height freestanding. Our features are not available in domestic models and meet the demanding requirements of the laboratory environment.

Not available outside North America.



**3.2 cu. ft. Value Freezer**  
Suitable for -25°C Applications  
Model: 3753

- Factory tested prior to shipment
- NRTL listed
- Environmentally-friendly, non CFC refrigerants
- Adjustable temperature control
- Easy-to-clean, durable finish
- Manual and auto defrost combination units
- Manual defrost refrigerators and freezers

**Warranty:** 13 months parts and labor

**Certifications:** CSA and or UL

**Alert:** These units are intended for basic laboratory use only. Storage of high-value, life science materials or clinical biologicals may require alarm, monitoring, back-up and temperature performance features. The need to add these features should be determined by your industrial hygienist or safety officer. These units are not suitable for corrosive environments.



**10 cu. ft. Value Combination Unit**  
Model: 3764A



**20.6 cu. ft. Value Freezer**  
Model: 3767A



**24.6 cu. ft. Value Combination Unit**  
Model: RCRF252A

Thermo Scientific Value Laboratory Refrigerators and Freezers (Not Available outside North America)

Refrigerators (Temperature Range 1° to 12°C)

Model No.	Capacity cu. ft.	Defrost	Voltage	Shelves	Key Lock	Interior Dimensions D x W x H (inches)	Exterior Dimensions D x W x H (inches)	Shipping Weight (lbs.)
3756	1.8	manual	115V / 60Hz	1 fixed, 2 door	yes	13 x 16.5 x 15.75**	21 x 19.75 x 20.75	55
RLR065A14*	5.5	manual	115V / 60Hz	3 fixed, 2 door	yes	15 x 20 x 26.5	23.5 x 23.75 x 33.75	110
3751-DB	5.6	manual	120V / 60Hz	4 adjustable	yes	16 x 20 x 27.5***	25.13 x 23.65 x 33.5	102
3754*	6.1	manual	115V / 60Hz	2 fixed, 3 door	no	15 x 21 x 26.5†	24 x 24 x 34	120
3766A	20.3	manual	115V / 60Hz	5 fixed	yes	20 x 27 x 57††	32 x 30.75 x 70	240
RLR4304A	43	manual	120V / 60Hz†††	8 adjustable	no	26.5 x 49 x 57 (fan is at 53")	32 x 53 x 76	645
RRC4304A	43	manual	120V / 60Hz†††	8 adjustable	no	26.5 x 49 x 57 (fan is at 53")	32 x 53 x 76	645

\* Features digital temperature display  
 \*\* 8" (for bottom 8"), 13" (for top 8")  
 \*\*\* 10" (for bottom 8"), 16" (for top 19.5")  
 † 8" (for bottom 5"), 15" (for top 21.5")  
 †† 17" (for bottom 20"), 20" (for top 40")  
 ††† NEMA 5-20P Plug

Freezers (Temperature Range -12° to -20°C)

Model No.	Capacity	Defrost	Voltage	Shelves	Key Lock	Interior Dimensions D x W x H (inches)	Exterior Dimensions D x W x H (inches)	Shipping Weight (lbs.)
3757	1.6	manual	115V / 60Hz	1	yes	13.5 x 15.5 x 14*	21 x 19.75 x 20.75	55
3753**	3.2	manual	120V / 60Hz	3 drawers	yes	16.5 x 17.75 x 24.25	23.5 x 23.5 x 33.5	102
3752-DB	5	manual	115V / 60Hz	2 fixed, 3 door	no	15 x 18 x 28.5	23.5 x 21 x 33.5	84
3767A	20.6	manual	115V / 60Hz	4 fixed	yes	20 x 27 x 57	31 x 32 x 71	200

\* 7" (for bottom 8"), 13.5" (for top 6")  
 \*\* -15 to -35°C, suitable for -25°C applications

Combination Refrigerator/Freezers (Temperature Range 1° to 12°C Ref., -12° to -20°C Frz.)\*

Model No.	Capacity cu. ft.	Defrost	Voltage	Shelves Interior/Door	Drawers	Interior Dimensions D x W x H (inches)	Exterior Dimensions D x W x H (inches)	Shipping Weight (lbs.)
3750	5.6	manual	120V / 60Hz	3 / 3	-	18 x 15 x 21 (ref.) 11 x 15 x 5 (frz.)	25.5 x 21.25 x 33.5	78
3771	8.8	automatic	115V / 60Hz	2 / 2	1	13 x 19.25 x 33 (ref.) 8.75 x 18.5 x 14.5 (frz.)	23.5 x 23.25 x 60.75	160
3764A	10	manual	120V / 60Hz	Ref: 3 / 4 Frz: 1	-	13.5 x 19.5 x 38.5 (ref.) 15 x 17.5 x 11.5 (frz.)	26.5 x 23.5 x 60	150
3772†	12	manual (frz.)	115V / 60Hz	Ref: 4	Frz: 3	13.5 x 19 x 35 (ref.) 13.5 x 17 x 28 (frz.)	23.75 x 23.75 x 79	251
RRCF192A	17.6	automatic	115V / 60Hz	2	2	17.5 x 25 x 312 (ref.) 15 x 23 x 15.25 (frz.)	32 x 28 x 66.5	225
3763A	18.2	automatic	120V / 60Hz	Ref: 4 half, 1 full Frz: 1 / 2	3	18 x 26 x 38 (ref.) 15 x 23 x 17.25 (frz.)	30.5 x 30 x 65.25	240
RRCF252A**	24.6	automatic	120V / 60Hz	Ref: 4 adj. half 1 full / 6 door Frz: 1 / 3	3	20 x 32 x 37 (ref.) 17 x 31 x 18.5 (frz.)	36 x 35.75 x 68	242

\* Key locks are not available on any combination units  
 \*\* Features digital temperature display  
 † Dual compressor system, automatic refrigerator and manual defrost freezer



5 cu. ft. Value Freezer  
Model: 3752-DB



5.6 cu. ft. Value Refrigerator  
Model: 3751-DB



© 2011 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

North America: USA/Canada +1 866 984 3766 (866-9-THERMO) [www.thermoscientific.com/cold](http://www.thermoscientific.com/cold)  
 Europe: Austria +43 1 801 40 0, Belgium +32 53 73 42 41, France +33 2 2803 2180, Germany national toll free 08001-536 376,  
 Germany international +49 6184 90 6940, Italy +39 02 95 05 93 54, Netherlands +31 76 579 55 55, Nordic/Baltic/CIS countries +358 9 329 10200,  
 Russia +7 812 703 42 15, Spain/Portugal +34 93 223 09 18, Switzerland +41 44 454 12 12, UK/Ireland +44 870 609 9203  
 Asia: Australia +61 45 453 9220, China +86 21 6865 4588 or +86 10 8419 3588, India toll free 1800 22 8374, India +91 22 6716 2200,  
 Japan +81 45 453 9220, New Zealand +64 9 980 6700, Other Asian countries +852 2885 4613 Countries not listed: +49 6184 90 6940  
 DSCSVALUE 0311

