

Catalog Number	Description	Ship Wt.	List Price
<b>4 LITRE BLENDERS</b>			
CB15E	4.0 Litre 3-Speed Commercial Blender with 3.75 Peak HP, plus Pulse, Electronic Membrane Keypad, Stainless Steel Container with 2 Handles & Rubberized Lid. (CE Approved; RoHS)	38.9 lbs.	\$2,485.00
CB15EB	Base Only. Same as CB15E, except with no Container or Lid. CE Approved. RoHS.	32 lbs.	\$1,682.00
LBC15E	Same as CB15E, except with Stainless Steel Lid. CE Approved. RoHS.	39 lbs.	\$2,772.00
CB15TE	4.0 Litre 3-Speed Commercial Blender with 3 minute Timer, 3.75 Peak HP, plus Pulse, Electronic Membrane Keypad, Stainless Steel Container with 2 Handles & Rubberized Lid. (CE Approved; RoHS)	38.9 lbs.	\$2,666.00
CB15VE	4.0 Litre Variable-Speed Commercial Blender with 3.75 Peak HP, Adjustable Speeds plus Pulse, Electronic Membrane Keypad, Stainless Steel Container with 2 Handles & Rubberized Lid. (CE Approved; RoHS)	38.9 lbs.	\$2,858.00
CB15VXE	4.0 Litre Variable-Speed Commercial Blender with 3.75 Peak HP, Adjustable Speeds, Built-in RPM Display, Electronic Membrane Keypad, Stainless Steel Container with 2 Handles & Rubberized Lid. (CE Approved; RoHS)	38.9 lbs.	\$2,911.00
CB15VXEB	Base Only. Same as CB15VXE, except with no Container or Lid. CE Approved. RoHS.	32 lbs.	\$1,969.00
<b>ACCESSORIES FOR 4 LITRE BLENDERS</b>			
CAC70	2 Handle 4.0 Litre Stainless Steel Container with Blade Assembly & Stainless Steel Lid	9 lbs.	\$1,285.00
CAC72	2 Handle 4.0 Litre Stainless Steel Container with Blade Assembly & Rubberized Lid	8 lbs.	\$1,216.00
CBL10	Polycarbonate Hot Blending Lid for 4.0 Litre Blender Polycarbonate Hot Blending Lid for 4.0 Litre Blender Polycarbonate Hot Blending Lid for 4.0 Litre Blender	8 lbs.	\$162.00
AD1	Adapter, 4.0 Litre, allows use of 1.0 Litre and smaller containers on 4.0 Litre Blenders	3.1 lbs.	\$399.00
2610C	4 Litre Cool Base Stainless Steel Container	7 lbs.	\$2,719.00
2610F	4 Litre Stainless Steel Foam Arrestor	4 lbs.	\$597.00
2610T	4 Litre Stainless Steel Temperature Control Coil	4 lbs.	\$1,088.00
<b>1.5 LITRE BLENDERS - TBB SERIES</b>			
TBB145E	1.5 Litre Commercial Blender Complete with BPA-free Copolyester Stackable Container, Toggle Switch Controls, 2 HP, 1500W 2-speed motor, 24,000 max RPM (CE Approved )	11.6 lbs.	\$660.00
TBB145EB	Base Only. Same as TBB145E, except with no Container or Lid. CE Approved.	9 lbs.	\$396.00
TBB160E	1.5 Litre Commercial Blender Complete with BPA-free Copolyester Stackable Container, Electronic Touchpad Controls with 60-second Countdown Timer, 2 HP, 1500W 2-speed motor, 24,000 max RPM (CE Approved )	11.6 lbs.	\$712.00
TBB160EB	Base Only. Same as TBB160E, except with no Container or Lid. CE Approved.	9 lbs.	\$416.00

Catalog Number	Description	Ship Wt.	List Price
TBB175E	1.5 Litre Commercial Blender Complete with BPA-free Copolyester Stackable Container, Electronic Touchpad Controls with Variable Speed Control Dial, 2 HP, 1500W 2-speed motor, 24,000 max RPM (CE Approved ) <b>Available end of 1st Quarter 2018 1.5 Litre Commercial Blender Complete with BPA-free Copolyester Stackable Container, Electronic Touchpad Controls with Variable Speed Control Dial, 2 HP, 1500W 2-speed motor, 24,000 max RPM (CE Approved ) Available end of 1st Quarter 2018</b>	11.6 lbs.	\$770.00
TBB175EB	Base Only. Same as TBB175E, except with no Container or Lid. CE Approved.		\$437.00
<b>ACCESSORIES FOR TBB SERIES 1.5 LITRE BLENDERS</b>			
CAC139	1.5 Litre BPA-free Stackable Copolyester Container Complete with Blade Assembly & Lid (lids not shown) (for use on TBB145, TBB160 & TBB175)		\$229.00
CAC140	Blending Assembly (for use on TBB145, TBB160 & TBB175)		\$161.00
CAC145	Lid with Measuring Cup (for use on TBB145, TBB160 & TBB175)		\$31.00
CAC146	Jar Pad 035091 (for use on TBB145, TBB160 & TBB175)		\$31.00
CAC147	Clutch Repair Kit (12 pieces of 035251) (for use on TBB145, TBB160 & TBB175)		\$380.00
CAC152	1.4 Litre Stainless Steel Container complete with Blade Assembly & Lid (for use on TBB145, TBB160 & TBB175)	4 lbs.	\$357.00

Catalog Number	Description	Ship Wt.	List Price
<b>2 LITRE BLENDERS - MX SERIES</b>			
MX1000XTXEE	2.0 Litre Commercial Blender Complete with BPA-free Copolyester Container. Heavy Duty 3.5 Peak HP High Speed Motor, 2 Speed Paddle Switch with Pulse, Replaceable Blending Assembly.	17 lbs.	\$970.00
MX1000XTEES	Same as MX1000XTXEE except with 2.0 Litre Stainless Container with Blending Assembly & Lid	19 lbs.	\$1,101.00
MX1050XTXEE	2.0 Litre Commercial Blender Complete with BPA-free Copolyester Container. Heavy Duty 3.5 Peak HP High Speed Motor, 2 Speed Electronic Control with Pulse, Replaceable Blending Assembly.	17 lbs.	\$992.00
MX1100XTXEE	2.0 Litre Commercial Blender Complete with BPA-free Copolyester Container. Heavy Duty 3.5 Peak HP High Speed Motor, 2 Speed Electronic Control with Pulse and Timer. Replaceable Blending Assembly.	17 lbs.	\$1,007.00
MX1100XTEES	Same as MX1100XTXEE except with 2.0 Litre Stainless Container with Blending Assembly & Lid	19 lbs.	\$1,154.00
MX1200XTXEE	2.0 Litre Variable-Speed Commercial Blender Complete with with BPA-free Stackable Copolyester Container, Heavy Duty 3.5 Peak HP High Speed Motor with Pulse. Replaceable Blending Assembly.	17 lbs.	\$1,026.00
MX1200XTEES	Same as MX1200XTXEE except with 2.0 Litre Stainless Container with Blending Assembly & Lid	19 lbs.	\$1,227.00
MX1200XTXSEE	2.0 Litre Variable-Speed Commercial Blender Complete with with BPA-free Stackable Copolyester Container, Heavy Duty 3.5 Peak HP High Speed Motor with Pulse, Large Sound Enclosure included. Replaceable Blending Assembly.	23 lbs.	\$1,167.00
MX1300XTXEE	2.0 Litre Programmable (4 Programs) Commercial Blender Complete with with BPA-free Stackable Copolyester Container, Electronic Control, Heavy Duty 3.5 Peak HP High Speed Motor. Replaceable Blending Assembly.	17 lbs.	\$1,167.00
MX1300XTEES	Same as MX1300XTXEE except with 2.0 Litre Stainless Container with Blending Assembly & Lid	19 lbs.	\$1,251.00
MX1500XTXSEE	2.0 Litre Programmable (4 Programs) Beverage Station Complete with Heavy Duty 3.5 Peak HP High Speed Motor, Electronic Control. BPA-free Stackable Copolyester Container. Large Sound Enclosure included. Replaceable Blending Assembly.	23 lbs.	\$1,312.00
<b>ACCESSORIES FOR MX SERIES 2 LITRE BLENDERS</b>			
CAC 90I	2.0 Litre Stainless Steel Container Complete with Blade Assembly & Lid (for use on MX blenders)	6 lbs.	\$491.00
CAC 95I	2.0 Litre BPA-free Copolyester Container Complete with Replaceable Blending Assembly & Lid (for use on MX blenders)	3.5 lbs.	\$257.00
SE1000	Large Sound Enclosure for MX series 2.0 Litre blenders.	5.5 lbs.	\$331.00
AD3	Adapter allows use of 1.0 Litre and smaller containers on 2.0 Litre MX Series Blenders	2.3 lbs.	\$451.00
SE500	Small Sound Enclosure for MX series 2.0 Litre blenders.	4 lbs.	\$380.00

Catalog Number	Description	Ship Wt.	List Price
<b>2 LITRE BLENDERS - HGB SERIES</b>			
HGB500	2.0 Litre Commercial Blender Complete with BPA-free Copolyester Container and Stainless Steel Housing, (NON-CE) <b>DISCONTINUED</b>	13.8 lbs.	
HGB50E	2.0 Litre Commercial Blender Complete with BPA-free Copolyester Container and Stainless Steel Housing, (CE Approved) <b>DISCONTINUED</b>	13.8 lbs.	
HGB550	2.0 Litre Commercial Blender Complete with Stainless Steel Container and Stainless Steel Housing, (NON-CE) <b>DISCONTINUED</b>	14.4 lbs.	
HGB55E	2.0 Litre Commercial Blender Complete with Stainless Steel Container and Stainless Steel Housing, (CE Approved) <b>DISCONTINUED</b>	14.4 lbs.	
<b>ACCESSORIES FOR HGB SERIES 2 LITRE BLENDERS</b>			
CAC15	Hot Blending Lid for 2.0 Litre Blenders	1.5 lbs.	\$587.00
CAC21	1.5 Litre BPA-free Copolyester Container Complete with Blade Assembly & Lid (for use on HGB500/HGB50E & HGB550/HGB55E)	2.4 lbs.	\$124.00
CAC31	2.0 Litre Stainless Steel Container Complete with Blade Assembly & Lid (for use on HGB500/HGB50E & HGB550/HGB55E)	5 lbs.	\$346.00
CAC59	2.0 Litre BPA-free Copolyester Container Complete with Blade Assembly & Lid (for use on HGB500/HGB50E & HGB550/HGB55E )	3 lbs.	\$217.00
SS810	2.0 Litre Stainless Steel Container Complete with Blade Assembly & Stainless Steel Lid (for use on HGB500/HGB50E & HGB550/HGB55E)	6 lbs.	\$430.00
SS715	1.0 Litre Stainless Steel Container Complete with Blade Assembly & 2 piece Vinyl Lid (for direct use on 2 Litre Blenders)	2.2 lbs.	\$215.00
AD2	Adapter, 2.0 Litre, allows use of 1.0 Litre and smaller containers on 2.0 Litre HGB Series Blenders	2.3 lbs.	\$418.00
<b>1 LITRE BLENDERS .</b>			
800S	Single Speed Blender with 1.0 Litre Stainless Steel Container and Baked Enamel Base (NON-CE)	8.4 lbs.	\$484.00
800ES	Single Speed Blender with 1.0 Litre Stainless Steel Container and Baked Enamel Base (CE Approved)	8.4 lbs.	\$518.00
800G	Single Speed Blender with 1.25 Litre Glass Container and Baked Enamel Base (NON-CE)	10.1 lbs.	\$391.00
800EG	Single Speed Blender with 1.25 Litre Glass Container and Baked Enamel Base (CE Approved)	10.1 lbs.	\$425.00
800BU	Single Speed Blender Base Only (NON-CE)	6.7 lbs.	\$323.00
800EBU	Single Speed Blender Base Only (CE Approved)	6.7 lbs.	\$346.00
8010S	2 Speed Blender with Timer, 1.0 Litre Stainless Steel Container and Baked Enamel Base (Non-CE)	9.4 lbs.	\$745.00
8010ES	2 Speed Blender with Timer, 1.0 Litre Stainless Steel Container and Baked Enamel Base (CE Approved)	9.4 lbs.	\$815.00
8010G	2 Speed Blender with Timer, 1.25 Litre Glass Container and Baked Enamel Base (NON-CE)	8.0 lbs.	\$654.00
8010EG	2 Speed Blender with Timer, 1.25 Litre Glass Container and Baked Enamel Base (CE Approved)	11.4 lbs.	\$726.00
8010BU	2 Speed Blender with Timer, Base Only (Non-CE)	8.0 lbs.	\$601.00
8010EB	2 Speed Blender with Timer, Base Only (CE App'd.)	8.0 lbs.	\$659.00
8011S	2 Speed Blender with 1.0 Litre Stainless Steel Container and Baked Enamel Base (NON-CE)	9.4 lbs.	\$686.00
8011ES	2 Speed Blender with 1.0 Litre Stainless Steel Container and Baked Enamel Base (CE Approved)	9.4 lbs.	\$726.00

Catalog Number	Description	Ship Wt.	List Price
8011G	2 Speed Blender with 1.25 Litre Glass Container and Baked Enamel Base (NON-CE)	11.4 lbs.	\$596.00
8011EG	2 Speed Blender with 1.25 Litre Glass Container and Baked Enamel Base (CE Approved)	11.4 lbs.	\$642.00
8011BU	2 Speed Blender, Base Only (NON-CE)	8.0 lbs.	\$548.00
8011EB	2 Speed Blender, Base Only (CE App'd.)	8.0 lbs.	\$593.00
8009L	2 Speed Blender with 1.25 Litre Stainless Steel Container and Baked Enamel Base with Jack for Digital Tachometer (NON-CE)	11.4 lbs.	\$946.00
LB20EG	Variable Speed Blender, Feedback, Maintains Torque between 500-20,000 RPM, 1.25 Litre Glass Container (CE Approved)	11.4 lbs.	\$991.00
LB20ES	Variable Speed Blender, Feedback, Maintains Torque between 500-20,000 RPM, 1.0 Litre Stainless Steel Container (CE Approved)	9.4 lbs.	\$1,083.00

Catalog Number	Description	Ship Wt.	List Price
<b>ACCESSORIES FOR 1 LITRE BLENDERS (also fit AD1, AD2, AD3 Adapters as noted)</b>			
CAC19	1.25 Litre BPA-free Copolyester Container Complete with Blade Assembly & Lid (for use on 800, 8010, 8011, AD1, AD2 & AD3 Adapters)	1.8 lbs.	\$138.00
CAC32	1.25 Litre Glass Container Complete with Blade Assembly & Lid (for use on 800, 8010, 8011, AD1, AD2 & AD3 Adapters)	2.2 lbs.	\$107.00
CAC33	1.0 Litre Stainless Steel Container Complete with Blade Assembly & Lid (for use on 800, 8010, 8011, AD1, AD2 & AD3 Adapters)	2.2 lbs.	\$174.00
CAC64	500ml. Polypropylene Container	1.1 lbs.	\$152.00
SS510A	Stainless Steel Cover with Aspirator for SS510 Series Container	1.0 lb.	\$265.00
SS510C	1.0 Litre Cool Base Container	3.0 lbs.	\$787.00
SS510F	1.0 Litre Foam Arrestor for SS510 Series Container	1.0 lb.	\$437.00
SS510S	1.0 Litre Stainless Steel Container with Sealed Cover	4.0 lbs.	\$1,509.00
SS510T	1.0 Litre Stainless Steel Temperature Control Coil for SS510 Series Container	1.0 lb.	\$1,293.00
SS515	1.0 Litre Stainless Steel Container with Stainless Steel Lid & Teflon Gaskets	4.0 lbs.	\$792.00
SS610	1.0 Litre Stainless Steel Container with Stainless Steel Lid	4.0 lbs.	\$243.00
SS610T	1.0 Litre Stainless Steel Temperature Control Coil for SS610 Series Container	1.0 lb.	\$774.00
SS610F	1.0 Litre Stainless Steel Foam Arrestor for SS610 Series Container	1.0 lb.	\$372.00
<b>MINI CONTAINERS (fit 1 litre units &amp; AD1, AD2, AD3 Adapters)</b>			
MC-1	Mini Container: 12-37 ml Capacity	1.0 lb.	\$287.00
MC-2	Mini Container: 37-110 ml Capacity	1.0 lb.	\$291.00
MC-3	Mini Container: 50-250- ml Capacity	1.0 lb.	\$331.00
<b>DRY BLENDING CONTAINERS (fit 1 litre units &amp; AD1, AD2, AD3 Adapters)</b>			
SS110	Stainless Steel Container/Plastisol lid: 10-75 Grams Dry Capacity	1.0 lb.	\$397.00
SS115	500 ml Capacity Stainless Steel Container with Lid	4.0 lb.	\$614.00
<b>HAZARDOUS LOCATION BLENDER &amp; ACCESSORIES</b>			
BHL240*	Hazardous Location Blender Base Only, 1.0 Litre Capacity (NON-CE) *Price does not include explosion-proof switch or special container required to retain Class 1, Group D hazardous location rating. These are available at additional cost. (See below.)	19 lbs.	\$7,175.00
BHLSW	Hazardous Location Switch for BHL240 Blender	4 lbs.	\$906.00
HL515	1.0 Litre Stainless Steel Container for use with BHL240 Blender	4 lbs.	\$723.00
<b>LAB MIXERS</b>			
WDM120E	Single Spindle 3 Speed Lab Mixer with Cup and 2 Agitators	15.5 lbs.	\$587.00
CAC20	Stainless Steel Cup (0.83 Litre) (for Lab Mixers)	0.8 lbs.	\$64.00
<b>ICE CRUSHER</b>			
IC20	Lab Ice Crusher with Funnel and Euro cord (GS/CE Mark)	24 lbs.	\$928.00
IC70E	Brushed Stainless Steel Commercial Ice Crusher, 2.4 Litre Bowl, 230V, 50Hz, CE & RoHS	8 lbs.	\$266.00
<b>CHOPPER / GRINDERS</b>			

Catalog Number	Description	Ship Wt.	List Price
WCG75E	Chopper Grinder - 0.75 Litre Capacity. 2 Separate Bowl Assemblies - 1 for Chopping, 1 for Grinding (CE Approved)Chopper Grinder - 0.75 Litre Capacity. 2 Separate Bowl Assemblies - 1 for Chopping, 1 for Grinding (CE Approved)Chopper Grinder - 0.75 Litre Capacity. 2 Separate Bowl Assemblies - 1 for Chopping, 1 for Grinding (CE Approved)	11.5 lbs.	\$311.00
WSG30E	Laboratory Particle Grinder with 3 grinding bowls, each with storage lid	4 lbs.	\$320.00
CAC103	Grind Bowl with lid - for WSG-30E	1 lb.	\$85.00
WSG60E	Laboratory Particle Grinder with 2 grinding bowls, each with storage lid. 1HP motor, 3 cup capacity	9 lbs.	\$361.00
CAC128	Grind Bowl with lid - for WSG- 60E	2 lbs.	\$111.00
<b>LABORATORY MICRONIZER</b>			
WFP14SE	14 cup Heavy Duty Laboratory Micronizer with LiquiLock™ Seal System, 3/4 Peak HP, Clear Bowl, Three SS Blades, 230V, 50/60Hz, CE & RoHS	20.25 lbs.	\$821.00
WFP16SE	16 Cup/3.75 Litre Heavy Duty Laboratory Micronizer with LiquiLock™ Seal System, 2 Peak HP, Clear Bowl, Sealed S Blade, Adjustable Slicing Disc, Reversible Shredding Disc, Sealed Whipping Disk, 230V, 50/60Hz, CE & RoHS	29.95 lbs.	\$981.00
<b>IMMERSION BLENDERS</b>			
WSB25XE	Quick Stik™ Immersion Blender Xpress - Cordless, with 9.6V Battery Pack and Wall Mount Recharging Station, 10"/25 cm Fixed Shaft Quick Stik™ Immersion Blender Xpress - Cordless, with 9.6V Battery Pack and Wall Mount Recharging Station, 10"/25 cm Fixed Shaft Quick Stik™ Immersion Blender Xpress - Cordless, with 9.6V Battery Pack and Wall Mount Recharging Station, 10"/25 cm Fixed Shaft	9 lbs.	\$631.00

Catalog Number	Description	Ship Wt.	List Price
WSB38XE	Immersion Blender - Cordless, 5-speed with 9.6V Battery Pack and Wall Mount Recharging Station, 7"/18 cm Removable Shaft	5 lbs.	\$541.00
WSB38XST	Replacement Shaft for WSB38XE Immersion Blender - 7"/18 cm	1 lb.	\$109.00
WSB38XBP	Replacement Battery Pack for WSB38XE Immersion Blender	1 lb.	\$159.00
WSB38XCS	Replacement Charger for WSB38XE Immersion Blender	2 lb.	\$79.00
WSB33XE	Quick Stik™ Immersion Blender - Light Duty, Fixed Shaft, 2 Speed Quick Stik™ Immersion Blender - Light Duty, Fixed Shaft, 2 Speed	14.2 lbs.	\$88.00
WSB40E	Medium Duty Quick Stik™ Immersion Blender with 10"/25 cm Fixed Shaft, 22 litre capacity, 350 watt, 2-speed, 1/2 HP heavy duty motor, Rubberized comfort grip , All purpose stainless steel blade.	5.5 lbs.	\$343.00
<b>IMMERSION BLENDERS cont'd.</b>			
WSB50CE	Immersion Blender with Heavy Duty, 12"/30.5 cm Removable Shaft, 37.9 litre capacity, 1 HP heavy duty motor - 18,000 RPM on Hi, 750 Watt variable speed motor - 9 settings, Continuous ON feature, Rubberized comfort grip and second handle, for safe and controlled operation, All purpose stainless steel blade, Completely sealed stainless steel shaft is easily removable and dishwasher-safe, 12 foot (365 cm) cord with 2 prong plug	21.3 lbs.	\$596.00
WSB55CE	Immersion Blender with Heavy Duty, Removable Shaft ( <i>Same as WSB50CE, except:</i> ) 14" / 35.6 cm removable shaft, 56.8 litre capacity Immersion Blender with Heavy Duty, Removable Shaft ( <i>Same as WSB50CE, except:</i> ) 14" / 35.6 cm removable shaft, 56.8 litre capacity	22.2 lbs.	\$619.00
WSB60CE	Immersion Blender with Heavy Duty, Removable Shaft ( <i>Same as WSB50E, except:</i> ) 16" / 40.6 cm removable shaft, 94.6 litre capacity Immersion Blender with Heavy Duty, Removable Shaft ( <i>Same as WSB50E, except:</i> ) 16" / 40.6 cm removable shaft, 94.6 litre capacity	22.8 lbs.	\$644.00
WSB65CE	Immersion Blender with Heavy Duty, Removable Shaft ( <i>Same as WSB50CE, except:</i> ) 18" / 45.7 cm removable shaft, 132.5 litre capacity Immersion Blender with Heavy Duty, Removable Shaft ( <i>Same as WSB50CE, except:</i> ) 18" / 45.7 cm removable shaft, 132.5 litre capacity	24.1 lbs.	\$671.00
WSB70CE	Immersion Blender with Heavy Duty, Removable Shaft ( <i>Same as WSB50CE, except:</i> ) 21" / 53.4 cm removable shaft, 189.3 litre capacity Immersion Blender with Heavy Duty, Removable Shaft ( <i>Same as WSB50CE, except:</i> ) 21" / 53.4 cm removable shaft, 189.3 litre capacity	25.2 lbs.	\$691.00
<b>Accessory SHAFTS for WSB Series Heavy Duty Immersion Blenders</b>			
WSB50ST	12" / 30.5 cm Stainless Steel Removable Shaft, 37.9 litre capacity, Interchangeable Removable Shaft for WSB Series Heavy Duty Immersion Blenders	16.3 lbs.	\$208.00
WSB55ST	(Same as WSB50ST, except:) 14" / 35.6 cm removable shaft, 56.8 litre capacity	19.7 lbs.	\$219.00
WSB60ST	(Same as WSB50ST, except:) 16" / 40.6 cm removable shaft, 94.6 litre capacity	14.5 lbs.	\$230.00
WSB65ST	(Same as WSB50ST, except:) 18" / 45.7 cm removable shaft, 132.5 litre capacity	22 lbs.	\$246.00
WSB70ST	(Same as WSB50ST, except:) 21" / 53.4 cm removable shaft, 189.3 litre capacity	25.2 lbs.	\$261.00
<b>Accessories for WSB Series Heavy Duty Immersion Blenders</b>			



Catalog Number	Description	Ship Wt.	List Price
WSBPPWE	Complete motor and whisk for commercial stick blenders, 230 volt, CE/RoHS approved.	27.35 lbs.	\$666.00
WSBPPCE	Motor head only for commercial stick blenders, 230 volt, CE/RoHS approved.	7 lbs.	\$350.00
WSB2W	10" / 25 cm Stainless Steel Whisk Attachment, for WSB Series Heavy Duty Immersion Blenders	21.2 lbs.	\$347.00
WSBBC	Bowl Clamp for hands-free operation of WSB Series Heavy Duty Immersion Blenders. (Shown in use.)	4 lbs.	\$173.00