

LIST PRICE COO

NIST TRACEABLE ISO 17025 CALIBRATION SERVICE FEES

CPT-100C	1 Temperature Calibration Point -78C to -1C (-140 to +31F)	\$107.05	SERVICE / US
CPT150C	1 Temperature Calibration Point 0 to +150C (+32 to +300F)	\$79.97	SERVICE / US
CPT300C	1 Temperature Calibration Point +151 to +300C (+301 to +572F)	\$107.05	SERVICE / US
CPT500C	1 Temperature Calibration Point +301 to +500C (+573 to +932F)	\$152.98	SERVICE / US
CAL SVC	FEE FOR TEMPS BETWEEN 300/500C & 572/932F <i>*one time fee per service</i>	\$50.00	SERVICE / US
CPT RH	1 Relative Humidity (Rh) Calibration Point	\$101.89	SERVICE / US
CPTHYD3_ISO	3 Calibration Points_ISO 17025	\$272.25	SERVICE / US
CPTHYD1	1 Calibration Point (MIDPOINT of scale)	\$103.84	SERVICE / US
CPTHYD2	2 Calibration Points (LOW & HIGH of scale)	\$155.75	SERVICE / US
CPTHYD3	3 Calibration Points (LOW/MID/HIGH of scale)	\$207.67	SERVICE / US

**thermo-hydrometers may require additional service codes for temperature calibration*

FLUORO-POLYMER COATING FOR THERMOMETERS

TFE-COAT 1	Up to 200C/400F LIG Thermometers	\$16.55	US
TFE-COAT 2	Up to 260C/500F LIG Thermometers	\$23.65	US

THERMOMETER CASES

87006	12" Leatherette, Felt-Lined Thermometer Case	\$41.22	CN
87008	18" Leatherette, Felt-Lined Thermometer Case	\$49.22	CN
87009	24" Leatherette, Felt-Lined Thermometer Case	\$58.53	CN
87011	27" Leatherette, Felt-Lined Thermometer Case	\$58.53	CN
87010	9 Place, Leatherette, Felt-Lined Thermometer Case	\$108.53	CN

THERMOMETER ARMORS

ACC6ARM	Thermometer Armor 6"	\$29.04	CN
ACC8ARM	Thermometer Armor 8"	\$29.04	CN
ACC12ARM	Thermometer Armor 12"	\$31.61	CN
ACC16ARM	Thermometer Armor 16"	\$34.29	CN
ACC18ARM	Thermometer Armor 18" DISCONTINUED		
ACC24ARM	Thermometer Armor 24" DISCONTINUED		

		<u>LIST PRICE</u>	<u>COO</u>
<u>THEMOMETER CUP CASES</u>			
ACC272CASE	Scoop Master, Plastic Case	\$80.17	DE
ACCWB	Wood Back, Brass Cup Case	\$59.47	DE
<u>THERMOMETER ACCESSORIES</u>			
ACCTM	Thermometer Magifier	\$81.36	US
ACCTH1000	Stainless Thermometer/pH Electrode/Test Tube Holder	\$105.37	DE
ACCTR100	Thermometer Rack, Holds 25 thermometers	\$99.21	CN
NR-100-25	Non-Roll Devise for Thermometers, 25/Pack	\$5.41	CN
<u>THERMOMETER BUFFERS</u>			
*bottles w/ buffer exclude grommet/cap for thermometer insertion (KITS ARE LISTED ON PRECISION THERM TAB)			
ACCBOT30BIO	30ml PETG plastic bottle w/30 of BIOSAFE Liquid	\$13.32	US
ACCBOT60BIO	60ml PETG plastic bottle w/60 of BIOSAFE Liquid	\$16.02	US
ACCBOT250BIO	250ml PETG plastic bottle w/60 of BIOSAFE Liquid	\$39.94	US
ACCBOT30SD	30ml PC plastic bottle w/30 of Sand	\$16.02	US
ACCBOT60SD	60ml PC plastic bottle w/60 of Sand	\$18.61	US
ACCBOT250SD	250ml PC plastic bottle w/60 of Sand	\$39.94	US
<u>HYDROMETER & GRADUATED CYLINDERS</u>			
UCY600H	600ml HYDROMETER CYLINDER, GLASS, 375 x 50mm	\$78.68	US
UCY1200H	1200ml HYDROMETER CYLINDER, GLASS, 450 x 65mm	\$99.66	US
UCY500BA	500ml GRADUATED CLASS "A", GLASS, 380 x 60mm	\$117.98	US
UCY500PP	500ml GRADUATED, POLYPROPYLENE, 360 x 60mm	\$44.25	US
DCY250PP	250ml GRADUATED, POLYPROPYLENE, 90 x 54mm	\$36.65	US
<u>HYDROMETER RACKS</u>			
ACCHR5	FOR SHORT HYDROMETERS	\$106.71	US
ACCHR7	FOR LONG HYDROMETERS	\$110.28	US
<u>HYDROMETER 8-PLACE CASE</u>			
GW3244CASE	CASE FOR 8 305mm HYDROMETERS	\$101.15	US

CPT Calibration service codes listed are for sourcing new products which require calibration.

<u>CAT. NO.</u>	<u>HTS Code 9025.11.4000</u>	<u>LIST PRICE</u>	<u>Haz(Hg)</u>	<u>COO</u>
<u>LIQUID-IN-GLASS THERMOMETERS</u>				
<u>N.I.S.T. & FACTORY CERTIFIED GLASS THERMOMETERS - MERCURY</u>				
ACC6113FC	-20/110C X 1.0C, 76mm Imm, Cert@ 0, 37, 56C	\$244.86	UN3506	DE
ACC611FC	-20/110C X 1.0C, Tol Imm, Cert@ 0, 37, 56C	\$244.86	UN3506	DE
ACC1003FC	-1/51C x 0.1C, Tol Imm, Cert@ 0, 10, 20, 30, 37, 40, 50C	\$538.92	UN3506	DE
ACC10033FC	-1/51C x 0.1C, 76mm Imm, Cert@ 0, 10, 20, 30, 37, 40, 50C	\$538.92	UN3506	DE
ACC1005FC	-1/101C x 0.1C, Tol Imm, Cert@ 0, 10, 20, 25, 37, 56, 70, 80, 90, 100C	\$818.55	UN3506	DE
ACC10053FC	-1/101C x 0.1C, 76mm Imm, Cert@ 0, 10, 20, 25, 37, 56, 70, 80, 90, 100C	\$818.55	UN3506	DE
ACC1007FC	-1/201C, x 0.2C, Tol Imm, Cert@ 0, 25, 37, 44.6, 56, 105, 121, 180, 200C	\$818.55	UN3506	DE
ACC10073FC	-1/201C, x 0.2C, 76mm Imm, Cert@ 0, 25, 37, 44.6, 56, 105, 121, 180, 200C	\$818.55	UN3506	DE
ACC10023FC	+30/124F x 0.2F, 76mm Imm, Cert@ 32, 40, 60, 100, 120F	\$421.87	UN3506	DE
ACC80135SCFC	+80/135C x 0.5C, Tol Imm, Cert@ 121C, Safety Coated	\$212.70	UN3506	DE
ACC80135SCFCA	+80/135C x 0.5C, Tol Imm, Cert@ 121C, Safety Coated, w/armor	\$228.98	UN3506	DE
ACC80135DSCFC	+80/135C x 0.5C, Tol Imm, Cert@ 120C, Safety Coated, For Dairies	\$240.44	UN3506	DE
ACC80135DSCFCA	+80/135C x 0.5C, Tol Imm, Cert@ 120C, Safety Coated, w/ armor, For Dairies	\$283.52	UN3506	DE
ACC10113FC	-1/110C x 0.2C, 76mm Imm, Cert@ 0, 4, 25, 37, 44.6, 56, 105C	\$539.07	UN3506	DE
ACC713FC	+34.5/46C, x 0.1C, 76mm Imm, Cert@ 37, 44.5C	\$203.65	UN3506	DE
ACC934FC	+24/38C x 0.05C, 95mm Imm, Cert@ 0, 25, 30, 37C	\$303.39	UN3506	DE
<u>N.I.S.T. & FACTORY CERTIFIED GLASS THERMOMETERS - SPIRIT</u>				
ACC1103BLSFC	-20/110C X 1.0C, 76mm Imm, Cert@ 0, 50, 100C	\$243.95		DE
ACC6113SFC	-20/110C x 0.5C, 76mm Imm, Cert@ -10, 0, 4, 25, 37, 56, 105C	\$441.91		DE
ACC6113SXFC	-20/110C x 1.0C, 76mm Imm, Cert@ 0, 50, 100C	\$244.86		DE
ACC6113SXYFC	-20/110C x 1.0C, 76mm Imm, Cert@ 0, 37, 56C	\$246.15		DE
ACC6113XXSFC	-20/110C x 0.5C, 76mm Imm, Spirit, CERT @ 0, 37, 56C	\$264.89		DE
ACC611SSCFC	-20/110C x 1.0C, Tol Imm, Cert@ 0, 37, 56C	\$244.86		DE
ACC6123SFC	0/300F x 2.0F, 76mm Imm, Cert@ 32, 100, 200, 300F	\$331.44		DE
ACC6133SFC	-20/150C x 1.0C, 76mm Imm, Cert @ 0, 60, 120	\$244.86		DE
ACC1003SFC	-1/51C x 0.1C, Tol Imm, Cert@ 0, 10, 20, 30, 37, 40, 50C	\$538.92		DE
ACC10033SFC	-1/51C x 0.1C, 76mm Imm, Cert@ 0, 10, 20, 30, 37, 40, 50C	\$538.92		DE
ACC1005SFC	-1/101C x 0.1C, Tol Imm, Cert@ 0, 10, 20, 25, 37, 44.5, 56, 70, 80, 90, 100C	\$818.55		DE
ACC10053SFC	-1/101C x 0.1C, 76mm Imm, Cert@ 0, 10, 20, 25, 37, 44.5, 56, 70, 80, 90, 100C	\$818.55		DE
ACC1007SFC	-1/201C, x 0.2C, Tol Imm, Cert@ 0, 25, 37, 44.6, 56, 105, 121, 180, 200C	\$818.55		DE
ACC10033BLSFC	-1/51C x 0.2C, 76mm Imm, Cert@ 0, 10, 20, 30, 37, 40, 50C	\$539.07		DE
ACC10613SFC	-1/61C x 0.1C, 76mm Imm, Cert@ 0, 37, 56C	\$264.89		DE
ACC1005BLSFC	-1/101C x 0.2C, Tol Imm, Cert@ 0, 10, 20, 25, 37, 56, 70, 80, 90, 100C	\$818.55		DE
ACC10053BLSFC	-1/101C x 0.2C, 76mm Imm, Cert@ 0, 10, 20, 25, 37, 56, 70, 80, 90, 100C	\$818.55		DE
ACC1007BLSFC	-1/201C, x 0.5C, Tol Imm, Cert@ 0, 25, 37, 44.5, 56, 105, 121, 180, 200C	\$818.55		DE
ACC10073BLSFC	-1/201C, x 0.5C, 76mm Imm, Cert@ 0, 25, 37, 44.5, 56, 105, 121, 180, 200C	\$818.55		DE
ACC603FC	-100/50C x 1.0C, Tol Imm, Cert@ -78, -20, 0, 4, 25, 37C	\$469.81		DE
ACC603XSFC	*-100/50C x 1.0C, Tol Imm, Cert@ -80, -70,-20, 0, 37C	\$469.81		DE
ACC6033FC	-100/50C x 1.0C, 76mm Imm, Cert@ -78, -20, 0, 4, 25, 37C	\$469.81		DE
ACC713BLSFC	+34.5/46C, x 0.1C, 76mm Imm, Cert@ 37, 44.5C	\$203.65		DE
ACC713TBLSSFC	+34.5/46C, x 0.1C, Toal Imm, w/ Safety Coat, Cert@ 37, 44.5C	\$221.21		DE
ACC10113BLSFC	-1/110C x 0.2C, 76mm Imm, Cert @ 0, 4, 25, 37, 44.6, 56, 105C	\$539.07		DE

<u>CAT. NO.</u>	<u>HTS Code 9025.11.4000</u>	<u>LIST PRICE</u>	<u>Haz(Hg)</u>	<u>COO</u>
ACC6103XXSFC	0/230F x 0.5F, 76mm, cert @ 32, 45, 161F	\$201.69		DE
<u>ASTM CERTIFIED TRACEABLE TO N.I.S.T. THERMOMETERS - MERCURY - Tol IMMERSION</u>				
ACC062CFC	-38/+2C x 0.1C, Cert@ -37, -30, -20, -10, 0C	\$422.63	UN3506	DE
ACC063CFC	-8/32C x 0.1C, Cert@ -7, 0, 10, 20, 30C	\$361.34	UN3506	DE
ACC064CFC	+25/55C x 0.1C, Cert@ 0, 25, 35, 45, 55C	\$406.22	UN3506	DE
ACC065CFC	+50/80C x 0.1C, Cert@ 0, 50, 60, 70, 80C	\$410.43	UN3506	DE
ACC066CFC	+75/105C x 0.1C, Cert@ 0, 75, 85, 95, 105C	\$441.69	UN3506	DE
ACC067CFC	+95/155C x 0.2C, Cert@ 0, 100, 110, 130, 150C	\$474.39	UN3506	DE
ACC068CFC	+145/205C x 0.2C, Cert@ 0, 150, 170, 190, 205C	\$668.03	UN3506	DE
ACC069CFC	+195/305C x 0.5C, Cert@ 0, 200, 235, 270, 305C	\$668.03	UN3506	DE
ACC070CFC	+295/405C x 0.5C, Cert@ 0, 300, 335, 370, 400C	\$781.20	UN3506	DE
<u>ASTM CERTIFIED TRACEABLE TO N.I.S.T. THERMOMETERS - MERCURY - Tol IMMERSION</u>				
ACC062FFC	-36/35F x 0.2F, Cert@ -35, -15, 0, 15, 32F	\$422.63	UN3506	DE
ACC063FFC	+18/89F x 0.2F, Cert@ 20, 32, 50, 70, 88F	\$361.34	UN3506	DE
ACC064FFC	+77/131F x 0.2F, Cert@ 32, 80, 95, 115, 130F	\$406.22	UN3506	DE
ACC065FFC	+122/176F x 0.2F, Cert@ 32, 125, 145, 160, 175F	\$410.43	UN3506	DE
ACC066FFC	+167/221F x 0.2F, Cert@ 32, 168, 185, 200, 220F	\$441.69	UN3506	DE
ACC067FFC	+203/311F x 0.5F, Cert@ 32, 205, 240, 275, 310F	\$474.39	UN3506	DE
ACC068FFC	+293/401F x 0.5F, Cert@ 32, 300, 340, 370, 400F	\$668.03	UN3506	DE
ACC069FFC	+383/581F x 1.0F, Cert @32, 400, 460, 520, 580F	\$668.03	UN3506	DE
ACC070FFC	+563/761F x 1.0F, Cert@ 32, 570, 640, 700, 760F	\$781.20	UN3506	DE
<u>ASTM CERTIFIED TRACEABLE TO N.I.S.T. THERMOMETERS-SET OF 9 THERMOMETERS</u>				
ACC6270CFC	ASTM 62C to 70C Tol IMMERSION, -38 to 405C w/ case	\$3,233.95	UN3506	DE
ACC6270FFC	ASTM 62F to 70F Tol IMMERSION, -36 to 761F w/ case	\$3,233.95	UN3506	DE
<u>ASTM CERTIFIED TRACEABLE TO N.I.S.T. THERMOMETERS - Tol IMMERSION BLUE SPIRIT</u>				
ACC062CBLSFC	0S62C, -38/+2C, Tol IMM, 400mm, BLUE SPIRIT, Cert@ -37,-15,-20,-10, 0C	\$573.40		DE
ACC062FBLSFC	0S62F, -36/+35F, Tol IMM, 400mm, BLUE SPIRIT, Cert @-35,-15, 0, 15, 32F	\$573.40		DE
ACC063CBLSFC	0S63C, -8/+32C, Tol IMM, 400mm, BLUE SPIRIT, Cert@ -7,0,10,20,30C	\$538.75		DE
ACC063FBLSFC	0S63F, +18/89F, Tol IMM, 400mm, BLUE SPIRIT, Cert@ 20,32,50,70 88F	\$538.75		DE
ACC064CBLSFC	0S64C, +25/55C, Tol IMM, 400mm, BLUE SPIRIT, Cert@ 0,25,35,45,55C	\$562.91		DE
ACC064FBLSFC	0S64F, +77/131F, Tol IMM, 400mm, BLUE SPIRIT, Cert@ 32,80,95,115,130F	\$562.91		DE
ACC065CBLSFC	0S65C, -50/80C, Tol IMM, 400mm, BLUE SPIRIT, Cert@ 0,50,60,70,80C	\$574.52		DE
ACC065FBLSFC	0S65F, +122/176F, Tol IMM, 400mm, BLUE SPIRIT Cert@ 32,125,145,160,175F	\$574.52		DE
ACC066CBLSFC	0S66C, +75/105C, Tol IMM, 400mm, BLUE SPIRIT, Cert@ 0,75,85,95,105C	\$592.23		DE
ACC066FBLSFC	0S66F, +167/221F, Tol IMM, 400mm, BLUE SPIRIT, Cert@32,168,185,200,220F	\$592.23		DE

<u>CAT. NO.</u>	<u>HTS Code 9025.11.4000</u>	<u>LIST PRICE</u>	<u>Haz(Hg)</u>	<u>COO</u>
ACC067CBLSFC	0S67C, +95/155C, Tol IMM, 400mm, BLUE SPIRIT, Cert@ 0,100,110,120,130C	\$592.23		DE
ACC067FBLSFC	0S67F, +203/311F, Tol IMM, 400mm, BLUE SPIRIT, Cert@ 32,205,240,275,310F	\$592.23		DE

ASTM CERTIFIED TRACEABLE TO N.I.S.T. THERMOMETERS-SET OF 6 THERMOMETERS - BLUE SPIRIT

ACC6267CBLSFC	ASTM 0S62C to 0S67C Tol IMMERSION, -38 to 155C w/ case	\$4,155.62		DE
ACC6267FBLSFC	ASTM 0S62F to 0S67F Tol IMMERSION, -36 to 311F w/ case	\$4,155.62		DE

"ACCUTHERM" BOTTLE-Enc CERTIFIED INSTRUMENT THERMOMETERS - RED SPIRIT

ACCULF010S	Ultra-Low Freezer, -90/+20C x 1.0C	\$86.57		DE
ACCF0101S	Freezer, -30/0C x 0.5C	\$71.89		CN, DE
ACCF100S	Freezer, -25/-5C x 0.1C	\$101.18		DE
ACCR0201S	Refrigerator, -5/+15C x 0.5C	\$71.89		CN, DE
ACCR20001S	Refrigerator, -2/+10C x 0.1C	\$101.12		DE
ACCBBS	Blood Bank Refrigerator, -5/+20C x 0.5C	\$101.12		DE
ACCB210S	Blood Bank Refrigerator, -2/+10C x 0.1C	\$130.46		DE
ACCI0301S	Incubator, +18/50C x 0.5C	\$71.89		CN, DE
ACCI300S	Incubator, +15/30C x 0.1C	\$101.12		DE
ACCI310S	Incubator, +25/45C x 0.1C	\$101.12		DE
ACCI370S	Incubator High Range, +30/110C x 0.5C	\$69.26		DE
ACCAM0401S	Ambient/Room, +10/30C x 0.5C	\$71.89		DE
ACCOV070S	Oven, +20/130C x 1.0C	\$71.89		DE
ACCOV700S	Oven, +50/115C x 0.5C	\$79.79		DE
ACCOV080S	Oven, +35/180C x 1.0C **new range	\$79.79		DE
ACCD301	Dry Bath, +18/60C x 0.5C	\$79.79		DE
ACCD40	Dry Bath, +50/110C x 0.5C	\$79.79		DE
ACC0040SIC	Cooler Thermometer -10/20C x 0.5C	\$63.86		DE

VACCINE THERMOMETER

ACC1482163	-30/25C x 0.5C Spirit Ref/Freezer Encased Thermometer w/CC	\$87.80		DE
------------	--	---------	--	----

"ACCUTHERM" BOTTLE-Enc CERTIFIED INSTRUMENT THERMOMETERS -BLUE SPIRIT

ACCF0101BLS	Freezer, -30/0C x 0.5C	\$63.86		CN, DE
ACCR0201BLS	Refrigerator, -5/+15C x 0.5C	\$63.86		CN, DE
ACCI0301BLS	Incubator, +18/50C x 0.5C	\$63.86		CN, DE

"ACCUTHERM" BOTTLE-Enc CERTIFIED INSTRUMENT THERMOMETERS -BLUE SPIRIT, SAFETY-COATED

ACCF0101BLST	Freezer, -30/0C x 0.5C	\$66.50		CN, DE
ACCR0201BLST	Refrigerator, -5/+15C x 0.5C	\$66.50		CN, DE
ACCI0301BLST	Incubator, +18/50C x 0.5C	\$66.50		CN, DE

<u>CAT. NO.</u>	<u>HTS Code 9025.11.4000</u>	<u>LIST PRICE</u>	<u>Haz(Hg)</u>	<u>COO</u>
<u>"ACCUTHERM" CERTIFIED INSTRUMENT THERMOMETERS, BLUE SPIRIT, REFILLS</u>				
ACCF0101BLS U	Freezer, -30/0C x 0.5C	\$39.92		CN, DE
ACCR0201BLS U	Refrigerator, -5/+15C x 0.5C	\$39.92		CN, DE
ACCI0301BLS U	Incubator, +18/50C x 0.5C	\$39.92		CN, DE
<u>"ACCUTHERM" 6ml VIAL-Enc CERTIFIED INSTRUMENT THERMOMETERS - RED SPIRIT IN PROTECTIVE PLASTIC TUBE W/MAGNETS</u>				
ACCF2401S	Freezer, -30/0C x 0.5C	\$66.50		CN, DE
ACCR2402S	Refrigerator, -5/+15C x 0.5C	\$66.50		CN, DE
ACCI2403S	Incubator, +18/50C x 0.5C	\$66.50		CN, DE
ACCULF2404S	Ultra-Low Freezer, -90/+20C x 1.0C	\$86.57		DE
<u>"ACCUTHERM" 6ml VIAL-Enc CERTIFIED INSTRUMENT THERMOMETERS - BLUE SPIRIT IN PROTECTIVE PLASTIC TUBE W/MAGNETS</u>				
ACCF2401BLS	Freezer, -30/0C x 0.5C	\$66.50		DE
ACCR2402BLS	Refrigerator, -5/+15C x 0.5C	\$66.50		DE
ACCI2403BLS	Incubator, +18/50C x 0.5C	\$64.57		DE
<u>"ACCUTHERM" BOTTLE-Enc CERTIFIED INSTRUMENT THERMOMETERS - MERCURY</u>				
ACCF0101	Freezer, -30/0C x 0.5C	\$74.57	UN3506	DE
ACCF100	Freezer, -25/+5C x 0.1C	\$103.77	UN3506	DE
ACCR0201	Refrigerator, -5/+15C x 0.5C	\$74.57	UN3506	DE
ACCR200	Refrigerator, -2/+10C x 0.1C	\$103.77	UN3506	DE
ACCI0301	Incubator, +18/50C x 0.5C	\$74.57	UN3506	DE
ACCI300	Incubator, +15/30C x 0.1C	\$103.77	UN3506	DE
ACCI310	Incubator, +25/45C x 0.1C	\$103.77	UN3506	DE
ACCOV070	Oven, +20/130C x 1.0C	\$74.57	UN3506	DE
ACCOV700	Oven, +95/115C x 0.1C	\$93.04	UN3506	DE
ACCOV080	Oven, +35/200C x 1.0C	\$79.86	UN3506	DE
<u>"ACCUTHERM" LABORATORY INSTRUMENT THERMOMETERS - SPIRIT</u>				
ACC6471S	Dry Bath/Block Heater 0/110C x 0.5C, 35mm Imm	\$69.26		DE
ACC647SSC	Dry Bath/Block Heater 0/110C x 0.5C, 35mm Imm, Safety Coated	\$82.58		DE
ACC715S	Dry Bath/Block Heater 25/57C x 0.5C, 35mm Imm	\$69.26		DE
ACC715SSC	Dry Bath/Block Heater 25/57C x 0.5C, 35mm Imm, Safety coated	\$79.86		DE
ACC641S	Incubator, 0/80C x 0.5C, 100mm Imm	\$50.54		DE
ACC641SSC	Incubator, 0/80C x 0.5C, 100mm Imm, Safety Coated	\$63.82		DE
ACC6371S	Incubator, 0/50C x 0.5C, 150mm Imm	\$57.25		DE
ACC6371SSC	Incubator, 0/50C x 0.5C, 150mm Imm, Safety Coated	\$69.26		DE
ACC6391S	Incubator, 0/70C x 1.0C, 150mm Imm	\$62.55		DE
ACC6391SSC	Incubator, 0/70C x 1.0C, 150mm Imm, Safety Coated	\$73.28		DE
ACC711S	Incubator, 20/100C x 0.5C, 76mm Imm	\$71.85		DE

<u>CAT. NO.</u>	<u>HTS Code 9025.11.4000</u>	<u>LIST PRICE</u>	<u>Haz(Hg)</u>	<u>COO</u>
ACC711SSC	Incubator, 20/100C x 0.5C, 76mm Imm, Safety Coated	\$82.58		DE
ACC2457S	Incubator, 25,37,& 56C x 0.5C, 76mm Imm	\$54.55		DE
ACC2457SSC	Incubator, 25,37 & 56C x 0.5C, 76mm Imm, Safety Coated	\$66.55		DE
ACC713BLS	Coliform Incubator, 34.5/46C x 0.1C, 76mm Imm	\$107.77		DE
ACC713BLSSC	Coliform Incubator, 34.5/46C x 0.1C, 76mm Imm, Safety Coated	\$122.50		DE
ACC713TBLS	Coliform Incubator, 34.5/46C x 0.1C, Tol Imm	\$107.76		DE
ACC713TBLSSC	Coliform Incubator, 34.5/46C x 0.1C, Tol Imm, Safety Coated	\$122.49		DE
ACC717S	Blood Bank Refrigerator -5/20C x 1.0C, Tol Imm	\$41.23		DE
ACC717SSC	Blood Bank Refrigerator -5/20C x 1.0C, Tol Imm, Safety Coated	\$54.55		DE
ACC535S1	Tissue Flootation Bath, -15/105C x 1.0C Tol Imm	\$31.93		DE
ACC535S1SC	Tissue Flootation Bath, -15/105C x 1.0C Tol Imm, Safety Coated	\$42.66		DE
ACC637S	Oven, 0/200C x 1.0C, 76mm Imm	\$67.84		DE
ACC637SSC	Oven, 0/200C x 1.0C, 76mm Imm Safety Coated	\$79.86		DE
ACC639S	Oven, 0/200C x 1.0C, 125mm Imm	\$67.84		DE
ACC639SSC	Oven, 0/200C x 1.0C, 125 Imm, Safety Coated	\$79.86		DE
ACC24130S	Oven, 20/130C x 1.0C, Tol Imm	\$59.95		DE
ACC24130SSC	Oven, 20/130C x 1.0C, Tol Imm, Safety Coated	\$71.85		DE
ACC7200S	Oven, 35/200C x 1.0C, Tol Imm	\$59.95		DE
ACC7200SSC	Oven, 35/200C x 1.0C, Tol Imm, Safety Coated	\$71.85		DE

"ACCUTHERM" LABORATORY INSTRUMENT THERMOMETERS - MERCURY

ACC645	Short Immersion, 12mm Imm, 75mm long	\$55.96	UN3506	DE
ACC645SC	Short Immersion, 12mm Imm, 75mm long, Safety Coated	\$70.25	UN3506	DE
ACC647	Dry Bath/Block Heater, 0/110C x 0.5C, 35mm Imm	\$71.85	UN3506	DE
ACC6471SC	Dry Bath/Block Heater, 0/110C x 0.5C, 35mm Imm, Safety Coated	\$87.07	UN3506	DE
ACC715	Dry Bath/Block Heater, 25/57C x 0.5C, 35mm Imm	\$71.85	UN3506	DE
ACC715SC	Dry Bath/Block Heater, 25/57C x 0.5C, 35mm Imm, Safety coated	\$84.31	UN3506	DE
ACC641	Incubator, 0/80C x 0.5C, 100mm Imm	\$58.02	UN3506	DE
ACC641SC	Incubator, 0/80C x 0.5C, 100mm Imm, Safety Coated	\$70.52	UN3506	DE
ACC6371	Incubator, 0/50C x 0.5C, 150mm Imm	\$60.86	UN3506	DE
ACC6371SC	Incubator, 0/50C x 0.5C, 150mm Imm, Safety Coated	\$73.15	UN3506	DE
ACC6391	Incubator, 0/70C x 0.5C, 150mm Imm	\$66.42	UN3506	DE
ACC6391SC	Incubator, 0/70C x 0.5C, 150mm Imm, Safety Coated	\$78.68	UN3506	DE
ACC711	Incubator, 20/100C x 0.5C, 76mm Imm	\$74.70	UN3506	DE
ACC711SC	Incubator, 20/100C x 0.5C, 76mm Imm, Safety Coated	\$87.07	UN3506	DE

"ACCUTHERM" LABORATORY INSTRUMENT THERMOMETERS - MERCURY

ACC934	Enzymology thermometer, 24/38C x 0.05C, 95mm Immersion	\$130.23	UN3506	DE
ACC713	Coliform Incubator, 34.5/46C x 0.1C, 76mm Imm	\$112.04	UN3506	DE
ACC713SC	Coliform Incubator, 34.5/56C x 0.1C, 76mm Imm, Safety Coated	\$125.08	UN3506	DE
ACC80135SC	Autoclave Maximum Registering, +80/135C x 0.5C, Safety Coated	\$105.37	UN3506	DE
ACC80135SCA	Autoclave Maximum Registering, +80/135C x 0.5C, Safety Coated, w/arm	\$117.27	UN3506	DE
ACC717	Blood Bank Refrigerator -5/20C x 1.0C, Tol Imm	\$44.59	UN3506	DE
ACC717SC	Blood Bank Refrigerator -5/20C x 1.0C, Tol Imm, Safety Coated	\$58.73	UN3506	DE
ACC629	Oven, 0/300C x 1.0C,125mm Imm	\$62.04	UN3506	DE
ACC631	Oven, 0/300C x 1.0C,150mm Imm	\$77.40	UN3506	DE

<u>CAT. NO.</u>	<u>HTS Code 9025.11.4000</u>	<u>LIST PRICE</u>	<u>Haz(Hg)</u>	<u>COO</u>
ACC633	Oven, 0/300C x 1.0C, 305mm Imm	\$103.77	UN3506	DE
ACC635	Oven, 0/200C x 1.0C, 305mm Imm	\$103.77	UN3506	DE
ACC635SC	Oven, 0/200C x 1.0C, 305mm Imm, Safety Coated	\$118.94	UN3506	DE
ACC637	Oven, 0/200C x 1.0C, 76mm Imm	\$66.43	UN3506	DE
ACC637SC	Oven, 0/200C x 1.0C, 76mm Imm Safety Coated	\$74.65	UN3506	DE
ACC639	Oven, 0/200C x 1.0C, 125mm Imm	\$71.85	UN3506	DE
ACC639SC	Oven, 0/200C x 1.0C, Tol Imm, Safety Coated	\$87.44	UN3506	DE
ACC7200	Oven, 35/200C x 1.0C, Tol Imm	\$63.32	UN3506	DE
ACC7200SC	Oven, 35/200C x 1.0C, Tol Imm, Safety Coated	\$75.90	UN3506	DE
ACC10113	-1 / 110C x 0.2C, 76mm Imm	\$111.78	UN3506	DE
ACC10113SC	-1 / 110C x 0.2C, 76mm Imm, Safety-Coated	\$127.96	UN3506	DE
ACCFISJOH1	Melting Point Thermometers, Fisher-Johns, 0/160C x 1.0C	\$63.96	UN3506	DE
ACCFISJOH2	Melting Point Thermometers, Fisher-Johns, 150/300C x 1.0C	\$77.55	UN3506	DE

SAMA PRECISION LABORATORY THERMOMETERS - MERCURY - 76mm/PARTIAL IMMERSION

ACC6093	-35/50C x 1.0C	\$43.17	UN3506	DE
ACC6113	-20/110C x 1.0C	\$43.17	UN3506	DE
ACC6113XX	-20/110C x 0.5C	\$69.26	UN3506	DE
ACC6133	-20/150C x 1.0C	\$43.17	UN3506	DE
ACC6153	-10/200C x 1.0C	\$43.17	UN3506	DE
ACC6173	-10/260C x 1.0C	\$57.51	UN3506	DE
ACC6233	-10/400C x 2.0C	\$66.43	UN3506	DE
ACC6083	-30/120F x 1.0F	\$43.17	UN3506	DE
ACC6103	0 /230F x 2.0F	\$43.17	UN3506	DE
ACC6103XX	0/230F x 0.5F	\$114.31	UN3506	DE
ACC61031X	0/230F x 1.0F	\$61.26	UN3506	DE
ACC6123	0/300F x 2.0F	\$43.17	UN3506	DE
ACC6143	+20/400F x 2.0F	\$43.17	UN3506	DE
ACC6163	+20/500F x 2.0F	\$60.34	UN3506	DE
ACC6203	+20/750F x 5.0F	\$69.26	UN3506	DE
ACC7043	+20/930F x 5.0F	\$149.10	UN3506	DE
ACC610113	-20/110C x 1.0C & 0/230F x 2.0F	\$46.00	UN3506	DE
ACC612133	-20/150C x 1.0C & 0/300F x 2.0F	\$48.98	UN3506	DE
ACC616173	-10/260C x 1.0C & +20/500F x 2.0F	\$60.34	UN3506	DE
ACC10033	-1/51C x 0.1C	\$117.72	UN3506	DE
ACC10053	-1/101C x 0.1C	\$153.07	UN3506	DE
ACC10073	-1/201C x 0.2C	\$169.95	UN3506	DE
ACC10023	+30/124F x 0.2F	\$117.72	UN3506	DE
ACC10023XX	+24/124F x 0.1F	\$177.34	UN3506	DE
ACC10043	+30/214F x 0.2F	\$137.10	UN3506	DE
ACC10063	+30/394F x 0.5F	\$152.22	UN3506	DE

SAMA PRECISION LABORATORY THERMOMETERS - MERCURY - Tol IMMERSION

ACC609	-35/50C x 1.0C	\$43.29	UN3506	DE
ACC611	-20/110C x 1.0C	\$43.29	UN3506	DE
ACC613	-20/150C x 1.0C	\$43.29	UN3506	DE

<u>CAT. NO.</u>	<u>HTS Code 9025.11.4000</u>	<u>LIST PRICE</u>	<u>Haz(Hg)</u>	<u>COO</u>
ACC615	-10/200C x 1.0C	\$43.29	UN3506	DE
ACC6151X	-10/220C x 0.1C	\$108.54	UN3506	DE
ACC617	-10/260C x 1.0C	\$60.34	UN3506	DE
ACC623	-10/400C x 2.0C	\$69.26	UN3506	DE

SAMA PRECISION LABORATORY THERMOMETERS - MERCURY - ToI IMMERSION

ACC608	-30/120F x 1.0F	\$43.29	UN3506	DE
ACC610	0/220F x 2.0F	\$43.29	UN3506	DE
ACC612	0/300F x 2.0F	\$43.29	UN3506	DE
ACC614	+20/400F x 2.0F	\$60.34	UN3506	DE
ACC616	+20/500F x 2.0F	\$60.34	UN3506	DE
ACC620	+20/750F x 5.0F	\$69.26	UN3506	DE
ACC1001	-35/50C x 0.1C	\$136.71	UN3506	DE
ACC1003	-1/51C x 0.1C	\$117.72	UN3506	DE
ACC1005	-1/101C x 0.1C	\$152.22	UN3506	DE
ACC1007	-1/201C x 0.2C	\$159.72	UN3506	DE
ACC1000	-30/120F x 0.2F	\$137.10	UN3506	DE
ACC1002	+30/124F x 0.2F	\$117.72	UN3506	DE
ACC1004	+30/214F x 0.2F	\$137.10	UN3506	DE
ACC1006	+30/394F x 0.5F	\$152.22	UN3506	DE

SAMA PRECISION FLUOROPOLYMER-COATED THERMOMETERS - MERCURY - PARTIAL IMMERSION

ACC6093SC	-35/50C x 1.0C, Safety Coated	\$58.73	UN3506	DE
ACC6113SC	-20/110C x 1.0C, Safety Coated	\$58.73	UN3506	DE
ACC6113XXSC	-20/110C x 0.5C, Safety Coated	\$87.59	UN3506	DE
ACC6133SC	-20/150C x 1.0C, Safety Coated	\$58.73	UN3506	DE
ACC6153SC	-10/200C x 1.0C, Safety Coated	\$58.73	UN3506	DE
ACC6173SC	-10/260C x 1.0C, Safety Coated	\$74.40	UN3506	DE
ACC6083SC	-30/120F x 1.0F, Safety Coated	\$58.73	UN3506	DE
ACC6103SC	0 /230F x 2.0F, Safety Coated	\$58.73	UN3506	DE
ACC6103XXSC	0 /230F x 0.5F, Safety Coated	\$87.59	UN3506	DE
ACC6123SC	0/300F x 2.0F, Safety Coated	\$58.73	UN3506	DE
ACC6143SC	+20/400F x 2.0F, Safety Coated	\$58.73	UN3506	DE
ACC6163SC	+20/500F x 2.0F, Safety Coated	\$74.40	UN3506	DE
ACC610113SC	-20/110C x 1.0C & 0/230F x 2.0F, Safety Coated	\$64.61	UN3506	DE
ACC612133SC	-20/150C x 1.0C & 0/300F x 2.0F, Safety Coated	\$64.61	UN3506	DE
ACC616173SC	-10/260C x 1.0C & +20/500F x 2.0F, Safety Coated	\$75.90	UN3506	DE
ACC10033SC	-1/51C x 0.1C, Safety Coated	\$132.85	UN3506	DE
ACC10053SC	-1/101C x 0.1C, Safety Coated	\$150.81	UN3506	DE
ACC10073SC	-1/201C x 0.2C, Safety Coated	\$168.10	UN3506	DE
ACC10023SC	+30/124F x 0.2F, Safety Coated	\$132.85	UN3506	DE
ACC10043SC	+30/214F x 0.2F, Safety Coated	\$150.81	UN3506	DE
ACC10063SC	+30/394F x 0.5F, Safety Coated	\$168.10	UN3506	DE

SAMA PRECISION FLUOROPOLYMER-COATED THERMOMETERS - MERCURY - ToI IMMERSION

ACC609SC	-35/50C x 1.0C, Safety Coated	\$58.73	UN3506	DE
----------	-------------------------------	---------	--------	----

<u>CAT. NO.</u>	<u>HTS Code 9025.11.4000</u>	<u>LIST PRICE</u>	<u>Haz(Hg)</u>	<u>COO</u>
ACC611SC	-20/110C x 1.0C, Safety Coated	\$58.73	UN3506	DE
ACC613SC	-20/150C x 1.0C, Safety Coated	\$58.73	UN3506	DE
ACC615SC	-10/200C x 1.0C, Safety Coated	\$58.73	UN3506	DE
ACC6151XSC	-10/220C x 0.5C, Safety Coated	\$125.33	UN3506	DE
ACC617SC	-10/260C x 1.0C, Safety Coated	\$76.03	UN3506	DE
ACC608SC	-30/120F x 1.0F, Safety Coated	\$58.73	UN3506	DE
ACC610SC	0 /230F x 2.0F, Safety Coated	\$58.73	UN3506	DE
ACC612SC	0/300F x 2.0F, Safety Coated	\$58.73	UN3506	DE
ACC614SC	+20/400F x 2.0F, Safety Coated	\$58.73	UN3506	DE
ACC616SC	+20/500F x 2.0F, Safety Coated	\$76.03	UN3506	DE
ACC1001SC	-35/50C x 0.1C, Safety Coated	\$151.18	UN3506	DE
ACC1003SC	-1/51C x 0.1C, Safety Coated	\$132.85	UN3506	DE
ACC1005SC	-1/101C x 0.1C, Safety Coated	\$150.81	UN3506	DE
ACC1007SC	-1/201C x 0.2C, Safety Coated	\$168.10	UN3506	DE
ACC1000SC	-30/120F x 0.2F, Safety Coated	\$147.91	UN3506	DE
ACC1002SC	+30/124F x 0.2F, Safety Coated	\$132.85	UN3506	DE
ACC1004SC	+30/214F x 0.2F, Safety Coated	\$150.81	UN3506	DE
ACC1006SC	+30/394F x 0.5F, Safety Coated	\$168.10	UN3506	DE

PRECISION RED-SPIRIT LABORATORY THERMOMETERS - PARTIAL IMMERSION

ACC6093S	-35/50C x 1.0C	\$41.49		DE
ACC6113S	-20/110C x 1.0C	\$41.49		DE
ACC6113XXS	-20/110C x 0.5C	\$54.01		DE
ACC6133S	-20/150C x 1.0C	\$41.49		DE
ACC6133XXS	-20/150C x 0.5C	\$63.82		DE
ACC6153S	-10/200C x 1.0C	\$41.49		DE
ACC6173S	-10/260C x 1.0C	\$41.49		DE
ACC6083S	-30/120F x 1.0F	\$41.49		DE
ACC6103S	0/230F x 2.0F	\$41.49		DE
ACC6103XXS	0/230F x 0.5F	\$100.83		DE
ACC6123S	0/300F x 2.0F	\$41.49		DE
ACC6123XS	+30/300F x 1.0F	\$93.05		DE
ACC6143S	+20/400F x 2.0F	\$41.49		DE
ACC6013	-50/+50C x 1.0C	\$44.99		DE
ACC6033	-100/+50C X 1.0C	\$79.87		DE
ACC10023S	+30/124F x 0.2F	\$123.79		DE
ACC10033S	-1/51C x 0.1C	\$123.79		DE
ACC10613S	-1/61C x 0.1C	\$127.80		DE
ACC10043S	+30/214F x 0.2F	\$137.10		DE
ACC10053S	-1/101C x 0.1C	\$179.69		DE
ACC10053BSCERT	-1/101C x 0.1C, *Blue Spirit, 630mm	\$169.44		DE
ACC10063S	+30/394F x 0.5F	\$146.41		DE
ACC10073BSCERT	-1/201C x 0.2C, *Blue Spirit, 630mm	\$222.91		DE
ACC10113BLS	-1/110C x 0.2C, Blue Spirit	\$124.81		DE

PRECISION RED-SPIRIT LABORATORY THERMOMETERS - ToI IMMERSION

<u>CAT. NO.</u>	<u>HTS Code 9025.11.4000</u>	<u>LIST PRICE</u>	<u>Haz(Hg)</u>	<u>COO</u>
ACC609S	-35/50C x 1.0C	\$41.49		DE
ACC611S	-20/110C x 1.0C	\$41.49		DE
ACC613S	-20/150C x 1.0C	\$41.49		DE
ACC615S	-10/200C x 1.0C	\$41.49		DE
ACC6151XS	-10/220C x 0.5C	\$108.54		DE
ACC608S	-30/120F x 1.0F	\$41.49		DE
ACC610S	0/230F x 2.0F	\$41.49		DE
facc6033	0/300F x 2.0F	\$41.49		DE
ACC601	-50/+50C x 1.0C	\$44.99		DE
ACC603	-100/+50C X 1.0C	\$79.87		DE
ACC605	-200/+30C x 1.0C *Blue Spirit, 350mm	\$211.06		DE
ACC1002S	+30/124F x 0.2F	\$123.79		DE
ACC1003S	-1/51C x 0.1C	\$123.79		DE
ACC1004S	+30/214F x 0.2F	\$137.10		DE
ACC1005S	-1/101C x 0.1C	\$186.34		DE
ACC1005BSCERT	-1/101C x 0.1C, *Blue Spirit, 630mm	\$169.44		DE
ACC1006S	+30/394F x 0.5F	\$154.57		DE
ACC1007BSCERT	-1/201C x 0.2C, *Blue Spirit, 630mm	\$191.06		DE

PRECISION RED-SPIRIT FLUOROPOLYMER-COATED THERMOMETERS - PARTIAL IMMERSION

ACC6093SSC	-35/50C x 1.0C, Safety Coated	\$56.60		DE
ACC6113SSC	-20/110C x 1.0C, Safety Coated	\$56.60		DE
ACC6113XXSSC	-20/110C x 0.5C, Safety Coated	\$67.75		DE
ACC6133SSC	-20/150C x 1.0C , Safety Coated	\$56.60		DE
ACC6153SSC	-10/200C x 1.0C, Safety Coated	\$56.60		DE
ACC6151XSSC	-10/220C x 0.5C, Safety Coated	\$125.22		DE
ACC6083SSC	-30/120F x 1.0F, Safety Coated	\$56.60		DE
ACC6103SSC	0/230F x 2.0F, Safety Coated	\$56.60		DE
ACC6103XXSSC	0/230F x 0.5F, Safety Coated	\$125.22		DE
ACC6123SSC	0/300F x 2.0F, Safety Coated	\$56.60		DE
ACC6013SC	-50/+50C x 1.0C, Safety Coated	\$59.34		DE
ACC6033SC	-100/+50C X 1.0C, Safety Coated	\$95.47		DE
ACC10113BLSSC	-1/110C x 0.2C, Safety Coated	\$125.22		DE

PRECISION RED-SPIRIT FLUOROPOLYMER-COATED THERMOMETERS - ToI IMMERSION

ACC609SSC	-35/50C x 1.0C, Safety Coated	\$56.60		DE
ACC611SSC	-20/110C x 1.0C, Safety Coated	\$56.60		DE
ACC613SSC	-20/150C x 1.0C , Safety Coated	\$56.60		DE
ACC608SSC	-30/120F x 1.0F, Safety Coated	\$56.60		DE
ACC610SSC	0/230F x 2.0F, Safety Coated	\$56.60		DE
ACC612SSC	0/300F x 2.0F, Safety Coated	\$56.60		DE
ACC601SC	-50/+50C x 1.0C, Safety Coated	\$59.34		DE
ACC603SC	-100/+50C x 1.0C, Safety Coated	\$95.47		DE
ACC605SC	-200/+30C x 1.0C, Safety Coated	\$150.54		DE

"ACCUSAFE" BIODEGRADABLE BLUE-SPIRIT LABORATORY THERMOMETER - PARTIAL IMMERSION

<u>CAT. NO.</u>	<u>HTS Code 9025.11.4000</u>	<u>LIST PRICE</u>	<u>Haz(Hg)</u>	<u>COO</u>
B50CW3BLS	-35/50C x 1.0C	\$18.35		DE, GB
ACC1103BLS	-20/110C x 1.0C	\$17.07		CN, DE, GB
ACC1503BLS	-20/150C x 1.0C	\$18.35		CN, DE, GB
ACC2603BLS	-10/260C x 1.0C	\$31.02		GB
B260CW3BLS16	-10/260C x 1.0C	\$31.02		GB
B330CW3BRN	-10/330C x 2.0C	\$32.56		GB
B230FW3BLS	0/230F x 2.0F	\$17.07		DE, GB
B3003BLS	0/300F x 2.0F	\$18.35		GB
B500FW3BLS	+20/500F x 2.0F	\$31.02		GB
ACC10033BLS	-1/51C x 0.2C	\$91.86		DE
ACC10053BLS	-1/101C x 0.2C	\$121.08		DE
ACC10073BLS	-1/201C x 0.5C	\$129.74		DE
ACC81102BLS	-10/110C x 1.0C	\$15.78		CN, GB
ACC81502BLS	-10/150C x 1.0C	\$15.78		CN, GB

"ACCUSAFE" BIODEGRADABLE BLUE-SPIRIT LABORATORY THERMOMETER - ToI IMMERSION

B50CWBL	-35/50C x 1.0C	\$18.35		DE, GB
ACC110BLS	-20/110C x 1.0C	\$17.07		CN, DE, GB
ACC150BLS	-20/150C x 1.0C	\$18.35		CN, DE, GB
ACC210BLS	-10/210C x 1.0C	\$31.02		DE, GB
B230FWBL	0/230F x 2.0F	\$17.07		DE, GB
B300FWBL	0/300F x 2.0F	\$18.35		DE, GB
B1050CWBL	-10/50C x 0.5C	\$39.42		DE, GB
ACC1003BLS	-1/51C x 0.2C	\$91.86		DE
ACC1005BLS	-1/101C x 0.2C	\$123.43		DE
ACC1007BLS	-1/201C x 0.5C	\$129.74		DE
ACC8110BLS	-10/110C x 1.0C	\$15.78		DE
ACC8150BLS	-10/150C x 1.0C	\$15.78		DE
ACC61050BLS	-10/+50C x 1.0C	\$14.62		DE

"ACCUSAFE" BIODEGRADABLE BLUE-SPIRIT FLUORPOLYMER-COATED - PARTIAL IMMERSION

B50CW3BLSSC	-35/50C x 1.0C	\$32.12		DE, GB
-------------	----------------	---------	--	-----------

<u>CAT. NO.</u>	<u>HTS Code 9025.11.4000</u>	<u>LIST PRICE</u>	<u>Haz(Hg)</u>	<u>COO</u>
ACC1103BLSSC	-20/110C x 1.0C	\$29.48		CN, DE, GB
ACC1503BLSSC	-20/150C x 1.0C	\$32.12		CN, DE, GB
ACC2603BLSSC	-10/260C x 1.0C	\$53.50		DE, GB
B260CW3BL16SC	-10/260C x 1.0C	\$57.72		DE, GB
B230FW3BLSSC	0/230F x 2.0F	\$29.48		GB
B3003BLSSC	0/300F x 2.0F	\$32.12		GB
B5003BLSSC	+20/500F x 2.0F	\$45.23		GB
B10503BLSC	-10/50C X 0.5C	\$51.82		GB
ACC10033BLSC	-1/51C x 0.2C	\$108.40		GB
ACC10053BLSC	-1/101C x 0.2C	\$139.63		DE
ACC10073BLSC	-1/201C x 0.5C	\$146.41		DE
B8TC1102BLS	-10/110C x 1.0C	\$25.98		DE, GB
B8TC1502BLS	-10/150C x 1.0C	\$25.98		DE, GB
<u>"ACCUSAFE" BIODEGRADABLE BLUE-SPIRIT FLUORPOLYMER-COATED - To1 IMMERSION</u>				
B50CWBLSSC	-35/50C x 1.0C	\$32.12		DE, GB
ACC110BLSSC	-20/110C x 1.0C	\$29.48		CN, DE, GB
ACC150BLSSC	-20/150C x 1.0C	\$32.12		CN, DE, GB
ACC210BLSSC	-10/210C x 1.0C	\$45.23		DE, GB
B230FWBLSSC	0/230F x 2.0F	\$29.48		DE, GB
B300FWBLSSC	0/300F x 2.0F	\$32.12		DE, GB
ACC1050BLSSC	-10/50C x 0.5C	\$51.82		GB
ACC1003BLSSC	-1/51C x 0.2C	\$108.40		GB
ACC1005BLSSC	-1/101C x 0.2C	\$139.63		DE
ACC1007BLSSC	-1/201C x 0.5C	\$146.41		DE
B8TC110BLS	-10/110C x 1.0C	\$25.98		DE, GB
B8TC150BLS	-10/150C x 1.0C	\$25.98		DE, GB
<u>"ACCUTHERM" 10/30 JOINTED THERMOMETERS - RED SPIRIT</u>				
ACCST03110S	-10/110C x 1.0C, 76mm Imm	\$79.35		DE
ACCST04110S	-10/110C x 1.0C, 100mm Imm	\$71.85		DE
ACCST01150S	-10/150C x 1.0C, 25mm Imm	\$82.31		DE
ACCST02150S	-10/150C x 1.0C, 51mm Imm	\$82.31		DE
ACCST03150S	-10/150C x 1.0C, 76mm Imm	\$85.42		DE

<u>CAT. NO.</u>	<u>HTS Code 9025.11.4000</u>	<u>LIST PRICE</u>	<u>Haz(Hg)</u>	<u>COO</u>
ACCST04150S	-10/150C x 1.0C, 100mm Imm	\$86.88		DE
ACCST05150S	-10/150C x 1.0C, 125mm Imm	\$81.15		DE
ACCST06150S	-10/150C x 1.0C, 150mm Imm	\$92.97		DE
ACCST01250S	-10/250C x 1.0C, 25mm Imm	\$85.42		DE
ACCST02250S	-10/250C x 1.0C, 50mm Imm	\$85.43		DE
ACCST03250S	-10/250C x 1.0C, 76mm Imm	\$92.97		DE
ACCST04250S	-10/250C x 1.0C, 100mm Imm	\$95.18		DE
ACCST05250S	-10/250C x 1.0C, 125mm Imm	\$99.18		DE
ACCST06250S	-10/250C x 1.0C, 150mm Imm	\$102.73		DE
ACCST07½250S	-10/250C x 1.0C, 190mm Imm	\$90.45		DE
ACCST08½250S	-10/250C x 1.0C, 215mm Imm	\$93.16		DE
ACCST09250S	-10/250C x 1.0C, 228mm Imm	\$98.46		DE
ACCST10250S	-10/250C x 1.0C, 250mm Imm	\$95.88		DE
ACCST11250S	-10/250C x 1.0C, 275mm Imm	\$103.48		DE
ACCST012250S	-10/250C x 1.0C, 305mm Imm	\$106.48		DE

"ACCUTHERM" 10/30 JOINTED THERMOMETERS - MERCURY FILLED

ACCST03110	-10/110C x 1.0C, 75mm Imm		UN3506	DE
ACCST03150	-10/150C x 1.0C, 75mm Imm	\$71.85	UN3506	DE
ACCST03250	-10/250C x 1.0C, 75mm Imm	\$89.98	UN3506	DE
ACCST03360	-10/360C x 1.0C, 75mm Imm	\$83.08	UN3506	DE
ACCST04110	-10/110C x 1.0C, 100mm Imm	\$74.57	UN3506	DE
ACCST04250	-10/250C x 1.0C, 100mm Imm	\$80.24	UN3506	DE
ACCST05360	-10/360C x 1.0C, 125mm Imm	\$88.52	UN3506	DE
ACCST06360	-10/360C x 1.0C, 150mm Imm	\$91.37	UN3506	DE

"ACCUTHERM" 10/18 JOINTED THERMOMETERS - RED SPIRIT

ACCST1801150S	-10/150C x 1.0C, 25mm Imm	\$74.57		DE
ACCST1802150S	-10/150C x 1.0C, 51mm Imm	\$74.57		DE
ACCST1803150S	-10/150C x 1.0C, 76mm Imm DISCONTINUED			DE
ACCST1801250S	-10/250C x 1.0C, 25mm Imm	\$79.87		DE
ACCST1802250S	-10/250C x 1.0C, 51mm Imm	\$79.87		DE
ACCST1803250S	-10/250C x 1.0C, 76mm Imm	\$79.87		DE

PYCNOMETER THERMOMETERS, 10/18 JOINTED-BLUE SPIRIT

ACCST011040S	+14/38C x 0.2C, 25mm Immersion	\$148.88		DE
ACCST0451040S	+14/38C x 0.2C, 45mm Immersion	\$127.80		DE

ADJUSTABLE BUSHING THERMOMETERS-MERCURY FILLED

ACCAB2506	-10/250C x 1.0C, UP TO 6" IMMERSION	\$71.85	UN3506	DE
ACCAB25012	-10/250C x 1.0C UP TO 12" IMMERSION	\$85.16	UN3506	DE
ACCAB3606	-10/360C x 1.0C, UP TO 6" IMMERSION	\$85.16	UN3506	DE

"ACCUTHERM" VARIABLE IMMERSION THERMOMETERS - SPIRIT

ACC20150615S	-20/150C x 1.0C, 6" to 15" Immersion	\$98.46		DE
ACC20150615SSC	-20/150C x 1.0C, 6" to 15" Immersion, Safety Coated	\$119.78		DE
ACC1025037S	-10/250C x 1.0C, 3" to 7" Immersion`	\$82.57		DE

<u>CAT. NO.</u>	<u>HTS Code 9025.11.4000</u>	<u>LIST PRICE</u>	<u>Haz(Hg)</u>	<u>COO</u>
ACC1025037SSC	*-10/250C x 1.0C, 3" to 7" Immersion, Safety Coated	\$101.18		DE
ACC10250612S	-10/250C x 1.0C, 6" to 12" Immersion	\$93.16		DE
ACC10250612SSC	*-10/250C x 1.0C, 6" to 12" Immersion, Safety Coated	\$111.78		DE
ACC102501224S	-10/250C x 1.0C, 12" to 24" Immersion	\$111.78		DE
ACC102501224SSC	-10/250C x 1.0C, 12" to 24" Immersion, Safety Coated	\$138.39		DE

"ACCUTHERM" DEEP IMMERSION THERMOMETERS - SPIRIT

ACC100506S	*-100/50C x 1.0C, 6" Immersion	\$100.79		DE
ACC1005012S	*-100/50C x 1.0C, 12" Immersion	\$111.13		DE
ACC201505S	-20/150C x 1.0C, 5" Immersion	\$63.82		DE
ACC201505SSC	-20/150C x 1.0C, 5" Immersion, Safety Coated	\$79.87		DE
ACC201506S	-20/150C x 1.0C, 6" Immersion	\$63.82		DE
ACC201506SSC	*-20/150C x 1.0C, 6" Immersion, Safety Coated	\$79.87		DE
ACC2015012S	-20/150C x 1.0C, 12" Immersion	\$71.85		DE
ACC2015012SSC	-20/150C x 1.0C, 12" Immersion, Safety Coated	\$85.16		DE
ACC102506S	-10/250C x 1.0C, 6" Immersion`	\$66.55		DE
ACC102506SSC	-10/250C x 1.0C, 6" Immersion, Safety Coated	\$82.57		DE
ACC1025012S	10/250C x 1.0C, 12" Immersion	\$74.57		DE
ACC1025012SSC	-10/250C x 1.0C, 12" Immersion, Safety Coated	\$90.45		DE
ACC1025018S	10/250C x 1.0C, 18" Immersion	\$79.87		DE
ACC1025018SSC	10/250C x 1.0C, 18" Immersion, Safety Coated	\$98.46		DE

"ACCUTHERM" BULB IMMERSION THERMOMETERS - SPIRIT

ACCSM15020S	*-10/150C x 1.0C, Bulb Immersion, Length 150mm	\$50.54		DE
ACCSM25020S	*-10/250C x 1.0C, Bulb Immersion, Length 225mm	\$55.95		DE
ACC5990BLS	30/100C x 0.5C, 30mm Immersion, Length 150mm	\$50.54		DE

"ACCUTHERM" VARIABLE IMMERSION THERMOMETERS - MERCURY

ACC102501224	-10/250C x 1.0C, 12" to 24" Immersion	\$111.78	UN3506	DE
ACC102501224SC	-10/250C x 1.0C, 12" to 24" Immersion, Safety Coated	\$138.39	UN3506	DE

"ACCUTHERM" DEEP IMMERSION THERMOMETERS - MERCURY

ACC102506	10/250C x 1.0C, 6" Immersion	\$53.25	UN3506	DE
ACC102506SC	10/250C x 1.0C, 6" Immersion, Safety Coated	\$69.26	UN3506	DE
ACC102509	10/250C x 1.0C, 9" Immersion	\$58.54	UN3506	DE
ACC102509SC	10/250C x 1.0C, 9" Immersion, Safety Coated	\$74.57	UN3506	DE
ACC1025012	10/250C x 1.0C, 12" Immersion	\$73.20	UN3506	DE
ACC1025012SC	10/250C x 1.0C, 12" Immersion, Safety Coated	\$79.87	UN3506	DE
ACC103606	10/360C x 1.0C, 6" Immersion	\$95.89	UN3506	DE
ACC1036012	10/360C x 1.0C, 12" Immersion	\$103.77	UN3506	DE

"ACCUTHERM" Safety THERMOMETER ADAPTERS - INLET

ACC87011	10/18 JOINT	\$45.62		DE
ACC87013	12/30 JOINT	\$45.62		DE

PRECISION SPIRIT-FILLED ARMORED THERMOMETERS, 76mm Immersion

ACC516SA	-30/120F x 1.0F, 12" ARMOR	\$74.57		DE
----------	----------------------------	---------	--	----

<u>CAT. NO.</u>	<u>HTS Code 9025.11.4000</u>	<u>LIST PRICE</u>	<u>Haz(Hg)</u>	<u>COO</u>
ACC517SA	-35/50C x 1.0C, 12" ARMOR	\$74.57		DE
ACC518SA	0/230 F x 2.0F, 12" ARMOR	\$71.85		DE
ACC519SA	-20/110C x 1.0C, 12" ARMOR	\$71.85		DE
ACC521SA	-20/150C x 1.0C, 12" ARMOR	\$71.85		DE
ACC522SA	0/300F x 2.0F, 12" ARMOR	\$71.85		DE
ACC524SA	+20/400F x 2.0F, 12" ARMOR	\$77.15		DE
ACC525SA	-10/200C x 1.0C, 12" ARMOR	\$74.57		DE
ACC526SA	+30/500F x 2.0F, 14" ARMOR	\$74.57		DE
ACC527SA	-10/260C x 1.0C, 16" ARMOR	\$80.24		DE

PRECISION SPIRIT-FILLED ARMORED THERMOMETER REFILL, 76mm Immersion

ACC516S	-30/120F x 1.0F, 12"	\$44.34		DE
ACC517S	-35/50C x 1.0C, 12"	\$44.34		DE
ACC518S	0/230 F x 2.0F, 12"	\$44.34		DE
ACC519S	-20/110C x 1.0C, 12"	\$44.34		DE
ACC521S	-20/150C x 1.0C, 12"	\$44.34		DE
ACC522S	0/300F x 2.0F, 12"	\$44.34		DE
ACC524S	+20/400F x 2.0F, 12"	\$47.03		DE
ACC525S	-10/200C x 1.0C, 12"	\$47.03		DE
ACC526S	+30/500F x 2.0F, 14"	\$47.03		DE
ACC527S	-10/260C x 1.0C, 16"	\$59.34		DE

PRECISION MERCURY-FILLED ARMORED THERMOMETERS, 76mm Immersion

ACC516A	-30/120F x 1.0F, 12" ARMOR	\$74.57	UN3506	DE
ACC518A	0/230F x 2.0F, 12" ARMOR	\$71.85	UN3506	DE
ACC520A	0/300F x 2.0F, 12" ARMOR	\$71.85	UN3506	DE
ACC524A	+20/400F x 2.0F, 14" ARMOR	\$77.15	UN3506	DE
ACC526A	+20/500F x 2.0F, 16" ARMOR	\$80.24	UN3506	DE
ACC517A	-35/50C x 1.0C, 12" ARMOR	\$74.57	UN3506	DE
ACC519A	-20/105C x 1.0C, 12" ARMOR	\$71.85	UN3506	DE
ACC521A	-20/150C x 1.0C, 12" ARMOR	\$71.85	UN3506	DE
ACC525A	-2/205C x 1.0C, 14" ARMOR	\$77.15	UN3506	DE
ACC527A	-2/260C x 1.0C, 14" ARMOR	\$80.24	UN3506	DE
ACC529A	-2/315C x 1.0C, 16" ARMOR	\$83.08	UN3506	DE
ACC531A	-2/400C x 1.0C, 16" ARMOR	\$88.64	UN3506	DE

PRECISION MERCURY-FILLED ARMORED THERMOMETER REFILLS 76mm Immersion

ACC516	-30/120F x 1.0F	\$47.03	UN3506	DE
ACC518	0/230F x 2.0F	\$44.34	UN3506	DE
ACC520	0/300F x 2.0F	\$44.34	UN3506	DE
ACC524	+20/400F x 2.0F	\$49.76	UN3506	DE
ACC517	-35/50C x 1.0C	\$47.03	UN3506	DE
ACC519	-20/105C x 1.0C	\$44.34	UN3506	DE
ACC521	-20/150C x 1.0C	\$44.34	UN3506	DE
ACC523	-2/150 x 0.5C	\$52.60	UN3506	DE
ACC525	-2/205C x 1.0C	\$49.87	UN3506	DE

<u>CAT. NO.</u>	<u>HTS Code 9025.11.4000</u>	<u>LIST PRICE</u>	<u>Haz(Hg)</u>	<u>COO</u>
ACC527	-2/260C x 1.0C	\$52.60	UN3506	DE
ACC529	-2/315C x 1.0C	\$55.33	UN3506	DE
ACC531	-2/400C x 1.0C	\$58.01	UN3506	DE

"ACCUTHERM" PRECISION POCKET GLASS THERMOMETERS - CLOSED CASE - SPIRIT

ACC532PXS	+30/120F x 1.0F, Tol Imm, 160mm length	\$39.94		DE
ACC532PS	-30/120F x 2.0F, Tol Imm, 160mm length	\$39.94		DE
ACC534PS	0/220F x 2.0F, Tol Imm, 160mm length	\$39.94		DE
ACC533PS	-35/50C x 1.0C, Tol Imm, 160mm length	\$39.94		DE
ACC535PS	-15/105C x 1.0C, Tol Imm, 160mm length	\$38.77		DE

"ACCUTHERM" PRECISION POCKET GLASS THERMOMETERS - OPEN FACE ARMOR CASE - SPIRIT

ACC532XSA	+30/120F x 1.0F, Tol Imm, 160mm length	\$45.25		DE
ACC532SA	-30/120F x 2.0F, Tol Imm, 160mm length	\$39.94		DE
ACC534SA	0/220F x 2.0F, Tol Imm, 160mm length	\$39.94		DE
ACC533SA	-35/50C x 1.0C, Tol Imm, 160mm length	\$39.94		DE
ACC535SA	-15/105C x 1.0C, Tol Imm, 160mm length	\$39.94		DE

"ACCUTHERM" PRECISION POCKET GLASS REFILL THERMOMETERS - SPIRIT

ACC532XS	+30/120F x 1.0F, Tol Imm, 150mm length	\$29.34		DE
ACC532S	-30/120F x 2.0F, Tol Imm, 150mm length	\$29.34		DE
ACC534S	0/220F x 2.0F, Tol Imm, 150mm length	\$29.34		DE
ACC533S	-35/50C x 1.0C, Tol Imm, 150mm length	\$29.34		DE
ACC535S	-15/105C x 1.0C, Tol Imm, 150mm length	\$29.34		DE

ASPHALT POCKET THERMOMETERS, MERCURY

ACC626BT	+100/450F x 5.0F, Tol Imm. 150mm length	\$42.66	UN3506	DE
ACC626P	+100/450F x 5.0F, Tol Imm. 150mm length, w/metal pocket case	\$50.54	UN3506	DE
ACC626A	+100/450F x 5.0F, Tol Imm. 150mm length, w/metal armor case	\$50.54	UN3506	DE

MAXIMUM REGISTERING THERMOMETERS IN POCKET CASE - MERCURY - Tol IMMERSION

ACC540P	-30/120F x 2.0F, 160mm length	\$83.75	UN3506	DE
ACC542P	0/220F x 2.0F, 160mm length	\$83.75	UN3506	DE
ACC544P	+50/300F x 2.0F, 160mm length	\$83.75	UN3506	DE
ACC545P	-20/150C x 1.0C, 160mm length	\$83.75	UN3506	DE
ACC546P	50/400F x 5F, 160mm length	\$83.75	UN3506	DE
ACC547P	0/200C x 2.0C, 150mm length	\$86.49	UN3506	DE
ACC548P	+50/500F x 5.0F, 165mm length	\$83.75	UN3506	DE
ACC549P	0/260C x 2.0C, 165mm length	\$82.57	UN3506	DE
ACC555P	-10/200C x 1.0C, 305mm length	\$115.79	UN3506	DE
ACC5050P	0/220F x 2.0F, 160mm length	\$83.75	UN3506	DE
ACC5054P	+50/300F x 2.0F, 160mm length	\$83.75	UN3506	DE
ACC5058P	+50/400F x 5.0F, 160mm length	\$83.75	UN3506	DE
ACC5062P	+200/500F x 5.0F, 160mm length	\$83.75	UN3506	DE

MAXIMUM REGISTERING THERMOMETERS IN OA CASE - MERCURY - Tol IMMERSION

ACC540A	-30/120F x 2.0F, 160mm length	\$83.75	UN3506	DE
---------	-------------------------------	---------	--------	----

<u>CAT. NO.</u>	<u>HTS Code 9025.11.4000</u>	<u>LIST PRICE</u>	<u>Haz(Hg)</u>	<u>COO</u>
ACC542A	0/220F x 2.0F, 160mm length	\$83.75	UN3506	DE
ACC544A	+50/300F x 2.0F, 160mm length	\$83.75	UN3506	DE
ACC545A	-20/150C x 1.0C, 160mm length	\$83.75	UN3506	DE
ACC547A	0/200C x 2.0C, 150mm length	\$83.75	UN3506	DE
ACC549A	0/260C x 2.0C, 165mm length	\$83.75	UN3506	DE
ACC555A	-10/200C x 1.0C, 305mm length	\$115.79	UN3506	DE
ACC5050A	0/220F x 2.0F, 160mm length	\$83.75	UN3506	DE
ACC5054A	+50/300F x 2.0F, 160mm length	\$83.75	UN3506	DE
ACC5058A	+50/400F x 5.0F, 160mm length	\$83.75	UN3506	DE
ACC5062A	+200/500F x 5.0F, 160mm length	\$83.75	UN3506	DE
S1141741SCA	-100/200C x 1.0C, 300mm Length, Safety Coated	\$127.80	UN3506	DE

MAXIMUM REGISTERING REFILL THERMOMETERS - MERCURY - Tol IMMERSION

ACC540	-30/120F x 2.0F, 150mm length	\$69.26	UN3506	DE
ACC542	0/220F x 2.0F, 150mm length	\$69.26	UN3506	DE
ACC544	+50/300F x 2.0F, 150mm length	\$69.26	UN3506	DE
ACC545	-20/150C x 1.0C, 150mm length	\$69.26	UN3506	DE
ACC547	0/200C x 2.0C, 150mm length	\$76.73	UN3506	DE
ACC549	0/260C x 2.0C, 140mm length	\$74.57	UN3506	DE
ACC555	-10/200C x 1.0C, 300mm length	\$102.47	UN3506	DE
ACC5050	0/220F x 2.0F, 160mm length	\$69.26	UN3506	DE
ACC5054	+50/300F x 2.0F, 160mm length	\$69.26	UN3506	DE
ACC5058	+50/400F x 5.0F, 160mm length	\$69.26	UN3506	DE
ACC5062	+200/500F x 5.0F, 160mm length	\$69.26	UN3506	DE
S1141741	0/200C x 1.0C, 300mm Length	\$106.48	UN3506	DE

FLOATING DAIRY THERMOMETER

ACC0610CS	0/220F x 2.0F, Tol imm, Red Spirit	\$45.16		US
-----------	------------------------------------	---------	--	----

MIN-MAX THERMOMETER

ACC710BLS	-40/50C x1.0C & -40/120F x 2.0F Push Button Reset Blue Spirit	\$47.95		CN
-----------	---	---------	--	----

REFRIGERATOR/FREEZER THERMOMETERS

BFZ80	-40/80F & -40/30C	\$18.62		CN
-------	-------------------	---------	--	----

WALL THERMOMETERS

ACC90121	WHITE ALUMINUM BACK -40/50C & -40/120F, (8 1/2")	\$16.69		US
----------	--	---------	--	----

<u>CAT. NO.</u>	<u>DESCRIPTION</u>	<u>HTS Code 9025.80.2000</u>	<u>LIST PRICE</u>	<u>Haz (Hg)</u>	<u>COO</u>
<u>PLASTIC HYDROMETERS, POLYCARBONATE</u>					
<u>SPECIFIC GRAVITY, PLASTIC HYDROMETERS</u>					
ACC8539PC	0.920/1.080 x 0.002		\$112.00		DE
ACC8540PC	1.020/1.220 x 0.002* new range		\$104.87		DE
ACC8541PC	1.200/1.420 x 0.002		\$104.87		DE
ACC8542PC	1.400/1.620 x 0.002		\$104.87		DE
<u>SALT BRINE PLASTIC HYDROMETERS</u>					
ACC8511PC	0/100% x 1.0%, Saturation of Salt		\$104.87		DE
ACC85112PC	0/26.4% x 0.2%, Percent by Weight		\$104.87		DE
<u>BRIX, SUGAR PLASTIC HYDROMETERS</u>					
ACC8523PC	0/35% x 0.5%		\$104.87		DE
ACC8595PC	0/50% x 0.5%		\$104.87		DE
ACC8515PC	0/15% x 0.2%		\$104.87		DE
ACC8530PC	15/30% x 0.2%		\$104.87		DE
ACC8575PC	60/75% x 0.2%		\$104.87		DE
ACC8590PC	75/90% x 0.2%		\$104.87		DE
<u>BAUME PLASTIC HYDROMETERS, HEAVIER THAN WATER</u>					
ACC8501PC	0/12° x 0.2°		\$104.87		DE
ACC8502PC	9/21° x 0.2°		\$104.87		DE
ACC8503PC	19/31° x 0.2°		\$104.87		DE
ACC8504PC	29/41° x 0.2°		\$104.87		DE
ACC8505PC	39/51° x 0.2°		\$104.87		DE
ACC8506PC	49/61° x 0.2°		\$104.87		DE
ACC8507PC	59/71° x 0.2°		\$104.87		DE
ACC8525PC	0/25° x 0.2°		\$104.87		DE
<u>ASTM / API HYDROMETERS - PLAIN - 330mm</u>					
GW001H	ASTM 1H, -1/11 x 0.1, 330mm		\$35.39		US
GW002H	ASTM 2H, 9/21 x 0.1, 330mm		\$35.39		US
GW003H	ASTM 3H, 19/31 x 0.1, 330mm		\$35.39		US
GW004H	ASTM 4H, 29/41 x 0.1, 330mm		\$35.39		US
GW005H	ASTM 5H, 39/51 x 0.1, 330mm		\$35.39		US
GW006H	ASTM 6H, 49/61 x 0.1, 330mm		\$35.39		US
GW007H	ASTM 7H, 59/71 x 0.1, 330mm		\$35.39		US
GW008H	ASTM 8H, 69/81 x 0.1, 330mm		\$35.39		US
GW009H	ASTM 9H, 79/91 x 0.1, 330mm		\$35.39		US
GW010H	ASTM 10H, 89/101 x 0.1, 330mm		\$35.39		US
GW011H	ASTM 11H, 37/49 x 0.1, 330mm		\$35.39		US
GW012H	ASTM 12H, 64/76 x 0.1, 330mm		\$35.39		US

<u>CAT. NO.</u>	<u>DESCRIPTION</u>	<u>HTS Code 9025.80.2000</u>	<u>LIST PRICE</u>	<u>Haz (Hg)</u>	<u>COO</u>
<u>ASTM / API HYDROMETERS - PLAIN - 163mm</u>					
GW021H	ASTM 21H, 0/6 x 0.1, 163mm		\$35.39		US
GW022H	ASTM 22H, 5/11 x 0.1, 163mm		\$35.39		US
GW023H	ASTM 23H, 10/16 x 0.1, 163mm		\$35.39		US
GW024H	ASTM 24H, 15/21 x 0.1, 163mm		\$35.39		US
GW025H	ASTM 25H, 20/36 x 0.1, 163mm		\$35.39		US
GW026H	ASTM 26H, 25/31 x 0.1, 163mm		\$35.39		US
GW027H	ASTM 27H, 30/36 x 0.1, 163mm		\$35.39		US
GW028H	ASTM 28H, 35/41 x 0.1, 163mm		\$35.39		US
GW029H	ASTM 29H, 40/46 x 0.1, 163mm		\$35.39		US
GW030H	ASTM 30H, 45/51 x 0.1, 163mm		\$35.39		US
GW031H	ASTM 31H, 50/56 x 0.1, 163mm		\$35.39		US
GW032H	ASTM 32H, 55/61 x 0.1, 163mm		\$35.39		US
GW033H	ASTM 33H, 60/66 x 0.1, 163mm		\$35.39		US
GW034H	ASTM 34H, 65/71 x 0.1, 163mm		\$35.39		US
GW035H	ASTM 35H, 70/76 x 0.1, 163mm		\$35.39		US
GW036H	ASTM 36H, 75/81 x 0.1, 163mm		\$35.39		US
GW037H	ASTM 37H, 80/86 x 0.1, 163mm		\$35.39		US
GW040H	ASTM 40H, 95/101 x 0.1, 163mm		\$35.39		US
<u>ASTM/API COMBINED W/Hg THERMOMETER IN BODY HYDROMETERS, 0/150F - 380MM</u>					
GW041HL	ASTM 41HL, 15/23 x 0.1 & 0/150F Hg, 380mm		\$70.76	UN3506	US
GW042HL	ASTM 42HL, 22/30 x 0.1 & 0/150F Hg, 380mm		\$70.76	UN3506	US
GW043HL	ASTM 43HL, 29/37 x 0.1 & 0/150F Hg, 380mm		\$70.76	UN3506	US
GW044HL	ASTM 44HL, 36/44 x 0.1 & 0/150F Hg, 380mm		\$70.76	UN3506	US
GW045HL	ASTM 45HL, 43/51 x 0.1 & 0/150F Hg, 380mm		\$70.76	UN3506	US
GW051HL	ASTM 51HL, -1/11 x 0.1 & 0/150F Hg, 380mm		\$70.76	UN3506	US
GW052HL	ASTM 52HL, 9/21 x 0.1 & 0/150F Hg, 380mm		\$70.76	UN3506	US
GW053HL	ASTM 53HL, 19/31 x 0.1 & 0/150F Hg, 380mm		\$70.76	UN3506	US
GW054HL	ASTM 54HL, 29/41 x 0.1 & 0/150F Hg, 380mm		\$70.76	UN3506	US
GW055HL	ASTM 55HL, 39/51 x 0.1 & 0/150F Hg, 380mm		\$70.76	UN3506	US
GW056HL	ASTM 56HL, 49/61 x 0.1 & 0/150F Hg, 380mm		\$70.76	UN3506	US
GW057HL	ASTM 57HL, 59/71 x 0.1 & 0/150F Hg, 380mm		\$70.76	UN3506	US
GW058HL	ASTM 58HL, 69/81 x 0.1 & 0/150F Hg, 380mm		\$70.76	UN3506	US
GW059HL	ASTM 59HL, 79/91 x 0.1 & 0/150F Hg, 380mm		\$70.76	UN3506	US
GW060HL	ASTM 60HL, 89/101 x 0.1 & 0/150F Hg, 380mm		\$70.76	UN3506	US
<u>ASTM/API COMBINED W/SPIRIT FILLED THERMOMETER IN BODY HYDROMETERS, 0/150F - 380MM</u>					
GW041HLS	ASTM 41HL, 15/23 x 0.1 & 0/150F Spirit, 380mm		\$70.76		US
GW042HLS	ASTM 42HL, 22/30 x 0.1 & 0/150F Spirit, 380mm		\$70.76		US
GW043HLS	ASTM 43HL, 29/37 x 0.1 & 0/150F Spirit, 380mm		\$70.76		US
GW044HLS	ASTM 44HL, 36/44 x 0.1 & 0/150F Spirit, 380mm		\$70.76		US
GW045HLS	ASTM 45HL, 43/51 x 0.1 & 0/150F Spirit, 380mm		\$70.76		US
GW051HLS	ASTM 51HL, -1/11 x 0.1 & 0/150F Spirit, 380mm		\$70.76		US
GW052HLS	ASTM 52HL, 9/21 x 0.1 & 0/150F Spirit, 380mm		\$70.76		US
GW053HLS	ASTM 53HL, 19/31 x 0.1 & 0/150F Spirit, 380mm		\$70.76		US
GW054HLS	ASTM 54HL, 29/41 x 0.1 & 0/150F Spirit, 380mm		\$70.76		US
GW055HLS	ASTM 55HL, 39/51 x 0.1 & 0/150F Spirit, 380mm		\$70.76		US
GW056HLS	ASTM 56HL, 49/61 x 0.1 & 0/150F Spirit, 380mm		\$70.76		US
GW057HLS	ASTM 57HL, 59/71 x 0.1 & 0/150F Spirit, 380mm		\$70.76		US
GW058HLS	ASTM 58HL, 69/81 x 0.1 & 0/150F Spirit, 380mm		\$70.76		US
GW059HLS	ASTM 59HL, 79/91 x 0.1 & 0/150F Spirit, 380mm		\$70.76		US

<u>CAT. NO.</u>	<u>DESCRIPTION</u>	<u>HTS Code</u> 9025.80.2000	<u>LIST PRICE</u>	<u>Haz (Hg)</u>	<u>COO</u>
<u>ASTM/API COMBINED W/Hg THERMOMETER IN BODY</u>					
<u>HYDROMETERS, 30/180F - 380MM</u>					
GW042HM	ASTM 42HM, 22/30 x 0.1 & 30/180F Hg, 380mm		\$70.76	UN3506	US
GW043HM	ASTM 43HM, 29/37 x 0.1 & 30/180F Hg, 380mm		\$70.76	UN3506	US
GW044HM	ASTM 44HM, 36/44 x 0.1 & 30/180F Hg, 380mm		\$70.76	UN3506	US
GW052HM	ASTM 52HM, 9/21 x 0.1 & 30/180F Hg, 380mm		\$70.76	UN3506	US
GW053HM	ASTM 53HM, 19/31 x 0.1 & 30/180F Hg, 380mm		\$70.76	UN3506	US
GW054HM	ASTM 54HM, 29/41 x 0.1 & 30/180F Hg, 380mm		\$70.76	UN3506	US
GW055HM	ASTM 55HM, 39/51 x 0.1 & 30/180F Hg, 380mm		\$70.76	UN3506	US
GW056HM	ASTM 56HM, 49/61 x 0.1 & 30/180F Hg, 380mm		\$70.76	UN3506	US
GW057HM	ASTM 57HM, 59/71 x 0.1 & 30/180F Hg, 380mm		\$70.76	UN3506	US
GW060HM	ASTM 60HM, 89/101 x 0.1 & 30/180F Hg, 380mm		\$70.76	UN3506	US
<u>ASTM/API COMBINED W/SPIRIT-FILLED THERMOMETER IN</u>					
<u>BODY HYDROMETERS, 30/180F - 380MM</u>					
GW041HMS	ASTM 41HM, 15/23 x 0.1 & 30/180F Spirit, 380mm		\$70.76		US
GW042HMS	ASTM 42HM, 22/30 x 0.1 & 30/180F Spirit, 380mm		\$70.76		US
GW043HMS	ASTM 43HM, 29/37 x 0.1 & 30/180F Spirit, 380mm		\$70.76		US
GW044HMS	ASTM 44HM, 36/44 x 0.1 & 30/180F Spirit, 380mm		\$70.76		US
GW045HMS	ASTM 45HM, 43/51 x 0.1 & 30/180F Spirit, 380mm		\$70.76		US
GW051HMS	ASTM 51HM, -1/11 x 0.1 & 30/180F Spirit, 380mm		\$70.76		US
GW052HMS	ASTM 52HM, 9/21 x 0.1 & 30/180F Spirit, 380mm		\$70.76		US
GW053HMS	ASTM 53HM, 19/31 x 0.1 & 30/180F Spirit, 380mm		\$70.76		US
GW054HMS	ASTM 54HM, 29/41 x 0.1 & 30/180F Spirit, 380mm		\$70.76		US
GW055HMS	ASTM 55HM, 39/51 x 0.1 & 30/180F Spirit, 380mm		\$70.76		US
GW056HMS	ASTM 56HM, 49/61 x 0.1 & 30/180F Spirit, 380mm		\$70.76		US
GW057HMS	ASTM 57HM, 59/71 x 0.1 & 30/180F Spirit, 380mm		\$70.76		US
GW058HMS	ASTM 58HM, 69/81 x 0.1 & 30/180F Spirit, 380mm		\$70.76		US
<u>ASTM/API COMBINED W/Hg THERMOMETER IN BODY</u>					
<u>HYDROMETERS, 60/220F - 380MM</u>					
GW043HH	ASTM 43HH, 29/37 x 0.1 & 60/220F Hg, 380mm		\$70.76	UN3506	US
GW045HH	ASTM 45HH, 43/51 x 0.1 & 60/220F Hg, 380mm		\$70.76	UN3506	US
GW051HH	ASTM 51HH, -1/11 x 0.1 & 60/220F Hg, 380mm		\$70.76	UN3506	US
GW052HH	ASTM 52HH, 9/21 x 0.1 & 60/220F Hg, 380mm		\$70.76	UN3506	US
GW053HH	ASTM 53HH, 19/31 x 0.1 & 60/220F Hg, 380mm		\$70.76	UN3506	US
GW054HH	ASTM 54HH, 29/41 x 0.1 & 60/220F Hg, 380mm		\$70.76	UN3506	US
GW055HH	ASTM 55HH, 39/51 x 0.1 & 60/220F Hg, 380mm		\$70.76	UN3506	US
GW056HH	ASTM 56HH, 49/61 x 0.1 & 60/220F Hg, 380mm		\$70.76	UN3506	US
GW057HH	ASTM 57HH, 59/71 x 0.1 & 60/220F Hg, 380mm		\$70.76	UN3506	US
GW058HH	ASTM 58HH, 69/81 x 0.1 & 60/220F Hg, 380mm		\$70.76	UN3506	US
<u>ASTM/API COMBINED W/Hg THERMOMETER IN STEM</u>					
<u>HYDROMETER, 30/220F - 380MM</u>					
GW071H	ASTM 71H, -1/11 x 0.1 & 30/220F Hg, 380mm		\$87.82	UN3506	US
GW072H	ASTM 72H, 9/21 x 0.1 & 30/220F Hg, 380mm		\$87.82	UN3506	US
GW073H	ASTM 73H, 19/31 x 0.1 & 30/220F Hg, 380mm		\$87.82	UN3506	US
GW074H	ASTM 74H, 29/41 x 0.1 & 30/220F Hg, 380mm		\$87.82	UN3506	US
<u>ASTM/API COMBINED W/ SPIRIT-FILLED THERMOMETER</u>					
<u>IN STEM HYDROMETER, 30/22F - 380MM</u>					
GW071HLS	ASTM 71H, -1/11 x 0.1 & 30/220F Spirit, 380mm		\$87.82		US
GW072HLS	ASTM 72H, 9/21 x 0.1 & 30/220F Spirit, 380mm		\$87.82		US
GW073HLS	ASTM 73H, 19/31 x 0.1 & 30/220F Spirit, 380mm		\$87.82		US
GW074HLS	ASTM 74H, 29/41 x 0.1 & 30/220F Spirit, 380mm		\$87.82		US

<u>CAT. NO.</u>	<u>DESCRIPTION</u>	<u>HTS Code 9025.80.2000</u>	<u>LIST PRICE</u>	<u>Haz (Hg)</u>	<u>COO</u>
<u>ASTM LPG COMBINED IN BODY HYDROMETER, HG, +30/90F - 380mm</u>					
GW101H	ASTM 101H, 0.500/0.650 x 0.001 & 30/90F Hg, 380mm		\$75.99	UN3506	US
<u>ASTM SPECIFIC GRAVITY, PLAIN HYDROMETERS - 330mm</u>					
GW082H	ASTM 82H - 0.650/0.700 x 0.0005, 330mm		\$42.01		US
GW083H	ASTM 83H - 0.700/0.750 x 0.0005, 330mm		\$42.01		US
GW084H	ASTM 84H - 0.750/0.800 x 0.0005, 330mm		\$42.01		US
GW085H	ASTM 85H - 0.800/0.850 x 0.0005, 330mm		\$42.01		US
GW086H	ASTM 86H - 0.850/0.900 x 0.0005, 330mm		\$42.01		US
GW087H	ASTM 87H - 0.900/0.950 x 0.0005, 330mm		\$42.01		US
GW088H	ASTM 88H - 0.950/1.000 x 0.0005, 330mm		\$42.01		US
GW089H	ASTM 89H - 1.000/1.050 x 0.0005, 330mm		\$42.01		US
GW090H	ASTM 90H - 1.050/1.100 x 0.0005, 330mm		\$42.01		US
GW111H	ASTM 111H - 1.000/1.050 x 0.0005, 330mm		\$42.01		US
GW112H	ASTM 112H - 1.050/1.100 x 0.0005, 330mm		\$42.01		US
GW113H	ASTM 113H - 1.100/1.150 x 0.0005, 330mm		\$42.01		US
GW114H	ASTM 114H - 1.150/1.200 x 0.0005, 330mm		\$42.01		US
GW115H	ASTM 115H - 1.200/1.250 x 0.0005, 330mm		\$42.01		US
GW116H	ASTM 116H - 1.250/1.300 x 0.0005, 330mm		\$42.01		US
GW117H	ASTM 117H - 1.300/1.350 x 0.0005, 330mm		\$42.01		US
GW118H	ASTM 118H - 1.350/1.400 x 0.0005, 330mm		\$42.01		US
GW119H	ASTM 119H - 1.400/1.450 x 0.0005, 330mm		\$42.01		US
GW120H	ASTM 120H - 1.450/1.500 x 0.0005, 330mm		\$42.01		US
<u>ASTM SPECIFIC GRAVITY HYDROMETERS, PLAIN FORM 260mm</u>					
GW102H	ASTM 102H - 0.650/0.700 x 0.001, 260mm				US
GW103H	ASTM 103H - 0.700/0.750 x 0.001, 260mm		\$42.01		US
GW104H	ASTM 104H - 0.750/0.800 x 0.001, 260mm		\$42.01		US
GW105H	ASTM 105H - 0.800/0.850 x 0.001, 260mm		\$42.01		US
GW106H	ASTM 106H - 0.850/0.900 x 0.001, 260mm		\$42.01		US
GW107H	ASTM 107H - 0.900/0.950 x 0.001, 260mm		\$42.01		US
GW108H	ASTM 108H - 0.950/1.000 x 0.001, 260mm		\$42.01		US
GW125H	ASTM 125H - 1.000/1.050 x 0.001, 260mm		\$42.01		US
GW126H	ASTM 126H - 1.050/1.100 x 0.001, 260mm		\$42.01		US
GW127H	ASTM 127H - 1.100/1.150 x 0.001, 260mm		\$42.01		US
GW128H	ASTM 128H - 1.150/1.200 x 0.001, 260mm		\$42.01		US
GW129H	ASTM 129H - 1.200/1.250 x 0.001, 260mm		\$42.01		US
GW130H	ASTM 130H - 1.250/1.300 x 0.001, 260mm		\$42.01		US
GW131H	ASTM 131H - 1.300/1.350 x 0.001, 260mm		\$42.01		US
GW132H	ASTM 132H - 1.350/1.400 x 0.001, 260mm		\$42.01		US
GW133H	ASTM 133H - 1.400/1.450 x 0.001, 260mm		\$42.01		US
GW134H	ASTM 134H - 1.450/1.500 x 0.001, 260mm		\$42.01		US
GW135H	ASTM 135H - 1.500/1.550 x 0.001, 260mm		\$42.01		US
GW136H	ASTM 136H - 1.550/1.600 x 0.001, 260mm		\$42.01		US
GW137H	ASTM 137H - 1.600/1.650 x 0.001, 260mm		\$42.01		US
GW138H	ASTM 138H - 1.650/1.700 x 0.001, 260mm		\$42.01		US
GW139H	ASTM 139H - 1.700/1.750 x 0.001, 260mm		\$42.01		US
GW140H	ASTM 140H - 1.750/1.800 x 0.001, 260mm		\$42.01		US
GW141H	ASTM 141H - 1.800/1.850 x 0.001, 260mm		\$42.01		US
GW144H	ASTM 144H - 1.950/2.000 x 0.001, 260mm		\$42.01		US
<u>ASTM SOIL ANALYSIS, PLAIN HYDROMETER, SPECIFIC GRAVITY 280mm</u>					
GW151H	ASTM 151H, 0.995/1.038 x 0.001, 280mm		\$65.55		US
<u>ASTM SOIL ANALYSIS, PLAIN HYDROMETER, GRAMS OF SOLID COLLOIDS 280mm</u>					
GW152H	ASTM 152H, -5/60 grams x 1, 280mm		\$62.87		US

<u>CAT. NO.</u>	<u>DESCRIPTION</u>	<u>HTS Code 9025.80.2000</u>	<u>LIST PRICE</u>	<u>Haz (Hg)</u>	<u>COO</u>
<u>ASTM DENSITY HYDROMETERS, COMBINED</u>					
<u>w/ HG-FILLED THERMOMETER - 20/65C - 330MM</u>					
GW300HL	ASTM 300H, 600/650 x 0.5 & -20/65C, 330mm		\$70.76	UN3506	US
GW301HL	ASTM 301H, 650/700 x 0.5 & -20/65C, 330mm		\$70.76	UN3506	US
GW302HL	ASTM302H, 700/750 x 0.5 & -20/65C, 330mm		\$70.76	UN3506	US
GW303HL	ASTM 303H, 750/800 x 0.5 & -20/65C, 330mm		\$70.76	UN3506	US
GW304HL	ASTM 304H, 800/850 x 0.5 & -20/65C, 330mm		\$70.76	UN3506	US
GW305HL	ASTM 305H, 850/900 x 0.5 & -20/65C, 330mm		\$70.76	UN3506	US
GW306HL	ASTM 306H, 900/950 x 0.5 & -20/65C, 330mm		\$70.76	UN3506	US
GW307HL	ASTM 307H, 950/1000 x 0.5 & -20/65C, 330mm		\$70.76	UN3506	US
GW308HL	ASTM 308H, 1000/1050 x 0.5 & -20/65C, 330mm		\$70.76	UN3506	US
GW309HL	ASTM 309H, 1050/1100 x 0.5 & -20/65C, 330mm		\$70.76	UN3506	US
GW310HL	ASTM 310H, 500/650 x 1.0 & 0/35C, 387mm		\$70.76	UN3506	US
<u>ASTM DENSITY HYDROMETERS, PLAIN - 330mm</u>					
GW311H	ASTM 311H, 600/650 x 0.5, 330mm		\$42.01		US
GW312H	ASTM 312H, 650/700 x 0.5, 330mm		\$42.01		US
GW313H	ASTM313H, 700/750 x 0.5, 330mm		\$42.01		US
GW314H	ASTM 314H, 750/800 x 0.5, 330mm		\$42.01		US
GW315H	ASTM 315H, 800/850 x 0.5, 330mm		\$42.01		US
GW316H	ASTM 316H, 850/900 x 0.5, 330mm		\$42.01		US
GW317H	ASTM 317H, 900/950 x 0.5, 330mm		\$42.01		US
GW318H	ASTM 318H, 950/1000 x 0.5, 330mm		\$42.01		US
GW319H	ASTM 319H, 1000/1050 x 0.5, 330mm		\$42.01		US
GW320H	ASTM 320H, 1050/1100 x 0.5, 330mm		\$42.01		US
<u>API, PLAIN HYDROMETERS - 6"</u>					
GW2500	10/45 x 1.0, 6"		\$39.34		US
GW2501	-1/11 x 0.2, 6"		\$39.34		US
GW2502	9/21 x 0.2, 6"		\$39.34		US
GW2503	19/31 x 0.2, 6"		\$39.34		US
GW2504	29/41 x 0.2, 6"		\$39.34		US
GW2505	39/51 x 0.2, 6"		\$39.34		US
GW2506	49/61 x 0.2, 6"		\$39.34		US
GW2507	59/71 x 0.2, 6"		\$39.34		US
GW2508	69/81 x 0.2, 6"		\$39.34		US
GW2509	79/91 x 0.2, 6"		\$39.34		US
<u>API, COMBINED HYDROMETER W/HG-FILLED THERMOMETER, 0/150F -7.5"</u>					
GW2500C	10/45 x 0.1 & 0/150F Hg, 7.5"		\$68.22	UN3506	US
GW2501C	-1/11 x 0.2 & 0/150F Hg, 7.5"		\$68.22	UN3506	US
GW2502C	9/21 x 0.2 & 0/150F Hg, 7.5"		\$68.22	UN3506	US
GW2503C	19/31 x 0.2 & 0/150F Hg, 7.5"		\$68.22	UN3506	US
GW2504C	29/41 x 0.2 & 0/150F Hg, 7.5"		\$68.22	UN3506	US
GW2505C	39/51 x 0.2 & 0/150F Hg, 7.5"		\$68.22	UN3506	US
GW2506C	49/61 x 0.2 & 0/150F Hg, 7.5"		\$68.22	UN3506	US
GW2507C	59/71 x 0.2 & 0/150F Hg, 7.5"		\$68.22	UN3506	US
GW2508C	69/81 x 0.2 & 0/150F Hg, 7.5"		\$68.22	UN3506	US
<u>API, COMBINED HYDROMETER W/SPIRIT-FILLED THERMOMETER, 0/150F - 7.5"</u>					
GW2500CS	10/45 x 0.1 & 0/150F Spirit, 7.5"		\$68.22		US
GW2501CS	-1/11 x 0.2 & 0/150F Spirit, 7.5"		\$68.22		US
GW2502CS	9/21 x 0.2 & 0/150F Spirit, 7.5"		\$68.22		US
GW2503CS	19/31 x 0.2 & 0/150F Spirit, 7.5"		\$68.22		US
GW2504CS	29/41 x 0.2 & 0/150F Spirit, 7.5"		\$68.22		US
GW2505CS	39/51 x 0.2 & 0/150F Spirit, 7.5"		\$68.22		US

<u>CAT. NO.</u>	<u>DESCRIPTION</u>	<u>HTS Code</u> 9025.80.2000	<u>LIST PRICE</u>	<u>Haz (Hg)</u>	<u>COO</u>
GW2506CS	49/61 x 0.2 & 0/150F Spirit, 7.5"		\$68.22		US
GW2507CS	59/71 x 0.2 & 0/150F Spirit, 7.5"		\$68.22		US
GW2508CS	69/81 x 0.2 & 0/150F Spirit, 7.5"		\$68.22		US
GW2509CS	79/91 x 0.2 & 0/150F Spirit, 7.5"		\$68.22		US

API, COMBINED HYDROMETER W/HG-FILLED THERMOMETER, 20/130F - 7.5"

GW3000	10/50 x 0.5 & 20/130F Hg, 7.5"		\$68.22	UN3506	US
GW3001	25/65 x 0.5 & 20/130F Hg, 7.5"		\$68.22	UN3506	US
GW3002	40/80 x 0.5 & 20/130F Hg, 7.5"		\$68.22	UN3506	US

API, COMBINED HYDROMETERS W/HG-FILLED THERMOMETER, 20/130F - 7.5"

GW5001	45/90 x 1.0 & 20/130F Hg, 7.5"		\$70.76	UN3506	US
GW5003	20/40 x 0.5 & 20/130F Hg, 7.5"		\$70.76	UN3506	US
GW5004	30/50 x 0.5 & 20/130F Hg, 7.5"		\$70.76	UN3506	US
GW5005	40/60 x 0.5 & 20/130F Hg, 7.5"		\$70.76	UN3506	US
GW5006	50/70 x 0.5 & 20/130F Hg, 7.5"		\$70.76	UN3506	US
GW5007	60/80 x 0.5 & 20/130F Hg, 7.5"		\$70.76	UN3506	US
GW5008	70/90 x 0.5 & 20/130F Hg, 7.5"		\$70.76	UN3506	US
GW5009	80/100 x 0.5 & 20/130F Hg, 7.5"		\$70.76	UN3506	US

API, COMBINED HYDROMETERS W/Hg THERMOMETER 0/150F - 15"

GW6000	10/45 x 0.5 & 0/150F Hg, 15"		\$70.76	UN3506	US
GW6000X	0/20 & 0/150F Hg, 15"		\$70.76	UN3506	US
GW6001	45/90 x 0.5 & 0/150F Hg, 15"		\$70.76	UN3506	US
GW6001X	30/90 & 0/150F Hg, 15"		\$70.76	UN3506	US
GW6002	10/30 x 0.2 & 0/150F Hg, 15"		\$70.76	UN3506	US
GW6003	20/40 x 0.2 & 0/150F Hg, 15"		\$70.76	UN3506	US
GW6004	30/50 x 0.2 & 0/150F Hg, 15"		\$70.76	UN3506	US
GW6005	40/60 x 0.2 & 0/150F Hg, 15"		\$70.76	UN3506	US
GW6006	50/70 x 0.2 & 0/150F Hg, 15"		\$70.76	UN3506	US
GW6007	60/80 x 0.2 & 0/150F Hg, 15"		\$70.76	UN3506	US
GW6007X	40/80 & 0/150F Hg, 15"		\$70.76	UN3506	US
GW6008	70/90 x 0.2 & 0/150F Hg, 15"		\$70.76	UN3506	US

API, COMBINED HYDROMETERS W/SPIRIT-FILLED THERMOMETER, 0/150F - 15"

GW6000S	10/45 x 0.5 & 0/150F Spirit, 15"		\$70.76		US
GW6001S	45/90 x 0.5 & 0/150F Spirit, 15"		\$70.76		US
GW6001XS	30/90 & 0/150F Spirit, 15"		\$70.76		US
GW6002S	10/30 x 0.2 & 0/150F Spirit, 15"		\$70.76		US
GW6003S	20/40 x 0.2 & 0/150F Spirit, 15"		\$70.76		US
GW6004S	30/50 x 0.2 & 0/150F Spirit, 15"		\$70.76		US
GW6005S	40/60 x 0.2 & 0/150F Spirit, 15"		\$70.76		US
GW6006S	50/70 x 0.2 & 0/150F Spirit, 15"		\$70.76		US

SPECIFIC GRAVITY HYDROMETERS, PLAIN - 150mm

GW2536	0.700/0.800 x 0.002, 150mm		\$32.83		US
GW2537	0.800/0.900 x 0.002, 150mm		\$32.83		US
GW2538	0.900/1.000 x 0.002, 150mm		\$32.83		US
GW2539	0.700/1.000 x 0.010, 150mm		\$32.83		US
GW2548	1.000/1.225 x 0.005, 150mm		\$32.83		US
GW2549	1.200/1.425 x 0.005, 150mm		\$32.83		US
GW2550	1.400/1.625 x 0.005, 150mm		\$32.83		US
GW2551	1.600/1.825 x 0.005, 150mm		\$32.83		US
GW2552	1.800/2.000 x 0.005, 150mm		\$32.83		US
GW2552A	1.000/2.000 x 0.020, 150mm		\$32.83		US
GW2554X	0.640/0.710 x 0.001, 150mm		\$42.01		US

<u>CAT. NO.</u>	<u>DESCRIPTION</u>	<u>HTS Code 9025.80.2000</u>	<u>LIST PRICE</u>	<u>Haz (Hg)</u>	<u>COO</u>
GW2555X	0.700/0.770 x 0.001, 150mm		\$42.01		US
GW2556X	0.760/0.830 x 0.001, 150mm		\$42.01		US
GW2557X	0.820/0.890 x 0.011, 150mm		\$42.01		US
GW2558X	0.880/0.950 x 0.001, 150mm		\$42.01		US
GW2559X	0.940/1.010 x 0.001, 150mm		\$42.01		US
GW2560X	1.000/1.070 x 0.001, 150mm		\$42.01		US
GW2561X	1.060/1.130 x 0.001, 150mm		\$42.01		US
GW2562X	1.120/1.190 x 0.001, 150mm		\$42.01		US
GW2563X	1.180/1.250 x 0.001, 150mm		\$42.01		US
GW2564X	1.240/1.310 x 0.001, 150mm		\$42.01		US
GW2565X	1.300/1.370 x 0.001, 150mm		\$42.01		US
GW2566X	1.360/1.430 x 0.001, 150mm		\$42.01		US
GW2567X	1.420/1.490 x 0.001, 150mm		\$42.01		US
GW2568X	1.480/1.550 x 0.001, 150mm		\$42.01		US
GW2569X	1.540/1.610 x 0.001, 150mm		\$42.01		US
GW2570X	1.600/1.670 x 0.001, 150mm		\$42.01		US
GW2571X	1.660/1.730 x 0.001, 150mm		\$42.01		US
GW2572X	1.720/1.790 x 0.001, 150mm		\$42.01		US
GW2573X	1.780/1.850 x 0.001, 150mm		\$42.01		US
GW2574X	1.840/1.920 x 0.001, 150mm		\$42.01		US

SPECIFIC GRAVITY HYDROMETERS, PLAIN - 300mm

GW2532	0.695/0.810 x 0.001, 300mm **NEW RANGE 0.700/0.810 x 0.001, 300mm		\$32.83		US
GW2533	0.795/0.910 x 0.001, 300mm **NEW RANGE 0.800/0.910 x 0.001, 300mm		\$32.83		US
GW2534	0.895/1.000 x 0.001, 300mm **NEW RANGE 0.900/1.010 x 0.001, 300mm		\$32.83		US
GW2535	0.650/1.000 x 0.005, 300mm		\$32.83		US
GW2540	1.000/1.220 x 0.002, 300mm		\$32.83		US
GW2541	1.200/1.420 x 0.002, 300mm		\$32.83		US
GW2542	1.400/1.620 x 0.002, 300mm		\$32.83		US
GW2543	1.600/1.820 x 0.002, 300mm		\$32.83		US
GW2544	1.800/2.000 x 0.002, 300mm		\$32.83		US
GW2545	1.000/2.000 x 0.010, 300mm		\$32.83		US
GW2546	1.000/1.600 x 0.005, 300mm		\$32.83		US
GW2547	1.500/2.000 x 0.005, 300mm		\$32.83		US

SPECIFIC GRAVITY HYDROMETERS, COMBINED FORM
w/ SPIRIT THERMOMETER, 20/120F - 380MM

GW2532C	0.695/0.810 x 0.001, 20/120F Spirit, 380mm		\$70.76		US
GW2533C	0.795/0.910 x 0.001, 20/120F Spirit, 380mm		\$70.76		US
GW2534C	0.895/1.000 x 0.001, 20/120F Spirit, 380mm **NEW RANGE 0.900/1.010 x 0.001,		\$70.76		US
GW2535C	0.650/1.000 x 0.010, 20/120F Spirit, 380mm		\$70.76		US
GW2540C	1.000/1.220 x 0.002, 20/120F Spirit, 380mm		\$70.76		US
GW2541C	1.200/1.420 x 0.002, 20/120F Spirit, 380mm		\$70.76		US
GW2542C	1.400/1.620 x 0.002, 20/120F Spirit, 380mm		\$70.76		US
GW2543C	1.600/1.820 x 0.002, 20/120F Spirit, 380mm		\$70.76		US
GW2544C	1.800/2.000 x 0.002, 20/120F Spirit, 380mm		\$70.76		US
GW2545C	1.000/2.000 x 0.010, 20/120F Spirit, 380mm		\$70.76		US
GW2546C	1.000/1.600 x 0.005, 20/120F Spirit, 380mm		\$70.76		US
GW2547C	1.500/2.000 x 0.005, 20/120F Spirit, 380mm		\$70.76		US

<u>CAT. NO.</u>	<u>DESCRIPTION</u>	<u>HTS Code 9025.80.2000</u>	<u>LIST PRICE</u>	<u>Haz (Hg)</u>	<u>COO</u>
<u>SPECIFIC GRAVITY HYDROMETERS, PLAIN - 330mm</u>					
GW2553	0.590/0.660 x 0.0005, 330mm		\$44.54		US
GW2554	0.640/0.710 x 0.0005, 330mm		\$44.54		US
GW2555	0.700/0.770 x 0.0005, 330mm		\$44.54		US
GW2556	0.760/0.830 x 0.0005, 330mm		\$44.54		US
GW2557	0.820/0.890 x 0.0005, 330mm		\$44.54		US
GW2558	0.880/0.950 x 0.0005, 330mm		\$44.54		US
GW2559	0.940/1.010 x 0.0005, 330mm		\$44.54		US
GW2560	1.000/1.070 x 0.0005, 330mm		\$44.54		US
GW2561	1.060/1.130 x 0.0005, 330mm		\$44.54		US
GW2562	1.120/1.190 x 0.0005, 330mm		\$44.54		US
GW2563	1.180/1.250 x 0.0005, 330mm		\$44.54		US
GW2564	1.240/1.310 x 0.0005, 330mm		\$44.54		US
GW2565	1.300/1.370 x 0.0005, 330mm		\$44.54		US
GW2566	1.360/1.430 x 0.0005, 330mm		\$44.54		US
GW2567	1/420/1.490 x 0.0005, 330mm		\$44.54		US
GW2568	1.480/1.550 x 0.0005, 330mm		\$44.54		US
GW2569	1.540/1.610 x 0.0005, 330mm		\$44.54		US
GW2570	1.600/1.670 x 0.0005, 330mm		\$44.54		US
GW2571	1.660/1.730 x 0.0005, 330mm		\$44.54		US
GW2572	1.720/1.790 x 0.0005, 330mm		\$44.54		US
GW2573	1.780/1.850 x 0.0005, 330mm		\$44.54		US
GW2574	1.840/1.920 x 0.0005, 330mm		\$44.54		US
<u>SPECIFIC GRAVITY HYDROMETER SET OF 8, PLAIN - 300mm, w/ case</u>					
GW3244SET	0.700/2.000 x 0.002		\$397.16		US
<u>BAUME HYDROMETERS, PLAIN - 300mm</u>					
GW2575	0/11 x 0.1, 300mm		\$32.83		US
GW2775	-1/11 x 0.1, 300mm		\$32.83		US
GW2576	9/21 x 0.1, 300mm		\$32.83		US
GW2577	19/31 x 0.1, 300mm		\$32.83		US
GW2578	29/41 x 0.1, 300mm		\$32.83		US
GW2579	39/51 x 0.1, 300mm		\$32.83		US
GW2580	49/61 x 0.1, 300mm		\$32.83		US
GW2581	59/71 x 0.1, 300mm		\$32.83		US
GW2582	0/35 x 0.5, 300mm		\$32.83		US
GW2583	35/70 x 0.5, 300mm		\$32.83		US
GW2584	0/50 x 1.0, 300mm		\$32.83		US
GW2585	0/70 x 1.0, 300mm		\$32.83		US
<u>BAUME HYDROMETERS, PLAIN - 150mm</u>					
GW2575X	0/11 x 0.1				US
GW2576X	9/21 x 0.1		\$32.83		US
<u>DOUBLE SCALE SPECIFIC GRAVITY/BAUME HYDROMETERS, PLAIN - 300MM</u>					
GW2532DS	0.700/0.810 x 0.001 SG & 70/42 x 0.2BE, 300mm		\$47.23		US
GW2533DS	0.800/0.910 x 0.001 SG & 45/25 x 0.2 BE, 300mm		\$47.23		US
GW2535DS	0.600/1.000 x 0.005 SG & 100/10 x 1.0 BE, 300mm		\$47.23		US
GW2540DS	1.000/1.200 x 0.002 SG & 0/26 x 0.2 BE, 300mm		\$47.23		US
GW2541DS	1.200/1.400 x 0.002 SG & 24/43 x 0.2 BE, 300mm		\$47.23		US
GW2542DS	1.400/1.620 x 0.002 SG & 42/55 x 0.2 BE, 300mm		\$47.23		US
GW2543DS	1.600/1.820 x 0.002 SG & 54/65 x 0.2 BE, 300mm		\$47.23		US
GW2544DS	1.800/2.000 x 0.002 SG & 64/72 x 0.2 BE, 300mm		\$47.23		US
GW2545DS	1.000/2.000 x 0.010 SG & 0/72 x 1.0 BE, 300mm		\$47.23		US
GW2546DS	1.000/1.600 x 0.005 SG & 0/55 x 0.5 BE, 300mm		\$47.23		US
GW2547DS	1.400/2.000 x 0.005 SG & 40/72 x 0.5 BE, 300mm		\$47.23		US

<u>CAT. NO.</u>	<u>DESCRIPTION</u>	<u>HTS Code 9025.80.2000</u>	<u>LIST PRICE</u>	<u>Haz (Hg)</u>	<u>COO</u>
<u>DOUBLE SCALE SPECIFIC GRAVITY/BAUME HYDROMETERS, PLAIN - 150MM</u>					
GW2536DS	0.700/0.800 x 0.002 SG & 70/43 x 0.5 BE, 150mm		\$47.23		US
GW2537DS	0.800/0.900 x 0.002 SG & 45/24 x 0.5 BE, 150mm		\$47.23		US
GW2538DS	0.900/1.000 x 0.002 SG & 25/10 x 0.5 BE , 150mm		\$47.23		US
GW2539DS	0.700/1.000 x 0.010 SG & 70/10 x 1.0 BE, 150mm		\$47.23		US
GW2545ADS	1.000/2.000 x 0.020 SG & 0/72 x 1.0 BE, 150mm		\$47.23		US
GW2548DS	1.000/1.200 x 0.005 SG & 0/26 x 0.5 BE, 150mm		\$47.23		US
GW2549DS	1.200/1.400 x 0.005 SG & 24/43 x 0.5 BE, 150mm		\$47.23		US
GW2549XDS	1.000/1.400 x 0.010 SG & 0/41 x 1.0 BE , 150mm		\$47.23		US
GW2550DS	1.400/1.620 x 0.005 SG & 42/55 x 0.5 BE, 150mm		\$47.23		US
GW2550XDS	1.400/2.000 x 0.010 SG & 40/72 x 1.0 BE, 150mm		\$47.23		US
GW2551DS	1.600/1.820 x 0.005 SG & 54/65 x 0.5 BE , 150mm		\$47.23		US
GW2552DS	1.800/2.000 x 0.005 SG & 64/72 x 0.5 BE, 150mm		\$47.23		US
<u>UNIVERSAL DOUBLE SCALE SPECIFIC GRAVITY/BAUME HYDROMETER, PLAIN</u>					
GW7200	0.700/2.000 x 0.050 S/G & 70/10 x 1.0 LIGHT AND 0/70 x 1.0 HEAVY BAUME		\$36.65		US
<u>BRIX (% SUGAR IN AQUEOUS SOLUTION), CAL 20C, PLAIN - 300mm</u>					
GW2740	-5/+5 x 0.1		\$35.39		US
GW2741	5/+15 x 0.1		\$35.39		US
GW2742	0/6 x 0.1		\$35.39		US
GW2743	6/12 x 0.1		\$35.39		US
GW2516	0/12 x 0.1		\$35.39		US
GW2517	9/21 x 0.1		\$35.39		US
GW2518	19/31 x 0.1		\$35.39		US
GW2519	29/41 x 0.1		\$35.39		US
GW2520	39/51 x 0.1		\$35.39		US
GW2521	49/61 x 0.1		\$35.39		US
GW2522	59/71 x 0.1		\$35.39		US
GW2522X	69/81 x 0.1		\$35.39		US
GW2523	0/35 x 0.5		\$35.39		US
GW2524	35/70 x 0.5		\$35.39		US
<u>BRIX (% SUGAR IN AQUEOUS SOLUTION), CAL 20C, COMBINED, 0/50C - 300MM</u>					
GW2740C	-5/+5 x 0.1		\$68.22		US
GW2741C	5/+15 x 0.1		\$68.22		US
GW2742C	0/6 x 0.1		\$68.22		US
GW2516C	0/12 x 0.1		\$68.22		US
GW2516CXX	0/12 x 0.1, CAL 27.5C		\$68.22		US
GW2517C	9/21 x 0.1		\$68.22		US
GW2518C	19/31 x 0.1		\$68.22		US
GW2519C	29/41 x 0.1		\$68.22		US
GW2520C	39/51 x 0.1		\$68.22		US
GW2522C	59/71 x 0.1		\$68.22		US
GW2522XC	69/81 x 0.1		\$68.22		US
GW2523C	0/35 x 0.5		\$68.22		US
GW2524C	35/70 x 0.5		\$68.22		US
<u>BALLING (% SUGAR IN AQUEOUS SOLUTION), CAL 20C, COMBINED 0/50C - 300MM</u>					
GW2512C	0 to 8.5 x 0.1		\$68.22		US
GW2513C	7.5 to 16 x 0.1		\$68.22		US
GW2514C	16 to 24 x 0.1		\$68.22		US
GW2515CB	0 to 22 x 0.2		\$68.22		US

<u>CAT. NO.</u>	<u>DESCRIPTION</u>	<u>HTS Code</u> 9025.80.2000	<u>LIST PRICE</u>	<u>Haz (Hg)</u>	<u>COO</u>
<u>SALT BRINE HYDROMETERS, SATURATED SALT IN AQUEOUS SOLUTION CAL 60F</u>					
GW2511	0/100% x 1.0, SODIUM CHLORIDE, PLAIN		\$32.83		US
GW2511X	0/26.5% x 0.5 % BY WEIGHT, PLAIN FORM		\$34.12		US
GW2590	0/120% x 1.0, CALCIUM CHLORIDE, PLAIN		\$32.83		US
<u>PROPYLENE GLYCOL HYDROMETER, DUAL SCALE, PLAIN FORM - 300MM</u>					
GW9000	0.990/ 1.050 x 0.002 SPECIFIC GRAVITY & 0/60% x 1.0%		\$36.65		US
<u>LACTOMETERS, COMBINED, +20/100F, Hg, CALIBRATED AT 60F - 330mm</u>					
GW9500	QUEVENNE, 14/42 x 1.0 QUEVENNE & 1.014/1.042 x 0.001 S/G		\$70.76	UN3506	US
GW9520	SPENCE (NY STATE), 0/120 x 1.0 SPENCE		\$70.76	UN3506	US
<u>LACTOMETERS, PLAIN FORM, CALIBRATED AT 60F - 300mm</u>					
GW9525	SPENCE (NY STATE), 0/120 x 1.0 SPENCE		\$42.01		US
GW9540	USDA, 26/37 x 0.2 LACTOMETER SCALE		\$52.45		US
<u>ALCOHOL HYDROMETERS, TRALLE & PROOF SCALES - 300mm</u>					
GW2510	0/100% x 1.0% TRALLE & 0/200% x 2% PROOF CAL. @ 60F, PLAIN FORM		\$32.83		US
GW2510CS	0/100% x 1.0% TRALLE & 0/200% x 2% PROOF CAL. @ 60F, COMBINED FORM -10/100F		\$68.22		US
<u>ALCOHOL HYDROMETERS, IRS SPECIFICATIONS, PROOF SCALES, CAL @ 60F - PLAIN FORM</u>					
GW2601	IRS SZ. F, 0/20% x 0.2%, 305mm length		\$99.18		DE
GW2602	IRS SZ. G, 20/40% x 0.2%, 305mm length		\$99.18		DE
GW2603	IRS SZ. H, 40/60% x 0.2%, 305mm length		\$99.18		DE
GW2604	IRS SZ. I, 60/80% x 0.2%, 305mm length		\$99.18		DE
GW2605	IRS SZ. K, 75/95% x 0.2%, 305mm length		\$99.18		DE
GW2606	IRS SZ. L, 90/110% x 0.2%, 305mm length		\$99.18		DE
GW2607	IRS SZ. M, 105/125% x 0.2%, 305mm length		\$99.18		DE
GW2608	IRS SZ. N, 125/145% x 0.2%, 305mm length		\$99.18		DE
GW2609	IRS SZ. P, 145/165% x 0.2%, 305mm length		\$99.18		DE
GW2610	IRS SZ. Q, 165/185% x 0.2%, 305mm length		\$99.18		DE
GW2611	IRS SZ. R, 185/205% x 0.2%, 305mm length		\$99.18		DE
GW2612	IRS SZ. A, 0/100% x 1.0%, 230mm length		\$99.18		DE
GW2613	IRS SZ. B, 80/120% x 0.5%, 230mm length		\$99.18		DE
GW2614	IRS SZ. C, 100/140% x 0.5%, 230mm length		\$99.18		DE
GW2615	IRS SZ. D, 130/170% x 0.5%, 230mm length		\$99.18		DE
GW2616	IRS SZ. E, 160/206% x 0.5%, 230mm length		\$99.18		DE
<u>THERMCO'S ASTM/API "GOLD BRAND" HYDROMETERS</u>					
<u>ASTM / API HYDROMETERS - PLAIN - 330mm</u>					
ACC001H	ASTM 1H, -1/11 x 0.1, 330mm		\$56.77		DE
ACC002H	ASTM 2H, 9/21 x 0.1, 330mm		\$56.77		DE
ACC003H	ASTM 3H, 19/31 x 0.1, 330mm		\$56.77		DE
ACC004H	ASTM 4H, 29/41 x 0.1, 330mm		\$56.77		DE
ACC005H	ASTM 5H, 39/51 x 0.1, 330mm		\$56.77		DE
ACC006H	ASTM 6H, 49/61 x 0.1, 330mm		\$56.77		DE
ACC007H	ASTM 7H, 59/71 x 0.1, 330mm		\$56.77		DE
ACC008H	ASTM 8H, 69/81 x 0.1, 330mm		\$56.77		DE
ACC009H	ASTM 9H, 79/91 x 0.1, 330mm		\$56.77		DE
ACC010H	ASTM 10H, 89/101 x 0.1, 330mm		\$56.77		DE
<u>ASTM / API HYDROMETERS - PLAIN - 163mm</u>					

<u>CAT. NO.</u>	<u>DESCRIPTION</u>	<u>HTS Code 9025.80.2000</u>	<u>LIST PRICE</u>	<u>Haz (Hg)</u>	<u>COO</u>
ACC021H	ASTM 21H, 0/6 x 0.1, 163mm		\$56.77		DE
ACC022H	ASTM 22H, 5/11 x 0.1, 163mm		\$56.77		DE
ACC023H	ASTM 23H, 10/16 x 0.1, 163mm		\$56.77		DE
ACC024H	ASTM 24H, 15/21 x 0.1, 163mm		\$56.77		DE
ACC025H	ASTM 25H, 20/36 x 0.1, 163mm		\$56.77		DE
ACC026H	ASTM 26H, 25/31 x 0.1, 163mm		\$56.77		DE
ACC027H	ASTM 27H, 30/36 x 0.1, 163mm		\$56.77		DE
ACC028H	ASTM 28H, 35/41 x 0.1, 163mm		\$56.77		DE
ACC029H	ASTM 29H, 40/46 x 0.1, 163mm		\$56.77		DE
ACC030H	ASTM 30H, 45/51 x 0.1, 163mm		\$56.77		DE
ACC031H	ASTM 31H, 50/56 x 0.1, 163mm		\$56.77		DE
ACC033H	ASTM 33H, 60/66 x 0.1, 163mm		\$56.77		DE
ACC034H	ASTM 34H, 65/71 x 0.1, 163mm		\$56.77		DE
ACC035H	ASTM 35H, 70/76 x 0.1, 163mm		\$56.77		DE
ACC039H	ASTM 39H, 90/96 x 0.1, 163mm		\$56.77		DE
ACC321H	ASTM 321H - 775/825 x 0.5, 325mm		\$70.79		DE
<u>ASTM/API COMBINED IN BODY HYDROMETERS - Hg therm. 0/150F - 380mm</u>					
ACC042HL	ASTM 42HL, 22/30 x 0.1 & 0/150F Hg, 380mm		\$99.88	UN3506	DE
ACC051HL	ASTM 51HL, -1/11 x 0.1 & 0/150F Hg, 380mm		\$99.88	UN3506	DE
ACC052HL	ASTM 52HL, 9/21 x 0.1 & 0/150F Hg, 380mm		\$99.88	UN3506	DE
ACC053HL	ASTM 53HL, 19/31 x 0.1 & 0/150F Hg, 380mm		\$99.88	UN3506	DE
ACC054HL	ASTM 54HL, 29/41 x 0.1 & 0/150F Hg, 380mm		\$99.88	UN3506	DE
ACC055HL	ASTM 55HL, 39/51 x 0.1 & 0/150F Hg, 380mm		\$99.88	UN3506	DE
ACC056HL	ASTM 56HL, 49/61 x 0.1 & 0/150F Hg, 380mm		\$99.88	UN3506	DE
ACC057HL	ASTM 57HL, 59/71 x 0.1 & 0/150F Hg, 380mm		\$99.88	UN3506	DE
ACC058HL	ASTM 58HL, 69/81 x 0.1 & 0/150F Hg, 380mm		\$99.88	UN3506	DE
ACC060HL	ASTM 60HL, 89/101 x 0.1 & 0/150F Hg, 380mm		\$99.88	UN3506	DE
<u>ASTM/API COMBINED IN BODY HYDROMETERS, Hg therm. +30/180F - 380mm</u>					
ACC041HM	ASTM 41HM, 15/23 x 0.1 & 30/180F Hg, 380mm		\$99.88	UN3506	DE
ACC042HM	ASTM 42HM, 22/30 x 0.1 & 30/180F Hg, 380mm		\$99.88	UN3506	DE
ACC043HM	ASTM 43HM, 29/37 x 0.1 & 30/180F Hg, 380mm		\$99.88	UN3506	DE
ACC044HM	ASTM 44HM, 36/44 x 0.1 & 30/180F Hg, 380mm		\$99.88	UN3506	DE
ACC045HM	ASTM 45HM, 43/51 x 0.1 & 30/180F Hg, 380mm		\$99.88	UN3506	DE
ACC051HM	ASTM 51HM, -1/11 x 0.1 & 30/180F Hg, 380mm		\$99.88	UN3506	DE
ACC052HM	ASTM 52HM, 9/21 x 0.1 & 30/180F Hg, 380mm		\$99.88	UN3506	DE
ACC053HM	ASTM 53HM, 19/31 x 0.1 & 30/180F Hg, 380mm		\$99.88	UN3506	DE
ACC054HM	ASTM 54HM, 29/41 x 0.1 & 30/180F Hg, 380mm		\$99.88	UN3506	DE
ACC055HM	ASTM 55HM, 39/51 x 0.1 & 30/180F Hg, 380mm		\$99.88	UN3506	DE
ACC056HM	ASTM 56HM, 49/61 x 0.1 & 30/180F Hg, 380mm		\$99.88	UN3506	DE
ACC057HM	ASTM 57HM, 59/71 x 0.1 & 30/180F Hg, 380mm		\$99.88	UN3506	DE
ACC058HM	ASTM 58HM, 69/81 x 0.1 & 30/180F Hg, 380mm		\$99.88	UN3506	DE
ACC059HM	ASTM 59HM, 79/91 x 0.1 & 30/180F Hg, 380mm		\$99.88	UN3506	DE
ACC060HM	ASTM 60HM, 89/101 x 0.1 & 30/180F Hg, 380mm		\$99.88	UN3506	DE
<u>ASTM/API COMBINED IN BODY HYDROMETERS, Hg therm. +60/220F - 380mm</u>					
ACC041HH	ASTM 41HH, 15/23 x 0.1 & 60/220F Hg, 380mm		\$99.88	UN3506	DE
ACC042HH	ASTM 42HH, 22/30 x 0.1 & 60/220F Hg, 380mm		\$99.88	UN3506	DE
ACC043HH	ASTM 43HH, 29/37 x 0.1 & 60/220F Hg, 380mm		\$99.88	UN3506	DE
ACC044HH	ASTM 44HH, 36/44 x 0.1 & 60/220F Hg, 380mm		\$99.88	UN3506	DE
ACC045HH	ASTM 45HH, 43/51 x 0.1 & 60/220F Hg, 380mm		\$99.88	UN3506	DE
ACC051HH	ASTM 51HH, -1/11 x 0.1 & 60/220F Hg, 380mm		\$99.88	UN3506	DE
ACC052HH	ASTM 52HH, 9/21 x 0.1 & 60/220F Hg, 380mm		\$99.88	UN3506	DE
ACC053HH	ASTM 53HH, 19/31 x 0.1 & 60/220F Hg, 380mm		\$99.88	UN3506	DE
ACC054HH	ASTM 54HH, 29/41 x 0.1 & 60/220F Hg, 380mm		\$99.88	UN3506	DE
ACC055HH	ASTM 55HH, 39/51 x 0.1 & 60/220F Hg, 380mm		\$99.88	UN3506	DE
ACC056HH	ASTM 56HH, 49/61 x 0.1 & 60/220F Hg, 380mm		\$99.88	UN3506	DE

<u>CAT. NO.</u>	<u>DESCRIPTION</u>	<u>HTS Code 9025.80.2000</u>	<u>LIST PRICE</u>	<u>Haz (Hg)</u>	<u>COO</u>
ACC057HH	ASTM 57HH, 59/71 x 0.1 & 60/220F Hg, 380mm		\$99.88	UN3506	DE
ACC058HH	ASTM 58HH, 69/81 x 0.1 & 60/220F Hg, 380mm		\$99.88	UN3506	DE
<u>ASTM/API COMBINED IN STEM HYDROMETER, Hg therm. +30/220F - 380mm</u>					
ACC072H	ASTM 72H, 9/21 x 0.1 & 30/220F Hg, 380mm		\$113.40	UN3506	DE
ACC073H	ASTM 73H, 19/31 x 0.1 & 30/220F Hg, 380mm		\$113.40	UN3506	DE
ACC074H	ASTM 74H, 29/41 x 0.1 & 30/220F Hg, 380mm		\$113.40	UN3506	DE
<u>API, COMBINED HYDROMETERS, 0/150F, Spirit - 15"</u>					
ACC6001S	45/90 x 1.0, 0/150F Spirit, 15"		\$116.37		DE
ACC6001XXS	30/90 x 0.5, 0/150F Spirit, 15"		\$116.37		DE
ACC6003S	20/40 x 0.5, 0/150F Spirit, 15"		\$116.37		DE
ACC6004S	30/50 x 0.5, 0/150F Spirit, 15"		\$116.37		DE
ACC6005S	40/60 x 0.5, 0/150F Spirit, 15"		\$116.37		DE
ACC6006S	50/70 x 0.5, 0/150F Spirit, 15"		\$116.37		DE
ACC6007S	60/80 x 0.5, 0/150F Spirit, 15"		\$116.37		DE
<u>RASCHER & BETZOLD A.S.B.C. Top Reading Pyrex® Brand GLASS HYDROMETERS</u>					
*standardization temperature 20°C / Hg / 0 to 30°C					
R164978	BEER/WORT, -2/6% X 0.1, 270MM length, 110ml of sample needed *		\$277.82	UN3506	DE
R164980	BEER/WORT, 0/8% X 0.1, 270MM length, 110ml of sample needed		\$277.82	UN3506	DE
R165000	BEER/WORT, 8/16% X 0.1, 270MM length, 110ml of sample needed		\$277.82	UN3506	DE
R165020	BEER/WORT, 16/24% X 0.1, 270MM length, 110ml of sample needed		\$277.82	UN3506	DE
R165030	BEER/WORT, 0/8% X 0.1, 330MM length, 240ml of sample needed		\$277.82	UN3506	DE
R165032	BEER/WORT, 8/16% X 0.1, 330MM length, 240ml of sample needed		\$277.82	UN3506	DE
R165034	BEER/WORT, 16/24% X 0.1, 330MM length, 240ml of sample needed		\$277.82	UN3506	DE
R165040	BEER/WORT, 0/14% X 0.1, 406MM length, 430ml of sample needed		\$277.82	UN3506	DE
R165060	BEER/WORT, 12/26% X 0.1, 406MM length, 430ml of sample needed		\$277.82	UN3506	DE
R165160	BEER/WORT, 2/8% X 0.1, 406MM length, 430ml of sample needed		\$277.82	UN3506	DE
R165161	BEER/WORT, 8/14% X 0.1, 406MM length, 430ml of sample needed		\$277.82	UN3506	DE
R165162	BEER/WORT, 10/16% X 0.1, 406MM length, 430ml of sample needed		\$277.82	UN3506	DE
R165168	BEER/WORT, 0/10% X 0.1, 432MM length, 430ml of sample needed		\$277.82	UN3506	DE
R165170	BEER/WORT, 10/20% X 0.1, 432MM length, 430ml of sample needed		\$277.82	UN3506	DE
R165078	BEER/WORT, -1/2% X 0.05, 356MM length, 370ml of sample needed		\$277.82	UN3506	DE
R165082	BEER/WORT, 2/5% X 0.05, 356MM length, 370ml of sample needed		\$277.82	UN3506	DE
R165089	BEER/WORT, 6/9% X 0.05, 356MM length, 370ml of sample needed		\$277.82	UN3506	DE
R165092	BEER/WORT, 9/12% X 0.05, 356MM length, 370ml of sample needed		\$277.82	UN3506	DE
R165096	BEER/WORT, 11/14% X 0.05, 356MM length, 370ml of sample needed		\$277.82	UN3506	DE
R165097	BEER/WORT, 12/15% X 0.05, 356MM length, 370ml of sample needed		\$277.82	UN3506	DE
R165100	BEER/WORT, 14/17% X 0.05, 356MM length, 370ml of sample needed		\$277.82	UN3506	DE

<u>CAT. NO.</u>	<u>DESCRIPTION</u>	<u>HTS Code</u>	<u>9025.19.8085</u>	<u>LIST PRICE</u>	<u>%</u>	<u>OUR PRICE</u>	<u>COO</u>
<u>THERMCO BRAND DIGITAL DATA LOGGERS</u>							
<u>USB TYPE</u>							
ACCDPICOLITE	Low Cost Micro-Therm USB Water Resistant Logger w/Software, -40 to +40C, 16K samples			\$40.50	25%	\$30.41	IL
ACCD00408	Micro-Therm USB Waterproof Logger w/Software, -40 to +80C, 32K samples			\$133.72	25%	\$100.28	IL
ACCD00408RH	Micro-Therm USB Waterproof Logger for Temperature & Humidity w/Software, -40 to 80C & 5 to 95% Rh, 32K samples			\$201.22	25%	\$150.90	IL
NISTUSB	NIST traceable ISO 17025 Calibration, 1 Temperature (SERVICE)			\$77.17	27%	\$56.19	US/ Service
<u>VFC PROBE TYPE</u>							
ACCRT8017	AccuTherm Dual buffered probe DDL, -50/70C, VFC compliant, 2 yr. NIST ISO Calibration -20C & 4C			\$186.85	39%	\$114.53	CN
ACCRTAMB	Ambient Room Temperature Sensor Probe for ACCRT8017 (accessory)			\$9.65	50%	\$4.84	CN
ACCRT8002	AccuTherm Single buffered probe DDL, -50/70C, VFC compliant, 2 yr. NIST ISO Calibration 4C			\$168.76	36%	\$108.48	CN
PBRT8002-STprobe	Straight probe for ACCRT8002 & ACCRT8017 (accessory)			\$24.13	50%	\$12.07	CN
ACCSL2021	AccuTherm SmartLOG 2021 DDL, -40/70C, w/ 1 wireless temperature sensor, email & SMS alerts, 2 year NIST Calibration 4C			\$211.14	25%	\$158.36	CN
ACCSLBLET	AccuTherm Temperature Sensor for SmartLOG DDL, -40/70C, w/ wireless temperature sensor only, 2 year NIST Calibration 4C			\$104.53	24%	\$79.19	CN
<u>ULTRA-LOW TYPE</u>							
ACCRT8019UL	AccuTherm Dual buffered probe DDL, ultra low -80 to +70C, 2 yr NIST Calibration -70C, VFC compliant			\$186.85	39%	\$114.53	CN
ACCCRY086	Ultra-Low Temperature Data Logger, -86 to +35C x 0.1C; does NOT include NIST Calibration Certificate			\$336.22	33%	\$225.32	US
ACCTIFC300	Windows Software/jump drive and cradle for ACCCRY086 (Required accessory)			\$174.21	33%	\$116.66	US
NISTCRYO	Certificate Traceable to N.I.S.T., 1 Temperature of ACCCRY086 (service)			\$165.92	35%	\$107.87	US/ Service
<u>HIGH TEMP/STERILIZATION TYPE</u>							
ACCHITEMP140	Autoclave Data Logger, -20 to 140C x 0.05C, Accuracy +/- 0.5C (formerly ACCHITEMP135)			\$606.23	33%	\$406.23	US
ACCHITEMP140FC	32,767 Readings; Does NOT include NIST Calibration Certificate Same as ACHITEMP140, except includes NIST Calibration @ 120.0C or 121.0C (does not include software for operation ACCTIFC400)			\$673.76	33%	\$451.44	US
ACCHITEMP140PT	High Temp Data Logger w/24" Flex Probe; Does NOT include NIST Calibration Certificate			\$606.23	33%	\$406.23	US
ACCHITEMP140-5	Autoclave Data Logger, -40 to 140C x 0.05C, w/5.25" probe, Accuracy +/- 0.5C; Does NOT include NIST Calibration Certificate			\$654.09	37%	\$412.91	US
ACCTIFC400	Windows Software for ACCHITEMP135/ACCHITEMP140 (Required accessory)			\$174.21	33%	\$116.66	US
ACCHITEMPBAT-2	Battery for ACCHITEMP135/ACCHITEMP140 (accessory)			\$94.51	33%	\$63.31	CN
NISTHITEMP	Certificate Traceable to N.I.S.T., 1 Temperature (service)			\$128.58	35%	\$83.64	US/ Service
<u>TEMPERATURE & RELATIVE HUMIDITY TYPE</u>							
HICSSN-22	Temperature & Humidity Data Logger, Digital Readout, 32K samples -30 to +70C (-35 to +176F) & 0 to 99% Rh, W/Audio & Visual Alarm			\$200.56	49%	\$102.25	CN

<u>CAT. NO.</u>	<u>DESCRIPTION</u>	<u>HTS Code</u>	<u>9025.19.8085</u>	<u>LIST PRICE</u>	<u>%</u>	<u>OUR PRICE</u>	<u>COO</u>
ACCDTA140	Thermo-Hygrometer Data Logger, w/Dew Point, Min-Max, Average, Values, Alarm, RS232 Interface, Wireless Connection of up to 8 transmitters. Does NOT include transmitter(s)			\$253.26	33%	\$169.68	DE
ACCDTA141	Transmitter temperature & Rh for the ACCDTA140			\$68.73	33%	\$46.06	DE
ACCDTA150RFV	Thermo-Hygrometer, Data Logger, Wireless w/vial & probe transmitter, NIST traceable certificate @ 4C			\$329.60	33%	\$220.85	DE
ACCDTA151	Transmitter w/ vial & probe for ACCDTA150RFV only			\$81.46	33%	\$54.57	DE

BERLINGER BRAND DIGITAL DATA LOGGERS

BERFREEZTAG2L	Vaccine Data Logger w/ Vial & Probe, Alarm Setting <-15C @ 60Min, 1.5m probe cable, fixed alarm setting Memory & Longer probe cable; VFC Compliant, 2 yr NIST Calibration @ -20C			\$175.10	21%	\$138.50	PH
BERFRIDGETAG2L	Vaccine Data Logger w/ Vial & Probe, Alarm Setting 2.0C to 8.0C, 1.5m probe cable, user configurable Memory & Longer probe cable; VFC Compliant, 2 yr NIST Calibration @ 4C			\$175.10	21%	\$138.50	PH
BERFT2ULF	Fridgetag ultra low DDL, -90/0C & -130/32F, non-buffered probe only, NIST Calibration @ -70C			\$261.17	27%	\$191.44	PH
BERVAXBKUP2L	Vaccine Data Logger w/ Vial & Probe, User Configurable Alarms, 1.5m probe cable, VFC compliant, 2 yr NIST Calibration -20 & 4C						PH

VFC@400/800 BRAND TEMPERATURE DATA LOGGERS

VFC PROBE TYPE

VFC400-VMKCS	Vaccine Monitoring Kit, includes VFC400 datalogger, probe, buffer, acrylic stand, wall mount bracket, ISO 17025 Certificate. **Excludes LTIHID, docking station			\$155.21	15%	\$131.96	HK
VFC400-SMKCS	Site Vaccine Monitoring Kit, includes (2) VFC400 dataloggers, probes w/buffer, acrylic stands, wall mount brackets, ISO 17025 Certificate. Includes 1 LTIHID, docking station			\$315.68	8%	\$290.11	HK
VFC400-SP	VFC400 w/ Smart Probe technology, hot swap calibrated probe. VFC compliant ISO 17025 certificate, requires LTIHID for operation			\$196.00	24%	\$148.05	HK
VFC800WIFI-SP	VFC800wifi w/ Smart Probe technology, hot swap calibrated probe. VFC compliant ISO 17025 certificate, requires LTO subscription			\$289.00	22%	\$224.00	HK
VFCSP-5	Smart probe, 5ft., calibrated w/ 60ml glycol, for VFC400®			\$69.00	19%	\$55.99	HK
VFCSP-10	Smart probe, 10ft., calibrated w/ 60ml glycol, for VFC800®			\$79.00	24%	\$59.99	HK

LOGTAG BRAND TEMPERATURE DATA LOGGERS

VFC PROBE TYPE

LTTRED3016R	VFC Approved Temperature Data Logger, 30 Day Summary, W/Probe & Vial & wall mount bracket LTBRACKET (does not include Docking Station LTIHID), 2 yr NIST Calibration @ 4C			\$117.27	25%	\$87.99	VN
LTTRED3016RUSB	VFC Approved Temperature Data Logger, 30 Day Summary, W/Probe & Vial, Docking Station LTIHID& wall mount bracket LTBRACKET; NIST Calibration @ 4C			\$180.04	26%	\$134.05	VN
LTUTRED16	LogTag DDL, -40/99C w/ buffed probe, USB 2.0, A-type plug, no interface cradle required, multi use, NIST Calibration @ 4C			\$131.98	25%	\$98.46	HK
LTUTRED30WIFI	VFC Approved Temperature Data Logger, 30 Day Summary, W/3M Probe, -40/90C & Vial, Real time WiFi enabled w/ LogTag Online & wall mount bracket LTBRACKETWIFI, NIST Calibration @ 4C			\$327.48	33%	\$218.32	HK
LTUTRED3016	VFC Datalogger, -40/99C, extra large display for easy viewing, 30 day summary, probe & vial, wall bracket, Micro USB 2.0			\$200.53	33%	\$133.68	HK

<u>CAT. NO.</u>	<u>DESCRIPTION</u>	<u>HTS Code</u>	<u>9025.19.8085</u>	<u>LIST PRICE</u>	<u>%</u>	<u>OUR PRICE</u>	<u>COO</u>
	NIST Calibration 4C						
LTUTRED3016-2	VFC Datalogger, -40/99C, extra large display for easy viewing, 30 day summary, dual Smart Probe sensors & vials, wall bracket, Micro USB 2.0, NIST Calibration -20C & 4C			\$314.15	29%	\$221.76	HK
UTRED30WIFI-SP	LogTag® Thermco WIFI DDL w/ Smart Probe technology, hot swap calibrated probe. Single probe. VFC compliant ISO 17025 certificate, requires LTO subscription			\$289.00	31%	\$197.99	HK
UTRED30WIFI-SP2	LogTag® Thermco WIFI DDL w/ Smart Probe technology, hot swap calibrated probe. Dual probe. VFC compliant ISO 17025 certificate, requires LTO subscription			\$343.00	25%	\$257.50	HK
<u>TEMPERATURE AND RELATIVE HUMIDITY TYPE</u>							
LTHAX08	Data Logger Temperature and Humidity; -40/85 C 0/100% RH. *Excludes Docking Station LTIHID			\$71.07	17%	\$59.13	HK
LTUHADO-16	Data Logger Temperature and Humidity w/ display; -30°C to +70°C and 0/100%RH. USB			\$101.97	13%	\$88.74	HK
<u>ULTRA LOW TYPE</u>							
LTUTREL3016	Ultra low Datalogger, -90/40C, extra large display for easy viewing, 30 day summary, w/ probe & buffer, wall bracket, Micro USB 2.0 NIST Calibration @ -70C			\$200.53	33%	\$133.68	HK
LTUTREL30WIFI	30 day WiFi Ultra low data logger, -90/40C & -130/210F, air probe w/ wall mount bracket. Real-time notification w/LogTag Online NIST Calibration @ -70C			\$272.10	28%	\$195.61	HK
LTUTREL16	LogTag ultra low DDL, -90/40C w/ aire probe, USB 2.0, A-type plug, no interface cradle required, multi use, NIST Calibration @ -70C			\$101.64	20%	\$80.87	HK
LTTREL3016	LogTag ultra low DDL, -90°C to +40°Cx 0.1C; w/ probe & wall bracket; does not include Docking Station LTIHID, NIST Calibration @ -70C			\$120.20	25%	\$90.19	HK
LTTREL3016USB	LogTag ultra low DDL, -90°C to +40°C x 0.1C; w/probe & wall bracket, includes Docking Station LTIHID, NIST Calibration @ -70C			\$182.38	25%	\$136.38	HK
LTTREL3016CVD	LogTag ultra low DDL, -90°C to +40°C (-130°F to +104°F) expanded range w/ probe and wall bracket, includes LTIHID interface cradle, NIST Calibration @ -70C			\$159.54	25%	\$120.21	HK
<u>LOW TYPE</u>							
<i>*excludes calibration</i>							
LTTREL8	Multi-Use Low Temperature Data Logger w/ external probe, -80 to +40C			\$55.62	21%	\$43.82	HK
LTTREL8	Dry Ice MULTI USE Low Temp Data Logger w/ internal probe, -80 to			\$83.80	13%	\$73.24	HK
LTSRIL8	Dry Ice SINGLE USE Low Temp Data Logger w/ internal probe, -80 to			\$66.95	21%	\$52.90	HK
LTTREL8	MULTI USE Data Logger, w/internal Sensor, -40 to 80C; does not inclu			\$32.96	28%	\$23.80	HK
LTUSRID16	Single Use USB PDF digital datalogger, -30/60C			\$29.87	25%	\$22.29	HK
<u>LOGTAG ACCESSORIES & SERVICES</u>							
LTIHID	Interface Docking Station for LogTag®Dataloggers (replaces LTIUSB) (Required accessory)			\$50.47	15%	\$43.00	HK
LTBRACKET	Wall Mount Bracket for LTTRED3016R & LTHAX08			\$6.18	17%	\$5.11	HK
LTBRACKETWIFI	Wall Mount Bracket for UTRED3016 & UTRED30WIFI			\$7.21	29%	\$5.11	HK
LTI-WM-WIFI	Wall mount Wi-Fi Interface Cradle LogTag®Dataloggers (retrofits any LogTag® datalogger w/ probe)			\$113.30	13%	\$98.85	HK
LTISOLOGTEMP	ISO17025/IEC Accredited Temperature Calibration for LogTag DDL			\$88.65	25%	\$66.49	US/ Service
LTISOLOGRH	ISO17025/IEC Accredited Temperature/Humidity Calibration for LogTag DDL			\$101.70	25%	\$76.25	US/ Service
LTCP100K-30	External Smart Probe sensor for dual channel UTRED30, 3.0 meter w/ buffer *for use w/ UTRED30 only (mfg after 2023)			\$60.77	13%	\$53.15	HK

<u>CAT. NO.</u>	<u>DESCRIPTION</u>	<u>HTS Code</u>	<u>9025.19.8085</u>	<u>LIST PRICE</u>	<u>%</u>	<u>OUR PRICE</u>	<u>COO</u>
<u>LASCAR BRAND DIGITAL DATA LOGGERS</u>							
EL-USB-2-LCD	Temperature, Humidity and Dew Point DDL w/ LCD screen. -35 to +80°C *excludes calibration			\$126.99	33%	\$84.66	GB
EL-WIFI-VACX	-40/60C, external buffered sensor, w/ calibration			\$353.80	35%	\$230.72	GB
EL-WIFI-TH	-20/60C & 0/100%, DDL *excludes calibration			\$206.63	19%	\$167.71	GB
EL-WIFI-21CFR-TH	-20/60C & 0/100%, DDL, 21 CFR COMPLIANT, *excludes calibration			\$257.59	29%	\$182.61	GB
EL-IOT-CO2	-20/60C, 0/100%RH, 0/40k PPM, AA Battery, air quality & CO2 DDL			\$219.19	20%	\$175.55	GB

<u>CAT. NO.</u>	<u>DESCRIPTION</u>	<u>LIST PRICE</u>	<u>COO</u>
	<u>HTS Code 84798999</u>		
	<u>THERMCO DRY ICE MACHINE</u>		
THDI800	Benchtop dy ice maker; 1lb. Blocks; simple lever operation	\$1,261.62	US
THDI800-010	Accessory; replacment membrane for THDI800	\$41.00	US

<u>CAT. NO.</u>	<u>DESCRIPTION</u>		<u>LIST PRICE</u>	<u>Haz (Hg)</u>	<u>COO</u>
<u>ASTM THERMOMETERS, MERCURY*</u>					
ACC001C	ASTM 1C, -20/150C x 1.0C, Hg, 76mm Imm	PARTIAL IMMERSION	\$42.27	UN3506	DE, FR
ACC001F	ASTM 1F, 0/302F x 2.0F, Hg, 76mm Imm	PARTIAL IMMERSION	\$42.27	UN3506	DE, FR
ACC002C	ASTM 2C, -5/300C x 1.0C, Hg, 76mm Imm	PARTIAL IMMERSION	\$54.45	UN3506	DE
ACC002F	ASTM 2F, 20/580F X 2.0F, Hg, 76mm Imm	PARTIAL IMMERSION	\$54.45	UN3506	DE, FR
ACC003C	ASTM 3C, -5/400C x 1.0C, Hg, 76mm Imm	PARTIAL IMMERSION	\$65.24	UN3506	DE
ACC003F	ASTM 3F, 20/760F x 2.0F, Hg, 76mm Imm	PARTIAL IMMERSION	\$65.24	UN3506	DE, FR
ACC004C	ASTM 4C, -1/105C x 0.5C, Hg, 152mm Imm	ACID HEAT	\$72.81	UN3506	DE
ACC004F	ASTM 4F, 30/220F x 1.0F, Hg, 152mm Imm	ACID HEAT	\$72.81	UN3506	DE
ACC005C	ASTM 5C, -38/50C x 1.0C, Hg, 108mm Imm	CLOUD & POUR	\$43.66	UN3506	DE, FR
ACC005F	ASTM 5F, -36/120F x 2.0F, Hg, 108mm Imm	CLOUD & POUR	\$43.66	UN3506	DE, FR
ACC006C	ASTM 6C, -80/20C x 1.0C, Non-Hg, 76mm Imm	LOW CLOUD & POUR*	\$48.71		DE
ACC006F	ASTM 6F, -112/70F x 2.0F, Non-Hg, 76mm Imm	LOW CLOUD & POUR*	\$48.71		DE, FR
ACC007C	ASTM 7C, -2/300C x 1.0C, Hg, Tol Imm	LOW DISTILLATION	\$60.99	UN3506	TW, DE
ACC007F	ASTM 7F, 30/580F x 2.0F, Hg, Tol Imm	LOW DISTILLATION	\$60.99	UN3506	DE, FR
ACC008C	ASTM 8C, -2/400C x 1.0C, Hg, Tol Imm	HIGH DISTILLATION	\$60.99	UN3506	TW, DE
ACC008F	ASTM 8F, 30/760F x 2.0F, Hg, Tol Imm	HIGH DISTILLATION	\$60.99	UN3506	DE, FR
ACC009C	ASTM 9C, -5/110C x 0.5C, Hg, 57mm Imm	PENSKY-MARTENS LOW	\$49.87	UN3506	TW, DE
ACC009F	ASTM 9F, 20/230F x 1.0F, Hg, 57mm Imm	PENSKY-MARTENS LOW	\$49.87	UN3506	TW, DE
ACC010C	ASTM 10C, 90/370C x 2.0C, Hg, 57mm Imm	PENSKY-MARTENS HIGH	\$60.99	UN3506	DE
ACC010F	ASTM 10F, 200/700F x 5.0F, Hg, 57mm Imm	PENSKY-MARTENS HIGH	\$60.99	UN3506	DE
ACC011C	ASTM 11C, -6/400C x 2.0C, Hg, 25mm Imm	OPEN FLASH	\$60.99	UN3506	TW, DE
ACC011F	ASTM 11F, 20/760F x 5.0F, Hg, 25mm Imm	OPEN FLASH	\$60.99	UN3506	TW, DE
ACC012C	ASTM 12C, -20/102C x 0.2C, Hg, Tol Imm	GRAVITY	\$69.62	UN3506	TW, DE
ACC012F	ASTM 12F, -5/215F x 0.5F, Hg, Tol Imm	GRAVITY	\$69.62	UN3506	TW, DE
ACC013C	ASTM 13C, 155/170C x 0.5C, Hg, Tol Imm	LOSS on HEAT	\$69.62	UN3506	DE, FR
ACC014C	ASTM 14C, 38/82C x 0.1C, Hg, 79mm Imm	PARRIFIN MELTING POINT	\$67.95	UN3506	DE
ACC014F	ASTM 14F, 100/180F x 0.2F Hg, 79mm Imm	PARRIFIN MELTING POINT	\$67.95	UN3506	DE, FR
ACC015C	ASTM 15C, -2/80C x 0.2C, Hg, Tol Imm	LOW SOFTENING POINT	\$64.43	UN3506	DE, FR
ACC015F	ASTM 15F, 30/180F x 0.5F, Hg, Tol Imm	LOW SOFTENING POINT	\$64.43	UN3506	DE, FR
ACC016C	ASTM 16C, 30/200C x 0.5C, Hg, Tol Imm	HIGH SOFTENING POINT	\$67.95	UN3506	DE, FR
ACC016F	ASTM 16F, 85/392F x 1.0F, Hg, Tol Imm	HIGH SOFTENING POINT	\$67.95	UN3506	DE
ACC017C	ASTM 17C, 19/27C x 0.1C, Hg, Tol Imm	SAYBOLT VISCOSITY	\$56.72	UN3506	DE, FR

<u>CAT. NO.</u>	<u>DESCRIPTION</u>		<u>LIST PRICE</u>	<u>Haz (Hg)</u>	<u>COO</u>
ACC017F	ASTM 17F, 66/80F x 0.2F, Hg, Tol Imm	SAYBOLT VISCOSITY	\$56.72	UN3506	DE
ACC018C	ASTM 18C, 34/42C x 0.1C, Hg, Tol Imm	SAYBOLT VISCOSITY	\$56.72	UN3506	DE, FR
ACC018F	ASTM 18F, 94/108F x 0.2F, Hg, Tol Imm	SAYBOLT VISCOSITY	\$56.72	UN3506	DE
ACC019C	ASTM 19C, 49/57C x 0.1C, Hg, Tol Imm	SAYBOLT VISCOSITY	\$56.72	UN3506	DE, FR
ACC019F	ASTM 19F, 120/134F x 0.2F, Hg, Tol Imm	SAYBOLT VISCOSITY	\$56.72	UN3506	DE, FR
ACC020C	ASTM 20C, 57/65C x 0.1C, Hg, Tol Imm	SAYBOLT VISCOSITY	\$56.72	UN3506	DE, FR
ACC020F	ASTM 20F, 134/148F x 0.2F, Hg, Tol Imm	SAYBOLT VISCOSITY	\$56.72	UN3506	DE
ACC021C	ASTM 21C, 79/87C x 0.1C, Hg, Tol Imm	SAYBOLT VISCOSITY	\$56.72	UN3506	DE, FR
ACC021F	ASTM 21F, 174F/188F x 0.2F, Hg, Tol Imm	SAYBOLT VISCOSITY	\$56.72	UN3506	DE
ACC022C	ASTM 22C, 95/103C x 0.1C, Hg, Tol Imm	SAYBOLT VISCOSITY	\$56.72	UN3506	DE, FR
ACC022F	ASTM 22F, 204/218F X 0.2F, Hg, Tol Imm	SAYBOLT VISCOSITY	\$56.72	UN3506	DE
ACC023C	ASTM 23C, 18/28C x 0.2C, Hg, 90mm Imm	ENGLER VISCOSITY	\$76.51	UN3506	DE, FR
ACC024C	ASTM 24C, 39/54C x 0.2C, Hg, 90mm Imm	ENGLER VISCOSITY	\$76.51	UN3506	DE, FR
ACC025C	ASTM 25C, 95/105C x 0.2C, Hg, 90mm Imm	ENGLER VISCOSITY	\$76.51	UN3506	DE, FR
ACC026C	ASTM 26C, 130/140C x 0.1C, Hg, Tol immersion	STABILITY TEST of NITROCELLULOSE	\$100.06	UN3506	DE, FR
ACC027C	ASTM 27C, 147/182C x 0.5C, Hg, 76mm Imm	TURPENTINE DISTILLATION	\$87.49	UN3506	DE
ACC028C	ASTM 28C, 36.6/39.4C x 0.05C, Hg, Tol Imm	KINEMATIC VISCOSITY	\$181.47	UN3506	DE
ACC028F	ASTM 28F, 97.5/102.5F x 0.1F, Hg, Tol Imm	KINEMATIC VISCOSITY	\$181.47	UN3506	DE, FR
ACC029C	ASTM 29C, 52.6/55.4C x 0.05C, Hg, Tol Imm	KINEMATIC VISCOSITY	\$181.47	UN3506	DE
ACC029F	ASTM 29F, 127.5/132.5F x 0.1F, Hg, Tol Imm	DISCONTINUED			
ACC030F	ASTM 30F, 207.5/212.5F x 0.1F, Hg, Tol Imm	KINEMATIC VISCOSITY	\$181.47	UN3506	DE
ACC031F	ASTM 31F, -30/120F x 1.0F, Hg, Tol Imm	REID VAPOR PRESSURE AIR CHAMBER	\$65.09	UN3506	DE
ACC033C	ASTM 33C, -38/42C x 0.2C, Hg, 50mm Imm	ANILINE POINT	\$98.39	UN3506	DE, FR
ACC033F	ASTM 33F, -36.5/107.5F x 0.2F, Hg, 50mm Imm	ANILINE POINT	\$98.39	UN3506	DE, FR
ACC034C	ASTM 34C, 25/105C x 0.2C, Hg, 50mm Imm	ANILINE POINT	\$80.03	UN3506	DE
ACC034F	ASTM 34F, 77/221F x 0.5F, Hg, 50mm Imm	ANILINE POINT	\$80.03	UN3506	DE
ACC035C	ASTM 35C, 90/170C x 0.2C 50mm Hg	ANILINE POINT	\$98.39	UN3506	DE, FR
ACC035F	ASTM 35F, 194/338f x 0.5F, Hg, 50mm Imm	ANILINE POINT	\$98.39	UN3506	DE, FR
ACC036C	ASTM36C, -2/68C x 0.2C, Hg, 45mm Imm	TITER TEST	\$80.20	UN3506	DE
ACC037C	ASTM 37C, -2/52C x 0.2C, Hg, 100mm Imm	SOLVENTS DISTILLATION	\$75.18	UN3506	DE, FR
ACC038C	ASTM 38C, 24/78C x 0.2C, Hg, 100mm Imm	SOLVENTS DISTILLATION	\$75.18	UN3506	DE, FR
ACC039C	ASTM 39C, 48/102C x 0.2, Hg, 100mm Imm	SOLVENTS DISTILLATION	\$81.26	UN3506	DE, FR
ACC040C	ASTM 40C, 72/126C x 0.2C, Hg, 101mm Imm	SOLVENTS DISTILLATION	\$81.26	UN3506	DE, FR
ACC041C	ASTM 41C, 98/152C x 0.2C, Hg, 100mm Imm	SOLVENTS DISTILLATION	\$81.26	UN3506	DE
ACC042C	ASTM 42C, 95/255C x 0.5C, Hg, 100mm Imm	SOLVENTS DISTILLATION	\$94.67	UN3506	DE
ACC044C	ASTM 44C, 18.6/21.4C x 0.05C, Hg, Tol Imm	KINEMATIC VISCOSITY	\$175.13	UN3506	DE
ACC044F	ASTM 44F, 66.5/71.5F x 0.1F, Hg, Tol Imm	KINEMATIC VISCOSITY	\$175.13	UN3506	DE
ACC045C	ASTM 45C, 23.6/26.4C x 0.05C, Hg, Tol Imm	KINEMATIC VISCOSITY	\$175.13	UN3506	DE

<u>CAT. NO.</u>	<u>DESCRIPTION</u>		<u>LIST PRICE</u>	<u>Haz (Hg)</u>	<u>COO</u>
ACC045F	ASTM 45F, 74.5/79.5F x 0.1F, Hg, Tol Imm	KINEMATIC VISCOSITY	\$175.13	UN3506	DE
ACC046C	ASTM 46C, 48.6/51.4C x 0.05C, Hg, Tol Imm	KINEMATIC VISCOSITY	\$175.13	UN3506	DE
ACC046F	ASTM 46F, 119.5/124.5F x 0.1F, Hg, Tol Imm	KINEMATIC VISCOSITY	\$175.13	UN3506	DE, FR
ACC047C	ASTM 47C, 58.6/61.4C x 0.05C, Hg, Tol Imm	KINEMATIC VISCOSITY	\$175.13	UN3506	DE
ACC047F	ASTM 47F, 137.5/142.5F x 0.1F, Hg, Tol Imm	KINEMATIC VISCOSITY	\$175.13	UN3506	DE, FR
ACC048C	ASTM 48C, 80.6/83.4C x 0.05C, Hg, Tol Imm	KINEMATIC VISCOSITY	\$175.13	UN3506	DE
ACC048F	ASTM 48F, 177.5/182.5F x 0.1F, Hg, Tol Imm	KINEMATIC VISCOSITY	\$175.13	UN3506	DE
ACC049C	ASTM 49C, 20/70C x 0.2C, Hg, 65mm Imm	STORMER VISCOSITY	\$71.66	UN3506	DE, FR
ACC050F	ASTM 50F, 54/101F x 0.1F, Hg, Tol Imm	GAS CALORIMETER, INLET	\$150.26	UN3506	DE
ACC051F	ASTM 51F, 69/116F x 0.1F, Hg, Tol Imm	GAS CALORIMETER, OUTLET	\$150.26	UN3506	DE
ACC053C	ASTM 53C, -0.6/10.4C x 0.1C, Hg, Tol Imm	BENZENE FREEZING POINT	\$84.45	UN3506	DE
ACC054C	ASTM 54C, 20/100.6C x 0.2C, Hg, Tol Imm	CONGEALING POINT	\$92.85	UN3506	DE
ACC054F	ASTM 54F, 68/213F x 0.5F, Hg, Tol Imm	CONGEALING POINT	\$92.85	UN3506	DE
ACC056C	ASTM 56C, 19/35C x 0.02C, Hg, Tol Imm	BOMB CALORIMETER	\$322.00	UN3506	DE
ACC056F	ASTM 56F, 66/95F x 0.05F, Hg, Tol Imm	BOMB CALORIMETER	\$269.58	UN3506	DE
ACC057C	ASTM 57C, -20/50C x 0.5F, Hg, 57mm Imm	TAG CLOSED TESTER, LOW RANGE	\$67.80	UN3506	DE, FR
ACC057F	ASTM 57F, -4/122F x 1.0F, Hg, 57mm Imm	TAG CLOSED TESTER, LOW RANGE	\$67.80	UN3506	DE
ACC058C	ASTM 58C, -34/49C x 0.5C, Hg, Tol Imm	TANK	\$49.33	UN3506	DE
ACC058F	ASTM 58F, -30/120F x 1.0F, Hg, Tol Imm	TANK	\$49.33	UN3506	DE, TW
ACC059C	ASTM 59C, -18/82C x 0.5C, Hg, Tol Imm	TANK	\$49.33	UN3506	DE
ACC059F	ASTM 59F, 0/180F x 1.0F, Hg, Tol Imm	TANK	\$49.33	UN3506	DE, TW
ACC060C	ASTM 60C, 77/260C x 1.0C, Hg, Tol Imm	TANK	\$70.72	UN3506	DE
ACC060F	ASTM 60F, 170/500F x 2.0F, Hg, Tol Imm	TANK	\$70.72	UN3506	DE
ACC061C	ASTM 61C, 32/127C x 0.2C, Hg, 79mm Imm	PETROATUM MELTING POINT	\$100.08	UN3506	DE
ACC061F	ASTM 61F, 90/260F x 0.5F, Hg, 79mm Imm	PETROATUM MELTING POINT	\$100.08	UN3506	DE, FR
ACC062C	ASTM 62C, -38/2C x 0.1C, Hg, Tol Imm	PRECISION	\$249.12	UN3506	DE, FR
ACC062F	ASTM 62F, -36/35F x 0.2F, Hg, Tol Imm	PRECISION	\$249.12	UN3506	DE
ACC063C	ASTM 63C, -8/32C x 0.1C, Hg, Tol Imm	PRECISION	\$234.55	UN3506	DE, FR
ACC063F	ASTM 63F, 18/89F x 0.2F, Hg, Tol Imm	PRECISION	\$234.55	UN3506	DE, FR
ACC064C	ASTM 64C, 25/55C x 0.1C, Hg, Tol Imm	PRECISION	\$234.55	UN3506	DE, FR
ACC064F	ASTM 64F, 77/131F x 0.2F, Hg, Tol Imm	PRECISION	\$234.55	UN3506	DE, FR
ACC065C	ASTM 65C, 50/80C x 0.1C, Hg, Tol Imm	PRECISION	\$249.12	UN3506	DE, FR
ACC065F	ASTM 65F, 122/176F x 0.2F, Hg, Tol Imm	PRECISION	\$249.12	UN3506	DE, FR
ACC066C	ASTM 66C, 75/105C x 0.1C, Hg, Tol Imm	PRECISION	\$262.34	UN3506	DE, FR
ACC066F	ASTM 66F, 167/221F x 0.2F, Hg, Tol Imm	PRECISION	\$262.34	UN3506	DE, FR
ACC067C	ASTM 67C, 95/155C x 0.2C, Hg, Tol Imm	PRECISION	\$271.21	UN3506	DE, FR
ACC067F	ASTM 67F, 203/311F x 0.5F, Hg, Tol Imm	PRECISION	\$271.21	UN3506	DE, FR
ACC068C	ASTM 68C, 145/205C x 0.2C, Hg, Tol Imm	PRECISION	\$281.55	UN3506	DE, FR

<u>CAT. NO.</u>	<u>DESCRIPTION</u>		<u>LIST PRICE</u>	<u>Haz (Hg)</u>	<u>COO</u>
ACC068F	ASTM 68F, 293/401F x 0.5F, Hg, Tol Imm	PRECISION	\$281.55	UN3506	DE
ACC069C	ASTM 69C, 195/305C x 0.5C, Hg, Tol Imm	PRECISION	\$321.38	UN3506	DE
ACC069F	ASTM 69F, 383/581F x 1.0F, Hg, Tol Imm	PRECISION	\$321.38	UN3506	DE
ACC070C	ASTM 70C, 295/405C x 0.5C, Hg, Tol Imm	PRECISION	\$357.25	UN3506	DE
ACC070F	ASTM 70F, 563/761F x 1.0F, Tol Imm	PRECISION	\$357.25	UN3506	DE
ACC071C	ASTM 71C, -37/21C x 0.5C, Hg, 76mm Imm	OIL in WAX	\$77.00	UN3506	DE, FR
ACC071F	ASTM 71F, -35/70 x 1.0F, Hg, 76mm Imm	OIL in WAX	\$77.00	UN3506	DE
ACC072C	ASTM 72C, -19.5/16.6C x 0.05C Hg, Tol Imm	DISCONTINUED			
ACC072F	ASTM 72F, -2.5/2.5F x 0.1F, Hg, Tol Imm	DISCONTINUED			
ACC075F	ASTM 75F, -35/35F x 0.5F, Hg, 100mm Imm	ANTIFREEZE FREEZING POINT	\$97.69	UN3506	DE
ACC077F	ASTM 77F, 245/265 x 0.5F, Hg, Tol Imm	SAYBOLT VISCOSITY	\$86.78	UN3506	DE
ACC078F	ASTM 78F, 295/315F x 0.5F, Hg, Tol Imm	SAYBOLT VISCOSITY	\$90.73	UN3506	DE
ACC079F	ASTM 79F, 345/365F x 0.5F, Hg, Tol Imm	SAYBOLT VISCOSITY	\$97.69	UN3506	DE
ACC080F	ASTM 80F, 395/415F x 0.5F, Hg, Tol Imm	SAYBOLT VISCOSITY	\$97.67	UN3506	DE
ACC081F	ASTM 81F, 445/465F x 0.5F, Hg, Tol Imm	SAYBOLT VISCOSITY	\$97.67	UN3506	DE
ACC082C	ASTM 82C, -15/105C x 1.0C, Hg, 30mm Imm	FUEL RATING, ENGINE	\$57.58	UN3506	DE, FR
ACC082F	ASTM 82F, 0/220F x 2.0F, Hg, 30mm Imm	FUEL RATING, ENGINE	\$57.58	UN3506	DE, FR
ACC083C	ASTM 83C, 15/70C x 0.5C, Hg, 40mm Imm	FUEL RATING, AIR	\$57.58	UN3506	DE, FR
ACC083F	ASTM 83F, 60/160F x 1.0F, Hg, 40mm Imm	FUEL RATING, AIR	\$57.58	UN3506	DE, FR
ACC084C	ASTM 84C, 25/80C x 1.0C, Hg, 249mm Imm	FUEL RATING, ORIFICE TANK	\$57.58	UN3506	DE
ACC084F	ASTM 84F, 75/175F x 1.0F, Hg, 249mm Imm	FUEL RATING, ORIFICE TANK	\$57.58	UN3506	DE
ACC085C	ASTM 85C, 40/150C x 1.0C, Hg, 181mm Imm	FUEL RATING, SURGE	\$57.58	UN3506	DE
ACC085F	ASTM 85F, 100/300F x 2.0F, Hg, 181mm Imm	FUEL RATING, SURGE	\$57.58	UN3506	DE
ACC086C	ASTM 86C, 95/175C x 1.0C, Hg, 35mm Imm	FUEL RATING, MIX	\$75.22	UN3506	DE, FR
ACC086F	ASTM 86F, 200/350F x 1.0F, Hg, 35mm Imm	FUEL RATING, MIX	\$75.22	UN3506	DE, FR
ACC087C	ASTM 87C, 150/205C x 1.0C, Hg, 40mm Imm	FUEL RATING, COOLANT	\$85.67	UN3506	DE
ACC087F	ASTM 87F, 300/400F x 1.0F, Hg, 40mm Imm	FUEL RATING, COOLANT	\$85.67	UN3506	DE
ACC088C	ASTM 88C, 10/200C x 1.0C, Hg, 57mm Imm	VEGETABLE OIL FLASH	\$90.67	UN3506	DE, FR
ACC088F	ASTM 88F, 50/392F x 2.0F, Hg, 57mm Imm	VEGETABLE OIL FLASH	\$90.67	UN3506	DE, FR
ACC089C	ASTM 89C, -20/10C x 0.1C, Hg, 76 Imm	SOLIDIFICATION POINT	\$101.93	UN3506	DE, FR
ACC090C	ASTM 90C, 0/30C x 0.1C, Hg, 76mm Imm	SOLIDIFICATION POINT	\$95.32	UN3506	DE, FR
ACC091C	ASTM 91C, 20/50C x 0.1C, Hg, 76mm Imm	SOLIDIFICATION POINT	\$95.32	UN3506	DE
ACC092C	ASTM 92C, 40/70C x 0.1C Hg, 76mm Imm	SOLIDIFICATION POINT	\$95.32	UN3506	DE, FR
ACC093C	ASTM 93C, 60/90C x 0.1C, Hg, 76mm Imm	SOLIDIFICATION POINT	\$105.63	UN3506	DE
ACC094C	ASTM 94C, 80/110C x 0.1C, Hg, 76mm Imm	SOLIDIFICATION POINT	\$105.63	UN3506	DE, FR
ACC095C	ASTM 95C, 100/130C x 0.1C, Hg, 76mm imm	SOLIDIFICATION POINT	\$109.02	UN3506	DE, FR
ACC096C	ASTM 96C, 120/150C x 0.1C, Hg, 76mm Imm	SOLIDIFICATION POINT	\$109.02	UN3506	DE, FR
ACC097C	ASTM 97C, -18/49C x 0.5C, Hg, Tol Imm	TANK	\$55.87	UN3506	DE, FR
ACC097F	ASTM 97F, 0/120F x 1.0F, Hg, Hook Top	TANK	\$55.87	UN3506	DE
ACC098C	ASTM 98C, 16/82C x 0.5C, Hg, Tol Imm	TANK	\$59.67	UN3506	DE

<u>CAT. NO.</u>	<u>DESCRIPTION</u>		<u>LIST PRICE</u>	<u>Haz (Hg)</u>	<u>COO</u>
ACC098F	ASTM 98F, 60/180F x 1.0F, Hg, Tol Imm	TANK	\$59.67	UN3506	DE
ACC100C	ASTM 100C, 145/205C x 0.2C, Hg, 76mm Imm	SOLIDIFICATION POINT	\$123.53	UN3506	DE
ACC101C	ASTM 101C, 195/305C x 0.2C, Hg, 76mm Imm	SOLIDIFICATION POINT	\$113.24	UN3506	DE
ACC102C	ASTM 102C, 123/177C x 0.2, Hg, 100mm Imm	SOLVENTS DISTILLATION	\$98.55	UN3506	DE
ACC103C	ASTM 103C, 148/202C x 0.2C, Hg, 100mm Imm	SOLVENTS DISTILLATION	\$101.63	UN3506	DE
ACC104C	ASTM 104C, 173/227C x 0.2C, Hg, 100mm Imm	SOLVENTS DISTILLATION	\$101.63	UN3506	DE
ACC105C	ASTM 105C, 198/252F x 1.0F, Hg, 40mm Imm	SOLVENTS DISTILLATION	\$107.97	UN3506	DE
ACC106C	ASTM 106C, 223/277C x 0.2C, Hg, 100mm Imm	SOLVENTS DISTILLATION	\$107.97	UN3506	DE
ACC107C	ASTM 107C, 248/302C x 0.2C, Hg, 100mm Imm	SOLVENTS DISTILLATION	\$107.97	UN3506	DE
ACC108F	ASTM 108F, 270/290F x 0.5F, Hg, Tol Imm	DISCONTINUED			
ACC109F	ASTM 109F, 320/340F x 0.5F, Hg, Tol Imm	DISCONTINUED			
ACC110F	ASTM 110F, 272.5/277.5F x 0.1F, Hg, Tol Imm	KINEMATIC VISCOSITY	\$205.84	UN3506	DE
ACC111C	ASTM 111C, 170/250C x 0.2C Hg, 100mm Imm	TAR ACIDS DISTILLATION	\$100.21	UN3506	DE
ACC112C	ASTM 112C, 4/6C x 0.02C, Hg, Tol Imm				DE,
ACC113C	ASTM 113C, -1/175C x 0.5C, Hg, Tol Imm	BITUMINOUS SOFTENING PT	\$75.06	UN3506	FR
ACC113F	ASTM 113F, 30/350F x 1.0F,Hg, Tol Imm	BITUMINOUS SOFTENING PT	\$75.06	UN3506	DE
ACC114C	ASTM 114C, -80/20C x 0.5C, Non-Hg, Tol Imm	AVIATION FUEL FREEZING PT ⁴	\$77.97		DE
ACC114F	ASTM 114F, -112/70F x 1.0F, Non-Hg, Tol Imm	AVIATION FUEL FREEZING PT ⁴	\$77.97		DE
ACC115C	ASTM 115U Beckmann 0/6	DISCONTINUED			
ACC116C	ASTM 116C, 18.9/25.1C x 0.01C, Hg, Tol Imm	BOMB CALORIMETER	\$338.06	UN3506	DE
ACC117C	ASTM 117C, 23.9/30.1C x 0.01C, Hg, Tol Imm	BOMB CALORIMETER	\$338.06	UN3506	DE
ACC118C	ASTM 118C, 28.6/31.4 x 0.05C, Hg, Tol Imm	KINEMATIC VISCOSITY	\$193.66	UN3506	DE
ACC118F	ASTM 118F, 83.5/88.5 x 0.1F, Hg, Tol Imm	KINEMATIC VISCOSITY	\$193.66	UN3506	DE
ACC119C	ASTM 119C, -38.3/-30C x 0.1C, Hg, 100mm Imm	DISCONTINUED			
ACC119F	ASTM 119F, -37/-22F x 0.2F, Hg, 100mm Imm	FREEZING POINT	\$239.97	UN3506	DE
ACC120C	ASTM 120C, 38.6/41.4C x 0.05C, Hg, Tol Imm	KINEMATIC VISCOSITY	\$190.87	UN3506	DE, FR
ACC121C	ASTM 121C 98.6/101.6C x 0.05C, Hg, Tol Imm	KINEMATIC VISCOSITY	\$190.87	UN3506	DE, FR
ACC123C	ASTM 123C, -35/-25C x 0.1C, Hg, Tol Imm	BROOKFIELD VISCOSITY	\$171.81	UN3506	DE
ACC124C	ASTM 124C -25/-15C x 0.1C, Hg, Tol Imm	BROOKFIELD VISCOSITY	\$171.81	UN3506	DE, FR
ACC125C	ASTM 125C -15/-5C x 0.1C, Hg, Tol Imm	BROOKFIELD VISCOSITY	\$171.81	UN3506	DE, FR
ACC126C	ASTM 126C, -27.4/-24.6C x 0.05C, Hg, Tol Imm	KINEMATIC VISCOSITY	\$237.18	UN3506	DE, FR
ACC126F	ASTM 126F, -17.5/-12.5 x 0.1F, Hg, Tol Imm	KINEMATIC VISCOSITY	\$237.18	UN3506	DE
ACC127C	ASTM 127C, -21.4/-18.6C x 0.05C, Hg, Tol Imm	KINEMATIC VISCOSITY	\$237.18	UN3506	DE, FR
ACC128C	ASTM 128C, -1.4/1.4C x 0.05C, Hg, Tol Imm	KINEMATIC VISCOSITY	\$203.18	UN3506	DE
ACC128F	ASTM 128F, 29.5/34.5 x 0.1F, Hg, Tol Imm	KINEMATIC VISCOSITY	\$203.18	UN3506	DE
ACC129C	ASTM 129C, 91.6/94.4C x 0.05C, Hg, Tol Imm	KINEMATIC VISCOSITY	\$186.79	UN3506	DE
ACC129F	ASTM 129F, 197.5/202.5F x 0.1F, Hg, Tol Imm	KINEMATIC VISCOSITY	\$186.79	UN3506	DE
ACC130C	ASTM 130C, -7/105C x 0.5C, Hg, Tol Imm	TANK	\$58.46	UN3506	DE
ACC130F	ASTM 130F 20/220F x 1.0F, Hg, Hook Top	TANK	\$53.59	UN3506	DE
ACC132C	ASTM 132C, 148.6/151.4 x 0.05C, Hg, Tol Imm	KINEMATIC VISCOSITY	\$199.48	UN3506	DE
ACC133C	ASTM 133C, -38/2C x 0.1C, Hg, 76mm Imm	PRECISION, 76mm IMMERSION	\$140.46	UN3506	DE, FR
ACC134C	ASTM 134C 144/156C x 0.2C, Hg, 100mm imm	SLUDGE	\$93.58	UN3506	DE
ACC135C	ASTM 135C, 38/93C x 1.0C, Hg, 40mm Imm	FUEL RATING, AIR-HIGH	\$85.12	UN3506	DE
ACC135F	ASTM 135F, 100/200F x 1.0F Hg, 40mm Imm	FUEL RATING, AIR-HIGH	\$85.12	UN3506	DE
ACC136C	ASTM 136C, -20/60C x 0.2C, Hg, Tol Imm	AVIATION FUEL DENSITY	\$105.63	UN3506	DE, FR
ACC136F	ASTM 136F, -5/140F x 0.5F, Hg, Tol Imm	AVIATION FUEL DENSITY	\$105.63	UN3506	DE
ACC137C	ASTM 137C 80/100C x 0.1C, Hg, 76mm Imm	OXIDATION CELL TEST	\$136.33	UN3506	DE

<u>CAT. NO.</u>	<u>DESCRIPTION</u>	<u>LIST PRICE</u>	<u>Haz (Hg)</u>	<u>COO</u>
<u>ASTM APPROVED NON-MERCURY, BLUE SPIRIT THERMOMETERS</u>				
<u>ASTM PRECISIONS</u>				
ACC062CBLS	0S62C, -38/+2C, BLUE SPIRIT	\$297.24		DE
ACC062FBLS	0S62F, -36/+35F, BLUE SPIRIT	\$297.24		DE
ACC063CBLS	0S63C, -8/+32C, BLUE SPIRIT	\$264.06		DE
ACC063FBLS	0S63F, +18/89F, BLUE SPIRIT	\$264.06		DE
ACC064CBLS	0S64C, +25/55C, BLUE SPIRIT	\$288.87		DE
ACC064FBLS	0S64F, +77/131F, BLUE SPIRIT	\$288.87		DE
ACC065CBLS	0S65C, +50/80C, BLUE SPIRIT	\$297.24		DE
ACC065FBLS	0S65F, +122/176F, BLUE SPIRIT	\$297.24		DE
ACC066CBLS	0S66C, +75/105C, BLUE SPIRIT	\$315.46		DE
ACC066FBLS	0S66F, +167/221F, BLUE SPIRIT	\$315.46		DE
ACC067CBLS	0S67C, +95/155C, BLUE SPIRIT	\$315.46		DE
ACC067FBLS	0S67F, +203/311F, BLUE SPIRIT	\$315.46		DE
<u>ASTM APPROVED NON-MERCURY, BLUE SPIRIT THERMOMETERS</u>				
<u>ASTM APPLICATIONS</u>				
ACC005CBLS	-38/50C x 1.0C	\$121.17		DE
ACC005FBLS	-36/+120F x 2.0F	\$121.17		DE
ACC012CBLS	-20/102C x 0.2C	\$208.85		DE
ACC012FBLS	-5/215F x 0.5F	\$208.85		DE
ACC015CBLS	-2/80C x 0.2C	\$159.81		DE
ACC015FBLS	+30/180F x 0.5F	\$159.81		DE
ACC018CBLS	+34/42C x 0.1C	\$143.50		DE
ACC018FBLS	+94/108F x 0.2F	\$143.50		DE
ACC022CBLS	+98/103C x 0.1C	\$159.81		DE
ACC022FBLS	+204/218F x 0.2F	\$159.81		DE
ACC056CBLS	+19/35C x 0.02C	\$621.95		DE
ACC056FBLS	+66/95F x 0.05F	\$621.95		DE
ACC058CBLS	-34/+49C x 0.5C	\$124.91		DE
ACC058FBLS	-30/+120F x 1.0F	\$124.91		DE
ACC091CBLS	+20/50C x 0.1C	\$172.42		DE
ACC116CBLS	+18.9/ 25.1C x 0.01C	\$690.95		DE
ACC117CBLS	+23.9/30.1C x 0.01C	\$699.08		DE
ACC120CBLS	+38.6/41.4C x 0.05	\$316.13		DE
ACC130CBLS	-7/+105C X 0.5C	\$123.85		DE
ACC130FBLS	+20/220F x 1.0F	\$123.85		DE
<u>ASTM LIKE NON-MERCURY, BLUE SPIRIT THERMOMETERS</u>				
<u>ASTM APPLICATIONS</u>				
ACC001CBLS	-20/150C x 1.0C	\$127.17		DE
ACC001FBLS	0/302F x 2.0F	\$127.17		DE
ACC009CS	-5/110C x 0.5C	\$104.04		DE
ACC009FS	+20/230F x 1.0F	\$104.04		DE
ACC013CBLS	+155/170C x 0.5	\$170.21		DE
ACC014CBLS	+38/82C x 0.1C	\$198.16		DE
ACC014FBLS	+100/180F x 0.2F	\$198.16		DE
ACC016CBLS	+30/200C x 0.5C	\$145.41		DE
ACC016FBLS	+85/392F 1.0F	\$145.41		DE
ACC017CBLS	+19/27C x 0.1C	\$129.59		DE
ACC017FBLS	+66/80F x 0.2F	\$129.59		DE
ACC019CBLS	+49/57C x 0.1C	\$141.96		DE
ACC019FBLS	+120//134F x 0.2F	\$141.96		DE

<u>CAT. NO.</u>	<u>DESCRIPTION</u>	<u>LIST PRICE</u>	<u>Haz (Hg)</u>	<u>COO</u>
ACC020CBLS	+57/65C x 0.1C	\$143.12		DE
ACC020FBLS	+134/148F x 0.2F	\$143.12		DE
ACC021CBLS	+79/87C x 0.1C	\$147.45		DE
ACC024CBLS	+39/54C x 0.2C	\$134.48		DE
ACC025CBLS	+95/105C x 0.2C	\$141.96		DE
ACC026CBLS	+130/140C x 0.1C	\$204.52		DE
ACC027CBLS	+147/182C x 0.5C	\$171.26		DE
ACC028CBLS	+36.6/39.4 x 0.05C	\$254.28		DE
ACC028FBLS	+97.5/102.5F x 0.1F	\$271.56		DE
ACC029CBLS	+52.6/55.4C x 0.05C	\$267.09		DE
ACC029FBLS	+127.5/132.5F x 0.1F	\$267.09		DE
ACC030FBLS	+207.5/212.5F x 0.1F	\$324.14		DE
ACC033CBLS	-38/42C x 0.2C	\$215.78		DE
ACC033FBLS	-36.5/107.5F x 0.5F	\$215.78		DE
ACC034CBLS	+25/105C x 0.2C	\$177.24		DE
ACC034FBLS	+77/221F x 0.5F	\$177.24		DE
ACC035FBLS	+194/338F x 0.5F	\$212.11		DE
ACC036CBLS	-2/68C x 0.2C	\$136.75		DE
ACC037CBLS	-2/52C x 0.2C	\$131.09		DE
ACC038CBLS	+24/78C x 0.2C	\$144.58		DE
ACC039CBLS	+48/102C x 0.2C	\$161.00		DE
ACC040CBLS	+72/126C x 0.2C	\$172.95		DE
ACC042CBLS	+95/255C x 0.5C	\$193.87		DE
ACC044CBLS	+18.6/21.4C x 0.05C	\$248.77		DE
ACC044FBLS	+66.5/71.5F x 0.1F	\$248.77		DE
ACC045CBLS	+23.6/26.4C x 0.05C	\$248.77		DE
ACC045FBLS	+74.5/79.5F x 0.1F	\$248.77		DE
ACC046CBLS	+48.6/51.4C x 0.05C	\$269.72		DE
ACC046FBLS	+119.5/124.5F x 0.1F	\$269.72		DE
ACC047CBLS	+58.6/61.4C x 0.05C	\$281.51		DE
ACC047FBLS	+137.5/142.5F x 0.1F	\$281.51		DE
ACC048CBLS	+80.6/83.4C x 0.05C	\$301.13		DE
ACC048FBLS	+177.5/182.5F x 0.1F	\$301.13		DE
ACC049CBLS	+20/70C x 0.2C	\$141.26		DE
ACC052CBLS	-10/5C x 0.1C	\$136.91		DE
ACC054CBLS	+20/100.6C x 0.2C	\$163.06		DE
ACC054FBLS	+68/213F x 0.5F	\$163.06		DE
ACC057CBLS	-20/50C x 0.5C	\$135.26		DE
ACC057FBLS	-4/122F x 1.0F	\$135.26		DE
ACC059CBLS	-8/82C x 0.5C	\$122.63		DE
ACC059FBLS	0/180F x 1.0F	\$168.99		DE
ACC061CBLS	+32/127C x 0.2C	\$168.99		DE
ACC061FBLS	+90/260F x 0.5F	\$157.13		DE
ACC071CBLS	-37/21XC x 0.5C	\$157.13		DE
ACC071FBLS	-35/70F x 1.0F	\$306.38		DE
ACC072CBLS	-19.4/16.6C x 0.05C	\$306.38		DE
ACC072FBLS	-2.5/2.5F x 0.1F	\$278.05		DE
ACC073CBLS	-41.4/-38.6C X 0.05C	\$278.05		DE
ACC073FBLS	-42.5/-37.5F X 0.1F	\$195.09		DE
ACC075FBLS	-35/35F x 0.5F	\$182.02		DE
ACC078FBLS	+245/265F x 0.5F	\$195.09		DE
ACC079FBLS	+295/315F x 0.5F	\$101.37		DE
ACC083CBLS	+15/70C x 1.0C	\$101.37		DE
ACC083FBLS	+60/160F x 1.0F	\$88.83		DE
ACC084CBLS	+25/80C x 1.0C	\$88.83		DE

<u>CAT. NO.</u>	<u>DESCRIPTION</u>	<u>LIST PRICE</u>	<u>Haz (Hg)</u>	<u>COO</u>
ACC084FBLS	+75/175F x 1.0F	\$107.80		DE
ACC085CBLS	+40/150C x 1.0C	\$107.80		DE
ACC085FBLS	+100/300F x 2.0F	\$138.72		DE
ACC088CBLS	+10/200C x 1.0C	\$138.72		DE
ACC088FBLS	+50/392F x 2.0F	\$215.78		DE
ACC089CBLS	-20/10C x 0.1C	\$166.47		DE
ACC090CBLS	0/30C x 0.1C	\$227.23		DE
ACC092CBLS	+40/70C x 0.1C	\$197.72		DE
ACC093CBLS	+60/90C x 0.1C	\$203.32		DE
ACC094CBLS	+80/110C x 0.1C	\$216.03		DE
ACC095CBLS	+100/130C x 0.1C	\$228.08		DE
ACC096CBLS	120/150C x 0.1C	\$102.73		DE
ACC097CBLS	-18/49C x 0.5C	\$102.73		DE
ACC097FBLS	0/120F x 1.0F	\$93.02		DE
ACC098CBLS	+16/82C x 0.5C	\$93.02		DE
ACC098FBLS	+60/180F x 1.0F	\$369.77		DE
ACC099CBLS	-50/5C x 0.2C	\$369.77		DE
ACC099FBLS	-58/41F x 0.5F	\$188.54		DE
ACC102CBLS	+123/177C x 0.2C	\$197.72		DE
ACC103CBLS	+148/202C x 0.2C	\$349.56		DE
ACC110CBLS	+133.6/136.4C x 0.05C	\$327.32		DE
ACC110FBLS	+272.5/277.5F x 0.1F	\$476.67		DE
ACC112CBLS	+4/6C x 0.02C	\$166.30		DE
ACC113CBLS	-1/175C x 0.5C	\$166.30		DE
ACC113FBLS	+30/350F x 1.0F	\$355.08		DE
ACC118CBLS	+28.6/31.4C x 0.05C	\$276.31		DE
ACC119CBLS	-38.3/-30.0C x 0.05C	\$316.13		DE
ACC120CBLS	+38.6/41.4C x 0.05C	\$319.35		DE
ACC121CBLS	+98.6/101.4C x 0.05C	\$249.63		DE
ACC122CBLS	-45/-35C x 0.1C	\$227.79		DE
ACC123CBLS	-35/-25C x 0.1C	\$208.77		DE
ACC124CBLS	-25/-15C x 0.1C	\$207.14		DE
ACC125CBLS	-15/-5C x 0.1C	\$322.44		DE
ACC126CBLS	-27.4/-24.6C x 0.05C	\$322.44		DE
ACC126FBLS	-17.5/-12.5F x 0.1F	\$379.04		DE
ACC127CBLS	-21.4/-18.6C x 0.05C	\$188.53		DE
ACC128CBLS	-1.4/+1.4C x 0.05C	\$188.53		DE
ACC128FBLS	+29.5/34.5F x 0.1F	\$288.46		DE
ACC129CBLS	+91.6/94.4F x 0.05C	\$288.46		DE
ACC129FBLS	+197.5/202.5F x 0.1F	\$289.68		DE
ACC132CBLS	+148.6/151.4C x 0.05C	\$284.27		DE
ACC133CBLS	-38/+2C x 0.1C	\$231.14		DE
ACC134CBLS	144/156C x 0.2C	\$136.40		DE
ACC135CBLS	+39/93C x 1.0C	\$122.17		DE
ACC135FBLS	+100/200F x 1.0F	\$197.04		DE
ACC136CBLS	-20/+60C x 0.2C	\$197.04		DE
ACC136FBLS	-5/140F x 0.5F	\$175.19		DE
ACC137CBLS	+80/100C x 0.1C	\$170.08		DE

ASTM TANK GAUGING THERMOMETERS

WOODBACK CUP CASE, MERCURY-FILLED THERMOMETERS

ACC272WB1	W/ASTM 58F -30/120F	\$103.65	UN3506	DE
ACC272WB2	W/ASTM 59F 0/180F	\$103.65	UN3506	DE
ACC272WB3	W/ASTM 130F 20/220F	\$103.65	UN3506	DE
ACC272WB4	W/ASTM 60F 170/500F	\$109.07	UN3506	DE

<u>CAT. NO.</u>	<u>DESCRIPTION</u>	<u>LIST PRICE</u>	<u>Haz (Hg)</u>	<u>COO</u>
ACC273WB1	W/ASTM 58C -34/+49C	\$103.65	UN3506	DE
ACC273WB2	W/ASTM 59C -18/82C	\$103.65	UN3506	DE
ACC273WB3	W/ASTM 130C -7/105C	\$103.65	UN3506	DE
ACC273WB4	W/ASTM 60C 77/260C	\$109.07	UN3506	DE
ACCWB	CASE ONLY	\$59.48		DE

WOODBACK CUP CASE: SPIRIT-FILLED THERMOMETERS

ACC272WB1S	W/ASTM Like 58F -30/120F	\$133.70		DE
ACC272WB2S	W/ASTM Like 59F 0/180F	\$133.70		DE
ACC272WB3S	W/ASTM Like 130F 20/220F	\$133.70		DE
ACC273WB1S	W/ASTM Like 58C -34/+49C	\$133.70		DE
ACC273WB2S	W/ASTM Like 59C -18/82C	\$133.70		DE
ACC273WB3S	W/ASTM Like 130C -7/105C	\$133.70		DE
ACCWB	CASE ONLY	\$59.48		DE

"SCOOP-MASTER" CUP CASE, MERCURY-FILLED THERMOMETERS

ACC272S1	W/ASTM 58F -30/120F	\$118.60	UN3506	DE
ACC272S2	W/ASTM 59F 0/180F	\$118.60	UN3506	DE
ACC272S3	W/ASTM 130F 20/220F	\$118.60	UN3506	DE
ACC273S1	W/ASTM 58C -34/+49	\$118.60	UN3506	DE
ACC273S2	W/ASTM 59C -18/82C	\$118.60	UN3506	DE
ACC273S3	W/ASTM 130C -7/105C	\$118.60	UN3506	DE
ACC272CASE	CASE ONLY	\$80.17		DE

"SCOOP-MASTER" CUP CASE, SPIRIT-FILLED THERMOMETERS

ACC272S1S	W/ASTM Like 58F -30/120F	\$164.91		DE
ACC272S2S	W/ASTM Like 59F 0/180F	\$164.91		DE
ACC272S3S	W/ASTM Like 130F 20/220F	\$164.91		DE
ACC273S1S	W/ASTM Like 58C -34/+49	\$164.91		DE
ACC273S2S	W/ASTM Like 59C -18/82C	\$164.91		DE
ACC273S3S	W/ASTM Like 130C -7/105C	\$164.91		DE
ACC272CASE	CASE ONLY	\$80.17		DE

ASTM CERTIFIED THERMOMETERS, Mercury*

CERTIFIED TEMPERATURES

ACC001CFC	ASTM 1C, -20/150C x 1.0C, Hg, 76mm Imm	-20, 0, 50, 100, 150C	\$378.95	UN3506
ACC001FFC	ASTM 1F, 0/302F x 2.0F, Hg, 76mm Imm	0, 32, 122, 212, 302F	\$378.95	UN3506
ACC002CFC	ASTM 2C, -5/300C x 1.0C, Hg, 76mm Imm	0, 75, 150, 225, 300C	\$376.31	UN3506
ACC002FFC	ASTM 2F, 20/580F X 2.0F, Hg, 76mm Imm	32, 150, 300, 450, 580F	\$376.31	UN3506
ACC003CFC	ASTM 3C, -5/400C x 1.0C, Hg, 76mm Imm	0, 100, 200, 300, 370C	\$381.06	UN3506
ACC003FFC	ASTM 3F, 20/760F x 2.0F, Hg, 76mm Imm	32, 200, 370, 540, 700F	\$381.06	UN3506
ACC004CFC	ASTM 4C, -1/105C x 0.5C, Hg, 152mm Imm	0, 50, 100C	\$253.62	UN3506
ACC004FFC	ASTM 4F, 30/220F x 1.0F, Hg, 152mm Imm	32, 122, 212F	\$253.62	UN3506
ACC005CFC	ASTM 5C, -38/50C x 1.0C, Hg, 108mm Imm	-35, 0, 50C	\$248.18	UN3506
ACC005FFC	ASTM 5F, -36/120F x 2.0F, Hg, 108mm Imm	-30, 32, 120F	\$248.18	UN3506
ACC006CFC	ASTM 6C, -80/20C x 1.0C, Non-Hg, 76mm Imm	-70, -35, 0, 20C	\$335.41	
ACC006FFC	ASTM 6F, -112/70F x 2.0F, Non-Hg, 76mm Imm	-94, -30, 32, 70F	\$335.41	
ACC007CFC	ASTM 7C, -2/300C x 1.0C, Hg, Tol Imm	0, 50, 100, 150, 200, 250, 300C	\$507.21	UN3506
ACC007FFC	ASTM 7F, 30/580F x 2.0F, Hg, Tol Imm	-32, 100, 200, 300, 400, 500, 570F	\$507.21	UN3506
ACC008CFC	ASTM 8C, -2/400C x 1.0C, Hg, Tol Imm	0, 100, 200, 300, 370C	\$376.31	UN3506
ACC008FFC	ASTM 8F, 30/760F x 2.0F, Hg, Tol Imm	32, 200, 370, 540, 700F	\$376.31	UN3506
ACC009CFC	ASTM 9C, -5/110C x 0.5C, Hg, 57mm Imm	0, 35, 70, 105C	\$305.37	UN3506
ACC009FFC	ASTM 9F, 20/230F x 1.0F, Hg, 57mm Imm	32, 100, 160, 220F	\$305.37	UN3506
ACC010CFC	ASTM 10C, 90/370C x 2.0C, Hg, 57mm Imm	100, 200, 300, 370C	\$340.83	UN3506
ACC010FFC	ASTM 10F, 200/700F x 5.0F, Hg, 57mm Imm	212, 390, 570, 700F	\$340.83	UN3506
ACC011CFC	ASTM 11C, -6/400C x 2.0C, Hg, 25mm Imm	0, 100, 200, 300, 370C	\$406.22	UN3506

<u>CAT. NO.</u>	<u>DESCRIPTION</u>		<u>LIST PRICE</u>	<u>Haz (Hg)</u>	<u>COO</u>
ACC011FFC	ASTM 11F, 20/760F x 5.0F, Hg, 25mm Imm	32, 200, 370, 540, 700F	\$406.22	UN3506	
ACC012CFC	ASTM 12C, -20/102C x 0.2C, Hg, Tol Imm	-20, -10, 0, 10, 20, 30, 40, 50, 60, 70, 80, 90	\$937.90	UN3506	
ACC012FFC	ASTM 12F, -5/215F x 0.5F, Hg, Tol Imm	-5, 15, 32, 60, 85, 110, 135, 160, 185, 210F	\$937.90	UN3506	
ACC013CFC	ASTM 13C, 155/170C x 0.5C, Hg, Tol Imm	155, 163, 170C	\$291.72	UN3506	
ACC014CFC	ASTM 14C, 38/82C x 0.1C, Hg, 79mm Imm	40, 50, 60, 70, 80C	\$384.52	UN3506	
ACC014FFC	ASTM 14F, 100/180F x 0.2F Hg, 79mm Imm	100, 120, 140, 160, 180F	\$384.52	UN3506	
ACC015CFC	ASTM 15C, -2/80C x 0.2C, Hg, Tol Imm	0, 20, 40, 60, 80C	\$381.73	UN3506	
ACC015FFC	ASTM 15F, 30/180F x 0.5F, Hg, Tol Imm	32, 70, 100, 140, 180F	\$381.73	UN3506	
ACC016CFC	ASTM 16C, 30/200C x 0.5C, Hg, Tol Imm	30, 60, 90, 120, 150, 180, 200C	\$515.29	UN3506	
ACC016FFC	ASTM 16F, 85/392F x 1.0F, Hg, Tol Imm	90, 140, 190, 240, 290, 340, 390F	\$515.29	UN3506	
ACC017CFC	ASTM 17C, 19/27C x 0.1C, Hg, Tol Imm	21, 25C	\$180.00	UN3506	
ACC017FFC	ASTM 17F, 66/80F x 0.2F, Hg, Tol Imm	70, 77F	\$180.00	UN3506	
ACC018CFC	ASTM 18C, 34/42C x 0.1C, Hg, Tol Imm	38, 41C	\$180.00	UN3506	
ACC018FFC	ASTM 18F, 94/108F x 0.2F, Hg, Tol Imm	100, 107F	\$180.00	UN3506	
ACC019CFC	ASTM 19C, 49/57C x 0.1C, Hg, Tol Imm	50, 54C	\$180.00	UN3506	
ACC019FFC	ASTM 19F, 120/134F x 0.2F, Hg, Tol Imm	122, 130F	\$180.00	UN3506	
ACC020CFC	ASTM 20C, 57/65C x 0.1C, Hg, Tol Imm	60, 64C	\$180.00	UN3506	
ACC020FFC	ASTM 20F, 134/148F x 0.2F, Hg, Tol Imm	140, 147F	\$180.00	UN3506	
ACC021CFC	ASTM 21C, 79/87C x 0.1C, Hg, Tol Imm	82, 86C	\$180.00	UN3506	
ACC021FFC	ASTM 21F, 174F/188F x 0.2F, Hg, Tol Imm	180, 187F	\$180.00	UN3506	
ACC022CFC	ASTM 22C, 95/103C x 0.1C, Hg, Tol Imm	99, 102C	\$180.00	UN3506	
ACC022FFC	ASTM 22F, 204/218F X 0.2F, Hg, Tol Imm	210, 212F	\$180.00	UN3506	
ACC023CFC	ASTM 23C, 18/28C x 0.2C, Hg, 90mm Imm	20, 25C	\$209.94	UN3506	
ACC024CFC	ASTM 24C, 39/54C x 0.2C, Hg, 90mm Imm	40, 50C	\$209.94	UN3506	
ACC025CFC	ASTM 25C, 95/105C x 0.2C, Hg, 90mm Imm	95, 100C	\$209.94	UN3506	
ACC026CFC	ASTM 26C, 130/140C x 0.1C, Hg, Tol immersion	130, 135, 140C	\$297.16	UN3506	
ACC027CFC	ASTM 27C, 147/182C x 0.5C, Hg, 76mm Imm	155, 165, 175C	\$286.30	UN3506	
ACC028CFC	ASTM 28C, 36.6/39.4C x 0.05C, Hg, Tol Imm	0, 37.8, 39C	\$378.95	UN3506	
ACC028FFC	ASTM 28F, 97.5/102.5F x 0.1F, Hg, Tol Imm	32, 100, 102F	\$378.95	UN3506	
ACC029CFC	ASTM 29C, 52.6/55.4C x 0.05C, Hg, Tol Imm	0, 54.4, 55C	\$378.95	UN3506	
ACC029FFC	ASTM 29F, 127.5/132.5F x 0.1F, Tol, Hg	32, 130, 132F	\$378.95	UN3506	
ACC030FFC	ASTM 30F, 207.5/212.5F x 0.1F, Hg, Tol Imm	32, 210, 212F	\$378.95	UN3506	
ACC031FFC	ASTM 31F, -30/120F x 1.0F, Hg, Tol Imm	-20, 32, 100F	\$267.24	UN3506	
ACC033CFC	ASTM 33C, -38/42C x 0.2C, Hg, 50mm Imm	-35, -20, 0, 20, 40C	\$439.03	UN3506	
ACC033FFC	ASTM 33F, -36.5/107.5F x 0.2F, Hg, 50mm Imm	-31, -4, 32, 68, 104F	\$439.03	UN3506	
ACC034CFC	ASTM 34C, 25/105C x 0.2C, Hg, 50mm Imm	25, 45, 65, 85, 100C	\$444.47	UN3506	
ACC034FFC	ASTM 34F, 77/221F x 0.5F, Hg, 50mm Imm	77, 113, 149, 185, 212F	\$444.47	UN3506	
ACC035CFC	ASTM 35C, 90/170C x 0.2C 50mm Hg	100, 120, 140, 160, 170C	\$411.78	UN3506	
ACC035FFC	ASTM 35F, 194/338f x 0.5F, Hg, 50mm Imm	212, 250, 285, 320, 338F	\$411.78	UN3506	
ACC036CFC	ASTM36C, -2/68C x 0.2C, Hg, 45mm Imm	0, 15, 30, 45, 65C	\$398.15	UN3506	
ACC037CFC	ASTM 37C, -2/52C x 0.2C, Hg, 100mm Imm	0, 15, 30, 50C	\$340.83	UN3506	
ACC038CFC	ASTM 38C, 24/78C x 0.2C, Hg, 100mm Imm	25, 40, 55, 75C	\$340.83	UN3506	
ACC039CFC	ASTM 39C, 48/102C x 0.2, Hg, 100mm Imm	50, 65, 80, 100C	\$340.83	UN3506	
ACC040CFC	ASTM 40C, 72/126C x 0.2C, Hg, 101mm Imm	75, 90, 105, 125C	\$334.08	UN3506	
ACC041CFC	ASTM 41C, 98/152C x 0.2C, Hg, 100mm Imm	100, 115, 130, 150C	\$346.26	UN3506	
ACC042CFC	ASTM 42C, 95/255C x 0.5C, Hg, 100mm Imm	100, 150, 200, 250C	\$346.26	UN3506	
ACC044CFC	ASTM 44C, 18.6/21.4C x 0.05C, Hg, Tol Imm	0, 20, 21C	\$373.53	UN3506	
ACC044FFC	ASTM 44F, 66.5/71.5F x 0.1F, Hg, Tol Imm	32, 68, 70F	\$373.53	UN3506	
ACC045CFC	ASTM 45C, 23.6/26.4C x 0.05C, Hg, Tol Imm	0, 25, 26C	\$373.53	UN3506	
ACC045FFC	ASTM 45F, 74.5/79.5 x 0.1F, Hg, Tol Imm	32, 77, 79F	\$373.53	UN3506	
ACC046CFC	ASTM 46C, 48.6/51.4C x 0.05C, Hg, Tol Imm	0, 50, 51C	\$373.53	UN3506	
ACC046FFC	ASTM 46F, 119.5/124.5F x 0.1F, Hg, Tol Imm	32, 122, 124F	\$373.53	UN3506	
ACC047CFC	ASTM 47C, 58.6/61/4C x 0.05C, Hg, Tol Imm	0, 60, 61C	\$373.53	UN3506	
ACC047FFC	ASTM 47F, 137.5/142.5F x 0.1F, Hg, Tol Imm	32, 140, 142F	\$373.53	UN3506	
ACC048CFC	ASTM 48C, 80.6/83.4C x 0.05C, Hg, Tol Imm	0, 82.2, 83C	\$373.53	UN3506	

<u>CAT. NO.</u>	<u>DESCRIPTION</u>	<u>LIST PRICE</u>	<u>Haz (Hg)</u>	<u>COO</u>
ACC048FFC	ASTM 48F, 177.5/182.5F x 0.1F, Hg, Tol Imm	32, 180, 182F	\$373.53	UN3506
ACC049CFC	ASTM 49C, 20/70C x 0.2C, Hg, 65mm Imm	20, 35, 50, 70C	\$329.99	UN3506
ACC050FFC	ASTM 50F, 54/101F x 0.1F, Hg, Tol Imm	55, 60, 65, 70, 75, 80F	\$534.47	UN3506
ACC051FFC	ASTM 51F, 69/116F x 0.1F, Hg, Tol Imm	70, 75, 80, 85, 90, 95F	\$534.47	UN3506
ACC053CFC	ASTM 53C, -0.6/10.4C x 0.1C, Hg, Tol Imm	0, 5, 10C	\$299.93	UN3506
ACC054CFC	ASTM 54C, 20/100.6C x 0.2C, Hg, Tol Imm	20, 50, 75, 100C	\$340.83	UN3506
ACC054FFC	ASTM 54F, 68/213F x 0.5F, Hg, Tol Imm	70, 120, 170, 210F	\$340.83	UN3506
ACC056CFC	ASTM 56C, 19/35C x 0.02C, Hg, Tol Imm	19, 21, 23, 25, 27, 29, 31, 33, 35C	\$861.67	UN3506
ACC056FFC	ASTM 56F, 66/95F x 0.05F, Hg, Tol Imm	66, 70, 74, 78, 82, 88, 92, 95F	\$796.15	UN3506
ACC057CFC	ASTM 57C, -20/50C x 0.5F, Hg, 57mm Imm	-20, 0, 25, 50C	\$332.62	UN3506
ACC057FFC	ASTM 57F, -4/122F x 1.0F, Hg, 57mm Imm	-3, 32, 77, 122F	\$332.62	UN3506
ACC058CFC	ASTM 58C, -34/49C x 0.5C, Hg, Tol Imm	-30, 0, 25, 45C	\$316.35	UN3506
ACC058FFC	ASTM 58F, -30/120F x 1.0F, Hg, Tol Imm	-20, 32, 80, 120F	\$316.35	UN3506
ACC059CFC	ASTM 59C, -18/82C x 0.5C, Hg, Tol Imm	0, 25, 55, 80C	\$302.70	UN3506
ACC059FFC	ASTM 59F, 0/180F x 1.0F, Hg, Tol Imm	32, 80, 130, 180F	\$302.70	UN3506
ACC060CFC	ASTM 60C, 77/260C x 1.0C, Hg, Tol Imm	100, 175, 255C	\$259.04	UN3506
ACC060FFC	ASTM 60F, 170/500F x 2.0F, Hg, Tol Imm	212, 350, 490F	\$259.04	UN3506
ACC061CFC	ASTM 61C, 32/127C x 0.2C, Hg, 79mm Imm	40, 60, 80, 100, 120C	\$409.00	UN3506
ACC061FFC	ASTM 61F, 90/260F x 0.5F, Hg, 79mm Imm	100, 150, 200, 250F	\$346.26	UN3506
ACC062CFC	ASTM 62C, -38/2C x 0.1C, Hg, Tol Imm	-37, -30, -20, -10, 0C	\$422.63	UN3506
ACC062FFC	ASTM 62F, -36/35F x 0.2F, Hg, Tol Imm	-35, -15, 0, 15, 32F	\$422.63	UN3506
ACC063CFC	ASTM 63C, -8/32C x 0.1C, Hg, Tol Imm	-7, 0, 10, 20, 30C	\$361.34	UN3506
ACC063FFC	ASTM 63F, 18/89F x 0.2F, Hg, Tol Imm	20, 32, 50, 70, 88F	\$361.34	UN3506
ACC064CFC	ASTM 64C, 25/55C x 0.1C, Hg, Tol Imm	0, 25, 35, 45, 55C	\$406.22	UN3506
ACC064FFC	ASTM 64F, 77/131F x 0.2F, Hg, Tol Imm	32, 80, 95, 115, 130F	\$406.22	UN3506
ACC065CFC	ASTM 65C, 50/80C x 0.1C, Hg, Tol Imm	0, 50, 60, 70, 80C	\$410.43	UN3506
ACC065FFC	ASTM 65F, 122/176F x 0.2F, Hg, Tol Imm	32, 125, 145, 160, 175F	\$410.43	UN3506
ACC066CFC	ASTM 66C, 75/105C x 0.1C, Hg, Tol Imm	0, 75, 85, 95, 105C	\$441.69	UN3506
ACC066FFC	ASTM 66F, 167/221F x 0.2F, Hg, Tol Imm	32, 168, 185, 200, 220F	\$441.69	UN3506
ACC067CFC	ASTM 67C, 95/155C x 0.2C, Hg, Tol Imm	0, 100, 110, 130, 150C	\$474.39	UN3506
ACC067FFC	ASTM 67F 203/311F x 0.5F, Hg, Tol Imm	32, 205, 240, 275, 310F	\$474.39	UN3506
ACC068CFC	ASTM 68C, 145/205C x 0.2C, Hg, Tol Imm	0, 150, 170, 190, 205C	\$668.03	UN3506
ACC068FFC	ASTM 68F, 293/401F x 0.5F, Hg, Tol Imm	32, 300, 340, 370, 400F	\$668.03	UN3506
ACC069CFC	ASTM 69C, 195/305C x 0.5C, Hg, Tol Imm	0, 200, 235, 270, 305C	\$668.03	UN3506
ACC069FFC	ASTM 69F, 383/581F x 1.0F, Hg, Tol Imm	32, 400, 460, 520, 580F	\$668.03	UN3506
ACC070CFC	ASTM 70C, 295/405C x 0.5C, Hg, Tol Imm	0, 300, 335, 370, 400C	\$781.20	UN3506
ACC070FFC	ASTM 70F, 563/761F x 1.0F, Tol Imm	32, 570, 640, 700, 760F	\$781.20	UN3506
ACC071CFC	ASTM 71C, -37/21C x 0.5C, Hg, 76mm Imm	-35, -18, 0, 20C	\$354.47	UN3506
ACC071FFC	ASTM 71F, -35/70 x 1.0F, Hg, 76mm Imm	-30, 0, 32, 70F	\$354.47	UN3506
ACC072CFC	ASTM 72C, -19.5/16.6C x 0.05C Hg, Tol Imm	-19, 17.8, 0C	\$414.42	UN3506
ACC072FFC	ASTM 72F, -2.5/2.5F x 0.1F, Hg, Tol Imm	-2, 0, 32F	\$414.42	UN3506
ACC075FFC	ASTM 75F, -35/35F x 0.5F, Hg, 100mm Imm	-35, 0, 32F	\$294.50	UN3506
ACC077FFC	ASTM 77F, 245/265 x 0.5F, Hg, Tol Imm	250, 260F	\$212.70	UN3506
ACC078FFC	ASTM 78F, 295/315F x 0.5F, Hg, Tol Imm	300, 310F	\$212.70	UN3506
ACC079FFC	ASTM 79F, 345/365F x 0.5F, Hg, Tol Imm	350, 360F	\$212.70	UN3506
ACC080FFC	ASTM 80F, 395/415F x 0.5F, Hg, Tol Imm	400, 410F	\$239.97	UN3506
ACC081FFC	ASTM 81F, 445/465F x 0.5F, Hg, Tol Imm	450, 460F	\$239.97	UN3506
ACC082CFC	ASTM 82C, -15/105C x 1.0C, Hg, 30mm Imm	0, 50, 100C	\$256.25	UN3506
ACC082FFC	ASTM 82F, 0/220F x 2.0F, Hg, 30mm Imm	32, 100, 200F	\$256.25	UN3506
ACC083CFC	ASTM 83C, 15/70C x 0.5C, Hg, 40mm Imm	25, 70C	\$190.87	UN3506
ACC083FFC	ASTM 83F, 60/160F x 1.0F, Hg, 40mm Imm	85, 135F	\$190.87	UN3506
ACC084CFC	ASTM 84C, 25/80C x 1.0C, Hg, 249mm Imm	30, 80C	\$190.87	UN3506
ACC084FFC	ASTM 84F, 75/175F x 1.0F, Hg, 249mm Imm	100, 150F	\$190.87	UN3506
ACC085CFC	ASTM 85C, 40/150C x 1.0C, Hg, 181mm Imm	50, 150C	\$190.87	UN3506
ACC085FFC	ASTM 85F, 100/300F x 2.0F, Hg, 181mm Imm	150, 250F	\$190.87	UN3506

<u>CAT. NO.</u>	<u>DESCRIPTION</u>		<u>LIST PRICE</u>	<u>Haz (Hg)</u>	<u>COO</u>
ACC086CFC	ASTM 86C, 95/175C x 1.0C, Hg, 35mm Imm	100, 175C	\$199.08	UN3506	
ACC086FFC	ASTM 86F, 200/350F x 1.0F, Hg, 35mm Imm	225, 325F	\$199.08	UN3506	
ACC087CFC	ASTM 87C, 150/205C x 1.0C, Hg, 40mm Imm	160, 200C	\$207.28	UN3506	
ACC087FFC	ASTM 87F, 300/400F x 1.0F, Hg, 40mm Imm	300, 400F	\$207.28	UN3506	
ACC088CFC	ASTM 88C, 10/200C x 1.0C, Hg, 57mm Imm	40, 100, 150, 200C	\$354.47	UN3506	
ACC088FFC	ASTM 88F, 50/392F x 2.0F, Hg, 57mm Imm	110, 212, 300, 392F	\$354.47	UN3506	
ACC089CFC	ASTM 89C, -20/10C x 0.1C, Hg, 76 Imm	-20, -10, 0, 10C	\$378.95	UN3506	
ACC090CFC	ASTM 90C, 0/30C x 0.1C, Hg, 76mm Imm	0, 10, 20, 30C	\$343.62	UN3506	
ACC091CFC	ASTM 91C, 20/50C x 0.1C, Hg, 76mm Imm	20, 30, 40, 50C	\$343.62	UN3506	
ACC092CFC	ASTM 92C, 40/70C x 0.1C Hg, 76mm Imm	40, 50, 60, 70C	\$343.62	UN3506	
ACC093CFC	ASTM 93C, 60/90C x 0.1C, Hg, 76mm Imm	60, 70, 80, 90C	\$354.47	UN3506	
ACC094CFC	ASTM 94C, 80/110C x 0.1C, Hg, 76mm Imm	80, 90, 100, 110C	\$354.47	UN3506	
ACC095CFC	ASTM 95C, 100/130C x 0.1C, Hg, 76mm Imm	100, 110, 120, 130C	\$354.47	UN3506	
ACC096CFC	ASTM 96C, 120/150C x 0.1C, Hg, 76mm Imm	120, 130, 140, 150C	\$354.47	UN3506	
ACC097CFC	ASTM 97C, -18/49C x 0.5C, Hg, Tol Imm	-15, 0, 20, 45C	\$332.62	UN3506	
ACC097FFC	ASTM 97F, 0/120F x 1.0F, Hg, Hook Top	0, 32, 70, 110F	\$332.62	UN3506	
ACC098CFC	ASTM 98C, 16/82C x 0.5C, Hg, Tol Imm	20, 40, 60, 80C	\$332.62	UN3506	
ACC098FFC	ASTM 98F, 60/180F x 1.0F, Hg, Tol Imm	60, 100, 140, 180F	\$332.62	UN3506	
ACC100CFC	ASTM 100C, 145/205C x 0.2C, Hg, 76mm Imm	145, 165, 185, 205C	\$384.52	UN3506	
ACC101CFC	ASTM 101C, 195/305C x 0.2C, Hg, 76mm Imm	200, 250, 300C	\$346.26	UN3506	
ACC102CFC	ASTM 102C, 123/177C x 0.2, Hg, 100mm Imm	125, 140, 155, 175C	\$354.47	UN3506	
ACC103CFC	ASTM 103C, 148/202C x 0.2C, Hg, 100mm Imm	150, 165, 180, 200C	\$359.90	UN3506	
ACC104CFC	ASTM 104C, 173/227C x 0.2C, Hg, 100mm Imm	175, 190, 205, 225C	\$359.90	UN3506	
ACC105CFC	ASTM 105C, 198/252F x 1.0F, Hg, 40mm Imm	200, 215, 230, 250C	\$545.34	UN3506	
ACC106CFC	ASTM 106C, 223/277C x 0.2C, Hg, 100mm Imm	225, 240, 255, 275C	\$545.34	UN3506	
ACC107CFC	ASTM 107C, 248/302C x 0.2C, Hg, 100mm Imm	250, 265, 280, 300C	\$545.34	UN3506	
ACC108FFC	ASTM 108F, 270/290F x 0.5F, Hg, Tol Imm	DISCONTINUED			
ACC109FFC	ASTM 109F, 320/340F x 0.5F, Hg, Tol Imm	DISCONTINUED			
ACC110FFC	ASTM 110F, 272.5/277.5F x 0.1F, Hg, Tol Imm	32, 275, 277F	\$403.58	UN3506	
ACC111CFC	ASTM 111C, 170/250C x 0.2C Hg, 100mm Imm	170, 200, 250C	\$327.21	UN3506	
ACC112CFC	ASTM 112C, 4/6C x 0.02C, Hg, Tol Imm	DISCONTINUED			
ACC113CFC	ASTM 113C, -1/175C x 0.5C, Hg, Tol Imm	0, 50, 100, 150, 175C	\$409.00	UN3506	
ACC113FFC	ASTM 113F, 30/350F x 1.0F,Hg, Tol Imm	32, 122, 212, 302, 347F	\$409.00	UN3506	
ACC114CFC	ASTM 114C, -80/20C x 0.5C, Non-Hg, Tol Imm	-75, -60, -40, 0C	\$377.62		
ACC114FFC	ASTM 114F, -112/70F x 1.0F, Non-Hg, Tol Imm	-103, -76, -40, 32F	\$378.95		
ACC116CFC	ASTM 116C, 18.9/25.1C x 0.01C, Hg, Tol Imm	19, 20, 21, 22, 23, 24, 25C	\$796.15	UN3506	
ACC117CFC	ASTM 117C, 23.9/30.1C x 0.01C, Hg, Tol Imm	24, 25, 26, 27, 28, 29, 30C	\$796.15	UN3506	
ACC118CFC	ASTM 118C, 28.6/31.4 x 0.05C, Hg, Tol Imm	0, 30, 31C	\$389.94	UN3506	
ACC118FFC	ASTM 118F, 83.5/88.5 x 0.1F, Hg, Tol Imm	32, 86, 88F	\$389.94	UN3506	
ACC119CFC	ASTM 119C, -38.3/-30C x 0.1C, Hg, 100mm Imm	DISCONTINUED			
ACC119FFC	ASTM 119F, -37/-22F x 0.2F, Hg, 100mm Imm	-36, -22, 32F	\$463.53	UN3506	
ACC120CFC	ASTM 120C, 38.6/41.4C x 0.05C, Hg, Tol Imm	0, 40, 41C	\$387.16	UN3506	
ACC121CFC	ASTM 121C 98.6/101.6C x 0.05C, Hg, Tol Imm	0, 100, 101C	\$387.16	UN3506	
ACC123CFC	ASTM 123C, -35/-25C x 0.1C, Hg, Tol Imm	-35, -30, -25C	\$403.58	UN3506	
ACC124CFC	ASTM 124C -25/-15C x 0.1C, Hg, Tol Imm	-25, -20, -15C	\$403.58	UN3506	
ACC125CFC	ASTM 125C -15/-5C x 0.1C, Hg, Tol Imm	-15, -10, -5C	\$403.58	UN3506	
ACC126CFC	ASTM 126C, -27.4/-24.6C x 0.05C, Hg, Tol Imm	-27, -26.1, 0C	\$460.73	UN3506	
ACC126FFC	ASTM 126F, -17.5/-12.5 x 0.1F, Hg, Tol Imm	-17, -15, 32F	\$460.73	UN3506	
ACC127CFC	ASTM 127C, -21.4/-18.6C x 0.05C, Hg, Tol Imm	-21, -20, 0C	\$460.73	UN3506	
ACC128CFC	ASTM 128C, -1.4/1.4C x 0.05C, Hg, Tol Imm	0, 1C	\$335.41	UN3506	
ACC128FFC	ASTM 128F, 29.5/34.5 x 0.1F, Hg, Tol Imm	32, 34F	\$335.41	UN3506	
ACC129CFC	ASTM 129C, 91.6/94.4C x 0.05C, Hg, Tol Imm	0, 93.3, 94C	\$384.52	UN3506	
ACC129FFC	ASTM 129F, 197.5/202.5F x 0.1F, Hg, Tol Imm	32, 200, 202F	\$384.52	UN3506	
ACC130CFC	ASTM 130C, -7/105C x 0.5C, Hg, Tol Imm	0, 35, 70, 105C	\$313.58	UN3506	
ACC130FFC	ASTM 130F 20/220F x 1.0F, Hg, Hook Top	32, 100, 160, 220F	\$313.58	UN3506	

<u>CAT. NO.</u>	<u>DESCRIPTION</u>	<u>LIST PRICE</u>	<u>Haz (Hg)</u>	<u>COO</u>
ACC132CFC	ASTM 132C, 148.6/151.4 x 0.05C, Hg, Tol Imm 0, 150, 151C	\$381.73	UN3506	
ACC133CFC	ASTM 133C, -38/2C x 0.1C, Hg, 76mm Imm -24, -18, -12, -6, 0C	\$602.65	UN3506	
ACC134CFC	ASTM 134C 144/156C x 0.2C, Hg, 100mm imm 145, 150, 155C	\$291.72	UN3506	
ACC135CFC	ASTM 135C, 38/93C x 1.0C, Hg, 40mm Imm 50, 90C	\$207.28	UN3506	
ACC135FFC	ASTM 135F, 100/200F x 1.0F Hg, 40mm Imm 125, 195F	\$207.28	UN3506	
ACC136CFC	ASTM 136C, -20/60C x 0.2C, Hg, Tol Imm -20, -10, 0, 10, 20, 30, 40, 50, 60C	\$708.92	UN3506	
ACC136FFC	ASTM 136F, -5/140F x 0.5F, Hg, Tol Imm -5, 15, 32, 60, 85, 110, 135F	\$564.37	UN3506	
ACC137CFC	ASTM 137C 80/100C x 0.1C, Hg, 76mm Imm 80, 90, 100C	\$332.62	UN3506	

ASTM CERTIFIED THERMOMETERS, SPIRIT FILLED

<u>CAT. NO.</u>	<u>DESCRIPTION</u>	<u>LIST PRICE</u>
ACC001CBLSF-20/150C x 1.0C	<u>CERTIFIED TEMPERATURES</u> -20, 0, 50, 100, 150C	\$493.00
ACC001FBLSF(0/302F x 2.0F	0, 32, 122, 212, 302F	\$493.00
ACC005CBLSF-38/50C x 1.0C	-35, 0, 50C	\$330.23
ACC005FBLSF(-36/+120F x 2.0F	-30, 32, 120F	\$330.23
ACC009CSFC -5/110C x 0.5C	0, 35, 70, 105C	\$330.23
ACC009FSFC +20/230F x 1.0F	32, 100, 160, 220F	\$330.23
ACC012CBLSF-20/102C x 0.2C	-20, -10, 0, 10, 20, 30, 40, 50, 60, 70C	\$1,172.21
ACC012FBLSF(-5/215F x 0.5F	-5, 15, 32, 60, 85, 110, 135, 160, 180F	\$1,172.21
ACC013CBLSF+155/170C x 0.5	155, 163, 170C	\$336.50
ACC014CBLSF+38/82C x 0.1C	40, 50, 60, 70, 80C	\$508.65
ACC014FBLSF(+100/180F x 0.2F	100, 120, 140, 160, 180F	\$508.65
ACC015CBLSF-2/80C x 0.2C	0, 20, 40, 60, 80C	\$530.57
ACC015FBLSF(+30/180F x 0.5F	32, 70, 100, 140, 180F	\$530.57
ACC016CBLSF+30/200C x 0.5C	30, 60, 90, 120, 150, 180, 200C	\$672.96
ACC016FBLSF(+85/392F 1.0F	90, 140, 190, 240, 290, 340, 390F	\$672.96
ACC017CBLSX+19/27C x 0.1C	21, 25C	\$281.70
ACC017FBLSX+66/80F x 0.2F	70, 77F	\$281.70
ACC018CBLSF+34/42C x 0.1C	38, 41C	\$281.70
ACC018FBLSF(+94/108F x 0.2F	100, 107F	\$281.70
ACC019CBLSF+49/57C x 0.1C	50, 54C	\$281.70
ACC019FBLSF(+120//134F x 0.2F	122, 130F	\$281.70
ACC020CBLSF+57/65C x 0.1C	60, 64C	\$281.70
ACC020FBLSF(+134/148F x 0.2F	140, 147F	\$281.70
ACC021CBLSF+79/87C x 0.1C	82, 86C	\$281.70
ACC021FBLSF(+174/188 x 0.2F	180, 187F	\$281.70
ACC022CBLSF+95/103C x 0.1C	99, 102C	\$281.70
ACC022FBLSF(+204/218F x 0.2F	210, 212F	\$281.70
ACC024CBLSF+39/54C x 0.2C	40, 50C	\$254.00
ACC025CBLSF+95/105C x 0.2C	95, 100C	\$254.00
ACC026CBLSF+130/140C x 0.1C	130, 135, 140C	\$356.14
ACC027CBLSF+147/182C x 0.5C	155, 165, 175C	\$356.14
ACC028CBLSF+36.6/39.4 x 0.05C	0, 37.8, 39C	\$439.92
ACC028FBLSF(+97.5/102.5F x 0.1F	32, 100, 102F	\$439.92
ACC029CBLSF+52.6/55.4C x 0.05C	0, 54.4, 55C	\$439.92
ACC029FBLSF(+127.5/132.5F x 0.1F	32, 130, 132F	\$439.92
ACC030FBLSF(+207.5/212.5F x 0.1F	32, 210, 212F	\$439.92
ACC033CBLSF-38/42C x 0.2C	-35, -20, 0, 20, 40C	\$572.80
ACC033FBLSF(-36.5/107.5F x 0.5F	-31, -4, 32, 68, 104F	\$572.80
ACC034CBLSF+25/105C x 0.2C	25, 45, 65, 85, 100C	\$585.32
ACC034FBLSF(+77/221F x 0.5F	77, 113, 149, 185, 212F	\$585.32
ACC035CBLSF+90/170C x 0.2C	100, 120, 140, 160, 170C	\$490.99
ACC035FBLSF(+194/338F x 0.5F	212, 250, 285, 320, 338F	\$490.99
ACC036CBLSF-2/68C x 0.2C	0, 15, 30, 45, 65C	\$535.24
ACC037CBLSF-2/52C x 0.2C	0, 15, 30, 50C	\$435.08
ACC038CBLSF+24/78C x 0.2C	25, 40, 55, 75F	\$453.86

<u>CAT. NO.</u>	<u>DESCRIPTION</u>	<u>LIST PRICE</u>	<u>Haz (Hg)</u>	<u>COO</u>
ACC039CBLSF+48/102C x 0.2C	50, 65, 80, 100C	\$453.86		
ACC040CBLSF+72/126C x 0.2C	75, 90, 105, 125C	\$449.18		
ACC042CBLSF+95/255C x 0.5C	100, 150, 200, 250C	\$488.29		
ACC044CBLSF+18.6/21.4C x 0.05C	0, 20, 21C	\$437.30		
ACC044FBLSF(+66.5/71.5F x 0.1F	32, 68, 70F	\$437.30		
ACC045CBLSF+23.6/26.4C x 0.05C	0, 25, 26C	\$437.30		
ACC045FBLSF(+74.5/79.5F x 0.1F	32, 77, 79F	\$437.30		
ACC046CBLSF+48.6/51.4C x 0.05C	0, 50, 51C	\$437.30		
ACC046FBLSF(+119.5/124.5F x 0.1F	32, 122, 124F	\$437.30		
ACC047CBLSF+58.6/61.4C x 0.05C	0, 60, 61C	\$437.30		
ACC047FBLSF(+137.5/142.5F x 0.1F	32, 140, 142F	\$437.30		
ACC048CBLSF+80.6/83.4C x 0.05C	0, 82.2, 83C	\$437.30		
ACC048FBLSF(+177.5/182.5F x 0.1F	32, 180, 182F	\$437.30		
ACC049CBLSF+20/70C x 0.2C	20, 35, 50, 70C	\$517.42		
ACC052CBLSF-10/5C x 0.1C	-10, 0, 5C	\$666.88		
ACC054CBLSF+20/100.6C x 0.2C	20, 50, 75, 100C	\$452.30		
ACC054FBLSF(+68/213F x 0.5F	70, 120, 170, 210F	\$452.30		
ACC056CBLSF+19/35C x 0.02C	19, 21, 23, 25, 27, 29, 31, 33, 35C	\$1,139.97		
ACC056FBLSF(+66/95F x 0.05F	66, 70, 74, 78, 82, 88, 92, 95F	\$1,139.97		
ACC057CBLSF-20/50C x 0.5C	-20, 0, 25, 50C	\$1,126.81		
ACC057FBLSF(-4/122F x 1.0F	-3, 32, 77, 122F	\$1,126.81		
ACC058CBLSF-34/+49C x 0.5C	-30, 0, 25, 45C	\$432.24		
ACC058FBLSF(-30/+120F x 1.0F	-20, 32, 80, 120F	\$432.24		
ACC059CSFC -18/+82C x 0.5C	0, 25, 55, 80C	\$432.24		
ACC059FSFC 0/180F x 1.0F	32, 80, 130, 180F	\$432.24		
ACC061CBLSF+32/127C x 0.2C	40, 60, 80, 100, 120C	\$535.24		
ACC061FBLSF(+90/260F x 0.5F	100, 150, 200, 250F	\$464.82		
ACC071CBLSF-37/21XC x 0.5C	-35, -18, 0, 20C	\$400.64		
ACC071FBLSF(-35/70F x 1.0F	-30, 0, 32, 70F	\$400.64		
ACC072CBLSF-19.4/16.6C x 0.05C	-19, 17.8, 0C	\$487.06		
ACC072FBLSF(-2.5/2.5F x 0.1F	-2, 0, 32F	\$487.06		
ACC073CBLSF-41.4/-38.6C X 0.05C	-41, -40, 0C	\$386.16		
ACC073FBLSF(-42.5/-37.5F X 0.1F	-42, -40, 32F	\$386.16		
ACC075FBLSF(-35/35F x 0.5F	-35, 0, 32F	\$345.66		
ACC078FBLSF(+245/265F x 0.5F	300, 310F	\$257.94		
ACC079FBLSF(+295/315F x 0.5F	350, 360F	\$257.94		
ACC083CBLSF+15/70C x 1.0C	25, 70C	\$338.05		
ACC083FBLSF(+60/160F x 1.0F	85, 135F	\$338.05		
ACC084CBLSF+25/80C x 1.0C	30, 80C	\$338.05		
ACC084FBLSF(+75/175F x 1.0F	100, 150F	\$338.05		
ACC085CBLSF+40/150C x 1.0C	50, 150C	\$338.05		
ACC085FBLSF(+100/300F x 2.0F	150, 250F	\$338.05		
ACC088CBLSF+10/200C x 1.0C	40, 100, 150, 200C	\$466.38		
ACC088FBLSF(+50/392F x 2.0F	110, 212, 300, 392F	\$466.38		
ACC089CBLSF-20/10C x 0.1C	-20, -10, 0, 10C	\$494.56		
ACC090CBLSF0/30C x 0.1C	0, 10, 20, 30C	\$444.46		
ACC091CBLSF+20/50C x 0.1C	20, 30, 40, 50C	\$437.22		
ACC092CBLSF+40/70C x 0.1C	40, 50, 60, 70C	\$371.83		
ACC093CBLSF+60/90C x 0.1C	60, 70, 80, 90C	\$400.64		
ACC094CBLSF+80/110C x 0.1C	80,90, 100, 110C	\$406.22		
ACC095CBLSF+100/130C x 0.1C	100, 110, 120, 130C	\$400.64		
ACC096CBLSF120/150C x 0.1C	120, 130, 140, 150C	\$400.64		
ACC097CBLSF-18/49C x 0.5C	-15, 0, 20, 45C	\$428.81		
ACC097FBLSF(0/120F x 1.0F	0, 32, 70, 110F	\$428.81		
ACC098CBLSF+16/82C x 0.5C	20, 40, 60, 80C	\$428.81		

<u>CAT. NO.</u>	<u>DESCRIPTION</u>	<u>LIST PRICE</u>	<u>Haz (Hg)</u>	<u>COO</u>
ACC098FBLSF(+60/180F x 1.0F	60, 100, 140, 180F	\$428.81		
ACC099CBLSF-50/5C x 0.2C	*-46, -32, -18, 0C	\$500.18		
ACC099FBLSF(-58/41F x 0.5F	*-50, -25, 0, 32F	\$500.18		
ACC102CBLSF+123/177C x 0.2C	125, 140, 155, 175C	\$401.96		
ACC103CBLSF+148/202C x 0.2C	150, 165, 180, 200C	\$411.12		
ACC110CBLSF+133.6/136.4C x 0.05C	0, 135, 138C	\$484.42		
ACC110FBLSF(+272.5/277.5F x 0.1F	32, 275, 277F	\$484.42		
ACC112CBLSF+4/6C x 0.02C	0, 4, 5, 6C	\$612.74		
ACC113CBLSF-1/175C x 0.5C	0, 50, 100, 150, 175C	\$524.30		
ACC113FBLSF(+30/350F x 1.0F	32, 122, 212, 302, 347F	\$524.30		
ACC116CBLSF+18.9/ 25.1C x 0.01C	19, 20, 21, 22, 23, 24, 25C	\$1,156.34		
ACC117CBLSF+23.9/30.1C x 0.01C	24, 25, 26, 27, 28, 29, 30C	\$1,156.34		
ACC118CBLSF+28.6/31.4C x 0.05C	0, 30, 31C	\$453.01		
ACC119CBLSF-38.3/-30.0C x 0.05C	-38, -30, 0C	\$536.80		
ACC120CBLSF+38.6/41.4C x 0.05	0, 40, 41C	\$513.34		
ACC121CBLSF+98.6/101.4C x 0.05C	0, 100, 101C	\$429.44		
ACC122CBLSF-45/-35C x 0.1C	-45, -40, -35C	\$438.46		
ACC123CBLSF-35/-25C x 0.1C	-35, -30, -25C	\$468.71		
ACC124CBLSF-25/-15C x 0.1C	-25, -20, -15C	\$468.71		
ACC125CBLSF-15/-5C x 0.1C	-15, -10, -5C	\$468.71		
ACC126CBLSF-27.4/-24.6C x 0.05C	-27, -26.1, 0C	\$471.34		
ACC126FBLSF(-17.5/-12.5F x 0.1F	-17, -15, 32F	\$471.34		
ACC127CBLSF-21.4/-18.6C x 0.05C	-21, -20, 0C	\$471.34		
ACC128CBLSF-1.4/+1.4C x 0.05C	0, 1C	\$387.54		
ACC128FBLSF(+29.5/34.5F x 0.1F	32, 34F	\$387.54		
ACC129CBLSF+91.6/94.4F x 0.05C	0, 93.3, 94C	\$434.68		
ACC129FBLSF(+197.5/202.5F x 0.1F	32, 200, 202F	\$434.68		
ACC130CBLSF-7/105C x 0.5C	0, 35, 70, 105C	\$416.31		
ACC130FBLSF(20/220F x 1.0F	32, 100, 160, 220F	\$416.31		
ACC132CBLSF+148.6/151.4C x 0.05C	0, 150, 151C	\$453.01		
ACC133CBLSF-38/+2C x 0.1C	-24, -18, -12, -6, 0C	\$670.34		
ACC134CBLSF144/156C x 0.2C	145, 150, 155C	\$453.01		
ACC135CBLSF+39/93C x 1.0C	50, 90C	\$235.68		
ACC135FBLSF(+100/200F x 1.0F	125, 195F	\$235.68		
ACC136CBLSF-20/+60C x 0.2C	-20, -10, 0, 10, 20, 30, 40, 50, 60C	\$792.11		
ACC136FBLSF(-5/140F x 0.5F	-5, 15, 32, 60, 85, 110, 135F	\$646.77		
ACC137CBLSF+80/100C x 0.1C	80, 90, 100C	\$379.69		

<u>Catalog No.</u>	<u>Range</u>	<u>Div</u>	<u>Filling</u>	<u>Imm</u>	<u>L</u>	<u>LIST PRICE</u>	<u>COO</u>	<u>Haz (Hg)</u>
<u>Partial Immersion, Spirit (Red & Blue)</u>							HTS# 9025.11. 4000	
B50CW3S	-35 to 50C	1.0C	RS	76mm	12"	\$16.21	DE, GB	
B110CW3S u	-20 to 110C	1.0C	RS	76mm	12"	\$14.93	DE, GB	
B150CW3S u	-20 to 150C	1.0C	RS	76mm	12"	\$16.21	DE, GB	
B200CW3S u	-10 to 200C	1.0C	RS	76mm	12"	\$18.92	DE, GB	
B260CW3BLS	-10 to 260C	1.0C	BS	76mm	12"	\$21.62	GB	
B330CW3BRN	-10 to 330C	2.0C	RS	76mm	12"	\$32.40	GB	
B120FW3S	-30 to 120F	2.0F	RS	76mm	12"	\$16.21	DE, GB	
B230FW3S	0 to 230F	2.0F	RS	76mm	12"	\$14.93	DE, GB	
B300FW3S	0 to 300F	2.0F	RS	76mm	12"	\$16.21	DE, GB	
B110C230FW3S	-20 to 110C	1.0C	RS	76mm	12"	\$21.62	GB	
	& 0 to 230F & 2.0F							
B150C300FW3S	-20 to 150C	1.0C	RS	76mm	12"	\$23.03	GB	
	& 0 to 300F & 2.0F							
ACC260C500FW3BS	-10 to 260C	1.0C	BS	76mm	12"	\$27.02	GB	
	& 20 to 500F & 2.0F							
<u>Tol Immersion, Spirit</u>								
B50CWBL S	-35 to 50C	1.0C	BS	Tol	12"	\$18.27	GB	
B110CWS	-20 to 110C	1.0C	RS	Tol	12"	\$14.93	DE, GB	
B150CWS	-20 to 150C	1.0C	RS	Tol	12"	\$16.21	DE, GB	
B200CWBL S	-10 to 200C	1.0C	BS	Tol	12"	\$18.92	DE, GB	
B120FWS	-30 to 120F	1.0F	RS	Tol	12"	\$16.21	DE, GB	
B230FWS	0 to 230F	2.0F	RS	Tol	12"	\$14.93	DE, GB	
B300FWS	0 to 300F	2.0F	RS	Tol	12"	\$16.21	DE, GB	
<u>Low Temperature, Partial and Tol Immersion, Spirit</u>								
B5050CW3S	-50 to 50C	1.0C	RS	76mm	12"	\$27.02	DE, GB	
B5050CWS	-50 to 50C	1.0C	RS	Tol	12"	\$27.02	DE, GB	
<u>Fluoro-Polymer Coated Thermometers, Spirit-Filled</u>								
<u>Partial Immersion, Red Spirit</u>								
B50CW3SSC	-35 to 50C	1.0C	RS	76mm	12"	\$24.30	DE, GB	
B110CW3SSC	-20 to 110C	1.0C	RS	76mm	12"	\$24.30	DE, GB	
B150CW3SSC	-20 to 150C	1.0C	RS	76mm	12"	\$24.30	DE, GB	
B230FW3SSC	0 to 230F	2.0F	RS	76mm	12"	\$24.30	DE, GB	
B300FW3SSC	0 to 300F	2.0F	RS	76mm	12"	\$24.30	DE, GB	
B110C230FW3SSC	-20 to 110C	1.0C	RS	76mm	12"	\$29.72	GB	
	&							
	& 0 to 230F	2.0F						
B150C300FW3SSC	-20 to 150C	1.0C	RS	76mm	12"	\$29.72	GB	
	&							
	& 0 to 300F	2.0F						
ACC260C500FW3BSSC	-10 to 260C	1.0C	BS	76mm	12"	\$37.82	GB	
	&							
	& 20 to 500F	2.0F						

<u>Catalog No.</u>	<u>Range</u>	<u>Div</u>	<u>Filling</u>	<u>Imm</u>	<u>L</u>	<u>LIST PRICE</u>	<u>COO</u>	<u>Haz (Hg)</u>
<u>Tol Immersion, Spirit, Safety Coated</u>								
B50CWBLSSC	-35 to 50C	1.0C	BS	Tol	12"	\$31.97	GB	
B110CWSSC	-20 to 110C	1.0C	RS	Tol	12"	\$24.30	DE, GB	
B150CWSSC	-20 to 150C	1.0C	RS	Tol	12"	\$24.30	DE, GB	
B200CWSSC	-10 to 200C	1.0C	RS	Tol	12"	\$29.72	DE, GB	
B230FWSSC	0 to 230F	2.0F	RS	Tol	12"	\$24.30	DE, GB	
B300FWSSC	0 to 300F	2.0F	RS	Tol	12"	\$24.30	DE, GB	
<u>Partial Immersion, Length 8" (200mm), Spirit, Safety Coated</u>								
B8TC1102BLS	-10 to 110C	1.0C	BS	50mm	8"	\$25.86	GB	
B8TC1502BLS	-10 to 150C	1.0C	BS	50mm	8"	\$25.86	GB	
B8TC2302BLS	0 to 230F	2.0F	BS	50mm	8"	\$21.62	GB	
B8TC3002BLS	0 to 300F	2.0F	BS	50mm	8"	\$21.62	GB	
<u>Tol Immersion, Length 8" (200mm), Spirit, Safety Coated</u>								
B8TC110BLS	-10 to 110C	1.0C	BS	Tol	8"	\$25.86	GB	
B8TC150BLS u	-10 to 150C	1.0C	BS	Tol	8"	\$25.86	GB	
B8TC230BLS	0 to 230F	2.0F	BS	Tol	8"	\$21.62	GB	
<u>Low Temperature, Partial and Tol Immersion, Spirit, Safety Coated</u>								
B5050CW3SSC	-50 to 50C	1.0C	RS	76mm	12"	\$37.82	DE, GB	
B5050CWSSC	-50 to 50C	1.0C	RS	Tol	12"	\$37.82	DE, GB	
<u>Mercury-Filled Thermometers</u>								
<u>Partial Immersion, Mercury Filled</u>								
B50CY3	-35 to 50C	1.0C	Hg	76mm	12"	\$21.62	DE, GB	UN3506
B110CY3	-20 to 110C	1.0C	Hg	76mm	12"	\$21.62	DE, GB	UN3506
B150CY3	-20 to 150C	1.0C	Hg	76mm	12"	\$21.62	DE, GB	UN3506
B260CY3	-10 to 260C	1.0C	Hg	76mm	16"	\$27.02	DE, GB	UN3506
B400CY3	-10 to 420C	1.0C	Hg	76mm	13"	\$29.72	DE, GB	UN3506
B120FY3	-30 to 120F	1.0F	Hg	76mm	12"	\$21.62	DE, GB	UN3506
B230FY3	0 to 230F	2.0F	Hg	76mm	12"	\$21.62	DE, GB	UN3506
B300FY3	0 to 300F	2.0F	Hg	76mm	12"	\$21.62	DE, GB	UN3506
B500FY3	20 to 500F	2.0F	Hg	Tol	16"	\$27.02	DE, GB	UN3506
B110C230FW3	-20 to 110C	1.0C	Hg	76mm	12"	\$27.02	GB	UN3506
		&						
	& 0 to 230F	2.0F						
B110C230FY3	-20 to 110C	1.0C	Hg	76mm	12"	\$27.02	DE, GB	UN3506
		&						
	& 0 to 230F	2.0F						
B150C300FW3	-20 to 150C	1.0C	Hg	76mm	12"	\$27.02	GB	UN3506
		&						
	& 0 to 300F	2.0F						
B260C500FW3	-10 to 260C	1.0C	Hg	76mm	12"	\$35.11	GB	UN3506
		&						
	& 20 to 500F	2.0F						
B400C750FY3	-10 to 400C	2.0C	Hg	76mm	12"	\$43.21	GB	UN3506
		&						
	& 20 to 750F	5.0F						

<u>Catalog No.</u>	<u>Range</u>	<u>Div</u>	<u>Filling</u>	<u>Imm</u>	<u>L</u>	<u>LIST PRICE</u>	<u>COO</u>	<u>Haz (Hg)</u>
<u>Tol Immersion, Mercury Filled</u>								
B50CY	-35 to 50C	1.0C	Hg	Tol	12"	\$21.62	DE, GB	UN3506
B110CY	-20 to 110C	1.0C	Hg	Tol	12"	\$21.62	DE, GB	UN3506
B150CY	-20 to 150C	1.0C	Hg	Tol	12"	\$21.62	DE, GB	UN3506
B260CY	-10 to 200C	1.0C	Hg	Tol	16"	\$27.02	DE, GB	UN3506
B400CY	-10 to 400C	2.0C	Hg	Tol	16"	\$29.72	DE, GB	UN3506
B120FY	-30 to 120F	1.0F	Hg	Tol	12"	\$21.62	DE, GB	UN3506
B230FY	0 to 230F	2.0F	Hg	Tol	12"	\$21.62	DE, GB	UN3506
B300FY	0 to 300F	2.0F	Hg	Tol	12"	\$21.62	DE, GB	UN3506
B500FY	20 to 500F	2.0F	Hg	Tol	16"	\$27.02	DE, GB	UN3506
B750FY	20 to 750F	5.0F	Hg	Tol	16"	\$29.72	DE, GB	UN3506
B110C230FY	-20 to 110C	1.0C	Hg	76mm	12"	\$27.02	DE, GB	UN3506
		&						
	& 0 to 230F	2.0F						
B150C300FY	-20 to 150C	1.0C	Hg	76mm	12"	\$27.02	DE, GB	UN3506
		&						
	& 0 to 300F	2.0F						
B400C750FY	-10 to 400C	2.0C	Hg	76mm	12"	\$43.21	GB	UN3506
		&						
	& 20 to 750F	5.0F						
<u>Fluoro-Polymer Coated Thermometers, Mercury-Filled</u>								
<u>Partial Immersion, Mercury Filled, Coated</u>								
B50CY3SC	-35 to 50C	1.0C	Hg	76mm	12"	\$32.40	DE, GB	UN3506
B110CY3SC	-20 to 110C	1.0C	Hg	76mm	12"	\$32.40	DE, GB	UN3506
B150CY3SC	-20 to 150C	1.0C	Hg	76mm	12"	\$32.40	DE, GB	UN3506
B260CY3SC	-10 to 260C	1.0C	Hg	76mm	16"	\$37.82	DE, GB	UN3506
B120FY3SC	-30 to 120F	1.0F	Hg	76mm	12"	\$32.40	DE, GB	UN3506
B230FY3SC	0 to 230F	2.0F	Hg	76mm	12"	\$32.40	DE, GB	UN3506
B300FY3SC	0 to 300F	2.0F	Hg	76mm	12"	\$32.40	DE, GB	UN3506
B500FY3SC	20 to 500F	2.0F	Hg	Tol	16"	\$37.82	DE, GB	UN3506
B110C230FW3SC	-20 to 110C	1.0C	Hg	76mm	12"	\$27.02	GB	UN3506
		&						
	& 0 to 230F	2.0F						
B110C230FY3SC	-20 to 110C	1.0C		76mm	12"	\$27.02	DE, GB	UN3506
		&						
	& 0 to 230F	2.0F	Hg					
B150C300FW3SC	-20 to 150C	1.0C	Hg	76mm	12"	\$27.02	GB	UN3506
		&						
	& 0 to 300F	2.0F						
B150C300FY3SC	-20 to 150C	1.0C	Hg	76mm	12"	\$27.02	DE, GB	UN3506
		&						
	& 0 to 300F	2.0F						
B260C500FW3SC	-10 to 260C	1.0C	Hg	76mm	12"	\$35.11	GB	UN3506
		&						
	& 20 to 500F	2.0F						
<u>Tol Immersion, Mercury Filled, Coated</u>								
B50CYSC	-35 to 50C	1.0C	Hg	Tol	12"	\$32.40	DE, GB	UN3506
B110CYSC	-20 to 110C	1.0C	Hg	Tol	12"	\$32.40	DE, GB	UN3506
B150CYSC	-20 to 150C	1.0C	Hg	Tol	12"	\$32.40	DE, GB	UN3506

<u>Catalog No.</u>	<u>Range</u>	<u>Div</u>	<u>Filling</u>	<u>Imm</u>	<u>L</u>	<u>LIST PRICE</u>	<u>COO</u>	<u>Haz (Hg)</u>
B260CYSC	-10 to 200C	1.0C	Hg	Tol	16"	\$37.82	DE, GB	UN3506
B120CYSC	-30 to 120F	1.0F	Hg	Tol	12"	\$32.40	GB	UN3506
B230FYSC	0 to 230F	2.0F	Hg	Tol	12"	\$32.40	DE, GB	UN3506
B300FYSC	0 to 300F	2.0F	Hg	Tol	12"	\$32.40	DE, GB	UN3506
B500FYSC	20 to 500F	2.0F	Hg	Tol	16"	\$37.82	DE, GB	UN3506
B110C230FWSC	-20 to 110C	1.0C	Hg	76mm	12"	\$37.82	GB	UN3506
	&							
	& 0 to 230F	2.0F						
B110C230FYSC	-20 to 110C	1.0C	Hg	76mm	12"	\$37.82	DE, GB	UN3506
	&							
	& 0 to 230F	2.0F						
B150C300FYSC	-20 to 150C	1.0C	Hg	76mm	12"	\$37.82	DE, GB	UN3506
	&							
	& 0 to 300F	2.0F						
B260C500FYSC	-10 to 260C	1.0C	Hg	76mm	12"	\$43.21	DE, GB	UN3506
	&							
	& 20 to 500F	2.0F						

Blue-Spirit Pocket Test Thermometers, Aluminum Cases

- Blue Spirit_Filled
- In Enc and Armor Aluminum Cases
- Aluminum Cases have Pocket Clips
- Length 6" (155mm)
- Permanent Markings
- Toughened Bulbs
- Tol Immersion
- Thermometers Replacable

<u>Catalog No.</u>	<u>Range</u>	<u>Div</u>		<u>Imm</u>	<u>Cs Typ</u>			
<u>w/ Full Pocket Case</u>								
BM2011PBLs	-5 to 50c	1.0C	Spirit	Tol	Enc	\$21.62	GB	
BM2015PBLs	-15 to 105C	1.0C	Spirit	Tol	Enc	\$21.62	GB	
BM2012PBLs	-30 to 120F	2.0F	Spirit	Tol	Enc	\$21.62	GB	
BM2020PBLs	0 to 220F	2.0F	Spirit	Tol	Enc	\$21.62	GB	

w/ Armor Case

BM2011ABLs	-5 to 50c	1.0C	Spirit	Tol	OA	\$21.62	GB	
BM2015ABLs	-15 to 105C	1.0C	Spirit	Tol	OA	\$21.62	GB	
BM2012ABLs	-30 to 120F	2.0F	Spirit	Tol	OA	\$21.62	GB	
BM2020ABLs	0 to 220F	2.0F	Spirit	Tol	OA	\$21.62	GB	

Pocket Test Thermometer Refill

B2011BLS	-5 to 50c	1.0C	Spirit	Tol	None	\$13.52	GB	
B2015BLS	-15 to 105C	1.0C	Spirit	Tol	None	\$13.52	GB	
B2012BLS u	-30 to 120F	2.0F	Spirit	Tol	None	\$13.52	GB	
B2020BLS	0 to 220F	2.0F	Spirit	Tol	None	\$13.52	GB	

Red-Spirit Pocket Test Thermometers, Aluminum Cases

- Red Spirit_Filled
- In Enc and Armor Aluminum Cases

<u>Catalog No.</u>	<u>Range</u>	<u>Div</u>	<u>Filling</u>	<u>Imm</u>	<u>L</u>	<u>LIST PRICE</u>	<u>COO</u>	<u>Haz (Hg)</u>
--------------------	--------------	------------	----------------	------------	----------	-------------------	------------	-----------------

- Aluminum Cases have Pocket Clips
- Length 6" (155mm)
- Permanent Markings
- Toughened Bulbs
- Tol Immersion
- Thermometers Replacable

<u>Catalog No.</u>	<u>Range</u>	<u>Div</u>	<u>Imm</u>	<u>Cs Typ</u>
--------------------	--------------	------------	------------	---------------

w/ Full Pocket Case

BM2011PS	-35 to 50C	1.0C	Spirit	Tol	Enc	\$21.62	GB
BM2019PS	-10 to 110C	1.0C	Spirit	Tol	Enc	\$21.62	GB
BM2015PS	-15 to 105C	1.0C	Spirit	Tol	Enc	\$21.62	GB
BM2012PS	+30 to 120F	1.0F	Spirit	Tol	Enc	\$21.62	GB
BM2014PS	-30 to 120F	2.0F	Spirit	Tol	Enc	\$21.62	GB
BM2020PS	0 to 220F	2.0F	Spirit	Tol	Enc	\$21.62	GB

w/ Armor Case

BM2011AS	-35 to 50C	1.0C	Spirit	Tol	Ar	\$21.62	GB
BM2019AS	-10 to 110C	1.0C	Spirit	Tol	Ar	\$21.62	GB
BM2015AS	-15 to 105C	1.0C	Spirit	Tol	Ar	\$21.62	GB
BM2012AS	+30 to 120F	1.0F	Spirit	Tol	Ar	\$21.62	GB
BM2014AS	-30 to 120F	2.0F	Spirit	Tol	Ar	\$21.62	GB
BM2020AS	0 to 220F	2.0F	Spirit	Tol	Ar	\$21.62	GB

Pocket Test Thermometer Refill

B2011S	-35 to 50C	1.0C	Spirit	Tol	None	\$12.22	GB
B2015S	-15 to 105C	1.0C	Spirit	Tol	None	\$12.22	GB
B2012S	+30 to 120F	1.0F	Spirit	Tol	None	\$12.22	GB
B2014S	-30 to 120F	2.0F	Spirit	Tol	None	\$12.22	GB

Pocket Test Thermometers, Clear-Plastic Case, Spirit-Filled

- Red Spirit_Filled
- In Clear Plastic Pocket Case
- Cases have Pocket Clips
- Length 6" (155mm)
- Permanent Markings
- Toughened Bulbs
- Tol Immersion

BCP2011PS	-35 to 50C	1.0C	Spirit	Tol	CP	\$21.62	GB
BCP2015PS	-15 to 105C	1.0C	Spirit	Tol	CP	\$21.62	GB
BCP2012PS	+30 to 120F	1.0F	Spirit	Tol	CP	\$21.62	GB
BCP2014PS	-30 to 120F	2.0F	Spirit	Tol	CP	\$21.62	GB
BCP2020PS	0 to 220F	2.0F	Spirit	Tol	CP	\$21.62	GB

PRYMARY BRAND

LABORATORY THERMOMETERS

By THERMCO PRODUCTS, INC

- General purpose laboratory glass thermometers
- Red Spirit liquid filling Suffix "S"

<u>Catalog No.</u>	<u>Range</u>	<u>Div</u>	<u>Filling</u>	<u>Imm</u>	<u>L</u>	<u>LIST PRICE</u>	<u>COO</u>	<u>Haz (Hg)</u>
--------------------	--------------	------------	----------------	------------	----------	-------------------	------------	-----------------

- Blue Spirit liquid filling Suffix "BLS"
- Uniquid Non-Roll devise attached
- Packaged in square plastic tubes for storage
- Verticle numbers for easy reading
- Permanent markings, numbers and lines
- Partial and Tol immersions
- Several Ranges and Lengths

<u>Catalog No.</u>	<u>Range</u>	<u>Div</u>	<u>Imm</u>	<u>Length</u>
--------------------	--------------	------------	------------	---------------

Partial Immersion 12" (305mm) and 8" (200mm)

CS110CW3S	-20 to 110C	1.0C	RS	76mm	12"	\$4.11	CN
CS110CW3BLS	-20 to 110C	1.0C	BS	76mm	12"	\$6.06	CN
CS150CW3S	-20 to 150C	1.0C	RS	76mm	12"	\$4.11	CN
CS150CW3BLS	-20 to 150C	1.0C	BS	76mm	12"	\$6.06	CN
CS230FW3S	0 to 230F	2.0F	RS	76mm	12"	\$4.11	CN
		1.0C					
	-20 to 110C	&					
CS150C300FW3S	& 0 to 300F	2.0F	RS	76mm	12"	\$26.64	CN

Tol Immersion 12" (305mm) and 8" (200mm)

CS110CWS	-20 to 110C	1.0C	RS	Tol	12"	\$4.11	CN
CS110CWBLs	-20 to 110C	1.0C	BS	Tol	12"	\$6.06	CN
CS150CWS	-20 to 150C	1.0C	RS	Tol	12"	\$4.11	CN
CS150CWBLs	-20 to 150C	1.0C	BS	Tol	12"	\$6.06	CN
CS230FWS	0 to 230F	2.0F	RS	Tol	12"	\$4.11	CN
		1.0C					
	-20 to 110C	&					
CS150C300FWS	& 0 to 300F	2.0F	RS	Tol	12"	\$26.64	CN

<u>CAT. NO.</u>	<u>DESCRIPTION</u>	<u>HTS Code</u> 9025.19.8080	<u>LIST PRICE</u>	<u>COO</u>
<u>1" POCKET BI-METAL DIAL THERMOMETERS - 5" STEMS w/Stem Cover, Lexan Lens</u>				
ACCP110C	-10/110C x 1.0C		\$23.40	CN
ACCP150C	0/150C x 1.0C		\$23.40	CN
ACCP250C	0/250C x 5.0C		\$23.40	CN
ACCP120F	-40/120F x 2.0F		\$23.40	CN
ACCP125F	+25/125F x 1.0F		\$23.40	CN
ACCP220F	0/220F x 2.0F		\$23.40	CN
ACCP550F	+50/550F x 5.0F		\$23.40	CN
<u>1" POCKET BI-METAL DIAL THERMOMETERS - 5" STEMS w/Stem Cover, Lexan Lens, NSF Appro</u>				
YP125FNSF	+25/125F x 1.0F		\$26.69	CN
YP220FNSF	0/220F x 2.0F		\$26.69	CN
<u>1 3/4" LABORATORY BI-METAL DIAL THERMOMETERS - 8" STEMS, Lexan Lens</u>				
ACCG100CX	-50/100C x 1.0C		\$29.59	CN
ACCG110C	-10/110C x 1.0C		\$29.59	CN
ACCG150C	0/150C x 1.0C		\$29.59	CN
ACCG250C	0/250C x 2.0C		\$29.59	CN
ACCG125F	+25/125F x 1.0F		\$29.59	CN
ACCG160F	-40/160F x 2.0F		\$29.59	CN
ACCG180F	0/180F x 2.0F		\$29.59	CN
ACCG220F	0/220F x 2.0F		\$29.59	CN
ACCG300F	+50/300F x 2.0F		\$29.59	CN
ACCG500F	+50/500F x 5.0F		\$29.59	CN
ACCG550F	+50/550F x 5.0F		\$29.59	CN
<u>1 3/4" LABORATORY BI-METAL DIAL THERMOMETERS - 5" STEMS, Lexan Lens</u>				
ACCG160F-5	-40/160F x 2.0F		\$29.59	CN
ACCG220F-5	0/220F x 2.0F		\$29.59	CN
<u>1 3/4" LABORATORY BI-METAL DIAL THERMOMETERS - 8" STEMS, Lexan Lens</u>				
ACCG160F-C	-40/160F x 2.0F & -40/70C x1.0C		\$45.85	CN

<u>CAT. NO.</u>	<u>DESCRIPTION</u>	<u>HTS Code 9025.19.4000 (digital thermometer) 9025.90.0000 (probes for digital thermometer)</u>	<u>LIST PRICE</u>	<u>COO</u>	<u>Haz (Hg)</u>
<u>THERMCO'S DIGITAL BLOOD BANK THERMOMETER for BLOOD BAGS</u>					
HCBB1000	Digital Blood Bank Thermometer, -40 to +40C x 0.1C, ISO CERT @ 0C		\$693.00	US	
<u>PRECISION PT100 DIGITAL STEM THERMOMETER</u>					
ACCD0215	V215 Pt100 Precision Digital Thermometer, -50/250C, 8" stem, 30 x 4.5mm, includes DaKKS calibration certificate		316.72	DE	
ACCD0315	V215 Pt100 Precision Digital Thermometer, -50/250C, 12" stem, 30 x 4.5mm, includes DaKKS calibration certificate		316.72	DE	
ACCD0216	USB Adapter for V215 & V315 Pt100 Precision Digital Thermometer, allows USB export of data		134.58	DE	
<u>PT100 PLATINUM HIGH PRECISION THERMOMETERS, RESOLUTION 0.001°C</u>					
<u>DUAL PROBE HAND HELD, W/RS232</u>					
ACCD795P	Pt100 High Precision Thermometer, Dual Probe, -200/850C & F X 0.001		\$1,612.26	DE	
ACCD790P	Pt100 High Precision Thermometer, Dual Probe, -200/1760C & F X 0.001, thermocouple & humidity w/ automatic probe adjustment and detection, bluetooth		\$1,685.99	DE	
ACCD1723	High Precision Immersion Probe, -200 to 450C & F, 300mm x 1.5mm		\$448.29	DE	
ACCD1772	High Precision Immersion Probe, -200 to 450C & F, 230mm x 3.0mm		\$448.29	DE	
ACCD1773	High Precision Immersion Probe, -200 to 450C & F, 150mm x 3.0mm		\$422.07	DE	
ACCD1774	High Precision Immersion Probe, -200 to 450C & F, 300mm x 3.0mm		\$448.29	DE	
ACCD1778	High Precision Immersion Probe, -200 to 450C & F, 150mm x 6.0mm		\$422.07	DE	
ACCD1779	High Precision Immersion Probe, -200 to 450C & F, 630mm x 6.0mm		\$532.18	DE	
ACCD1784	High Precision Immersion Probe, -200 to 450C & F, 300mm x 6.0mm		\$491.55	DE	
ACCD1793	High Precision Immersion Probe, -200 to 450C & F, 375mm x 6.0mm		\$491.55	DE	
ACCD0081PD	Windows Software Smartgraph for On-Line Documentation		\$258.24	DE	
ACCD0007PD	Service Case for Thermometer, Probe & Accessories		\$146.81	DE	
ACCD10795	Power Supply		\$86.39	DE	
<u>PT100 PLATINUM HIGH PRECISION THERMOMETERS, RESOLUTION 0.01°C</u>					
<u>SINGLE & DUAL PROBE HAND HELD, W/RS232</u>					
ACCD700PEX	Pt100 High Precision Thermometer, Single Probe, -200/850C & F, Resolution 0.1C Intrinsically Safe/Explosion Proof		\$1,165.29	DE	
ACCD750P	Pt100 High Precision Thermometer, Single Probe, -200/850C & F, Resolution 0.01C		\$1,009.30	DE	
ACCD755P	Pt100 High Precision Thermometer, Dual Probe, -200/850C & F, Resolution 0.01C		\$1,123.35	DE	
ACCD750PEX	Pt100 High Precision Thermometer, Single Probe, -200/850C & F, Resolution 0.01C Intrinsically Safe/Explosion Proof		\$1,684.36	DE	
ACCD755PEX	Pt100 High Precision Thermometer, Dual Probe, -200/850C & F, Resolution 0.01C Intrinsically Safe/Explosion Proof		\$1,852.12	DE	
ACCD755LOG	Pt100 High Precision Thermometer & data logger, Dual Probe, -200/850C & F, Resolution 0.01C		\$1,169.24	DE	
ACCD700P	Pt100 High Precision Thermometer, Single Probe, -200/850C & F, Resolution 0.1C		\$482.96	DE	
ACCD1018	Pt100 Immersion Probe, 150mm x 3mm, -200/450C & -330/840F		\$393.25	DE	
ACCD1019	Pt100 Immersion Probe, 300mm x 3mm, -200/450C & -330/840F		\$416.84	DE	
ACCD1023	Pt100 Immersion Probe, 300mm x 1.5mm, -200/450C & -330/840F		\$424.69	DE	
ACCD1110	Pt100 Immersion Micro Probe, 1" x 0.03", -200/420C & -330/788F		\$326.39	DE	
ACCD1078	Pt100 Immersion Probe, 300mm x 6mm, -200/450C & -330/840F		\$427.32	DE	
ACCD1001	Pt100 Temperature & Humidity Probe, -30/80C & 0/100% Rh, 120mm x 20mm ABS		\$427.32	DE	
ACCD1001PT	Pt100 Immersion Probe, -50/350C, 150mm x 3.0mm		\$172.59	DE	
ACCD1002	Pt100 Immersion Probe, 300mm x 3mm, -50/350C & -58/660F		\$222.85	DE	
ACCD1005	Pt100 Immersion Probe, 500mm x 3mm, -50/350C & -58/660F		\$222.85	DE	
ACCD1007	Pt100 Insertion Probe, 300mm x 4mm, -50/350C & -58/660F		\$222.85	DE	
ACCD1023S	Pt100 Probe, 50mmx1mm		\$326.39	DE	
ACCD1075	Pt100 Surface Probe, w/Handle & Buffer Block, 150mm x 4mm, -40/300C & -40/572F		\$280.52	DE	
ACCD1055	Pt100 Air Probe, 250mm x 4mm, -50/250C		\$231.46	DE	
ACCD1059	Pt100 Surface Probe, Self-Adhesive, Bendable, 35mm x 13mm x 2mm, -20/250C & -4/480F		\$280.52	DE	
ACCD1082	Pt100 Immersion Probe for Tanks		\$523.02	DE	
ACCD1083	Pt100 Probe, M8thread, -100/450C, 50x30mm		\$300.18	DE	

<u>CAT. NO.</u>	<u>DESCRIPTION</u>	<u>HTS Code 9025.19.4000 (digital thermometer)</u>	<u>LIST PRICE</u>	<u>COO</u>	<u>Haz (Hg)</u>
ACCD0020PT-1.8	RTD Probe, Cable 1.8 Meters, -50/180C, 50x3mm		\$234.18	DE	
ACC50900002	PC-Adapter cable for RS232 interface		\$117.99	DE	
ACCD0070	USB cable and power supply		\$60.85	DE	
ACCD0081PD	Windows Software Smartgraph for On-Line Documentation		\$258.24	DE	
ACCD0007PD	Service Case for Thermometer, Probe & Accessories		\$146.81	DE	
ACC50900017	Protective Splash Cover for P700 Series, w/ strap		\$87.83	DE	
ACCD0092	Rubber Boot		\$31.47	DE	
ACCD0074	Extension Cable		\$66.86	DE	
ACCD0072	Probe Connector for 4000, 700, 795 Series		\$15.92	DE	

PT100 PLATINUM HIGH PRECISION THERMOMETERS,
 Grade "A" Pasteurized Milk Ordinance

ACCD750PD12C	PMO High Precision Digital, -200/850 C&F w/ 12" probe, Calibrated @ 0, 63, 72, 83, 89C		\$1,800.27	DE	
ACCD750PD12F	PMO High Precision Digital, -328/1760 C&F w/ 12" probe, Calibrated @ 32, 145, 161.6, 180, 191F		\$1,800.27	DE	
ACCD750PD6C	PMO High Precision Digital, -200/850 C&F w/ 6" probe, Calibrated @ 0, 63, 72, 83, 89C		\$1,713.35	DE	
ACCD750PD6F	PMO High Precision Digital, -328/1760 C&F w/ 6" probe, Calibrated @ 32, 145, 161.6, 180, 191F		\$1,713.35	DE	

"SUPRA"-PRECISION

PT100 PLATINUM HIGH PRECISION THERMOMETERS, RESOLUTION 0.01°C
 HAND HELD, MODEL CT376DIG W/16,000 DATA LOGGER

CT375DIG	RTD, -100/400C & -148/752F X 1.01C/0.02F (includes CTTPR04)		\$485.00	TW	
CT376DIG	Pt100 Platinum Thermometer w/16,000 Temp Data Logger		\$554.87	TW	
CTTPR04	Pt 100 Immersion Probe		\$170.14	TW	
CTTPR05	Pt 100 Air Probe		\$186.27	TW	
CTTPR06	Pt 100 Insertion/Piercing Probe		\$170.14	TW	
CT1776RS	RS-232 Adapter		\$56.76	TW	
CT1776AC	AC Adapter		\$45.87	TW	

GENERAL PURPOSE HAND-HELD Pt100 PLATINUM DIGITAL THERMOMETER,
 WATER PROOF, CERTIFIED

ACCD370P	Pt100 Platinum Thermometer, -100/300C & F x 0.1C & F		\$430.73	TW	
ACCD372P	Pt100 Platinum Thermometer, Dual Probe, -100/300C & F x 0.1C & F		\$503.61	TW	
ACCD370PFC	ACCD370P Certified at -20, 0, 4, 25, 37, 56, 95, 121C		\$648.06	TW	
ACCD370PUFC	ACCD370P Certified at -80, -20, 0, 4, 25, 37C		\$648.06	TW	
ACCD370PDFC	ACCD370P Certified at -20, 0, 4, 23, 32		\$575.17	TW	
ACCD370PWFC	ACCD370P Certified at -20, 0, 4, 25, 37, 44.5, 105, 121, 180C		\$648.06	TW	
CT370PT	Pt100 Probe for ACCD370P thermometers		\$116.14	TW	

GENERAL PURPOSE HAND-HELD Pt100 PLATINUM DIGITAL THERMOMETER, W/RS232

ACCD4000PW	IP 65 Water Proof Pt100 Platinum Digital Thermometer w/o Probes		\$272.77	DE	
ACCD4000P	Pt100 Digital 1-Channel, (Pt100-4 wire), -200/850C, USB, w/o Probe		\$229.71	DE	
ACCD4005P	Pt100 Digital 2-Channel, (Pt100-4 wire), -200/850C, USB, w/o Probe		\$273.15	DE	
ACCD3001	Pt100 Immersion Probe, -50/350C & F, 150mm x 3mm		\$143.14	DE	
ACCD3001W	Pt100 Immersion Probe, -50/350C & F, 150mm x 3mm; Waterproof for ACCD400PW		\$196.61	DE	
ACCD3002	Pt100 Immersion Probe, -50/350C & F, 300mm x 3mm		\$171.46	DE	
ACCD3006	Pt100 Insertion Probe, -50/350C & F, 150mm x 3mm		\$152.59	DE	
ACCD0016	Carrying/Storage Case		\$59.38	DE	

HIGH PRECISION (+/- 0.2C) THERMOCOUPLE DIGITAL THERMOMETER, K,J,L,N,T, W/RS232

ACCD705P	Pt100 Precision Thermometer, Dual Probe, -200/850C & F, Resolution 0.1C &		\$576.50	DE	
ACCD710P	Thermomcouple Thermometer, Single Channel, -200C to 1760C & F, Resolution 0.1C		\$514.33	DE	
ACCD715P	Thermomcouple Thermometer, Dual Channel, -200C to 1760C & F, Resolution 0.1C		\$518.42	DE	
ACCD715PEX	Thermomcouple Thermometer, Dual Channel, -200C to 1760C & F, Resolution 0.1C Intrinsically Safe		\$951.90	DE	

<u>CAT. NO.</u>	<u>DESCRIPTION</u>	<u>HTS Code 9025.19.4000 (digital thermometer)</u>	<u>LIST PRICE</u>	<u>COO</u>	<u>Haz (Hg)</u>
ACCD700P	Pt100 & Thermocouple Thermometer, Single Channel PT100 -200/850C & F and Thermocouple -200/1760C & F, Resolution 0.1C		\$482.96	DE	
ACCD1016	Thermocouple K Immersion Probe, -40 to 400C & F, 120mm x 3.5mm		\$113.24	DE	
ACCD1010	Thermocouple K Immersion Probe, -100 to 1100C & F, 300MM x 3mm		\$156.64	DE	
ACCD1007K	Thermocouple K Immersion Probe, -100 to 1100C & F, 500MM x 3mm		\$222.85	DE	
ACCD1020	Thermocouple K Surface Probe w/90 degree bend, -40 to 900C & F		\$167.09	DE	
ACCD0048	Thermocouple Probe Adapter for Mini Connector		\$56.76	DE	
ACCD0088-6	Thermocouple w/steel braided cable, 6000mm Cable, -40/480C		\$83.90	DE	
ACC50900002	PC-Adapter cable for RS232 interface		\$117.99	DE	
ACCD0035PD	USB Cable & Adapter for RS232 interface		\$108.01	DE	
ACCD10795	USB Power Pack for Above Thermometers		\$86.39	DE	
ACCD0070	Power Supply		\$60.85	DE	
ACCD0081PD	Windows Software Smartgraph for On-Line Documentation		\$258.24	DE	
ACCD0007PD	Service Case for Thermometer, Probe & Accessories		\$146.81	DE	
ACC50900017	Protective Splash Cover		\$86.39	DE	
ACC99050005	Heat-Conductive Paste 20g		\$26.07	DE	

PRECISION THERMOCOUPLE THERMOMETER w/NIST TRACEABLE CERTIFICATE AND RS232

No certificate for individual probe purchase

ACCD4010K	1 Channel (Type K), -200/1760C, USB, w/ Probe		\$201.15	DE	
ACCD4015K	2 Channel (Type K), -200/1760C, USB, w/ Probe		\$235.91	DE	
ACCD0013P	K Thermocouple Thermometer, -100C to 1370C & F, +/- 0.5C & 1.0F		\$102.25	DE	
ACCD0507	K Thermocouple Immersion Probe w/handle, -40 to 400C & F, 120mm x 3mm		\$70.28	DE	
ACCD0501	K Thermocouple Immersion Probe w/handle, -100 to 1100C & F, 300mm x 1.5mm		\$128.33	DE	
ACCD0503	K Thermocouple Immersion Probe w/handle, -100 to 1100C & F, 300mm x 3.0mm		\$137.78	DE	
ACCD0511	K Thermocouple Insertion Probe w/handle, -100 to 1100C & F, 300mm x 4.0mm		\$156.64	DE	
ACCD0062	K Thermocouple beaded wire probe, teflon insulated, -40 to 250C & F		\$40.50	DE	
ACCD0016	Carring/Storage Case for ACCD410T		\$59.38	DE	

PRECISION THERMOCOUPLE THERMOMETER, 2 & 4 CHANNEL DATA LOGGER W/
NIST TRACEABLE CERTIFICATE AND RS232

CT305	K-Type -200/1370C, 328/2498F x 0.1C&F, w/ logger, PC Interface		\$253.26	TW	
CT306	2 Channel Thermocouple Data Logger, -200 to 1370C (-328 to 2498F)		\$332.17	TW	
CT309	4 Channel Thermocouple Data Logger, -200 to 1370C (-328 to 2498F)		\$415.79	TW	
CT307K	Single Channel Thermocouple -200 to 1370C (-328 to 2498F)		\$137.78	TW	
CTK01	Thermocouple Beaded Wire Probe, -50 to 200C		\$23.59	TW	
CTK02	Thermocouple Immersion Probe, -50 to 700C		\$68.19	TW	
CTK04	Thermocouple Insertion Probe, -50 to 600C		\$82.32	TW	
CTK05	Thermocouple, Surface Probe, -50 to 400C		\$109.33	TW	
CTK06	Thermocouple, Air/Gas Probe, -50 to 800C		\$109.33	TW	

GENERAL PURPOSE THERMOCOUPLE THERMOMETER, K-TYPE SINGLE & DOUBLE PROBE MODELS
W/°C,°F,°KMIN,MAX,AVE, HOLD TI/T2, MEMORY, HI/LOW ALARM, REAL TIME CLOCK

ACC950K	SINGLE CHANNEL TYPE K THERMOCOUPLE THERMOMETER		\$161.24	CN	
ACCD952K	TYPE K GENERAL PURPOSE PROBE W/HANDLE		\$99.24	CN	
ACC954K	TYPE K PENETRATION PROBE W/HANDLE		\$86.11	CN	
ACC958K	TYPE K SURFACE PROBE W/HANDLE		\$97.38	CN	

GENERAL PURPOSE DIGITAL THERMOMETER W/PROBE, MAX HOLD FUNCTIONS

ACCD300P	DIGITAL THERMOMETER, C & F, -40 to 200C & -40 to 400F X 0.1C & F		\$140.28	DE	
----------	--	--	----------	----	--

GENERAL LABORATORY DIGITAL THERMOMETERS W/ NIST TRACEABLE CERTIFICATE

ACCD101DIG	Digital Thermometer w/4" Insertion Probe, -40/200C		\$72.37	CN	
ACCD102DIG	Digital Thermometer w/Probe, -50 to 70C & F x 0.1 C & F, Min-Max, Alarm, Stand		\$69.48	CN	
ACCD200DIG	Gen. Purpose Digital Thermometer w/4" Probe, -40/200C & -40/400F		\$70.92	DE	
ACC300DIG	1 1/2" Pocket Dial Thermometer, -58/302F & -50/150C, 5" Stem		\$37.50	CN	

<u>CAT. NO.</u>	DESCRIPTION	HTS Code 9025.19.4000 (digital thermometer)	<u>LIST PRICE</u>	<u>COO</u>	<u>Haz (Hg)</u>
ACC9291DIG	2" Pocket Dial Thermometer, -58/302F & -50/150C, 5" Stem		\$56.27	DE	
ACC304MCF1	Waterproof & NSF approved , NIST Compliant, 1½" Dial, 2½" Stem, -50/200C & -58/392F x 0.1, Min-Max		\$36.18	CN	
ACC310DIG	Digital Lab Thermometer, -58/536F & -50/280C, 8" Stem		\$37.50	CN	
ACC310SUP	Digital Lab Thermometer, -58/536F & -50/280C, 8" Stem accuracy 0.2C & 0.5C		\$58.86	CN	
ACC310DLS	Digital Lab Thermometer, Long Probe, -58/536F & -50/280C, 12" Stem		\$49.94	CN	
ACC310DLSUP	Digital Lab Thermometer, Long Probe, -58/536F & -50/280C, 12" Stem, accuracy 0.5C		\$64.76	CN	
ACC3105DIG	Digital Lab Thermometer, -58/302F & -50/150C, 5" Stem		\$37.50	CN	
ACC3105SUP	Digital Lab Thermometer, -58/302F & -50/150C, 5" Stem, accuracy 0.2C & 0.5C		\$58.48	CN	
ACC3108DIG	Digital Lab Thermometer, -58/392F & -50/200C, 8" Stem		\$37.50	CN	
ACC3108SUP	Digital Lab Thermometer, -58/302F & -50/150C, 8" Stem, accuracy 0.2C & 0.5C		\$58.48	CN	
ACC310DMM	Digital Lab Thermometer, -50/200C, 8" Stem w Min-Max Recall		\$45.87	CN	
ACC310MMSUP	Digital Lab Thermometer, -58/392F & -50/200C, 8" Stem w Min-Max, accuracy 0.2C & 0.5C		\$58.60	CN	
ACC320DIG	"T" Barrel Head Thermometer, -58/302F & -50/150C, 5" Stem		\$37.50	CN	
ACC330DIG	"Non-Roll" Thermometer, -58/302F & -50/150C, 5" Stem		\$37.50	CN	
ACC350DIG	Waterproof, Rotary Head Thermometer, -58/302 & -50/150C		\$37.50	CN	
ACC370DIG	Deluxe Waterproof Pocket Thermometer. 31/2" Probe, -50/300C & -58/572F		\$37.50	CN	
ACC370SUP	Deluxe Waterproof Pocket Thermometer, -50/300C & -58/572F, accuracy 0.2C & 0.5C		\$58.86	CN	
ACC400DIG	Waterproof/Shockproof Digital Thermometer, -50/300C & -58/572F, 5" Stem		\$37.50	CN	
ACC400SUP	Waterproof/Shockproof Digital Thermometer, -50/300C & -58/572F, 5" Stem, accuracy 0.5C		\$58.48	CN	
ACC410PW	Digital Ultra Low Freezer Thermometer, -100/200C x 0.1C		\$99.34	CN	
ACC420DIG	Waterproof/Shockproof Digital Thermometer, -50/150C & -58/302F, 5" Stem		\$37.50	CN	
ACC420SUP	Waterproof/Shockproof Digital Thermometer, -50/300C & -58/572F, 5" Stem, accuracy 0.5C		\$58.48	CN	
ACCD1017MM	MIN/MAX -40/70C & -40/158F x 0.1C&F, NIST Compliant w/cert		\$55.57	DE	
ACCD0552DIG	"Flip-Probe" Digital Pocket Thermometer, -50/350C & -58/662F, 4-1/2" Probe		\$110.50	DE	
ACCD0553DIG	"Flip-Probe" Digital Pocket Thermometer, -40/250C & -40/482F, 2-3/4" Probe		\$61.09	DE	
ACC500DIG	"Flip Probe" Digital Thermometer, -50/300C & -58/572F		\$41.70	CN	
ACC550DIG	N.I.S.T. Pocket Thermometer w/Probe, 31/2" Cable, -58/572F & -50/300C		\$58.48	CN	
ACC600DIG	Waterproof Pocket Digital Thermometer, -58/302F & -50/150C, 4 1/2" Probe		\$37.50	CN	
ACC600SUP	Waterproof Pocket Digital Thermometer, -58/302F & -50/150C, 4 1/2" Probe, accuracy 0.5C		\$58.48	CN	
ACC655DIG	Waterproof Pocket Digital Thermometer, -58/572F & -50/300C, 4 1/2" Probe		\$37.50	CN	
ACC655SUP	Waterproof Pocket Digital Thermometer, -58/572F & -50/300C, 4 1/2" Probe, accuracy 0.5C		\$58.48	CN	
ACC6558DIG	Waterproof Pocket Digital Thermometer, -58/572F & -50/300C, 8" Probe		\$51.28	CN	
ACC6558SUP	Waterproof Pocket Digital Thermometer, -58/572F & -50/300C, 8" Probe, accuracy 0.5C		\$64.76	CN	
ACC610DIG	Waterproof, T-Handle Thermometer, -58/302F & -50/150C, 4 1/2" Probe		\$58.24	CN	
ACC610SUP	Waterproof, T-Handle Thermometer, -58/302F & -50/150C, 4 1/2" Probe, accuracy 0.5C		\$58.48	CN	
ACC800DIG	In-Door/Out-Door Digital Thermometer -10/50C INT, -50/70C EXT, Min-Max		\$40.50	CN	
ACC801DIG	Digital Thermometer w/Probe, -50/70C & F, Min-Max, Alarm		\$54.13	CN	
ACC840DIG	In-Door/Out-Door Digital Thermometer, Large Digit -0/50C&F & -50/70C&F, Min-Max,		\$51.28	CN	
ACC850DIG	Triple Display Digital Thermometer, Current, Min, Max -50 to 70C & F		\$56.76	CN	
ACC895WDIG	Large Digit Triple Screen Digital Thermometer, -50/70C & F, Display Current Temp, probe only		\$58.48	CN	
ACC9242DIG	Pocket Digital Thermometer w/4" Probe, -50/150C & -58/302F, MIN-MAX		\$54.13	DE	
ACC9270DIG	Digital Thermometer w/ Timer, Alarm, Min-Max, Stainless Steel Probe		\$58.48	DE	
ACC9276DIG	Remote Probe Thermometer, -50/300C & -58/572F		\$47.33	DE	
301025C	Window Thermometer, 3 Readings, Current, Minimum, Maximum in F		\$40.50	CN	
301025F	Window Thermometer, 3 Readings, Current, Minimum, Maximum in C		\$40.50	CN	

DIGITAL MINIMUM / MAXIMUM THERMOMETERS

ACCD1017MM	MIN/MAX -40/70C & -40/158F x 0.1C&F, NIST Compliant w/cert	\$55.57	DE
ACCD1039MM	MIN/MAX -20/70C & -4/158F x 0.1C&F, NIST Compliant w/cert	\$55.57	DE

TRIPLE DISPLAY LARGE DIGIT ENVIRONMENTAL CHAMBER THERMOMETERS

CURRENT, MIN & MAX DISPLAYED, CERTIFIED, w/ALARM & Probe in 30ml vial

ACC895FRE	For Freezers, -50 to 70C, Certified at -20C	\$72.37	CN
ACC895REF	For Refrigerators, -50 to 70C, Certified at +4C	\$72.37	CN

<u>CAT. NO.</u>	<u>DESCRIPTION</u>	<u>HTS Code 9025.19.4000 (digital thermometer)</u>	<u>LIST PRICE</u>	<u>COO</u>	<u>Haz (Hg)</u>
ACC895INC	For Incubators, -50 to 70C, Certified at +37C		\$72.37	CN	
ACC895WB	For Block Heaters, -50 to 70C, Certified at +37C		\$72.37	CN	
ACC895AMB	For Ambient/rooms -55 to 70C, Certified at +22C		\$72.37	CN	
<u>w/ PROBE IN 6ml VIAL</u>					
ACC895FRE	For Freezers, -50 to 70C, Certified at -20C		\$72.37	CN	
ACC895REFV	For Refrigerators, -50 to 70C, Certified at +4C		\$72.37	CN	
ACC895INCV	For Incubators, -50 to 70C, Certified at +37C		\$72.37	CN	
<u>w/ PROBE IN 1.5ml MICRO-CENTRIFUGE TUBE</u>					
ACC895BLK	For Dry-Bath/Block Heaters, Certified +37C		\$72.37	CN	
<u>TRIPLE DISPLAY LARGE DIGIT W/WIRELESS SENSOR DIGITAL THERMOMETER</u>					
<u>CURRENT, MIN & MAX DISPLAYED, CERTIFIED, w/ALARM NO PROBE CABLE REQUIRED</u>					
ACC895WRFT	For Refrigerators/Freezers/Ambient, -50/70C, Certified at +4C & -20C		\$74.91	CN	
ACC895WRFTCOV	Accessory - protective plastic cover for ACC895WRFT		\$5.86	CN	
<u>TRIPLE DISPLAY LARGE DIGIT ENVIRONMENTAL CHAMBER THERMOMETERS</u>					
ACC895ULF	For -80/70C ultra low freezer or refrigerator, air probe, no buffer, Certified at -70C		\$68.93	CN	
<u>DIGITAL WORKSTATION THERMOMETERS, CERTIFIED</u>					
WS0013FC	Digital Workstation Thermometer, w/bare wire probe, Certified at 37C		\$110.10	DE	
SW0013FC	Slide Warmer Thermometer, Certified at 40C & 54C		\$129.13	DE	
<u>ENVIRONMENTAL CHAMBER DIGITAL THERMOMETER FOR ULTRA-LOW FREEZERS</u>					
<i>periodic monitoring</i>					
ACCULF307	-100/160C x 0.1C, C&F Digital w/60ml sand filled bottle, Cert @ -78C		\$207.24	HK	
ACCULF307FC	-100/160C x 0.1C, C&F Digital w/60ml sand filled bottle, Cert @ -80, 0, 20C		\$318.66	HK	
<u>TRIPLE DIGIT VACCINE THERMOMETERS</u>					
<u>MIN & MAX DISPLAYED, CERTIFIED, w/ALARM & TIME-DATE STAMP</u>					
ACC821FRE	Time/Date Stamp Digital Thermometer w/Audible & Visual Alarm for Freezers,		\$80.63	CN	
ACC821REFV	Time/Date Stamp Digital Thermometer w/Audible & Visual Alarm for Refrigerators		\$80.63	CN	
ACC821INC	Time/Date Stamp Digital Thermometer w/Audible & Visual Alarm for Incubators		\$80.63	CN	
ACC821DIG	Thermometer & Probe Only		\$66.72	CN	
<u>DUAL PROBE TRIPLE DIGIT VACCINE THERMOMETERS</u>					
<u>MIN & MAX DISPLAYED, CERTIFIED, w/ALARM & TIME-DATE STAMP</u>					
ACC821RFV2	Time/Date Stamp Digital Thermometer w/Audible & Visual Alarm for Freezers & Refrigerators		\$89.02	CN	
ACC821INC2	Time/Date Stamp Digital Thermometer w/Audible & Visual Alarm for Incubators		\$89.02	CN	
ACC821DIG2	Thermometer & Probe Only		\$80.63	CN	
<u>DUAL DISPLAY LARGE DIGIT VACCINE THERMOMETERS</u>					
<u>MIN & MAX DISPLAYED, CERTIFIED, w/ALARM</u>					
ACC8100DIG	Dual Display Refrigerator/Freezer, -40 to 70C x 0.1C Cert at -20C & +4C		\$75.10	CN	
ACC8002CL	Dual Display Refrigerator/Freezer -40 to 70C x 0.1C Cert at -20C & +4C W/2 VIALS		\$75.10	CN	
<u>DUAL DISPLAY LARGE DIGIT VACCINE THERMOMETERS</u>					
<u>MIN & MAX DISPLAYED, CERTIFIED, w/ALARM & TIME-DATE STAMP</u>					
ACC8100MATB	Dual Display Freezer/Refrigerator, -40 to 70C x 0.1C Cert at -20C & +4C		\$80.63	CN	
<u>TRIPLE DISPLAY DIGITAL ENVIRONMENTAL CHAMBER THERMOMETERS</u>					
<u>CURRENT, MIN & MAX DISPLAYED, CERTIFIED, w/ALARM</u>					
ACC850DIG	Thermometer only		\$56.76	CN	
ACC850FRE	Thermometer w/ Probe and 5ml Liquid Vial for Freezers		\$64.76	CN	
ACC850REF	Thermometer w/ Probe and 5ml Liquid Vial for Refrigerators		\$64.76	CN	

<u>CAT. NO.</u>	<u>DESCRIPTION</u>	<u>HTS Code 9025.19.4000 (digital thermometer)</u>	<u>LIST PRICE</u>	<u>COO</u>	<u>Haz (Hg)</u>
ACC850INC	Thermometer w/ Probe and 5ml Liquid Vial for Incubators		\$64.76	CN	
ACC850AMB	Thermometer w/ Probe and 5ml Liquid Vial for Ambient/Room		\$64.76	CN	
ACC850BLK	Thermometer w/ Probe and 5ml Liquid Vial for Block Heaters/Dry Baths		\$64.76	CN	
ACC850WB	Thermometer w/ Probe and 30ml Plastic Bottle of Sand for Water Baths		\$64.76	CN	

DUAL CHANNEL FREEZER-REFRIGERTOR THERMOMETER
CERTIFIED, MIN & MAX, w/ALARM

ACC8002CL	DUAL DISPLAY FREEZER/REFRIGERATOR -40 to 70C x 0.1C Cert at -20C & +4C		\$75.10	CN	
ACC8002CLV	DUAL PROBE w/ DUAL PROBE & VIAL -40 to 70C x 0.1C Cert at -20C & +4C		\$77.87	CN	

ENVIRONMENTAL CHAMBER DIGITAL THERMOMETERS w/ALARM, MIN/MAX, CERTIFIED

ACC801FRE	For Freezers, -50 to 70C, Certified at -20C		\$69.48	CN	
ACC801REF	For Refrigerators, -50 to 70C, Certified at +4C		\$69.48	CN	
ACC801INC	For Incubators, -50 to 70C, Certified at +37C		\$69.48	CN	
ACC801BLC	For Block Heaters, -50 to 200C, Certified at +37C		\$69.48	CN	
ACC801AMB	For Ambient/Rooms, -50 to 70C, Certified at +22C		\$69.48	CN	
ACC801WB	For Water Baths, -50 to 200C, Certified at +37C		\$69.48	CN	
ACC801OV	For Ovens, -50 to 200C, Certified at +100C		\$69.48	CN	
ACC801DIG	Thermometer & Probe Only		\$52.69	CN	

ENVIRONMENTAL CHAMBER DIGITAL THERMOMETERS w/ALARM
MIN-MAX, CERTIFIED, for REFRIGERATORS FREEZERS, INCUBATORS, OVENS

ACC900FRE	For Freezers		\$75.10	CN	
ACC900REF	For Refrigerators		\$75.10	CN	
ACC900INC	For Incubators		\$75.10	CN	
ACC900OV	For Ovens		\$75.10	CN	
ACC900DIB	Thermometer & Probe Only		\$55.57	CN	

ENVIRONMENTAL CHAMBER DIGITAL THERMOMETERS
MIN-MAX, CERTIFIED, for REFRIGERATORS FREEZERS, INCUBATORS, DRY BLOCKS, ROOMS
THERMOMETERS W/LIQUID FILLING, ACCURACY $\pm 1.0C$

ACC304FREL	For Freezers		\$58.48	CN	
ACC304REFL	For Refrigerators		\$58.48	CN	
ACC304AMBL	For Ambient/Rooms		\$58.48	CN	
ACC304INCL	For Incubators		\$58.48	CN	
ACC304DBLL	For Dry Blocks (Low)		\$58.48	CN	
ACC304DBHL	For Dry Blocks (High)		\$58.48	CN	

"SUPRA"

ACC304FRELS	For Freezers		\$72.37	CN	
ACC304REFLS	For Refrigerators		\$72.37	CN	
ACC304INCLS	For Incubators		\$72.37	CN	
ACC304AMBLS	For Ambient/Rooms		\$72.37	CN	
ACC304DBLLS	For Dry Blocks		\$72.37	CN	

THERMOMETERS W/SAND FILLING, ACCURACY $\pm 1.0C$

ACC304FRES	For Freezers		\$58.48	CN	
ACC304REFS	For Refrigerators		\$58.48	CN	
ACC304AMBS	For Ambient/Rooms		\$58.48	CN	
ACC304INCS	For Incubators		\$58.48	CN	

"SUPRA"

ACC304FRESS	For Freezers		\$72.37	CN	
ACC304REFSS	For Refrigerators		\$72.37	CN	
ACC304AMBSS	For Ambient/Rooms		\$72.37	CN	

<u>CAT. NO.</u>	<u>DESCRIPTION</u>	<u>HTS Code 9025.19.4000 (digital thermometer)</u>	<u>LIST PRICE</u>	<u>COO</u>	<u>Haz (Hg)</u>
ACC304INCSS	For Incubators		\$72.37	CN	
<u>WINDOW THERMOMETER</u>					
<u>WATER RESISTANT, °C & °F, CERTIFIED</u>					
301025C	Window Thermometer, 3 Readings, Current, Minimum, Maximum in F		\$40.50	CN	
301025F	Window Thermometer, 3 Readings, Current, Minimum, Maximum in C		\$40.50	CN	
<u>PRECISION INFRA-RED NON-CONTACT THERMOMETERS, NOT FOR HUMAN USE</u>					
311125IR	FLASH PEN, -33 TO 220C X 0.1, °C & °F		\$56.12	DE	
ACCD0414IR	MINIFLASH 2, -33 TO 220C X 0.1, °C & °F MIN, MAX, HOLD, LOCK 1.3: 1 RES.		\$63.97	DE	
ACCD0385IR	Economy IR Thermometer, w/laser, -35 to 365C & F x 0.1.		\$109.84	DE	
ACCD0556IR	FLASH 3, W/LASER SIGHTING -33 TO 220C X 0.1, °C & °F MIN, MAX, HOLD, LOCK 1.3: 1 RES.		\$97.01	DE	
ACCD0413IR	DUAL-TEMP, INFRA-RED AND THERMOCOUPLE PROBE, -55 to 220C & F, MIN, MAX, HOLD, LOCK		\$171.46	DE	
ACCD0503IR	SCAN-TEMP 410, -33 to 500C & F, LAZER SIGHTING, MAX, HOLD, 11:1,		\$186.27	DE	
ACCD0503IRK	SCAN-TEMP 440, -33 to 500C & F, LAZER SIGHTING, MAX, MIN, HOLD, LOCK, DIF, AVE, 11:1,		\$214.71	DE	
ACCD0485IR	HI-LOW ALARM, W/THERMOCOUPLE PROBE CONNECTOR, W/CASE SCAN-TEMP 485, -50 to 800C & F, ACCURACY ±2.0C, RES. 20:1, DOUBLE LAZER SIGHTING MAX, HOLD, LOCK FUNCTIONS, W/CASE		\$272.77	DE	
TCT8530	IR THERMOMETER, 12:1, -50/530C & -58/986F		\$102.64	CN	
TCT8580	IR THERMOMETER, 16:1, -50/580C & -58/1076F		\$221.40	CN	
TCT8850H	IR THERMOMETER, 16:1, -50/850C & -58/1562F		\$110.77	CN	
TCT8010H	IR THERMOMETER, 50:1, -50/1000C & -58/1832F		\$259.27	CN	
TCT8011H	IR THERMOMETER, 50:1, -50/1100C & -58/2012F		\$264.65	CN	
<u>PRECISION INFRA-RED NON-CONTACT THERMOMETERS, APPLICABLE FOR HUMAN USE & SURFACE</u>					
ACCD0478IR	Medi-Scan Forehead & Surface InfraRed Thermometer, 32/42.5C & 93.2/108.5F		\$100.62	CN	
ACCD04900IR	AccuTherm Care4U Forehead & Surface InfraRed Thermometer, 32/42.5C & 89.6/108.5F		\$66.56	CN	
<u>THERMO-HYGROMETERS, DEW POINT, WET BULB</u>					
CT317TH	Handheld High Precision Hygrometer/Thermometer		\$309.22	TW	
CT31TH	Hygrometer/Thermometer, High Precision		\$174.23	TW	
CTRP32	Extension Cable for CT31Th		\$39.21	TW	
CT315TH	Compact Hygrometer/Thermometer, High Precision		\$187.72	TW	
CT310TH	Hygrometer/Thermometer, High Precision		\$268.72	TW	
CT33CELL	Calibration Cell @ 33% Rh		\$66.21	TW	
CT75CELL	Calibration Cell @ 75% Rh		\$66.21	TW	
TCT908TH	Economy Digital Hygrometer/Thermometer w/clock		\$56.76	DE	
30500002	Economy Digital Hygrometer/Thermometer		\$47.33	DE	
ACCD5013TH	Indoor/Outdoor Hygrometer/Thermometer		\$89.14	DE	
ACCDTH03	Large Digit Hygrometer/Thermometer		\$68.19	DE	
ACCD560TH	Thermo-Hygrometer, Dew Point, Wet Bulb, Min-Max-Hold, 0/100%Rh, -40/70C		\$165.43	DE	
ACCDTA100	Thermo-Hygrometer, Dew Point, Wet Bulb, Min-Max-Hold, 0/100%Rh, -40/70C		\$100.15	DE	
ACCDTA100FC	Thermo-Hygrometer, Dew Point, Wet Bulb, Min-Max-Hold, 0/100%Rh, -40/70C		\$325.35	DE	
30515	Triple Screen Digital Thermo-Hygrometer displays Current Temperature & Rh Min/Max Temperature and Rh		\$66.72	DE	
ACC811CE	In & out Temperature, Humidity & Clock		\$68.19	CN	
ACC9216DIG	Large Digit, Indoor/Outdoor Thermo-Hygrometer, -50/70C & F, 20/100% Rh, Min-Max		\$68.19	CN	
ACCHT5063	Analog Thermo-Hygrometer, brass 6", glass, Certified		\$354.18	CN	
ACCHT5063A	Analog Thermo-Hygrometer, economy, Certified		\$124.26	CN	
TT8020TH	Indoor/Outdoor Hygrometer/Thermometer/Clock, -50/70C & -58/158F, 10/99%Rh		\$60.47	CN	
<u>SLING PSYCHROMETERS</u>					
BPS100FS	Sling Psychrometer w/ Spirit-Filled Thermometers, 0° to 120°F, Slide-Rule Calculator		\$97.38	GB	

<u>CAT. NO.</u>	<u>DESCRIPTION</u>	<u>HTS Code 9025.19.4000 (digital thermometer)</u>	<u>LIST PRICE</u>	<u>COO</u>	<u>Haz (Hg)</u>
BPS100CS	Sling Psychrometer w/ Spirit-Filled Thermometers, -5° to 50°C, Slide-Rule Calculator		\$97.38	GB	
BPS100F	Sling Psychrometer w/ Mercury-Filled Thermometers, 0° to 120°F, Slide-Rule Calculator		\$97.38	GB	UN3506
BPS100C	Sling Psychrometer w/ Mercury-Filled Thermometers, -5° to 50°C, Slide-Rule Calculator		\$97.38	GB	UN3506
<u>DIGITAL TACHOMETERS</u>					
ACCPTC900	MULTIFUNCTIONAL, PHOTO, CONTACT & SURFACE TACHOMETER		\$333.48	CN	
ACCDCT801	CONTACT (TOUCH) & PHOTO TACHOMETER		\$164.78	CN	
ACCDLC2234C	PHOTO-LASER (TOUCHLESS) TACHOMETER		\$137.78	CN	
<u>BAROMETERS AND WEATHER STATIONS</u>					
ACCD300632	ANALOG BAROMETER W/TEMPERATURE AND HUMIDITY		\$104.98	DE	

<u>CAT. NO.</u>	<u>DESCRIPTION</u>	<u>LIST PRICE</u>
EXPEDITE FEE	3-5 BUSINESS DAY PROCESSING TURN AROUND, RECALIBRATION SVC ONLY, *excludes transit time **some restrictions may apply	\$145.00
GENLAB SVC FEE	DIAGNOSTICS, REPAIR, LABOR, MINOR PARTS, FAILED SVC	\$40.00
<u>GLASS AND GENERAL LABORATORY DIGITAL THERMOMETERS</u>		
RECERT 1	RECERTIFICATION 1 TEMPERATURE POINT	\$65.00
RECERT 2	RECERTIFICATION 2 to 4 TEMPERATURE POINTS	\$107.00
RECERT 3	RECERTIFICATION 5 to 9 TEMPERATURE POINTS	\$129.75
RECERT 4	RECERTIFICATION 10 to 15 TEMPERATURE POINTS	\$151.50
CAL SVC	FEE FOR TEMPS BETWEEN 300/500C & 572/932F	\$50.00
	<i>*one time fee per service</i>	
RECERT B	RECERTIFICATION 1 TEMPERATURE POINT BOTTLES and DIGITALS (fully accredited cert; default temp)	\$43.25
<u>HIGH PRECISION Pt100 PLATINUM THERMOMETERS 1 PROBE</u>		
RECERT DIG 1	RECERTIFICATION 1 TEMPERATURE POINT	\$96.75
RECERT DIG 2	RECERTIFICATION 2 TEMPERATURE POINTS	\$151.50
RECERT DIG 3	RECERTIFICATION 3 TEMPERATURE POINTS	\$162.75
RECERT DIG 4	RECERTIFICATION 4 TEMPERATURE POINTS	\$173.00
RECERT DIG 5	RECERTIFICATION 5 TEMPERATURE POINTS	\$194.50
RECERT DIG 6	RECERTIFICATION 6 to 8 TEMPERATURE POINTS	\$215.25
CAL SVC	FEE FOR TEMPS BETWEEN 300/500C & 572/932F	\$50.00
	<i>*one time fee per service</i>	
<u>HIGH PRECISION Pt100 PLATINUM THERMOMETERS 2 PROBE</u>		
RECERT DIG 1/2	RECERTIFICATION 1 TEMPERATURE POINT	\$151.50
RECERT DIG 2/2	RECERTIFICATION 2 TEMPERATURE POINTS	\$215.50
RECERT DIG 3/2	RECERTIFICATION 3 TEMPERATURE POINTS	\$226.75
RECERT DIG 4/2	RECERTIFICATION 4 TEMPERATURE POINTS	\$237.00
RECERT DIG 5/2	RECERTIFICATION 5 TEMPERATURE POINTS	\$258.50
RECERT DIG 6/2	RECERTIFICATION 6 to 8 TEMPERATURE POINTS	\$279.00
CAL SVC	FEE FOR TEMPS BETWEEN 300/500C & 572/932F	\$50.00
	<i>*one time fee per service</i>	

<u>CAT. NO.</u>	<u>DESCRIPTION</u>	<u>LIST PRICE</u>
<u>INFRARED THERMOMETERS</u>		
RECERT IR 1	RECERTIFICATION 1 TEMPERATURE POINT	\$65.00
RECERT IR 2-4	RECERTIFICATION 2 to 4 TEMPERATURE POINTS	\$107.00
RECERT IR 5-6	RECERTIFICATION 5 to 6 TEMPERATURE POINTS	\$129.75
<u>HYDROMETERS:</u>		
REHYD 1	RECERTIFICATION 1 POINT	\$81.25
REHYD 2	RECERTIFICATION 2 POINTS	\$109.00
REHYD 3	RECERTIFICATION 3 POINTS	\$136.00
<u>Rh / HYGROMETERS W/ THERMOMETER 1 TEMPERATURE</u>		
RERh1	RECERTIFICATION Rh 1 POINT & 1 Temperature Point	\$81.25
RERh 3	RECERTIFICATION Rh 3 POINTS & 1 Temperature Point	\$126.75
<u>RECERTIFICATION OF VFC DATALOGGER</u>		
REVFCLOG 1	RECERTIFICATION OF 1 POINT VFC DATALOGGER	\$50.00
REVFCLOG 2	RECERTIFICATION OF 2 POINT VFC DATALOGGER	\$65.00
<u>RECERTIFICATION OF DATA LOGGER</u>		
RELOGGER 1	RECERTIFICATION 1 POINT	\$86.50
RELOGGER 2	RECERTIFICATION 2 POINTS	\$104.00
RELOGGER 3	RECERTIFICATION 3 POINTS	\$126.75
RELOGGER 4	RECERTIFICATION 4 POINTS	\$147.25
RELOGGER 5	RECERTIFICATION 5 POINTS	\$213.25
<u>RECERTIFICATION OF TACHOMETERS</u>		
RECAL	RECERTIFICATION OF 6 RPM	\$173.00

Cat#	
RECERT 1	RECERTIFICATION 1 TEMPERATURE POINT
RECERT 2	RECERTIFICATION 2 to 4 TEMPERATURE POINTS
RECERT 3	RECERTIFICATION 5 to 9 TEMPERATURE POINTS
RECERT 4	RECERTIFICATION 10 to 15 TEMPERATURE POINTS
RECERT B	RECERTIFICATION 1 TEMPERATURE POINT BOTTLES and DIGITALS (fully accredited cert; default temp)
RECERT DIG 1	RECERTIFICATION 1 TEMPERATURE POINTS
RECERT DIG 2	RECERTIFICATION 2 TEMPERATURE POINTS
RECERT DIG 3	RECERTIFICATION 3 TEMPERATURE POINTS
RECERT DIG 4	RECERTIFICATION 4 TEMPERATURE POINTS
RECERT DIG 5	RECERTIFICATION 5 TEMPERATURE POINTS
RECERT DIG 6	RECERTIFICATION 6 to 8 TEMPERATURE POINTS
RECERT DIG 1/2	RECERTIFICATION 1 TEMPERATURE POINTS
RECERT DIG 2/2	RECERTIFICATION 2 TEMPERATURE POINTS
RECERT DIG 3/2	RECERTIFICATION 3 TEMPERATURE POINTS
RECERT DIG 4/2	RECERTIFICATION 4 TEMPERATURE POINTS
RECERT DIG 5/2	RECERTIFICATION 5 TEMPERATURE POINTS
RECERT DIG 6/2	RECERTIFICATION 6 to 8 TEMPERATURE POINTS
RECERT IR 1	RECERTIFICATION 1 TEMPERATURE POINT
RECERT IR 2-4	RECERTIFICATION 2 to 4 TEMPERATURE POINTS
REHYD 1	RECERTIFICATION 1 POINT
REHYD 2	RECERTIFICATION 2 POINTS
REHYD 3	RECERTIFICATION 3 POINTS
RERh 1	RECERTIFICATION Rh 1 POINT
RERh 3	RECERTIFICATION Rh 3 POINT
RELOGGER 1	RECERTIFICATION 1 POINT
RELOGGER 2	RECERTIFICATION 2 POINTS
RELOGGER 3	RECERTIFICATION 3 POINTS
ACC052F	GAS CALORIMETER, INLET: DELETE/OBSOLETE
UCY600H	change list price from \$37.57 to \$63.65
ACCG500F-5	+50/500F x 5.0F / DELETE / OBSOLETE
ACCD700P	Pt100 High Precision Thermometer, Single Probe, -200/850C & F, Resolution 0.1C
ACCD700P	Pt100 & Thermocouple Thermometer, Single Channel PT100 -200/850C & F and Thermocouple -200/1760C & F, Resolution 0.1C
ACCD1001PT	Pt100 Immersion Probe, -350/350C, 150mm x 23.0mm
ACC112C	4/6C X 0.02C, Hg / DELETE / OBSOLETE
ACC112CBLS	4/6C X 0.02C, Spirit DELETE / OBSOLETE
TT8020TH	Indoor/Outdoor Hygrometer/Thermometer/Clock, -50/70C & -58/158F, 10/99%Rh
ACCD0330TH	Hygrometer/Thermometer/Dew Point (DISCONTINUED / REPLACE w/ CT310TH)
ACC549P	0/260C x 2.0C, 165mm length
ACC549A	0/260C x 2.0C, 165mm length
ACC549	0/260C x 2.0C, 140mm length
ACC370SUP	Deluxe Waterproof Pocket Thermometer, -50/300C & -58/572F, Cal @ 0.2C & 0.5C
ACC082FBLS	0/220F x 2.0F; discontinued no spirit product available
ACC082FBLSFC	0/220F x 2.0F; discontinued; no spirit product available
ACC082CBLS	-15/105C x 1.0C; discontinued; no spirit product available
ACC082CBLSFC	-15/105C x 1.0C; discontinued; no spirit product available
ACC086FBLS	+200/350F x 2.0F; discontinued; no spirit product available
ACC086FBLSFC	+200/350F x 2.0F; discontinued; no spirit product available
ACC086CBLS	+95/175C x 1.0C; discontinued; no spirit product available
ACC086CBLSFC	+95/175C x 1.0C; discontinued; no spirit product available
ACC600SUP	DISCONTINUED / REPLACE w/ ACC370SUP
ACC655SUP	DISCONTINUED / REPLACE w/ ACC370SUP
ACC310MM	Digital Lab Thermometer, -50/200C, 8" Stem w Min-Max Recall (CHANGE TO ACC310DMM)
ACC1007BSCERT	-1/201C x 0.2C, *Blue Spirit, 630mm

Cat#	
ACC10073BSCERT	-1/201C x 0.2C, *Blue Spirit, 630mm
ACC1005BSCERT	-1/101C x 0.1C, *Blue Spirit, 630mm
ACC10053BSCERT	-1/101C x 0.1C, *Blue Spirit, 630mm
ACCD650PD12C	PMO High Precision Digital, -200/850 C&F w/ 12" probe, Calibrated @ 0, 63, 72, 83, 89C
ACCD650PD12F	PMO High Precision Digital, -328/1760 C&F w/ 12" probe, Calibrated @ 32, 145, 161.6, 180, 191F
ACCD650PD6C	PMO High Precision Digital, -200/850 C&F w/ 6" probe, Calibrated @ 0, 63, 72, 83, 89C
ACCD650PD6F	PMO High Precision Digital, -328/1760 C&F w/ 6" probe, Calibrated @ 32, 145, 161.6, 180, 191F
ACCD0070	USB cable and power supply
BERFREEZTAG2L	Vaccine Data Logger w/ Vial & Probe, Alarm Setting -15.0C w/112 Day
ACCD4000P	Pt100 Digital 1-Channel, (Pt100-4 wire), -200/850C, USB, w/o Probe
ACCD4005P	Pt100 Digital 2-Channel, (Pt100-4 wire), -200/850C, USB, w/o Probe
ACCD4010K	1 Channel (Type K), -200/1760C, USB, w/o Probe
ACCD4015K	2 Channel (Type K), -200/1760C, USB, w/o Probe
ACC8538PC	0.920/1.080 x 0.002 (change part no to ACC8539PC)
ACC713TBLSSFC	+34.5/46C, x 0.1C, 76mm Imm, w/ Safety Coat, Cert@ 37, 44.5C
87011	27" Leatherette, Felt-Lined Thermometer Case
BERFRIDGETAG3	Vaccine Data Logger w/ Vial & 1.5m Probe, Wireless real-time SMS msg certified for FRIDGE (4C) & FREEZER (-20C), 2 yr NIST Calibration
FRIDGETAG2USA	Vaccine Data Logger w/ Vial & Probe, User Sets Alarm Setting 2.0C to 8.0C w/112 Day, Memory & 1.5m probe cable; 2 yr NIST Calibration
ACC038CFC	ASTM 38C, certified 25, 40, 55, 75C
LTTRED307	VFC Approved Temperature Data Logger, 30 Day Summary, W/Probe & Vial (does not include Docking Station LTIUSB); DISCONTINUED / REPLACE w/ LTTRED3016R
LTTRED307USB	VFC Approved Temperature Data Logger, 30 Day Summary, W/Probe & Vial and Docking Station DISCONTINUED / REPLACE w/ LTTRED3016RUSB
LTTRED3016R	VFC Approved Temperature Data Logger, 30 Day Summary, W/Probe & Vial (does not include Docking Station LTIUSB)
LTTRED3016RUSB	VFC Approved Temperature Data Logger, 30 Day Summary, W/Probe & Vial and Docking Station LTIUSB
SW0013FC	Slide Warmer Thermometer, Certified at 40C & 54C
	MID YEAR REVISION PUBLISHED AUG 2017
ACC099CFC	certified -46, -32, -18, 0C
ACC099FFC	certified -50, -25, 0, 32F
ACC0610CS	0/220F x 2.0F, Tol imm, Red Spirit
ACCD1055	Pt100 Air Probe, 250mm x 4mm, -50/250C
ACCD0072	Probe Connector for 4000, 700, 795 Series
ACCHT5063	Analog Thermo-Hygrometer, Certified
HIR90TTW	2 channel Wireless/WiFi Data Logger, Digital Readout, Memory 65,000
BERCLMdocLR	Multi-Use Data Logger, w/ Alarm, -25/55C
LASGFX2V	GFX Dual Probe Data Logger -34/+40C & -34/+104F
LASPROBE	Thermocouple Probe w/SS Tip & Glycol Vial
LASNIST	1 Temperature Point NIST-ISO17025/ILAC Certification
ACCD0552DIG	"Flip-Probe" Digital Pocket Thermometer, -50/350C & -58/662F, 4-1/2" Probe
ACC6267CBLSFC	ASTM 0S62C to 0S67C Tol IMMERSION, -38 to 155C w/ case
ACC6267FBLSFC	ASTM 0S62F to 0S67F Tol IMMERSION, -36 to 311F w/ case
ACC074FBLS	discontinued; no spirit product available
ACC134CBLS	144/156C x 0.2C
ACC134CBLSFC	144/156C x 0.2C, CERT 145, 150, 155C
ACCTIFC300	Windows Software/jump drive and cradle for ACCCRYO86
ACCRT8017	Monitoring and logging up to 24 million readings w/ 2 channel & WiFi capability
Flash Alert	WiFi Card for ACCRT8017
ACCTIFC300	Windows Software/jump drive and cradle for ACCCRYO86
CPHYD1	1 CERTIFIED POINT (remove spaces in part no.)
CPHYD2	2 CERTIFIED POINTS (remove spaces in part no.)

Cat#	
CPTHYD3	3 CERTIFIED POINTS (remove spaces in part no.)
ACCD600PEX	Pt100 High Precision Thermometer, Single Probe, -200/850C & F, Resolution 0.1C Intrinsically Safe/Explosion Proof; DISCONTINUED / REPLACE w/ ACCD700PEX
ACCD650P	Pt100 High Precision Thermometer, Single Probe, -200/850C & F, Resolution 0.01C DISCONTINUED / REPLACE w/
ACCD655P	Pt100 High Precision Thermometer, Dual Probe, -200/850C & F, Resolution 0.01C DISCONTINUED / REPLACE w/
ACCD650PEX	Pt100 High Precision Thermometer, Single Probe, -200/850C & F, Resolution 0.01C Intrinsically Safe/Explosion Proof DISCONTINUED / REPLACE w/ ACCD750PEX
ACCD655PEX	Pt100 High Precision Thermometer, Dual Probe, -200/850C & F, Resolution 0.01C Intrinsically Safe/Explosion Proof DISCONTINUED / REPLACE w/ ACCD755PEX
ACC410PW	Digital Ultra Low Freezer Thermometer, -100/200C x 0.1C
B260CW3BLS	change length from 12" to 16"
ACC605P	Pt100 Precision Thermometer, Dual Probe, -200/850C & F, DISCONTINUED / REPLACE w/ ACCD705P
ACCD610P	Thermomcouple Thermometer, Single Channel, -200C to 1760C & F, DISCONTINUED / REPLACE w/ ACCD710P
ACCD615P	Thermomcouple Thermometer, Dual Channel, -200C to 1760C & F, DISCONTINUED / REPLACE w/ ACCD715P
ACCD610PEX	Thermomcouple Thermometer, Single Channel, -200C to 1760C & F, DISCONTINUED / REPLACE w/ ACCD710PEX
ACC0486IR	CHANGE PART NO. TO ACCD0486IR
ACC647SC	Dry Bath/Block Heater, 0/110C x 0.5C, 35mm Imm, Safety Coated (CHANGE PART NO. TO ACC6471SC)
ACC646SC	Short Immersion, 12mm Imm, 75mm long, Safety Coated (CHANGE PART NO. TO ACC645SC)
GW3000	change discription of thermometer from spirit to Hg
GW3001	change discription of thermometer from spirit to Hg
GW3002	change discription of thermometer from spirit to Hg
GW3003	change discription of thermometer from spirit to Hg
GW5000	change discription of thermometer from spirit to Hg
GW5001	change discription of thermometer from spirit to Hg
GW5002	change discription of thermometer from spirit to Hg
GW5003	change discription of thermometer from spirit to Hg
GW5004	change discription of thermometer from spirit to Hg
GW5005	change discription of thermometer from spirit to Hg
GW5006	change discription of thermometer from spirit to Hg
GW5007	change discription of thermometer from spirit to Hg
GW5008	change discription of thermometer from spirit to Hg
GW5009	change discription of thermometer from spirit to Hg
GW2500C	change discription of thermometer from spirit to Hg
GW2501C	change discription of thermometer from spirit to Hg
GW2502C	change discription of thermometer from spirit to Hg
GW2503C	change discription of thermometer from spirit to Hg
GW2504C	change discription of thermometer from spirit to Hg
GW2505C	change discription of thermometer from spirit to Hg
GW2506C	change discription of thermometer from spirit to Hg
GW2507C	change discription of thermometer from spirit to Hg
GW2508C	change discription of thermometer from spirit to Hg
GW2509C	change discription of thermometer from spirit to Hg
GW2509AC	change discription of thermometer from spirit to Hg
GW2532C	change description, define thermometer as spirit
GW2533C	change description, define thermometer as spirit
GW2534C	change description, define thermometer as spirit
GW2535C	change description, define thermometer as spirit
GW2540C	change description, define thermometer as spirit
GW2541C	change description, define thermometer as spirit
GW2542C	change description, define thermometer as spirit
GW2543C	change description, define thermometer as spirit
GW2544C	change description, define thermometer as spirit
GW2545C	change description, define thermometer as spirit
GW2546C	change description, define thermometer as spirit
GW2547C	change description, define thermometer as spirit
ACCD1079	Pt100 Immersion Probe, 300mm x 6mm, -200/650C & -330/1202F (DISC/REPLACE w/ ACCD1092)

Cat#	
ACCD1092	Pt100 Immersion Probe, 400mm x 6mm, -200/660C & -330/1220F, w/out handle, w/ 1500mm cable
Flash Alert	WiFi Card for ACCRT8017 / (DISCONTINUE/REPLACE w/ ACCFA1)
ACCFA1	Flash Alert WiFi Card for ACCRT8017
ACCD400PW	IP 65 Water Proof Pt100 Platinum Digital Thermometer w/o Probes (CHANGE PART NO TO ACCD4000PW)
ACC10073BLSFC	-1/201C, x 0.2C, 76mm Imm, Cert@ 0, 25, 37, 44.6, 56, 105, 121, 180, 200C change div & cert temp from 44.6 to 44.5
ACC0623FC	-38/+2C x 0.1, 76mm Imm, Cert@ -24, -18, -12, -6, 0C (REPLACE w/ ACC133CFC)
ACC6558DIGSUP	Waterproof Pocket Digital Thermometer, -58/572F & -50/300C, 8" Probe, Cal @ 0.5C (CHANGE PART NO. TO ACC
ACCD755LOG	Pt100 High Precision Thermometer & data logger, Dual Probe, -200/850C & F, Resolution 0.01C
FRIDGETAG2USA	Vaccine Data Logger w/ Vial & Probe, User Sets Alarm Setting -30°C to 50°C, in 0.1°C increments 1.5m probe cable; 3 yr NIST Calibration (specs changed)
ACC067FBLSFC	0S67F, +206/311F, Tol IMM, 400mm, BLUE SPIRIT, Cert@ 32,205,240,275,310F (specs changed)
TFA144002	Freezer/Fridge Thermometer, -30/20C, 6", Hook Top (DISCONTINUED)
TFA144005	Round Freezer/Fridge Thermometer, -40/30C, 2" Dial, Hook Top (DISCONTINUED)
ACCRTAMB	Ambient Room Temperature Sensor Probe forr ACCRT8017
ACCFA1	Flash Alert WiFi Card for ACCRT8017 (RENAME part)
ACC067FBLS	0S67F, +203/311F, BLUE SPIRIT (specs changed from 206/311F to 203/311F)
ACC048FBLS	+177.5/182.5F x 0.41F (specs changed from 0.41F div to 0.1F div)
ACC048FBLSFC	+177.5/182.5F x 0.41F (specs changed from 0.41F div to 0.1F div)
ACC524A	+20/400F x 2.0F, 16" ARMOR (change to 14" armor)
PBRT8002-STprobe	Straight probe for ACCRT8002 & ACCRT8017 (accessory)
BERFRIDGETAG2V	REPLACED BY BERFRIDGETAG2L
BERFREEZTAG2V	REPLACED BY BERFREEZTAG2L
ACC400	13" Galileo Thermometer (DISCONTINUE)
BRIX COMBINED FOR	BRIX combined form length change from 300mm to 375mm
ACC713TBLSSFC	+34.5/46C, x 0.1C, Toal Imm, w/ Safety Coat, Cert@ 37, 44.5C (change from 76mm imm)
ACC097FBLSFC	0/120F x 1.0F (changed div from 0.5F)
ACCCHITEMP135-5	Autoclve Data Logger, -40 to 140C x 0.05C, w/5.25" probe, Accuracy +/- 0.5C change from 150C
ACCCHITEMP140-5	Autoclve Data Logger, -40 to 140C x 0.05C, w/5.25" probe, Accuracy +/- 0.5C change from 150C
ACCD650PD12C	PMO High Precision Digital, -200/850 C&F w/ 12" probe, Calibrated @ 0, 63, 72, 83, 89C (change to 750P)
ACCD650PD12F	PMO High Precision Digital, -328/1760 C&F w/ 12" probe, Calibrated @ 32, 145, 161.6, 180, 191F (change to 750P)
ACCD650PD6C	PMO High Precision Digital, -200/850 C&F w/ 6" probe, Calibrated @ 0, 63, 72, 83, 89C (change to 750P)
ACCD650PD6F	PMO High Precision Digital, -328/1760 C&F w/ 6" probe, Calibrated @ 32, 145, 161.6, 180, 191F (change to 750P)
B260C500FW3S	-10 to 260C X 1.0c - Replace w/ACC260C500FW3BS
ACC310SUP	Digital Lab Thermometer, -58/536F & -50/280C, 8" Stem , Cal @ 0.2C & 0.5C accuracy
ACC310DLSSUP	Digital Lab Thermometer, Long Probe, -58/536F & -50/280C, 12" Stem, Cal @ 0.5C accuracy
ACC3105SUP	Digital Lab Thermometer, -58/302F & -50/150C, 5" Stem, Cal @ 0.2C & 0.5C accuracy
ACC3108SUP	Digital Lab Thermometer, -58/302F & -50/150C, 8" Stem, Cal @ 0.2C & 0.5C accuracy
ACC310MMSUP	Digital Lab Thermometer, -58/392F & -50/200C, 8" Stem w Min-Max, Cal @ 0.2C & 0.5C accuracy
ACC370SUP	Deluxe Waterproof Pocket Thermometer, -50/300C & -58/572F, Cal @ 0.2C & 0.5C accuracy
ACC400SUP	Waterproof/Shockproof Digital Thermometer, -50/300C & -58/572F, 5" Stem, Cal @ 0.5C accuracy
ACC420SUP	Waterproof/Shockproof Digital Thermometer, -50/300C & -58/572F, 5" Stem, Cal @ 0.5C accuracy
ACC600SUP	Waterproof Pocket Digital Thermometer, -58/302F & -50/150C, 4 1/2" Probe, CAL@ 0.5C accuracy
ACC655SUP	Waterproof Pocket Digital Thermometer, -58/572F & -50/300C, 4 1/2" Probe, Cal @ 0.5C accuracy
ACC6558SUP	Waterproof Pocket Digital Thermometer, -58/572F & -50/300C, 8" Probe, Cal @ 0.5C accuracy
ACC610SUP	Waterproof, T-Handle Thermometer, -58/302F & -50/150C, 4 1/2" Probe, Cal @ 0.5C accuracy
B110CWBLK	discontinued; no black spirit available
B150CWBLK	discontinued; no black spirit available
B120FWBLK	discontinued; no black spirit available
ACC006C/ACC006C	note product is non-mercury, *Toluene or other suitable liquid colored red w/ a permanent dye shall be used as tl
ACC006F/ACC006F	note product is non-mercury, *Toluene or other suitable liquid colored red w/ a permanent dye shall be used as tl

Cat#		
ACC114C/ACC114C	note product is non-mercury, *Toluene or other suitable liquid colored red w/ a permanent dye shall be used as tl	
LTTREL3016R	Ultra-Low Temperature Data Logger, -86 to +35C x 0.1C; does not include Docking Station LTIUSB	
LTTREL3016RUSB	Ultra-Low Temperature Data Logger, -86 to +35C x 0.1C includes Docking Station LTIUSB	
ACC9604	+50/450F x5.0F; discontinued no replacement	
ACCD0330TH	Thermo-Hygrometer, Dew Point, & Wet Bulb Temp, -40/70C Accuracy 0.5C, 1/99RH Accuracy +/-3% (I	
ACCD1044WS	WIRELESS WEATHER STATION, TEMP, HUMIDITY, BAROMETER & TIME (discontinue/no replacement)	
ACCD033147	TRANSMITTER FOR ACCD1044WS (discontinue/no replacement)	
ACCD0412IR	REPLACED w/ ACCD0414IR	
ACC0530IR	PRO-SCAN 5230, HIGH ACCURACY ±0.75C, -35 to 900C & THERMOCOUPLE -64 to 1370C & F	
ACC043C	KINEMATIC VISCOSITY; *contains mercury thalium; no longer produced DISCONTINUED	
ACC043F	KINEMATIC VISCOSITY; *contains mercury thalium; no longer produced DISCONTINUED	
ACC074C	KINEMATIC VISCOSITY; *contains mercury thalium; no longer produced DISCONTINUED	
ACC074F	KINEMATIC VISCOSITY; *contains mercury thalium; no longer produced DISCONTINUED	
ACC076F	ANTIFREEZE FREEZING POINT; *contains mercury thalium; no longer produced DISCONTINUED	
ACC099F	WEATHERING TEST; *contains mercury thalium; no longer produced DISCONTINUED	
ACC099C	WEATHERING TEST; *contains mercury thalium; no longer produced DISCONTINUED	
ACC122C	BROOKFIELD VISCOSITY; *contains mercury thalium; no longer produced DISCONTINUED	
87007	16" Leatherette, Felt-Lined Thermometer Case; DISCONTINUED	\$41.22
ACCTH2000	DISCONTINUE	
ACCD790P	Pt100 High Precision Thermometer, Dual Probe, -200/1760C & F X 0.001, thermocouple & humidity w/ automatic probe adjustment and detection, bluetooth	
LTTREL3016	Ultra Low Cryotemp data logger; does not include Docking Station LTIUSB	
LTTREL3016USB	Ultra Low Cryotemp data logger; includes Docking Station LTIUSB	
LTIWIFI	Wi-Fi Interface Docking Station for LogTag@Dataloggers	
ACC522SA	CHANGE FROM 30 TO 0/300F x 2.0F, 12" ARMOR	
ACC526SA	30/500F x 2.0F, 14" ARMOR	
ACC091CBL5FC	ASTM S91C, 20/50 x 0.1C, 76mm Imm, CERT @ 20, 30, 40, 50C	
NR-100-10	Non-Roll Devise for Thermometers, 10/Pack (DISCONTINUE); SOLD AS 25/PACK ONLY	
CS150C300FWS	-20/110C x 1.0C & 0/300F x 2.0F, Tol immersion, red spirit, 12" (305mm)	
CS150C300FW3S	-20/110C x 1.0C & 0/300F x 2.0F, partial immersion, red spirit, 12" (305mm)	
ACCST1220150S	-10/150C x 1.0C, 20mm Imm; DISCONTINUED NO LONGER CARRYING 7/12 GROUND JT	
ACCST1240150S	-10/150C x 1.0C, 40mm Imm; DISCONTINUED NO LONGER CARRYING 7/12 GROUND JT	
ACCST1220250S	-10/250C x 1.0C, 20mm Imm; DISCONTINUED NO LONGER CARRYING 7/12 GROUND JT	
ACCST1240250S	-10/250C x 1.0C, 40mm Imm; DISCONTINUED NO LONGER CARRYING 7/12 GROUND JT	
ACC960K2	DUAL CHANNEL TYPE K THERMOCOUPLE THERMOMETER; DISCONTINUE replacement ACCD4015K	
ACCD0485IR	SCAN-TEMP 485,-50 to 800C & F, ACCURACY ±2.0C, RES. 20:1, DOUBLE LAZER SIGHTING MAX, HOLD, LOCK FUNCTIONS, W/CASE (previously ACCD485IR)	
GW2575X	0/11 x 0.1, baume, 150mm	
GW2576X	9/21 x 0.1, baume, 150mm	
PT310F	0 to 300F x 2.0F (DUPLICATED; WRONG PRICE)	
ACCHT5063	Analog Thermo-Hygrometer, 6" brass, glass, Certified (insert face type)	
ACCHT5063A	Analog Thermo-Hygrometer, economy, Certified (insert face type)	
ACCD4010K	1 Channel (Type K), -200/1760C, USB, *change from w/o probe to w/ probe	
ACCD4015K	2 Channel (Type K), -200/1760C, USB, w/ Probe *change from w/o probe to w/ probe	
ACC800MAEFV	VACCINE FREEZER THERMOMETER, -50 to 70C x 0.1C, Certified at -20C (DISCONTINUE)	
ACC800MAERV	VACCINE REFRIGERATOR THERMOMETER, -50 to 70C x 0.1C, Certified at +4C (DISCONTINUE)	
ACC800MAEIV	VACCINE INCUBATOR THERMOMETER, -50 to 70C x 0.1C, Certified at +37C (DISCONTINUE)	
ACC800MAE	THERMOMETER & Probe Only (DISCONTINUE)	

Cat#	
RELOGGER 4	RECERTIFICATION 4 POINTS
RELOGGER 5	RECERTIFICATION 5 POINTS
S1141741SCA	-100/200C x 1.0C, 300mm Length, Safety Coated (change range previously 0/200C)
LTUTRED30WIFI	VFC Approved Temperature Data Logger, 30 Day Summary, W/3M Probe & Vial, Real time WiFi enabled & wall mount bracket LTBRACKETWIFI
LTBRACKETWIFI	Wall Mount Bracket for UTRED3016 & UTRED30WIFI
ACC8010V	For Ovens, -50 to 200C, Certified at +100C (previously +80C)
LTIUSB	USB Interface Docking Station for LogTag@Dataloggers (DISCONTINUED)
LTIHID	Desktop USB Interface Docking Station for LogTag@Dataloggers (replaces LTIUSB)
ACCOV080S	Oven, +35/200C x 1.0C; NEW RANGE Oven, +35/180C x 1.0C
ACCD0069	K Thermomcouple beaded wire probe, teflon insulated, -40 to 250C & F (DISCONTINUED)
ACCD0062	K Thermomcouple beaded wire probe, teflon insulated, -40 to 250C & F (REPLACEMENT FOR ACCD0069)
B400CY3	change specifications from -10/400 x 2.0C, 12" to -10/420 x 1.0C, 13"
ACCF100S	Freezer, -25/+5C x 0.1C; NEW RANGE Freezer, -25/-5C x 0.1C
ACC6113XXSFC	-20/110C x 0.5C, 76mm Imm, Spirit, CERT @ 0, 37, 56C
ACC895WRFT	For Refrigerators/Freezers/Ambient, -50/70C, Certified at +4C & -20C
ACC895WRFTCOV	Accessory - protective plastic cover for ACC895WRFT
CPT RH	Certification of 1 Rh point
RECAL	RECERTIFICATION OF 6 RPM
PT312F	0 to 250F x 2.0F - DISCONTINUED
PT314F	+50 to 750F x 10.0F - DISCONTINUED
PT312C	-20 to 120C x 1.0C - DISCONTINUED
PT311CCOMP	+20 to 180C x 2.0C - DISCONTINUED
PT309F	-50 to 250F x 2.0F (DISCONTINUED)
PT311F	+70 to 370F x 2.0f (DISCONTINUED)
ACC073CFC	DISCONTINUE - No Hg Thallium
THDI800	Benchtop dy ice maker; 1lb. Blocks; simple lever operation
ACC304MCF1	Waterproof & NSF approved , NIST Compliant, 1½" Dial, 2½ Stem, -50/200C & -58/392F x 0.1, Min-Max
LTI-WM-WIFI	Wall mount Wi-Fi Interface Cradle LogTag@Dataloggers (retrofits any LogTag@ datalogger w/ probe)
ACC031F	2020 PRICE INCREASE
ACC089C	2020 PRICE INCREASE
ACC096C	2020 PRICE INCREASE
ACC10073S	2020 PRICE INCREASE
ACC1025012	2020 PRICE INCREASE
ACC1025012SC	2020 PRICE INCREASE
ACC10250618S	2020 PRICE INCREASE
ACC10250618SSC	2020 PRICE INCREASE
ACC107C	2020 PRICE INCREASE
ACC5050A	2020 PRICE INCREASE
ACC519SA	2020 PRICE INCREASE
ACC525A	2020 PRICE INCREASE
ACC527A	2020 PRICE INCREASE
ACC6103XX	2020 PRICE INCREASE
ACC6103XXS	2020 PRICE INCREASE
ACC610XX	2020 PRICE INCREASE
ACC703	2020 PRICE INCREASE
ACCCRPENR	2020 PRICE INCREASE
ACCHR5	2020 PRICE INCREASE
ACCHR7	2020 PRICE INCREASE
ACCST11250S	2020 PRICE INCREASE
GW2601	2020 PRICE INCREASE
GW2602	2020 PRICE INCREASE
GW2603	2020 PRICE INCREASE

Cat#	
GW2604	2020 PRICE INCREASE
GW2605	2020 PRICE INCREASE
GW2606	2020 PRICE INCREASE
GW2607	2020 PRICE INCREASE
GW2608	2020 PRICE INCREASE
GW2609	2020 PRICE INCREASE
GW2610	2020 PRICE INCREASE
GW2611	2020 PRICE INCREASE
GW2612	2020 PRICE INCREASE
GW2613	2020 PRICE INCREASE
GW2614	2020 PRICE INCREASE
GW2615	2020 PRICE INCREASE
GW2616	2020 PRICE INCREASE
ACC605	2020 PRICE INCREASE
ACC1006S	+30/394F x 0.5F * 2020 PRICE INCREASE
ACC3108DIG	Digital Lab Thermometer, -58/392F & -50/200C, 8" Stem (CHANGE SPECS FROM -58/302F & -50/150C)
R164978	BEER/WORT, -2/6% X 0.1, 270MM length, 110ml of sample needed
R164980	BEER/WORT, 0/8% X 0.1, 270MM length, 110ml of sample needed
R165000	BEER/WORT, 8/16% X 0.1, 270MM length, 110ml of sample needed
R165020	BEER/WORT, 16/24% X 0.1, 270MM length, 110ml of sample needed
R165030	BEER/WORT, 0/8% X 0.1, 330MM length, 240ml of sample needed
R165032	BEER/WORT, 8/16% X 0.1, 330MM length, 240ml of sample needed
R165034	BEER/WORT, 16/24% X 0.1, 330MM length, 240ml of sample needed
R165040	BEER/WORT, 0/14% X 0.1, 406MM length, 430ml of sample needed
R165060	BEER/WORT, 12/26% X 0.1, 406MM length, 430ml of sample needed
R165160	BEER/WORT, 2/8% X 0.1, 406MM length, 430ml of sample needed
R165161	BEER/WORT, 8/14% X 0.1, 406MM length, 430ml of sample needed
R165162	BEER/WORT, 10/16% X 0.1, 406MM length, 430ml of sample needed
R165168	BEER/WORT, 0/10% X 0.1, 432MM length, 430ml of sample needed
R165170	BEER/WORT, 10/20% X 0.1, 432MM length, 430ml of sample needed
R165078	BEER/WORT, -1/2% X 0.05, 356MM length, 370ml of sample needed
R165082	BEER/WORT, 2/5% X 0.05, 356MM length, 370ml of sample needed
R165089	BEER/WORT, 6/9% X 0.05, 356MM length, 370ml of sample needed
R165092	BEER/WORT, 9/12% X 0.05, 356MM length, 370ml of sample needed
R165096	BEER/WORT, 11/14% X 0.05, 356MM length, 370ml of sample needed
R165097	BEER/WORT, 12/15% X 0.05, 356MM length, 370ml of sample needed
R165100	BEER/WORT, 14/17% X 0.05, 356MM length, 370ml of sample needed
B50CWBL5	-35/50C x 1.0C (multiple listings in different tabs)
B330CW3BRN	-10/330C x 2.0C (multiple listings in different tabs)
ACC120CBLS	+38.6/41.4C x 0.05C (multiple listings in different tabs)
ACC022FBLS	+204/218F x 0.2F (multiple listings in different tabs)
ACC022CBLS	+95/103C x 0.1C (multiple listings in different tabs)
ACC015FBLS	+30/180F x 0.5F
ACC015CBLS	-2/80C x 0.2C
ACC801DIG	Digital Thermometer w/Probe, -50/70C & F, Min-Max, Alarm (change range from 200C to 70C)
ACCG500F-C	+50/500F x 5.0F % 20/260C x 2.0C
B50CWBLSSC	-30/50C X 1.0C w/ SAFTEY COAT (multiple listings in different tabs)
ACCG250F	+50/250F x 2.0F (DISCONTINUE; replace w/ ACCG300F)
BERFZTAG60	Monitoring Frost-Sensitive Materials, 60 Minutes Below 0C (DISCONTINUED)
ACC10613SFC	-1/61C x 0.1C, 76mm Imm, Cert@ 0, 37, 56C (Change description from 1/61C to -1/61C)
ACCD0478IR	Medi-Scan Forehead & Surface InfraRed Thermometer, 32/42.5C & 93.2/108.5F

Cat#		
ACCD04900IR	AccuTherm Care4U Forehead & Surface InfraRed Thermometer, 32/42.5C & 89.6/108.5F	
ACCD0973IR	Wall Mount Fever Alarm Non-Contact Forehead InfraRed Thermometer, 0/50C & 32/122F	
ACCD4000P	Pt100 Platinum Digital Thermometer w/o Probes (change from ACCD400P)	
ACCD0973IR	Wall Mount Fever Alarm Non-Contact Forehead InfraRed Thermometer, 0/50C & 32/122F (PRICE CHAN	
REVFCLOG 1	RECERTIFICATION OF 1 POINT VFC DATALOGGER (NEW CATEGORY)	
REVFCLOG 2	RECERTIFICATION OF 1 POINT VFC DATALOGGER (NEW CATEGORY)	
ACC6103XXSFC	0/230F x 0.5F, 76mm, cert @ 32, 45, 161F (new part number)	\$174.83
ACCCR87C	Temperature Circular Chart/Digital Recorder, Celsius, 120V (DISCONTINUED)	\$544.57
ACCCR87F	Temperature Circular Chart/Digital Recorder, Fahrenheit, 120V (DISCONTINUED)	\$544.57
ACCCR87C220	Temperature Circular Chart/Digital Recorder, Celsius, 220V (DISCONTINUED)	\$576.82
ACCCR87F220	Temperature Circular Chart/Digital Recorder, Fahrenheit, 220V (DISCONTINUED)	\$576.82
ACCC800CH	Temperature/Humidity Circular Chart/Digital Recorder, 115V (DISCONTINUED)	\$888.31
ACCCR8716	Range -40 to 0C, Time 7 Days (DISCONTINUED)	\$50.28
ACCCR8717	Range -30 to 10C, Time 7 Days (DISCONTINUED)	\$50.28
ACCCR8718	Range +10 to 50C, Time 7 Days (DISCONTINUED)	\$50.28
ACCCR8722	Range -40 to 0C, Time 31 Days (DISCONTINUED)	\$50.28
ACCCR8723	Range -30 to 10C, Time 31 Days (DISCONTINUED)	\$50.28
ACCCR8724	Range +10 to 50C, Time 31 Days (DISCONTINUED)	\$50.28
E3050C7	RANGE -30 to 50C & 0 to 100%Rh, TIME 7 Days (DISCONTINUED)	\$60.04
E3010C7	RANGE -30 to 10C & 0 to 100%Rh, TIME 7 Days (DISCONTINUED)	\$60.04
E545C7	RANGE +5 to 45C & 0 to 100%Rh, TIME 7 Days (DISCONTINUED)	\$60.04
E3050C31	RANGE -30 to 50C & 0 to 100%Rh, TIME 31 Days (DISCONTINUED)	\$60.04
E3010C31	RANGE -30 to 10C & 0 to 100%Rh, TIME 31 Days (DISCONTINUED)	\$60.04
ACCCRPEN	PEN, BLACK (DISCONTINUED)	\$18.48
ACCCRPENB	PEN, BLUE (DISCONTINUED)	\$18.48
ACCCRPENR	PEN, RED (DISCONTINUED)	\$21.04
ACCRHTEMP101	Dual Channel Temperature / Humidity Logger, -40 to +80C&F and 0 to 95% Rh, 21K readings (DISCONT	
PT313F	0 to 500F x 5.0F (DISCONTINUED)	
PT315F	0 to 150F x 1.0F (DISCONTINUED)	
PT313C	-20 to 250C x 2.0C (DISCONTINUED)	
PT314C	+10 to 400C x 5.0C (DISCONTINUED)	
PT315C	-15 to 65C x 1.0C (DISCONTINUED)	
PT571F	+50 to 320F x 1.0F (1-3/4" dia.) (DISCONTINUED)	
PT310F	0 to 300F x 2.0F (DISCONTINUED)	
PT309C	-50 to 110C x 2.0C (DISCONTINUED)	
PT310C	-10 to 150C x 2.0C (DISCONTINUED)	
PT311C	+20 to 180C x 2.0C (DISCONTINUED)	
PT309FCOMP	-50 to 250F x 2.0F (DISCONTINUED)	
PT310FCOMP	0 to 300F x 2.0F (DISCONTINUED)	
PT311FCOMP	+70 to 370F x 2.0F (DISCONTINUED)	
PT309CCOMP	-50 to 110C x 2.0C (DISCONTINUED)	
PT310CCOMP	-10 to 150C x 2.0C (DISCONTINUED)	
ACCGG500FC	DISCONTINUED USE (ACCG500F) 10/260C & 50/500F, 1 3/4" Dial, 8" Stem	
ACCG250F	+50/250F x 2.0F (DISCONTINUED I413+I414	
ACC57089	PHYSICS 1000 Pt100 HAND HELD DIGITAL THERMOMETER (DISCONTINUED)	
ACC59954	Pt100 Immersion Probe, -90 to 350C, 250mm x 3mm (DISCONTINUED)	
ACC59966	Pt100 Insertion Probe, -50 to 200C, 250mm x 3mm (DISCONTINUED)	

Cat#	
ACC59270	Pt100 Glass-Type ASTM Immersion Probe, -50 to 310C, 250mm x 6mm (DISCONTINUED)
ACC59265	Pt100 ASTM immersion Probe, -90 to 410C, 400mm x 3mm (DISCONTINUED)
ACC58306	"Physics-View Software (DISCONTINUED)
ACC703	-10/500C x 2.0C (DISCONTINUED)
LTREL8	Multi-Use Low Temperature Data Logger w/ external probe, -80 to +40C; does not include Docking Station
LTTRIL8	Dry Ice Low Temperature Data Logger w/ internal probe, -80 to +40C; does not include Docking Station
LTSRIL8	Dry Ice SINGLE USE Low Temp Data Logger w/ internal probe, -80 to +40C; does not include Docking Station
ACCD0974IR	Wall Mount Fever Alarm Non-Contact Forehead InfraRed Thermometer, 0/50C & 32/122F
ACCD0973IR	Wall Mount Fever Alarm Non-Contact Forehead InfraRed Thermometer, 0/50C & 32/122F
BPS120	Replacement Wicks
BERFRIDGETAG3	Vaccine Data Logger w/ Vial & 1.5m Probe, Wireless real-time SMS msg (DISCONTINUED) certified for FRIDGE (4C) & FREEZER (-20C), 2 yr NIST Calibration <i>*Fridgetag@3 product sales sunset date 12/28/2020 due to 3G network being disbanded</i>
LTSRIL8	Dry Ice SINGLE USE Low Temp Data Logger w/ internal probe, *duplicate entries / match pricing
ACC1119DIG	Pocket Digital Thermometer w/4" Probe, -40/200C & -40/400F (DISCONTINUED)
ACCD952K	TYPE K GENERAL PURPOSE PROBE W/HANDLE *Change part no. from ACC952K to ACCD952K
ACCP160F	-40/160F x 2.0F (DISCONTINUED)
ACCL250C	0/250C x 2.0C (DISCONTINUED)
ACCL160FC	-40/160F x 2.0F & -40/70C x 1.0C (DISCONTINUED)
ACCL160F	-40/160F x 2.0F (DISCONTINUED)
ACCL125FC	+25/125F x 1.0F & 0/50C x 0.5C (DISCONTINUED)
ACCL125F	+25/125F x 1.0F (DISCONTINUED)
ACCL110C	-10/110C x 1.0C (DISCONTINUED)
CAL SVC	FEE FOR TEMPS BETWEEN 300/500C & 572/932F (NEW CATEGORY)
ACC14ARM	Thermometer Armor 14" (DISCONTINUED)
ACC345DIG	Rotary Head Thermometer, -58/302F & -50/150C, 5" Stem (DISCONTINUED)
B2020S	0 to 220F x 2.0F (DISCONTINUED)
GW2553X	0.590/0.660 x 0.001, 150mm (DISCONTINUED)
ACC809DIG	For Refrigerator or Indoor/Outdoor Temperatures -10 to 50C x 0.1C (DISCONTINUED)
ACC809SUP	For Refrigerator or Indoor/Outdoor Temperatures -10 to 50C x 0.1C, CERTIFIED to 0.5C (DISCONTINUED)
ACCD0974IR	Wall Mount Fever Alarm Non-Contact Forehead InfraRed Thermometer, 0/50C & 32/122F (DISCONTINUED)
ACCG40C	-100/40C x 1.0C (DISCONTINUED)
ACC094CBLS	+80/110C x 0.1C
ACC094CBLSFC	+80/110C x 0.1C, CERT @ 80,90, 100, 110C
ACC605	-200/+30C x 1.0C *Blue Spirit, 350mm (update specs)
ACCL100CX	-50/100C x 1.0C (DISCONTINUED / REPLACE w/ ACCG100CX)
ACCL150C	0/150C x 1.0C (DISCONTINUED / REPLACE w/ ACCG150C)
ACCF1	Flash Alert WiFi Card for ACCRT8017 (accessory) (DISCONTINUED)
ACC1119DIG	Pocket Digital Thermometer w/4" Probe, -40/200C & -40/400F (DISCONTINUED)
ACCG220F-C	0/220F x 2.0F & -10/100C x 1.0C (DISCONTINUED)
BERQTAG	Multi-Limit Temperature Indicator Logger, w/Alarm, -25/+50C (DISCONTINUED SPECIAL ORDER)
BERQTAGUSB	Multi-Limit Temperature Indicator Logger, w/Alarm, -25/+50C & USB Interface (DISCONTINUED SPECIAL ORDER)
BERQTAGCRT	Multi-Limit Temperature Indicator Logger, w/Alarm, -25/+50C, Certified Room Temp (DISCONTINUED)
ACCRT8017	Monitoring and logging up to 24 million readings w/ 2 channel & WiFi capability (change description)
ACCRT8002	Monitoring and logging up to 24 million readings (change description)
ACCRT8019UL	Dual buffered probe DDL, ultra low -80 to +70C, 2 yr NIST Calibration -70C, VFC compliant (NEW PRODUCT)
ACCD100LOG	Temperature Data Logger, Digital Readout, Memory 60,000, (DISCONTINUED) -30 to +70C, W/Audio & Visual Alarm
ACCD110LOG	Temperature & Humidity Data Logger, Digital Readout, Memory 60,000, (DISCONTINUED) -30 to +70C & 0 to 99% Rh, W/Audio & Visual Alarm

Cat#		
ACCD110LS	Windows Software for LOG100 & LOG110 (DISCONTINUED)	
ACCD103LOG	External Temperature Probe, 9' (3 meter) Cable (DISCONTINUED)	
ACCD108LOG	External Temperature Probe, 25' (8 meter) Cable (DISCONTINUED)	
ACCD067LOG	Spare Battery (DISCONTINUED)	
ACCD0088	Lockable Wall Mount (DISCONTINUED)	
ACCTEMP101A	DTemperature Recorder, -40 to +80C x 0.1, Accuracy +/- 0.5C, 32,676 Readings (DISCONTINUED)	
ACCTC4000	DDual Channel Thermocouple /Ambient Logger Thermocouple (DISCONTINUED) Types J,K,T,E,S,B,N, 16,383 Readings / Channel	
MGNISTRH	Certificate Traceable to N.I.S.T., Temperature & Humidity	
ACCTIFC200	Windows Software for General Purpose Data Loggers	
ACC895WDIG	CHANGE PART NO. FROM ACC895DIG TO ACC895WDIG	
ACC895ULF	For -80/70C ultra low freezer or refrigerator, air probe, no buffr, Certified at -70C (NEW PRODUCT)	
THID800	price increase from supplier	
B50CWBLSSC	duplicate part number on multiple tabs; corrected pricing to reflect same on all	
<u>PRECISION THERMOCOUPLE THERMOMETER w/NIST TRACEABLE CERTIFICATE AND RS232</u>		
	<ADD> No certificate for individual probe purchase	
ACC340DIG	DRotary Head Thermometer, -58/302F & -50/150C, 8" Stem (DISCONTINUED)	
ACC99050005	Heat-Conductive Paste 20g	
LTIWIFI	DWi-Fi Interface Docking Station for LogTag®Dataloggers (accessory) DISCONTINUED	
ACCD1001PT	Pt100 Immersion Probe, -350/350C, 150mm x 23.0mm (change range from -350/350C to -50/350C ANI	
ACC24130	DOven, 20/130C x 1.0C, Tol Imm (DISCONTINUED)	
ACC24130SC	DOven, 20/130C x 1.0C, Tol Imm, Safety Coated (DISCONTINUED)	
ACC6391S	Incubator, 0/70C x 1.0C, 150mm Imm (CHANGE SPEC FROM 0.5C TO 1.0C DIV)	
ACC6391SSC	Incubator, 0/70C x 1.0C, 150mm Imm, Safety Coat (CHANGE SPEC FROM 0.5C TO 1.0C DIV)	
ACCSL2021	AccuTherm SmartLOG DDL, -40/70C, w/ wireless temperature sensor, email & SMS alerts, 2 year NIST C	
ACCSLBLET	AccuTherm Temperature Sensor for SmartLOG DDL, -40/70C, w/ wireless temperature sensor only, 2 y	
303018IT	DWIRELESS WEATHER STATION, TEMP, HUMIDITY, BAROMETER & TIME (DISCONTINUED)	
ACC077FBLS	D-65/5F x 0.5F (DISCONTINUE)	
ACC077FBLSFC	D-65/5F X 0.5F CERT @ 250, 260F (DISCONTINUE)	
ACC102504	D10/250C x 1.0C, 4" Immersion (DISCONTINUE)	\$46.01
ACC201501221S	D-20/150C x 1.0C, 12" to 21" Immersion (DISCONTINUE)	
ACC2015018S	D-20/150C x 1.0C, 18" Immersion (DISCONTINUE)	
ACC2015018SSC	D-20/150C x 1.0C, 18" Immersion, Safety Coated (DISCONTINUE)	
ACC201501221SSD*	D*-20/150C x 1.0C, 12" to 21" Immersion, Safety Coated (DISCONTINUE)	
ACC201503-7	D-20/150C x 1.0C, 3" to 7" Immersion (DISCONTINUE)	
ACC201503-7SC	D-20/150C x 1.0C, 3" to 7" Immersion, Safety Coated (DISCONTINUE)	
ACC2457	DIncubator, 25,37,& 56C x 0.5C, 76mm Imm (DISCONTINUE)	
ACC2457SC	DIncubator, 25,37 & 56C x 0.5C, 76mm Imm, Safety Coated (DISCONTINUE)	
ACC6000	D10/45 x 0.5 & 0/150F Hg, 15" (DISCONTINUE)	
ACC6001	D45/90 x 0.5 & 0/150F Hg, 15" (DISCONTINUE)	
ACC6002	D10/30 x 0.2 & 0/150F Hg, 15" (DISCONTINUE)	
ACC6003	D20/40 x 0.2 & 0/150F Hg, 15" (DISCONTINUE)	
ACC6004	D30/50 x 0.2 & 0/150F Hg, 15" (DISCONTINUE)	
ACC6005	D40/60 x 0.2 & 0/150F Hg, 15" (DISCONTINUE)	
ACC6006	D50/70 x 0.2 & 0/150F Hg, 15" (DISCONTINUE)	
ACC6007	D60/80 x 0.2 & 0/150F Hg, 15" (DISCONTINUE)	
ACC6008	D70/90 x 0.2 & 0/150F Hg, 15" (DISCONTINUE)	
ACC6009	D80/100 x 0.2 & 0/150F Hg, 15" (DISCONTINUE)	
ACC6110S	D-10/110C x 1.0C C518 DISCONTINUE	
ACC9213TH	DLarge Digit, Indoor/Outdoor Thermo-Hygrometer, -50/70C & F, 20/100% Rh, Min-Max DISCONTINUE	
ACC620233	D-10/400C x 2.0C & +20/750F x 5.0F DISCONTINUED	

Cat#	
B430BLS	D+20/430F x 2.0F DISCONTINUED
ACCST05110S	D-10/110C x 1.0C, 125mm Imm DISCONTINUE
ACCST06110S	D-10/110C x 1.0C, 150mm Imm DISCONTINUE
ACCST04150	D-10/150C x 1.0C, 100mm Imm DISCONTINUE
ACCST05150	D-10/150C x 1.0C, 125mm Imm DISCONTINUE
ACCST05250	D-10/250C x 1.0C, 125mm Imm DISCONTINUE
ACCST06150	D-10/150C x 1.0C, 150mm Imm DISCONTINUE
ACCAB36012	D-10/360C x 1.0C UP TO 12" IMMERSION DISCONTINUE
ACC102504	D10/250C x 1.0C, 4" Immersion DISCONTINUE
ACC102504SC	D10/250C x 1.0C, 4" Immersion, Safety Coated DISCONTINUE
ACC526	D+20/500F x 2.0F (DISCONTINUE - SUITABLE SUB ACC6163)
ACC532AXS	CHANGE PART NUMBER TO ACC532XSA
ACC532AS	CHANGE PART NUMBER TO ACC532SA
ACC534AS	CHANGE PART NUMBER TO ACC534SA
ACC533AS	CHANGE PART NUMBER TO ACC533SA
ACC535AS	CHANGE PART NUMBER TO ACC535SA
ACC534A	D-30/120F x 2.0F, Tol Imm, 160mm length DISCONTINUE
ACC5051P	D-15/105C X 1.0C, 160m length, w/ pocket armor DISCONTINUE
ACC5051A	D-15/105C X 1.0C, 160m length, w/ armor DISCONTINUE
ACC5051	D-15/105C x 1.0C, 160m length, refill DISCONTINUE
ACC062C4	DASTM Like 62C, -38/+2C x 0.1C, w/ auxillary scale, cert @ 0C (DISCONTINUE; DIGITAL REPLACEMENT .
ACC063C4	DASTM Like 63C, -8/32C x 0.1C, w/ auxillary scale, cert @ 0C (DISCONTINUE; DIGITAL REPLACEMENT A
ACC064C4	DASTM Like 64C, +25/55C x 0.1C (DISCONTINUE; DIGITAL REPLACEMENT ACCD0215)
ACC065C4	DASTM Like 65C, +50/80C x 0.1C (DISCONTINUE; DIGITAL REPLACEMENT ACCD0215)
ACC066C4	DATM Like 66C, +75/105C x 0.1C (DISCONTINUE; DIGITAL REPLACEMENT ACCD0215)
ACC067C4	DASTM Like 67C, +95/155C x 0.2C (DISCONTINUE; DIGITAL REPLACEMENT ACCD0215)
ACC068C4	DASTM Like 68C, +145/205C x 0.2C (DISCONTINUE; DIGITAL REPLACEMENT ACCD0215)
ACC069C4	DASTM Like 69C, +195/305C x 0.5C (DISCONTINUE; DIGITAL REPLACEMENT ACCD0215)
ACC070C4	DASTM Like 70C, +295/405C x 0.5C (DISCONTINUE; DIGITAL REPLACEMENT ACCD0215)
ACC062F4	D-36/35F x 0.2F (DISCONTINUE; DIGITAL REPLACEMENT ACCD0215)
ACC063F4	D+18/89F x 0.2F (DISCONTINUE; DIGITAL REPLACEMENT ACCD0215)
ACC064F4	D+77/131F x 0.2F (DISCONTINUE; DIGITAL REPLACEMENT ACCD0215)
ACC065F4	D+122/176F x 0.2F (DISCONTINUE; DIGITAL REPLACEMENT ACCD0215)
ACC066F4	D+167/221F x 0.2F (DISCONTINUE; DIGITAL REPLACEMENT ACCD0215)
ACC067F4	D+203/311F x 0.5F (DISCONTINUE; DIGITAL REPLACEMENT ACCD0215)
ACC068F4	D+293/401F x 0.5F (DISCONTINUE; DIGITAL REPLACEMENT ACCD0215)
ACC069F4	D+383/581F x 1.0F (DISCONTINUE; DIGITAL REPLACEMENT ACCD0215)
ACC070F4	D+563/761F x 1.0F (DISCONTINUE; DIGITAL REPLACEMENT ACCD0215)
ACC021FBLS	D+174/188 x 0.2F DISCONTINUE
ACC035CBLS	D+90/170C x 0.2C DISCONTINUE
ACC108FBLS	D+270/290F x 0.5F DISCONTINUE
ACC109FBLS	D+320/340F x 0.5F DISCONTINUE
ACC052CFC	DASTM 52C, -10/5C x 0.1C, Hg, Tol Imm, CERTIFIED DISCONTINUE
ACC109FBLSFC	D+320/340F x 0.5F DISCONTINUE
ACC052C	DASTM 52C, -10/5C x 0.1C, Hg, Tol Imm DISCONTINUE
ACC059CS	D-18/+82C x 0.5C (DISCONTINUE; REPLACEMENT ACC059CBLS)
ACC059FS	D0/180F x 1.0F (DISCONTINUE; REPLACEMENT ACC059FBLS)
ACCD0020	DK Surface Probe 90° Angle, w/handle, -40/900C & F, 130mm x 8mm DISCONTINUE
ACCD0058	DK Thermocouple beaded wire probe, fiberglass insulated, -40 to 400C & F DISCONTINUE
B950K	DDigital Min/Max Thermocouple K, -50/1300 C&F DISCONTINUE
ACCD303410	DRECEIVER & TRANSMITTER -40 to 65C (-40 to 149F) DISCONTINUE

Cat#	
ACCD303163	DTRANSMITTER ONLY DISCONTINUE
ACCD05011R	DSCAN-TEMP 385, -35 TO 365C & F, RES. 8:1, LAZER SIGHTING, MAX FUNCTION DISCONTINUE SCAN-TEMP 485,-50 to 800C & F, ACCURACY ±2.0C, RES. 20:1, DOUBLE LAZER SIGHTING
ACCD0486IR	DDISCONTINUE
ACC9237TH	DPocket Thermo-Hygrometer, 0/50C & F, +20/100%Rh, Accuracy +/-1C & 5%Rh DISCONTINUE
ACCD0343H	DNon-Contact Moisture Meter, Digital w/Min-Max, Penetration Depth 40mm DISCONTINUE
303018IT	DWIRELESS WEATHER STATION, TEMP, HUMIDITY, BAROMETER & TIME DISCONTINUE
ACCD452019	DThermohygrometer 0/100% Rh & -15/55C DISCONTINUE
ACCD0020PT-1.2M	DPt100 Probe, Cable 1.2 Meters, -50/350C, DISCONTINUE
ACCD1091	DPt100 Probe, w/ out handle, -200/450C, 400 x 6mm, 1500mm cable DISCONTINUE
ACCD1009	DPt100 Temperature & Humidity Probe, -30/100C & 0/100% Rh, 230mm x 12mm Alu. DISCONTINUE
ACCD1006	DPt100 Insertion Probe, 150mm x 4mm, -50/350C & -58/660F DISCONTINUE
ACCD0215	V315 Pt100 Precision Digital Thermometer, -50/250C, 8" stem, 30 x 4.5mm NEW PRODUCT includes 3 point calibration DaKKS Certificate, 0, 75, 150C
ACCD0315	V315 Pt100 Precision Digital Thermometer, -50/250C, 12" stem, 30 x 4.5mm NEW PRODUCT includes 3 point calibration DaKKS Certificate, 0, 75, 150C
ACCD0216	USB Adapter for V215 & V315 Pt100 Precision Digital Thermometer, NEW PRODUCT allows USB export of data to LogConnect App
ACCP70C	D-40/70C x 1.0C DISCONTINUE
YP160FNSF	D-40/160F x 2.0F DISCONTINUE
ACCL550F	D+50/550F x 5.0F DISCONTINUE
ACCL220F	D0/220F x 2.0F DISCONTINUE
B150C300FY3	D-20 to 150C DISCONTINUE
B110C230FW	D-20 to 110C & 0 TO 230F DISCONTINUE
B260C500FY	D -10 to 260C & 20 TO 500F DISCONTINUE
BM2040PS	D-50 to 120F SPIRITI POCKET DISCONTINUED
BM2040AS	D-50 to 120F SPIRIT ARMOR DISCONTINUED
B2040S	D-50 to 120F SPIRIT REFILL DISCONTINUE
CS8110CW2BLS	D-10 to 110C DISCONTINUE
CS8150CW2BLS	D-10 to 150C DISCONTINUE
CS300FW3S	D0 to 300F DISCONTINUE
CS8110CWBLS	D-10 to 110C DISCONTINUE
CS8150CWBLS	D-10 to 150C DISCONTINUE
CS5110CWS	D-10 to 110C DISCONTINUE
CS300FWS	D0 to 300F DISCONTINUE
B200CWS	CHANGE PART NUMBER TO B200CWBLs, BLUE SPIRIT
<u>Partial Immersion, Length 8" (200mm), Spirit</u>	
B81102BLS	D-10 to 110C DISCONTINUE
B81502BLS	D-10 to 150C DISCONTINUE
B82302BLS	D0 to 230F DISCONTINUE
B83002BLS	D0 to 300F DISCONTINUE
<u>Tol Immersion, Length 8" (200mm), Spirit</u>	
B8110BLS	D-10 to 110C DISCONTINUE
B8150BLS	D-10 to 150C DISCONTINUE
B8230BLS	D0 to 230F DISCONTINUE
B8300BLS	D0 to 300F DISCONTINUE
B260CW3BLSSC	D -10 to 260C DISCONTINUE
ACC10073SFC	D-1/201C, x 0.2C, 76mm Imm, Cert@ 0, 25, 37, 44.6, 56, 105, 121, 180, 200C DISCONTINUE SUB ACC1007
ACC10073S	D-1/201C x 0.2C DISCONTINUE SUBSTITUTE ACC10073BLS
GW2743C	D6/12 x 0.1 DISCONTINUE
ACC6000	D10/45 x 0.5 & 0/150F Hg, 15" DISCONTINUE

Cat#	
ACC6001	D45/90 x 0.5 & 0/150F Hg, 15" DISCONTINUE
ACC6002	D10/30 x 0.2 & 0/150F Hg, 15" DISCONTINUE
ACC6003	D20/40 x 0.2 & 0/150F Hg, 15" DISCONTINUE
ACC6004	D30/50 x 0.2 & 0/150F Hg, 15" DISCONTINUE
ACC6005	D40/60 x 0.2 & 0/150F Hg, 15" DISCONTINUE
ACC6006	D50/70 x 0.2 & 0/150F Hg, 15" DISCONTINUE
ACC6007	D60/80 x 0.2 & 0/150F Hg, 15" DISCONTINUE
ACC6008	D70/90 x 0.2 & 0/150F Hg, 15" DISCONTINUE
ACC6009	D80/100 x 0.2 & 0/150F Hg, 15" DISCONTINUE
GW6000XS	D0/20 & 0/150F Spirit, 15" DISCONTINUE
GW6007S	D60/80 x 0.2 & 0/150F Spirit, 15" DISCONTINUE
GW6009S	D80/100 x 0.2 & 0/150F Spirit, 15" DISCONTINUE
ACC062C4FC	D-38/+2C x 0.1C, Cert@ -37, -30, -20, -10, 0C DISCONTINUE
ACC063C4FC	D-8/32C x 0.1C, Cert@ -7, 0, 10, 20, 30C DISCONTINUE
ACC064C4FC	D+25/55C x 0.1C, Cert@ 0, 25, 35, 45, 55C DISCONTINUE
ACC065C4FC	D+50/80C x 0.1C, Cert@ 0, 50, 60, 70, 80C DISCONTINUE
ACC066C4FC	D+75/105C x 0.1C, Cert@ 0, 75, 85, 95, 105C DISCONTINUE
ACC067C4FC	D+95/155C x 0.2C, Cert@ 0, 100, 110, 130, 150C DISCONTINUE
ACC068C4FC	D+145/205C x 0.2C, Cert@ 0, 150, 170, 190, 205C DISCONTINUE
ACC069C4FC	D+195/305C x 0.5C, Cert@ 0, 200, 235, 270, 305C DISCONTINUE
ACC070C4FC	D+295/405C x 0.5C, Cert@ 0, 300, 335, 370, 400C DISCONTINUE
ACC062F4FC	D-36/35F x 0.2F, Cert@ -35, -15, 0, 15, 32F DISCONTINUE
ACC063F4FC	D+18/89F x 0.2F, Cert@ 20, 32, 50, 70, 88F DISCONTINUE
ACC064F4FC	D+77/131F x 0.2F, Cert@ 32, 80, 95, 115, 130F DISCONTINUE
ACC065F4FC	D+122/176F x 0.2F, Cert@ 32, 125, 145, 160, 175F DISCONTINUE
ACC066F4FC	D+167/221F x 0.2F, Cert@ 32, 168, 185, 200, 220F DISCONTINUE
ACC067F4FC	D+203/311F x 0.5F, Cert@ 32, 205, 240, 275, 310F DISCONTINUE
ACC068F4FC	D+293/401F x 0.5F, Cert@ 32, 300, 340, 370, 400F DISCONTINUE
ACC069F4FC	D+383/581F x 1.0F, Cert@ 32, 400, 460, 520, 580F DISCONTINUE
ACC070F4FC	D+563/761F x 1.0F, Cert@ 32, 570, 640, 700, 760F DISCONTINUE
ACC6270C4FC	DASTM 62C to 70C PARTIAL IMMERSION (104mm), -38 to 405C w/ case DISCONTINUE
ACC6270F4FC	DASTM 62F to 70F PARTIAL IMMERSION (104mm), -36 to 761F w/ case DISCONTINUE
ACC532P	D+30/120F x 1.0F, Tol Imm, 160mm length DISCONTINUE
ACC534P	D-30/120F x 2.0F, Tol Imm, 160mm length DISCONTINUE
ACC536P	D0/220F x 2.0F, Tol Imm, 160mm length DISCONTINUE
ACC533P	D-35/50C x 1.0C, Tol Imm, 160mm length DISCONTINUE
ACC535P	D-15/105C x 1.0C, Tol Imm, 160mm length DISCONTINUE
ACC539P	D-10-110C x 1.0C, Tol Imm, 160mm length DISCONTINUE
ACC532A	D+30/120F x 1.0F, Tol Imm, 160mm length DISCONTINUE
ACC534A	D-30/120F x 2.0F, Tol Imm, 160mm length DISCONTINUE
ACC536A	D0/220F x 2.0F, Tol Imm, 160mm length DISCONTINUE
ACC533A	D-35/50C x 1.0C, Tol Imm, 160mm length DISCONTINUE
ACC535A	D-15/105C x 1.0C, Tol Imm, 160mm length DISCONTINUE
ACC532	D+30/120F x 1.0F, Tol Imm, 150mm length DISCONTINUE
ACC534	D-30/120F x 2.0F, Tol Imm, 150mm length DISCONTINUE
ACC536	D0/220F x 2.0F, Tol Imm, 150mm length DISCONTINUE
ACC533	D-35/50C x 1.0C, Tol Imm, 150mm length DISCONTINUE
ACC535	D-15/105C x 1.0C, Tol Imm, 150mm length DISCONTINUE
ACC539	D-10/110C x 1.0C, Tol Imm, 150mm length DISCONTINUE
ACCD144003/TFAD	Freezer/Fridge Thermometer, 8", Hook Top w/ Freeze Zone indicator DISCONTINUE
ACC3303BRNA	D-10/320C x 2.0C, 12" ARMOR DISCONTINUE
ACC3303BRN	D-10/320C x 2.0C, 12" DISCONTINUE
TCT8011T	DIR THERMOMETER, 12:1, -50/1100C & -58/2012F DISCONTINUE

Cat#	
ACCD1800IRK	HI-TEMP IR +100 to 1800C w thermocouple -64 to 1370C, Hold, Min-Max, Dif, Averaging, DLock DISCONTINUE
ACCD2400IRK	HI-TEMP IR +100 to 2400C w thermocouple -64 to 1370C, Hold, Min-Max, Dif, Averaging, DLock DISCONTINUE
ACCD1003	DThermocouple K Surface Probe w/handle, -100 to 1100C & F DISCONTINUE
ACC108FBL SFC	D270/290 F X 0.5F, CERT 275/285F DISCONTINUE
ACC8540PC	1.020/1.220 x 0.002* new range eff. 12/2020
THDI800-010	Accessory; replacment membrane for THDI800
ACC059FBLS	0/180F X 1.0F ADD back to price file
ACC059FBLSFC,	ASTM Like 59F, 0/180F Blue Spirit, CERT @ 32, 80, 130, 180F
ACC10250618S	D-10/250C x 1.0C, 6" to 18" Immersion : DISCONTINUED
ACC10250618SSC	D*-10/250C x 1.0C, 6" to 18" Immersion, Safety Coated : DISCONTINUED
ACC102501224S	D-10/250C x 1.0C, 12" to 24" Immersion : DISCONTINUED
ACC102501224SSD	D-10/250C x 1.0C, 12" to 24" Immersion, Safety Coated : DISCONTINUED
ACC631350	DOven, 0/350C x 1.0C, 125mm Imm DISCONTINUED
ACC1007S	D-1/201C x 0.2C Red spirit DISCONTINUED
ACC610XXSC	D0/230F x 0.5F, Safety Coated DISCONTINUED
B1050CW3BLS	D-10/50C x 0.5C DISCONTINUED
ACCAB	DREPLACEMENT BUSHING DISCONTINUED
ACC8545PC	D30/45% x 0.2% DISCONTINUED
ACC8560PC	D45/60% x 0.2% DISCONTINUED
GW038H	DASTM 38H, 85/91 x 0.1, 163mm DISCONTINUED
GW039H	DASTM 39H, 90/96 x 0.1, 163mm DISCONTINUED
GW060HLS	DASTM 60HL, 89/101 x 0.1 & 0/150F Spirit, 380mm DISCONTINUED
GW142H	DASTM 142H - 1.850/1.900 x 0.001, 260mm DISCONTINUED
GW143H	DASTM 143H - 1.900/1.950 x 0.001, 260mm DISCONTINUED
GW293H	DASTM 293H, 5.83/6.24 x 0.005, 330mm DISCONTINUED
GW294H	DASTM 294H, 6.24/6.66 x 0.005, 330mm DISCONTINUED
GW295H	DASTM 295H, 6.66/7.08 x 0.005, 330mm DISCONTINUED
GW296H	DASTM 296H, 7.08/7.50 x 0.005, 330mm DISCONTINUED
GW297H	DASTM 297H, 7.50/7.91 x 0.005, 330mm DISCONTINUED
GW298H	DASTM 298H, 7.91/8.23 x 0.005, 330mm DISCONTINUED
GW2509C	D79/91 x 0.2 & 0/150F Hg, 7.5" DISCONTINUE
GW2509AC	D89/101 & 0/150F Hg, 7.5" DISCONTINUE
GW2509ACS	D89/101 & 0/150F Spirit, 7.5" DISCONTINUE
GW3003	D80/120 x 0.5 & 20/130F Hg, 7.5" DISCONTINUE
GW6009	D80/100 x 0.2 & 0/150F Hg, 15" DISCONTINUE
GW6007XS	D40/80 & 0/150F Spirit, 15" DISCONTINUE
GW6008S	D70/90 x 0.2 & 0/150F Spirit, 15" DISCONTINUE
ACC059HL	DASTM 59HL, 79/91 x 0.1 & 0/150F Hg, 380mm DISCONTINUE
GW2521C	D49/61 x 0.1 DISCONTINUE
GW2515C	D0 to 14 x 0.1 DISCONTINUE
GW5000	D10/45 x 1.0 & 20/130F Hg, 7.5" DISCONTINUE
GW5002	D10/30 x 0.5 & 20/130F Hg, 7.5" DISCONTINUE
GW041HM	DASTM 41HM, 15/23 x 0.1 & 30/180F Hg, 380mm DISCONTINUE
GW045HM	DASTM 45HM, 43/51 x 0.1 & 30/180F Hg, 380mm DISCONTINUE
GW051HM	DASTM 51HM, -1/11 x 0.1 & 30/180F Hg, 380mm DISCONTINUE
GW058HM	DASTM 58HM, 69/81 x 0.1 & 30/180F Hg, 380mm DISCONTINUE
GW059HM	DASTM 59HM, 79/91 x 0.1 & 30/180F Hg, 380mm DISCONTINUE
GW059HMS	DASTM 59HM, 79/91 x 0.1 & 30/180F Spirit, 380mm DISCONTINUE
GW060HMS	DASTM 60HM, 89/101 x 0.1 & 30/180F Spirit, 380mm DISCONTINUE

Cat#		
GW041HH	DASTM 41HH, 15/23 x 0.1 & 60/220F Hg, 380mm	DISCONTINUE
GW042HH	DASTM 42HH, 22/30 x 0.1 & 60/220F Hg, 380mm	DISCONTINUE
GW044HH	DASTM 44HH, 36/44 x 0.1 & 60/220F Hg, 380mm	DISCONTINUE
GW059HH	DASTM 59HH, 79/91 x 0.1 & 60/220F Hg, 380mm	DISCONTINUE
GW060HH	DASTM 60HH, 89/101 x 0.1 & 60/220F Hg, 380mm	DISCONTINUE
GW041HHS	DASTM 41HH, 15/23 x 0.1 & 60/220F Spirit, 380mm	DISCONTINUE
GW042HHS	DASTM 42HH, 22/30 x 0.1 & 60/220F Spirit, 380mm	DISCONTINUE
GW043HHS	DASTM 43HH, 29/37 x 0.1 & 60/220F Spirit, 380mm	DISCONTINUE
GW044HHS	DASTM 44HH, 36/44 x 0.1 & 60/220F Spirit, 380mm	DISCONTINUE
GW045HHS	DASTM 45HH, 43/51 x 0.1 & 60/220F Spirit, 380mm	DISCONTINUE
GW051HHS	DASTM 51HH, -1/11 x 0.1 & 60/220F Spirit, 380mm	DISCONTINUE
GW052HHS	DASTM 52HH, 9/21 x 0.1 & 60/220F Spirit, 380mm	DISCONTINUE
GW053HHS	DASTM 53HH, 19/31 x 0.1 & 60/220F Spirit, 380mm	DISCONTINUE
GW054HHS	DASTM 54HH, 29/41 x 0.1 & 60/220F Spirit, 380mm	DISCONTINUE
GW055HHS	DASTM 55HH, 39/51 x 0.1 & 60/220F Spirit, 380mm	DISCONTINUE
GW056HHS	DASTM 56HH, 49/61 x 0.1 & 60/220F Spirit, 380mm	DISCONTINUE
GW057HHS	DASTM 57HH, 59/71 x 0.1 & 60/220F Spirit, 380mm	DISCONTINUE
GW058HHS	DASTM 58HH, 69/81 x 0.1 & 60/220F Spirit, 380mm	DISCONTINUE
GW059HHS	DASTM 59HH, 79/91 x 0.1 & 60/220F Spirit, 380mm	DISCONTINUE
GW060HHS	DASTM 60HH, 89/101 x 0.1 & 60/220F Spirit, 380mm	DISCONTINUE
ACC6000S	D10/45 x 1.0, 0/150F Spirit, 15"	DISCONTINUE
ACC6002S	D10/30 x 0.5, 0/150F Spirit, 15"	DISCONTINUE
ACC6008S	D70/90 x 0.5, 0/150F Spirit, 15"	DISCONTINUE
ACC6009S	D80/100 x 0.5, 0/150F Spirit, 15"	DISCONTINUE
ACC032H	DASTM 32H, 55/61 x 0.1, 163mm	DISCONTINUE
ACC036H	DASTM 36H, 75/81 x 0.1, 163mm	DISCONTINUE
ACC037H	DASTM 37H, 80/86 x 0.1, 163mm	DISCONTINUE
ACC038H	DASTM 38H, 85/91 x 0.1, 163mm	DISCONTINUE
ACC040H	DASTM 40H, 95/101 x 0.1, 163mm	DISCONTINUE
ACC041HL	DASTM 41HL, 15/23 x 0.1 & 0/150F Hg, 380mm	DISCONTINUE
ACC043HL	DASTM 43HL, 29/37 x 0.1 & 0/150F Hg, 380mm	DISCONTINUE
ACC044HL	DASTM 44HL, 36/44 x 0.1 & 0/150F Hg, 380mm	DISCONTINUE
ACC045HL	DASTM 45HL, 43/51 x 0.1 & 0/150F Hg, 380mm	DISCONTINUE
ACC059HH	DASTM 59HH, 79/91 x 0.1 & 60/220F Hg, 380mm	DISCONTINUE
ACC060HH	DASTM 60HH, 89/101 x 0.1 & 60/220F Hg, 380mm	DISCONTINUE
ACC071H	DASTM 71H, -1/11 x 0.1 & 30/220F Hg, 380mm	DISCONTINUE
LTUSRID16	Single Use USB PDF digital datalogger, -30/60C **new product	\$29.00
VFC400-VMKCS	Vaccine Monitoring Kit, includes VFC400 datalogger, probe, buffer, acrylic stand, wall mount bracket, ISO 17025 Certificate. **Excludes LTIHID, docking station *NEW PRODUCT	\$147.00
VFC400-SMKCS	Site Vaccine Monitoring Kit, includes (2) VFC400 dataloggers, probes w/buffer, acrylic stands, wall mount brackets, ISO 17025 Certificate. Includes 1 LTIHID, docking station *NEW PRODUCT	\$299.00
LTUHADO-16	Data Logger Temperature and Humidity w/ display; -30°C to +70°C and 0/100%RH. USB **NE	\$102.41
ACC110C	DASTM 110C, 133.6/136.4 x 0.05C, Hg, Tol Imm	DISCONTINUE - unable to meet specifications
ACC110CFC	DASTM 110C, 133.6/136.4 x 0.05C, Hg, Tol Imm, CERT @ 0, 135, 138C	DISCONTINUE - unable to meet spe
B8TC300BLS	E0 to 300F	DISCONTINUE
B930FY	E20 to 930	DISCONTINUE
ACCD1091	EPt100 Probe, w/ out handle, -200/450C, 400 x 6mm, 1500mm cable	DISCONTINUE
ACCD1092	EPt100 Immersion Probe, 400mm x 6mm, -200/660C & -330/1220F, w/out handle, w/ 1500mm cable	DISCON
ACCD50900017	IProtective Splash Cover for P700 Series, w/ strap **UPDATE PART NUMBER AND DESCRIPTION	
ACCD710PEX	IThermomcouple Thermometer, Single Channel, -200C to 1760C & F, Resolution 0.1C	DISCONTINUE

Cat#	
ACCD1037	LThermocouple K Insertion Probe, -100 to 1100C & F, 300MM x 4mm DISCONTINUE
ACCD0014	DK Surface Probe Straight, w/handle, -40/900C & F, 130mm x 8mm DISCONTINUE
TCT307K	LThermocouple probe DISCONTINUE
ACC8500V	LThermometer w/ Probe and 100ml Glass Bottle of Sand for Ovens DISCONTINUE
ACC0490IR	SCAN-TEMP 490, -60 to 1000C, W/THERMOCOUPLE INPUT -64 to 1370C, 50:1 DISCONTINUE
ACC0510IR	PRO-SCAN 510 HIGH ACCURACY, ±1.0C, REPRODUCABILITY ±0.7C, -32 to 530C & F, DISCONTINUE
ACC0520IR	PRO-SCAN 520, SAME AS PRO-SCAN 510 BUT W/40:1 RES AND THERMOCOUPLE INPUT DISCONTINUE
ACC0896	SCAN-TEMP RH 896, INFRA-RED, DEW POINT HUMIDITY, & THERMOCOUPLE, IR -50 to 500C, DISCONTINUE
ACCDTA122	LThermo-Hygrometer, Dew Point, Wet Bulb, Min-Max, 0/99%Rh, -40/70C DISCONTINUE
ACCD0890DP	LDifferential Pressure Meter w/USB Interface, Range 0 to 140 mbar x 0.01, Acc. +/- 0.2% FS DISCONTINUE
ACC9201AV	LAnemometer, Air Flow, w/ Min-Max, Averaging, Dual Readings, Range 0.40 to 30 m/s DISCONTINUE
ACCD0328SL	LSound Level Meter, Range 32 to 130 dB x 0.1 dB, Accuracy 1.4dB DISCONTINUE
ACCD0322L	High Precision Sound Level Meter, w/ Data Logger Function 32,000 Value Memory, w/Software, LMin-Max, Calibration Adjustable, Thread for Tripod, Carrying Case DISCONTINUE
ACCD0326CAL	LHigh Precision Sound Meter Calibrator, 94 dB / 114 dB DISCONTINUE
ACCD037LM	LLight Meter, Range 0.00 to 40,000 lx or 0.000 to 4000 fc, Accuracy +/- 3% + 5 digits. DISCONTINUE
ACCBR1000	LHIGH PRECISION BAROMETER DISCONTINUE
BPS108FS	LReplacement Tube, Spirit-Filled, 0° to 120°F DISCONTINUE
BPS108CS	LReplacement Tube, Spirit-Filled, -5° to 50°C DISCONTINUE
ACCD0342H	LMoisture Mini Meter, Digital 6 to 42% (wood), Temp. 0 to 40C & F DISCONTINUE
ACCD0340H	LMoisture Mini Meter, Analog 6 to 44% (wood), Accuracy +/-1% DISCONTINUE
ACC5024TH	LIn & out Temperature, Humidity w/Min-Max DISCONTINUE
ACC811C	LLarge Digit, Indoor/Outdoor Thermo-Hygrometer, -50/70C & F, 20/100% Rh, Min-Max DISCONTINUE
ACCDTA120T	LWireless Transmitter for Temperature & Humidity DISCONTINUE
ACCD0080	Windows Software **UPDATE PART NUMBER & DESCRIPTION ACCD0081PD, DE-Graph Windows Software for P470/ P600/P700/T900-series
ACCD3001-20	LPt100 Immersion Probe, w/20' Cable -50/350C & F, 150mm x 3mm DISCONTINUE
ACC955K	LTYPE K BEADED WIRE PROBE, FIBERGLASS INSULATED DISCONTINUE
ACCG300F-C	L+60/300F x 2.0F & +16/140C x 1.0C DISCONTINUE
ACCGG125F	L+25/125F x 1.0F DISCONTINUE
ACCL220FC	L0/220F x 2.0F & -10/100C x 1.0C DISCONTINUE
ACCL550FC	L+50/550F x 5.0F & +10/290C DISCONTINUE
	3-5 PROCESSING TURN AROUND, RECALIBRATION, *excludes transit time
EXPEDITE FEE (new)	
ACC1103BLSFC	-20/110C X 1.0C, 76mm Imm, Cert@ 0, 50, 100C NEW PRODUCT
B2019S	D-10/110Cx1.0C pocket (DISCONTINUE)
ACC023CBLS	D +18/28C x 0.2C DISCONTINUE
ACC023CBLSFC	D +18/28C x 0.2C, CERTIFIED DISCONTINUE
ACC0400C	D Melting Point Thermometers, 0 to 400C x 1.0C, 30mm Immersion DISCONTINUE
ACC041CBLS	D +98/152C x 0.2C DISCONTINUE
ACC041CBLSFC	D +98/152C x 0.2CM CERTIFIED DISCONTINUE
ACC044CBLS	D +18.6/21.4C x 0.05C DISCONTINUE
ACC044CBLSFC	D +18.6/21.4C x 0.05C, CERTIFIED DISCONTINUE
ACC520SA	D0/300F x 2.0F, 12" ARMOR DISCONTINUE
ACC520S	D0/300F x 2.0F, 12" DISCONTINUE
ACC537	DDISCONTINUE
ACC610XXS	D0/230F x 0.5F, DISCONTINUE
ACC610XXSSC	D0/230F x 0.5F, Safety Coated, DISCONTINUE

Cat#

ACC613XXS U	
ACC7033	D-10/500C x 2.0C DISCONTINUE
ACCST04360	D-10/360C x 1.0C, 100mm Imm DISCONTINUE
ACCST06250	D-10/250C x 1.0C, 150mm Imm DISCONTINUE
ACCD0341H	EMoisture Meter, Digital 0 to 100%, 8 Calibration Scales for 150 Wood Species DISCONTINUE
B10050CWS	D-100/50 C Tol Red Spirit, 1C, 305mm DISCONTINUED
B10050CW3S	D-100/50 C Partial Red Spirit, 1C, 305mm DISCONTINUED
B10050CWSSC	D-100/50 C Tol Red Spirit, 1C, safety coat, 305mm DISCONTINUED
B10050CW3SSC	D-100/50 C Partial Red Spirit, 1C, safety coat, 305mm DISCONTINUED
EL-USB-2-LCD	Temperature, Humidity and Dew Point DDL w/ LCD screen. -35 to +80°C *excludes calibration NEW PRODUCT
EL-WIFI-VACX	-40/60C, external buffered sensor*excludes calibration NEW PRODUCT
EL-WIFI-TH	-20/60C & 0/100%, DDL *EXCLUDES CALIBRATION NEW PRODUCT
EL-WIFI-21CFR-TH	-20/60C & 0/100%, DDL, 21 CFR COMPLIANT*EXCLUDES CALIBRATION NEW PRODUCT
ACC1005SFC	-1/101C x 0.1C, Tol Imm, Cert@ 0, 10, 20, 25, 37, 44.5, 56, 70, 80, 90, 100C ADDED TEMP POINT
ACC10053SFC	-1/101C x 0.1C, 76mm Imm, Cert@ 0, 10, 20, 25, 37, 44.5, 56, 70, 80, 90, 100C ADDED TEMP POINT
ACCGG40C	d-74/40C & -100/+100F x 2.0C & 2.0F DISCONTINUE
ACCGG110C	d-10/110C x 1.0C DISCONTINUE
ACCGG150C	d0/150C x 1.0C DISCONTINUE
ACCGG250C	d0/250C x 2.0C DISCONTINUE
ACCGG180F	d0/180F x 2F DISCONTINUE
ACCGG220F	d0/220F x 2.0F DISCONTINUE
ACCGG300F	d+50/300F x 2.0F DISCONTINUE
ACCGG500F	d+50/550F x 5.0F DISCONTINUE
ACC016CBLS	+30/200C X 1.0C CHANGE SPECIFICATIONS TO 30/200C X 0.5C
ACC016CBLSFC	+30/200C X 1.0C CHANGE SPECIFICATIONS TO 30/200C X 0.5C
GW2532	
GW2533	
GW2534	
GW2534c	
ACC118CBLS	+28.6/31.4C x 0.05C *revised pricing - price increase (\$221.92 / \$110.43)
EL-IOT-CO2	-20/60C, 0/100%RH, 0/40k PPM, AA Battery, air quality & CO2 DDL
VFC400-SP	VFC400 w/ Smart Probe technology, hot swap calibrated probe. VFC compliant ISO 17025 certificate, requires LTIHID for operation
VFC800WIFI-SP	VFC800wifi w/ Smart Probe technology, hot swap calibrated probe. VFC compliant ISO 17025 certificate, requires LTO subscription
UTRED30WIFI-SP	LogTag® Thermco WIFI DDL w/ Smart Probe technology, hot swap calibrated probe. Single probe. VFC compliant ISO 17025 certificate, requires LTO subscription
UTRED30WIFI-SP2	LogTag® Thermco WIFI DDL w/ Smart Probe technology, hot swap calibrated probe. Dual probe. VFC compliant ISO 17025 certificate, requires LTO subscription
VFCSP-5	Smart probe, 5ft., calibrated w/ 60ml glycol, for VFC400®
VFCSP-10	Smart probe, 10ft., calibrated w/ 60ml glycol, for VFC800®