

Item	Description	List Price	Comments	Replacement Information			
				Item	Part #	Replacement Description	List Price
<b>METERS- PORTABLE pH</b>							
475660	CHEK-MITE pH-10 pH TESTER	Discontinued	Replaced w/ EcoSense Model	EcoSense pH10	606pH10	EcoSense® pH10 pH/temperature pen	\$115.00
475661	CHEK-MITE pH-15 pH TESTER	Discontinued	Replaced w/ EcoSense Model	EcoSense pH10	606pH10	EcoSense® pH10 pH/temperature pen	\$115.00
475662	CHEK-MITE pH-20 pH TESTER	Discontinued	Replaced w/ EcoSense Model	EcoSense pH10	606pH10	EcoSense® pH10 pH/temperature pen	\$115.00
475663	CHEK-MITE pH 25 pH TESTER	Discontinued	Obsolete, no replacement				
475875	MODEL M 313P pH/MV/°C METER ONLY	Discontinued	Replaced w/ SI Analytics model	Handylab pH 11	285203940	handylab pH11 BNC, individual meter	\$417.00
475870	MODEL M 313P pH/MV/°C METER W/ ELECTRODE	Discontinued	Replaced w/ SI Analytics model	Handylab pH 11 kit	285203960	handylab pH11 BNC, case set, complete, ready to use	\$493.00
475876	MODEL M 314P pH/MV/°C METER ONLY	Discontinued	Replaced w/ SI Analytics model	Handylab pH 12	285203990	handylab pH12 BNC, individual meter	\$493.00
475871	MODEL M 314P pH/MV/°C METER W/ELECTRODE	Discontinued	Replaced w/ SI Analytics model	Handylab pH 11 kit	285204020	handylab pH12 BNC, case set, complete, ready to use	\$744.00
475872	MODEL M 315P pH/MV/°C METER W/ELECTRODE	Discontinued	Obsolete, no replacement				
475877	MODEL M 315P pH/MV/°C METER ONLY	Discontinued	Obsolete, no replacement				
<b>METERS- PORTABLE CONDUCTIVITY/TDS</b>							
475873	MODEL M 316P CONDUCTIVITY METER W/ ELECTRODE	Discontinued	Replaced w/ SI Analytics model	Handylab LF 12 Kit	285203810	handylab LF12, case set, complete, ready to use	\$1,040.00
475878	MODEL M 316P CONDUCTIVITY METER ONLY	Discontinued	Replaced w/ SI Analytics model	Handylab LF 12	285203840	handylab LF12, individual meter	\$755.00
<b>SCHOLAR BENCHTOP pH METERS</b>							
475894	SCHOLAR 425 pH METER ONLY (120V)	Discontinued	Replaced w/ SI Analytics model	Lab 850	285204050	Measuring parameters pH, mV, temp., micro	\$574.00
475895	SCHOLAR 425 pH METER WITH REFILLABLE ELECTRODE (120V)	Discontinued	Replaced w/ SI Analytics model	Lab 850 kit	285204070	Measuring parameters pH, mV, temp., micro	\$689.00
475898	SCHOLAR 425 pH METER WITH GEL FILLED ELECTRODE (120V)	Discontinued	Replaced w/ SI Analytics model	Lab 850 kit	285204080	Measuring parameters pH, mV, temp., micro	\$629.00
475896	SCHOLAR 425 pH METER ONLY (220V)	Discontinued	Replaced w/ SI Analytics model	Lab 850	285204050	Measuring parameters pH, mV, temp., micro	\$574.00
475897	SCHOLAR 425 pH METER WITH REFILLABLE ELECTRODE (220V)	Discontinued	Replaced w/ SI Analytics model	Lab 850 kit	285204070	Measuring parameters pH, mV, temp., micro	\$689.00
475899	SCHOLAR 425 pH METER WITH GEL FILLED ELECTRODE (220V)	Discontinued	Replaced w/ SI Analytics model	Lab 850 kit	285204080	Measuring parameters pH, mV, temp., micro	\$629.00
<b>METERS- BENCHTOP PH UNIVERSAL POWER SUPPLY</b>							
475730	M 530P pH METER W/ARM, 120V	Discontinued	Replaced w/ SI Analytics model	Lab 850	285204050	Measuring parameters pH, mV, temp., micro	\$574.00
475731	M 530P pH METER W/ARM & pH ELECTRODE, 120V	Discontinued	Replaced w/ SI Analytics model	Lab 850 kit	285204070	Measuring parameters pH, mV, temp., micro	\$689.00
475734	M 540P pH METER W/ARM, 120V	Discontinued	Replaced w/ SI Analytics model	Lab 850	285204050	Measuring parameters pH, mV, temp., micro	\$574.00
475735	M 540P pH METER W/ARM & pH ELECTRODE, 120V	Discontinued	Replaced w/ SI Analytics model	Lab 850 kit	285204070	Measuring parameters pH, mV, temp., micro	\$689.00
475771	M 545P pH METER W/ARM, 120V	Discontinued	Replaced w/ SI Analytics model	Lab 860	285204100	Meter only, measuring parameters pH, mV, t	\$657.00
475772	M 545P pH METER W/ARM & pH ELECTRODE, 120V	Discontinued	Replaced w/ SI Analytics model	Lab 860 kit	285204120	Measuring parameters pH, mV, temp., micro	\$1,050.00
475738	M 555P pH/ION METER W/ARM, 120V	Discontinued	Replaced w/ SI Analytics model	ProLab 3000	285204250	Meter only, sensor + user identification. QVG	\$1,698.00
475739	M 555P, pH/ION METER W/ARM & HP pH ELEC, 120V	Discontinued	Replaced w/ SI Analytics model	ProLab 3000 kit	285203860	Sensor + user identification. QVGA colour dis	\$1,998.00
<b>METERS- BENCH TOP CONDUCTIVITY - 120V</b>							
475742	M 541P CONDUCTIVITY METER W/ARM, 120V	Discontinued	Replaced w/ SI Analytics model	Lab 960	285204170	Meter only, measuring ranges 0.000 µS/cm...	\$875.00
475743	M 541P CONDUCTIVITY METER W/ARM & COND PROBE, 120V	Discontinued	Replaced w/ SI Analytics model	Lab 960 Kit	285204180	Measuring ranges 0.000 µS/cm...500 mS/cm,	\$1,258.00
<b>CONDUCTIVITY PROBES- For use with Meters specified</b>							
476501	CARBON CONDUCTIVITY, 0-1000MS, EPOXY BODY, 36" CABLE (Model M542)	Discontinued	Obsolete, no replacement				
476501P	CONDUCTIVITY PROBE (Model, M541P, M316P)	Discontinued	Replaced w/ SI Analytics model	LF 413 T	285106172	4 pole cell, plastic shaft, 1.5 m cable with 8 p	\$509.00
476502	GLASS CONDUCTIVITY, LOW RANGE 0-500µS (Model M542)	Discontinued	Obsolete, no replacement				
476503	GLASS CONDUCTIVITY, HIGH RANGE 0-500MS (Model M542)	Discontinued	Obsolete, no replacement				
<b>COMBINATION ELECTRODES</b>							
476043	SCHOLAR ELECTRODE - GEL FILLED	\$120.00	Same product & p/n. Rebranded to SI Analytics.				
476044	SCHOLAR ELECTRODE - REFILLABLE	\$116.00	Same product & p/n. Rebranded to SI Analytics.				
476086	GENERAL PURPOSE COMBO, REFILLABLE, BNC	\$224.00	Same product & p/n. Rebranded to SI Analytics.				
476146	HIGH PERFORMANCE, COMBO, REFILLABLE, BNC	\$445.00	Same product & p/n. Rebranded to SI Analytics.				
476156	SEMI-MICRO, COMBO, REFILLABLE, BNC	\$323.00	Same product & p/n. Rebranded to SI Analytics.				
476286	FLAT SURFACE, COMBO, REFILLABLE, BNC	\$310.00	Same product & p/n. Rebranded to SI Analytics.				
476296	RUGGED BULB COMBO, BNC	\$315.00	Same product & p/n. Rebranded to SI Analytics.				
476306P	DEEP VESSEL ELECTRODE, EPOXY BODY	\$232.00	Same product & p/n. Rebranded to SI Analytics.				
476316	FAST-FLOW COMBO, REFILLABLE, BNC	Discontinued	Obsolete, no replacement				
476326	GENERAL PURPOSE GEL COMBO, BNC	\$170.00	Same product & p/n. Rebranded to SI Analytics.				
476346	SEMI MICRO COMBO, BNC	\$229.00	Same product & p/n. Rebranded to SI Analytics.				
476436	3-IN-1 COMBO, REFILLABLE, BNC & RCA (ATC)	\$268.00	Same product & p/n. Rebranded to SI Analytics.				
476436P	PH-E-T 3-IN-1 COMB, ELECTRODE W/ RJ, BNC+BANANA PLUG, EPOXY BODY	\$416.00	Same product & p/n. Rebranded to SI Analytics.				

476456	SPEAR COMBO, REFILLABLE, BNC	\$439.00	Same product & p/n. Rebranded to SI Analytics.				
476466	HIGH PERFORMANCE 3-IN-1 GEL COMBO, BNC & RCA	\$252.00	Same product & p/n. Rebranded to SI Analytics.				
476466P	HIGH PERFORMANCE 3-IN-1 GEL COMBO, BNC & BANANNA	\$252.00	Same product & p/n. Rebranded to SI Analytics.				
476476	SPEAR GEL COMBO, BNC	\$386.00	Same product & p/n. Rebranded to SI Analytics.				
476566	PREMIUM GEL COMBO, GLASS	\$350.00	Same product & p/n. Rebranded to SI Analytics.				
<b>REPLACEABLE JUNCTIONS</b>							
477269	REPLACEABLE JUNCTIONS, CERAMIC, FOR USE WITH ALL ELECTRODES (3 PAC	\$89.00	Same product & p/n. Rebranded to SI Analytics.				
<b>pH ELECTRODES (REFERENCE ELECTRODE REQUIRED)</b>							
476026	HALF CELL pH, BNC	\$209.00	Same product & p/n. Rebranded to SI Analytics.				
<b>REFERENCE ELECTRODES (pH ELECTRODE REQUIRED)</b>							
476002	GENERAL PURPOSE, CALOMEL, PIN TIP	\$258.00	Same product & p/n. Rebranded to SI Analytics.				
476416P	DOUBLE JUNCTION, REFERENCE,	\$298.00	Same product & p/n. Rebranded to SI Analytics.				
<b>METALLIC ELECTRODES (REDOX)</b>							
476516	REDOX COMBO, REFILLABLE, BNC	\$390.00	Same product & p/n. Rebranded to SI Analytics.				
476526	SILVER BILLET, 2MM PIN	\$345.00	Same product & p/n. Rebranded to SI Analytics.				
<b>PORTABLE METER - pH ELECTRODES (FOR USE WITH MODEL REFERENCED)</b>							
473040	REPLACEMENT pH ELECTRODE FOR CHEK-MITE pH-10	Discontinued	Obsolete, no replacement, Refer to EcoSense pH10				
473051	REPLACEMENT pH ELECTRODE FOR CHEK-MITE pH-15	Discontinued	Obsolete, no replacement, Refer to EcoSense pH10				
473049	REPLACEMENT pH ELECTRODE FOR CHEK-MITE pH-20	Discontinued	Obsolete, no replacement, Refer to EcoSense pH10				
476613	GEN PURPOSE, GEL FILL, 6' CABLE (M 313, M 314, M 315)	\$296.00	Same product & p/n. Rebranded to SI Analytics.				
476613P	PINNACLE GEN PURPOSE, GEL FILL, 6' CABLE (M 313P, M 314P, M 315P)	\$296.00	Same product & p/n. Rebranded to SI Analytics.				
476614	PLASTIC, REFILLABLE, IP67 RATING (MODEL M 313, M 314, M 315)	\$315.00	Same product & p/n. Rebranded to SI Analytics.				
476625	GENERAL PURPOSE, 15' CABLE (MODEL M 313, M 314, M 315)	\$347.00	Same product & p/n. Rebranded to SI Analytics.				
<b>ION SELECTIVE ELECTRODES</b>							
476126	CHLORIDE SOLID STATE, HALF-CELL, BNC	Discontinued	Replaced w/ WTW model	CI 500 BNC	106610	Chloride selective electrode with 3.3 ft (1 m)	\$797.00
479030	CYANIDE SOLID STATE, HALF-CELL, BNC	Discontinued	Replaced w/ WTW model	CN 500 BNC	106612	Cyanide selective electrode with 3.3 ft (1 m) c	\$797.00
479031	IODIDE SOLID STATE, HALF-CELL, BNC	Discontinued	Replaced w/ WTW model	I 500 BNC	106618	Iodide selective electrode with 3.3 ft (1 m) ca	\$831.00
476042	FLUORIDE SOLID STATE, HALF-CELL, BNC	Discontinued	Replaced w/ WTW model	F 500 BNC	106616	Fluoride selective electrode with 3.3 ft (1 m) c	\$797.00
476129	SILVER SULFIDE SOLID STATE, HALF-CELL, BNC	Discontinued	Replaced w/ WTW model	Ag/S 500 BNC	106600	Silver ion/sulfide selective electrode 3.3 ft (1	\$854.00
476130	AMMONIA, THREE MEMBRANE CAPS AND FILL SOLUTION, BNC	Discontinued	Replaced w/ WTW model	NH 500/2	106395	Ammonia gas-sensitive combination electrodo	\$824.00
476134	NITRATE PVC, TWO MEMBRANE CAPS AND FILL SOLUTION, BNC	Discontinued	Replaced w/ WTW model	NO 500 BNC	106624	Nitrate selective electrode with 3.3 ft (1 m) c	\$816.00
476132	POTASSIUM PVC, TWO MEMBRANE CAPS AND FILL SOLUTION, BNC	Discontinued	Replaced w/ WTW model	K 500 BNC	106620	Potassium ion selective electrode with 3.3 ft (	\$816.00
476123	SODIUM COMBINATION, BNC	Discontinued	Replaced w/ WTW model	DX 223 NA	106375	Sodium ion-sensitive glass electrode with plu	\$387.00
479022	CHLORIDE COMBINATION, BNC	Discontinued	Replaced w/ WTW model	CI 800 BNC	106660	Combined Chloride selective electrode, with	\$1,025.00
476137	NITRATE COMBINATION, BNC	Discontinued	Replaced w/ WTW model	NO 800 BNC	106674	Combined Nitrate selective electrode, with 3	\$1,043.00
479028	NITRATE COMBINATION, BNC	Discontinued	Obsolete, no replacement				
479000	REFERENCE ELECTRODE, BANANA PLUG	Discontinued	Replaced w/ WTW model	R 503/D	106571	Reference electrode with electrolyte bridge, f	\$434.00
<b>pH REAGENTS AND STANDARDS</b>							
478551	4.00 BUFFER RED 1BT/CS 500ML	\$18.00	Same product & p/n. Rebranded to SI Analytics.				
478540	4.00 BUFFER RED 2BT/CS 500ML	\$39.00	Same product & p/n. Rebranded to SI Analytics.				
478552	7.00 BUFFER YELLOW 1BT/CS 500ML	\$18.00	Same product & p/n. Rebranded to SI Analytics.				
478570	7.00 BUFFER YELLOW 2BT/CS 500ML	\$39.00	Same product & p/n. Rebranded to SI Analytics.				
478553	10.01 BUFFER BLUE 1BT/CS 500ML	\$18.00	Same product & p/n. Rebranded to SI Analytics.				
478510	10.01 BUFFER BLUE 2BT/CS 500ML	\$39.00	Same product & p/n. Rebranded to SI Analytics.				
478554	RAINBOW BUFFER PACK, 1X7, 1X4, 1X10, 3BT/CS 500ML	\$54.00	Same product & p/n. Rebranded to SI Analytics.				
478574	RAINBOW BUFFER PACK, 2X7, 2X4, 2X10, 6BT/CS 500ML	\$106.00	Same product & p/n. Rebranded to SI Analytics.				
477000	SATURATED KCl 6BT/CS 125ML	\$77.00	Same product & p/n. Rebranded to SI Analytics.				
477006	3M KCl 6BT/CS 125ML	\$144.00	Same product & p/n. Rebranded to SI Analytics.				
477007	CASE (2) 3M KCl SOLUTION, 125 ML BTLS	\$42.00	Same product & p/n. Rebranded to SI Analytics.				
473651	4.00 BUFFER SACHET	\$21.00	Same product & p/n. Rebranded to SI Analytics.				
473650	7.00 BUFFER SACHET	\$21.00	Same product & p/n. Rebranded to SI Analytics.				
473652	10.01 BUFFER SACHET	\$21.00	Same product & p/n. Rebranded to SI Analytics.				
473676	MULTI-PACK BUFFER SACHETS	\$69.00	Same product & p/n. Rebranded to SI Analytics.				

<b>CONDUCTIVITY STANDARDS</b>						
473623	CONDUCTIVITY STANDARD; 1413µS 500ML	\$71.00	Same product & p/n. Rebranded to SI Analytics.			
473624	CONDUCTIVITY STANDARD; 12.88MS 500ML	\$71.00	Same product & p/n. Rebranded to SI Analytics.			
473737	CONDUCTIVITY STD. SACHETS 12.88MS	\$29.00	Same product & p/n. Rebranded to SI Analytics.			
473738	CONDUCTIVITY STD. SACHETS 1413µS	\$29.00	Same product & p/n. Rebranded to SI Analytics.			
473740	CONDUCTIVITY STD. SACHETS 84µS	\$29.00	Same product & p/n. Rebranded to SI Analytics.			
473741	CONDUCTIVITY STD. SACHETS 80MS	\$29.00	Same product & p/n. Rebranded to SI Analytics.			
<b>DISSOLVED OXYGEN SOLUTION</b>						
473739	ZERO OXYGEN STD. SACHETS	\$29.00	Same product & p/n. Rebranded to SI Analytics.			
<b>ISE REAGENTS AND STANDARDS</b>						
478172	DILUENT, FLUORIDE ANALYSIS, 3.785 LITERS	\$145.00	Obsolete; Reffer to WTW products			
478153	MOLARITY STANDARD, AMMONIUM SULFATE, 2.0M, 500ML	Discontinued	Obsolete			
478161	MOLARITY STANDARD, AMMONIUM 0.100M, 500ML	Discontinued	Obsolete			
478162	MOLARITY STANDARD, CALCIUM 0.100M, 500ML	\$71.00	Same product & p/n. Rebranded to SI Analytics.			
478168	MOLARITY STANDARD, CHLORIDE 0.100M, 500ML	\$71.00	Same product & p/n. Rebranded to SI Analytics.			
478165	MOLARITY STANDARD, FLUORIDE 0.100M, 500ML	\$71.00	Same product & p/n. Rebranded to SI Analytics.			
478163	MOLARITY STANDARD, NITRATE 0.100M, 500ML	Discontinued	Obsolete			
478154	MOLARITY STANDARD, POTASSIUM NITRATE, 1.0M (10%), 500ML	\$71.00	Same product & p/n. Rebranded to SI Analytics.			
478160	PPM STANDARD, AMMONIA AS N, 1000PPM, 500ML	\$71.00	Same product & p/n. Rebranded to SI Analytics.			
478171	PPM STANDARD, CHLORIDE, 100PPM, 500ML	\$71.00	Same product & p/n. Rebranded to SI Analytics.			
478170	PPM STANDARD, FLUORIDE, 100PPM, 500ML	\$71.00	Same product & p/n. Rebranded to SI Analytics.			
478169	PPM STANDARD, SODIUM, 100PPM, 500ML	Discontinued	Obsolete			
476148	0.1M NaHCO3 STANDARD, 475ML	Discontinued	Obsolete			
476161	1000 PPM SODIUM STD., 475ML	Discontinued	Obsolete			
479013	ISA/FK ionic strength adapter NaNO3,5 mol/l, 2 bottles, 250 mL each	\$100.00	Replacement; Reffer to WTW product 140111	ISA/FK	140110	ISA sample conditioning solution for for solid
479014	ISA/Ca ionic strength adapter KCl, 4 mol/l, 250 mL	\$92.00	Replacement; Reffer to WTW product 140121	ISA/Ca	140120	ISA sample conditioning solution for calcium
479015	ISA/Ag/FI for photographic fixing baths	\$228.00	Replacement; Reffer to WTW product 140141	ISA/Ag/FI	140140	ISA for determination of silver in fixer solution
479016	TISAB/NO3 ionic strength buffer	\$228.00	Replacement; Reffer to WTW product 150121	TISAB/NO3	150120	TISAB total ionic strength buffer solution for r
479017	MZ/NH3/CN NaOH 10 mol/l	\$92.00	Replacement; Reffer to WTW product 150131	MZ/NH3/CN	150130	Alkaline reagent, NaOH, 10 mol/l, 1 bottle with
<b>ISE REPLACEMENT PARTS</b>						
477693	MEMBRANE CAP, AMMONIA, W/FILL SOLUTION	\$172.00	Same product & p/n. Rebranded to SI Analytics.			
477394	FILL SOLUTION, AMMONIA	\$89.00	Same product & p/n. Rebranded to SI Analytics.			
<b>ACCESSORIES - BENCH METERS</b>						
470493P	PINNACLE THERMAL PRINTER	\$791.00	Obsolete			
477758	SOFTWARE INTERFACE CABLE (MODEL 340/345/440/442/443I/445)	\$150.00				
478930P	PINNACLE SERIES METER ELECTRODE ARM	\$139.00	Same product & p/n. Rebranded to SI Analytics.			
478939	POWER SUPPLY FOR PINNACLE SERIES METERS	\$63.00	Same product & p/n. Rebranded to SI Analytics.			
470445P	COMPUTER COMMUNICATION CABLE FOR PINNACLE SERIES BENCHTOP MET	\$70.00	Same product & p/n. Rebranded to SI Analytics.			
470455P	COMPUTER COMMUNICATION CABLE FOR PINNACLE SERIES HANDHELD MET	\$48.00	Same product & p/n. Rebranded to SI Analytics.			
470136P	REPLACEMENT PAPER FOR PINNACLE PRINTER	\$11.00	Same product & p/n. Rebranded to SI Analytics.			
471226	AC ADAPTER 300/400 SERIES BT METERS	\$96.00				
<b>ACCESSORIES - OTHER METERS</b>						
470530	REPLACEMENT BATTERY FOR CHEK-MITE SERIES	Discontinued	Obsolete, no replacement			
<b>AUTOMATIC TEMPERATURE COMPENSATORS</b>						
477968	AUTOMATIC TEMP COMPENSATOR 30K	\$234.00				
477968P	PINNACLE STAINLESS STEEL	\$234.00	Same product & p/n. Rebranded to SI Analytics.			
477969	STAINLESS STEEL, PT1000, FOR USE WITH MODELS 350, 355, 450 AND 455	\$430.00				
<b>ADAPTERS</b>						
478381	ADAPTER, BNC METER TO U.S. STANDARD ELECTRODE	\$131.00				
478383	ADAPTER, U.S. STANDARD METER TO BNC COMBINATION ELECTRODE	\$131.00				


## Portables

Part #	Item	Parameter	Replacement Description	List Price
285203940	handylab pH11 BNC	pH	handylab pH11 BNC, individual meter	\$417
285203960	handylab pH11 BNC/25pH	pH	handylab pH11 BNC, case set, complete, ready to use with pH combination electrode BlueLine 25pH, calibration solutions and plastic beakers	\$493
285203990	handylab pH12 BNC	pH	handylab pH12 BNC, individual meter	\$493
285204020	handylab pH12 BNC/26pH	pH	handylab pH12 BNC, case set, complete, ready to use with pH combination electrode BlueLine 26pH, calibration solutions and plastic beakers	\$744
285203810	handylab LF12/413T	Conductivity	handylab LF12, case set, complete, ready to use with 4-pole conductivity cell LF 413T, calibration solutions and plastic beaker	\$1,040
285203840	handylab LF12	Conductivity	handylab LF12, individual meter	\$755



## Bench Top Units/Sets

Part #	Item	Parameter	Replacement Description	List Price
285204050	Lab 850 BNC	pH	Measuring parameters pH, mV, temp., microprocessor, BNC connection. Including cover Z 880, stand and power supply Z 850 and buffers.	\$574
285204070	Lab 850 BNC Set/BL15pH	pH	Measuring parameters pH, mV, temp., microprocessor, BNC connection. Including electrode BlueLine 15pH, cover Z 880, stand and power supply Z 850 and buffers.	\$689
285204080	Lab 850 BNC Set/BL25pH	pH	Measuring parameters pH, mV, temp., microprocessor, BNC connection. Including electrode BlueLine 25pH, cover Z 880, stand and power supply Z 850 and buffers.	\$629
285204100	Lab 860 BNC meter only	pH	Meter only, measuring parameters pH, mV, temp., microprocessor, RS 232 C and USB (slave) interface, data storage for 800 data, GLP conform, BNC connection.	\$657
285204120	Lab 860 BNC Set/BL15pH	pH	Measuring parameters pH, mV, temp., micropr., RS 232 C and USB (slave) interface, storage for 800 data, GLP conform, BNC connection. Including cover Z 880, stand, power supply Z 850, pH-temp. comb. electrode BlueLine 15 pH, calibr. sol. (DIN).	\$1,050
285204170	Lab 960 meter only	Conductivity	Meter only, measuring ranges 0.000 µS/cm...500 mS/cm, salinity, total dissolved solids (TDS), temp., RS 232 C and USB (slave) interface, microprocessor, data storage for 800 data, GLP conform.	\$875
285204180	Lab 960 Set/LF413T	Conductivity	Measuring ranges 0.000 µS/cm...500 mS/cm, salinity, total dissolved solids (TDS), temp., RS 232 C and USB (slave) interface, micropr., storage for 800 data, GLP conform. Including cover Z 880, stand, power supply Z 850, cond. cell LF 413 T.	\$1,258
285204250	ProLab 3000 BNC meter only	pH/ISE 2-Channel	Meter only, sensor + user identification. QVGA colour display. Menustructured handling. Recorder function.. Dual pH, mV, Temp, ISE. RS232-, USB Host + USB Slave interfaces. BNC connection. Incl. Z850 power supply.	\$1,698
285203860	ProLab 3000 BNC Set/IL A170	pH/ISE 2-Channel	Sensor + user identification. QVGA colour display. Menustructured handling. Recorder function.. Dual pH, mV, Temp, ISE. RS232-, USB Host + USB Slave interfaces. BNC connection. Incl. Z880, Z865, Z850, ILpHT A170MF BNC-ID, DIN-Buffer.	\$1,998



## Sensors

Part #	Item	Parameter	Replacement Description	List Price
285129155	BlueLine 15 pH	pH	Glass shaft, platinum diaphragm, electrolyte KCl 3 mol/l, Ag/AgCl-reference system, temp.-sensor NTC 30 kOhm, cone membrane, A-glass, 1 m fixed cable with BNC- +4-mm banana plug, length 120 mm, 12 mm Ø, -5...+100° C, 0...14 pH	\$244
285129258	BlueLine 25 pH	pH	Plastic shaft, fibre-diaphragm, gel electrolyte, Ag/AgCl-reference system, cylinder membrane, A glass, 1 m fixed cable with BNC plug, length 120 mm, 12 mm Ø, -5...+80 °C, 0...14 pH	\$129
285129266	BlueLine 26 pH	pH	Plastic shaft, fibre-diaphragm, gel electrolyte, Ag/AgCl-reference system, temp.-sensor NTC 30 kOhm, cylinder membrane, A glass, 1 m fixed cable with BNC- +4-mm banana plug, length 120 mm, 12 mm Ø, -5...+80 °C, 0...14 pH	\$209
285106172	LF 413T, Cond.	Conductivity	4 pole cell, plastic shaft, 1.5 m cable with 8 pole plug, sensor material graphite, cell constant 0.475 cm-1, temp.-sensor NTC 30 kOhm, length 120 mm, 15.3 mm Ø, -5...+80 °C	\$509





## Chek-Mite Replacement Chart

	Replacement Model
<b>Specs</b>	<b>606pH10</b>
<b>Temperature</b>	
Range	0.0 to 100.0 °C (32.0 to 212 °F)
Resolution	0.1 °C (0.2 °F)
Accuracy	±0.3 °C (±0.6 °F)
<b>pH</b>	
Range	0.00 to 14.00 units
Resolution	0.01 units
Accuracy	±0.1 unit within 10 °C of calibration ±0.2 unit within 20 °C ±0.1 unit within 10 °C of calibration ±0.2 unit within 20 °C
<b>Additional Technical Specifications</b>	
Operating range	0.0 to 50.0 °C (32.0 to 122.0 °F)
Water resistance	IP67 waterproof case
Weight with batteries	105 grams (3.7 ounces)
ATC Probe	Thermistor, 10kΩ / 25°C
Battery	Four LR44 alkalines (included)
Battery life	>200 hrs. (low batt indicator)
# of Calibration Points	1, 2, or 3-point calibration
pH temp compensation	Auto 0.0 to 100.0 °C (32.0 to 212 °F)
pH buffer recognition	USA (4.01, 7.00, & 10.01) or NIST (4.01, 6.86, & 9.18)
pH calibration temp	0.0 to 60 °C (32.0 to 140.0 °F)
Memory	50 sets (pH, temperature, date & time)
Warranty	1 yr Instrument, 6 mo Electrode
Auto Read	"Hold" feature locks readings on display
Other	Graphic display with on-screen instr. User-replaceable single- or double-junction electrodes
What's in the box?	606pH10A instrument, batteries, manual.



Discontinued Models		
Chek-Mite pH-10	Chek-Mite pH-15	Chek-Mite pH-20
n/a	n/a	Unavailable
n/a	n/a	0.1 °C (0.2 °F)
n/a	n/a	Unavailable
0.00 to 14.00 units	0.00 to 14.00 units	0.00 to 14.00 units
0.01 units	0.01 units	0.01 units
±0.2	±0.04	±0.04
n/a	n/a	0.0 to 50.0 °C (32.0 to 122.0 °F)
IP66	IP66	IP66
115 grams (approx. 5 ozs.)	115 grams (approx. 5 ozs.)	115 grams (approx. 5 ozs.)
None	None	Unavailable
Two CR-2032 (included)	Two CR-2032 (included)	Two CR-2032 (included)
>100 hrs.	>100 hrs.	>100 hrs.
2 of 3	2 of 3	2 of 3
Set @ 25°C	Set @ 25°C	0.0 to 50.0 °C (32.0 to 122.0 °F)
pH 4, 7, 10.01	pH 4, 7, 10.01	pH 4, 7, 10.01
Set @ 25°F	Set @ 25°F	Unavailable
None	None	None
6 Months	6 Months	6 Months
"Hold" feature locks readings on display	"Hold" feature locks readings on display	"Hold" feature locks readings on display
LCD readout	LCD readout	LCD readout
User replaceable electrodes	User replaceable electrodes	User replaceable electrodes
Chek-Mite 10, battery, manual.	Chek-Mite 15, battery, manual.	Chek-Mite 20, battery, manual.